



***GE Medical Systems***

---

# **Technical Publications**

**2233657–100**

**Revision 0**

**LOGIQ™ 500**

**CE**<sub>0459</sub>

**Advanced Reference Manual**

**Volume 1**

Copyright © 1999 By General Electric Co.

**Operating Documentation**



## *Regulatory Requirement*

This product complies with regulatory requirements of the following European Directive 93/42/EEC concerning medical devices



This manual is a reference for the LOGIQ™ 500 MD MR3 Plus. It applies to all versions of 4.2 software for the LOGIQ™ 500.



***GE Medical Systems***

---

*GE Medical Systems: Telex 3797371  
P.O. Box 414, Milwaukee, Wisconsin 53201 U.S.A.  
(Asia, Pacific, Latin America, North America)*

*GE Medical Systems — Europe: Tel: +49 (0) 212 28 02 208  
Beethovenstraße 239, Postfach 11 05 60, D-42655 Solingen  
GERMANY*



## Revision History

| REV | DATE           | REASON FOR CHANGE |
|-----|----------------|-------------------|
| 0   | April 16, 1999 | Initial Release   |

### LIST OF EFFECTIVE PAGES

| PAGE<br>NUMBER                  | REVISION<br>NUMBER | PAGE<br>NUMBER                  | REVISION<br>NUMBER |
|---------------------------------|--------------------|---------------------------------|--------------------|
| Volume 1 Title Page             | 0                  | Volume 2 Title Page             | 0                  |
| Revision History A and B        | 0                  | Revision History A and B        | 0                  |
| Regulatory Requirements 1 and 2 | 0                  | Regulatory Requirements 1 and 2 | 0                  |
| Table of Contents 1 thru 16     | 0                  | Table of Contents 1 thru 16     | 0                  |
| 1–1 thru 1–14                   | 0                  | 13–1 thru 13–36                 | 0                  |
| 2–1 thru 2–22                   | 0                  | 14–1 thru 14–106                | 0                  |
| 3–1 thru 3–50                   | 0                  | 15–1 thru 15–48                 | 0                  |
| 4–1 thru 4–8                    | 0                  | 16–1 thru 16–60                 | 0                  |
| 5–1 thru 5–88                   | 0                  | 17–1 thru 17–120                | 0                  |
| 6–1 thru 6–26                   | 0                  | 18–1 thru 18–40                 | 0                  |
| 7–1 thru 7–26                   | 0                  | 19–1 thru 19–50                 | 0                  |
| 8–1 thru 8–12                   | 0                  | 20–1 thru 20–18                 | 0                  |
| 9–1 thru 9–72                   | 0                  | Index 1 thru 24                 | 0                  |
| 10–1 thru 10–104                | 0                  |                                 |                    |
| 11–1 thru 11–20                 | 0                  |                                 |                    |
| 12–1 and 12–10                  | 0                  |                                 |                    |
| Index 1 thru 24                 | 0                  |                                 |                    |

Please verify that you are using the latest revision of this document. Information pertaining to this document is maintained on GPC (GE Medical Systems Global Product Configuration). If you need to know the latest revision, contact your distributor, local GE Sales Representative or in the USA call the GE Ultrasound Clinical Answer Center at 1-800-682-5327 or 414-524-5698.



## *Revision History*

---

This page left blank intentionally.





# Regulatory Requirements

This product complies with the regulatory requirements of the following:

- Council Directive 93/42/EEC concerning medical devices:  
the  <sup>0459</sup> label affixed to the product testifies compliance to the Directive.  
The location of the CE marking is shown on 2–21 of this manual.  
European registered place of business:  
*GE Medical Systems Europe*  
Quality Assurance Manager  
BP 34  
F 78533 BUC CEDEX France  
Tel: +33 (0)1 30 70 40 40
- Medical Device Good Manufacturing Practice Manual issued by the FDA (Food and Drug Administration, Department of Health, USA).
- Underwriters' Laboratories, Inc. (UL), an independent testing laboratory.
- Canadian Standards Association (CSA).
- International Electrotechnical Commission (IEC), international standards organizations, when applicable.



**For USA  
Only**

Caution: United States law restricts this device to sale or use by or on the order of a physician.

- *General Electric Medical Systems* is ISO 9001 and EN 46001 certified.
- The original document was written in English.



## *Regulatory Requirements*

---

**NOTE:** This equipment generates, uses and can radiate radio frequency energy. The equipment may cause radio frequency interference to other medical and non-medical devices and radio communications. To provide reasonable protection against such interference, this product complies with emissions limits for a Group 1, Class A Medical Devices Directive as stated in EN 60601-1-2. However, there is no guarantee that interference will not occur in a particular installation.

**NOTE:** If this equipment is found to cause interference (which may be determined by turning the equipment on and off), the user (or qualified service personnel) should attempt to correct the problem by one or more of the following measure(s):

- reorient or relocate the affected device(s)
- increase the separation between the equipment and the affected device
- power the equipment from a source different from that of the affected device
- consult the point of purchase or service representative for further suggestions

**NOTE:** The manufacturer is not responsible for any interference caused by using other than recommended interconnect cables or by unauthorized changes or modifications to this equipment. Unauthorized changes or modifications could void the users' authority to operate the equipment.

**NOTE:** To comply with the regulations on electromagnetic interference for a Class A FCC Device, all interconnect cables to peripheral devices must be shielded and properly grounded. Use of cables not properly shielded and grounded may result in the equipment causing radio frequency interference in violation of the FCC regulations.

**NOTE:** Do not use devices which intentionally transmit RF Signals (cellular phones, transceivers, or radio controlled products) in the vicinity of the equipment as it may cause performance outside the published specifications. Keep the power to these type devices turned off when near this equipment.

The medical staff in charge of this equipment is required to instruct technicians, patients, and other people who may be around this equipment to fully comply with the above requirement.



# Table of Contents

## VOLUME 1

### Front Matter

|                         |                     |
|-------------------------|---------------------|
| Title Page              |                     |
| Revision History        | A                   |
| Regulatory Requirements | Regulatory Req 1    |
| Table of Contents       | Table of Contents 1 |

### Chapter 1—Introduction

|                                          |             |
|------------------------------------------|-------------|
| <b>System Overview</b>                   | <b>1-2</b>  |
| Attention                                | 1-2         |
| Documentation                            | 1-2         |
| Introduction                             | 1-3         |
| Physical Principle Used                  | 1-4         |
| Prescription Device                      | 1-4         |
| General Indications for Use              | 1-5         |
| Contraindications                        | 1-6         |
| LOGIQ™ 500's Features                    | 1-6         |
| <b>Who To Contact</b>                    | <b>1-7</b>  |
| Contacting GE Medical Systems—Ultrasound | 1-7         |
| Manufacturer                             | 1-10        |
| <b>How This Book is Organized</b>        | <b>1-11</b> |
| Manual Content                           | 1-11        |
| Manual Format                            | 1-13        |

### Chapter 2—Safety

|                                |            |
|--------------------------------|------------|
| <b>Safety Precautions</b>      | <b>2-2</b> |
| Precaution Levels              | 2-2        |
| Hazard Symbols                 | 2-3        |
| Patient Safety                 | 2-4        |
| Equipment and Personnel Safety | 2-5        |
| Device Labels                  | 2-7        |
| Acoustic Output                | 2-14       |
| Warning Label Locations        | 2-16       |



## Table of Contents

---

### Chapter 3—Preparing the System for Use

|                                                 |             |
|-------------------------------------------------|-------------|
| <b>Site Requirements</b> .....                  | <b>3-2</b>  |
| Introduction .....                              | 3-2         |
| Before the system arrives .....                 | 3-3         |
| Environmental Requirements .....                | 3-3         |
| <b>Console Overview</b> .....                   | <b>3-4</b>  |
| Console graphics .....                          | 3-4         |
| Peripheral/Accessory Connection .....           | 3-7         |
| <b>System Positioning/Transporting</b> .....    | <b>3-10</b> |
| Moving the System .....                         | 3-10        |
| Transporting the System .....                   | 3-13        |
| Wheels .....                                    | 3-14        |
| <b>Powering On the System</b> .....             | <b>3-15</b> |
| Connecting and Using the System .....           | 3-15        |
| <b>Adjusting the Display Monitor</b> .....      | <b>3-22</b> |
| Rotate, tilt, raise and lower the monitor ..... | 3-22        |
| Brightness and Contrast .....                   | 3-23        |
| Speakers .....                                  | 3-26        |
| <b>Probes</b> .....                             | <b>3-27</b> |
| Introduction .....                              | 3-27        |
| Selecting a probe .....                         | 3-27        |
| Connecting the Probe .....                      | 3-27        |
| Cable Handling .....                            | 3-28        |
| Activating the Probe .....                      | 3-29        |
| Deactivating the Probe .....                    | 3-30        |
| Disconnecting the Probe .....                   | 3-30        |
| Transporting Probes .....                       | 3-30        |
| Storing the Probe .....                         | 3-30        |
| <b>Operator Controls</b> .....                  | <b>3-31</b> |
| Control Panel Map .....                         | 3-31        |
| Key Illumination .....                          | 3-32        |
| Keyboard .....                                  | 3-33        |
| Soft Menu Control Panel .....                   | 3-34        |
| Mode, Display and Record .....                  | 3-46        |
| Measurement and Annotation .....                | 3-48        |
| VCR Controls .....                              | 3-50        |



---

## Table of Contents

---

### Chapter 4—Preparing for an Exam

|                                                |            |
|------------------------------------------------|------------|
| <b>Beginning an Exam</b> .....                 | <b>4-2</b> |
| Introduction .....                             | 4-2        |
| Beginning a New Patient .....                  | 4-3        |
| ID/Name .....                                  | 4-6        |
| <b>Exam Application Preset Selection</b> ..... | <b>4-7</b> |
| Introduction .....                             | 4-7        |
| Selecting a probe .....                        | 4-7        |

### Chapter 5—Modes

|                                                           |             |
|-----------------------------------------------------------|-------------|
| <b>B-Mode</b> .....                                       | <b>5-3</b>  |
| Introduction .....                                        | 5-3         |
| B-Mode Key Operation .....                                | 5-3         |
| Reading the B-Mode Display .....                          | 5-4         |
| Optimizing the Image .....                                | 5-10        |
| <b>Adding Color</b> .....                                 | <b>5-25</b> |
| Gray Scale Color .....                                    | 5-25        |
| Color Flow Mode .....                                     | 5-25        |
| Activating Color Flow .....                               | 5-26        |
| Reading the Color Flow Display .....                      | 5-27        |
| Optimizing the Color Flow Image .....                     | 5-28        |
| Power Doppler Imaging (option) .....                      | 5-41        |
| <b>Doppler</b> .....                                      | <b>5-42</b> |
| Introduction .....                                        | 5-42        |
| Pulsed Wave Doppler .....                                 | 5-42        |
| Continuous Wave Doppler .....                             | 5-44        |
| Reading the Doppler Display .....                         | 5-45        |
| Activating Doppler Mode .....                             | 5-47        |
| Doppler Optimization .....                                | 5-48        |
| <b>M-Mode</b> .....                                       | <b>5-62</b> |
| Introduction .....                                        | 5-62        |
| Reading the M-Mode or Doppler Spectrum Only Display ..... | 5-62        |
| Reading the Dual Doppler Spectrum Only Display .....      | 5-63        |
| Optimizing the Timeline .....                             | 5-64        |



## Table of Contents

---

|                                                  |             |
|--------------------------------------------------|-------------|
| <b>3Dview Mode (Option)</b> .....                | <b>5-70</b> |
| Overview .....                                   | 5-70        |
| Functionality in 3D-Mode .....                   | 5-71        |
| Functionality while MIP Image is Displayed ..... | 5-71        |
| Mode Changes in 3D-Mode .....                    | 5-71        |
| CFM Map .....                                    | 5-72        |
| Activating 3D-Mode .....                         | 5-72        |
| Creating a MIP Image Rendering .....             | 5-75        |
| 3D Option Techniques .....                       | 5-77        |
| <b>3D-Surface Mode (Option)</b> .....            | <b>5-79</b> |
| Overview .....                                   | 5-79        |
| Standard Procedure .....                         | 5-80        |
| <b>Mixed Mode Display Formats</b> .....          | <b>5-84</b> |
| Display Formats .....                            | 5-84        |

## Chapter 6—Scanning/Display Functions

|                                                |            |
|------------------------------------------------|------------|
| <b>Zooming an Image</b> .....                  | <b>6-2</b> |
| Introduction .....                             | 6-2        |
| Zoom Methods .....                             | 6-2        |
| Zooming an M-Mode Image .....                  | 6-4        |
| Multi-Image Zoom .....                         | 6-5        |
| <b>Freezing an Image</b> .....                 | <b>6-6</b> |
| Introduction .....                             | 6-6        |
| Freezing an Image (Freeze Key) .....           | 6-7        |
| Freezing an Image (Foot Switch option) .....   | 6-7        |
| <b>Using Cine</b> .....                        | <b>6-8</b> |
| Introduction .....                             | 6-8        |
| Accessing Cine .....                           | 6-9        |
| Using Cine Loop .....                          | 6-10       |
| Cine Loop Speed .....                          | 6-12       |
| Multipl CINE .....                             | 6-12       |
| Side Change .....                              | 6-13       |
| CINE Gauge .....                               | 6-13       |
| CINE Capture (option—color images only) .....  | 6-13       |
| Capture Frame (option—color images only) ..... | 6-14       |
| Exiting Cine .....                             | 6-14       |
| Helpful Hints .....                            | 6-14       |
| ECG/Cine Gauge/Image Tracking .....            | 6-14       |



## Table of Contents

|                                   |             |
|-----------------------------------|-------------|
| <b>Annotating an Image .....</b>  | <b>6–15</b> |
| Introduction .....                | 6–15        |
| Annotation Library .....          | 6–17        |
| Adding Comments to an Image ..... | 6–19        |
| Special Annotation Keys .....     | 6–20        |
| Editing Annotations .....         | 6–21        |
| Body Patterns .....               | 6–22        |

## Chapter 7—General Measurements and Calculations

|                                        |            |
|----------------------------------------|------------|
| <b>Introduction .....</b>              | <b>7–2</b> |
| Overview .....                         | 7–2        |
| Measurement Controls .....             | 7–3        |
| Cursors .....                          | 7–4        |
| General Mode Measurements Method ..... | 7–4        |
| Measurement Key .....                  | 7–5        |
| General Instructions .....             | 7–6        |
| Erasing Measurements .....             | 7–6        |
| <b>Mode Measurements .....</b>         | <b>7–7</b> |
| B-Mode Measurements .....              | 7–7        |
| CFM B-Mode Measurements .....          | 7–12       |
| Doppler Mode Measurements .....        | 7–14       |
| M-Mode Measurements .....              | 7–23       |

## Chapter 8—Abdomen and Small Parts

|                                      |            |
|--------------------------------------|------------|
| <b>General Calculations .....</b>    | <b>8–2</b> |
| Overview .....                       | 8–2        |
| Volume .....                         | 8–2        |
| Angle .....                          | 8–3        |
| Stenosis Ratio (% stenosis) .....    | 8–3        |
| S/D Ratio, RI, A/B Ratio or PI ..... | 8–3        |
| Heart Rate .....                     | 8–3        |
| Trace Auto .....                     | 8–3        |
| Transf Calcs .....                   | 8–3        |
| Max PG .....                         | 8–3        |
| Mean PG .....                        | 8–3        |
| Cardiac Output (CO) .....            | 8–4        |
| Stroke Volume Ratio (SV) .....       | 8–5        |
| Heart Rate (HR) .....                | 8–6        |
| Flow Volume (FV) .....               | 8–7        |
| Trace Auto .....                     | 8–8        |
| Flow Volume Output (FVO) .....       | 8–8        |
| Helpful hints .....                  | 8–8        |
| Hip Dysplasia Measurement .....      | 8–9        |
| General Calculation Formulas .....   | 8–11       |



## Table of Contents

---

### Chapter 9—OB/GYN (Basic OB software option)

|                                                        |             |
|--------------------------------------------------------|-------------|
| <b>Exam Preparation</b> .....                          | <b>9-3</b>  |
| Overview .....                                         | 9-3         |
| <b>Fetal Doppler</b> .....                             | <b>9-4</b>  |
| Doppler Mode for Fetal Exams .....                     | 9-4         |
| <b>Acoustic Output</b> .....                           | <b>9-5</b>  |
| Considerations .....                                   | 9-5         |
| <b>OB Measurements and Formulas</b> .....              | <b>9-6</b>  |
| Introduction .....                                     | 9-6         |
| OB Format Selection .....                              | 9-6         |
| OB Measurement Soft Menus and Formulas .....           | 9-7         |
| Helpful Hints .....                                    | 9-23        |
| <b>OB Summary Reports</b> .....                        | <b>9-24</b> |
| Starting an Exam .....                                 | 9-24        |
| OB Report Page Layout .....                            | 9-25        |
| Editing the Report .....                               | 9-31        |
| Recording Summary Reports .....                        | 9-32        |
| <b>Anatomical Survey</b> .....                         | <b>9-33</b> |
| Overview .....                                         | 9-33        |
| Editing .....                                          | 9-34        |
| User Programmed Features .....                         | 9-34        |
| <b>OB Graphs</b> .....                                 | <b>9-35</b> |
| Overview .....                                         | 9-35        |
| OB Graph Selection .....                               | 9-36        |
| <b>Fetal Trend Management (software option)</b> .....  | <b>9-39</b> |
| Overview .....                                         | 9-39        |
| Storing Patient Information .....                      | 9-39        |
| Growth Trending .....                                  | 9-44        |
| List ID Management .....                               | 9-45        |
| Data List Management .....                             | 9-51        |
| <b>OB—Multigestational (software option)</b> .....     | <b>9-53</b> |
| Overview .....                                         | 9-53        |
| Patient Entry Menu .....                               | 9-53        |
| Distinguishing Each Fetus .....                        | 9-54        |
| Measurements/Calculations .....                        | 9-54        |
| Change the Number of Fetuses .....                     | 9-55        |
| Report Page Layout .....                               | 9-56        |
| OB Graph .....                                         | 9-57        |
| Fetal Trend Management (Multigestational Option) ..... | 9-59        |





---

## Table of Contents

|                                                 |              |
|-------------------------------------------------|--------------|
| <b>Data Management Center (DMC)</b> .....       | <b>9-61</b>  |
| Overview .....                                  | 9-61         |
| Operational Setup .....                         | 9-61         |
| Transferring OB Data .....                      | 9-62         |
| Error Messages .....                            | 9-63         |
| Patient Data Input .....                        | 9-65         |
| <b>GYN Measurements</b> .....                   | <b>9-67</b>  |
| B-Mode .....                                    | 9-67         |
| Doppler Mode .....                              | 9-68         |
| <b>GYN Summary Report</b> .....                 | <b>9-70</b>  |
| GYN Report Layout .....                         | 9-70         |
| GYN Calculation Formulas .....                  | 9-71         |
| <b>Chapter 10—Cardiology (software option)</b>  |              |
| <b>Introduction</b> .....                       | <b>10-3</b>  |
| Overview .....                                  | 10-3         |
| Report Pages .....                              | 10-5         |
| BSA Calculation Methods .....                   | 10-6         |
| Heart Rate (HR) .....                           | 10-6         |
| <b>LV Measurement Methods</b> .....             | <b>10-7</b>  |
| Cubed Method .....                              | 10-7         |
| Teichholz Method .....                          | 10-9         |
| Bullet Method .....                             | 10-11        |
| Modified Simpson's Rule Method .....            | 10-13        |
| Single Plane Ellipsoid Method .....             | 10-15        |
| Bi Plane Ellipsoid Methods .....                | 10-17        |
| <b>Additional Cardiology Calculations</b> ..... | <b>10-19</b> |
| Volume .....                                    | 10-19        |
| Angle .....                                     | 10-22        |
| % stenosis (stenosis ratio) .....               | 10-23        |
| PHT (Pressure Half Time) .....                  | 10-24        |
| MVA (Mitral Valve Area) .....                   | 10-24        |
| ET (Ejection Time) .....                        | 10-25        |
| Max PG .....                                    | 10-26        |
| Mean PG .....                                   | 10-27        |
| Trace Auto .....                                | 10-28        |
| S/D (D/S) Ratio, RI, A/B Ratio or PI .....      | 10-29        |
| Heart Rate (HR) .....                           | 10-30        |
| Transf Calcs .....                              | 10-30        |



## Table of Contents

---

|                                                                |                  |
|----------------------------------------------------------------|------------------|
| <b>ECG Option .....</b>                                        | <b>10-31</b>     |
| Overview .....                                                 | 10-31            |
| Physio Sweep Speed .....                                       | 10-32            |
| ECG Sub-Menu .....                                             | 10-32            |
| ECG Lead Placement .....                                       | 10-33            |
| ECG Sync Mark Display .....                                    | 10-33            |
| ECG Sub-Menu Page 1 .....                                      | 10-34            |
| ECG Gain Pages 2 and 3 .....                                   | 10-36            |
| ECG/Cine Gauge/Image Tracking .....                            | 10-36            |
| <br><b>Advanced Cardiac Calculations (AMCAL option) .....</b>  | <br><b>10-37</b> |
| Overview .....                                                 | 10-37            |
| Measurement Sequences .....                                    | 10-38            |
| Re-measurement .....                                           | 10-39            |
| Automatic Determination of Systole and Diastole .....          | 10-41            |
| Auto Trace Measurements .....                                  | 10-42            |
| Continuous M-Mode Measurements .....                           | 10-43            |
| Advanced Cardiac Calculations Measurement Menus .....          | 10-44            |
| Cardiac Measurements .....                                     | 10-63            |
| Customizing Measurement Sequences .....                        | 10-63            |
| Auto Sequence Programming .....                                | 10-66            |
| Advanced Cardiac Specification Tables .....                    | 10-67            |
| LV Calculation Formulas (Cubed Method) .....                   | 10-67            |
| LV Calculation Formulas (Teichholz Method) .....               | 10-68            |
| LV Calculation Formulas (Bullet Method) .....                  | 10-69            |
| LV Calculation Formulas (LV SP-DISC Method) .....              | 10-69            |
| LV Calculation Formulas (LV BP-DISC Method) .....              | 10-71            |
| LV Calculation Formulas (Modified Simpson's Rule Method) ..... | 10-73            |
| LV Calculation Formulas (Single Plane Ellipsoid Method) .....  | 10-73            |
| LV Calculation Formulas (Bi Plane Ellipsoid Method) .....      | 10-74            |
| LV Calculation Formulas (Gibson Method) .....                  | 10-75            |
| B-Mode Analysis – Parasternal Long Axis .....                  | 10-77            |
| B-Mode Analysis – Parasternal Short Axis (PSAX-AV) .....       | 10-78            |
| B-Mode Analysis – Parasternal Short Axis (PSAX-MV) .....       | 10-79            |
| B-Mode Analysis – Parasternal Short Axis (PSAX-PAP) .....      | 10-80            |
| B-Mode Analysis – Apical 4 Chamber (AP-4CH) .....              | 10-81            |
| B-Mode Analysis – Apical 2 Chamber (AP-2CH) .....              | 10-83            |
| M-Mode Analysis – Left/Right Ventricle (M-LV/RV) .....         | 10-84            |
| M-Mode Analysis – Mitral Valve (M-MV) .....                    | 10-85            |
| M-Mode Analysis – Aortic Valve (M-AV) .....                    | 10-86            |
| M-Mode Analysis – Pulmonic Valve (M-PV) .....                  | 10-87            |
| M-Mode Analysis – Tricuspid Valve (M-TV) .....                 | 10-88            |
| Doppler Analysis – Mitral Valve (D-MV) .....                   | 10-89            |
| Doppler Analysis – Aortic Valve (D-AV) .....                   | 10-91            |
| Doppler Analysis – Pulmonic Valve (D-PV) .....                 | 10-93            |
| Doppler Analysis – Tricuspid Valve (D-TV) .....                | 10-95            |
| Advanced Cardiac Reports .....                                 | 10-97            |



---

## Table of Contents

---

### Chapter 11—Vascular (software option)

|                                                  |             |
|--------------------------------------------------|-------------|
| <b>Exam Preparation</b> .....                    | <b>11-2</b> |
| Introduction .....                               | 11-2        |
| General Guidelines .....                         | 11-2        |
| <b>Measurements</b> .....                        | <b>11-3</b> |
| Carotid Artery Measurements .....                | 11-3        |
| Heart Rate (HR) .....                            | 11-5        |
| Trace Auto .....                                 | 11-5        |
| Transf Calcs .....                               | 11-5        |
| Helpful Hints .....                              | 11-5        |
| <b>Vascular Summary Report</b> .....             | <b>11-6</b> |
| Introduction .....                               | 11-6        |
| Displaying the Summary Report .....              | 11-6        |
| Editing the Summary Report .....                 | 11-7        |
| Recording Summary Reports .....                  | 11-7        |
| Vascular Calculation Formulas .....              | 11-8        |
| <b>Advanced Vascular (software option)</b> ..... | <b>11-9</b> |
| Overview .....                                   | 11-9        |
| Menu Selections .....                            | 11-9        |
| Report Page Layout .....                         | 11-13       |
| Venous Comments .....                            | 11-16       |
| Vascular Calculation Formulas .....              | 11-18       |

### Chapter 12—Urology

|                                                    |                |
|----------------------------------------------------|----------------|
| <b>Urology Basic Calculations</b> .....            | <b>12-2</b>    |
| Overview .....                                     | 12-2           |
| Presumed Circle Area Ratio (PCAR) .....            | 12-3           |
| Stepper Volume (STVOL) .....                       | 12-3           |
| <b>Urology Calculation (software option)</b> ..... | <b>12-4</b>    |
| Urology Summary Report .....                       | 12-4           |
| Stepper Volume Calculation .....                   | 12-7           |
| <b>Index</b> .....                                 | <b>Index 1</b> |



## Table of Contents

---

### VOLUME 2

#### Front Matter

|                         |                     |
|-------------------------|---------------------|
| Title Page              |                     |
| Revision History        | A                   |
| Regulatory Requirements | Regulatory Req 1    |
| Table of Contents       | Table of Contents 1 |

#### Chapter 13—Recording Images

|                             |             |
|-----------------------------|-------------|
| <b>Recording Images</b>     | <b>13–2</b> |
| Image Memory                | 13–2        |
| Recall                      | 13–4        |
| Helpful hints               | 13–4        |
| Peripheral Devices          | 13–5        |
| Video Signal Specifications | 13–22       |
| Maintenance                 | 13–22       |
| MOD Image Archive (option)  | 13–23       |

#### Chapter 14—Customizing Your System

|                                                        |             |
|--------------------------------------------------------|-------------|
| <b>Time Adjustment</b>                                 | <b>14–3</b> |
| Overview                                               | 14–3        |
| Time Adjustment                                        | 14–3        |
| <b>Preset Parameters</b>                               | <b>14–5</b> |
| Overview                                               | 14–5        |
| <b>Custom Display</b>                                  | <b>14–7</b> |
| Overview                                               | 14–7        |
| Parameter Menu Command Lines                           | 14–8        |
| Changing a Parameter Value                             | 14–10       |
| Custom Display Contents                                | 14–12       |
| Page 1 of 18 (Imaging Parameter 1 – Probe Dependent 1) | 14–13       |
| Page 2 of 18 (Imaging Parameter 2 – Probe Dependent 2) | 14–15       |
| Page 3 of 18 (Imaging Parameter 3 – Probe Dependent 3) | 14–18       |
| Page 3 of 18 (Imaging Parameter 3 – Probe Dependent 3) | 14–19       |
| Page 4 of 18 (Imaging Parameter 4 – Probe Dependent 4) | 14–20       |
| Page 5 of 18 (Imaging Parameter 5 – Probe Dependent 5) | 14–22       |
| Page 6 of 18 (Imaging Parameter 6 – Probe Dependent 6) | 14–22       |
| Page 7 of 18 (Imaging Parameter 7 – Probe Dependent 7) | 14–25       |
| Page 8 of 18 (Imaging Parameter 8 – Probe Dependent 8) | 14–26       |
| Page 9 of 18 (Imaging Parameter 9 – Probe Dependent 9) | 14–29       |
| Page 10 of 18 (Imaging Parameter 10)                   | 14–30       |
| Page 11 of 18 (Imaging Parameter 11)                   | 14–32       |
| Page 12 of 18 (Imaging Parameter 12)                   | 14–34       |
| Page 13 of 18 (Imaging Parameter 13)                   | 14–36       |
| Page 14 of 18 (Imaging Parameter 14)                   | 14–39       |



## Table of Contents

|                                                             |              |
|-------------------------------------------------------------|--------------|
| Page 15 of 18 (Imaging Parameter 15) .....                  | 14-40        |
| Page 16 of 18 (Imaging Parameter 16) .....                  | 14-42        |
| Page 17 of 18 (Imaging Parameter 17) .....                  | 14-43        |
| Page 18 of 18 (Imaging Parameter 18) .....                  | 14-45        |
| <b>System Parameters .....</b>                              | <b>14-47</b> |
| Overview .....                                              | 14-47        |
| System Parameters Contents .....                            | 14-48        |
| Page 1 of 7 (System Setup) .....                            | 14-49        |
| Page 2 of 7 (System Setup) .....                            | 14-51        |
| Page 3 of 7 (System Setup) .....                            | 14-54        |
| Page 4 of 7 (System Setup – Body Pattern) .....             | 14-55        |
| Page 5 of 7 (System Setup – Recording) .....                | 14-59        |
| Page 6 of 7 (System Setup – User ID and Password) .....     | 14-61        |
| Page 7 of 7 (DICOM) .....                                   | 14-63        |
| <b>Preset Program .....</b>                                 | <b>14-64</b> |
| Overview .....                                              | 14-64        |
| Preset Program Contents .....                               | 14-65        |
| Page 1 of 10 (Application) .....                            | 14-66        |
| Page 2 of 10 (Application) .....                            | 14-67        |
| Page 3 of 10 (Application – Measurement) .....              | 14-68        |
| Page 4 of 10 (Application – Measurement) .....              | 14-71        |
| Page 4 of 10 (Application – Measurement) .....              | 14-72        |
| Page 5 & 6 of 10 (Application – Measurement Sub-Menu) ..... | 14-73        |
| Page 7 of 10 (Application – Measurement Submenu) .....      | 14-74        |
| Page 8 & 9 of 10 (Application – Annotation Library) .....   | 14-90        |
| Page 10 of 10 (Application – Patient Information) .....     | 14-90        |
| <b>Save Values .....</b>                                    | <b>14-91</b> |
| Overview .....                                              | 14-91        |
| Saving Scan Values .....                                    | 14-91        |
| <b>Exam Applications Presets .....</b>                      | <b>14-93</b> |
| Overview .....                                              | 14-93        |
| Defining a User Preset .....                                | 14-93        |
| Naming a User Preset .....                                  | 14-96        |
| Deleting User Presets and Names .....                       | 14-97        |
| Recall Preset .....                                         | 14-97        |
| <b>User Define Function .....</b>                           | <b>14-98</b> |
| Overview .....                                              | 14-98        |
| Programming the User Define Function .....                  | 14-98        |
| User Define Key Program Example .....                       | 14-100       |
| User Define Names & Lock/Unlock .....                       | 14-101       |
| Deleting User Define Functions .....                        | 14-103       |
| Saving User Define Functions .....                          | 14-103       |



## Table of Contents

---

|                                |               |
|--------------------------------|---------------|
| <b>User Data Back-up .....</b> | <b>14–104</b> |
| Overview .....                 | 14–104        |
| Saving Presets .....           | 14–104        |
| Loading Presets .....          | 14–105        |

## Chapter 15—Probes and Biopsy

|                                                           |              |
|-----------------------------------------------------------|--------------|
| <b>Probe Overview .....</b>                               | <b>15–2</b>  |
| Ergonomics .....                                          | 15–2         |
| Labeling .....                                            | 15–3         |
| Applications .....                                        | 15–6         |
| Specifications .....                                      | 15–7         |
| Probe Usage .....                                         | 15–10        |
| Care and Maintenance .....                                | 15–10        |
| Probe Safety .....                                        | 15–11        |
| Probe handling and infection control .....                | 15–13        |
| Coupling gels .....                                       | 15–19        |
| Planned Maintenance .....                                 | 15–19        |
| <b>Probe Discussion .....</b>                             | <b>15–20</b> |
| Introduction .....                                        | 15–20        |
| Curved Array (Convex) Probes .....                        | 15–21        |
| Linear Array Probes .....                                 | 15–24        |
| Sector Probes .....                                       | 15–26        |
| CWD Probes .....                                          | 15–29        |
| <b>Biopsy Special Concerns .....</b>                      | <b>15–30</b> |
| Precautions Concerning the Use of Biopsy Procedures ..... | 15–30        |
| <b>Preparing for a Biopsy .....</b>                       | <b>15–31</b> |
| Displaying the Guidezone .....                            | 15–31        |
| Determining Needle Length .....                           | 15–32        |
| Needle Guide Type Preset Selection .....                  | 15–33        |
| Preparing the Biopsy Guide Attachment .....               | 15–36        |
| E721 Probe Biopsy Guide .....                             | 15–44        |
| Biopsy Probes .....                                       | 15–48        |

## Chapter 16—User Maintenance

|                                                |             |
|------------------------------------------------|-------------|
| <b>System Data .....</b>                       | <b>16–3</b> |
| Features/Specifications .....                  | 16–3        |
| LOGIQ™ 500 Clinical Measurement Accuracy ..... | 16–5        |
| LOGIQ™ 500 Clinical Calculation Accuracy ..... | 16–6        |
| <b>Warranties .....</b>                        | <b>16–7</b> |
| Scope and Duration of Warranties .....         | 16–7        |
| Warranty Exclusions .....                      | 16–9        |
| Exclusive Warranty Remedies .....              | 16–10       |



## Table of Contents

|                                           |              |
|-------------------------------------------|--------------|
| <b>System Care and Maintenance .....</b>  | <b>16-11</b> |
| Overview .....                            | 16-11        |
| Inspecting the System .....               | 16-11        |
| Weekly Maintenance .....                  | 16-12        |
| Other Maintenance .....                   | 16-16        |
| <b>Troubleshooting .....</b>              | <b>16-19</b> |
| Introduction .....                        | 16-19        |
| Trouble images .....                      | 16-19        |
| Loose cables .....                        | 16-25        |
| Display Messages .....                    | 16-26        |
| System Error Message Description .....    | 16-27        |
| Operation Error Message Description ..... | 16-28        |
| Operation Guide Message Description ..... | 16-32        |
| Warning Message Description .....         | 16-33        |
| <b>Operator Diagnostics .....</b>         | <b>16-34</b> |
| Introduction .....                        | 16-34        |
| Probe Selection .....                     | 16-34        |
| Accessing Diagnostics .....               | 16-35        |
| System Test 1 (reduced) .....             | 16-36        |
| Test Pattern Black & White .....          | 16-37        |
| Test Pattern Color .....                  | 16-37        |
| Test Pattern Graphics .....               | 16-37        |
| <b>Assistance .....</b>                   | <b>16-38</b> |
| Clinical Questions .....                  | 16-38        |
| Service Questions .....                   | 16-38        |
| Literature .....                          | 16-38        |
| Accessories .....                         | 16-38        |
| Supplies/Accessories .....                | 16-39        |
| <b>Quality Assurance .....</b>            | <b>16-44</b> |
| Introduction .....                        | 16-44        |
| Typical Tests to Perform .....            | 16-45        |
| Baselines .....                           | 16-48        |
| Periodic Checks .....                     | 16-48        |
| Results .....                             | 16-49        |
| System Setup .....                        | 16-50        |
| Test Procedures .....                     | 16-50        |
| Setting up a Record Keeping System .....  | 16-59        |



## Table of Contents

---

### Chapter 17—Acoustic Output

|                                                             |              |
|-------------------------------------------------------------|--------------|
| <b>Bioeffects</b> .....                                     | <b>17-2</b>  |
| Concerns Surrounding the Use of Diagnostic Ultrasound ..... | 17-2         |
| Operator Awareness and Actions to Minimize Bioeffect .....  | 17-5         |
| Implementing ALARA Methods .....                            | 17-9         |
| Clinical instructions for fetal use .....                   | 17-10        |
| Efficacy of Fetal Doppler .....                             | 17-13        |
| Variance studies for fetal Doppler measurements .....       | 17-16        |
| Training and User Assistance .....                          | 17-17        |
| <b>FDA Acoustic Output Tables</b> .....                     | <b>17-18</b> |
| Maximum output summary .....                                | 17-18        |
| Maximum Thermal Indices .....                               | 17-35        |
| Key to Tables .....                                         | 17-81        |
| Measurement Precision and Uncertainty .....                 | 17-82        |
| Acoustic Output Display Operation and Accuracy .....        | 17-82        |
| Endnotes .....                                              | 17-83        |
| <b>IEC Acoustic Output Tables</b> .....                     | <b>17-85</b> |
| Key to Tables .....                                         | 17-85        |

### Chapter 18—OB Tables

|                         |      |
|-------------------------|------|
| List of OB Tables ..... | 18-1 |
|-------------------------|------|

### Chapter 19—Video Cassette Recorder

|                                         |             |
|-----------------------------------------|-------------|
| <b>VCR Operating Instructions</b> ..... | <b>19-2</b> |
| Operating Manuals .....                 | 19-2        |
| Recording .....                         | 19-2        |
| Cassette tapes .....                    | 19-2        |
| VTR—PB function .....                   | 19-2        |
| Introduction of VCR Features .....      | 19-3        |
| Safety .....                            | 19-5        |
| Setting Up the VCR .....                | 19-7        |
| Starting the VCR .....                  | 19-14       |
| Recording/Playback/Image Search .....   | 19-28       |
| Troubleshooting .....                   | 19-44       |





## Table of Contents

---

### Chapter 20—DICOM

|                                   |                |
|-----------------------------------|----------------|
| <b>Storage/Print Option .....</b> | <b>20-2</b>    |
| Overview .....                    | 20-2           |
| DICOM Presets .....               | 20-2           |
| Menu Access .....                 | 20-2           |
| Patient Entry Menu .....          | 20-3           |
| Network Configuration .....       | 20-3           |
| Printer Setup .....               | 20-6           |
| Host Verification .....           | 20-8           |
| Image Transfer .....              | 20-10          |
| <b>Worklist Option .....</b>      | <b>20-15</b>   |
| Overview .....                    | 20-15          |
| New Patient Selection .....       | 20-15          |
| Data Transferred .....            | 20-15          |
| Worklist (Schedule Menu) .....    | 20-16          |
| <b>Index .....</b>                | <b>Index 1</b> |



## *Table of Contents*

---

This page left blank intentionally.



# System Overview

## Attention

This manual contains enough information to operate the system safely. Advanced equipment training may be provided by a factory trained Applications Specialist for the agreed upon time period.

Read and understand all instructions in this manual before attempting to use the LOGIQ™ 500 system.

Keep this manual with the equipment at all times. Periodically review the procedures for operation and safety precautions.

## Documentation

LOGIQ™ 500 Documentation consists of three manuals:



- The Quick Start Guide (TRANSLATED) provides a step-by-step description of the basic features and operation of the LOGIQ™ 500. It is intended to be used in conjunction with the Basic User Manual in order to provide the information necessary to operate the system safely.



- The Basic User Manual (TRANSLATED) provides information needed by the user to operate the system safely. It describes basic functions of the system, safety features, operating modes, basic measurements/calculations, probes, user care and maintenance.



- The Advanced Reference Manual (ENGLISH ONLY) is intended for the trained, professional user. It contains all the information found in the Quick Start Guide and Basic User Manual, as well as information on options, advanced customization techniques and data tables.

The LOGIQ™ 500 manuals are written for users who are familiar with basic ultrasound principals and techniques. They do not include sonography training or clinical procedures.

## **Introduction**

The LOGIQ™ 500 digital Ultrasound System is a high performance ultrasound imaging system, intended for general purpose applications.

The system provides image generation in B-Mode, M-Mode, Pulsed, CW and Color Flow Doppler and Color M-Mode with all transducer types. Digital architecture allows maximum flexibility of all scanning modes and transducer types, throughout the full spectrum of operating frequencies.

All transducers are precise solid state array devices, allowing electronically controlled imaging with Phased Array Sector, Convex, Micro-convex and Steered Linear probes. Use of solid state digital designs allows a wide variety of scan parameters to be optimized including focusing, scan control, spatial resolution, temporal resolution and contrast resolution. The result is consistent generation of finely detailed anatomical resolution with excellent dynamic contrast tissue range and penetration.

LOGIQ™ 500 also features newly integrated specialized processing for Flow Data acquisition. Doppler information is displayed with low noise and clean spectral content to optimize measurements of important flow parameters. Selected probes can operate in Multifrequency Mode in order to Optimize Resolution in B-Mode and Sensitivity to flow in Doppler and Color Flow Modes.

The system display processor is highly versatile to produce the optimal set of imaging parameters and display formats without compromising important diagnostic information. Comprehensive graphical displays allow rapid and easy placement of Doppler sample volumes. In Color Flow Mode, combined B-Mode and Color Flow images can be steered independently so that optimal positioning is available in both modes.

Versatile, yet easy to use, the LOGIQ™ 500 system combines a wide variety of state-of-the-art operator features without complicating operation. The operator can customize all set-up parameters for a given mode, probe or clinical application. Operator controls have been placed in a logical clinical format with both hard controls and menu-driven soft control components. Three simultaneous probe connections allow rapid switching electronically between probes without delaying the examination.

The LOGIQ™ 500 System provides a total imaging solution for today's diverse ultrasound department needs, with investment security through reliable upgrades, application enhancements, and complete product support from GE.



## System Overview

---

### Physical Principle Used

The transmission and reception of mechanical high frequency waves through a transducer associated with a computer that creates the image in a digital memory, are used for the creation of medical ultrasound images. The spreading of mechanical ultrasound waves produces echoes when body changes density. In the case of human tissue, these echoes are created when the signal goes from an adipose tissue (fat) region to a muscular tissue region, among others. The echoes are returned through the same transducer that converts them back into electrical signals. These signals are highly amplified, processed by filters with several frequency and time response options, and finally scanned and stored in a digital memory. Once in the memory, the image can be displayed in real-time on a monitor. Several analog and digital circuits transform the electrical high frequency signals into a flow of digital signals, allowing the composition of the image in the memory. All the signal reception and transmission parameters are controlled by the main computer. Through the selection of these parameters by the operator, the system modifies the transmission and reception features allowing a wide range of uses, from obstetrics to peripheral vascular examinations. As its design is based on solid state components, the system is free from variations over time and requires very little maintenance. All the transducers are accurate solid state devices, allowing control of creation of images from convex, micro-convex and linear transducers. The use of a solid state design allows a wide range of sweep parameters that can be optimized resulting in a consistent creation of fine anatomical details with excellent penetration and dynamic contrast band in the tissue. The system features a sophisticated design, providing multiple functions of diagnostic and function setup keys. This makes the system user-friendly and easy to use.

### Prescription Device



**For USA  
Only**

Caution: United States law restricts this device to sale or use by or on the order of a physician.

## General Indications for Use

The LOGIQ™ 500 is a general purpose ultrasound imaging system intended for use in the dynamic evaluation of soft tissue and vascular diseases in the following areas:

- Head
- Neck
- Chest
- Abdomen
- Pelvis
- Male reproductive organs
- Female reproductive organs
- Limbs/Extremities
- Pregnant uterus
- Cardiac

## Indications for Fetal Doppler use

The LOGIQ™ 500 system can be used for fetal examination in Pulsed Wave Doppler, Continuous Wave Doppler, Color Flow Doppler, and Color M-Mode for the diagnosis of:

- Structural fetal cardiac anomalies for high-risk patients.
- Intrauterine growth retardation (IUGR) for high-risk patients with one or more of the following known or suspected conditions:
  - Multiple pregnancy
  - Maternal hypertension
  - Hydrops
  - Diabetes
  - Lupus
  - Placenta abnormality



## System Overview

---

### Contraindications

The system is **NOT** intended for use in the following areas:

Ophthalmic use (or any use causing the acoustic beam to pass through the eye).



Pulsed Wave Doppler, Continuous Wave Doppler, Color Flow Doppler, and Color M-Mode are not intended for routine fetal examination or screening nor are they intended for fetal examination in a low-risk population. The use of Doppler, even at minimal output levels, in fetal examination must be adjunctive with conventional fetal echocardiography and other clinical diagnostic methods, for high risk patients only.

### LOGIQ™ 500's Features

The LOGIQ™ 500 digital ultrasound system offers the following enhanced features:

#### Improved operator interface and system ergonomics

The LOGIQ™ 500 has been designed to streamline users' workflow, especially by:

- Creating intuitive user controls and prompts
- Grouping controls by mode or functionality
- Making the controls easy to recognize by touch

Assures users that with little effort and minimum time they can produce a complete exam with consistently high quality images. The sonographer can comfortably have full reach of all controls making the system easy to learn in order to perform a quality exam on any patient.

#### Improved sensitivity and resolution in each imaging mode

Benefits the user with improved acquisition and presentation of images and biometric information.

# Who To Contact

## Contacting GE Medical Systems—Ultrasound

For additional information or assistance, please contact your local distributor or the appropriate support resource listed below:

### USA

GE Medical Systems  
Ultrasound Service Engineering  
4855 W. Electric Avenue  
Milwaukee, WI 53219

TEL: (1) 800-437-1171  
FAX: (1) 414-647-4090

Customer Answer Center

TEL: (1) 800-682-5327  
(1) 414-524-5698

### CANADA

GE Medical Systems  
Ultrasound Service Engineering  
4855 W. Electric Avenue  
Milwaukee, WI 53219

TEL: (1) 800-664-0732

Customer Answer Center

TEL: (1) 414-524-5698

### LATIN & SOUTH AMERICA

GE Medical Systems  
Ultrasound Service Engineering  
4855 W. Electric Avenue  
Milwaukee, WI 53219

TEL: (1) 305-735-2304

Customer Answer Center

TEL: (1) 414-524-5698

### EUROPE

GE Ultraschall  
Deutschland GmbH & Co. KG  
Beethovenstraße 239  
Postfach 11 05 60  
D-42655 Solingen

TEL: 0130 81 6370 toll free  
TEL: (49)(0) 212.28.02.208  
FAX: (49)(0) 212.28.02.28

### ASIA

GE Medical Systems Asia  
Asia Support Center  
67-4 Takakura cho, Hachiouji-shi  
Tokyo, 192  
JAPAN

TEL: (81) 426-56-0033  
FAX: (81) 426-56-0053





## Who To Contact

---

### Contacting GE Medical Systems—Ultrasound (cont'd)

#### **ARGENTINA**

GEME S.A  
Miranda 5237  
Buenos Aires – 1407

TEL: (1) 639-1619  
FAX: (1) 567-2678

#### **AUSTRIA**

GE GesmbH Medical Systems Austria  
Prinz Eugen Strasse 8/8  
A-1040 WIEN

TEL: 0660 8459 toll free  
FAX: +43 1 505 38 74  
TLX: 136314

#### **BELGIUM**

GE Medical Systems Benelux  
Gulkenrodestraat 3  
B-2160 WOMMELGEM

TEL: 0 800 11733 toll free  
FAX: +32 0 3 320 12 59  
TLX: 72722

#### **BRAZIL**

GE Sistemas Médicos  
Av Nove de Julho 5229  
01407-907 São Paulo SP

TEL: 0800-122345  
FAX: (011) 3067-8298

#### **DENMARK**

GE Medical Systems Danmark  
Skovlytoften 4  
DK-2840 HOLTE

TEL: +45 45 51 00 55  
FAX: +45 42 42 59 89

#### **FRANCE**

GE Medical Systems  
738 rue Yves Carmen  
F-92658 BOULOGNE CEDEX

TEL: 05 49 33 71 toll free  
FAX: +33 1 46 10 01 20

#### **GERMANY**

GE Ultraschall  
Deutschland GmbH & Co. KG  
Beethovenstraße 239  
Postfach 11 05 60  
D-42655 Solingen

TEL: 0130 81 6370 toll free  
TEL: (49)(0) 212.28.02.208  
FAX: (49)(0) 212.28.02.28

#### **GREECE**

GE Medical Systems Hellas  
41, Nikolaou Plastira Street  
G-171 21 NEA SMYRNI

TEL: +30 1 93 24 582  
FAX: +30 1 93 58 414

#### **ITALY**

GE Medical Systems Italia  
Via Monte Albenza 9  
I-20052 MONZA

TEL: 1678 744 73 toll free  
FAX: +39 39 73 37 86  
TLX: 3333 28

#### **LUXEMBOURG**

TEL: 0800 2603 toll free



## Who To Contact

---

### Contacting GE Medical Systems—Ultrasound (cont'd)

#### **MEXICO**

GE Sistemas Médicos de Mexico S.A. de C.V.  
Rio Lerma #302, 1º y 2º Pisos  
Colonia Cuauhtémoc  
06500—México, D.F.

TEL: (5) 228-9600  
FAX: (5) 211-4631

#### **NETHERLANDS**

GE Medical Systems Nederland B.V.  
Atoomweg 512  
NL-3542 AB UTRECHT

TEL: 06 022 3797 toll free  
FAX: +31 304 11702

#### **POLAND**

GE Medical Systems Polska  
Krzywickiego 34  
P-02-078 WARSZAWA

TEL: +48 2 625 59 62  
FAX: +48 2 615 59 66

#### **PORTUGAL**

GE Medical Systems Portuguesa S.A.  
Rua Sa da Bandeira, 585  
Apartado 4094  
P-4002 PORTO CODEX

TEL: 05 05 33 7313 toll free  
FAX: +351 2 2084494  
TLX: 22804

#### **RUSSIA**

GE VNIIEM  
Mantulinskaya Ul. 5A  
123100 MOSCOW

TEL: +7 095 956 7037  
FAX: +7 502 220 32 59  
TLX: 613020 GEMED SU

#### **SPAIN**

GE Medical Systems España  
Hiero 1 Arturo Gimeno  
Poligono Industrial I  
E-28850 TORREJON DE ARDOZ

TEL: 900 95 3349 toll free  
FAX: +34 1 675 3364  
TLX: 22384 A/B GEMDE

#### **SWEDEN**

GE Medical Systems  
PO-BOX 1243  
S-16428 KISTA

TEL: 020 795 433 toll free  
FAX: +46 87 51 30 90  
TLX: 12228 CGRSWES

#### **SWITZERLAND**

GE Medical Systems (Schweiz) AG  
Sternmattweg 1  
CH-6010 KRIENS

TEL: 155 5306 toll free  
FAX: +41 41 421859



## *Who To Contact*

---

### **Contacting GE Medical Systems—Ultrasound (cont'd)**

#### **TURKEY**

GE Medical Systems Turkiye A.S.  
Mevluk Pehliran Sodak  
Yilmaz Han, No 24 Kat 1  
Gayretteppe  
ISTANBUL

TEL: +90 212 75 5552  
FAX: +90 212 211 2571

#### **UNITED KINGDOM**

IGE Medical Systems  
Coolidge House  
352 Buckingham Avenue  
SLOUGH  
Berkshire SL1 4ER

TEL: 0800 89 7905 toll free  
FAX: +44 753 696067

#### **OTHER COUNTRIES**

NO TOLL FREE

TEL: international code +  
33 1 39 20 0007

### **Manufacturer**

#### **GE YOKAGAWA MEDICAL SYSTEMS**

67-4 Takakura cho, Hachioji-shi  
Tokyo, 192  
JAPAN

# How This Book is Organized

## Manual Content

The LOGIQ™ 500 Advanced Reference Manual is organized to provide the information needed to start scanning right away. Detailed information is also provided for more time-intensive studies.

- **Getting started.** These sections give an overview of the system to help the operator start scanning as soon as possible.
  - *Introduction.* Information concerning indications/contraindications for use, who to contact and how this documentation is organized.
  - *Safety.* Important information concerning the safe operation of the LOGIQ™ 500 system.
  - *Preparing the System for Use.* How to prepare the system for use and a map of the control layout.
  - *Preparing for an Exam.* How to enter patient information, select an exam category and application preset.
- **Image optimization.** These sections detail how to improve image, trace, or spectral information.
  - *Modes.* How to adjust and optimize B-Mode, Color Flow, Doppler and M-Mode imaging.
  - *Scanning and Display Functions.* Information concerning Zoom, Freeze, Cine and Annotation functions.



## How This Book is Organized

---

### Manual Content (cont'd)

- **Measurements and Reports.** Shows how to do general and exam category specific measurements and calculations.
  - *General Measurements and Calculations.* Emphasis on basic measurements for each mode.
  - *Exam Categories.*
    - *Abdomen and Small Parts.*
    - *OB/GYN.*
    - *Cardiology.*
    - *Vascular.*
    - *Urology.*
- **Recording Images.** Explains the use of image archive and peripheral options.
- **Customizing your system.** Shows how to customize the system for your particular institution, clinic, or exam type.
- **Probes and Biopsy.** Provides intended uses, specifications, care and maintenance, and biopsy capability instructions for each probe.
- **User Maintenance.** Provides information concerning system specifications, error messages, user diagnostics, quality assurance, system care and assistance.
- **Data Tables.** Provides necessary data for reference.
  - *Acoustic Output.*
  - *OB Tables.*
- **Imaging/Recording Devices.** Provides information concerning interfacing with imaging/recording devices alone or on a network.
  - *VCR Operating Instructions.*
  - *DICOM.*



---

## How This Book is Organized

---

### Manual Format

Information has been arranged and provided to help find information easily and quickly.

### Finding information

**Tables of Contents**      Locate topics in the main table of contents.


**Tabs**      Chapter tabs are provided.

**Headers/Footers**      The section name and page number appear on the outer corners of every page.

**References**      See also page references that are noted.

**Index**      Meant for frequent and easy reference. Extensive tool that presents ideas, topics, terms, titles, headings, and cross references. Also, use it to find all entries of a like topic throughout the manual.

**2-Column Layout**      The right column contains text; the left column contains headers and graphics to highlight the text.




### Graphics


Graphics provide a visual guide to the text when possible.



Push the top or bottom of a rocker switch to get the desired result.



Turn rotary knobs to the left (counterclockwise) and right (clockwise).




Press a key to activate a function or change a parameter.




## How This Book is Organized

---

### Finding information (cont'd)



Move **TGC** slidepots to the left and right.



Move the **Trackball** around with the palm of a hand or fingertips.

### Notes

Notes are set in *italics*.

### Icons

Various icons highlight safety issues.



Indicates precautions or prudent use recommendations that should be used in the operation of the ultrasound system.

### References

References to other chapters appear in *italics*.



### Hints

Scanning hints help save time.

# Safety Precautions

## Precaution Levels

### Icon Description

Various levels of safety precautions may be found on the equipment and different levels of concern are identified by one of the following flag words which precede the precautionary statement.

#### **DANGER**



Indicates that a specific hazard is known to exist which through inappropriate conditions or actions will cause:

- Severe or fatal personal injury
- Substantial property damage.

#### **WARNING**



Indicates that a specific hazard is known to exist which through inappropriate conditions or actions may cause:

- Severe personal injury
- Substantial property damage.

#### **CAUTION**



Indicates that a potential hazard may exist which through inappropriate conditions or actions will or can cause:

- Minor injury
- Property damage.



Indicates precautions or prudent use recommendations that should be used in the operation of the ultrasound system, specifically:

- Use of the ultrasound system as a prescription device, under the order of a physician
- Maintaining an optimum system environment
- Using this Manual
- Notes to emphasize or clarify a point.



## Hazard Symbols

### Icon Description

Potential hazards are indicated by the following icons:

| Icon                                                                                                               | Potential Hazard                                                                                                                                                                                                                                                                                       | Usage                                                                                                                                                | Source            |
|--------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|
|  <b>Biological Hazard</b>         | <ul style="list-style-type: none"> <li>• Patient/user infection due to contaminated equipment.</li> </ul>                                                                                                                                                                                              | <ul style="list-style-type: none"> <li>• Cleaning and care instructions</li> <li>• Sheath and glove guidelines</li> </ul>                            | ISO 7000 No. 0659 |
|  <b>Electrical Hazard</b>         | <ul style="list-style-type: none"> <li>• Electrical micro-shock to patient, e.g., ventricular fibrillation initiated.</li> <li>• Electrical micro-shock to patient/user.</li> </ul>                                                                                                                    | <ul style="list-style-type: none"> <li>• Probes</li> <li>• ECG</li> <li>• Connections to back panel</li> </ul>                                       |                   |
|  <b>Moving Hazard</b>             | <ul style="list-style-type: none"> <li>• Console, accessories or optional storage devices fall on patient, user, or others.</li> <li>• Collision with persons or objects results in injury while maneuvering or during system transport.</li> <li>• Injury to user from moving the console.</li> </ul> | <ul style="list-style-type: none"> <li>• Moving</li> <li>• Using brakes</li> <li>• Transporting</li> </ul>                                           |                   |
|  <b>Acoustic Output Hazard</b>  | <ul style="list-style-type: none"> <li>• Patient injury or tissue damage from ultrasound radiation.</li> </ul>                                                                                                                                                                                         | <ul style="list-style-type: none"> <li>• ALARA, the use of acoustic output following the <u>as low as reasonably achievable</u> principle</li> </ul> |                   |
|  <b>Explosion Hazard</b>        | <ul style="list-style-type: none"> <li>• Risk of explosion if used in the presence of flammable anesthetics.</li> </ul>                                                                                                                                                                                | <ul style="list-style-type: none"> <li>• Flammable anesthetic</li> </ul>                                                                             |                   |
|  <b>Smoke &amp; Fire Hazard</b> | <ul style="list-style-type: none"> <li>• Patient/user injury or adverse reaction from fire or smoke.</li> <li>• Patient/user injury from explosion and fire.</li> </ul>                                                                                                                                | <ul style="list-style-type: none"> <li>• Replacing fuses</li> <li>• Outlet guidelines</li> </ul>                                                     |                   |
|  <b>Non-Ionizing Radiation</b>  | <ul style="list-style-type: none"> <li>• Console failure, erratic operation or output error due to RF interference.</li> </ul>                                                                                                                                                                         | <ul style="list-style-type: none"> <li>• RF</li> </ul>                                                                                               | IEC 878 No. 03-04 |

Table 2–1. Potential Hazards



## Safety Precautions

---

### Important Safety Considerations

The following topic headings (*Patient Safety*, and *Equipment and Personnel Safety*) are intended to make the equipment user aware of particular hazards associated with the use of this equipment and the extent to which injury can occur if precautions are not observed. Additional precautions may be provided throughout the manual. The equipment user is obligated to be familiar with these concerns and avoid conditions that could result in injury.

### Patient Safety

#### Related Hazards

#### WARNING



The concerns listed can seriously affect the safety of patients undergoing a diagnostic ultrasound examination.

#### Patient identification

Always include proper identification with all patient data and verify the accuracy of the patient's name or ID numbers when entering such data. Make sure correct patient ID is provided on all recorded data and hard copy prints. Identification errors could result in an incorrect diagnosis.

#### Diagnostic information

Equipment malfunction or incorrect settings can result in measurement errors or failure to detect details within the image. The equipment user must become thoroughly familiar with the equipment operation in order to optimize its performance and recognize possible malfunctions. Applications training is available through the local GE representative. Added confidence in the equipment operation can be gained by establishing a quality assurance program.

#### Mechanical hazards

Damaged probes or improper use and manipulation of intracavitary probes can result in injury or increased risk of infection. Inspect probes often for sharp, pointed, or rough surface damage that could cause injury or tear protective barriers. Never use excessive force when manipulating intracavitary probes. Become familiar with all instructions and precautions provided with special purpose probes.



#### Electrical Hazard

A damaged probe can also increase the risk of electric shock if conductive solutions come in contact with internal live parts. Inspect probes often for cracks or openings in the housing and holes in and around the acoustic lens or other damage that could allow liquid entry. Become familiar with the probe's use and care precautions outlined in *Probes and Biopsy*.



## Safety Precautions

---

### Related Hazards (cont'd)



#### Acoustic Output Hazard

Ultrasound energy, even at diagnostic levels, is capable of damaging sensitive tissues if adequate precautions are not followed. The wrong combination of equipment settings, probe positioning, and tissue type can result in injury. Please become thoroughly familiar with equipment controls that affect acoustic output levels as well as the output display.

Follow the principle of as low as reasonably achievable (ALARA) when scanning a patient. During each ultrasound examination, the clinical user is expected to weigh the medical benefit of the diagnostic information obtained against the risk of potential harmful effects. Once an optimal image is achieved the need for increasing acoustic output or prolonging the exposure cannot be justified.

#### Training

It is recommended that all users receive proper training in applications before performing them in a clinical setting. Please contact the local GE representative for training assistance. ALARA training is provided by GE Application Specialists.

## Equipment and Personnel Safety

### Related Hazards

#### WARNING



This equipment contains dangerous voltages that are capable of serious injury or death.

There are no user serviceable components inside the console. Refer all servicing to qualified service personnel only.

#### DANGER



The concerns listed below can seriously affect the safety of equipment and personnel during a diagnostic ultrasound examination.



#### Explosion Hazard

Risk of explosion if used in the presence of flammable anesthetics.



## Safety Precautions

---

### Related Hazards (cont'd)



#### Electrical Hazard

To avoid injury:

- Do not remove protective covers. No user serviceable parts are inside. Refer servicing to qualified service personnel.
- To assure adequate grounding, connect the attachment plug to a reliable (hospital grade) grounding outlet (having equalization conductor ).
- Do not place liquids on or above the console. Spilled liquid may contact live parts and increase the risk of shock.



#### Smoke & Fire Hazard

The system must be supplied from an adequately rated electrical circuit. The capacity of the supply circuit must be as specified on 3–3.



#### Biological Hazard

For patient and personnel safety, beware of biological hazards while performing invasive procedures. To avoid the risk of disease transmission:

- Use protective barriers (gloves and probe sheaths) whenever possible. Follow sterile procedures when appropriate.
- Thoroughly clean probes and reusable accessories after each patient examination and disinfect or sterilize as needed. Refer to *Probes and Biopsy* for probe use and care instructions.
- Follow all infection control policies established by your office, department or institution as they apply to personnel and equipment.

### CAUTION



Devices containing latex may cause severe allergic reaction in latex sensitive individuals. USA customers should refer to the FDA's March 29, 1991 Medical Alert on latex products.

## Device Labels

### Label Icon Description

The following table describes the purpose and location of safety labels and other important information provided on the equipment.

| Label/Icon                                                                          | Purpose/Meaning                                                                                                                                                                                                  | Location                           |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|
| Identification and Rating Plate                                                     | <ul style="list-style-type: none"> <li>Manufacturer's name and address</li> <li>Date of manufacture</li> <li>Model and serial numbers</li> <li>Electrical ratings (Volts, Amps, phase, and frequency)</li> </ul> | Rear of console near power inlet   |
| Type/Class Label                                                                    | Used to indicate the degree of safety or protection.                                                                                                                                                             |                                    |
| IP Code (IPX1)                                                                      | Indicates the degree of protection provided by the enclosure per IEC 529. IPX1 indicates drip proof.                                                                                                             | Foot Switch                        |
|   | Equipment Type BF (man in the box symbol) IEC 878-02-03 indicates B Type equipment having a floating applied part.                                                                                               | Probe connectors and PCG connector |
|  | Equipment Type CF (heart in the box symbol) IEC 878-02-05 indicate equipment having a floating applied part having a degree of protection suitable for direct cardiac contact.                                   | ECG connector and surgical probes  |
| Device Listing/Certification Labels                                                 | Laboratory logo or labels denoting conformance with industry safety standards such as UL or IEC.                                                                                                                 | Rear of console                    |
| "DANGER – Risk of explosion used in..."                                             | The system is not designed for use with flammable anesthetic gases.                                                                                                                                              |                                    |
|  | "CAUTION" The equilateral triangle is usually used in combination with other symbols to advise or warn the user.                                                                                                 | Various                            |

Table 2–2. Label Icons

## Safety Precautions

### Label Icon Description (cont'd)

| Label/Icon                                                                          | Purpose/Meaning                                                                                                                                                                                                                                           | Location                                 |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------|
|    | "ATTENTION – Consult accompanying documents" is intended to alert the user to refer to the operator manual or other instructions when complete information cannot be provided on the label.                                                               | Various                                  |
|    | "CAUTION – Dangerous voltage" (the lightning flash with arrowhead) is used to indicate electric shock hazards.                                                                                                                                            |                                          |
|    | "Mains OFF" Indicates the power off position of the mains power switch.                                                                                                                                                                                   | Rear of system, adjacent to mains switch |
|    | "Mains ON" Indicates the power on position of the mains power switch.                                                                                                                                                                                     |                                          |
|   | "ON" Indicates the power on position of the power switch.<br><b>CAUTION:</b> This Power Switch <b>DOES NOT ISOLATE</b> Mains Supply.                                                                                                                      | Adjacent to On-Off/Standby Switch        |
|  | "Off/Standby" Indicates the power off/standby position of the power switch.<br><b>CAUTION:</b> This Power Switch <b>DOES NOT ISOLATE</b> Mains Supply.                                                                                                    |                                          |
|  | "Equipotentiality" Indicates the terminal to be used for connecting equipotential conductors when interconnecting (grounding) with other equipment.<br><b>CAUTION:</b> This is only for " <b>FUNCTIONAL GROUNDING</b> ", NOT " <b>PROTECTIVE EARTH</b> ". | Rear of console                          |

Table 2–2. Label Icons (cont'd)



## Safety Precautions

---

### Classifications

Type of protection against electric shock  
Class I Equipment (\*1)

Degree of protection against electric shock  
Type BF Equipment (\*2) (Except ECG)  
Type CF Equipment (\*3) (ECG Only)  
Ordinary Equipment  
Continuous Operation

#### \*1. Class I EQUIPMENT

EQUIPMENT in which protection against electric shock does not rely on BASIC INSULATION only, but includes an earth ground. This additional safety precaution prevents exposed metal parts from becoming LIVE in the event of an insulation failure.

#### \*2. Type BF EQUIPMENT

TYPE B EQUIPMENT with an F-TYPE APPLIED PART  
TYPE B EQUIPMENT: EQUIPMENT providing a specified degree of protection against electric shock, with particular regard to allowable LEAKAGE CURRENT.

|                         | Normal Mode           | Single fault condition |
|-------------------------|-----------------------|------------------------|
| Patient leakage current | Less than 100 $\mu$ A | Less than 500 $\mu$ A  |

#### \*3. Type CF EQUIPMENT

EQUIPMENT providing a degree of protection higher than that for TYPE BF EQUIPMENT against electric shock particularly regarding allowable LEAKAGE CURRENTS, and having an F-TYPE APPLIED PART.

|                         | Normal Mode          | Single fault condition |
|-------------------------|----------------------|------------------------|
| Patient leakage current | Less than 10 $\mu$ A | Less than 50 $\mu$ A   |



## Safety Precautions

---

### Classifications (cont'd)

#### \*4. EMC (Electromagnetic Compatibility)

##### 4.1 EMC Performance

All types of electronic equipment may characteristically cause electromagnetic interference with other equipment, either transmitted through air or connecting cables. The term EMC (Electromagnetic Compatibility) indicates the capability of equipment to curb electromagnetic influence from other equipment and at the same time not affect other equipment with similar electromagnetic radiation from itself.

This product is designed to fully comply with the EN60601–1–2 (IEC60601–1–2) in medical electric equipment EMC regulations.

Proper installation following the service manual is required in order to achieve the full EMC performance of the product.

The product must be installed as stipulated in 4.2, Notice upon Installation of Product.

In case of issues related to EMC, please call your service personnel.

#### CAUTION



Do not use the following devices near this equipment. Use of these devices near this equipment could cause this equipment to malfunction.

- Devices which intrinsically transmit radio waves such as: Cellular phone, radio transceiver, mobile radio transmitter, radio-controlled toys, etc.

Keep power to these devices turned off when near this equipment.

Medical staff in charge of this equipment is required to instruct technicians, patients and other people who may be around this equipment to fully comply with the above regulation.



## **Classifications (cont'd)**

### **4.2 Notice upon Installation of Product**

1. Use either power supply cords provided by GE Medical Systems or ones designated by GE Medical Systems. Products equipped with a power source plug should be plugged into the fixed power socket which has the protective grounding conductor. Never use any adaptor or converter to connect with a power source plug (i.e. three-prong-to-two-prong converter).
2. Locate the equipment as far away as possible from other electronic equipment.
3. Be sure to use only the cables provided by or designated by GE Medical Systems. Connect these cables following the installation procedures (i.e. wire power cables separately from signal cables).
4. Lay out the main equipment and other peripherals following the installation procedures described in the Option Installation manuals.

### **4.3 General Notice**

1. Designation of Peripheral Equipment Connectable to This Product.

The equipment indicated on 13–5 can be hooked up to the product without compromising its EMC performance.

Avoid using equipment not designated in the list. Failure to comply with this instruction may result in poor EMC performance of the product.

2. Notice against User Modification

The user should never modify this product. User modifications may cause degradation in EMC performance.

Modification of the product includes changes in:

- a. Cables (length, material, wiring, etc.)
  - b. System installation/layout
  - c. System configuration/components
  - d. Securing system parts (cover open/close, cover screwing)
3. Operate the system with all covers closed. If a cover is opened for some reason, be sure to shut it before starting/resuming operation.

Operating the system with any cover open may affect EMC performance.

## Safety Precautions

### Classifications (cont'd)

#### \*5. Patient Environmental Devices




Figure 2-1. Patient Environmental Devices

**Classifications (cont'd)****5.1 Acceptable Devices**

The devices shown in Figure 2–1 are specified to be suitable for use within the PATIENT ENVIRONMENT.

**CAUTION**

**DO NOT** connect any probes or accessories without approval by GE.

Please refer to Peripheral Devices on 13–5 and Assistance on 16–38 for more details. The devices listed have been tested and verified to be compatible with the LOGIQ™ 500 system.

**5.2 Unapproved Devices****CAUTION**

The user takes All Responsibility for connecting unapproved devices.

If devices are connected without the approval of GE, the warranty will be **INVALID**.

Any device connected to the LOGIQ™ 500 must conform to one or more of the requirements listed below:

1. IEC 50, IEC 65, IEC 335, IEC 348, IEC 414, IEC 820, IEC 950, IEC 1010–1, ISO 7767, ISO 8185, ISO 8359 or IEC 60601–1.
2. The devices shall be connected to PROTECTIVE EARTH (GROUND).



## Safety Precautions

---

### Acoustic Output

#### Controls Affecting Output

The potential for producing mechanical or thermal bioeffects is influenced by the controls listed below (refer to 2–15).

**Direct.** The **Acoustic Output** control has the most significant effect on Acoustic Output.

**Indirect.** Indirect effects may occur when adjusting the controls listed on 2–15.

Always observe the output display for possible effects.

#### Best practices while scanning



##### Hints

- Raise the Acoustic Output **only** after attempting image optimization with controls that have **no** affect on Acoustic Output, such as **Gain** and **TGC**.



*NOTE: Refer to the Optimization sections of the Modes chapter for a complete discussion of each control.*

##### WARNING



Be sure to have read and understood control explanations for each Mode intended to be used before attempting to adjust the Acoustic Output control or any control that can affect Acoustic Output.



##### Acoustic Output Hazard

Use the minimum necessary output to get the best diagnostic image or measurement during an examination. Begin the exam with the probe that provides an optimum focal depth and penetration.

## Best practices while scanning (cont'd)

### Controls

| Mode    | Control                           | Affect                 | Default Setting                                                                                                                                                        |
|---------|-----------------------------------|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| All     | Acoustic Output                   | Direct.<br>Significant | The middle setting is a factory preset determined to be a reasonable setting for all exams. Use presets to set the output preferred by scan mode and exam combination. |
| B       | Focus Comb                        | Indirect. Minor        | Off.                                                                                                                                                                   |
| B/M/CFD | Depth (FOV)                       | Indirect. Minor        | Probe-dependent operator preset.                                                                                                                                       |
| B/M     | Focal Zone Position and Number    | Indirect. Minor        | Probe-dependent system preset.                                                                                                                                         |
| B/CFD   | Scan Area                         | Indirect.              | Off.                                                                                                                                                                   |
| Doppler | M/D Cursor                        | Indirect. Minor        | Off.                                                                                                                                                                   |
| Doppler | Doppler Sample Volume Gate Length | Indirect. Minor        | Application-dependent system preset.                                                                                                                                   |
| Doppler | Velocity Scale                    | Indirect. Minor        | Application-dependent operator preset.                                                                                                                                 |
| PWD/CFD | Freq. Hi/Low                      | Indirect. Minor        | Low.                                                                                                                                                                   |
| CFD     | Scan Area                         | Indirect. Minor        | Off.                                                                                                                                                                   |
| B/M/CFD | Zoom                              | Indirect. Minor        | Off.                                                                                                                                                                   |

Table 2–3. Controls Affecting Acoustic Output

### Acoustic Output Default Levels

In order to assure that an exam does not start at a high output level, the LOGIQ™ 500 initiates scanning at a reduced default output level. This reduced level is preset programmable and depends upon exam category, application preset and probe selected. It takes effect when the system is powered on or **New Patient** is selected.



## Safety Precautions

### Warning Label Locations

#### Overview

LOGIQ™ 500 warning labels are provided in seven different languages (English, Japanese, German, French, Italian, Portuguese and Spanish).

#### Monitor Labels

A temporary label is placed on the monitor face to warn not to move the monitor support arm without the monitor attached. This label is removed after installation of the monitor.




Figure 2–2. Temporary Warning Label and Location



Safety Precautions

Monitor Labels (cont'd)

Two caution labels are found on the back of the monitor. One warns to only move the console with the monitor in its lowest position; the second warns not to push the console from the side.





Figure 2–3. Lowering Monitor Caution Labels and Locations



## Safety Precautions

### Monitor Labels (cont'd)

One caution label is found on the top of the monitor.

| <b>⚠ CAUTION 注意 ACHTUNG ATTENTION ATTENZIONE ATENÇÃO ATENCIÓN</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Possible injury. Placing objects on top of the monitor may cause the monitor to tilt with the falling objects resulting in injury to the operator. Do not place any objects on the monitor.</p> <p>怪我の可能性あり。モニタ上に物を載せるとモニタが傾き、物が落下して怪我をすることがあります。モニタの上に物を載せないこと。</p> <p>Eventuelle Verletzungsgefahr. Auf den Monitor gelegte Gegenstände können ihn neigen und herunterfallende Gegenstände evtl. Verletzungen verursachen. Keinerlei Gegenstände auf den Monitor legen.</p> <p>Risque de blessure possible : le poids d'objets posés sur le moniteur peut le faire s'incliner. Afin d'éviter tout risque de blessure occasionné par la chute de ces objets, ne rien poser sur le moniteur.</p> | <p>Rischio di lesione. Eventuali oggetti riposti sul monitor, ne possono provocare l'inclinamento. Per evitare il rischio di caduta di tali oggetti con conseguenti ferite, non riporre alcun oggetto sul monitor.</p> <p>Existe a possibilidade de ocorrer acidentes. Caso algum objeto seja colocado sobre o monitor, este poderá inclinar-se e o objeto poderá cair, provocando algum acidente. Não coloque nenhum objeto sobre o monitor.</p> <p>Existen posibilidades de lesiones. Al colocar objetos sobre el monitor, puede inclinarse el monitor, caerse los objetos y causar lesiones. No colocar objetos sobre el monitor.</p> |




Figure 2-4. Top of Monitor Caution and Location



## Console Labels

Labels found on the back and side of the console will either be translated to the twelve languages or be specific to the region.

## Defibrillator Caution




Figure 2–5. Defibrillator Label Location

## Safety Precautions

---

### Ground Point




Figure 2–6. Signal Ground Point Location and Label

### CAUTION



This is only for **“FUNCTIONAL GROUNDING”**, NOT **“PROTECTIVE EARTH”**.

## Regulatory Labels (European Systems)

| <b>CAUTION 注意 ACHTUNG ATTENTION ATTENZIONE ATENÇÃO ATENCIÓN</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |  |                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |  |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| <p>Do not use the following devices near this equipment. Cellular phone, radio transceiver, mobile radio transmitter, radio controlled toy, etc. Use of these devices near this equipment could cause this equipment to perform outside the published specifications. Keep power to these devices turned off when near this equipment.</p> <p>本機の近くでは以下の機器を使用しないでください。本機が仕様から外れた動作をすることがあります。携帯電話、トランシーバー、携帯無線機、ラジコンのおもちゃなど。本機の近くではこれらの機器の電源スイッチを切ってください。</p> <p>In der Nähe dieses Systems folgende Geräte nicht benutzen: Funktelefone, Radiomöbiler, mobile Radiosender, ferngesteuertes Spielzeug, usw. Der Gebrauch dieser Geräte könnte während des Betriebes Fehlfunktionen verursachen. Schalten Sie die Stromversorgung dieser Geräte in der Nähe der Ausrüstung ab.</p> <p>Près de ne pas utiliser les appareils suivants à proximité de cet équipement: téléphones portables, émetteurs-récepteurs, appareils radio-commandés, etc. Ce type d'appareil peut faire dévier les performances du système en dehors des spécifications annoncées. Mettre ces appareils hors tension lorsqu'ils se trouvent à proximité de cet équipement.</p> |  |                                                                                                                                                                                                                                                                                                                                                                                | <p>Non utilizzare mai i seguenti dispositivi nelle vicinanze di questa apparecchiatura: telefoni cellulari, ricevatrasmettitori radio, trasmettitori radiomobili, giocattoli telecomandati, ecc. L'utilizzo di tali dispositivi potrebbe provocare il difunzionamento di questa apparecchiatura. Tenere spenti tali dispositivi nelle vicinanze di questa apparecchiatura.</p> <p>Favor não utilizar os seguintes aparelhos perto deste equipamento. Telefone móvel, rádio trans-receptor, rádio transmissor móvel, brinquedos de controle remoto, etc. A utilização desses aparelhos perto deste equipamento pode fazer com que não funcione segundo as especificações publicadas.</p> <p>Favor manter esses aparelhos desligados quando estiverem perto deste equipamento.</p> <p>No utilice los siguientes aparatos cerca de esta equipo. Teléfono móvil, transmisor-receptor, transmisor de radio móvil, juguete radio controlado, etc. La utilización de estos aparatos cerca de esta equipo puede hacer que el equipo no funcione según las especificaciones publicadas.</p> <p>Mantenga estos aparatos apagados cuando esté cerca de este equipo.</p> |                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |  |
| <b>WARNING 警告 WARNUNG AVERTISSEMENT AVVERTENZE ADVERTÊNCIA ADVERTENCIA</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |  |                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |  |
| <p>Possible shock hazard. Do not remove Covers or Panels. Refer servicing to qualified personnel.</p> <p>感電の危険あり。カバー（またはパネル）を取り除かないこと。修理作業はサービスセンターに依頼すること。</p> <p>Gefahr eines elektrischen Schlags! Abdeckungen oder Verkleidungen nicht entfernen. Wartungsarbeiten nur durch qualifiziertes Fachpersonal durchführen lassen.</p> <p>Risque de décharge électrique. Les capots et panneaux ne doivent être retirés que par un personnel qualifié.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |  | <p>Rischio di scosse elettriche. I pannelli di chiusura devono essere rimossi solo da personale qualificato.</p> <p>Possível choque. Não remova tampas ou painéis. Solicite assistência técnica por pessoal qualificado.</p> <p>Peligro de descargas eléctricas. No saque las cubiertas ni los paneles. El servicio debe realizarlo el personal cualificado.</p>               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | <p>For continued protection against fire or shock hazard. Replace only with same type and rating of Fuse.</p> <p>火災の危険あり。表示されたヒューズと交換のこと。</p> <p>Zum dauerhaften Geräteschutz und zur Vermeidung elektrischer Schläge nur Sicherungen des gleichen Type und Nennwertes verwenden.</p> <p>Pour éviter tout risque de décharge électrique ou d'incendie, les fusibles de remplacement doivent respecter le type et le calibre spécifiés pour la tension réseau choisie.</p> |  |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |  | <p>Sostituire i fusibili utilizzando soltanto fusibili dello stesso tipo, dimensione e amperaggio.</p> <p>Para proteção contínua contra fogo ou choque. Substitua somente com fusível do mesmo tipo e capacidade.</p> <p>Como protección continua contra descargas eléctricas e incendios, reemplace el fusible sólo por otro que sea del mismo tipo y del mismo amperaje.</p> |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | <p>Grounding reliability can only be achieved when this equipment is connected to a receptacle marked "Hospital Only" or "Hospital Grade".</p> <p>この機器は「医用コンセント（Hマーク）」に接続したときにのみ保護接地の信頼性が保証されます。</p>                                                                                                                                                                                                                                                                      |  |
| <b>CAUTION</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |  |                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |  |
| <p>United States law restricts this device to sale or use by or on the order of a physician.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |  |                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |  |
| <b>DANGER 危険 GEFAHR DANGER PERICOLO PERIGO PELIGRO</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |  |                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |  |
| <p>Possible explosion hazard if used in the presence of flammable anesthetics.</p> <p>爆発の危険あり。引火性麻酔剤のある場所では使用しないこと。</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |  | <p>Explosionsgefahr! Nicht in Gegenwart brennbarer Narkosegase verwenden.</p> <p>Risque d'explosion. Ne pas utiliser en présence d'anesthésiques inflammables.</p>                                                                                                                                                                                                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | <p>Rischio di esplosione se il sistema è usato in presenza di gas anestetici infiammabili.</p> <p>Possível explosão se usado na presença de anestésicos inflamáveis.</p>                                                                                                                                                                                                                                                                                                  |  |




Figure 2-7. Regulatory Label Location (European)



Safety Precautions

Regulatory Labels (American Systems)

**CAUTION 注意 ACHTUNG ATTENTION ATTENZIONE ATENÇÃO ATENCIÓN**

Do not use the following devices near this equipment. Cellular phone, radio transceiver, mobile radio transmitter, radio controlled toy, etc. Use of these devices near this equipment could cause this equipment to perform outside the published specifications. Keep power to these devices turned off when near this equipment.

本機の近くでは以下の機器を使用しないでください。本機が仕掛けから外れた動作をすることがあります。携帯電話、トランシーバー、携帯無線機、ラジオコンのおもちゃなど。本機の近くではこれらの機器の電源スイッチを切ってください。

In der Nähe dieses Systems folgende Geräte nicht benutzen: Funktelefone, Radioempfänger, mobile Radiosender, ferngesteuertes Spielzeug, usw. Der Gebrauch dieser Geräte könnte während des Betriebes Fehlfunktionen verursachen. Schalten Sie die Stromversorgung dieser Geräte in der Nähe der Ausrüstung ab.

Prière de ne pas utiliser les appareils suivants à proximité de cet équipement: téléphones portables, émetteurs-récepteurs, appareils radio-commandés, etc. Ce type d'appareil peut faire dévier les performances du système en dehors des spécifications annoncées. Mettre ces appareils hors tension lorsqu'ils se trouvent à proximité de cet équipement.

Non utilizzare mai i seguenti dispositivi nelle vicinanze di questa apparecchiatura: telefoni cellulari, ricevatrasmettitori radio, trasmettitori radiomobili, giocattoli telecomandati, ecc. L'utilizzo di tali dispositivi potrebbe provocare il disfunzionamento di questa apparecchiatura. Tenere spenti tali dispositivi nelle vicinanze di questa apparecchiatura.

Favor não utilizar os seguintes aparelhos perto deste equipamento. Telefone móvel, rádio trans-receptor, rádio transmissor móvel, brinquedos de controle remoto, etc. A utilização desses aparelhos perto deste equipamento pode fazer com que não funcione segundo as especificações publicadas. Favor manter esses aparelhos desligados quando estiverem perto deste equipamento.

No utilices los siguientes aparatos cerca de este equipo. Teléfono móvil, transmisor-receptor, transmisor de radio móvil, juguete radio controlado, etc. La utilización de estos aparatos cerca de este equipo puede hacer que el equipo no funcione según las especificaciones publicadas. Mantenga estos aparatos apagados cuando esté cerca de este equipo.

**WARNING 警告 WARNUNG AVERTISSEMENT AVVERTENZE ADVERTÊNCIA ADVERTENCIA**

Possible shock hazard. Do not remove Covers or Panels. Refer servicing to qualified personnel.

感電の危険あり。カバー（または背板）を取り除かないこと。修理点検はサービスセンターに依頼すること。

Gefahr eines elektrischen Schlags! Abdeckungen oder Verkleidungen nicht entfernen. Wartungsarbeiten nur durch qualifiziertes Personal durchführen lassen.

Risque de décharge électrique. Les capots et panneaux ne doivent être retirés que par un personnel qualifié.

Rischio di scosse elettriche. I pannelli di chiusura devono essere rimossi solo da personale qualificato.

Possível choque. Não remova tampas ou painéis. Solicite assistência técnica por pessoal qualificado.

Peligro de descargas eléctricas. No saque las cubiertas ni los paneles. El servicio debe realizarlo el personal cualificado.

For continued protection against fire or shock hazard. Replace only with same type and rating of Fuse.

火災の危険あり。表示されたヒューズと交換のこと。

Zum dauerhaften Geräteschutz und zur Vermeidung elektrischer Schläge nur Sicherungen des gleichen Typs und Nennwertes verwenden.

Pour éviter tout risque de décharge électrique ou d'incendie, les fusibles de remplacement doivent respecter le type et le calibre spécifiés pour la tension réseau choisie.

Sostituire i fusibili utilizzando soltanto fusibili dello stesso tipo, dimensione e amperaggio.

Para proteção contínua contra fogo ou choque. Substitua somente com fusível do mesmo tipo e capacidade.

Como protección continua contra descargas eléctricas e incendios, reemplace el fusible sólo por otro que sea del mismo tipo y del mismo amperaje.

Grounding reliability can only be achieved when this equipment is connected to a receptacle marked "Hospital Only" or "Hospital Grade".

この機器は「医用コンセント（Hマーク）」に接続したときのみ保護接地の信頼性が保証されます。

**CAUTION**

United States law restricts this device to sale or use by or on the order of a physician.

**DANGER 危険 GEFAHR DANGER PERICOLO PERIGO PELIGRO**

Possible explosion hazard if used in the presence of flammable anesthetics.

爆発の危険あり。引火性麻酔剤のある場所では使用しないこと。

Explosionsgefahr! Nicht in Gegenwart brennbarer Narkosegase verwenden.

Risque d'explosion. Ne pas utiliser en présence d'anesthésiques inflammables.

Rischio di esplosione se il sistema è usato in presenza di gas anestetici infiammabili.

Possível explosão se usado na presença de anestésicos inflamáveis.

Riesgo de explosión. No emplear en presencia de anestésicos inflamables.




Figure 2–8. Regulatory Label Location (Americas)



# Site Requirements

## Introduction



Only qualified physicians or sonographers should perform ultrasound scanning on human subjects for medical diagnostic reasons. Request training, if needed.

Do not attempt to install the system alone. General Electric, Affiliate, or Distributor Field Engineers and Application Specialists will install and setup the system. Refer to Who To Contact on 1-7.

The LOGIQ™ 500 does not contain any operator serviceable internal components. Ensure that unauthorized personnel do not tamper with the unit.

Perform regular preventive maintenance. Refer to 16-11 for maintenance instructions.

Maintain a clean environment. Turn off the system circuit breaker before cleaning the unit. Refer to 16-12 for cleaning instructions.

Never set liquids on the unit to ensure that liquid does not drip into the control panel or unit.

## Before the system arrives

### NOTICE

This medical equipment is approved, in terms of the prevention of radio wave interference, to be used in hospitals, clinics and other institutions which are environmentally qualified. The use of this equipment in an inappropriate environment may cause some electronic interference to radios and televisions around the equipment. This equipment can be used in residential areas only under the supervision of physicians or qualified technicians.

Ensure that the following is provided for the new system:

- A separate power outlet with a 20 amp circuit breaker for 120 VAC (USA) or 7.5 amp circuit breaker for 220–240 VAC (Europe, Latin America).
- Take precautions to ensure that the console is protected from electromagnetic interference.

Precautions include:

- Operate the console at least 15 feet away from motors, typewriters, elevators, and other sources of strong electromagnetic radiation.
- Operation in an enclosed area (wood, plaster or concrete walls, floors and ceilings) helps prevent electromagnetic interference.
- Special shielding may be required if the console is to be operated in the vicinity of radio broadcast equipment.

## Environmental Requirements

The system should be operated, stored, or transported within the parameters outlined below. Either its operational environment must be constantly maintained or the unit must be turned off.


|             | Operational                 | Storage                      | Transport (<16 hrs.)          |
|-------------|-----------------------------|------------------------------|-------------------------------|
| Temperature | 10° - 40° C<br>50° - 104° F | -10° - 60° C<br>14° - 140° F | -40° - 60° C<br>-40° - 140° F |
| Humidity    | 30-85%<br>non-condensing    | 30-90%<br>non-condensing     | 30-90%<br>non-condensing      |
| Pressure    | 700-1060hPa                 | 700-1060hPa                  | 700-1060hPa                   |

Table 3–1. System Environmental Requirements

## Console Overview

### Console graphics


The following are illustrations of the console:



- 1 Optional Probe Cable Holder (wire holder standard for Americas' systems)
- 2 Probe and Gel Bottle Holder (Removable for Cleaning)
- 3 Optional Cable Clipper for Probe Cable

Figure 3-1. LOGIQ™ 500 System (right and left side views)

## Console graphics (cont'd)



- 1 Task Light Switch
- 2 Task Light
- 3 Optional B/W Video Page Printer
- 4 Air Filter (front and back of system)
- 5 Power Supply Air Filter
- 6 Optional Probe Cable Holder (wire holder standard for Americas' systems)
- 7 VCR Microphone
- 8 Release Button—to raise and lower video monitor
- 9 Optional Physiological Input Panel
- 10 Swivel Lock on Video Monitor Arm
- 11 Hinged Peripheral Cable Access Panel Door
- 12 Cable Access Channel

Figure 3-2. LOGIQ™ 500 System (front and back views)





## Console Overview

---


### External disk drive

The external MOD drive is a 3.5 inch 128MB/230MB MOD drive located below the keyboard.

It can be used to perform software upgrades, image archiving (option) and service diagnostics.

### Storage areas

Several convenient storage areas are provided within the console as shown by the shaded areas in Figure 3–3. Use them to store gel, options, probe cables, accessories, etc.



1 Storage

Figure 3–3. Storage Areas

## Peripheral/Accessory Connection

### Peripheral/Accessory Connector Panel

LOGIQ™ 500 peripherals and accessories can be properly connected using the rear connector panel located behind the rear door. Only the B/W Page Printer (UP-890) can be connected to the front accessory panel.

Located on the rear panel are video input and output, audio input and output, camera expose, foot switch, power and control connections for VCR, printer, MIC and service tools. See Figure 3–5.

#### CAUTION



Each outer (case) ground line of peripheral/accessory connectors are **Earth Grounded**.

Signal ground lines are **Not Isolated**, except the Service Port.

All of the signal lines (including the signal ground) of the Service Port are **Isolated**.



Figure 3–4. Service Port

#### CAUTION



Use only approved probes, peripherals or accessories.



Console Overview

Peripheral/Accessory Connector Panel (cont'd)





Figure 3-5. Peripheral/Accessory Connector Panel


## Foot Switch (option)



Use only the recommended optional multi-functional Foot Switch in parallel with or as an alternative to the **Freeze** and **Record 1** controls to:

- Freeze a real-time image (left switch).
- Send an image to the hard copy device (right switch).

The Foot Switch connection is located at the back of the console on the left-hand side of the back panel.



- 1 Accessory Panel for Foot Switch Connector
- 2 Store Foot Switch Here
- 3 Console Air Filter Screen
- 4 Power Supply Air Filter

Figure 3–6. Foot Switch Storage and Connectors

## Other Peripherals/Accessories (options)

Refer to the *Recording Images* chapter of this manual for more information.



# System Positioning/Transporting

## Moving the System


When moving or transporting the system, follow the precautions below to ensure the maximum safety for personnel, the system, and other equipment.

Before moving the system:

1. Turn the system power switch off.
2. Unplug the power cord.
3. All cables from off-board peripheral devices (IIE camera, external printer, VCR, etc.) must be disconnected from the console.
4. Ensure that no loose items are left on the console.
5. Loop the cord around the handle on the back of the system or wrap the cord in a bundle and store it behind the rear storage panel.

To prevent damage to the Power Cord, **DO NOT** pull excessively on the cord or make sharp bends while wrapping.

## Moving the System (cont'd)



- 1 Handle
- 2 Storage Area Behind Door

Figure 3–7. Location of Storage Area

6. Connect all probes to be used while off site. Ensure that probe cables are out of the way from the wheels and not protruding beyond the console.



**NOTE:** *If more than three (3) probes are intended to be used, store the additional probes securely in the front storage area.*

7. Store all other probes in their original cases or in soft cloth or foam to prevent damage.
8. Store sufficient gel, optical disks, and other essential accessories in the provided space.
9. Adjust the monitor to its lowest position possible. Ensure that the monitor arm is locked in place.
10. Unlock the front wheels.



## System Positioning/Transporting

---

### Moving the System (cont'd)

When moving the system:

1. Always use the rear handle grips to move the system.
2. Take extra care when moving the system long distances and on inclines. Ask for help if necessary.  
Avoid ramps that are steeper than ten degrees to avoid tipping over the system.



**NOTE:** *Wheel chair ramps are usually less than five degrees.*

Utilize additional care and personnel when moving on steep incline ( $>5^\circ$ ) or loading into a vehicle for transport.



**NOTE:** ***DO NOT*** attempt to move the console using any cables or fixtures, such as the probe connectors.

3. Use the brake, located on the bottom of the system in the front, when necessary.
4. Do not let the system strike walls or door frames.
5. Use extra care when crossing door or elevator thresholds.
6. Once the destination is reached, lock the wheels.

### CAUTION



The system weighs approximately 180 kg (397 lbs). To avoid possible injury and equipment damage:

- Be sure the pathway is clear.
- Limit movement to a slow careful walk.
- Use two or more persons to move the system on inclines or long distances.

## Transporting the System

Use extra care when transporting the system using vehicles. In addition to the instructions used when moving the system (refer to 3–10), do the following:

1. Only use vehicles that are designed for transport of the LOGIQ™ 500 system.
2. Load and unload the system to a vehicle parked on a level surface.
3. Ensure that the transporting vehicle can handle the weight of the system plus the passengers.
4. Ensure that the load capacity of the lift (a minimum of 180 kg [397 lbs] is recommended) is capable of handling the weight of the system.
5. Ensure that the lift is in good working order.
6. Secure the system while it is on the lift so that it cannot roll. Use either wood chocks, restraining straps, or other similar types of constraints. Do not attempt to hold it in place by hand.



### WARNING



**NOTE:** Strap the system below its handle so that the system does not break loose.

Never ride on the lift with the system. A person's weight coupled with the weight of the system may exceed the load capacity of the lift.

7. Employ two to three persons to load and unload safely from a vehicle.
8. Load the unit aboard the vehicle carefully and over its center of gravity. Keep the unit still and upright.



**NOTE:** Do not lay the unit down.

9. Ensure that the system is firmly secured while inside the vehicle. Any movement, coupled with the weight of the system, could cause it to break loose.
10. Secure system with straps or as directed otherwise to prevent motion during transport.
11. Prevent vibration damage by driving cautiously. Avoid unpaved roads, excessive speeds, and erratic stops or starts.





## System Positioning/Transporting

---

### Wheels



Examine the wheels frequently for any obvious defects that could cause them to break or bind.

#### Front wheels

The front wheels swivel, pivot, and lock.

#### Back wheels

The back wheels swivel and pivot but do not lock.



*NOTE: For the USA version console, the back wheels do not pivot.*

### Setting the lock

To engage the wheel lock:

Press down on the lock pedal (located at the front of the wheels). The pedal remains depressed.

To release the lock:

Press down on the pedal again. The pedal returns to its normal position.

# Powering On the System

## Connecting and Using the System

To connect the system to the electrical supply:

1. Ensure that the wall outlet is of the appropriate type.
2. Ensure that the power switch is turned off.
3. Unwrap the power cable. Make sure to allow sufficient slack in the cable so that the plug is not pulled out of the wall if the system is moved slightly.
4. Push the power plug securely into the wall outlet.


### WARNING




To avoid risk of fire, the system power must be supplied from a separate, properly rated outlet. See *Site Requirements, Before the system arrives* on 3–3 for rating information.

Under no circumstances should the AC power plug be altered, changed, or adapted to a configuration rated less than specified. Never use an extension cord or adapter plug.

To help assure grounding reliability, connect to a “hospital grade” or “hospital only” grounded power outlet.



120 VAC, 1350 VA  
Plug and Outlet Configuration  
(USA)



220–240 VAC, 1350 VA  
Plug and Outlet Configuration  
(Europe)

Figure 3–8. Example Plug and Outlet Configurations



## Powering On the System

### Acclimation Time

After being transported, the unit requires one hour for each 2.5° increment its temperature is below 10° C or above 40° C.

|       |     |     |     |     |     |    |    |    |    |    |    |
|-------|-----|-----|-----|-----|-----|----|----|----|----|----|----|
| ° C   | 60  | 55  | 50  | 45  | 40  | 35 | 30 | 25 | 20 | 15 | 10 |
| ° F   | 140 | 131 | 122 | 113 | 104 | 95 | 86 | 77 | 68 | 59 | 50 |
| hours | 8   | 6   | 4   | 2   | 0   | 0  | 0  | 0  | 0  | 0  | 0  |

|       |    |    |    |     |     |     |     |     |     |     |
|-------|----|----|----|-----|-----|-----|-----|-----|-----|-----|
| ° C   | 5  | 0  | -5 | -10 | -15 | -20 | -25 | -30 | -35 | -40 |
| ° F   | 41 | 32 | 23 | 14  | 5   | -4  | -13 | -22 | -31 | -40 |
| hours | 2  | 4  | 6  | 8   | 10  | 12  | 14  | 16  | 18  | 20  |

Table 3–2. System Acclimation Time Chart

### Power On/Off Standby

#### CAUTION



Press the top portion of the **Power On/Off Standby** switch to turn the power on. The circuit breaker, on the rear of the unit, must also be in the on position (see 3–21).




Figure 3–9. Power Switch Location

## Power Up Sequence

- The monitor and console power indicator light up.
- The system is initialized. During this time:
  - Two beeps sound during the sequence.
  - All lighted buttons on the keyboard light.
  - System diagnostics run. Its status is reflected on the monitor by the graphics in *Figure 3–10*.




Figure 3–10. Power Up Graphic Sequence



## Powering On the System

---

### Power Up Sequence (cont'd)



**NOTE:** If errors occur, an error message appears at the bottom of the screen. See *User Maintenance, Troubleshooting* for more information.



#### Hints

If problems occur, freeze the image and take a picture for reference. This will help if there is a need to call for service.

- Probes are initialized for immediate operation.



**NOTE:** If no probes are connected, the system goes into standby mode.

- Peripheral devices are activated on power up.

After initialization has been completed, the new patient entry menu will be displayed.

### Password Protection

A maximum of eight personal IDs and associated passwords can be preset in the LOGIQ™ 500 (see 14–61).

If IDs and passwords have been entered and the Password Ask parameter is on, the following occurs in the power up sequence:

- The message USERID:\_\_\_\_\_ appears in the middle of the start up sequence.
- Type in a User ID and press **Return**.
- Keyboard lights go out.
- The message PASSWORD:\_\_\_\_\_ is highlighted.
- Type in the password that corresponds to the User ID and press **Return**.

Password Protection (cont'd)




Figure 3–11. System Startup Screen with Password Ask on

If the correct User ID/Password is entered, the system continues with the power up sequence.

If an incorrect User ID/Password is entered, the system allows four additional attempts. After the fifth incorrect entry, the message "Login incorrect" appears. The power up sequence will not continue. The system must be turned off, then on, to begin again.



*NOTE: If User ID/Passwords are registered and the Password Ask function is off, press **Return** at the ID\_\_\_\_\_ prompt. The system continues the power up sequence.*



## Powering On the System

---

### Power Off

When switching off the system:

- Move the **ON, OFF/STAND-BY** switch to the OFF position. The message “WARNING: NOW STARTING THE POWER OFF PROCESS” appears at the bottom of the display.
- The LOGIQ™ 500 takes a few seconds to save current scan parameter data in the temporary files to the hard drive before turning the power off.

During this time a message flashes on the screen:


***“Do not pull Power Cable.  
Do not turn off Breaker.”***

Pulling the power cable or turning off the circuit breaker, while the hard drive is working, may corrupt the system operating software on the hard drive.

If the system has not turned off five minutes after pressing the power switch off, listen for hard drive activity. If there is no hard disk drive activity, the circuit breaker on the bottom of the power supply can be used to turn off the system. Do NOT turn off the circuit breaker while the hard disk is working.

- Disconnect the probes.  
Clean or sanitize all probes as necessary. Store them in their shipping cases to avoid damage.
- If daily maintenance is to be performed, turn off the circuit breaker in the back of the system.

## Circuit breaker



The Circuit Breaker is located on the back of the console, at the bottom of the system. **On** supplies main power to all internal systems. **Off** removes main power from all internal systems. The circuit breaker automatically shuts off power to the system in case of a power overload.

If a power overload occurs:

1. Turn off all peripheral devices.
2. Turn off the Main Power Switch to the console.
3. Reactivate the Circuit Breaker switch.


The Circuit Breaker switch should stay in the **On** position; **DO NOT** hold the switch in the **On** position. If the Circuit Breaker switch remains **On**, follow the *Power On* procedure previously described.



**NOTE:** If the Circuit Breaker switch does **not** remain in the **On** position or trips again:

1. Disconnect the Power Cable.
2. Call Service immediately.

**DO NOT** attempt to use the system.



1 Circuit Breaker

Figure 3–12. Location of Circuit Breaker



## Adjusting the Display Monitor

### Rotate, tilt, raise and lower the monitor

The monitor position can be adjusted for easy viewing.

- The monitor can be rotated around its central pivot point.
- The monitor can be tilted for the optimum viewing angle.
- The monitor arm can swing forward or backwards.
- The monitor arm can be raised or lowered for the best viewing height.




Figure 3–13. Display Monitor Movement

### CAUTION



Movement of the monitor swing arm or height adjustment requires the release of the locking mechanism. After an adjustment is made, ensure that the mechanism is locked to prevent unexpected motion.



When moving the LOGIQ™ 500 system, lower the monitor to its lowest possible position to improve stability.

## Brightness and Contrast

Adjusting the monitor's contrast and brightness is one of the most important factors for proper image quality. If these controls are set incorrectly, the Gain, TGC, Dynamic Range and even Acoustic Output may have to be changed more often than necessary to compensate.

The proper setup displays a complete gray scale. The lowest level of black should just disappear into the background and the highest white should be bright, but not saturated.

To adjust the Brightness and Contrast:

1. Select the Set Up Top Menu. From the Set Up Sub-Menus, select Utility. The Utility Menu appears as shown in Figure 3-14.




Figure 3-14. Utility Menu



## Brightness and Contrast (cont'd)

2. Enter menu number "11" (Monitor Adjustment) in the Utility Menu and press **Return**. The Monitor Adjustment Menu appears as shown in Figure 3–15.




Figure 3–15. Monitor Adjustment Menu

3. Enter the number of a gray scale map for adjustment according to the brightness of the exam room and the exam category and press **Return**.
  - 1: Dark Room**—for a dimly lit room with a color monitor.
  - 2: Bright Room**—for a well lit room with a color monitor.
  - 3: Cardiac**—for a Cardiac Exam with a color monitor.
  - 4: Dark Room (B/W)**—for a dimly lit room with a black and white monitor.
  - 5: Bright Room (B/W)**—for a well lit room with a black and white monitor.
  - 6: Cardiac (B/W)**—for a Cardiac Exam with a black and white monitor.

## Brightness and Contrast (cont'd)

The adjustment maps' gray scale bar appears on the monitor as shown in Figure 3-16.




Figure 3-16. Brightness and Contrast

4. Access the Brightness/Contrast controls by pushing on the top center of the control panel door, located under the display screen. Push **Brightness** and **Contrast** slidepots to the maximum (to the right).
5. Decrease the **Brightness** by degrees until "B" in the gray scale bar is not visible.
6. Decrease the Contrast by degrees until "A" in the gray scale bar is not visible.
7. Push **Brightness** slidepots to the maximum (to the right) again.
8. Decrease the **Brightness** by degrees until "C" in the gray scale bar is not visible.
9. Press the Utility rocker switch to return to the previous scan mode.



## Adjusting the Monitor

---

### Brightness and Contrast (cont'd)

Generally speaking, do not change the controls once they have been set. Once set, the display then becomes the reference for the hard copy device(s).



*NOTE: After readjusting the monitor's Contrast and Brightness, readjust all preset and peripheral settings.*



*NOTE: Monitor degaussing (demagnetizing) is done automatically when the system is turned on.*

### Speakers

Stereo audio, provided by speakers located on the sides of the Display Monitor, is standard for:

- Audio Doppler operation (left side is blood flow away/right side is blood flow toward)
- Audio playback of videotaped scan sessions
- Audio error notification




Figure 3–17. Display Monitor Speakers



# Probes

## Introduction

Only use approved probes.

All imaging probes can be plugged into any of the three probe ports.

**For more  
information**

Refer to the *Probes* chapter.

## Selecting a probe



### Hints

- Always start out with a probe that provides optimum focal depths and penetration for the patient size and application.
- Begin the scan session using the default Acoustic Output setting for the probe and application.

*NOTE: Selecting a new probe unfreezes the image.*

## Connecting the Probe

Probes can be connected at any time, regardless of whether the console is powered on or off. To ensure that the ports are not active, place the system in the image freeze condition.

To connect a probe:

1. Place the probe's carrying case on a stable surface and open the case.
2. Carefully remove the probe and unwrap the probe cord.
3. **DO NOT** allow the probe head to hang free. Impact to the probe head could result in irreparable damage.
4. Turn the connector locking handle counterclockwise.
5. Align the connector with the probe port and carefully push into place.

**CAUTION**





## Connecting the Probe (cont'd)




Figure 3-18. Probe Connector Panel

6. Turn the connector locking handle clockwise to secure the probe connector.
7. Carefully position the probe cord so it is free to move and is not resting on the floor.




Figure 3-19. Connecting a Probe

## Cable Handling

Take the following precautions with probe cables:

- Keep free from wheels
- Do not bend the cable acutely
- Avoid crossing cables between probes.



Activating the Probe

To activate the probe, press the **Probe Select** key that corresponds to the probe port to which the desired probe is connected. The **Single CWD** key is used to activate the dedicated continuous wave Doppler (CWD) probe.



Figure 3–20. Probe Select Keys

Key Light Indicators:

- Dark—no probe attached to probe port.
- Half Bright—probe attached to probe port but not active.
- Brightly Lit—probe attached to probe port and is active.

The probe activates in the currently selected operating mode. The probe’s default settings for the mode and selected application are used automatically.

CWD Mode, split Crystal or Doppler-Only CW probes are options available on the LOGIQ™ 500.

Probe Name Menu

If the Display Probe Name preset in the Setup/System Parameters page one is ON, by pressing a half bright or brightly lit key causes the probe names for each port to be displayed in the soft menu. The desired probe can then be selected from the soft menu or by pressing the **Probe Select** keys. Press the **Mode Top Menu** key to exit the probe name display.

| B           | Preset      | Set Up | ECG |
|-------------|-------------|--------|-----|
| C364<br>(1) | L739<br>(2) | ???(3) |     |
|             |             |        |     |

Figure 3–21. Probe Name Menu





## Probes

---

### Deactivating the Probe

When deactivating the probe, the probe is automatically placed in standby mode.

To deactivate a probe:



- Press the **Freeze** key.
- Gently wipe the excess gel from the face of the probe.
- Carefully slide the probe around the right side of the keyboard, toward the probe holder.
- Ensure that the probe is placed gently in the probe holder.

### Disconnecting the Probe

Probes can be disconnected at any time. However, the probe should not be selected as the active probe.

- Move the probe locking handle counterclockwise.
- Pull the probe and connector straight out of the probe port.
- Carefully slide the probe and connector away from the probe port and around the right side of the keyboard.
- Ensure the cable is free.
- Be sure that the probe head is clean before placing the probe in its storage box.

### Transporting Probes

Secure the probe in its holder for moving short distances. When transporting a probe a long distance, store it in its carrying case.

### Storing the Probe

It is recommended that all probes be stored in the carrying case provided.

- First place the probe connector into the carrying case.
- Carefully wind the cable into the carrying case.
- Carefully place the probe head into the carrying case. **DO NOT** use excessive force or impact the probe head.

# Operator Controls

## Control Panel Map

Controls are grouped together by function for ease of use. See the callouts for this illustration on the following page.




Figure 3–22. Control Panel



### Control Panel Map (cont'd)

1. Alphanumeric Keyboard (see 3–31).
2. Patient Information (see 4–3).
3. Probe Selection (see 3–27).
4. Soft Menu (see 3–34).
5. Doppler CFM (see 5–28 and 5–48).
6. User Define (see 14–98).
7. Mode Display Record (see 3–46).
8. Measurements/Annotations (see 3–48).
9. TGC and Acoustic Output (see 5–11).
10. VCR (see 3–50).

### Key Illumination

The front panel keys are all illuminated according to their availability.

- **Disabled** or unavailable selections have the key illumination turned off.
- **Enabled** or available selections have the key illumination half lighted (back lit).
- **Active** or cancelable selections have the key illumination fully lighted.



## Keyboard


The standard alphanumeric keyboard has some special functions.



The Cursor Home key brings the alphanumeric cursor to the very upper left corner of the available field.




The Escape key is used to exit or cancel specified functions or modes.




Control is used in conjunction with other keys to activate special keyboard functions.



Back Space is used to delete previous characters while annotating.



Return is used to move to the next line of annotation.



Tab is used to move forward or backwards through the text one word at a time, or eight characters at a time.



Caps Lock locks all alpha characters in the upper case mode.




Red Shift is used to activate the special characters highlighted in red on the keys to the right side of the keyboard. See 6–20 for details.

Red



Blue Shift activates the VCR controls on the keyboard for the approved Sony SVO-9500MD. When **Blue Shift** is activated and the VCR is in play mode, the left/right pointers control searching reverse or forward. The up/down pointers control shifting backwards/forwards one frame at a time while the VCR is paused.

Blue



If the VCR is not in play mode, the left/right pointers will cause the VCR to rewind or fast forward.

## User Define Function

Alphanumeric keys (A~Z and 1~0) can also be used for the User Define Key function. Refer to 14–98 for details.



## Soft Menu Control Panel

Additional functionality, not available as a control or key on the front panel, can be found via the Soft Menus. Different soft menus appear depending on the mode, special function or calculation package selected.

The Soft Menu consists of **Top Menu Select** keys and **Sub-Menu Select** rocker switches.




Figure 3–23. Soft Menu Control Panel

The Soft Menu Display is in the center of the Soft Menu Controls. The display is divided into twelve sections, four top menu categories and eight sub-menu categories.




Figure 3–24. Top Menu Page Display



# Top Menu Organization

Top Menu Select



The Top Menu Select key cycles through the top level menu page selections. The far left side top menu is the default selection and its sub-menus are automatically displayed.

The top menu groups are divided into two pages. The first top menu page displays:

|              |        |        |     |
|--------------|--------|--------|-----|
| Mode Default | Preset | Set Up | ECG |
|--------------|--------|--------|-----|

Figure 3–25. Mode Default Top Menu (page one)

- Mode Default Menu is the current highest priority active mode. The five possibilities are “B”, “M”, “PWD”, “CWD” and “CFM”.
- Preset is for user programmable application and factory exam parameter preset selections.
- Set Up enables the system customization sub-menus.
- ECG (option) is used to adjust the ECG waveform and ECG synchronized scanning.

The second top menu page shows:

|         |       |          |      |
|---------|-------|----------|------|
| Archive | DICOM | Auto Seq | CINE |
|---------|-------|----------|------|

Figure 3–26. Mode Default Top Menu (page two)

- Archive is used with the Image Archive option for storing images to MOD.
- DICOM (option) displays sub-menus for DICOM setup.
- Auto Sequence displays the user programmed measurement sequences.
- CINE is used when the image is frozen for review of the accumulated image data.



## Operator Controls

---

### Top Menu Organization (cont'd)



The four top menu select keys select and highlight the corresponding top menu and display its associated sub-menus.

Four keys on the front panel automatically disable all top menus and display specific sub-menu selections.

- **Body Pattern** displays the available graphic selections in the sub-menu area.
- **Comment** displays the available selections from the comment library.
- **Measurement** displays the available measurements suitable for the current exam category, image format or individual user sequences.
- **Image Recall** displays information pertaining to images temporarily stored in system memory.



*NOTE: To return functions to the Top-Menus, press **Clear**.*

### Sub-Menu Organization


Different Sub-Menus are displayed according to the front panel key pressed or Top Menu selected.

Each selection or parameter in the Sub-Menus relates directly to the Top Menu that is fully illuminated.

*Sub-Menu Select*



The Sub-Menu Select rocker switch cycles to the next (▶) or previous (◀) sub-menu page.



The Sub-Menu rocker switches allow for the increase/decrease of a parameter value (i.e. Dynamic Range) or enabling/disabling of a parameter (i.e. Image Softener).



## Sub-Menu Displays

### B-Mode Top Menu

| B             | Preset    | Set Up       | ECG        |
|---------------|-----------|--------------|------------|
| Dynamic Range | Gray Map  | Focus Number | ATO Create |
| ▲▼<br>48      | ▲▼<br>B-2 | ▲▼<br>2      | ATO On/Off |

| B             | Preset       | Set Up        | ECG     |
|---------------|--------------|---------------|---------|
| Frame Average | Imaging Freq | Image Softner |         |
| ▲▼<br>OFF     | ▲▼<br>3 MHz  | Color         | 3D Mode |

| B            | Preset | Set Up    | ECG        |
|--------------|--------|-----------|------------|
| Focus Positn |        |           | Tag Positn |
| ▲▼           |        | Color Tag | ▲▼         |

| B           | Preset       | Set Up   | ECG          |
|-------------|--------------|----------|--------------|
| Biopsy Zone | Image Rotatn | Rejectn  | Edge Enhance |
| ▲▼<br>SGL   | ▲▼<br>0DEG   | ▲▼<br>40 | ▲▼<br>MID    |

Figure 3–27. B-Mode Sub-Menu (pages 1 through 4)





## Operator Controls

### CFM Top Menu

| CFM     | Preset     | Set Up     | ECG        |
|---------|------------|------------|------------|
| CFM Map | Slant Scan | Diag. Mode | MTI Filter |
| ▼▲ VT-1 | ▼▲ 0       | ▼▲ Map     | ▼▲ LOW     |

| CFM           | Preset       | Set Up          | ECG     |
|---------------|--------------|-----------------|---------|
| Frame Average | Penet.       | Display Thrshld | Capture |
| ▼▲ LOW        | High Resoltn | ▼▲ 53           | ▼▲ 0.5  |

| CFM         | Preset         | Set Up      | ECG        |
|-------------|----------------|-------------|------------|
| Packet Size | Spatial Filter | W.E. Cancel | Tag Positn |
| ▼▲ MID      | ▼▲ OFF         | Color Tag   | ▼▲         |

| CFM    | Preset        | Set Up       | ECG     |
|--------|---------------|--------------|---------|
| ACE    | Noise Blanker | Persist ence | MR-Flow |
| ▼▲ OFF | ▼▲ OFF        | ▼▲ OFF       | 3D Mode |

Figure 3–28. CFM Sub-Menu (pages 1 through 4)

(Sub-Menu Page 4 is the CFM/PDI Enhancement Option)



PWD Top Menu

| PWD           | Preset     | Set Up      | ECG         |
|---------------|------------|-------------|-------------|
| Dynamic Range | Slant Scan | Wall Filter | S.V. Length |
| <br>30        | <br>0      | <br>20.0    | <br>5       |

| PWD         | Preset | Set Up    | ECG        |
|-------------|--------|-----------|------------|
| Sweep Speed | Penet. | Color     | Tag Positn |
| <br>MID     |        | Color Tag | <br>       |

| PWD  | Preset  | Set Up        | ECG        |
|------|---------|---------------|------------|
| HPRF | Rejectn | CFM/PWD Ratio | CFM Shrink |
|      | <br>40  | <br>1/2       |            |

| PWD           | Preset    | Set Up       | ECG |
|---------------|-----------|--------------|-----|
| Realtim Trace | Calc Dir. | Trace Method |     |
| <br>OFF       | <br>Compo | <br>PEAK     |     |

Figure 3–29. PWD Sub-Menu (pages 1 through 4)



## Operator Controls

---

### CWD Top Menu

| CWD           | Preset     | Set Up      | ECG |
|---------------|------------|-------------|-----|
| Dynamic Range | Slant Scan | Wall Filter |     |
| ▲▼<br>20      | ▲▼<br>0    | ▲▼<br>20.0  |     |

| CWD         | Preset | Set Up    | ECG        |
|-------------|--------|-----------|------------|
| Sweep Speed |        | Color     | Tag Positn |
| ▲▼<br>MID   |        | Color Tag | ▲▼         |

| CWD | Preset | Set Up   | ECG |
|-----|--------|----------|-----|
|     |        | Rejectn  |     |
|     |        | ▲▼<br>40 |     |

| CWD           | Preset      | Set Up       | ECG |
|---------------|-------------|--------------|-----|
| Realtim Trace | Calc Dir.   | Trace Method |     |
| ▲▼<br>OFF     | ▲▼<br>Compo | ▲▼<br>PEAK   |     |

Figure 3–30. CWD Sub-Menu (pages 1 through 4)

### M-Mode Top Menu

| M             | Preset    | Set Up   | ECG          |
|---------------|-----------|----------|--------------|
| Dynamic Range | Gray Map  | Rejectn  | Edge Enhance |
| ▲▼<br>48      | ▲▼<br>M-2 | ▲▼<br>40 | ▲▼<br>MID    |

| M           | Preset | Set Up    | ECG        |
|-------------|--------|-----------|------------|
| Sweep Speed |        | Color     | Tag Positn |
| ▲▼<br>MID   |        | Color Tag | ▲▼         |

Figure 3–31. M-Mode Sub-Menu (pages 1 and 2)

### Preset Top Menu

| B             | Preset        | Set Up        | ECG           |
|---------------|---------------|---------------|---------------|
| User Preset–1 | User Preset–3 | User Preset–5 | User Preset–7 |
| User Preset–2 | User Preset–4 | User Preset–6 | User Preset–8 |

| B                | Preset           | Set Up           | ECG              |
|------------------|------------------|------------------|------------------|
| Factory Preset–1 | Factory Preset–3 | Factory Preset–5 | Factory Preset–7 |
| Factory Preset–2 | Factory Preset–4 | Factory Preset–6 | Factory Preset–8 |

| B             | Preset | Set Up | ECG |
|---------------|--------|--------|-----|
| Recall Preset |        |        |     |
|               |        |        |     |

Figure 3–32. Preset Sub-Menu (pages 1 through 3)



## Operator Controls

---

### Set Up Top Menu

| B              | Preset         | Set Up      | ECG         |
|----------------|----------------|-------------|-------------|
| Custom Display | Preset Program | Save Values | User Define |
| System Paramtr |                | Utility     | Diag.       |

Figure 3–33. Set Up Sub-Menu (page 1 of 1)

### ECG Top Menu

| B      | Preset        | Set Up   | ECG      |
|--------|---------------|----------|----------|
| Single | Sync. Selectn | R Delay  | PCG Wave |
| Dual   | Ref. Scan     | ECG Wave | Aux Wave |

| B                 | Preset           | Set Up            | ECG              |
|-------------------|------------------|-------------------|------------------|
| ECG Gain<br>▲▼ 10 | ECG Positn<br>▲▼ | PCG Gain<br>▲▼ 10 | PCG Positn<br>▲▼ |

| B                 | Preset           | Set Up | ECG |
|-------------------|------------------|--------|-----|
| Aux Gain<br>▲▼ 10 | Aux Positn<br>▲▼ |        |     |

Figure 3–34. ECG Sub-Menu (pages 1 through 3)

### Image Archive Option Top Menu

| Archive        | DICOM       | Auto Seq    | CINE     |
|----------------|-------------|-------------|----------|
| Patient Search | Store Image |             | MO Eject |
| Media Search   |             | DEFF Format |          |

Figure 3–35. Image Archive Sub-Menu (page 1 of 1)



### Cine Top Menu

| Archive     | DICOM       | Auto Seq   | CINE        |
|-------------|-------------|------------|-------------|
| Start Frame | Review Loop | Loop Speed | Side Change |
| End Frame   | Multpl CINE | 1/1        | Cine Gauge  |

| Archive       | DICOM | Auto Seq | CINE |
|---------------|-------|----------|------|
| Cine Capture  |       |          |      |
| Capture Frame |       |          |      |

Figure 3–36. Cine Sub-Menu (pages 1 and 2)

### Auto Sequence Top Menu

| Archive | DICOM | Auto Seq | CINE  |
|---------|-------|----------|-------|
| Seq 1   | Seq 3 | Seq 5    | Seq 7 |
| Seq 2   | Seq 4 | Seq 6    | Seq 8 |

Figure 3–37. Auto Sequence Sub-Menu (page 1 of 1)

### Probe Name Menu

| B        | Preset   | Set Up   | ECG |
|----------|----------|----------|-----|
| C364 (1) | L739 (2) | ???? (3) |     |
|          |          |          |     |

Figure 3–38. Probe Name Menu



## Operator Controls

---

### Body Pattern

|                |                |                |                |
|----------------|----------------|----------------|----------------|
|                |                |                |                |
| Pattern name-1 | Pattern name-3 | Pattern name-5 | Pattern name-7 |
| Pattern name-2 | Pattern name-4 | Pattern name-6 | Pattern name-8 |

Figure 3–39. Body Pattern Sub-Menu

### Comment

|              |              |              |              |
|--------------|--------------|--------------|--------------|
|              |              |              |              |
| Annotation 1 | Annotation 3 | Annotation 5 | Annotation 7 |
| Annotation 2 | Annotation 4 | Annotation 6 | Annotation 8 |

Figure 3–40. Comment Library Sub-Menu

### Measurement (GYN calculation menu)

|         |         |         |            |
|---------|---------|---------|------------|
|         |         |         |            |
| Lt Ov-L | Lt Ov-H | Lt Ov-W |            |
| Rt Ov-L | Rt Ov-H | Rt Ov-W | GYN Report |

Figure 3–41. Typical Measurement Sub-Menu

## Image Recall

|              |               |               |                |
|--------------|---------------|---------------|----------------|
|              |               |               |                |
| B<br>08:48   | BM<br>08:51   | BD<br>08:54   | B<br>08:58     |
| C,B<br>08:50 | C,BM<br>08:53 | C,BD<br>08:57 | C,B/B<br>08:59 |

Figure 3–42. Image Recall Sub-Menu



**NOTE:** The number of images that can be saved for recall will depend on the availability of the extended memory option.

## Advanced Cardiac Measurement Option

|             |               |              |                 |
|-------------|---------------|--------------|-----------------|
|             |               |              |                 |
| LV<br>CUBED | LV<br>BULLET  | LV<br>SP–ELP | LV<br>GIBSON    |
| LV<br>TEICH | LV<br>SIMPSON | LV<br>BP–ELP | CARD<br>Report  |
|             |               |              |                 |
| PLAX        | PSAX<br>–MV   | AP–4CH       |                 |
| PSAX<br>–AV | PSAX<br>–PAP  | AP–2CH       | CARD<br>Report  |
|             |               |              |                 |
| D–MV        | D–PV          |              | TRACE<br>AUTO   |
| D–AV        | D–TV          |              | Transf<br>CALCS |

Transf CALCS is only available with Realtime Doppler Calculation option installed

|         |      |      |                |
|---------|------|------|----------------|
|         |      |      |                |
| M–LV/RV | M–AV | M–TV |                |
| M–MV    | M–PV |      | CARD<br>Report |

Figure 3–43. Advanced Cardiac Measurement Option Sub-Menu (pages 1 through 4)






## Operator Controls

---


### Mode, Display and Record

This group of controls provides various functions relating to the display mode, display orientation, image recording/saving, freeze, gain and cine scroll.



The Mode Controls select the desired display mode or combinations of display modes.

During dual display modes the **L** and **R** keys activate the Left or Right displayed image.



The Depth knob controls the image display depth.



The Reverse key toggles the left/right orientation of the scan image.



The Image Memory key stores the current frozen displayed image in system image memory.

Maximum number of B-Mode images is 8.  
Maximum number of Timeline images is 4.

All images are erased when the New Patient key is pressed or there is a loss of system power.



The Image Recall key displays a menu of the images stored in memory. After pressing **Image Recall**, select the desired image from the Soft Menu for display.



Mode, Display and Record (cont'd)

Record 2



Record 1 and Record 2 are used to activate/print the designated recording device (i.e. video page printer, multi-image camera, image archive option).


Record 1



Freeze

The Freeze key is used to stop the acquisition of ultrasound data and freeze the image in system memory.

Pressing **Freeze** a second time continues live image data acquisition.



Controls the gain of the displayed echoes during B-Mode scanning. Controls the gain of the displayed timeline echoes during B/M- and M-Mode. During B/M-Mode, B-Mode gain can be controlled by the Doppler/CFM Angle Control.

The gain value displayed on the monitor is:

|          |             |
|----------|-------------|
| B-Mode   | B-Mode Gain |
| M-Mode   | M-Mode Gain |
| B/M-Mode | M-Mode Gain |

When the image is frozen it controls scrolling forwards and backwards through the cine loop images in temporary memory.



## Operator Controls

---

### Measurement and Annotation

This group of controls performs various functions related to making measurements, annotating and adjusting the image information.



The Comment key enables the image text editor and displays the annotation library soft menu. For more details, refer to 6–15.




The Body Pattern key enables the body pattern Soft Menu and displays the default pattern on the screen. For more details, refer to 6–22.



The M/D Cursor key enables Trackball control of the M-Mode or Doppler cursor line (not angle correction) or the CFM window in real-time Color Flow Mode.



The Scan Area key enables Trackball control of the B-Mode image area size and position in B/W, the CFM window size and position in Color mode and the Zoomed area size and position in Zoom.



The Measurement key is used in all types of basic measurements. It displays the measurement soft menu for the current exam category.



## Measurement and Annotation (cont'd)

*Ellipse*



The Ellipse rocker switch is used to activate the ellipse measurement function after the first distance measurement has been set. It also toggles which cursor is the movable cursor during the ellipse adjustment.

*Set*



The Set key is used for various functions, but is generally used to fix or finish an operation (i.e. to fix a measurement cursor or exit scan area size/position in B-Mode and zoom).

*Clear*



The Clear key is generally used to erase or exit functions such as annotations/comments, measurements and zoom. This will return the system to the basic mode top menu.


*Zoom*



Press the **Zoom** key to activate the zoom function. For more details, refer to 6-2.



This knob is used to control the zoom size with the zoom function in real-time or rotate the probe position indicator in the body pattern function.



The Trackball is used with almost every key function in this group. Trackball control depends on the last key function pressed.

## Operator Controls

---

### VCR Controls

These controls are used in conjunction with the approved optional video cassette recorder.



Press the **External Video** key to allow for the display of an external video signal (i.e. from a VCR) on the system monitor. Pressing **Play** automatically selects external video.

Press the key a second time to return to LOGIQ™ 500 video display.



Press this key to enable or disable the on board microphone to allow for audio to be added to the taping of system video.



Pressing the **Stop/Play** key either places the VCR in the play mode or stops the play or record mode.



Pressing the **Pause/Record** key either places the VCR in the record mode or causes the VCR to pause if it is already in the record or play mode.



*NOTE: See Keyboard on 3–33 for additional VCR controls.*



# Beginning an Exam

## Introduction

Begin an exam by entering new patient information.

The operator should enter as much information as possible, such as:

- Exam category
- Patient name
- Patient ID
- Comments
- Exam Information

The patient's name and ID number is retained with each patient's image and transferred with each image during archiving or hard copy printing.

## Beginning a New Patient



The **New Patient** key should be pressed at the beginning of each patient study. Pressing this key automatically erases all patient data, annotations, measurements, calculations and summary report pages. A preset parameter can be set to prompt the user if they wish to erase all patient data or not. The system also resets every parameter (i.e. Zoom, 3D Mode, Triplex Mode, etc.).



**NOTE:** *Oper ID and Ref MD will not be erased.*

Enter patient data with the alphanumeric keyboard.

The **Soft Menu** automatically defaults to the last exam preset top menu when **New Patient** is pressed. An exam applications preset can be selected from the sub-menu to provide a suitable starting point for system scan parameters.



**NOTE:** *After a specific time interval the **Soft Menu** defaults to the "Mode" **Top Menu**.*

### CAUTION



To avoid patient identification errors, always verify the identification with the patient. Make sure the correct patient identification appears on all screens and hard copy prints.

The Patient Entry Menu appears on the display monitor.

```

Do you really want to change patient? (y/n)      : 
[ PATIENT ENTRY MENU ]
EXAM CATEGORY SELECTION : 1
1: RAD/ABDOMEN          2: OBSTETRICS          3: GYNECOLOGY
4: CARDIOLOGY           5: VASCULAR            6: UROLOGY
7: SMALL PARTS
PT NAME :
PT ID   :
NOTE    :
Oper ID:

EXAM INFORMATION
AGE :   yrs                DOB :   /   /

Ref MD:
COMMENT :

EXIT
    
```

Figure 4–1. Patient Entry Menu



## **Beginning a New Patient (cont'd)**

The first Data Entry field is presented in reversed display.

- **Trackball** to the exam category selection. Select from 7 examination categories: Radiology/Abdominal, Obstetrics, Gynecology, Cardiology, Vascular, Urology and Small Parts. Input the appropriate number. The patient information input menu changes. Information pertinent to the selected exam category appears in the menu.

To navigate through the Patient Entry Menu, use the **Return** key or the **Trackball** to move the cursor.

- Input the Patient Name (29 characters maximum).
- Input Patient ID Number (14 characters maximum).
- Input any desired Note (30 characters maximum).
- Input the desired Operator ID (four digits maximum).



*NOTE: Patient Name, Patient ID, Notes, Oper ID, Ref MD, and Comments are common to all exam category menus.*

*The Patient Name and Patient ID will appear on the image screen. All other information is automatically entered into the appropriate exam report page.*

*Information in the Exam Category patient entry menus is considered necessary for that type of exam. Fill in all information possible.*



## Beginning a New Patient (cont'd)

- The display units of measure for items such as weight or height can be selected from the Set Up/Preset program menu page 10. Choose the priority and unit of measure on this Preset Menu page.
- Input Ref MD (16 characters maximum).
- Input comment field (2 lines of 50 characters each).
- When all patient data entries have been completed, highlight **Exit** and press **Return** or the **New Patient** key.



*NOTE: If patient information needs to be edited or the exam category changed, use the **ID Name** key. Pressing **ID Name** allows for modification of the Patient Entry Menus without erasing accumulated patient images, measurements, annotations, calculations and summary reports.*



*NOTE: Patient age entry information (years, months, weeks, days) is selected in Set Up/Preset Program page 10, "Display Unit Age". Patient height is "Display Unit Height" and Patient weight is "Display Unit Weight". Choose the unit values that are to appear on the patient entry menu.*

## OB/GYN Exam Categories

BBT is a pregnancy origin data selection choice that appears in the Tokyo University, Osaka University and European OB formats only. LMP, EDC and GA are the only choices in the USA version.

If the Multigestational option is installed, "Fetus Number:" will appear to the right of Ref. MD:. If there is more than one fetus, enter the correct number of fetuses. If the Multigestational option is not installed, each fetus must be scanned separately. Re-enter the patient's information for each fetus.


It is possible to read patient data from a PC (personal computer) into the LOGIQ™ 500 Patient Entry Menu. See 9–65 for details.



## Beginning an Exam

---

### ID/Name



The **ID/Name** key functions the same as the **New Patient** key except it does **NOT** erase previous patient data (i.e. measurements, calculations, etc.) Use the **ID/Name** key to enter or replace patient data without changing the current status of the system. One common reason might be to change the exam category.

Pressing **ID/Name** enables the Patient Entry Menu. Edit the menu as described in *Beginning a New Patient*.

No other function can be started until ID/Name is completed. To complete the ID/Name process, press **ID/Name** a second time or move the highlighted cursor to **Exit** and press **Return**.



### Hints

If power is lost during the ID/Name function, any data that was added or modified will not be saved.

Patient information/setup is saved to the system hard drive at power off only if the System Parameters (page 1) preset “Power On Status” is set to *Keep-Latest*. The setup is then returned to its power down state when power is turned on.

The following rules apply when filling in the New Patient menu:

- Press **Caps Lock** to type uppercase letters. Press **Caps Lock** again to type lowercase letters.
- Press **Back Space** to erase characters and correct errors.
- To change information, press **Return** or use the **Trackball** to move to the field, then type over the existing information with correct information.
- Press **Return** to move to the next field.
- Use the **Trackball** to move the reversed cursor to the desired item.
- When pressing **Return** at the last data entry field, the system returns to real-time scanning.
- Standard keyboard keys repeat when held down.
- To start over, press **New Patient**.

Remember, user and factory-defined presets are dependent upon the exam category selected when filling in the New Patient menu.

# Exam Application Preset Selection

## Introduction

Select the exam category application preset, that best describes the desired study to be performed, from the factory default preset selections displayed on the Soft Menu.

A typical example for OB is shown in Figure 4–2.

|               |                |  |  |
|---------------|----------------|--|--|
|               |                |  |  |
| General<br>OB | Fetal<br>Heart |  |  |
| Tech<br>Diff  | Newborn        |  |  |

Figure 4–2. OB Preset Selection

Use these parameters as a starting point for the exam.

## Selecting a probe

- Always start out with a probe that provides optimum focal depths and penetration for the patient size and application.
- Begin the scan session using the default Acoustic Output setting for the probe and application.

The “Change Probe and Category” preset, found on System Parameters page 1, allows the exam category to automatically change depending on the probe selected.

Some probes are used in more than one exam category. If the default for “Change Probe and Category” is set to *Simultaneous* and the wrong exam category is automatically activated, press **ID/Name** to select the correct exam category.



## *Exam Applic. Preset Selection*

---

This page left blank intentionally.



# B-Mode

## Introduction

The LOGIQ™ 500 offers a wide variety of display formats. Each format shows the operator valuable information relating to patient data and system scan parameters.

B-Mode, also referred to as 2-D Mode, displays a two dimensional gray scale image of the structures within the probe's field of view. The anatomy imaged can be evaluated, measurements and calculations performed as well as making an initial survey prior to entering other modes.

## B-Mode Key Operation

Pressing the **B**-Mode key while in a combination of other modes will result in the following:

| Current Mode    | Resultant Mode |
|-----------------|----------------|
| B/PD            | B              |
| B/M             | B              |
| B/CFM and PD    | B/CFM          |
| B/CFM and M/CFM | B/CFM          |
| B/CFM           | B              |
| M/CFM           | B/CFM          |

Table 5–1. B-Mode Key Operation



## B-Mode

### Reading the B-Mode Display

Figure 5–1 shows the basic B-Mode display information.




Figure 5–1. B-Mode Display Format

## Reading the B-Mode Display (cont'd)

|               |              |      |           |
|---------------|--------------|------|-----------|
| Patient Name  | ID           | Date | Time      |
| Hospital Name | LMP=xx/xx/xx | P100 | MHz Probe |

| B-Mode Display                                                                                                | Description, Format, Values                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|---------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Patient Name                                                                                                  | A maximum of 29 alphanumeric characters. Input at Patient Entry Menu.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| ID                                                                                                            | Patient identification number. A maximum of 14 alphanumeric characters. Input at Patient Entry Menu.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Date                                                                                                          | Today's date according to the system settings. Can be preset to display as MM/DD/YR, DD/MM/YR, or YR/MM/DD.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Time                                                                                                          | Displays the current time during normal operation. Can be preset to a 12 or 24 hour clock. Displays the frame acquisition time when in Cine Mode.                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Hospital Name                                                                                                 | A maximum of 29 alphanumeric characters. Input at Set Up/System Parameters, page 1.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| LMP or GA                                                                                                     | In GYN Exam Category: LMP=xx/xx/xx<br>In OB Exam Category: GA (xxx)=##W#D<br>The xxx after the GA shows how gestational age was calculated:<br>LMP, EDD or OPE (operator entered) – These are input in the Patient Entry Menu.                                                                                                                                                                                                                                                                                                                                                                       |
| Acoustic Output Percentage                                                                                    | The percentage is preceded by a "P" and then a number from 0–100.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Frequency of Operation (Imaging Freq)                                                                         | Frequency of Operation (for B- and M-Modes). The image frequency value is displayed in MHz. The Frequency of Operation value coincides with the "B Image Frequency" and "M Image Frequency" presets located on Set Up/Custom Display pages 2 and 3, respectively. The image frequency can range between Normal, Resolution, Penetration, Further Penetration, Harmonic 1 and Harmonic 2. The table on 5–19 shows the correlation between the image frequency optimization setting and the frequency of operation value. When Harmonic is active, the letter "T" is shown before the frequency value. |
| Probe                                                                                                         | Probe name or designation of the active probe.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Scan Orientation<br>GE or  | GE or the marker is used for scan orientation. This should coincide with the probe orientation marking on the probe body. This marker can be turned off in Set Up/Custom Display page 10.                                                                                                                                                                                                                                                                                                                                                                                                            |

Table 5–2. B-Mode Display Explanation




| Graphic Display         | Description, Format, Values                                                                                                                                                                                                                           |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| B Gray Scale Map Number | The B gray scale map number is displayed above the gray scale bar on only B Mode. The B gray scale map number disappears when B-Color is active. When ATO map is active, an asterisk (*) is displayed on the right side of the gray scale map number. |
| Gray Scale              | Shows the B-Mode Gray Scale assignment.                                                                                                                                                                                                               |
| Color Scale             | In Color Mode, the right half shows the Color Scale assignment.                                                                                                                                                                                       |

Table 5–3. B-Mode Display Explanation




## B-Mode

### Reading the B-Mode Display (cont'd)



| Graphic Display | Description, Format, Values                                                                                                                                                                              |
|-----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| FR/Cine         | Shows display acquisition frame rate.<br>###Hz (Real-time) Can be turned off by a parameter in Set Up/Custom Display page 10.<br>In Cine Mode it shows the display playback frame number. CN### (Frozen) |
| Depth           | Shows the display in cm. ##cm<br>Dual Mode: ##cm##                                                                                                                                                       |
| DR              | Dynamic Range shows the range over which echo intensities are converted to gray scale. Displayed in dB. DR##<br>Dual Mode: ##DR##                                                                        |
| Gain            | Displays the overall B-Mode or M-Mode Receive Gain. G##<br>Dual Mode: ##G##                                                                                                                              |
| Body Pattern    | Shows the pattern selected for scan orientation.                                                                                                                                                         |
| VCR Counter     | Indicates the video tape location shown as: ##.##.##                                                                                                                                                     |
| ▶▶              | VCR status indication shows the current operational status of the VCR.                                                                                                                                   |
| TV Counter      | A two digit indication of the TV Frame count during VCR playback.                                                                                                                                        |

Table 5-4. B-Mode Display Explanation



| Graphic Display   | Description, Format, Values                                                                                                                                                                                  |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Measurements      | Lines of image measurement data are displayed in this area.<br>The format and value depends on the type of measurement.<br>The bottom line fills first and scrolls up as additional information is obtained. |
| Operator Messages | System generated messages are displayed on this line (i.e. error messages)                                                                                                                                   |
| HR                | Heart Rate is displayed here in beats per minute. ###BPM (Requires ECG input)                                                                                                                                |
| A Output          | Acoustic Output display for MI (Mechanical Index) or TI (Thermal Index). The larger one is displayed. If both are less than 0.4, the message "TI=<0.4" is displayed.                                         |

Table 5-5. B-Mode Display Explanation



## Acoustic Output Display

The acoustic output display in the lower right corner of the screen is intended to keep the user informed of the potential bioeffects that can result when controls are changed.

The Acoustic Output display shows MI (Mechanical Index) while in B-Mode. When other modes are selected TIS, TIB or TIC are displayed. See Figure 5-2.

LOGIQ Display increments are 0.2 for values of one or less and 1.0 for values greater than one. Display accuracy for MI is -15% to +38%. Display accuracy for TI is -37% to +95%. For all index values, if the calculated value is less than 0.4, the display will be "<0.4".




Figure 5-2. Acoustic Output Display Selection



## *B-Mode*

---

### **Adjusting the Acoustic Output**

The initial acoustic output may be preset within the range of 0% and 100% in Set Up/Custom Display page 1. Acoustic power settings are available for B/M and D/CFM/PDI with the following presets:

- Acoustic Power B/M [%]
- Acoustic Power D/CFM/PDI [%]

While scanning, rotate the dial clockwise to increase the acoustic output, or counterclockwise to decrease it.



Reading the B-Mode Display (cont'd)




Figure 5-3. B-Mode Display Format

Scale  
Marker

(0 cm)  
(1 cm)  
  
(5 cm)

Scale markers are presented along the right side of the display as large marks every 5cm and small marks every 1cm.



The TGC curve displays the relative position of the TGC slide pots compared to their depth.

This display can be turned on or off in the Setup/Custom Display Menu page 10.

10.0

The number displayed next to the probe orientation symbol is the zoom or scroll depth starting point for the image displayed.



## B-Mode

### Optimizing the Image

#### Control Layout




Figure 5–4. Controls

Mode and Display controls are on the lower right side, while Acoustic Output and TGC are on the lower left side.

## Controls

|                                                                                                                                        | Acoustic Output                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | TGC                                                                                                                                                                                           |
|----------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                                                                     | Allows for varying the power output of the transmitted sound wave.<br><br>Follow the principle of “as low as reasonably achievable” (ALARA) when scanning patients. Always use the least amount of acoustic output possible to obtain an optimal image.                                                                                                                                                                                                                                                                                                          | Amplifies returning signals to correct for the attenuation caused by tissue at increasing depths. The individual eight slide pots correspond to the maximum display depths for each probe.    |
| <b>Accessing/Changing</b>                                                                                                              | <i>NOTE: Always optimize gain before increasing the acoustic output.</i><br><br>Clockwise rotation increases the acoustic output level. Counterclockwise rotation decreases the acoustic output level.                                                                                                                                                                                                                                                                                                                                                           | <i>NOTE: Active slide pots backlight.</i>                                                                                                                                                     |
| <b>Benefits</b>                                                                                                                        | Optimizes the image quality thereby minimizing exposure time to the patient while maximizing penetration and echo return.                                                                                                                                                                                                                                                                                                                                                                                                                                        | Balances the image so that the density of echoes is the same throughout the image.                                                                                                            |
| <b>Values</b>                                                                                                                          | The mid-range acoustic output setting is factory preset to be a reasonable setting for each type of exam category.<br><br>Acoustic output levels are returned to the factory preset value when changing the following: <ul style="list-style-type: none"> <li>• Probe</li> <li>• Exam category</li> <li>• New patient</li> </ul>                                                                                                                                                                                                                                 | When display depth is changed, TGC is rescaled to affect the new depth range. Each TGC pot is proportionally scaled across the display depth range as shown in Figure 5-5 on 5-12.            |
| <b>Affects on Other Controls</b>                                                                                                       | Acoustic output and gain interact to affect acoustic output exposure. To a point, gain can be increased to compensate for a reduction in acoustic output. This minimizes acoustic output exposure.                                                                                                                                                                                                                                                                                                                                                               | In Zoom, the TGC adjusts automatically if the “TGC Depth Remap” preset, on Custom Display page 10, is set to ON. Only the TGC pots that affect the zoomed region of interest will be lighted. |
| <b>Bioeffects</b><br> <b>Acoustic Output Hazard</b> | <p>Increasing acoustic output increases Mechanical and Thermal Index values; decreasing acoustic output decreases Mechanical and Thermal Index values.</p> <p>However, the improvement in image quality should allow the patient exposure time to be decreased.</p> <p><b>CAUTION</b>  Observe the effect of acoustic output adjustment on the TI/MI display. Increasing the index values to levels greater than 1.0 should take into account the risk/benefit potential.</p> |                                                                                                                                                                                               |



Controls (cont'd)




Figure 5-5. Time Gain Compensation

## Controls (cont'd)

|                                                                                                                                                | Depth                                                                                                                                                                                                                                                                                                 | B/M Gain                                                                                                                                                                                                                                                                                                                                                    |
|------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                                                                             | Controls the distance over which the B-Mode images anatomy. Display depth may be changed according to the anatomical size or to the region of interest. Minimum/maximum values available depend on the probe. Select from 4, 5, 6, 7, 8, 10, 12, 14, 16, 18, 20, 22, and 24 cm for the desired depth. | Increases or decreases the amount of echo information displayed in an image. It has the effect of brightening or darkening all displayed echoes at any depth.                                                                                                                                                                                               |
| <b>Accessing/<br/>Changing</b>                                                                                                                 | To reduce depth (look at a shallower image), turn the <b>Depth</b> rotary encoder clockwise.<br><br>To increase depth (look at a deeper image), turn the <b>Depth</b> rotary encoder counterclockwise.                                                                                                | <i>NOTE: It is not possible to change the gain on a frozen image.</i><br><br><i>Changing the gain while in another mode does not affect the B-Mode image gain.</i><br><br>The Angle knob controls B-Mode gain in B/M and CFM-B/M-Modes.                                                                                                                     |
| <b>Benefits</b>                                                                                                                                | Depth adjusts the field of view. Increase the field of view to look at larger or deeper structures; decrease the field of view to look at structures near the skin line.                                                                                                                              | Allows for the balance of echo contrast so that cystic structures appear echofree and reflecting tissue fills in.                                                                                                                                                                                                                                           |
| <b>Values</b>                                                                                                                                  | Depth increments vary by probe. Values may be preset for each probe. Depth displays on the monitor in centimeters.                                                                                                                                                                                    | Gain displays on the monitor in dB. Gain increments are available every 2 dB within the range of 0 to 98 dB, depending on the selected probe.                                                                                                                                                                                                               |
| <b>Affects on<br/>Other Controls</b>                                                                                                           | After adjusting the depth, it may be necessary to adjust the TGC, focus, frame rate and edge enhance.<br><br>Changing depth while scanning clears Cine memory, but not Timeline Replay. A bar is displayed in the timeline at which point it is not possible to scroll past this change.              | After adjusting acoustic output, there may be a need to adjust the gain. Generally speaking, if acoustic output increases, the gain may need to decrease; a decrease in acoustic output may require an increase in gain.<br><br>Gain and TGC interact together. Gain changes overall echo amplification while TGC changes amplification at specific depths. |
| <b>Bioeffects</b><br> <b>Acoustic<br/>Output<br/>Hazard</b> | Increasing the depth tends to decrease the MI and TI because the frame rate slows down.                                                                                                                                                                                                               | Gain has no affect on Acoustic Output. However, with increased Gain, the output level can usually be reduced to produce equivalent image quality.                                                                                                                                                                                                           |
|  <b>Hints</b>                                               | Make sure enough space is left below the anatomy of interest to demonstrate shadowing or enhancement.                                                                                                                                                                                                 |                                                                                                                                                                                                                                                                                                                                                             |



## B-Mode

### Controls (cont'd)

|                                                                                                                                        | Scan Area Size                                                                                                                                                                                                                                                                                                                                         | Scan Area Position                                                                                                                                                                                    |
|----------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                                                                     | Widen or narrow the size of the sector to maximize the image's field of view for convex and sector probes. Adjusts color window for linear probes.                                                                                                                                                                                                                                                                                      | The reduced sector angle can be steered to get more information without moving the probe. Colored markers around the image act as guides while steering the reduced sector angle.                                                                                                        |
| <b>Accessing/Changing</b>                                                                                                              | <p>To narrow the angle:</p> <ul style="list-style-type: none"> <li>Press <b>Scan Area</b>.</li> <li>Move the <b>Trackball</b> toward the left.</li> </ul> <p>To widen the angle:</p> <ul style="list-style-type: none"> <li>Press <b>Scan Area</b> (if not already activated).</li> <li>Move the <b>Trackball</b> toward the right.</li> </ul> <p><i>NOTE: Pressing <b>Scan Area</b> toggles between size and steering control.</i></p> | <p>To assign trackball control to the scan area position mode and steer the angle:</p> <ul style="list-style-type: none"> <li>Press <b>Scan Area</b> (adjust scan area size).</li> <li>Press <b>Scan Area</b> again.</li> <li>Move the <b>Trackball</b> to the right or left.</li> </ul> |
| <b>Benefits</b>                                                                                                                        | Increase the sector width to see a wide anatomical structure, decrease the sector width to have a faster frame rate.                                                                                                                                                                                                                                                                                                                    | Allows the movement of a reduced sector angle laterally. Beneficial in GYN.                                                                                                                                                                                                              |
| <b>Values</b>                                                                                                                          | <p>Sector/Convex Probe: Range from 10° to full "B" width.</p> <p>Linear Probe: Range from 10mm to full "B" width.</p>                                                                                                                                                                                                                                                                                                                   | Continuous steer within the full maximum probe width.                                                                                                                                                                                                                                    |
| <b>Affects on Other Controls</b>                                                                                                       | Increasing sector size decreases the frame rate; decreasing sector size increases the frame rate.                                                                                                                                                                                                                                                                                                                                       | In scan area position mode, press <b>Scan Area</b> to return trackball control to the scan area sizing mode.                                                                                                                                                                             |
| <b>Bioeffects</b><br> <b>Acoustic Output Hazard</b> | <p>As the sector width narrows, the TI tends to increase since the target is being hit more often and the MI may decrease since the peak power is reduced.</p> <p>Observe the output display for possible effects.</p>                                                                                                                                                                                                                  |                                                                                                                                                                                                                                                                                          |





## B-Mode

### Controls (cont'd)

|                                  | Reverse                                                                                                                                                                                                                                  | Display Format (Dual)                                                                                                                                                                                                                                                                                                                                            |
|----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Controls the horizontal orientation of the image on the screen. It toggles the left/right orientations of the image display.                                                                                                             | Each scan mode or combinations of modes can be displayed in dual format (side by side images). The <b>Left Image</b> and <b>Right Image</b> keys allow for the dual format to be displayed.                                                                                                                                                                                                                                                        |
| <b>Accessing/ Changing</b>       | Press the <b>Reverse</b> key to toggle image reverse on or off.                                                                                                                                                                          | To initiate the dual format display, press the <b>Left Image</b> key. The current display mode will be reduced to the left half of the image area.<br><br>Press the <b>Right Image</b> key to activate the same image format on the right half of the image area.<br><br><i>NOTE: When the <b>Left Image</b> key and <b>Right Image</b> key are pressed simultaneously, while in single B-Mode imaging, the dual live image mode is activated.</i> |
| <b>Benefits</b>                  | The <b>Reverse</b> key allows for changing the orientation of the image display without physically rotating the probe 180°.                                                                                                              | Allows for viewing two images side-by-side.                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Values</b>                    | On or Off.<br><br>The GE logo at the top of the sector wedge corresponds to the orientation mark on the probe body. Pressing the <b>Reverse</b> key flips the image left/right and switches the GE logo to correspond to the probe mark. | System parameter displays relocated for each format.                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Affects on Other Controls</b> |                                                                                                                                                                                                                                          | Other controls only affect the active image.                                                                                                                                                                                                                                                                                                                                                                                                       |



## B-Mode

### B-Mode Sub-Menu Page 1

| B                                                                                                                                                                         | Preset                                                                                                                                                                     | Set Up                                                                                                                                                                   | ECG        |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| Dynamic Range                                                                                                                                                             | Gray Map                                                                                                                                                                   | Focus Number                                                                                                                                                             | ATO Create |
| <br> 48 | <br> B-2 | <br> 2 | ATO On/Off |

Figure 5–6. B-Mode Sub-Menu Page 1

|                                  | Dynamic Range                                                                                                                                                                                                                                                                                                                                                                                                                                 | Gray Scale Mapping (Gray Map)                                                                                                                                                                                                                                                                              |
|----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Controls how echo intensities are converted to shades of gray.                                                                                                                                                                                                                                                                                                                                                                                | Determines how the echo intensity levels received are presented as shades of gray.                                                                                                                                                                                                                         |
| <b>Benefits</b>                  | Useful for optimizing tissue texture to differentiate between echo levels that are close together.<br><br>Dynamic range should be adjusted so that the highest amplitude edges appear as white while lowest levels (such as blood) are just visible.                                                                                                                                                                                          | Displays the received echo levels with different weights on specific levels of gray. For example, a certain gray map may enhance mid level echoes over a wider range of grays verses high or low level echoes.<br><br>Allows for better differentiation between echo levels through gray levels displayed. |
| <b>Values</b>                    | The settings cycle in 6 dB steps from 30 dB to 90 dB. Use presets to set the output preferred to any scan mode and exam combination. The default for B Dynamic Range is set in the Set Up/Custom Display menu page 2. Dynamic range levels are returned to the preset value when changing the following: <ul style="list-style-type: none"><li>• Application</li><li>• Exam category</li><li>• New patient</li></ul> Pre-processing function. | There are 16 selections of Gray Scale Mapping.<br><br>See Figure 14–6 on 14–16 for a graph of the Gray Scale Map.<br><br>Post-processing function.                                                                                                                                                         |
| <b>Affects on Other Controls</b> | Operates only in real-time, not in freeze, Cine, Timeline replay, or VCR playback.                                                                                                                                                                                                                                                                                                                                                            |                                                                                                                                                                                                                                                                                                            |

## B-Mode Sub-Menu Page 1 (cont'd)

|                                                                                                                                        | Focus Number                                                                                                                                                                                             | ATO Create (software option)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|----------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                                                                     | Changes the number of focal zones so that the beam can be tightened or expanded for a specific area. A graphic caret corresponding to the focal zone position(s) appears on the right edge of the image. | ATO, Automatic Tissue Optimization, optimizes the image based upon a specified Region Of Interest (ROI) or anatomy within the display. ATO gray scale map is available only in B-Mode.<br><br><b>ACCESSING/CHANGING:</b><br>Once ATO Create has been selected, a ROI box cursor is displayed on the image. Use the trackball for changing the size and location of the ROI box cursor. Press <b>Scan Area</b> to toggle between sizing and moving the cursor.<br><br>Press <b>Set</b> . The system automatically optimizes the entire image and the highlighted sub-menu selection changes to ATO ON/OFF. |
| <b>Benefits</b>                                                                                                                        | Optimizes the image by increasing the resolution for a specific area.                                                                                                                                    | Improves productivity and patient throughput with one touch image optimization. Improves contrast of the displayed image data.                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Values</b>                                                                                                                          | Choose from 1, 2, 3 or 4 focal zones.<br>Pre-processing function.                                                                                                                                        | On or off.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Affects on Other Controls</b>                                                                                                       | Changing the focus number affects frame rate. The greater number of focus points, the slower the frame rate.                                                                                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Bioeffects</b><br> <b>Acoustic Output Hazard</b> | Adding focal zones tends to increase TI, although it may decrease MI.<br><br>Observe the output display for possible effects.                                                                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |

|                    | ATO On/Off (software option)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | The ATO On/Off selection is highlighted when the ATO map is active on the image. ATO ON/OFF can be toggled on and off to display that same ATO map again. When ATO map is active, the asterisk (*) is displayed on the right side of the gray scale map number. The B gray scale map number is displayed under the gray scale bar.<br><br>ATO gray scale map is available until changing the following: <ul style="list-style-type: none"> <li>• Probe</li> <li>• New Patient</li> <li>• Exam Category</li> <li>• User Preset</li> <li>• Gray Scale Map</li> <li>• Power Off</li> <li>• New ATO map creation.</li> </ul> Only one ATO gray map is available at a specific time. |
| <b>Benefits</b>    | Ability to access the previously created ATO map without any time delay.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Values</b>      | On or off.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |



## B-Mode

### B-Mode Sub-Menu Page 2

| B             | Preset          | Set Up        | ECG     |
|---------------|-----------------|---------------|---------|
| Frame Average | Imaging Freq    | Image Softner |         |
| ▲<br>▼<br>OFF | ▲<br>▼<br>3 MHz | Color         | 3D Mode |

Figure 5–7. B-Mode Sub-Menu Page 2

|                                  | Frame Averaging                                                                                                                                                                                                                                                                                  | Imaging Frequency                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Averages previous frames of image data with the current frame. Frame averaging uses more data points to make up 1 image. This has the effect of presenting a smoother, softer image.<br><i>NOTE: Not available on Cine replay.</i>                                                               | In B-Mode, imaging frequency can be changed to allow more echoes to pass through.                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Benefits</b>                  | Helps to average out brief, sudden changes in echo intensity information. Could help to filter out low intensity noise.                                                                                                                                                                          | Allows for an increase in echo intensity without changing gain, TGC or acoustic output with lower frequencies. Allows for increased resolution with higher frequencies.                                                                                                                                                                                                                                                                                                       |
| <b>Values</b>                    | Frame average values are off, 1, 2, 3, 4, 5 and 6. Frame average values are returned to the preset value when changing the following: <ul style="list-style-type: none"><li>• Probe</li><li>• Exam category</li><li>• New Patient</li><li>• Application type.</li></ul> Pre-processing function. | Values shown are probe dependent.<br>Table 5–6 shows the correlation between the frequency values displayed and the imaging frequency mode preference (penetration, normal, resolution, further penetration, harmonic 1 and harmonic 2).<br>To set the default image frequency for B-Mode and M-Mode, refer to Custom Display pages 2 and 3, respectively.<br>Further penetration mode is available on the C358 probe only.<br>Refer to 5–20 for details on Tissue Harmonics. |
| <b>Affects on Other Controls</b> |                                                                                                                                                                                                                                                                                                  | Acoustic output, gain and TGC may be able to be reduced due to increased penetration.                                                                                                                                                                                                                                                                                                                                                                                         |

## B-Mode Sub-Menu Page 2 (cont'd)

| Probe      | Frequency of Operation depends on Image Frequency Sub-Menu Selection |        |             |                     |            |            |
|------------|----------------------------------------------------------------------|--------|-------------|---------------------|------------|------------|
|            | Resolution                                                           | Normal | Penetration | Further Penetration | Harmonic 1 | Harmonic 2 |
| B510       | 6.5MHz                                                               | 5MHz   | 5MHz        | N/A                 | N/A        | N/A        |
| S317       | 4MHz                                                                 | 3.5MHz | 2.5MHz      | N/A                 | T5.0       | T4.7       |
| S222       | 3.5MHz                                                               | 3MHz   | 2MHz        | N/A                 | T3.5       | T3.2       |
| S611       | 6MHz                                                                 | 5MHz   | 4MHz        | N/A                 | N/A        | N/A        |
| S220 (W)   | 3MHz                                                                 | 2MHz   | 2MHz        | N/A                 | N/A        | N/A        |
| S316 (UC)  | 3.6MHz                                                               | 3MHz   | 3MHz        | N/A                 | N/A        | N/A        |
| P509       | 6.5MHz                                                               | 5.5MHz | 5MHz        | N/A                 | N/A        | N/A        |
| E721(MTZ)  | 7MHz                                                                 | 6MHz   | 5.5MHz      | N/A                 | N/A        | N/A        |
| C721       | 7MHz                                                                 | 6MHz   | 6MHz        | N/A                 | N/A        | N/A        |
| C551 (CAE) | 6MHz                                                                 | 5MHz   | 4.4MHz      | N/A                 | N/A        | N/A        |
| C364 (CBF) | 4MHz                                                                 | 3.5MHz | 3MHz        | N/A                 | N/A        | N/A        |
| C386       | 5MHz                                                                 | 3.8MHz | 3MHz        | N/A                 | N/A        | N/A        |
| C358       | 5MHz                                                                 | 4MHz   | 3MHz        | 2MHz                | T5.0       | T4.5       |
| 739L       | 9MHz                                                                 | 7MHz   | 6MHz        | N/A                 | N/A        | N/A        |
| LA39       | 11MHz                                                                | 9.6MHz | 8.2MHz      | N/A                 | N/A        | N/A        |
| 546L       | 6.6MHz                                                               | 5.2MHz | 4.6MHz      | N/A                 | N/A        | N/A        |
| LD         | 4MHz                                                                 | 3.5MHz | 3MHz        | N/A                 | N/A        | N/A        |
| L764 (LH)  | 7MHz                                                                 | 7MHz   | 7MHz        | N/A                 | N/A        | N/A        |
| I739       | 8MHz                                                                 | 7MHz   | 6MHz        | N/A                 | N/A        | N/A        |
| T739       | 8MHz                                                                 | 7MHz   | 6MHz        | N/A                 | N/A        | N/A        |
| CWD2       | N/A                                                                  | N/A    | N/A         | N/A                 | N/A        | N/A        |
| CWD5       | N/A                                                                  | N/A    | N/A         | N/A                 | N/A        | N/A        |
| ERB7C      | 7MHz                                                                 | 6MHz   | 5.5MHz      | N/A                 | N/A        | N/A        |
| ERB7L      | 7MHz                                                                 | 6MHz   | 5.5MHz      | N/A                 | N/A        | N/A        |

Table 5-6. Frequency of Operation



## B-Mode

---

### B-Mode Sub-Menu Page 2 (cont'd)

|                    | Image Softener                                                     | Color                                                                                                                                                      |
|--------------------|--------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | B-Mode images may be adjusted for the amount of smoothing applied. | Color allows for enabling B-Mode image colorization.                                                                                                       |
| <b>Benefits</b>    | Provides a smoother homogenous image display.                      | Displaying the gray scale as shades of color may allow for improved differentiation between echo levels.                                                   |
| <b>Values</b>      | On or Off.<br>Pre-processing function.                             | On or Off.<br>The B-Mode color presented is a selection in the Set Up/Custom Display menu page 11. There are 6 color choices.<br>Post-processing function. |

|                    | 3D Mode (option)                                                  |
|--------------------|-------------------------------------------------------------------|
| <b>Description</b> | Allows for 3D Mode option activation from the B-Mode menu.        |
| <b>Benefits</b>    | Allows for quick on/off of 3D Mode function from B-Mode Sub-Menu. |
| <b>Values</b>      | On or Off.                                                        |

### THI (Tissue Harmonic Imaging) (option)

**Description**—Tissue Harmonic Imaging utilizes the natural harmonic characteristic inherent in non-linear sound propagation properties in the human body. This function does not permit the use of any contrast agent. It is only available on the S222, S317 and C358 probes.

**Benefits**—Harmonic imaging diminishes low frequency high amplitude noise and often improves imaging technically difficult patients. Tissue Harmonic imaging is beneficial to contrast resolution and cystic clearing.

**Accessing/Changing**—Access Tissue Harmonic frequency using the Imaging Freq sub-menu selection. See 5–18 for details. To easily recognize a tissue harmonic frequency, the letter “T” precedes the frequency value on the monitor and soft menu.

## B-Mode Sub-Menu Page 3

| B                                                                                                    | Preset | Set Up       | ECG                                                                                                |
|------------------------------------------------------------------------------------------------------|--------|--------------|----------------------------------------------------------------------------------------------------|
| Focus<br>Positn<br> |        |              | Tag<br>Positn<br> |
|                                                                                                      |        | Color<br>Tag |                                                                                                    |

Figure 5–8. B-Mode Sub-Menu Page 3

|                                                                                                                                        | Focus Position                                                                                                                                                                                                                                                                        | Color Tag                                                                                                                                                          |
|----------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                                                                     | Changes the depth at which the selected number of focal zones are optimized. All graphic carets representing focal points will move with a change in focal position.                                                                                                                  | Enables the colorization of a specific gray scale level range. This causes the specified gray levels to be displayed as a predetermined color in the B-Mode image. |
| <b>Benefits</b>                                                                                                                        | Optimizes the image by centering the focal point(s) to the depth of the area of interest.                                                                                                                                                                                             | Allows for the quick recognition of specific gray levels by colorization.                                                                                          |
| <b>Values</b>                                                                                                                          | Relative to depth of the display (FOV). Focal point indicators vary with position change up (shallow) or down (deep).<br>Pre-processing function.                                                                                                                                     | On or Off.<br><br>The color tag size presented is a selection in the Set Up/Custom Display menu page 11.<br>Post-processing function.                              |
| <b>Affects on Other Controls</b>                                                                                                       | Changing the focus position affects frame rate. The deeper the focus, the slower the frame rate.                                                                                                                                                                                      |                                                                                                                                                                    |
| <b>Bioeffects</b><br> <b>Acoustic Output Hazard</b> | Moving the focal zone(s) may affect acoustic output requirements because of concentrating on a specific area of interest. Acoustic output may be decreased with the front panel control if the focal zone is properly placed.<br><br>Observe the output display for possible effects. |                                                                                                                                                                    |

|                    | Tag Position                                                                                                                                                                                                                                                                       |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Allows for the movement of the specified color tag range throughout the gray scale displayed.                                                                                                                                                                                      |
| <b>Benefits</b>    | Allows for setting the color tag to the desired level of gray scale colorization.                                                                                                                                                                                                  |
| <b>Values</b>      | If the Color Tag selection is not enabled when Tag Positn is selected, Color Tag is automatically activated.<br><br>Move the color tag up or down the gray scale. Color and size of the tag are determined in the Set Up/Custom Display menu page 11.<br>Post-processing function. |



## B-Mode

### B-Mode Sub-Menu Page 4

| B                                                                                     | Preset                                                                                 | Set Up                                                                               | ECG                                                                                   |
|---------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Biopsy<br>Zone                                                                        | Image<br>Rotatn                                                                        | Rejectn                                                                              | Edge<br>Enhance                                                                       |
|  SGL |  0DEG |  40 |  MID |

Figure 5–9. B-Mode Sub-Menu Page 4

|                       | Biopsy Zone                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | Image Rotate                                                                                                                                                   |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>    | <p>This selection enables any electronic biopsy guidezone(s) available for the active probe.</p> <p>If the single fixed angle (SGL) was selected as the needle guide type in Set Up/Custom Display page 1, the fixed zone angle is displayed.</p> <p>If the needle guide type selected is multiple angle (MBX), the Biopsy Zone rocker switch can be used to cycle through the angle selections MBX1, MBX2, MBX3 and OFF.</p>                                                                                                                                                                    | <p>Rotates the single real-time or zoomed B-Mode image in 90° increments.</p> <p>When B-Mode is rotated 180°, the M-Mode (if active) is also rotated 180°.</p> |
| <b>Benefits</b>       | <p>The electronic biopsy guidezone shows the expected needle path during insertion.</p> <p><b>DANGER</b>  Failure to match the guidezone displayed to the guide may cause the needle to track a path outside the zone.</p> <p>It is extremely important that when using the adjustable angle biopsy guides, the angle displayed on the screen matches the angle set on the guide, otherwise the needle will not follow the displayed guidezone which could result in repeated biopsies or patient injury.</p> | <p>Orient image display for easy reference.</p>                                                                                                                |
| <b>Safety Concern</b> | <p><b>DANGER</b>  When the biopsy guidezone is displayed, the message:</p> <p>“Confirm BX type of Bracket”</p> <p>is displayed at the bottom of the screen with the angle selected.</p> <p>Ensure that the (M)BX type selected for each probe is the same as the angle selected on the biopsy guide.</p> <p><i>NOTE: For specifics on availability of biopsy guide angles, see 15–33.</i></p>                                                                                                                 |                                                                                                                                                                |



## B-Mode Sub-Menu Page 4 (cont'd)

|                                  | Biopsy Lines                                                                                                                                                                                                                                                                                                                                                                                  | Image Rotate                                                                                                                                                                                                                                             |
|----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Values</b>                    | On (SGL, MBX1, MBX2, MBX3, TV0, TR5, 1~8, GBX) or Off.<br><i>NOTE: Press the <b>Measurement</b> key once to display the integrated biopsy depth cursor and center line while the guidezone is present. Use the <b>Trackball</b> to measure the needle length. Needle length is from the top of needle barrel to the target.</i><br>See Figure 15–12 on 15–32 for measuring the needle length. | The image rotates in 90° increments in a clockwise and counterclockwise direction.<br>Post-processing function.                                                                                                                                          |
| <b>Affects on Other Controls</b> |                                                                                                                                                                                                                                                                                                                                                                                               | <b>CAUTION</b>  When reading a rotated image, be careful to observe the probe orientation to avoid possible confusion over scan direction or right/left image reversal. |

|                                  | Rejection (Rejectn)                                                                                                             | B Edge Enhance                                                                                                                                                                                                                                                                                                                                                        |
|----------------------------------|---------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Allows for the elimination of low level echoes from the display. This is generally used to clear noise out of vessels or cysts. | B edge enhance brings out subtle tissue differences and boundaries by enhancing the gray scale differences corresponding to the edges of structures.                                                                                                                                                                                                                  |
| <b>Benefits</b>                  | Allows for the elimination from the display image of low level echoes caused by noise.                                          | Edge enhance modifies the B-Mode image by accentuating the interfaces between organs or vessels.                                                                                                                                                                                                                                                                      |
| <b>Values</b>                    | Off and 2 through 40 in increments of 2.<br>Post-processing function.                                                           | The four selections are Off, Low, Mid and High. The default for B Edge Enhance is set in the Set Up/Custom Display menu page 2. Edge enhance levels are returned to the preset value when changes are made to the following: <ul style="list-style-type: none"> <li>• Application</li> <li>• Exam category</li> <li>• New patient</li> </ul> Pre-processing function. |
| <b>Affects on Other Controls</b> | Rejection affects real-time imaging, frozen, Cine or VCR playback images.                                                       | Operates in real-time only, not in Freeze, Cine, or VCR playback.                                                                                                                                                                                                                                                                                                     |



## B-Mode

---

### B-Mode Optimization

| B-Mode Optimization |                                                                                                                                                                                                                                                                                                           |                                |                                                                                                                                                                                                                                                                                                                                                    |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Adjustments for...  | Do the Following...                                                                                                                                                                                                                                                                                       | Adjustments for...             | Do the Following...                                                                                                                                                                                                                                                                                                                                |
| Image too grainy    | <ol style="list-style-type: none"><li>1. Increase <b>Dynamic Range</b>.</li><li>2. Increase <b>Frame Average</b>.</li><li>3. Decrease <b>Edge Enhance</b>.</li><li>4. Change <b>Gray Map</b>.</li></ol>                                                                                                   | Image too soft                 | <ol style="list-style-type: none"><li>1. Decrease <b>Dynamic Range</b>.</li><li>2. Increase <b>Edge Enhance</b>.</li><li>3. Decrease <b>Frame Average</b>.</li><li>4. Change <b>Gray Map</b>.</li></ol>                                                                                                                                            |
| Image too noisy     | <ol style="list-style-type: none"><li>1. Decrease <b>B/M GAIN</b>.</li><li>2. Decrease <b>Dynamic Range</b>.</li><li>3. Increase <b>Frame Average</b>.</li><li>4. Increase <b>Edge Enhance</b>.</li></ol>                                                                                                 | Improve Uniformity             | <ol style="list-style-type: none"><li>1. Increase the number of focal zones.</li><li>2. Decrease <b>SCAN AREA SIZE</b>.</li><li>3. Adjust <b>TGC</b> to compensate for attenuation.</li></ol>                                                                                                                                                      |
| Cystic Imaging      | <ol style="list-style-type: none"><li>1. Decrease the <b>B/M GAIN</b>.</li><li>2. Decrease <b>Dynamic Range</b>.</li><li>3. Use <b>SCAN AREA SIZE</b> to reduce image width.</li><li>4. Change <b>Focus Number</b> to increase number of focal zones.</li><li>5. Optimize focal zone placement.</li></ol> | Technically Difficult Patients | <ol style="list-style-type: none"><li>1. Select the proper probe or change <b>Imaging Freq.</b> for the exam (larger patient, lower frequency).</li><li>2. Increase <b>ACOUSTIC OUTPUT</b>, if necessary.</li><li>3. Maintain a lower <b>Dynamic Range</b> (45 to 48).</li><li>4. Decrease <b>SCAN AREA SIZE</b> for faster frame rates.</li></ol> |



# Adding Color

## Gray Scale Color

Each mode selection (except CFM) can display the shades of gray as shades of a color. The system simply maps out the gray levels as hues of a color scale.

The gray scale bar changes color and the image display changes accordingly.

The color displayed is selected in the Set Up/Custom Display menu. The preset “Color Map” is chosen separately for B-Mode, M-Mode and Doppler displays. Details on the color map selection are found in *Customizing Your System*.

Colorized gray scale is an on/off selection found in each mode Top Menu.

In each mode (B, M, PD and CWD), the color selection is found on the Sub-Menu page two.

The **Color** rocker switch toggles colorized gray scale on or off.

## Color Flow Mode

Color Flow Mapping (CFM) may also be referred to as Color Doppler or Color Doppler Imaging. It is used to display color coded blood flow information super imposed on the B-Mode, M-Mode or Doppler displays.



## Adding Color

---

### Activating Color Flow



To activate Color Flow Mapping, press the **CFM** key (color flow). The key backlights while in Color Flow Mode. The color flow image appears after a short pause.

Color Flow Mapping can be used with B-Mode, M-Mode or Doppler Imaging.

In B-Mode, the Color Flow information overlays the B-Mode image.

The **Scan Area** key and **Trackball** now control the size and position of the color flow window.

### M-Mode

In M-Mode, the color flow information overlays the B-Mode image. The color displayed in the M-Mode is that which is seen along the timeline cursor.

The Color Map selected for B-Mode applies to M-Mode. The size and position of the color flow window in B-Mode determines the color information displayed in M-Mode.

### Exiting Color Flow



To exit Color, press **CFM** (color flow).

## Reading the Color Flow Display

The following information displays on the Color Flow Mode image:




Figure 5–10. Color Flow Mode Graphic Display

| Graphic Display        | Description, Format, Values                                                                                                                                                                                          |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Velocity Scale         | Raises/lowers the velocity scale on the color bar. Displayed in cm/m per second or Hz.                                                                                                                               |
| Color Scale            | Displays the selected velocity, velocity & variance, power, or velocity & power map.                                                                                                                                 |
| CFM Area Cursor        | Color displays only in this sector. Press the <b>Scan Area</b> key to toggle between cursor size and position control. Size this window using the <b>Trackball</b> . Position this window via the <b>Trackball</b> . |
| MTI Filter             | Displays on the Color Flow bar as a black area surrounding the baseline.                                                                                                                                             |
| Color Threshold Marker | The color display threshold based on the B-Mode gray scale level.                                                                                                                                                    |
| Units (V or F)         | Indicates the unit of measure. Velocity (V) in meters per second or Frequency (F) in KHz.                                                                                                                            |
| Color Gain (CG)        | Displays the Color Receive Gain value from 0 to 31.<br>CG## Dual Mode: ##CG##                                                                                                                                        |

Figure 5–11. Color Flow Mode Display Format

## Adding Color

### Optimizing the Color Flow Image

#### Control Layout




Figure 5-12. Color Doppler Controls

Color Flow Mapping is basically Doppler velocity/direction information mapped as a color on top of the gray scale B-Mode (or M-Mode) Image.

Keyboard controls used to affect the Doppler spectrum also affect CFM Mode. However, some Soft Menu controls affect only the color processing and display.

#### Common Controls

##### Description

After initial transmission/reception, CFM processing is generally separate from gray scale processing.

Acoustic output affects the transmit power for both B-Mode and Color Flow Mapping signals.

TGC, Depth, B/M Gain and Reverse function only on the B-Mode image displayed.

The **M/D Cursor** key is not relevant to CFM scanning; it only affects the Doppler/M-Mode.

Dual format keys (L-R) work the same as in dual B-Mode, but display two B-Mode and CFM images on the left and right side of the screen.




*NOTE: If "L" and "R" keys are pressed simultaneously, both left and right B-Mode images will be active. Pressing CFM will then provide a live B-Mode image on the left and a live B-Mode with CFM on the right.*



## Adding Color

### Select



If the color flow mode is enabled while in Doppler, the **Select** key switches the Doppler/CFM controls between Doppler and CFM.

### Controls

|                                  | Color Doppler Gain                                                                                                                                                                                                                           | CFM/Spectrum Invert                                                                                                                                                                                                                                           |
|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Amplifies the overall strength of the echoes being processed in the Color Doppler section.                                                                                                                                                   | Allows for viewing blood flow from a different perspective, i.e., red away (negative velocities) and blue toward (positive velocities). Color reverses a real-time or frozen image.<br><br><i>NOTE: Invert reverses the forward/reverse color assignment.</i> |
| <b>Accessing/Changing</b>        | Gain values change depending on the probe; gain is not associated with a particular position of the knob.<br><br>To increase gain, turn the <b>Gain</b> knob clockwise.<br><br>To decrease gain, turn the <b>Gain</b> knob counterclockwise. | To reverse the color assignment, press <b>CFM/Spect Invert</b> .                                                                                                                                                                                              |
| <b>Benefits</b>                  | Allows control of the amount of color within a vessel.                                                                                                                                                                                       | Allows viewing of blood flow according to personal preference, or display proper color without changing probe orientation.                                                                                                                                    |
| <b>Values</b>                    | The settings range is 32 levels, but values are not displayed on the image monitor. Usually set the dial in the middle, depending on the vessel.                                                                                             | Invert and non-invert.                                                                                                                                                                                                                                        |
| <b>Affects on Other Controls</b> | The system automatically adjusts the gain based on acoustic output.                                                                                                                                                                          |                                                                                                                                                                                                                                                               |
| <b>Bioeffects</b>                | Color Doppler gain has no affect on acoustic output. However, with increased color Doppler gain, the output level can usually be reduced to produce an equivalent spectrum image quality.                                                    | When “ <i>Spec. Inv. in Triplex</i> ” preset on Custom Display page 14 is set to Simultaneous, the Doppler spectrum is also inverted.                                                                                                                         |



## Adding Color

### Controls (cont'd)

|                                                                                                                                        | Velocity Scale                                                                                                                                                                                                                                          | Color Flow Baseline Shift                                                                                                                                                                                                                                                                                           |
|----------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                                                                     | Increases or decreases the Velocity Scale represented by the color bar.                                                                                                                                                                                 | Minimizes aliasing by reallocating the forward/reverse color velocity scale assignment. For example, it allows a greater portion of the color scale to be assigned to forward flow than to reverse flow.<br><br>Baseline Shift adjusts the point in the color spectrum where the color velocity is at zero (black). |
| <b>Accessing/Changing</b>                                                                                                              | To increase velocity scale, press the top of the <b>Velocity Scale</b> rocker switch until reaching the desired scale.<br><br>To decrease velocity scale, press the bottom of the <b>Velocity Scale</b> rocker switch until reaching the desired scale. | To adjust the baseline shift, press the top or bottom of the <b>Baseline Shift</b> rocker switch.<br><br>Baseline does not wrap. If the baseline is positioned at the top of the color bar, shifting the baseline again will cause the system to beep.                                                              |
| <b>Benefits</b>                                                                                                                        | If blood is flowing at a higher velocity, increase the velocity scale in order to avoid aliasing in the color display.                                                                                                                                  | Minimize aliasing during color flow imaging.                                                                                                                                                                                                                                                                        |
| <b>Values</b>                                                                                                                          | Velocity scale is displayed in cm per second.                                                                                                                                                                                                           | Color flow baseline can be shifted to one of 7 proportionately-spaced positions on the color bar.<br><br>Zero velocity follows the baseline. The total velocity range remains the same.                                                                                                                             |
| <b>Affects on Other Controls</b>                                                                                                       | Velocity scale changes frame rate, acoustic output, and wall filter.<br><br>When adjusting velocity scale, Cine memory is cleared.                                                                                                                      |                                                                                                                                                                                                                                                                                                                     |
| <b>Bioeffects</b><br> <b>Acoustic Output Hazard</b> | Changing the velocity scale tends to change the MI.<br><br>Observe the output display for possible effects.                                                                                                                                             |                                                                                                                                                                                                                                                                                                                     |



## Controls (cont'd)

|                                                                                                                                        | Color Flow Window Size (Scan Area)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|----------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                                                                     | <p>Color Flow data is placed on top of the gray scale image in a specific area. This is represented by a sector wedge outline with convex probes or a rectangle with linear probes.</p> <p>The size and position of this color area can be easily adjusted.</p> <p>The Color Flow window is a solid line when position may be adjusted by the Trackball. The Color Flow window is a dotted line when the size (CFM area) may be adjusted by the Trackball.</p>                                                                                                                                                                                                                                                                                                                                 |
| <b>Accessing/Changing</b>                                                                                                              | <p>While in the CFM Mode, press <b>Scan Area</b> once. This assigns trackball control to color window size adjustment.</p> <p>Use the <b>Trackball</b> to adjust the color window size.</p> <ul style="list-style-type: none"> <li>• Move the <b>Trackball</b> left to close the color window and right to open it.</li> <li>• Move the <b>Trackball</b> down to increase and up to decrease the color window vertical size.</li> </ul> <p>Press <b>Scan Area</b> a second time to assign trackball control to positioning the re-sized window.</p> <p><i>NOTE: Pressing <b>Scan Area</b> toggles between trackball control of window size or position.</i></p> <p>Pressing <b>Set</b> fixes the color window size and position. Trackball control is returned to its previous assignment.</p> |
| <b>Values</b>                                                                                                                          | <p><b>Sector/Convex Probe.</b> Ranges from 10° to full “B” width.</p> <p><b>Linear Probe.</b> Ranges from 10mm to full “B” width.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Bioeffects</b><br> <b>Acoustic Output Hazard</b> | <p>Narrowing the image size may increase the frame rate, thereby increasing the TI. This change also tends to decrease the MI.</p> <p>Observe the output display for possible effects.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |



## Adding Color

### CFM Mode Sub-Menu Page 1

| CFM        | Preset     | Set Up     | ECG        |
|------------|------------|------------|------------|
| CFM Map    | Slant Scan | Diag. Mode | MTI Filter |
| ▲▼<br>VT-1 | ▲▼<br>0    | ▲▼<br>Map  | ▲▼<br>LOW  |

Figure 5–13. CFM Mode Sub-Menu Page 1

|                    | Color Flow Maps                                                                                                                                                                                                                                                                           | Slant Scan                                                                                                  |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Allows for the selection of how Doppler velocities are mapped as color over the gray scale. The Color Bar represents the selected map.                                                                                                                                                    | Slant scan in CFM Imaging is used to control the position of the CFM window for <b>LINEAR PROBES ONLY</b> . |
| <b>Values</b>      | The 6 color map selections available are assigned in the Set Up/Custom Display page 15. The parameters CFM Color Map 1–6 can be assigned as Velocity maps (V1 thru V8), Velocity/Turbulence maps (VT1 thru VT8), Turbulence maps (T1 and T2) or Advanced Velocity Mode Maps (A1 thru A4). | + (left), 0 (center) and – (right). See Figure 5–14.<br>Pre-processing function.                            |
| <b>Benefits</b>    |                                                                                                                                                                                                                                                                                           | Provides a color flow window suitable for linear probe operation.                                           |




Figure 5–14. Slant Scan Window Selections

CFM Mode Sub-Menu Page 1 (cont'd)

|                                  | Diag Mode                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | Color Flow MTI Filter                                                                              |
|----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|
| <b>Description</b>               | <p>Allows for the selection of the best color flow display format. The 4 methods are:</p> <div>  <p><b>Survey Mode (Svy)</b>—The Survey Mode is the default Diag Mode display and is common for each of the 2 programs. The Survey Mode displays the B/W B-Mode image and has an adjustable, full depth CF wedge that can be steered with the <b>Trackball</b>.</p> </div> <div>  <p><b>Survey/Detail Mode (SvyDtl)</b>—The Survey/Detail Mode displays a high resolution color flow window that is half the depth of the B-Mode image. The CF wedge is adjustable horizontally and vertically and can be steered within the full B-Mode image with the <b>Trackball</b>.</p> </div> <div>  <p><b>Map Mode (Map)</b>—In the Map Mode, the B-Mode image and the color flow are the same size. The Map Mode offers higher CF frame rates and a full size adjustable CF wedge.</p> </div> <div>  <p><b>Map/Detail Mode (MapDtl)</b>—The Detail Mode shows a high resolution, high frame rate color flow window that is half of the B-Mode depth. The B-Mode and the CF are the same adjustable angle.</p> </div> | <p>Filters out low flow velocity color. This minimizes motion artifacts caused from breathing.</p> |
| <b>Benefits</b>                  | Gives the operator a choice of CFM displays to optimize scanning.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | Filters out the lowest velocities (colors).                                                        |
| <b>Values</b>                    | <p>4 sequences are available in the Set Up/Custom Display menu page 15:</p> <ul style="list-style-type: none"> <li>• Survey&gt;Map&gt;Map/Detail</li> <li>• Survey&gt;Survey/Detail</li> <li>• Survey&gt;Survey/Detail&gt;Map</li> <li>• Survey/Detail&gt;Map/Detail</li> </ul> <p>Pre-processing function.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | <p>Values are Low, M1, M2, M3, M4 and High.</p> <p>Pre-processing function.</p>                    |
| <b>Affects on Other Controls</b> | Display frame rate could be improved by minimized gray scale B-Mode image size or CFM window size.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |                                                                                                    |



## Adding Color

### CFM Mode Sub-Menu Page 2

| CFM           | Preset       | Set Up          | ECG           |
|---------------|--------------|-----------------|---------------|
| Frame Average | Penet.       | Display Thrshld | Capture       |
| ▲<br>▼<br>LOW | High Resoltn | ▲<br>▼<br>53    | ▲<br>▼<br>0.5 |

Figure 5–15. CFM Mode Sub-Menu Page 2

|                                  | Frame Average                                                                                                                                                          | Penetration (Penet.)                                                                                                                                                                                                       |
|----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Averages color information from previous frames with the current frame.                                                                                                | Penetration can be increased by lowering the operating frequency of the active probe.<br><br>For reference, Color Flow Doppler frequencies (with penetration On or Off) are shown according to probe in Table 5–7 on 5–43. |
| <b>Benefits</b>                  | Averages out 1 time events like noise spikes.                                                                                                                          | Allows for a slight increase in penetration without changing gain or acoustic output.                                                                                                                                      |
| <b>Values</b>                    | Off, low, medium and high.<br>Pre-processing function.                                                                                                                 | On or Off.<br>Pre-processing function.                                                                                                                                                                                     |
| <b>Affects on Other Controls</b> | Trade off between frame rate and color quality. As the color quality increases, the frame rate decreases and as the frame rate increases, the color quality decreases. | Acoustic output, gain and TGC may be able to be reduced due to increased penetration.                                                                                                                                      |
| <b>Bioeffects</b>                |                                                                                                                                                                        | Decreasing the probe transmitting frequency, without decreasing acoustic output, tends to increase Mechanical and Thermal Index values. However, improved penetration should allow for decreasing patient exposure time.   |

CFM Mode Sub-Menu Page 2 (cont'd)

|                                  | High Resoltn (Resolution)                                                                                                                            | Color Flow Display Threshold                                                                                                                                                                                                                                                                                                                                                    |
|----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Provides a quick way to maximize resolution for the CFM display.                                                                                     | Color flow display threshold is the gray scale level at which the overlay of color information stops.<br><br>If the display threshold is set at a gray shade percentage of 34, any color information available will be overlaid below 34% of the gray scale map. Color information will not be placed on top of any part of the gray scale greater than 34% of maximum (white). |
| <b>Benefits</b>                  | Optimizes CFM color display by changing system parameters to maximize color line density.                                                            | Can limit color flow overlay to low level echoes inside vessel walls. Helps minimize color "bleeding" outside vessel walls.                                                                                                                                                                                                                                                     |
| <b>Values</b>                    | When High Resoltn is highlighted in the sub-menu display, the function is on. When it is not highlighted, it is off.<br><br>Pre-processing function. | Post processing function. The color threshold can be set on real-time or frozen images.<br><br>The settings cycle through the following values:<br><br>0% to 100% of the gray scale in 3% increments.                                                                                                                                                                           |
| <b>Affects on Other Controls</b> | High Resolution will cause the scan frame rate to decrease.                                                                                          |                                                                                                                                                                                                                                                                                                                                                                                 |

|                    | Color Capture (option)                                                                                                                                                                                                                                                                                     |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Displays the highest mean velocity (average velocity) detected over a specific time interval.<br><br><i>NOTE: If a velocity map is selected, the image displays accumulated peak mean velocities. If a variance map is selected, the image displays variance for the accumulated peak mean velocities.</i> |
| <b>Benefits</b>    | Useful in making slide presentations. Also refer to the 6–13 for post processing Cine Capture and Capture Frame functions.                                                                                                                                                                                 |
| <b>Values</b>      | Off, 0.5 seconds, 1 second, 2 seconds and ECG.<br><br>Pre-processing function.                                                                                                                                                                                                                             |



## Adding Color

### CFM Mode Sub-Menu Page 3

| CFM         | Preset         | Set Up      | ECG        |
|-------------|----------------|-------------|------------|
| Packet Size | Spatial Filter | W.E. Cancel | Tag Positn |
| ▲▼<br>MID   | ▲▼<br>OFF      | Color Tag   | ▲▼         |

Figure 5–16. CFM Mode Sub-Menu Page 3

|                                                                                                                                        | Packet Size                                                                                                                                                                                                                                                                                  | Spatial Filter                                                                                                                                              |
|----------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                                                                     | Controls the number of samples gathered for a single color flow vector.                                                                                                                                                                                                                      | Smooths color information so that it is less grainy. It averages the color information out over time.                                                       |
| <b>Benefits</b>                                                                                                                        | Allows for the improvement of the color sensitivity.                                                                                                                                                                                                                                         | Smooths color information.                                                                                                                                  |
| <b>Values</b>                                                                                                                          | Soft Menu values are SMALL, MID and LARGE.                                                                                                                                                                                                                                                   | Off, Low, Mid and High.                                                                                                                                     |
| <b>Affects on Other Controls</b>                                                                                                       | When the packet size is decreased, the frame rate is increased at the expense of CFM image quality. When the packet size is increased, CFM image quality is improved at the expense of frame rate.<br><i>Note: This function does not work in Color M-Mode or B/CFM with Pulsed Doppler.</i> | Affects frame rate. The slower the spatial filter averaging, the faster the frame rate; the higher the spatial filter averaging, the slower the frame rate. |
| <b>Bioeffects</b><br> <b>Acoustic Output Hazard</b> | Changing packet size may change the TI and/or MI.<br>Observe the output display for possible effects.                                                                                                                                                                                        | Changing Spatial Average tends to decrease MI.<br>Observe the output display for possible effects.                                                          |

CFM Mode Sub-Menu Page 3 (cont'd)

|                                  | W. E. (Wall Echo) Cancel                                                                               | Color Flow Velocity Tag (Color Tag)                                                                                                                                                                                                                                        |
|----------------------------------|--------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Eliminates the low velocity echoes caused by the motion of vessel walls.                               | Allows for the assignment of a single color to a range of velocities.<br><br>Real-time or frozen images may be tagged. Velocity tags can be incorporated into velocity/variance color maps.<br><br><i>NOTE: The tag assignment appears on the color bar for reference.</i> |
| <b>Benefits</b>                  | Eliminate low velocity echoes caused by vessel wall motion to help clean up the color Doppler display. | Provides the ability to emphasize specific blood flow velocities so that they stand out to the eye.                                                                                                                                                                        |
| <b>Values</b>                    | On or Off.<br>Pre-processing function.                                                                 | On/Off.<br>Post-processing function.<br><br><i>NOTE: The "CFM Tag Width" is a parameter that can be preset in the Set Up/Custom Display menu page 15.</i>                                                                                                                  |
| <b>Affects on Other Controls</b> |                                                                                                        |                                                                                                                                                                                                                                                                            |

|                    | Tag Position (Positn)                                                                                                                                                                                                                                                                 |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Allows for the movement of the specified color tag range throughout the gray scale displayed.                                                                                                                                                                                         |
| <b>Benefits</b>    | Allows for setting the color tag to the desired level of gray scale colorization.                                                                                                                                                                                                     |
| <b>Values</b>      | If the Color Tag selection is not enabled when Tag Positn is selected, Color Tag is automatically activated.<br><br>Move the color tag up or down the gray scale. Color and size of the tag is determined in the Set Up/Custom Display menu page 15.<br><br>Post-processing function. |



## CFM Mode Sub-Menu Page 4

| CFM    | Preset        | Set Up      | ECG     |
|--------|---------------|-------------|---------|
| ACE    | Noise Blanker | Persistence | MR-Flow |
| ▲<br>▼ | ▲<br>▼        | ▲<br>▼      | 3D Mode |
| OFF    | OFF           | OFF         |         |

Figure 5–17. CFM Sub-Menu Page 4

An option is available to enhance CFM/PDI signal processing. It requires the necessary hardware (Color Processing Board) and software (CFM/PDI Enhancement). The following choices are added to the CFM Mode Soft-Menu:

- Adaptive Color Enhancement (ACE)
- Noise Blanker
- Persistence

|                                  | ACE (Adaptive Color Enhancement)                                                                                                                                                                                                                                                                                                                | Noise Blanker                                                                                                                                  |
|----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Designed to reduce “color artifact noise”. This noise can appear as a flash of color synchronized with the heart beat while scanning the abdomen. The noise may also appear to fill the color window when changing window size/position or when moving the probe. Noise due to valve motion in cardiac studies may be considered clutter noise. | Noise Blanker is intended to reduce the random noise generated due to the increase in Color Gain.                                              |
| <b>Benefits</b>                  | Reduces the “color artifact noise” as described in the description section.                                                                                                                                                                                                                                                                     | Noise Rejection increases color sensitivity with the ability to use higher color gains. Very effective when observing low velocity blood flow. |
| <b>Values</b>                    | On or Off.                                                                                                                                                                                                                                                                                                                                      | OFF    No Noise Blanking<br>LOW    Minimum Noise Blanking<br>MID    Medium Noise Blanking<br>HIGH    Maximum Noise Blanking                    |
| <b>Affects on Other Controls</b> | May have a tendency to reduce frame rate due to increased CFM processor functions.                                                                                                                                                                                                                                                              | Allows use of increased Color Gain or equivalent results with less gain.                                                                       |



CFM Mode Sub-Menu Page 4 (cont'd)

|                                  | Persistence                                                                                                                                                                                                                                                                                                                            | MR-Flow                                                                                                                                                                                                                                                                                                               |
|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Retains the largest pixel color value until a larger value is detected or preset time has elapsed.<br><br>When a higher color value for a pixel is detected, that value is retained and the persistence time is reset for that pixel. When the persistence time elapses for that pixel, color is deleted and it returns to gray scale. | Changes the number of focal zones in the CFM area to 3. A graphic caret corresponding to the focal zone position(s) appears on the right edge of the image. Allows for MR-Flow (ACE-2 option) to be turned on or off from the Color Flow Menu. Available only for LA39, C551 (CAE), C358, E721, 739L and 546L probes. |
| <b>Benefits</b>                  | Present a continuous, high sensitivity color image.                                                                                                                                                                                                                                                                                    | MR-Flow optimizes the image by increasing the resolution for the color flow window. Set to 3 focal zones to be concentrated in the CFM area.<br><br>Allows for quick on/off of MR-Flow function from CFM Sub-Menu. Requires ACE-2 option.                                                                             |
| <b>Values</b>                    | OFF No Persistence<br>SHORT Short Persistence time<br>MID Medium Persistence time<br>HIGH Long Persistence time                                                                                                                                                                                                                        | On and Off.                                                                                                                                                                                                                                                                                                           |
| <b>Affects on Other Controls</b> | The activation of the Persistence function will over ride any previous selections made for CFM Capture or CFM B Frame Average. When Persistence is deactivated CFM Capture and CFM B Frame Average will return to normal operation.                                                                                                    |                                                                                                                                                                                                                                                                                                                       |

|                    | 3D Mode                                                        |
|--------------------|----------------------------------------------------------------|
| <b>Description</b> | Allows for 3D Mode activation from the Color Flow menu.        |
| <b>Benefits</b>    | Allows for quick on/off of 3D Mode function from CFM Sub-Menu. |
| <b>Values</b>      | On and Off.                                                    |



## Adding Color

### CFM Optimization

| Color Flow Mode Optimization |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |                                  |                                                                                                                                                                                                                                                                                                                                                          |
|------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Adjustments for...           | Do the Following...                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | Adjustments for...               | Do the Following...                                                                                                                                                                                                                                                                                                                                      |
| Decrease Motion Artifact     | <ol style="list-style-type: none"><li>1. Increase the <b>VELOCITY SCALE</b>.</li><li>2. Increase <b>W.E. Cancel</b>.</li></ol>                                                                                                                                                                                                                                                                                                                                                                        | Eliminate Aliasing               | <ol style="list-style-type: none"><li>1. Increase <b>VELOCITY SCALE</b>.</li><li>2. Lower <b>BASELINE SHIFT</b>.</li></ol>                                                                                                                                                                                                                               |
| Improve Sensitivity          | <ol style="list-style-type: none"><li>1. Increase <b>GAIN</b>.</li><li>2. Decrease <b>VELOCITY SCALE</b>.</li><li>3. Increase <b>ACOUSTIC OUTPUT</b>.</li><li>4. Switch <b>Penet.</b> on and <b>High Resoltn</b> on.</li><li>5. Decrease <b>W.E. Cancel</b>.</li><li>6. Increase <b>Frame Average</b> (this allows gain to increase).</li><li>7. Increase <b>Packet Size</b>.</li><li>8. Reduce <b>SCAN AREA SIZE</b> to smallest size reasonable.</li><li>9. Optimize focal zone position.</li></ol> | Improve Frame Rate/Flow Dynamics | <ol style="list-style-type: none"><li>1. Decrease color window using <b>SCAN AREA SIZE</b> width.</li><li>2. Exit <b>COLOR FLOW</b> to reduce B-Mode image width (use <b>SCAN AREA SIZE</b>).</li><li>3. Decrease <b>DEPTH</b>.</li><li>4. Increase <b>VELOCITY SCALE</b> (if flow state allows it).</li><li>5. Decrease <b>Frame Average</b>.</li></ol> |
| Decrease Color Flow Bleeding | <ol style="list-style-type: none"><li>1. Decrease <b>GAIN</b>.</li><li>2. Increase <b>VELOCITY SCALE</b>.</li><li>3. Decrease <b>Display Thrshld.</b></li></ol>                                                                                                                                                                                                                                                                                                                                       |                                  |                                                                                                                                                                                                                                                                                                                                                          |

## Power Doppler Imaging (option)

|                                                                                                  | Power Doppler Imaging (option)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                               | Power Doppler Imaging, also known as PDI, (option) is a color flow mapping technique used to map the strength of the Doppler signal rather than the frequency shift or velocity of the signal. PDI is a different way of processing the power information from the color flow transmissions by using different color maps and not displaying directional information. Using this technique, the system plots color based on the number of reflectors that are moving, regardless of velocity.                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Accessing/Changing</b>                                                                        | While scanning in the CFM Mode, enable the Blue Shift key (key will be lighted) and press P.<br><br><i>NOTE: For USA systems, the User Define Key 6 has been locked and programmed to toggle the PDI mode on/off. User Define Key 6 can be reprogrammed for other functions if it is unlocked.</i><br><br>Although remaining in the CFM Mode Top Menu, the system will shift to the PDI presets and color maps.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Values</b>                                                                                    | Only power maps are available in the PDI mode.<br>P1 thru P4 are non-topo (power) maps.<br>P5 and P6 are Directional PDI Maps (option).<br>P7 and P8 are topo maps.<br><br>Changes between non-topo and topo maps can only be done in real-time, not on a frozen image.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Affects on Other Controls</b>                                                                 | Only changes between P1 thru P4, P5 and P6, or P7 and P8 are post-processing functions and do not affect image acquisition controls. PDI parameters need to be changed in real-time. Some controls that affect CFM also affect PDI.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Optimizing the PDI Option</b>                                                                 | Custom Display pages 8 and 16 have parameters specific to PDI.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|  <b>Hints</b> | <p>The PDI Velocity choice is a trade-off between low flow detectability and wall motion.</p> <p>The PDI B High Resolution choice is a trade-off between frame rate and color resolution.</p> <p>The PDI B Frame Average choice is a trade-off between a sharper noisy appearance or a smoother less noisy appearance to the image.</p> <p>The activation of the Persistence function will override any previous selections made for PDI B Frame Average. When Persistence is deactivated PDI B Frame Average will return to normal operation.</p> <p>The PDI B Wall Filter choice is made depending on the frequency of the probe.</p> <p><i>NOTE: Remember that parameters found on Custom Display page 8 are probe dependent. Set up these parameters for all probes that will be used for PDI.</i></p> <p><i>NOTE: Parameters on Custom Display page 16 are effective for all probes.</i></p> |



# Doppler

## Introduction

Time zero (the start of the trace) appears on the left side of the graph. As time progresses, the trace moves to the right. The baseline of the graph (representing zero velocity, zero frequency shift or no detected flow) appears as a solid line running horizontally across the display. By convention, movement toward the probe is positive frequencies or velocities which appear above the baseline. Movement away from the probe is negative frequencies or velocities which appear below the baseline.

Typically, blood flow is not uniform but is composed of a mix of blood cells moving at different velocities and in different directions. Thus, the display is composed of a spectrum of gray scale values. Strong signals are displayed as bright shades of gray while weak signals are displayed as darker shades.

Information about the Doppler display is automatically written on the screen and updated when scanning parameters are changed.

The following can be used:

- Pulsed Wave (PW) Doppler
- Continuous Wave (CW) Doppler
- CW Non-Imaging Doppler

## Pulsed Wave Doppler

PW Doppler is typically used for displaying the speed, direction, and spectral content of blood flow at selected anatomical sites.

PW Doppler operates in two different modes: conventional PW and High Pulse Repetition Frequency (HPRF).

PW Doppler can be combined with B-Mode for rapidly selecting the anatomical site for PW Doppler examinations. This is described under B/M-Mode operation. The site where PW Doppler data is derived appears graphically on the B-Mode image (Sample Volume Gate). The site can be moved anywhere within the B-Mode image.



## Frequencies Used

For reference, Doppler frequencies (with penetration On or Off) are shown according to probe in the table below.

| PROBE    | FREQUENCY (MHz) |     |
|----------|-----------------|-----|
|          | ON              | OFF |
| B510     | 4.0             | 5.0 |
| C358     | 2.5             | 3.3 |
| C364/CBF | 2.5             | 3.3 |
| C386     | 2.5             | 3.3 |
| C551/CAE | 4.0             | 5.0 |
| C721     | 5.0             | 6.6 |
| E721/MTZ | 5.0             | 6.6 |
| ERB7     | 5.0             | 6.6 |
| I739     | 5.0             | 6.6 |
| 546L     | 4.0             | 5.0 |
| 739L     | 5.0             | 6.6 |
| LA39     | 5.0             | 6.6 |
| L764/LH  | 5.0             | 6.6 |
| P509     | 4.0             | 5.0 |
| S316/UC  | 2.5             | 2.8 |
| S317     | 2.5             | 2.8 |
| S220/W   | 2.2             | 2.5 |
| S222     | 2.0             | 2.5 |
| S611     | 4.0             | 5.0 |
| T739     | 5.0             | 6.6 |

Table 5–7. Doppler/CFM Frequency

## High PRF

High Pulse Repetition Frequency (HPRF) is a special operating mode of PW Doppler. In conventional mode, a single energy pulse is used to obtain signals for each line in the PW spectrum. In HPRF mode, multiple energy pulses are used. This allows higher velocities to be detected without aliasing artifacts. HPRF mode is used when detected velocities exceed the processing capabilities of the currently selected PW Doppler scale or when the selected anatomical site is too deep for the selected PW Doppler scale.



## Doppler

---

### High PRF (cont'd)

HPRF is enabled when operating in PW Doppler Mode and the velocity scale is increased beyond a set value. The selection on page 3 of the Doppler Sub-Menu can be used to override this automatic switching to HPRF mode.

When HPRF is active, 3 sample volume gates appear along the Doppler Mode cursor. Doppler information can be received from any of the 3 sample volume gates. The Doppler signals from all 3 gates are added together and displayed in 1 spectrum.



*NOTE: Ensure that only 1 gate overlays a blood vessel at a time. Otherwise, signals from more than 1 flow area are superimposed.*

Refer to 5–59 on how to enable/disable the HPRF function.

### Continuous Wave Doppler

Allows examination of blood flow data all along the Doppler Mode cursor rather than from any specific depth. Gather samples along the entire Doppler beam for rapid scanning of heart. Range gated CW allows information to be gathered at higher velocity.



*NOTE: If split crystal is activated, CW is automatically selected.*

There are two CW Doppler operating modes: Steerable and Non-Imaging.

**Steerable**—Allows viewing of the B-Mode image to position the Doppler cursor to the area of interest while viewing the Doppler spectrum (shown below the B-Mode image) and listening to the Doppler audio signal.

Works with sector probes only.

**Non-Imaging**—Provides only Doppler Spectrum and Audio for ascending/descending aortic arch, other hard-to-get-to spaces or higher velocities.

Requires a single CWD probe and connection.

Works with CWD probes only (CWD2, CWD5).



Reading the Doppler Display



*NOTE: Timeline formats (M-Mode and Doppler) available for display are enabled on the Set Up/Custom Display screen menus.*

*First, choose the Side/Side or Top/Bottom display style on Custom Display page 12, "Timeline Format". The two styles cannot be mixed.*

*Second, make the Doppler format Enable/Disable selections on Custom Display page 13.*

The following information displays on the Doppler Mode image:




Figure 5–18. 1/2 B-Mode plus 1/2 D-Mode Display

The following additional parameters are displayed when in Doppler Mode.



| Graphic Display                | Description, Format, Values                                                                                                                                                                                                                                                                                                                                                                       |
|--------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Doppler Frequency (D)          | The value is displayed in MHz (MegaHertz).                                                                                                                                                                                                                                                                                                                                                        |
| Pulse Repetition Frequency (R) | Displays the frequency in KHz (KiloHertz).                                                                                                                                                                                                                                                                                                                                                        |
| Wall Filter (WV)<br>(WF)       | The velocity value is displayed in m/sec (meters per second) or the frequency value in Hz.<br>Dual Mode: WV###<br>###WV###                                                                                                                                                                                                                                                                        |
| Sample Volume (Length) (SV)    | Displayed in mm.<br>Dual Mode: SV##<br>##SV##                                                                                                                                                                                                                                                                                                                                                     |
| $\theta$ (Angle)               | Indicates the Doppler Angle in $\Theta$ (degrees) between the Doppler Mode cursor and angle correction indicator.<br>Dual Mode: $\theta$ ##<br>## $\theta$ ##                                                                                                                                                                                                                                     |
| D.G                            | Doppler signal gain is displayed as a two digit number.                                                                                                                                                                                                                                                                                                                                           |
| HPRF                           | Displayed to indicate that the system is in the high pulse repetition frequency mode.                                                                                                                                                                                                                                                                                                             |
| SPEC INV                       | Displayed to indicate that the spectrum display is inverted from the prescribed norm.                                                                                                                                                                                                                                                                                                             |
| D. Velocity                    | The Doppler velocity or frequency maximum value is displayed in meters per second when “ <i>Display Spectrm Ref. Value</i> ” preset, found on Set Up/Custom Display page 14, is set to “Off”. The forward velocity (frequency) is displayed above the spectrum while the reverse velocity (frequency) is displayed below. When the spectrum is inverted, the velocity scale readings also invert. |
| Time Scale                     | Three speeds:<br>Slow (4 second sweep). Major marker = 1 second; minor marker = 0.5 second.<br>Medium (2 second sweep). Major marker = 1 second; minor marker = 0.5 second.<br>Fast (1 second sweep). Major marker = 1 second; minor marker = 0.5 second.                                                                                                                                         |
| Spectrum Ref. Value            | Doppler velocity or frequency value is displayed as a reference scale when “ <i>Display Spectrm Ref. Value</i> ” preset, found on Set Up/Custom Display page 14, is set to “On”. The scale size varies depending on velocity or frequency.                                                                                                                                                        |

### Table 5–8. Doppler Mode Display Format



## Activating Doppler Mode

### Activating PW Doppler Mode



To activate PW Doppler Mode, press the **PD** key. The key backlights while in PW Doppler Mode. The Doppler sample line and cursor are displayed.

Press **PD** a second time and the Doppler spectrum displays along with the B-Mode image.



*NOTE: Doppler display formats available depend on the presettable parameters found in the Set Up/Custom Display menus.*

Position and size the sample volume to get a spectrum display. Use the system's audio to listen for when the sample volume is positioned over an area of flow.



To exit PW Doppler Mode, press **B**-Mode.



*NOTE: The Set Up/Preset Program allows the types of B/Doppler display formats to be chosen. Each time the **PD** (or **CWD**) key is pressed, it cycles to the next display format enabled in the Set Up/Preset Program.*

### Activating CW Doppler Mode



To activate CW Doppler Mode, press the **CWD** key, as necessary. The key backlights while in CW Doppler Mode.

The Doppler spectrum displays along with the B-Mode image (Steerable CW Doppler, **NOT** Non-Imaging CW Doppler).



To exit CW Doppler Mode, press **B**-Mode.



## Doppler

---

### Doppler Optimization



*NOTE: This section discusses the Pulsed Doppler controls and menu selections. If Continuous Wave Doppler is being used the keyboard controls are the same. Selections found in the CWD soft menu function the same as a like selection in the PWD soft menu.*

### Controls

Doppler Mode controls discussed in this section are shown below:




Figure 5–19. Doppler/CFM Controls

These functions are specific to pulsed or continuous wave Doppler when that mode is active, as well as color flow.

If the color flow mode is enabled while in Doppler, the **Select** key switches the Doppler/CFM controls between Doppler and CFM.



---

## B-Mode Controls

### Description

Several Front Panel controls affect the B-Mode portion of the display and not the echoes in the Doppler spectrum.

These are:

- TGC
- Depth
- B/M Gain
- Scan Area Size
- Scan Area Position
- Reverse
- Dual Format Keys (L/R)



*NOTE: If the scan area size is reduced and the position changed, the Doppler cursor will follow the position change to stay within the displayed scan area.*

One Front Panel control affects the B-Mode portion of the display as well as the echoes in the Doppler spectrum.

It is:

- Acoustic Output

See the *B-Mode* section of this chapter for details on B-Mode controls.



## Doppler

### Controls

|                                      | M/D Cursor                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | Audio Volume                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                   | Assigns trackball control to the B-Mode Doppler cursor. It also displays and erases the theta angle correction cursor.                                                                                                                                                                                                                                                                                                                                                                            | Controls Doppler audio output.<br><i>NOTE: The system automatically unwraps aliased audio.</i>                                                                                                                                                                                                                                                                                          |
| <b>Accessing/<br/>Changing</b>       | Pressing <b>PD</b> will automatically assign trackball control to the B-Mode Doppler cursor.<br>To terminate trackball control of the M/D cursor, press <b>SET</b> .<br>To reactivate trackball control of the Doppler cursor, press <b>M/D Cursor</b> .<br><i>NOTE: Pressing the <b>M/D Cursor</b> key while the trackball is actively controlling the B-Mode Doppler cursor will display the theta angle correction cursor. (Refer to Theta Angle Correction on 5–51 for more information.)</i> | Pressing <b>PD</b> automatically activates Doppler Audio.<br>The volume is initially set to the Audio Volume % chosen in the Set Up/Custom Display page 14.<br><b>CAUTION</b>  Doppler audio sounds change rapidly, often abruptly. Increase the volume in small steps to avoid startling the patient. |
| <b>Benefits</b>                      | Allows for repositioning of the Doppler cursors after the trackball has been assigned a different function.                                                                                                                                                                                                                                                                                                                                                                                       | An audio representation of the flow within a vessel can be used before activating the spectral trace or CWD.                                                                                                                                                                                                                                                                            |
| <b>Values</b>                        | On or Off.<br>The Doppler cursor displays in the center of the screen or changes to the active color if already present.                                                                                                                                                                                                                                                                                                                                                                          | The volume increases/decreases logarithmically, with positive flow through the right channel and negative flow through the left channel.                                                                                                                                                                                                                                                |
| <b>Affects on<br/>Other Controls</b> | Terminates trackball control from its current function and assigns it to the Doppler cursor.<br>When angle correction is active, the M/D cursor returns angle correction to zero degrees.                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                                                                                                                                                                                                                                         |



## Controls (cont'd)

|                                | Doppler Spectral Gain                                                                                                                                                                           | Theta Angle Correction                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>             | Amplifies the overall strength of the echoes being processed in the Doppler spectral trace.                                                                                                     | Estimates the flow velocity in a direction at an angle to the Doppler vector by computing the angle between the Doppler vector and the flow to be measured.<br><br>Estimated flow velocity equals the velocity toward the probe divided by the COS $\theta$ .                                                                                                                                                                                                                                                                                                |
| <b>Accessing/<br/>Changing</b> | Gain (DG#) values are shown on the left side of the display. The initial gain value is selected in the Set Up/Custom Display menu page 4.                                                       | Flow towards the probe is mapped above the baseline of the spectral display and left speaker.<br><br>To adjust the angle clockwise relative to the probe face, turn $\theta$ <b>Angle</b> clockwise.<br><br>To adjust the angle counterclockwise relative to the probe face, turn $\theta$ <b>Angle</b> counterclockwise.<br><br>Press <b>M/D Cursor</b> to quickly return the angle correction to zero degrees and erase the angle correction cursor from the display.<br><br>The initial angle value is selected in the Set Up/Custom Display menu page 4. |
| <b>Benefits</b>                | Fills in or cleans out spectrum information.                                                                                                                                                    | Optimizes the accuracy of the flow velocity. This is especially useful in vascular studies where it is necessary to measure velocity.                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Values</b>                  | 0 to 32 in increments of 2.<br><br>Usually set in the middle of the range.                                                                                                                      | For optimum velocity measurements, the angle of incidence should be less than 20° for cardiac applications. The usual setting is between 45° and 65° for vascular applications.<br><br>1° increments from 0° to 80°.                                                                                                                                                                                                                                                                                                                                         |
| <b>Bioeffects</b>              | Doppler spectral gain has no affect on acoustic output. However, with increased Doppler spectral gain, the output level can usually be reduced to produce an equivalent spectrum image quality. |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |




Figure 5–20. Phased Array Doppler Display with Angle Correct



## Doppler

### Controls (cont'd)

|                                      | CFM/Spectrum Invert                                                                                                                                                                                                                                                                                                                                                                      | Baseline Shift                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|--------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                   | Vertically inverts the spectrum trace without affecting the baseline position.                                                                                                                                                                                                                                                                                                           | Changes the spectrum baseline to accommodate higher velocity blood flow. Minimizes aliasing by displaying a greater range of forward flow than reverse flow. With baseline shift, higher velocities in one direction can be displayed without clipping off the peaks.                                                                                                                                                                                                                          |
| <b>Accessing/<br/>Changing</b>       | Forward flow (F) and Reverse flow (R) signs appear on the velocity scale which reverse when the spectrum is inverted. Positive velocities go to the bottom of the spectrum and to the right speaker. A spectral invert graphic appears on the display to indicate that the spectrum has been inverted.<br><br><i>NOTE: Inverting the spectrum does not affect the baseline position.</i> | Baseline shift adjusts the point in the spectrum where the velocity trace is at zero. The default baseline can be chosen in the Set Up/Custom Display menu page 15.<br><br>The baseline is displayed as a solid line running across the spectrum. The baseline is raised and lowered in equal increments, depending on the current Doppler scale factor.<br><br>The control does not wrap when the maximum baseline shift (in either direction) has been reached. This is a non-repeating key. |
| <b>Benefits</b>                      | If the blood flow is still moving in the same direction after changing the probe angle, the Doppler information may be reversed. It is easy to invert the spectrum instead of reversing the probe orientation.                                                                                                                                                                           | Unwraps the alias. Rearranges the velocity scale display without changing the velocity scale. Readjusts the positive and negative velocities limit without changing the total velocity range.                                                                                                                                                                                                                                                                                                  |
| <b>Values</b>                        |                                                                                                                                                                                                                                                                                                                                                                                          | –75%, –50%, –25%, 0, +25%, +50%, and +75%, with zero being the center of the display, +100% being the top edge of the display and –100% being the bottom edge of the display.                                                                                                                                                                                                                                                                                                                  |
| <b>Affects on<br/>Other Controls</b> | The polarity of the audio Doppler is also reversed.<br><br>When “Spec. Inv. in Triplex” preset on Custom Display page 14 is set to Simultaneous, the CFM spectrum is also inverted.                                                                                                                                                                                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |



## Doppler

### Controls (cont'd)

|                                                                                                                                        | B Pause                                                                                                                                                                                                                                                                                                                                                                                                                                                 | Velocity Scale                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                                                                     | Freezes the B-Mode image while keeping the Doppler spectrum display active.                                                                                                                                                                                                                                                                                                                                                                             | Adjusts the velocity scale to accommodate faster/slower blood flow velocities. Adjust the size of the presentation to show a low velocity scale for low flow and a higher velocity scale for high flow.<br><br>Velocity scale determines PRF. If the velocity goes off the scale (outside the limits of pulsed Doppler), high PRF activates and samples up to 3 gates.<br><br>The display updates velocity scale parameters after adjusting the velocity scale.                                                   |
| <b>Accessing/Changing</b>                                                                                                              | To pause the B-Mode image, press <b>B Pause</b> . The key backlights fully.<br>To cancel B-Mode image freezing, press <b>B Pause</b> . The key semi-backlights.                                                                                                                                                                                                                                                                                         | <i>NOTE: A bar is inserted in the display with each velocity scale change that cannot be measured across velocity scale changes.</i>                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Benefits</b>                                                                                                                        | The spectrum tracks the Doppler Sample Volume Gate in real-time.                                                                                                                                                                                                                                                                                                                                                                                        | Blood flow information is not cut off due to the effect of aliasing.                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Values</b>                                                                                                                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                         | Velocity range is in meters per second or KHz, depending upon preset velocity scale units. Velocity range is dependent on the probe frequency.                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Affects on Other Controls</b>                                                                                                       | B-Mode pause increases the Doppler display quality.<br><i>NOTE: When pressing <b>Freeze</b> while in B Pause and then unfreezing the image, a preset selection will determine if the B Pause key is released or the image remains in B Pause. This is determined by the Set Up/Custom Display menu page 12 "Auto B Melt at Unfreeze".</i><br><br><i>The echo level measurement is not available in B Pause. It is only available on a frozen image.</i> | When raising the velocity scale, the spectrum waveform decreases in size; when lowering the velocity scale, the spectrum waveform increases in size. Changes in the spectrum are relative to changes in the velocity scale, that is, it sizes accordingly.<br><br>When adjusting the velocity scale, Cine memory is cleared.<br><br>Velocity scale units change with sample volume size, depth, frame rate and Doppler wall filter.<br><br><i>NOTE: The velocity scale changes when angle correct is changed.</i> |
| <b>Bioeffects</b><br> <b>Acoustic Output Hazard</b> |                                                                                                                                                                                                                                                                                                                                                                                                                                                         | Adjusting this control may cause minor changes in acoustic output. Observe the output display for possible effects.                                                                                                                                                                                                                                                                                                                                                                                               |



## Doppler

### PW Doppler Mode Sub-Menu Page 1

| PWD                                                                                  | Preset                                                                              | Set Up                                                                                 | ECG                                                                                 |
|--------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Dynamic Range                                                                        | Slant Scan                                                                          | Wall Filter                                                                            | S.V. Length                                                                         |
|  30 |  0 |  20.0 |  5 |

Figure 5–21. PW Doppler Mode Sub-Menu Page 1

|                                  | Dynamic Range                                                                                                                                                                                                                                                                                                                                                               | Slant Scan                                                                                                                                                                                                                                                                                                                                                                                                |
|----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Controls how echo intensities are converted to shades of gray, thereby creating a range of gray scale that can be adjusted. (Also found in the CWD soft menu.)                                                                                                                                                                                                              | Slant Scan in Doppler imaging is used to control the position of the Doppler cursor for <b>LINEAR PROBES ONLY</b> . (Also found in the CWD soft menu.)<br><br>When Slant Scan is selected while in a B-Mode display, the B-Mode is not effected. Only the Doppler cursor will show at a slant.<br><br>The Doppler angle remains constant in relation to the anatomy when slant scan direction is changed. |
| <b>Benefits</b>                  | Useful for optimizing tissue texture to differentiate between echo levels that are close together. Dynamic range should be adjusted so that the highest amplitude edges appear as white while lowest levels (such as blood) are just visible.                                                                                                                               | Provides a Doppler cursor suitable for linear probe operation.                                                                                                                                                                                                                                                                                                                                            |
| <b>Values</b>                    | The settings cycle in 6 dB steps from 18 dB to 48 dB. Set the default value for PWD dynamic range in the Set Up/Custom Display menu page 4.<br><br>Dynamic range levels are returned to the preset value when changing the following: <ul style="list-style-type: none"><li>• Application</li><li>• Exam category</li><li>• New patient.</li></ul> Pre-processing function. | + (Left), 0 (Center) and – (Right). See Figure 5–22 and Figure 5–23.<br><br>Pre-processing function.                                                                                                                                                                                                                                                                                                      |
| <b>Affects on Other Controls</b> | Operates only in real-time, not in freeze, Cine, Timeline replay, or VCR playback.                                                                                                                                                                                                                                                                                          |                                                                                                                                                                                                                                                                                                                                                                                                           |



**PW Doppler Mode Sub-Menu Page 1 (cont'd)**




Figure 5-22. Slant Scan Cursor Selections




Figure 5-23. B-Mode Display with Slant Selection



## Doppler

### PW Doppler Mode Sub-Menu Page 1 (cont'd)

|                                                                                                                                        | Wall Filter                                                                                                                                                                                                                                                                                                                                   | Sample Volume (SV) Length                                                                                                                                                                                                                               |
|----------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                                                                     | Removes the low level, low frequency Doppler signal caused by movement of the vessel walls. (Also found in the CWD soft menu.)                                                                                                                                                                                                                | Changes the size of the sample volume gate length.<br><i>NOTE: Adjustments to the sample volume gate size are made from the center point of the sample volume position. A larger gate will increase the spectral broadening affect to the spectrum.</i> |
| <b>Benefits</b>                                                                                                                        | Eliminates unnecessary information. Filters out low level noise above and below the baseline so it cannot be seen or heard on the spectrum.                                                                                                                                                                                                   | A smaller gate produces accurate sampling results because it is more sensitive. The gate can be enlarged if there are problems hearing the Doppler audio or for sampling large chambers.                                                                |
| <b>Values</b>                                                                                                                          | The value can be displayed as velocity or frequency (Setup/Custom Display page 13). Wall Filter cutoff values are selected as cm/s in Setup/Custom Display page 4.<br><br>1.0, 2.0, 3.0, 5.0, 7.0, 10, 15 or 20 cm/sec. The values could go as high as 50 cm/s depending on the probe and the velocity scale.<br><br>Pre-processing function. | 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14 and 16 mm.<br><br>Pre-processing function.                                                                                                                                                                        |
| <b>Affects on Other Controls</b>                                                                                                       | Wall filter changes with velocity scale.                                                                                                                                                                                                                                                                                                      |                                                                                                                                                                                                                                                         |
| <b>Bioeffects</b><br> <b>Acoustic Output Hazard</b> |                                                                                                                                                                                                                                                                                                                                               | Adjusting this control may cause minor changes in acoustic output. Observe the output display for possible effects.                                                                                                                                     |

## PW Doppler Mode Sub-Menu Page 2

| PWD           | Preset | Set Up    | ECG        |
|---------------|--------|-----------|------------|
| Sweep Speed   | Penet. | Color     | Tag Positn |
| ▲<br>▼<br>MID |        | Color Tag | ▲<br>▼     |

Figure 5–24. PW Doppler Mode Sub-Menu Page 2

|                                  | Sweep Speed                                                                                                                                                                      | Penetration (Penet.)                                                                                                                                                                                                     |
|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Changes the speed at which the timeline is updated. (Also found in the CWD soft menu.)<br><i>NOTE: Time or distance measurements are not allowed across sweep speed changes.</i> | Penetration can be increased by lowering the operating frequency of the active probe.<br>For reference, Doppler frequencies (with penetration On or Off) are shown according to probe in Table 5–7 on 5–43.              |
| <b>Benefits</b>                  | Speed up or slow down the spectrum for analysis to record more or fewer occurrences over time.                                                                                   | Allows for a slight increase in penetration without changing gain or acoustic output.                                                                                                                                    |
| <b>Values</b>                    | Slow (16 sec), Mid (8 sec), Fast (4 sec) and Very-Fast (2sec).<br>Very-Fast is a function in Doppler mode only.<br>Pre-processing function.                                      | On or Off.<br>Pre-processing function.                                                                                                                                                                                   |
| <b>Affects on Other Controls</b> | If more cardiac cycles are seen, the spectrum appears smaller; if less are seen, the spectrum appears larger or more spread out.                                                 | Acoustic output, gain and TGC may be able to be reduced due to increased penetration.                                                                                                                                    |
| <b>Bioeffects</b>                |                                                                                                                                                                                  | Decreasing the probe transmitting frequency, without decreasing acoustic output, tends to increase Mechanical and Thermal Index values. However, improved penetration should allow for decreasing patient exposure time. |



## Doppler

---

### PW Doppler Mode Sub-Menu Page 2 (cont'd)

|                    | Color                                                                                                                               | Color Tag                                                                                                                                                                                                                         |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Color allows for enabling pulsed wave Doppler Mode image colorization. (Also found in the CWD soft menu.)                           | Color tag enables the colorization of a specific gray scale level range. This causes the specified gray levels to be displayed as a predetermined color in the pulsed wave Doppler Mode image. (Also found in the CWD soft menu.) |
| <b>Benefits</b>    | Displaying the gray scale as shades of color may allow for improved differentiation between echo levels.                            | Allows for the quick recognition of specific gray levels by colorization.                                                                                                                                                         |
| <b>Values</b>      | On or Off.<br>The PW Doppler color presented is a selection in the Set Up/Custom Display menu page 13.<br>Post-processing function. | On or Off.<br>The color tag size presented is determined by the Set Up/Custom Display menu page 13.<br>Post-processing function.                                                                                                  |

|                    | Tag Position (Positn)                                                                                                                                                                                                                                                           |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Tag position allows for the movement of the specified color tag range throughout the gray scale displayed. (Also found in the CWD soft menu.)<br><br>If the Color Tag selection is not enabled when <b>Tag Positn</b> is selected, <b>Color Tag</b> is automatically activated. |
| <b>Benefits</b>    | Allows for setting the color tag to the desired level of gray scale colorization.                                                                                                                                                                                               |
| <b>Values</b>      | Move the color tag up or down the gray scale. Color and size of the tag is determined in the Set Up/Custom Display menu page 13.<br>Post-processing function.                                                                                                                   |

## PW Doppler Mode Sub-Menu Page 3

| PWD  | Preset   | Set Up        | ECG        |
|------|----------|---------------|------------|
| HPRF | Rejectn  | CFM/PWD Ratio | CFM Shrink |
|      | ▲▼<br>40 | ▲▼<br>1/2     |            |

Figure 5–25. PW Doppler Mode Sub-Menu Page 3

|                                  | HPRF                                                                                                                         | Rejection (Rejectn)                                                                                                                                                          |
|----------------------------------|------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | HPRF allows the operator to enable or disable the High Pulse Repetition Frequency function while in PWD Mode. Refer to 5–43. | Rejection allows for the elimination of low level echoes from the display. This is generally used to clear noise out of vessels or cysts. (Also found in the CWD soft menu.) |
| <b>Benefits</b>                  | Allows for higher frequencies to be detected with aliasing artifacts.                                                        | Allows for the elimination from the display image of low level echoes caused by noise.                                                                                       |
| <b>Values</b>                    | HPRF is highlighted when enabled (on).<br>HPRF is not highlighted when disabled (off).                                       | Off and 4 through 40 in increments of 4.<br>Post-processing function.                                                                                                        |
| <b>Affects on Other Controls</b> | Frame rate (image update) may decrease when HPRF is active.                                                                  | Rejection affects real-time imaging, frozen, Cine, or VCR playback images.                                                                                                   |

|                                  | CFM/PWD Ratio                                                                                                                                                   | CFM Shrink                                                                                                                                                                         |
|----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | CFM/PWD Ratio is active in triplex mode. It is used to set the velocity ratio between PWD and CFM.                                                              | Reduces the CFM window to the size specified in Set Up/Custom Display page 4.<br><br>CFM Shrink can be turned on/off in Setup/Custom Display page 4.                               |
| <b>Benefits</b>                  | Used to optimize frame rate in Triplex mode. Without changing CFM or PDI velocity scale, PWD velocity scale can be changed to decrease aliasing.                | Used to optimize frame rate in Triplex mode.                                                                                                                                       |
| <b>Values</b>                    | Choose from 1/1, 1/2 or 1/4.<br>1/1 CFM and PWD Velocity the same.<br>1/2 PWD Velocity twice the CFM Velocity.<br>1/4 PWD Velocity four times the CFM Velocity. | When On is selected, CFM Window size is set by the value selected in Set Up/Custom Display page 4.<br><br>When Off is selected, CFM window size is adjusted by scan area function. |
| <b>Affects on Other Controls</b> | Affects display frame rate.                                                                                                                                     | Affects display frame rate.                                                                                                                                                        |



## Doppler

### PW Doppler Mode Sub-Menu Page 4

| PWD              | Preset          | Set Up          | ECG |
|------------------|-----------------|-----------------|-----|
| Realtim<br>Trace | Calc<br>Dir.    | Trace<br>Method |     |
| ▲<br>▼<br>OFF    | ▲<br>▼<br>Compo | ▲<br>▼<br>PEAK  |     |

Figure 5–26. PW Doppler Mode Sub-Menu Page 4

|                                                                                                  | Realtime Trace                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | Calc Dir.                                                                                                                                                                                   |
|--------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                               | Automatically traces (real-time) the parameter selected in the Set Up/Custom Display Menu page 13. The choices are Peak, Floor, Mean, and Mode. (Also found in the CWD soft menu.)                                                                                                                                                                                                                                                                                                           | Choose which part of the Doppler Trace will be used for automatic measurements/calculations. (Also found in the CWD soft menu.)                                                             |
| <b>Benefits</b>                                                                                  | No need to manually trace a selected value. Measurements and calculations will be automatic.                                                                                                                                                                                                                                                                                                                                                                                                 | Easy access to change which portion of the Doppler Trace is used for automatic measurements/calculations.                                                                                   |
| <b>Values</b>                                                                                    | <p>OFF RealTime Trace is turned off.</p> <p>ON RealTime Trace is active.</p> <p><i>NOTE: If the Realtime Doppler Calculation option is not installed, the CALC selection will not be available.</i></p> <p>CALC RealTime Trace is active and measurements/calculations are displayed.</p> <p>Trace and calculations are done internally, even if RealTime Trace is turned off. After the image is frozen, change the RealTime Trace to ON or CALC to display the trace and calculations.</p> | <p>Compo Composite—Uses both forward and reverse flow data.</p> <p>Fored Forward—Trace only data towards the transducer.</p> <p>Revrs Reverse—Trace only data away from the transducer.</p> |
|  <b>Hints</b> | <p>The best Doppler data is collected when parallel to flow, with orientation also parallel to the anatomic target; whereas, the best B-Mode image data is collected perpendicular to the anatomic target. Therefore, there usually is not both an ideal B-Mode image and ideal Doppler data simultaneously.</p>                                                                                                                                                                             |                                                                                                                                                                                             |

## PW Doppler Mode Sub-Menu Page 4 (cont'd)

|                    | Trace Method                                                                                                                                                                                                                                    |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Choose which method is used to trace the Doppler waveform in realtime. (Also found in the CWD soft menu.)                                                                                                                                       |
| <b>Benefits</b>    | Front panel access through the VFD to change trace method.                                                                                                                                                                                      |
| <b>Values</b>      | Select the default parameter to be traced in the Set Up/Custom Display menu page 13.<br>Peak    Peak or highest values<br>Floor    Bottom or lowest values and peak<br>Mean    Average values and peak<br>Mode    Brightest echo trace and peak |

## Doppler Mode Optimization

| Doppler Mode Optimization |                                                                                                                                                                                                                                                                                                                                                                                                                                                     |                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Adjustments for...        | Do the Following...                                                                                                                                                                                                                                                                                                                                                                                                                                 | Adjustments for...         | Do the Following...                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Increase Sensitivity      | <ol style="list-style-type: none"> <li>1. Increase <b>GAIN</b>.</li> <li>2. Increase <b>ACOUSTIC OUTPUT</b>.</li> <li>3. Decrease <b>VELOCITY SCALE</b> (if flow state allows it).</li> <li>4. Increase <b>S.V. Length</b>.</li> <li>5. Use lower frequency probe or use lower Doppler frequency on some applications.</li> <li>6. Activate <b>B-PAUSE</b> to freeze B-Mode image.</li> </ol> <i>Scan angle is critical to Doppler sensitivity.</i> | Improve Display Aesthetics | <ol style="list-style-type: none"> <li>1. Activate <b>B PAUSE</b> to freeze B-Mode image.</li> <li>2. Decrease <b>GAIN</b>.</li> <li>3. Increase <b>ACOUSTIC OUTPUT</b>.</li> <li>4. Use smaller sample volume size, whenever possible.</li> <li>5. Increase/decrease <b>Dynamic Range</b>.</li> <li>6. Adjust <b>BASELINE SHIFT</b> and <b>VELOCITY SCALE</b> to adjust spectral size.</li> <li>7. Lower Doppler <b>Sweep Speed</b>.</li> </ol> |
| Increase Spectral Clarity | <ol style="list-style-type: none"> <li>1. Activate <b>B PAUSE</b> to freeze B-Mode image.</li> <li>2. Increase <b>ACOUSTIC OUTPUT</b>.</li> <li>3. Decrease <b>S.V. Length</b>.</li> <li>4. Decrease <b>GAIN</b>.</li> <li>5. Decrease <b>Dynamic Range</b>.</li> </ol>                                                                                                                                                                             |                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                  |



# M-Mode

## Introduction

M-Mode is used to determine patterns of motion for objects within the ultrasound beam. The most common use is for viewing motion patterns of the heart.

Be sure to read and understand Acoustic Output considerations for each mode (refer to *Safety* chapter) before adjusting the Acoustic Output control or any control affecting acoustic output.

## Reading the M-Mode or Doppler Spectrum Only Display

The following information is displayed on the M-Mode or Doppler mode spectrum image:




Figure 5-27. M-Mode or Doppler Spectrum Only Display Format





Reading the M-Mode or Doppler Spectrum Only Display (cont'd)



*NOTE: Timeline formats (M-Mode and Doppler) available for display are enabled on the Set Up/Custom Display screen menus.*

*First, choose the Side/Side or Top/Bottom display style on Custom Display page 12, "Timeline Format". The two styles cannot be mixed.*

*Second, make the M-Mode format Enable/Disable selections on Custom Display page 12.*

Reading the Dual Doppler Spectrum Only Display




Figure 5-28. Dual D-Mode Display Format



*NOTE: Dual M-Mode format is not currently available.*



## M-Mode

---

### Optimizing the Timeline

#### Common Controls

##### Description

Since M-Mode is basically a single B-Mode scan vector displayed over time, basic controls that affect the B-Mode display also affect the M-Mode display.

Acoustic output, TGC, depth and B/M Gain affect both the M-Mode and B-Mode displays.

Scan area size, scan area position and reverse affect B-Mode only.



If the scan area size is reduced and the position changed, the M-Mode cursor will follow the position change to stay within the displayed scan area.

The Dual Format Keys (L–R) work the same as in dual B-Mode, but display both B-Mode and M-Mode on the left and right side of the screen.

#### Accessing/Changing

See the *B-Mode* section of this chapter for details on these controls.

## Controls

|                                                                                                                                        | M/D Cursor                                                                                                                                                                                                                              | Zoom (M-Mode)                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>                                                                                                                     | Entering M-Mode automatically assigns trackball control to the M/D Cursor.                                                                                                                                                              | Acoustic zoom can be accomplished while in M-Mode. The B-Mode reference image does not zoom but the M-Mode display will enlarge. Display (freeze) zoom is not available with M-Mode. See Figure 5–29.                                                                                                                                                                            |
| <b>Accessing/Changing</b>                                                                                                              | To enable trackball control of the M/D cursor, press the <b>M/D Cursor</b> key. The M-Mode cursor line changes color.<br><br>To terminate trackball control of the M/D cursor, press <b>Set</b> . The M-Mode cursor line changes color. | To access zoom while in B/M-Mode, press the <b>Zoom</b> key.<br><br>Rotate the <b>Zoom Size</b> control to adjust the size of the zoom area cursors.<br><br>Use the <b>Trackball</b> to move the zoom area vertically along the M-Mode cursor line.<br><br>Press <b>Set</b> to fix the zoom function and continue scanning.<br><br>Press <b>Clear</b> to exit the zoom function. |
| <b>Benefits</b>                                                                                                                        | Allows for repositioning of the M-Mode cursors after the trackball has been assigned a different function.                                                                                                                              | Allows for expanding a region of interest over the entire M-Mode timeline display.                                                                                                                                                                                                                                                                                               |
| <b>Values</b>                                                                                                                          | On or Off.<br><br>The M-Mode cursor displays in the center of the screen or changes to the active color if already present.                                                                                                             | Acoustic zoom has magnification values of 1, 1.2, 1.5, 2, 2.5, 3.0 and 4.0. The default value is set in the Set Up/Custom Display menu page 1.                                                                                                                                                                                                                                   |
| <b>Affects on Other Controls</b>                                                                                                       | Terminates trackball control from its current function and assigns it to the M-Mode cursor.                                                                                                                                             | Changes the transmit focal point position when moved vertically.                                                                                                                                                                                                                                                                                                                 |
| <b>Bioeffects</b><br> <b>Acoustic Output Hazard</b> |                                                                                                                                                                                                                                         | Moving the focal zone affects acoustic output by bringing the near field closer and by adding focal zones. The amount of increase varies depending on the probe and its frequency. Adding focal zones tends to increase TI, although it may decrease MI.<br><br>Observe the output display for possible effects.                                                                 |




Figure 5–29. Zoom (M-Mode) Display



## M-Mode

### M-Mode Sub-Menu Page 1

| M                                                                                                                                                                         | Preset                                                                                                                                                                     | Set Up                                                                                                                                                                    | ECG                                                                                                                                                                        |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Dynamic Range                                                                                                                                                             | Gray Map                                                                                                                                                                   | Rejectn                                                                                                                                                                   | Edge Enhance                                                                                                                                                               |
| <br> 48 | <br> M-2 | <br> 40 | <br> MID |

Figure 5–30. M-Mode Sub-Menu Page 1

|                                  | Dynamic Range                                                                                                                                                                                                                                                                                                                                            | Gray Scale Mapping (Gray Map)                                                                                                                                                                                                                                                                              |
|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Controls how echo intensities are converted to shades of gray, thereby creating a range of gray scale that can be adjusted. Adjustments to M-Mode's dynamic range affects the M-Mode timeline <u>only</u> .                                                                                                                                              | Gray scale mapping determines how the echo intensity levels received are presented as shades of gray.                                                                                                                                                                                                      |
| <b>Benefits</b>                  | Dynamic range is useful for optimizing tissue texture to differentiate between echo levels that are close together. Dynamic range should be adjusted so that the highest amplitude edges appear as white while lowest levels (such as blood) are just visible.                                                                                           | Displays the received echo levels with different weights on specific levels of gray. For example, a certain gray map may enhance mid level echoes over a wider range of grays verses high or low level echoes.<br><br>Allows for better differentiation between echo levels through gray levels displayed. |
| <b>Values</b>                    | The settings cycle in 6 dB steps from 30 dB to 90 dB. The default for M dynamic range is set in the Set Up/Custom Display menu page 3. Dynamic range levels are returned to the preset value when changing the following: <ul style="list-style-type: none"><li>• Probe</li><li>• Exam category</li><li>• New patient</li></ul> Pre-processing function. | There are 6 selections of gray scale mapping.<br><br>Post-processing function.                                                                                                                                                                                                                             |
| <b>Affects on Other Controls</b> | Operates only in real-time, not in freeze, Cine, Timeline replay, or VCR playback.                                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                            |



**M-Mode Sub-Menu Page 1 (cont'd)**

|                                  | Rejection (Rejectn)                                                                  | Edge Enhance                                                                                                                                                                                                                                                                                                                                                      |
|----------------------------------|--------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Deletes low level echoes.                                                            | Brings out subtle tissue differences and boundaries by enhancing the gray scale differences corresponding to the edges of structures. Adjustments to M-Mode's edge enhancement affects the M-Mode timeline only.                                                                                                                                                  |
| <b>Benefits</b>                  | Higher rejection values will remove weak, low level echoes in a displayed image.     | Modifies the M-Mode image by accentuating the interfaces between organs or vessels.                                                                                                                                                                                                                                                                               |
| <b>Values</b>                    | Off and 2 through 40 in increments of 2.<br>Rejection is a post-processing function. | The four selections are Off, Low, Mid and High. The default for M Edge Enhance is set in the Set Up/Custom Display Menu page 3. Edge enhance levels are returned to the preset value when changes are made to the following: <ul style="list-style-type: none"><li>• Application</li><li>• Exam category</li><li>• New patient</li></ul> Pre-processing function. |
| <b>Affects on Other Controls</b> | Functions in real-time, as well as in freeze, Cine, or VTR playback.                 | Operates in real-time only, not in freeze, Cine, or VCR playback.                                                                                                                                                                                                                                                                                                 |



## M-Mode

### M-Mode Sub-Menu Page 2

| M             | Preset | Set Up    | ECG        |
|---------------|--------|-----------|------------|
| Sweep Speed   |        | Color     | Tag Positn |
| ▲<br>▼<br>MID |        | Color Tag | ▲<br>▼     |

Figure 5–31. M-Mode Sub-Menu Page 2

|                                  | Sweep Speed                                                                                                                                                                                                                                        | Color                                                                                                                           |
|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Changes the speed at which the timeline updates across the display.<br><i>NOTE: Time or distance measurements are not allowed across sweep speed changes.</i>                                                                                      | Allows for enabling M-Mode image colorization.                                                                                  |
| <b>Benefits</b>                  | Speed up or slow down the timeline in order to view more or fewer occurrences over a period of time.<br>A fast speed shows less cycles but better transition definition.<br>A slow speed shows more cycles but less definition during transitions. | Displaying the gray scale as shades of color may allow for improved differentiation between echo levels.                        |
| <b>Values</b>                    | Slow (16 sec), Mid (8 sec) and Fast (4 sec).<br>Pre-processing function.                                                                                                                                                                           | On or Off.<br>The M-Mode color presented is a selection in the Set Up/Custom Display menu page 12.<br>Post processing function. |
| <b>Affects on Other Controls</b> | If more events are seen, the timeline appears smaller; if less events are seen, the timeline appears larger or more spread out.                                                                                                                    |                                                                                                                                 |

### M-Mode Sub-Menu Page 2

|                    | Color Tag                                                                                                                                                          | Tag Position                                                                                                                                                  |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Enables the colorization of a specific gray scale level range. This causes the specified gray levels to be displayed as a predetermined color in the M-Mode image. | Allows for the movement of the specified color tag range throughout the gray scale displayed.                                                                 |
| <b>Benefits</b>    | Allows for the quick recognition of specific gray levels by colorization.                                                                                          | Allows for setting the color tag to the desired level of gray scale colorization.                                                                             |
| <b>Values</b>      | On or Off.<br>The color tag size presented is determined by the Set Up/Custom Display menu page 12.<br>Post-processing function.                                   | Move the color tag up or down the gray scale. Color and size of the tag is determined in the Set Up/Custom Display menu page 12.<br>Post processing function. |



M-Mode Optimization

| M-Mode Optimization |                                                                              |                                   |                             |
|---------------------|------------------------------------------------------------------------------|-----------------------------------|-----------------------------|
| Adjustments for...  | Do the Following...                                                          | Adjustments for...                | Do the Following...         |
| Improve M-Mode      | 1. Increase/decrease <b>B/M GAIN</b> .<br>2. Increase <b>Dynamic Range</b> . | Increase Size of Area of Interest | 1. Use M-Mode <b>ZOOM</b> . |



## 3DvieW Mode (Option)

### Overview

A 3D image can be rendered or constructed from the cine data collected during a normal scan by using a process called MIP (Maximum/Minimum Intensity Projection).

When the scan data is properly collected MIP images can be used to display vascular structure using the Advanced Velocity Maps. MIP imaging is not a real-time function, but a post-processing function.

#### CAUTION



This 3DvieW option is designed for qualitative information only and has no provision for any measurement or computation ability.

No specific reconstruction resolution or accuracy is provided or claimed for this option.

#### CAUTION



Since positional information of the original scan is not available, all measurement functions cannot be performed on a MIP image. Information concerning the front or back of the MIP image may be uncertain.

Caution should be taken when viewing the MIP image due to the fact that the algorithm used in the MIP process may produce a skewness to the MIP image displayed.



*NOTE: Consider all of the above cautions when using MIP images to evaluate a patient. MIP images should not be used for quantitative information but a qualitative supplement to all other diagnostic data collected for that patient.*





## Functionality in 3D-Mode

Because of the lack of positional information for the displayed MIP image, the following functions are disabled while in 3D-Mode:

- Measurement Function
- Image Memory/Image Recall
- Biopsy Guideline Display
- Image Reverse
- Zoom Reference
- Dual Image Display
- ECG
- CWD Probes
- User Diagnostics
- 3D-Surface

## Functionality while MIP Image is Displayed

Additional disabled functions while 3D is displayed are:

- Cine Capture
- Cine Gauge Display
- Scale Marker Display

## Mode Changes in 3D-Mode

Some scan modes are not allowed while collecting images in 3D-Mode. The following are the only modes that changes are allowed:

- B-Mode > CFM-Mode
- CFM-Mode > B-Mode
- B-Mode > PDI-Mode
- PDI-Mode > B-Mode



*NOTE: 3D images cannot be rendered in PDI-Mode.*



## 3DvieW Mode

---

### CFM Map

CFM Map choices (1–6) used in 3D Mode are selected in the Set Up/Custom Display page 15.

The symmetrical maps A1 thru A4 for the Advanced velocity mode (option) are available for 3D Mode.

### Activating 3D-Mode



*NOTE: The 3D Option presets are located on Set Up/Custom Display pages 9 and 18.*

3D-Mode must be activated prior to collecting the desired scan data. To activate 3D-Mode:

2. While in normal B-Mode, display Sub-Menu page 2. Press the **Sub-Menu Select** rocker switch, if necessary, to display page 2.
3. Select **3D Mode** so that the menu item is highlighted. The message 3DvieW is displayed in the lower left portion of the screen.
4. Scan as desired in B- or CFM-Mode.



*NOTE: If 3D-Mode is enabled with the image left/right orientation reversed, the system will not automatically switch to normal orientation. The image stays in the reverse presentation because the reverse function is not available.*

### Activating 3D-Mode (cont'd)

To avoid confusion with the rendered 3D image, ensure that:

- The orientation mark on the probe and the displayed image are properly aligned.




Figure 5–32. Proper Probe Orientation



## 3Dview Mode

### Activating 3D-Mode (cont'd)

#### CAUTION



When acquiring the image data, pull the probe from Start to Finish. **DO NOT PUSH THE PROBE.**

Pushing the probe may cause the rendered image to appear confusing or improper.



*NOTE: 3D images cannot be rendered in PDI-Mode.*

3D Mode select can also be found in the CFM, PDI and Cine Sub-Menus.

#### CAUTION



When 3D-Mode is activated, Image memory will be cleared and the message displayed is:

**“Are you sure to erase all the image memory data? (y/n)”**

When 3D-Mode is activated, one half of cine memory is allocated for use by 3D-Mode.

### Cine Top Menu

Page 1 of 3

| Archive     | DICOM       | Auto Seq   | CINE        |
|-------------|-------------|------------|-------------|
| Start Frame | Review Loop | Loop Speed | Side Change |
| End Frame   | Multpl CINE | ▲▼ 1/1     | Cine Gauge  |

Page 2 of 3

| Archive       | DICOM       | Auto Seq     | CINE    |
|---------------|-------------|--------------|---------|
| CINE Capture  | Swing Speed | Aspect Ratio |         |
| Capture Frame | ▲▼ 1/2      | ▲▼ 1/3       | 3D Mode |

Page 3 of 3

| Archive     | DICOM     | Auto Seq      | CINE         |
|-------------|-----------|---------------|--------------|
| Start Frame | Render 3D | Project Tech. | View Directn |
| End Frame   | Review 3D | ▲▼ Min        | ▲▼ Normal    |

Figure 5–33. Cine Sub-Menus

## Creating a MIP Image Rendering

1. Perform a standard 2D scan in B- or CFM-Mode. The images to be used for 3D construction will be stored in Cine memory during scanning.
2. Press **Freeze** directly after the desired anatomy is scanned. Remember that Cine memory is limited. Activate Cine Gauge by rotating the Cine Scroll knob.
3. Display the Cine Sub-Menu page 3 by pressing the **Sub-Menu Select** rocker switch. Rotate the Cine Scroll knob to the position of the desired Start Frame, as shown in Figure 5–34. Select **Start Frame** on the Cine Sub-Menu page 3.




Figure 5–34. Start Frame

4. Rotate the Cine Scroll knob to the position of the desired End Frame as shown in Figure 5–35. Select **End Frame** on the Cine Sub-Menu page 3.




Figure 5–35. End Frame

### Creating a MIP Image Rendering (cont'd)

The images in Cine Memory between the start and end frames are the data that will be used for rendering the MIP image. See Figure 5–36.




Figure 5–36. Images used for MIP

5. Select **Render 3D** from the Cine Sub-Menu page 3.

Depending on the parameters set in Custom Display pages 9 and 18, the system will take some time to use the selected data to construct a set of composite images to be displayed as MIP Image.

6. When the rendering or MIP Image construction is complete, select **Review 3D** from the Cine Sub-Menu page 3 to view the MIP image.



## 3D Option Techniques

### Max/Normal

This technique takes the maximum intensity of each frame on one projection line. No threshold is used.

*Example: Data in the first frame along a projection line is compared to data in the second frame along that line. The maximum is then compared to data in the third frame. The maximum is compared to data in the fourth frame, etc.*

This is available in B- or CFM-Modes.

### Min

This technique operates the same as Max/Normal. The difference being the minimum intensity of each frame on one projection line is used. No threshold is used.

This is available only in B-Mode.

### Inverse

This technique operates the same as Min. The value of the result is inverted. No threshold is used.

This is available only in B-Mode.

### Gradient

This technique is a weighted projection method. A Threshold value (Gradient Threshold) is used.

### CFM-Mode (Gradient)

It is assumed that the anatomy of interest will have a higher intensity level than the CFM Gradient Threshold.

For each frame the object region is filled with a fixed value according to its frame number. The near frame takes a higher fixed value, while the far frame takes a lower fixed value. The value is changed by gradation.

After this process, a maximum intensity projection is made. The result is an image with the near side of the object having higher intensity and the far side having lower intensity.



## 3DviewW Mode

---

### B-Mode (Gradient)

Almost opposite from CFM processing. It is assumed that the anatomy of interest will have a lower intensity level than the B Gradient Threshold.

The near frame takes a lower fixed value, while the far frame takes a higher fixed value. The value is changed by gradation.

After this process, a minimum intensity projection is made. The result is an image with the near side of the object having lower intensity and the far side having higher intensity.

### Shade

This technique is basically a maximum intensity projection. However, this technique uses two thresholds:

- Shade Object Threshold      The value is higher than the set threshold.
- Shade Boundary Threshold      The value, to be compared to Object Threshold result, is lower than the Boundary Threshold.

Comparison will be made frame by frame after the above thresholds have been met. If the two thresholds are satisfied, no comparison will be made with the rest of the frame.

Using the Max/Normal method, lower intensity echoes in the near frame will be masked by higher intensity echoes in the far frame. The Shade Method makes it possible to show lower intensity echoes in the near frame despite the fact that there may be a higher echo projection in the far frame.

This is available in CFM-Mode only.

### 3D Presets

Presets in Set Up/Custom Display pages 9 and 18 allow for choices in display and rendering technique.

All Parameters are not available in the soft menu.

Ensure that these presets are adjusted for the desired result.





# 3D-Surface Mode

## (Option)

### Overview

This option creates a projected image different from 3DvieW. It uses images stored in Cine. The calculation used for the projection is based on the transmission principle of light using Volume-Rendering.

This function is available only in single B-Mode.

Refer to 3DvieW Mode section as needed.

#### CAUTION



This 3D-Surface option is designed for qualitative information only and has no provision for any measurement or computation ability.

No specific reconstruction resolution or accuracy is provided or claimed for this option.

#### CAUTION



Since positional information of the original scan is not available, all measurement functions cannot be performed on a 3D-Surface image.



*NOTE: Consider all of the above cautions when using 3D-Surface images to evaluate a patient. 3D-Surface images should not be used for quantitative information but a qualitative supplement to all other diagnostic data collected for that patient.*



## 3D–Surface Mode

---

### Standard Procedure

1. To access 3D-Surface Mode, the system should be active. Press 3D SURFACE from the Cine Sub-Menu page 4. While in 3D-Surface Mode, measurements cannot be performed. While 3D-Surface Mode is active, “3D-Surface” is displayed on the bottom of the screen.

| Archive     | DICOM           | Auto Seq                                                                               | CINE         |
|-------------|-----------------|----------------------------------------------------------------------------------------|--------------|
| Start Frame | Surface ROI     | Opacity                                                                                | Zoom Surface |
| End Frame   | Display Surface |  AUTO | 3D Surface   |

Figure 5–37. Cine Sub-Menu page 4

2. Scan good images for 3D-Surface and press **Freeze**. Cine is active.

## Standard Procedure (cont'd)

To avoid confusion with the 3D-Surface image, ensure that:

- The orientation mark on the probe and the displayed image are properly aligned.




Figure 5–38. Proper Probe Orientation

### CAUTION



When acquiring the image data, pull the probe from Start to Finish. **DO NOT PUSH THE PROBE.**

Pushing the probe may cause the 3D-Surface image to appear confusing or improper.



### Standard Procedure (cont'd)

3. From the Cine Sub-Menu, set the START FRAME and END FRAME for the desired images.

At this time, 3D-Surface utilizes all of Cine memory (When 3DvieW is active, half of the memory is used for 3Dview and the other half is for Cine memory. With 3D-Surface, 3DvieW is not active, so all of Cine memory is available.)

4. Select the Opacity setting. The available selections are:
  - AUTO-1 The standard setting.
  - AUTO-2 Lower contrast image displayed than AUTO-1.
  - AUTO-3 Higher contrast image displayed than AUTO-1.
5. Press SURFACE ROI. The ROI box cursor is displayed on the B-Mode image and the Cine Loop starts automatically.
6. After reviewing the selected CINE images, use the **Trackball** and **Scan Area** key to adjust the box cursor to the desired size and position for the 3D-Surface image.
7. To cancel the ROI box without rendering, press SURFACE ROI. The SURFACE ROI sub-menu selection is not highlighted, the ROI box cursor disappears and the Cine Loop stops.
8. Pressing **Set**, the following message is displayed:

**In progress, please wait.**

After a while, a temporary 3D-Surface image is displayed.
9. Depending on the probe movement used to collect the images, adjust the length/width scale with the **Ellipse** rocker switch.
10. Press **Set**, **Clear** or SURFACE ROI to fix the 3D-Surface Image.



## Standard Procedure (cont'd)



### Hints

- To zoom the 3D-Surface image, press ZOOM SURFACE.
- To return to the Cine image, select DISPLAY SURFACE. The 3D-Surface image disappears.
- When the “Not Enough Memory” message is displayed on the bottom of the screen, change the size of ROI to small until this message disappears and then press **SET**.
- To exit 3D-Surface Mode, select 3D SURFACE. Measurements will be available.



*In 3D Mode, the 3D Surface image is not active.*

*3D-Surface image is not available based on a MIP (3DvieW) image.*



# Mixed Mode Display Formats

## Display Formats

Some possible display formats are:




Figure 5–39. Top/Bottom B Mid Preset

## Display Formats (cont'd)




Figure 5-40. Top/Bottom B Large Preset




Figure 5-41. Top/Bottom B Small Preset



## Mixed Mode Display Formats

### Display Formats (cont'd)




Figure 5-42. Dual B-Mode Display Format




Figure 5-43. Zoom with Reference Mode Display Format



## Display Formats (cont'd)




Figure 5–44. Dual Format



## *Mixed Mode Display Formats*

---

This page left blank intentionally.



# Zooming an Image

## Introduction

Zoom is used to magnify a specified region of interest (ROI). The system adjusts all imaging parameters accordingly. Decreasing the size of the ROI increases the magnification factor.

Cine and VCR playback images must be frozen to be zoomed.


## Zoom Methods

The LOGIQ™ 500 offers two types of zoom capabilities, Acoustic (real-time) Zoom and Display (freeze) Zoom.

### Acoustic Zoom

Acoustic Zoom is accomplished while scanning live (real-time).

While scanning, press the **Zoom** key to activate the zoom function.



Rotate the **Zoom Size** control to the desired magnification. The magnification choices available are:

1.2, 1.5 (default), 2.0, 2.5, 3.0 and 4.0

The default value can be set in Set Up/Custom Display page 1 "Zoom Factor".

Use the **Trackball** to position the zoom region of interest.

Press **Clear** or **Zoom** to cancel the zoom function.




## Zooming an Image

### Display Zoom

Display (freeze) Zoom is accomplished after the image is frozen. This applies to the current image, a Cine image or a VCR playback image.

The magnification factor for Display Zoom is fixed at 2.0.

Press **Freeze** to stop image acquisition.



Press **Zoom** to activate the zoom function.

Use the **Trackball** to pan around the image area.

Press **Clear** or **Zoom** to cancel the zoom function.



### Hints

*NOTE: If an acoustic zoom image is frozen, display zoom will not function.*

The zoomed image can be panned around the B-Mode image by moving the **Trackball**.

Measurements can be performed on zoomed images.

Changing the depth does not affect the zoomed image, unless the zoomed area is no longer contained within the new depth.

TGC pots within the zoom ROI are active (and these pots' LEDs light up). When zoom is exited, the 8 TGC pots are reportioned to the displayed depth scale.




## Zooming an Image

---

### Zooming an M-Mode Image

If the **Zoom** key is pressed while in B/M-Mode, only the M-Mode image will be zoomed.

Press **Zoom** to display a magnified M-Mode image and a M-Mode zoom marker in the B-Mode image.



Adjust the **Zoom Size** control to the desired magnification factor. The six steps are the same as B-Mode.

Use the up/down movement of the **Trackball** to move the Zoom Markers vertically along the M-Mode cursor.

Use the side to side movement of the **Trackball** to move the M cursor and Zoom cursor throughout the image.

Press **Clear** or **Zoom** to cancel the zoom function.



**NOTE:** Display Zoom does not function on a frozen B- and M-Mode Timeline Image.

## Multi-Image Zoom

When using the Zoom function in a multiple image display format, a few basic rules apply.

- Display Zoom (Freeze) can be performed on any or all of the multiple images.




Figure 6–1. Dual Format Display Zoom

- A frozen acoustic zoom and active acoustic zoom can be displayed simultaneously.



Figure 6–2. Dual Format Display, Frozen and Active

- Acoustic Zoom and Display Zoom cannot be mixed on the same display.




Figure 6–3. Dual Format Display/Acoustic Zoom Unavailable



# Freezing an Image

## Introduction

Freezing a real-time image stops all acquisition of information into system memory.

This allows for measurements, annotations, printing or storage of images into temporary memory.

VCR playback images can also be frozen for similar reasons.

## Post processing of the image

The following post processing image parameters can be changed to affect the appearance of a frozen image:

- Display Zoom (x2 only)
- Reject
- Color (B Color or CFM Color Maps)
- Gray Scale Map Selection
- Image Rotation
- Color Tag



## Freezing an Image

---

### Freezing an Image (Freeze Key)

To freeze an image:

- Press **Freeze**. The key backlights.



*NOTE: If both B- and M-Modes are active, the B-Mode and timeline trace stops immediately. Use the Cine Scroll Control to start CINE review.*

To reactivate the image:

- Press **Freeze** again. Deactivating Freeze restarts the B-Mode and timeline after a black and white bar indicating discontinuity is inserted in the timeline (M-Mode Display).




*NOTE: Deactivating Freeze erases all measurements and calculations from the display (but not from the report page).*

*Selecting a new probe unfreezes the image.*

### Freezing an Image (Foot Switch option)

Foot Switch



Toggle **Freeze** on and off by pressing the foot switch. Pressing the left foot switch backlights the **Freeze** key on the console.





# Using Cine

## Introduction

Cine is useful for focusing on images during a specific part of the heart cycle or to find an image before the patient moved or breathed.

Cine images are constantly being stored by the system. Cine memory capacity is 31 (0–30) frames with 128 scan lines. The Cine Memory stores the most recent data available for playback or manual review via Cine. The amount of time represented depends on the system's frame rate, scan mode, image size and other parameters.

Timeline data is continually stored at four times the display width of timeline data (and updates the corresponding B-Mode images).

View Cine as a continuous loop via Cine Loop or manually review Cine images frame by frame via the Cine Scroll Control.

Data in Cine is available until new data is acquired. Cine is stored in the system's memory and can be transferred to image memory, video page printer, multi-image camera or optional optical disk (MOD).

An Extended Cine Memory Option is available to increase Cine storage capabilities.

## Cine memory

Cine memory is erased when changing the following:

- Probe
- Scan Mode
- Depth
- Display format (zoom, dual, rotate)
- Timeline Sweep Speed (D/M-Modes)
- Changing the PWD Velocity Scale (spectrum only)
- Changing the Color Flow Velocity Scale

## Cine functionality

Post Processing functions can be performed while in Cine such as:

- Measurements and calculations
- Color Baseline shift
- B Color
- Color Flow Velocity Tag
- Color Flow Display Threshold
- CFM Spectrum Invert
- Zoom
- Rejection
- Gray Scale Maps
- Edit annotations

## Accessing Cine

To access and manually review Cine:

1. Press **Freeze**.

*NOTE: One click of the **Cine Scroll** knob erases displayed measurements.*

2. Rotate the **B/M Gain/Cine Scroll** knob to activate Cine.
3. Rotate the **Cine Scroll** dial left (backward) and right (forward) to move through the images in Cine memory.
4. The current frame on the Cine gauge moves and the Cine frame number is displayed on the lower left side of the screen, above the depth.




Figure 6–4. Cine Gauge Display



## Accessing Cine (cont'd)



*NOTE: Cine frame number 0 is the most current image. The higher the Cine frame number, the older the image.*

*Depth, dynamic range, and gain parameters are valid for Cine frame number zero only.*

## Using Cine Loop

Specify which sections of Cine memory to playback by creating a Cine Loop.

To create a Cine Loop:

1. Press **Freeze**. The key backlights.
2. Rotate the **Cine Scroll** dial to the desired starting position in Cine memory. The Cine Sub-Menu is displayed on the Soft Menu display.

| Archive     | DICOM       | Auto Seq                                                                                | CINE        |
|-------------|-------------|-----------------------------------------------------------------------------------------|-------------|
| Start Frame | Review Loop | Loop Speed                                                                              | Side Change |
| End Frame   | Multpl CINE |  1/1 | Cine Gauge  |

Figure 6–5. Cine Sub-Menu

3. Select **START FRAME** on the Cine Sub-Menu.




Figure 6–6. Cine Gauge Frame Start

4. Rotate the **Cine Scroll** dial to the desired ending position in Cine memory.
5. Select **END FRAME** on the Cine Sub-Menu.

## Using Cine Loop (cont'd)




Figure 6–7. Cine Gauge Frame End

6. Select **REVIEW LOOP** on the Cine Sub-Menu.

Cine Review Loop playback begins right away at the selected speed.




Figure 6–8. Cine Loop Operation



*NOTE: Cine Loop is not available during Timeline Review.*

To deactivate Cine Loop:

- Select **REVIEW LOOP** on the Cine Sub-Menu or turn the Cine Scroll knob. This returns to manual Cine review.

To reactivate Cine Loop:

- Select **REVIEW LOOP** on the Cine Sub-Menu.



## Cine Loop Speed

To adjust the speed of viewing Cine Loop playback, access **LOOP SPEED** from the Cine Sub-Menu.

| Archive     | DICOM       | Auto Seq                                                                              | CINE        |
|-------------|-------------|---------------------------------------------------------------------------------------|-------------|
| Start Frame | Review Loop | Loop Speed                                                                            | Side Change |
| End Frame   | Multpl CINE |  1/1 | Cine Gauge  |

Figure 6–9. Cine Sub-Menu (Loop Speed)



*NOTE: Cine cannot be viewed faster than real-time.*

To increase playback speed, press the top of the **LOOP SPEED** rocker switch.

Press the bottom of the **LOOP SPEED** rocker switch to slow down playback speed.

Each press of the rocker switch cycles to the next speed setting (1/1, 1/2, 1/4 or 1/8).

## Multpl CINE

The Multpl CINE Soft-Menu selection is used in Dual B-Mode Cine Operation.

With Multpl CINE on (menu highlighted), Cine review of the two B-Mode images can be done simultaneously.

With Multpl CINE off, Cine review is done separately using the **L/R** keys to designate the active Cine image.

Review with Multpl CINE on/off can be accomplished with the **Cine** knob or review loop.

## Side Change

This menu selection is used with B/M (D)-Mode imaging.

Side Change is used to toggle between B image Cine scroll and Timeline (M/D) Cine scroll.

| Archive     | DICOM       | Auto Seq                                   | CINE        |
|-------------|-------------|--------------------------------------------|-------------|
| Start Frame | Review Loop | Loop Speed                                 | Side Change |
| End Frame   | Multpl CINE | <div> <div>▲</div> <div>▼</div> </div> 1/1 | Cine Gauge  |

Figure 6–10. Cine Sub-Menu

## CINE Gauge

Used to toggle the Cine gauge graphic display on or off.

## CINE Capture (option—color images only)

CINE Capture is a post processing version of the Capture selection in the CFM/PDI sub-menu.

| Archive       | DICOM | Auto Seq | CINE |
|---------------|-------|----------|------|
| Cine Capture  |       |          |      |
| Capture Frame |       |          |      |

Figure 6–11. Cine Sub-Menu

Selecting **CINE Capture** will search through all images between the start frame and end frame and display each peak or the highest velocity. Adjust the start frame and end frame points to limit the image frames used in the process.



## Capture Frame (option—color images only)

Capture frame can be used to eliminate specific image frames from the Cine Capture process.

Use the **Cine Scroll** knob to display an image frame to be eliminated, and press the **Capture Frame** rocker switch. A small mark will be displayed on the Cine gauge. Do this for all images that are not to be included in the Cine Capture process. When all frames to be eliminated are marked, press the **Cine Capture** rocker switch to display the peak velocity detected color image.

## Exiting Cine

To exit Cine, press **Freeze**.

## Helpful Hints



### Hints

The following hints can help when freezing an image:

- Color, multi-format or M/D images take up more memory than normal B-Mode. Therefore, less frames of information are available for Cine storage.

## ECG/Cine Gauge/Image Tracking

When an ECG waveform is displayed with the cine gauge, an arrow pointer appears above the ECG waveform.

As the cine gauge marker is moved with the Cine Scroll control, the arrow above the ECG waveform moves to indicate where on the ECG cycle the displayed image was taken.

# Annotating an Image

## Introduction

The annotation keyboard is always active. Upon starting a new patient or power up, the underscore cursor appears in the mode's home position. Annotation can commence after using the Trackball to specify where the comment should start.




Pressing the **Comment** key assigns the trackball function to controlling the cursor and displays the annotation library in the soft menu.

The Comment key enables the image text editor and displays the annotation library soft menu. After the **Comment** key is pressed, text can be added through the Soft Menu comment library or by typing comments from the alphanumeric keyboard. Turning off the cursor can be done by pressing the **Set** key. Pressing the **Clear** key erases all comments.



*NOTE: If the "Use A/N Keys for User Define" preset in Set Up/System Parameters page 3 is ON, the Comment key must be pressed first to activate the annotation function.*



1 Underscore Cursor or Block

Figure 6–12. Annotation Cursor in Home Position

Annotations are input in type-over, not insert, mode. Be careful not to write over text when editing.





## Annotating an Image

---

### Introduction (cont'd)



*NOTE: The comment function will work with any report page. When the cursor is in a field designated for comments, the comment key will illuminate and the annotation library will appear in the Soft Menu. All comment, edit and annotation functions are then available.*


All annotations are permanently retained with the image. However, annotations are erased at power down or when **Clear** or **New Patient** are pressed.

In addition, the display's home position can be changed (preferred annotation area) for each display so that all subsequent annotations begin in the same spot.



**Cursor Home:** Established in Set Up/Preset Program page 1, the cursor returns to the upper left part of the screen or the user specified position.


A new cursor home position is established by placing the cursor in the desired position and pressing **Ctrl** and **M**. The new cursor position is placed in the Preset Program values.



To end the Comment/Library annotation function, press **Clear**. Comments will be erased.


"Comment Clear Key Function" is a parameter found in the Set Up/Preset Program page 1. This parameter can be set to erase all comments on the screen or just the comment line the cursor is on. With the line selection, the second time **Clear** is pressed, all comments will be erased.

"Measurement Clear Operation" is a parameter found in the Set Up/Preset Program page 3. This parameter can be set to have only measurements erased with the **Clear** key or measurements and comments.




The **Trackball** and **Keyboard Arrow** keys are used to move the cursor to the desired position on the image.

When the blinking cursor is in the desired position, comments may be typed in or a selection can be made from the annotation library.




The **Set** key is used to end the Comment/Library function.

## Introduction (cont'd)



The **Tab** key will move the cursor to the right every eight characters or to the next word depending on the preset parameter. See 6–21 for more details.



**Shift** and **Tab** moves the cursor in the same manner but to the left.

## Annotation Library

To reduce the amount of time spent annotating an image, store often-used annotations in the Annotation (Comment) Library. These scripts can be up to 20 characters in length. As many as 24 scripts can be saved for each user application preset within each exam category.

## Entering/Editing the Library

Access the Annotation Library by selecting Set Up/Preset Program. Pages 7 and 8 of the Preset Program menus will be the 24 annotation library selections.

Use the **Trackball** to move the arrow cursor to the desired Annotation Library location number.



*NOTE: Remember the first 8 selections will appear on sub-menu page 1. Annotation Library 9 to 16 will be on sub-menu page 2. Annotation Library 17 to 24 will be on sub-menu page 3.*

Press **Set**. The 20 character space for that library location will be in reverse video. Add or edit the desired script.

Select the next library location and press **Set**.

Continue until all additions or edits are complete.

To save all entries and edits, **Trackball** to the SAVE selection and press **Set**.

The system prompt will read “overwrite existing data? ‘Y’ or ‘N’”. Press ‘Y’ to rewrite and save data. Press ‘N’ to program a new annotation library.

To avoid saving unwanted changes (all current changes), select RESET and press **SET**.

## Annotating an Image

---

### Displaying annotation scripts

To review the Comment Library scripts:

Press **Comment**. The Library menu for the designated exam category appears.

Use the **Sub-Menu Select** rocker switch to cycle through the three pages of scripts available.

|                 |                 |                 |                 |
|-----------------|-----------------|-----------------|-----------------|
|                 |                 |                 |                 |
| Annotation<br>1 | Annotation<br>3 | Annotation<br>5 | Annotation<br>7 |
| Annotation<br>2 | Annotation<br>4 | Annotation<br>6 | Annotation<br>8 |

Figure 6–13. Comment Library Sub-Menu Display

The Library Sub-Menu display shows the first 14 characters of each script.



### Hints

Print a hard copy of these codes from the Set Up/Preset Program Menus if there are several operators in the department.

There is a space following all programmed annotations.  
Example: PANCREAS\_

Type in the 8 most used scripts first, so that they appear on the first Sub-Menu page.



## Adding Comments to an Image

To annotate an image:

- Type comments where the cursor is currently located (the display's home position) or use the **Trackball/Arrow** keys to place the annotation cursor in the desired location before typing.
- Press **Return** to move to the next line.



*NOTE: Annotations wrap to the next line when they are within one character of the right margin.*

The word wrap starts one line below the start of that annotation.

Annotations appear on all prints, photos, and VCR recordings.



Figure 6–14. Next Line Word Wrap

If the cursor appears at the right edge of the lowest line, or a word cannot be completed in the lower right corner, word wrap cannot be executed.

To annotate an image using the Library:

- Press **Comment** and move the annotation cursor to the desired start position using the **Trackball/Arrow** keys.
- Use the **Sub-Menu Select** rocker switch, if necessary, to display the desired script.
- Press the top/bottom of the appropriate rocker switch to write the library script starting at the position of the cursor.
- The same word wrap principles apply for library scripts as typed comments.

## Special Annotation Keys

Some special annotation symbols can be used by activating the **Blue Shift** or **Red Shift** keys. **Red Shift** does not function while **Blue Shift** is active. The Blue Shift and Red Shift keys act as a lock function (similar to the Caps Lock key).

Activating **Blue Shift** will cause the **arrow**, **female** and **male** symbols to be printed on the screen during the comment function when the keys shown in Figure 6–15 are pressed.




Figure 6–15. Blue and Red Shift Keys

The **Red Shift** key enables the special symbols shown in red on the keyboard.

- The red symbols shown in the lower right portion of a key will print when the **Red Shift** is active.
- The red symbols shown in the upper half of a key will print if the normal shift key is held down while the **Red Shift** is active.
- Red symbols can be used in any language, but can only be used with the proper designated letters.



## Editing Annotations

On screen annotations can be revised. Revision can be accomplished by adding or deleting text, or completely removing all annotations by pressing **Clear**.

### Editing while annotating

Backspace over any error(s) made. Blank spaces take the place of the letter(s) that was there. Continue typing the annotation after backspacing over all incorrect letters.



*NOTE: Text does not adjust automatically.*

To delete previous character(s):

- Press **Backspace** as many times as necessary to make the deletion.
- Retype the annotation from the point where backspacing was stopped.
- Position the cursor and type over existing text.

To move through the text eight characters at a time:

- Press **Tab** to move to the right (Preset Keyboard Tab = Normal).
- Press **Shift** and **Tab** to move to the left.

To move through the text a word at a time:

- Press **Tab** to move to the right (Preset Keyboard Tab = Word)
- Press **Shift** and **Tab** to move to the left.



*NOTE: The Tab selection is found on the Set Up/System Parameter menu page 1.*

## Body Patterns

An additional way to annotate the image display is with body patterns. Body patterns are a simple graphic of a portion of the anatomy that is frequently scanned.

This body pattern is generally displayed in the lower left corner of the screen. Its placement can vary with the format of the display.

Along with the body part graphic is a marker that illustrates the probe position. This marker can be placed with the **Trackball** and rotated with the **Zoom Size/Rotation** control.

The body pattern and probe marker can serve as a reference for patient and probe positioning when images are archived.




Figure 6–16. Body Pattern with Probe Marker



Names of available body patterns can be displayed in the Soft-Menu by pressing the **Body Pattern** key.

Body Pattern packages are displayed according to exam category and preset.

The 8 body pattern packages may be customized to accommodate user preference. The 16 individual body patterns in the eight packages can be changed.

## Body Patterns (cont'd)

Use the sub-menu rocker switches to select the desired pattern to be displayed. The Set Up/System Parameters menu page 4 allows for the choice of displaying the pattern only during freeze or at all times.

Each Body Pattern Package (1–8) can be customized and the package for each preset selected in the Set Up/System Parameters menus.

The Body Pattern Package (1–8) is selected from the Set Up/Preset Program page 1.

Figure 6–17 shows the body patterns available to be preset and the order in which they appear in the selection cycle.




Figure 6–17. Body Patterns





## Annotating an Image

### Body Patterns (cont'd)




Figure 6–17. Body Patterns (cont'd)

## Body Patterns (cont'd)




Figure 6–17. Body Patterns (cont'd)



**NOTE:** Refer to 14–66 for details on selecting the body pattern package to be used.



## *Annotating an Image*

---

This page left blank intentionally.



# Introduction

## Overview

Measurements and calculations derived from ultrasound images are intended to supplement other clinical procedures available to the attending physician. The accuracy of measurements is not only determined by system accuracy, but also by the use of proper medical protocols by the user. When appropriate, be sure to note any protocols associated with a particular measurement or calculation. Formulas and databases used within the system software that are associated with specific investigators are so noted. Be sure to refer to the original article describing the investigator's recommended clinical procedures.

Measurements can be made in all modes and image formats, including real-time, freeze, CINE or VCR playback. The selected exam category defines the default calculations displayed on the Softkey Sub-Menu.

Measurement values update on the display until the measurement has been completed. 18 lines of measurements can be displayed on the screen at one time.

The measurement cursor is not automatically fixed when the **Record** key is pressed. Record all desired measurements displayed on the screen before erasing them.

## Measurement Controls




Figure 7–1. Locating Measurement Controls

|                    |                                                                                                                                                                                                                     |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Measurement</b> | The Measurement key is used to enable all types of basic measurements. It displays the measurement and calculation soft menu for the current exam category.                                                         |
| <b>Ellipse</b>     | The Ellipse rocker switch is used to activate the ellipse measurement function after the first distance measurement has been set. It also toggles which cursor is the movable cursor during the ellipse adjustment. |
| <b>Trackball</b>   | Used to move the measurement cursors.                                                                                                                                                                               |
| <b>Set</b>         | Fixes the cursor for measurements or completes the measurement sequence.                                                                                                                                            |
| <b>Clear</b>       | Erases measuring cursor and measurement data from the display during a measurement sequence. Clears all cursors and measurements from the display when not performing a measurement sequence.                       |

## Cursors

During the measurement process, the fixed and active cursors can either be an “X” or “+” symbol. The “*Active Measurement Cursor*” preset, located on System Parameters page 3, allows for the choice of the Measurement cursor.

Once the measurement sequence is complete, the cursor symbol changes to one of the eight shown below.



Figure 7-2. Measurement Completed Cursors

## General Mode Measurements Method

When performing measurements, it is possible to choose the display of and location of the caliper for subsequent measurements while in B, M or Doppler modes. These presets (“*B Auto Caliper Control*”, “*M Auto Caliper Control*” and “*D Auto Caliper Control*”) are located on Preset Program page 4.



## Measurement Key

The following table indicates the types of generic measurements available when the **Measurement** key is pressed and no specific calculation is chosen. The type of measurement depends on the current scan mode and the number of times the **Measurement** key is pressed.

| Key Pressed<br>(while frozen) | MODE                       |                                     |                                         |                            |
|-------------------------------|----------------------------|-------------------------------------|-----------------------------------------|----------------------------|
|                               | B                          | Spectral Doppler<br>(with B)        | M<br>(with B)                           | CFM<br>(with B)            |
| Once                          | Distance<br>Ellipse/Circle | Peak Velocity                       | Tissue Depth<br>(Distance)              | Distance<br>Ellipse/Circle |
| Twice                         | Trace/Circle               | TAMAX/TAMIN/<br>TAMEAN/TAMODE       | Time Interval                           | Trace/Circle               |
| Three Times                   | Gray Scale<br>Echo Level   | 2 Velocities<br>Slope/Time Interval | Depth Difference<br>Slope/Time Interval | Point Velocity             |
| Four Times                    | Arrow Mark *               | Time Interval                       | Arrow Mark *                            | Gray Scale<br>Echo Level   |
| Five Times                    |                            | Arrow Mark *                        |                                         | Arrow Mark *               |

\* Assuming the “Arrow Mark Display” preset is set to *On*.

Table 7–1. General Measurements by Mode

The “Arrow Mark Display” preset, found on Preset Program page 4, allows the display of an arrow mark when the Measurement key is pressed a certain number of times.



### General Instructions

Prior to making measurements, **Freeze** should be pressed to stop the acquisition of the image data.

For some types of measurements, press the **Measurement** key to toggle between active cursors for fine adjustment prior to completing the measurements.

Pressing **Clear** prior to completing the measurement sequence erases the active measuring marker and the current data measured.

Pressing **Clear** after the sequence is complete erases all data that has been measured to this point, but not data entered on to report pages.

Any measurement can be repeated by selecting that measurement again from the Sub-Menu.

### Erasing Measurements

The following actions erase measurements from the system's memory:

- Rotating the **Cine Scroll** knob erases measurements.
- Pressing **New Patient** erases all measurements and calculations on the display and clears the report pages.
- Unfreezing the image erases measurements from the display only, not from the report pages.
- Adding a new measurement that exceeds the maximum number of allowable measurements erases the first or oldest measurement.
- Pressing **Clear** erases all completed measurements and calculations on the display. Measurements and calculations, however, remain on the report pages.

"Measurement Clear Operation" is a parameter found in the Set Up/Preset Program page 3. This parameter can be set to have only measurements erased with the **Clear** key or measurements and comments.





# Mode Measurements

## B-Mode Measurements


Three basic measurements can be made in B-Mode.

- Distance
- Circumference/Area (Ellipse or Trace Methods)
- Gray Scale Echo Level

## Distance Measurement

Distance Measurements are typically made in the B-Mode portion of the image. To make a distance measurement:

Press the **Measurement** key once to display a “×” cursor.



Use the **Trackball** to move the “×” cursor to the measurement start point.

Press **Set** to fix the measurement start point and display a second “×” cursor.

Use the **Trackball** to move the second “×” cursor to the measurement end point.



## Mode Measurements

---

### Distance Measurement (cont'd)



#### Hints

Press **Measurement** to toggle between which cursor is activate.




Press **Clear** once to erase the second measuring marker and the current data measured and start the measurement again.

Press **Set** to complete the measurement and fix the Distance value displayed on the bottom part of the screen.

Pressing **Clear** after the sequence is complete erases all data that has been measured to this point, but not data entered on to report pages.

### Circumference/Area (Ellipse) Measurement

An ellipse can also be used to measure the circumference and area of an organ. To measure with an ellipse:



(continued)


Press **Measurement** to set up the measurement mode. The “×” cursor appears on the display.

Use the **Trackball** to move the “×” cursor to either end of the major axis of the area to measure.

Press **Set** to fix the start point cursor.

Use the **Trackball** to move the second cursor to the major axis measurement end point.

## Circumference/Area (Ellipse) Measurement (cont'd)



Press the top of the **Ellipse** rocker switch. An ellipse having an initial circle shape appears.

Use the **Trackball** to position the ellipse, as necessary, and to size the measured axis.

Press the top of the **Ellipse** rocker switch to increase the size of the minor axis.

Press the bottom of the **Ellipse** rocker switch to decrease the size of the minor axis.



## Hints

Press **Measurement** to toggle between which cursor is activate.

Press **Clear** once to erase the ellipse and the current data measured. The original cursor is displayed to restart the measurement.

Press **Set** to complete the ellipse measurement and display the circumference and area.

Pressing **Clear** a second time exits the measurement function without completing the measurement.







## Mode Measurements

---

### Circumference/Area (Trace) Measurement

To trace the circumference of a portion of the anatomy and calculate its area:



Press **Measurement** twice to display a dot “⊙” cursor on the screen. The display on the bottom of the screen shows the circumference in cm.

Use the **Trackball** to move the dot “⊙” cursor to the measurement start point.

Press **Set** to change the dot “⊙” cursor to a “×” cursor.

Use the **Trackball** to trace the measurement area. The circumference displayed on the bottom of the screen will change with the tracing.

*NOTE: Press **Clear** to erase the measuring marker and the current data measured and redo the measurement.*

Press **Set** to fix the Circumference value displayed on the bottom of the screen and calculate the area.





*NOTE: When using the trace method, the **Zoom Size/Rotation** knob can be used to edit the trace line. Turn it clockwise or counterclockwise to erase the line (bit by bit) back from its current point. This erase function is limited to a small portion of the most recent trace.*




## Echo Level Measurement

To make an echo level measurement:

Press **Measurement** three times to enable the echo level function. A box “□” cursor appears.



Use the **Trackball** to move the box “□” cursor over the measurement area.



Press **Set** to fix the Echo value displayed.



### Hints

- The size of the box cursor can be changed by using the Set Up/Preset Program menu page 3.
- The box is actual size and will appear to get larger or smaller with changes in depth (scale).
- The echo level measurement is only available on a frozen image, not on a B-paused image.



## *Mode Measurements*

---

### **CFM B-Mode Measurements**

While in the Color Flow and B-Mode combination, 4 basic measurements can be taken.

- Distance
- Circumference/Area (Ellipse or Trace Methods)
- Gray Scale Echo Level
- Velocity Point

#### **Distance**

Distance is measured the same as B/W B-Mode. Refer to 7–7.

#### **Trace**

Circumference and area are measured the same as in B/W B-Mode. Refer to 7–8.

#### **Gray Scale Echo Level**


Gray scale echo level is measured the same as in B/W B-Mode. Refer to 7–11.

## Velocity point


This measurement provides an estimate of the average (or mean) velocity at one point in the CFM B-Mode display. This is achieved by calculating the numerical value equivalent to the color at the selected point in the display and adjusting it to the selected flow angle. This procedure cannot be used to determine maximum, minimum or other specific velocity estimates.

To measure the velocity at a certain point in the color display:



Press **Measurement** three times. A “×” cursor appears.



Use the **Trackball** to move the cursor to the desired portion of the display.



Use the **Zoom Size/Rotation** control to adjust to the angle of blood flow.



Press **Set** to complete the measurement and to display the velocity and angle values.

## CAUTION



Potential for misdiagnosis. The CFM B-Mode velocity point measurement only provides an average velocity estimate and should not be used as a substitute for PW Doppler velocity measurement. Use PW Doppler Mode to obtain measurements from specific portions of the Doppler spectrum.

## Mode Measurements

---

### Doppler Mode Measurements

4 basic measurements can be made in Doppler Mode or the Doppler Spectrum portion of the display.

- Peak Velocity
- TAMAX/TAMIN/TAMEAN/TAMODE (Manual or Auto Trace)
- Two Velocities with the Slope and Time Interval between them.
- Time Interval

#### Peak Velocity

To measure peak velocity:



Press **Measurement**. The “×” cursor with a vertical dotted line appears.

Use the **Trackball** to move the cursor to the desired point of measurement in the Doppler Spectrum.

Press **Set** to fix the velocity cursor. The velocity measurement in m/s or cm/s will be displayed at the bottom of the screen.

The “*Velocity Unit*” preset, located on Preset Program page 3, allows the unit of display for all velocity measurements to be selected.

The “*Display MaxPG with Velocity*” preset, located on Preset Program page 4, allows MaxPG to be displayed with the Peak Velocity measurement.






### TAMAX/TAMIN/TAMEAN/TAMODE (Manual or Auto Trace)

The value measured depends upon “D Realtime Trace Method” selection in the Set Up/Custom Display page 13. The 4 selections available are: Peak (TAMAX), Floor (TAMIN), Mean (TAMEAN) and Mode (TAMODE).

To trace TAMAX/TAMIN/TAMEAN/TAMODE:



Press **Measurement** twice. A dot “⊙” cursor appears if the default method is manual. A “×” cursor with a vertical dotted line appears if the trace default method is auto.

Use the **Trackball** to move the cursor to the trace start point in the Doppler spectrum.

Press **Set** to fix the start point.

If D Realtime Auto Trace (on or calc) is the default method or Trace Auto is on, use the **Trackball** to position the vertical cursor at the second value (end point).


Press **Set**. The system automatically fixes both cursors and traces the maximum value between the two points. That value is displayed at the bottom of the screen.



## Mode Measurements

---

### TAMAX/TAMIN/TAMEAN/TAMODE (Manual or Auto Trace) (cont'd)



If D Realtime Auto Trace (off) is the default method, use the **Trackball** to trace the maximum values of the desired portion of the spectrum.

Press **Set** to complete the measurement.




**NOTE:** When using the manual trace method, the **Zoom Size/Rotation** knob can be used to edit the trace line.


Turn it clockwise or counterclockwise to erase the line (bit by bit) back from its current point. This erase function is limited to a small portion of the most recent trace.

## Slope/Time Interval


Two velocity values, the time interval (sec) and slope ( $m/s^2$ ) can be measured by:




Press **Measurement** three times. A "×" cursor with vertical dotted lines appears.




Use the **Trackball** to move the cursor to the measurement start point.



Press **Set** to fix the start point.



Use the **Trackball** to move the end cursor to the measurement end point.



Press **Set** to complete the measurement.

The two peak end point velocities, time interval and slope are displayed.




## Mode Measurements

---


### Time Interval

To measure a horizontal time interval:

Press **Measurement** four times. A “×” cursor with a vertical dotted line appears when the cursor is in the Doppler spectrum portion of the display.




Use the **Trackball** to move the cursor to the measurement start point.



Press **Set** to fix the first cursor.



Use the **Trackball** to move the end cursor to the measurement end point.



Press **Set** to complete the measurement.

The time interval between the two cursors is displayed at the bottom of the screen.



## Realtime Doppler Calculations (option)

The ultrasound system has the ability to automatically measure, calculate and display specific parameters using the Doppler Auto Trace in realtime. This applies to all exam categories except GYN and Cardiology with AMCAL option.

### Soft Menu Changes

Page 4 of the PWD soft menu contains 3 selections that affect Realtime Doppler Calculations. These are:

| PWD           | Preset    | Set Up       | ECG |
|---------------|-----------|--------------|-----|
| Realtim Trace | Calc Dir. | Trace Method |     |
| ▲▼ OFF        | ▲▼ Combo  | ▲▼ PEAK      |     |

Figure 7–3. PWD Sub-Menu Page 4

- Realtime Trace: On, Off or Calc
- Calc Dir.: Forwd, Revrs or Combo
- Trace Method: Peak, Mean, Mode, Floor

Refer to 5–60 for details.



*NOTE: Trace method selection is important in the determination of values used to calculate measurement selections.*

Example: “Mean” is the Trace Method selection.

- Enable Trace Auto and perform Doppler scan.
- TAMEAN and TAMAX (peak) are traced automatically.
- Freeze the image and select “Transf. CALCs”.
- Values are transferred to measurement area.
- Select PI and TAMAX value is used in calculations.
- Select other measurements like FV and TAMEAN is used in the calculations.

## Soft Menu Changes (cont'd)

Each measurement menu in all exam categories (except GYN and Cardiology with AMCAL option) has a selection for:

|              |              |               |                 |
|--------------|--------------|---------------|-----------------|
|              |              |               |                 |
| S/D<br>Ratio | A/B<br>Ratio | HEART<br>RATE | TRACE<br>AUTO   |
| PI           | RI           |               | Transf<br>CALCS |

Transf CALCS is only available with Realtime Doppler Calculation option installed

Figure 7–4. General Calculations Sub-Menu Page 2

- Trace Auto Turn Auto Trace on/off
- Transf. Calcs: Transfers automatic calculations currently displayed to the Measurement results area.

*NOTE: Only OB or Vascular PI, RI and S/D (D/S) ratio are transferred to the OB, Vascular or Carotid report pages.*



## Setup

The realtime measurement/calculation display is set with factory defaults for each exam category. This display arrangement, content and calculation method can be customized in Set Up/ Preset Program menu page 4 and Set Up/Custom Display menu page 13.

See 14-71 for details on setting:

- Diastole Velocity for PI Vmin, Vd
- Diastole Velocity for RI Vmin, Vd
- Diastole Velocity for S/D (D/S) Vmin, Vd
- Averaging Number for D Realtime Calc. 1–5
- RI Calculation to bidirectional flow: Off, On

See 14–88 for details on customizing the appearance and content of the realtime calculation display window.



## Mode Measurements

---

### Realtime Calc Window

This display appears in the upper left corner of the monitor. It consists of 8 lines of information. If some values are not detected, “????” will be displayed.

Example:

- “MEAN” is selected as the trace method. TAMEAN and TAMAX are traced automatically.
- If TAMAX, TAMEAN and TAMOD (mode) are selected to be displayed in the Realtime Doppler Calculation window, “????” will be displayed as the value for TAMOD because it was not traced.



*NOTE: It is important to remember the trace method selected in the soft menu or Set Up/Custom Display menu page 13.*

### Display Modes

Enter the Doppler Realtime Calc mode by selecting Realtime Trace = Calc from the PWD soft menu.

Ensure Realtime Trace = On in the measurement menu.

The Auto Trace display is available only in :

Display Format = Single D Format, Single B/D Format

Scan Mode = PD, B/PD, BCFM/PD or CWD Mode

Sweep Speed = Slow, Mid, Fast

Exit the Doppler Realtime Calc mode by turning Trace Auto OFF in a measurement menu or selecting Realtime Trace = Off in the PWD soft menu.

### Doppler Trace Indicators

Each time end diastole is detected in the Realtime Doppler Trace, a white or colored vertical line is placed on the display. This indicates the area before and after the heart rate that is used to calculate the displayed measurements/calculations.

Gray markers are used to indicate areas that are past or beyond the current displayed calculations. This area can be accessed after pressing **Freeze** and using the Cine function. If end diastole cannot be detected, no markers are displayed.



## Mode Measurements

---

### Transferring Calcs

Selecting **Transf. CALCS** copies data from the Realtime Doppler Calculation window to the normal calculation result window at the bottom of the screen.

The following is an example of transferring calculations in the OB exam category:

1. Select OB PI, OB RI and OB S/D (D/S) from the Doppler Realtime Calculation : Submenu selection in the Set Up/Preset Program menu page 7.
2. Turn on the Realtime Doppler Calculation selection in the Softmenu.
3. The selected measurements will be displayed in the upper left portion of the screen.
4. Press **Freeze**.
5. Choose "Select Location" from the Measurement Softmenu.
6. Select "Transf. CALCS" to transfer the calculations to the selected location on the report page.

The Transf. CALCS function is not available unless values are displayed in the Doppler Realtime Display window.



*NOTE: The Advanced Vascular Calculation Sub-Menu has an option called Manual Calc. When Manual Calc is activated, Auto Calc measurements will not be transferred to the Report Page or result window.*





## M-Mode Measurements

Basic measurements that can be taken in the M-Mode portion of the display are:


- Tissue Depth (Distance)
- Time Interval
- Depth Difference with Time Interval and Slope

### Tissue depth


Tissue depth measurement in M-Mode functions the same as the distance measurement in B-Mode. It measures only vertical distance between points.

Scan patient with M-Mode timeline displayed.


Press **Measurement** once. A “×” cursor with a vertical dotted line appears.




Use the **Trackball** to move the “×” cursor to the most anterior point to be measured.



Press **Set** to fix the start point.



Use the **Trackball** to move the second point to the most posterior point to be measured.



Press **Set** to complete the measurement.

The vertical distance between the two points is displayed at the bottom of the screen.




## Mode Measurements

---

### Time Interval



To measure a horizontal time interval:




Press **Measurement** twice. A “X” cursor with a vertical dotted line appears when the cursor is in the M-Mode portion of the display.




Use the **Trackball** to move the cursor to the measurement start point.



Press **Set** to fix the first cursor.



Use the **Trackball** to move the end cursor to the measurement end point.




Press **Set** to complete the measurement.

The time interval between the 2 cursors is displayed at the bottom of the screen.





## Depth Difference with Time and Slope


To measure depth difference:



Press **Measurement** 3 times. With the cursor in the M-Mode timeline area, a “×” cursor with a vertical dotted line is displayed.




Use the **Trackball** to position the measurement start point.



Press **Set** to fix the start point.



Use the **Trackball** to position the measurement end point.



Press **Set** to complete the measurement.

The depth difference, time interval and slope between the 2 end points is displayed at the bottom of the screen.



## *Mode Measurements*

---

This page left blank intentionally.



# General Calculations

## Overview

The General Calculations Sub-Menu can be invoked from the Rad/Abdomen exam category. The calculations available are:

|        |            |  |  |
|--------|------------|--|--|
|        |            |  |  |
| VOLUME | %<br>STENO |  |  |
| ANGLE  |            |  |  |

|              |              |               |                 |
|--------------|--------------|---------------|-----------------|
|              |              |               |                 |
| S/D<br>Ratio | A/B<br>Ratio | HEART<br>RATE | TRACE<br>AUTO   |
| PI           | RI           |               | Transf<br>CALCS |

Transf CALCS is only available with Realtime Doppler Calculation option installed

|         |    |               |               |
|---------|----|---------------|---------------|
|         |    |               |               |
| Max PG  | CO | HEART<br>RATE | TRACE<br>AUTO |
| Mean PG | SV | FV            | FVO           |

Figure 8–1. General Calculations Sub-Menus

## Volume

Refer to the procedure found on 10–19.



---

## General Calculations

---

### Angle

Refer to the procedure found on 10–22.

### Stenosis Ratio (% stenosis)

Refer to the procedure found on 10–23.

### S/D Ratio, RI, A/B Ratio or PI

Refer to the procedure found on 10–29.

### Heart Rate

Refer to the procedure found on 10–30.

### Trace Auto

Use the Trace Auto Sub-Menu selection to turn on/off the Doppler Auto Trace function. See 10–28 for more details.

### Transf Calcs

This selection is only available if the Realtime Doppler Calculation option is installed and Realtime Doppler Calculations are displayed in the upper left corner of the screen.

Use the Trans Calcs selection to turn on/off the transfer of Automatic Doppler Calculations displayed on the screen to the report page. See 7–19 for details on Realtime Doppler operation.

### Max PG

Refer to the procedure found on 10–26.

### Mean PG

Refer to the procedure found on 10–27.



## General Calculations

---

### Cardiac Output (CO)

To measure CO (Cardiac Output), a velocity measurement is taken in the Doppler Spectrum. An FCA (Flow Cross Sectional Area) is measured on the vessel in B-Mode. These two measurements are used to calculate SV (Stroke Volume). Finally, a HR (Heart Rate) measurement is taken in the Doppler Spectrum. SV and HR are then used to calculate cardiac output.

|         |    |            |            |
|---------|----|------------|------------|
|         |    |            |            |
| Max PG  | CO | HEART RATE | TRACE AUTO |
| Mean PG | SV | FV         | FVO        |

Figure 8–2. General Calculations Sub-Menu (CO)

Select **CO** from the Calculation Sub-Menu. The trace cursor (horizontal dotted line) appears in the Doppler spectrum.

See Stroke Volume Ratio on 8–5 for a step-by-step procedure on SV. After the stroke volume measurements have been completed, a vertical line cursor appears in the Doppler spectrum. Measure the Heart Rate. Refer to the procedure found on 10–30.

CO (Cardiac Output) is computed from the SV and HR values and is displayed.

### Automatic CO calculation

If the following measurements have previously been made in any order, CO (Cardiac Output) automatically calculates when it is selected from the Calculation Sub-Menu:

- Velocity in Doppler Mode
- Functional Cross Sectional Area in B-Mode (Circumference/Area)
- Heart Rate in Doppler Mode



### Stroke Volume Ratio (SV)


To measure stroke volume, a velocity measurement is taken in the Doppler spectrum. An FCA (Flow Cross Sectional Area) is measured on the vessel in B-Mode. The velocity and FCA values are then used to calculate SV (Stroke Volume).

|         |    |            |            |
|---------|----|------------|------------|
|         |    |            |            |
| Max PG  | CO | HEART RATE | TRACE AUTO |
| Mean PG | SV | FV         | FVO        |

Figure 8–3. General Calculations Sub-Menu (SV)

Select **SV** from the Calculation Sub-Menu. The trace cursor (horizontal dotted line) appears in the Doppler spectrum.

Move the cursor to the measurement start point dependent on the trace method chosen (peak, floor, mean or mode).




Press **Set**. The first cursor is fixed and a second is active.

If Trace Auto is enabled, use the **Trackball** to move the second cursor to the velocity measurement end point.

Press **Set**. The waveform maximum velocities are automatically traced and a velocity value is displayed.

If Trace Auto is disabled, use the **Trackball** to trace the waveform.

Press **Set**. The trace waveform is fixed and the velocity value is displayed.



(continued)






## General Calculations

---

### Stroke Volume Ratio (SV) (cont'd)

A cross-point cursor appears in the B-Mode image to measure FCA.



Use the **Trackball** to move the cursor to a point on the vessel wall.

Press **Set** to fix the start point cursor. Use the Ellipse or Trace method to measure the circumference and area of the vessel as described on 7–8 or 7–10, respectively.

Press **Set** to complete the circumference/area measurement. The FCA (Flow Cross Sectional Area) is displayed.

The SV (Stroke Volume) is also computed and displayed from velocity and FCA values.

### Automatic SV calculation

If velocity was previously measured while in Doppler Mode and the cross sectional area of a vessel was measured while in B-Mode, stroke volume will be automatically calculated from these previously measured values when SV is selected from the General Calculation Sub-Menu.

### Heart Rate (HR)

Refer to the procedure found on 10–30.



## Flow Volume (FV)

Flow Volume estimates the volume of blood that flows through a vessel per unit time. It is derived from a vessel's cross-sectional diameter obtained from the B-Mode portion of the image and the mean velocity of flow in the vessel obtained from the Doppler portion of the image. It is measured in milliliters. When the FV measurement is made, FVO is automatically calculated.

To measure flow volume:

|         |    |            |            |
|---------|----|------------|------------|
|         |    |            |            |
| Max PG  | CO | HEART RATE | TRACE AUTO |
| Mean PG | SV | FV         | FVO        |

Figure 8–4. General Calculations Sub Menu (FV)

Select **FV** from the Calculation Sub Menu. Place the dotted horizontal line cursor at each of the timebase on the Doppler spectrum.

If Trace Auto is selected, the waveform is automatically traced.

If Trace Auto is not selected, manually trace the desired portion of the waveform.

The cursor moves to the B-Mode area.

Use the Ellipse or Trace method to measure the circumference and area of the vessel as described on 7–8 or 7–10, respectively.

The flow volume (FV) is calculated and displayed in milliliters. The flow volume output (FVO) is also calculated and displayed in milliliters/minute.



## General Calculations

---

### Trace Auto

Use the Trace Auto Sub-Menu selection to turn on/off the Doppler Auto Trace function. See 10–28 for more details.

### Flow Volume Output (FVO)

This measurement is used to measure the flow volume output in a vessel on the Doppler spectrum. It is measured in milliliters/minute. When the FVO measurement is made, FV is automatically calculated.

Refer to the FV measurement on 8–7 for details.

### Helpful hints



#### Hints

The following hints can help when taking a measurement:

- If **PRINT** is pressed while making a measurement, the system completes the measurement and sends the data to the report page (unless the VCR is assigned to the **PRINT** key).
- Prior to making measurements, use the Cine function, if necessary, to display the best image.
- Measurements can continue to be made until all measurement/calculation cells are filled. The cells are displayed at the bottom center of the display.
- Once all measurement/calculation cells are filled on the display, any further measurements will cause the top (first) cell to be erased and the new measurement added last (“first in, first out”).
- If there is a series of measurements/calculations that are routinely taken in a specific sequence, the auto sequence feature of the system can be used to program these measurements/calculations to automatically appear on the display in the desired order. Refer to 14–81 for details on programming an auto sequence.

## Hip Dysplasia Measurement

A HIP Dysplasia measurement is available. The HIP calculation assists in assessing the development of the infant hip. In this calculation three straight lines are superimposed on the image and aligned with the anatomical features. The 2 angles are computed, displayed and can be used by the physician in making a diagnosis.

The three lines are (Source: R GRAF, journal of Pediatric Orthopedics, 4: 735–740(1984)).

1. The inclination line connects the osseous convexity to labrum acetabulare.
2. The Acetabulum roof line connects the lower edge of the osilium to the osseous convexity.
3. The Baseline connects the osseous acetabulum convexity to the point where the joint capsule and the perichondrium unite with the iliac bone.

The Angle  $\alpha$  (Alpha) is the supplement of the angle between 1 and 3. It characterizes the osseous convexity. The angle  $\beta$  (Beta) is the angle between lines 1 and 2. It characterizes the bone supplementing additional roofing by the cartilaginous convexity.



*NOTE: There are two methods (Cranial–left and Caudal–left) for the Hip Dysplasia Measurement Orientation. This preset, “Hip Orientation”, is located on Preset Program page 4.*




Figure 3–44. Hip Dysplasia



## **HIP Dysplasia Measurement (cont'd)**

To activate the Hip Dysplasia Measurement:

1. To add the HIP Measurement, select "*Measurement Menu B*" in Preset Program page 7. If necessary, select NEXT CALC until HIP DYSPLASIA is displayed. Select HIP DYSPLASIA.
2. Press **Measurement**. Select the HIP measurement on B-Mode Measurement soft-menu. The system prompts the user in the measurement process.
3. To set the 3 lines for the Hip measurement, use the **Trackball**, **Zoom Size/Rotation** knob and **Set** key. The hip measurements ( $\alpha$  and  $\beta$ ) are displayed on the bottom of the screen.

## General Calculation Formulas

| Calc Mnemonic | Calc Name                                                       | Input Measurements                                                                            | Formula                                                                                             |
|---------------|-----------------------------------------------------------------|-----------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|
| MaxPG         | Maximum Pressure Gradient                                       | two Doppler blood flow peak velocities                                                        | $\text{MaxPG}[\text{mmHg}] = 4 \times (v_1^2 - v_2^2)$                                              |
| MeanPG        | Mean Pressure Gradient                                          | flow velocities from one time marker to another time marker in a Doppler display              | $\text{MeanPG}[\text{mmHg}] = \frac{4 \times \sum_{i=1}^n (V_i^2)}{n}$                              |
| SV            | Stroke Volume                                                   | flow velocities from one time marker to another time marker in a Doppler display and one area | $\text{SV}[\text{ml}] = (\sum \{V_i\} \text{ from } t_1 \text{ to } t_2) \times a_1$                |
| CO            | Cardiac Output                                                  | two distances and one 2 beat time interval                                                    | $\text{CO}[\text{l/min}] = \text{SV} \times \text{HR} / 1000$                                       |
| FV            | Flow Volume                                                     | flow velocities from one time marker to another time marker in a Doppler display and one area | $\text{FV}[\text{ml}] = (\sum \{V_i\} \text{ from } t_1 \text{ to } t_2) \times a_1 [\text{cm}^2]$  |
| FVO           | Flow Volume Output                                              | flow velocities from one time marker to another time marker in a Doppler display and one area | $\text{FVO}[\text{ml/min}] = \text{FV}[\text{ml}] / (t_2 - t_1) \times 60000$                       |
| % Stenosis    | Stenosis Ratio                                                  | two areas (by ellipse, trace, circle or distance)                                             | $\% \text{ Stenosis} = [1 - (A_{\text{residual}} / A_{\text{lumen}})] \times 100$                   |
| PI            | Pulsatility Index                                               | two Doppler blood flow peak velocities and TAMAX                                              | $\text{PI} = (V_{\text{max}} - V_{\text{diastole}}) / \text{TAMAX}$<br>*NOTE*                       |
| RI            | Resistivity Index                                               | two Doppler blood flow peak velocities                                                        | $\text{RI} = (V_{\text{max}} - V_{\text{diastole}}) / V_{\text{max}}$<br>*NOTE*                     |
| HR            | Heart Rate (beats/minute)                                       | one 2 beat time interval                                                                      | $\text{HR}[\text{BPM}] = 120[\text{sec}] / 2 \text{ beat time}[\text{sec}]$                         |
| A/B Ratio     | Velocities Ratio                                                | two Doppler blood flow peak velocities                                                        | $A/B = V_1 / V_2$                                                                                   |
| TAMAX         | Time Averaged Maximum Velocity (Trace Method is Peak or manual) | two time marks in a Doppler display                                                           | $\text{TAMAX} = \sum \{V_i\} \text{ from } t_1 \text{ to } t_2 / (t_2 - t_1) [\text{cm/s or m/s}]$  |
| TAMIN         | Time Averaged Minimum Velocity (Trace method is Floor)          | two time marks in a Doppler display                                                           | $\text{TAMIN} = \sum \{V_i\} \text{ from } t_1 \text{ to } t_2 / (t_2 - t_1) [\text{cm/s or m/s}]$  |
| TAMEAN        | Time Averaged Mean Velocity (Trace method is Mean)              | two time marks in a Doppler display                                                           | $\text{TAMEAN} = \sum \{V_i\} \text{ from } t_1 \text{ to } t_2 / (t_2 - t_1) [\text{cm/s or m/s}]$ |
| TAMODE        | Time Averaged Mode Velocity (Trace method is Mode)              | two time marks in a Doppler display                                                           | $\text{TAMODE} = \sum \{V_i\} \text{ from } t_1 \text{ to } t_2 / (t_2 - t_1) [\text{cm/s or m/s}]$ |

\*NOTE\*:  $V_{\text{diastole}} = V_{\text{min}}$  or  $V_{\text{end-diastole}}$  (depends on preset selection)

Table 8–1. General Calculation Formulas



## *General Calculations*

---

This page left blank intentionally.



# Exam Preparation

## Overview

Prior to an ultrasound examination, the patient should be informed of the clinical indication, specific benefits, potential risks, and alternatives, if any. In addition, if the patient requests information about the exposure time and intensity, it should be provided. Patient access to educational materials regarding ultrasound is strongly encouraged to supplement the information communicated directly to the patient. Furthermore, these examinations should be conducted in a manner and take place in a setting which assures patient dignity and privacy.

- Prior material knowledge and approval of the presence of nonessential personnel with the number of such personnel kept to a minimum.
- An intent to share with the parents, either during the examination or shortly thereafter, the information derived.
- An offer of choice about viewing the fetus.
- An offer of choice about learning the sex of the fetus, if such information becomes available.

Ultrasound examinations performed solely to satisfy the family's desire to know the fetal sex, to view the fetus or to obtain a picture of the fetus should be discouraged.





# Fetal Doppler

## Doppler Mode for Fetal Exams

### Indications for Fetal Doppler use

The LOGIQ™ 500 system can be used for fetal examination in Pulsed Wave Doppler, Continuous Wave Doppler, Color Flow Doppler and Color M-Mode for the diagnosis of:

- Structural fetal cardiac anomalies (18 weeks gestation to term) for high-risk patients.
- Intrauterine growth retardation (IUGR) (28 weeks gestation to term) for high-risk patients with one or more of the following known or suspected conditions:
  - Multiple pregnancy
  - Maternal hypertension
  - Hydrops
  - Diabetes
  - Lupus
  - Placenta abnormality.

### Contraindications

Pulsed Wave Doppler, Continuous Wave Doppler, Color Flow Doppler and Color M-Mode are not intended for routine fetal examination or screening nor are they intended for fetal examination in a low-risk population. The use of Doppler, even at minimal output levels, in fetal examination must be adjunctive with conventional fetal echocardiography and other clinical diagnostic methods, for high risk patients only.

# Acoustic Output

## Considerations

### General warning

#### CAUTION



The LOGIQ™ 500 system is a multi-use device which is capable of exceeding FDA Pre-enactment acoustic output (spatial peak-temporal average) intensity limits for fetal applications.

### Prudent use

It is prudent to conduct an exam with the minimum amount and duration of acoustic output necessary to optimize the image's diagnostic value.

### Concerns surrounding fetal exposure

Always be aware of the acoustic output level by observing the Acoustic Output Display. In addition, become thoroughly familiar with the Acoustic Output Display and equipment controls affecting output.

### Training

It is recommended that all users receive proper training in fetal Doppler applications before performing them in a clinical setting. Please contact a local sales representative for training assistance.



# OB Measurements and Formulas

## Introduction

Measurements and calculations derived from ultrasound images are intended to supplement other clinical procedures available to the attending physician. The accuracy of measurements is not only determined by the system accuracy, but also by use of proper medical protocols by the user. When appropriate, be sure to note any protocols associated with a particular measurement or calculation. Formulas and databases used within the system software that are associated with specific investigators are so noted. Be sure to refer to the original article describing the investigator's recommended clinical procedures.

If obstetrics was selected as the exam category on the Patient Entry Menu, the OB Menu automatically appears in the Soft-Menu when the **Measurement** key is pressed.

## OB Format Selection

Four types of OB measurements and report pages can be selected from the Set Up/System Parameter Menu page 2.

The user may choose to measure and report items based on Tokyo University, Osaka University, USA or European methods.

## OB Measurement Soft Menus and Formulas

The factory default measurement menus are shown in Figure 9–1 through Figure 9–4. The user should be familiar with how to do basic distance and circumference measurements, as described in the General Measurements and Calculations chapter.

### Tokyo University Method

|     |     |            |           |
|-----|-----|------------|-----------|
|     |     |            |           |
| GS  | BPD | LV         | EFBW      |
| CRL | FL  | APTD & TTD | OB Report |

Page 1

|      |       |       |           |
|------|-------|-------|-----------|
|      |       |       |           |
| CTAR | USER1 | USER3 | OB Graph  |
| EF   | USER2 | USER4 | OB Report |

Page 2

|               |           |           |              |
|---------------|-----------|-----------|--------------|
|               |           |           |              |
| Select Locatn | S/D Ratio | A/B Ratio | TRACE AUTO   |
| 1             | RI        | PI        | Transf CALCS |

Page 3

|      |       |                |           |
|------|-------|----------------|-----------|
|      |       |                |           |
| Vmax | USER5 | HEART RATE     | OB Graph  |
| PLI  |       | ANATOMY SURVEY | OB Report |

Page 4

Transf CALCS is only available with Realtime Doppler Calculation option installed

Figure 9–1. Tokyo University Sub-Menus



**NOTE:** Tokyo EFBW selection is made in the Set Up/System Parameters menu page 2.



## OB Measurements and Formulas

### Tokyo University Method (cont'd)

| Calc Mnemonic | Calc Name                            | Input Measurements                                                                     | Formula                                                                           | Author Reference                                  |
|---------------|--------------------------------------|----------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|---------------------------------------------------|
| HR            | Heart Rate (beats/minute)            | one 2 beat time (measured manually or automatically)                                   | $HR[BPM] = \frac{120[sec]}{2beat\ time[sec]}$                                     | n/a                                               |
| GS            | Gestational Sac                      | one distance                                                                           | Curve data is available in the Advanced Reference Manual                          | Tokyo University Method 1986, 6 by Univ. of Tokyo |
| CRL           | Crown Rump Length                    | one distance                                                                           |                                                                                   |                                                   |
| FL            | Femur Length                         | one distance                                                                           |                                                                                   |                                                   |
| BPD           | Biparietal Diameter                  | one distance                                                                           |                                                                                   |                                                   |
| LV            | Length of Vertebra                   | one distance                                                                           |                                                                                   |                                                   |
| APTD          | Antero posterior Trunk Diameter      | one distance                                                                           | APTD = input                                                                      | Y.Chiba Fetal Diagn Ther 1990 : 5 : 175 – 188     |
| TTD           | Transverse Trunk Diameter            | one distance                                                                           | TTD = input                                                                       |                                                   |
| EFBW          | Estimated Fetal Body Weight          | Average of BPD, Average of APTD, Average of TTD and Average of FL                      | $EFBW = 1.07 \times BPD^3 + 3.42 \times APTD \times TTD \times FL$                |                                                   |
| CTAR          | Cardio–Thoracic Area Ratio           | four distances or two areas                                                            | $CTAR = \frac{(d1 \times d2)}{(d3 \times d4)} \times 100$ or $a1 / a2 \times 100$ |                                                   |
| PLI           | Preload Index                        | two Doppler blood flow velocities                                                      | $PLI =  v1 / v2 $                                                                 |                                                   |
| EF            | Ejection Fraction                    | two distances on M-Mode (End-diastolic dimension and End-systolic dimension on M-Mode) | $EF = (1 - Ds^3 / Dd^3)$                                                          | n/a                                               |
| Vmax          | Maximum Velocity in descending aorta | one Doppler blood flow peak velocities                                                 | Vmax = input                                                                      | n/a                                               |
| S/D Ratio     | Systolic / Diastolic Ratio           | two Doppler blood flow peak velocities                                                 | $S/D = V_{systole} / V_{diastole}$ *NOTE*                                         | n/a                                               |
| A/B Ratio     | Velocities Ratio                     | two Doppler blood flow peak velocities                                                 | $A/B = V_1 / V_2$                                                                 | n/a                                               |
| PI            | Pulsatility Index                    | two Doppler blood flow peak velocities and TA_MAX                                      | $PI = (V_{max} - V_{diastole}) / TA\_MAX$ *NOTE*                                  | n/a                                               |
| RI            | Resistivity Index                    | two Doppler blood flow peak velocities                                                 | $RI = (V_{max} - V_{diastole}) / V_{max}$ *NOTE*                                  | n/a                                               |

\*NOTE\*:  $V_{diastole} = V_{min}$  or  $V_{end-diastole}$  (depends on preset selection)

Table 9–1. OB Calculation Formulas—Tokyo University Factory Defaults



*OB Measurements and Formulas*

**Osaka University Method**

|     |     |       |           |
|-----|-----|-------|-----------|
|     |     |       |           |
| CRL | FTA | HL    | EFBW      |
| BPD | FL  | USER1 | OB Report |

Page 1

|      |       |       |           |
|------|-------|-------|-----------|
|      |       |       |           |
| CTAR | USER2 | USER4 | OB Graph  |
| EF   | USER3 | USER5 | OB Report |

Page 2

|                                                                                     |           |           |              |
|-------------------------------------------------------------------------------------|-----------|-----------|--------------|
|                                                                                     |           |           |              |
| Select Locatn                                                                       | S/D Ratio | A/B Ratio | TRACE AUTO   |
|  1 | RI        | PI        | Transf CALCS |

Page 3

|      |  |                |           |
|------|--|----------------|-----------|
|      |  |                |           |
| Vmax |  | HEART RATE     | OB Graph  |
| PLI  |  | ANATOMY SURVEY | OB Report |

Page 4

Transf CALCS is only available with Realtime Doppler Calculation option installed

Figure 9–2. Osaka University Sub-Menus



## OB Measurements and Formulas

### Osaka University Method (cont'd)

| Calc Mnemonic | Calc Name                            | Input Measurements                                                                     | Formula                                                                                                                                                                                 | Author Reference                                  |
|---------------|--------------------------------------|----------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------|
| HR            | Heart Rate (beats/minute)            | one 2 beat time (measured manually or automatically)                                   | $HR[BPM] = 120[sec] / 2beat\ time[sec]$                                                                                                                                                 | n/a                                               |
| CRL           | Crown Rump Length                    | one distance                                                                           | Curve data is available in the Advanced Reference Manual                                                                                                                                | Osaka University Method 1989, 3 by Univ. of Osaka |
| FL            | Femur Length                         | one distance                                                                           |                                                                                                                                                                                         |                                                   |
| BPD           | Biparietal Diameter                  | one distance                                                                           |                                                                                                                                                                                         |                                                   |
| HL            | Humerus Length                       | one distance                                                                           |                                                                                                                                                                                         |                                                   |
| FTA           | Fetal Trunk Cross Sectional Area     | one area                                                                               | $APTD = input$                                                                                                                                                                          |                                                   |
| EFBW          | Estimated Fetal Body Weight          | Average of BPD, Average of FTA and Average of FL                                       | $EFBW = 1.25647 \times BPD^3 + 3.50655 \times FTA \times FL + 6.3$ (<5000g)<br>$IUGR = 1.229 \times BPD^3 + 3.063 \times FTA \times FL - 24.6$<br>Curve data is available in Appendix F |                                                   |
| CTAR          | Cardio–Thoracic Area Ratio           | four distances or two areas                                                            | $CTAR = (d1 \times d2) / (d3 \times d4) \times 100$<br>or $a1 / a2 \times 100$                                                                                                          | Y.Chiba Fetal Diagn Ther 1990 : 5 : 175 – 188     |
| PLI           | Preload Index                        | two Doppler blood flow velocities                                                      | $PLI =  v1 / v2 $                                                                                                                                                                       | Y.Chiba Fetal Diagn Ther 1990 : 5 : 168 – 174     |
| EF            | Ejection Fraction                    | two distances on M–Mode (End–diastolic dimension and End–systolic dimension on M–Mode) | $EF = (1 - Ds^3 / Dd^3)$                                                                                                                                                                | n/a                                               |
| Vmax          | Maximum Velocity in descending aorta | one Doppler blood flow peak velocities                                                 | $Vmax = input$                                                                                                                                                                          | n/a                                               |
| S/D Ratio     | Systolic / Diastolic Ratio           | two Doppler blood flow peak velocities                                                 | $S/D = Vsystole / Vdiastole$<br><i>*NOTE*</i>                                                                                                                                           | n/a                                               |
| A/B Ratio     | Velocities Ratio                     | two Doppler blood flow peak velocities                                                 | $A/B = V_1 / V_2$                                                                                                                                                                       | n/a                                               |
| PI            | Pulsatility Index                    | two Doppler blood flow peak velocities and TA_MAX                                      | $PI = (Vmax - V_{diastole}) / TA\_MAX$<br><i>*NOTE*</i>                                                                                                                                 | n/a                                               |
| RI            | Resistivity Index                    | two Doppler blood flow peak velocities                                                 | $RI = (Vmax - V_{diastole}) / Vmax$<br><i>*NOTE*</i>                                                                                                                                    | n/a                                               |

*\*NOTE\**:  $V_{diastole} = V_{min}$  or  $V_{end-diastole}$  (depends on preset selection)

Table 9–2. OB Calculation Formulas—Osaka University Factory Defaults

## OB Measurements and Formulas

### USA Method

|     |     |    |           |
|-----|-----|----|-----------|
|     |     |    |           |
| GS  | BPD | AC | USER1     |
| CRL | FL  | HC | OB Report |

Page 1

|       |       |                |           |
|-------|-------|----------------|-----------|
|       |       |                |           |
| USER2 | USER4 | AFI            | OB Graph  |
| USER3 | USER5 | ANATOMY SURVEY | OB Report |

Page 2

|               |           |           |              |
|---------------|-----------|-----------|--------------|
|               |           |           |              |
| Select Locatn | S/D Ratio | A/B Ratio | TRACE AUTO   |
| 1             | RI        | PI        | Transf CALCS |

Page 3

|  |  |            |           |
|--|--|------------|-----------|
|  |  |            |           |
|  |  | HEART RATE | OB Graph  |
|  |  |            | OB Report |

Page 4

Transf CALCS is only available with Realtime Doppler Calculation option installed

Figure 9–3. USA Sub-Menus

| Calc Mnemonic | Calc Name                 | Input Measurements                                       | Formula                                                                 | Author Reference                            |
|---------------|---------------------------|----------------------------------------------------------|-------------------------------------------------------------------------|---------------------------------------------|
| HR            | Heart Rate (beats/minute) | one 2 beat time (Measured manually or automatically)     | HR[BPM]= 120[sec]/2 beat time[sec]                                      | n/a                                         |
| GS            | Gestational Sac           | three distances                                          | $GS=(d1+d2+d3)/3+3.6225$                                                | Hellman, A/OG, 103: 789,1969                |
| CRL           | Crown Rump Length         | one distance                                             | $CRL=1.684969+0.315646xd1+0.049306xd1^2+0.004057xd1^3+0.000120456xd1^4$ | Hadlock, Radiology, 182:501, 1992           |
| FL            | Femur Length              | one distance                                             | $FL = -3.91 + 0.427 (GA) - 0.0034 (GA)^2$                               | Hadlock et al, Radiology, 152:497–501, 1984 |
| BPD           | Biparietal Diameter       | one distance                                             | $BPD = -3.08+0.41 (GA) - 0.000061 (GA)^3$                               |                                             |
| AC            | Abdominal Circumference   | circumference by trace, ellipse, circle or two distances | $AC = -13.3 + 1.61 (GA) - 0.00998 (GA)^2$                               |                                             |

Table 9–3. OB Calculation Formulas Part 1—USA Version Factory Defaults





## OB Measurements and Formulas

### USA Method (cont'd)

| Calc Mnemonic | Calc Name                    | Input Measurements                                             | Formula                                                                                                           | Author Reference                            |
|---------------|------------------------------|----------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------|---------------------------------------------|
| HC            | Head Circumference           | circumference by trace, ellipse, circle or two distances       | $HC = -11.48 + 1.56(GA) - 0.0002548(GA)^3$                                                                        | Hadlock et al, Radiology, 152:497-501, 1984 |
| OFD           | Occipito Frontal Diameter    | HC by ellipse or two distance                                  | Longest axis from HC measurement                                                                                  | n/a                                         |
| S/D Ratio     | Systolic/Diastolic Ratio     | two Doppler blood flow peak velocities                         | $S/D = V_{\text{systole}}/V_{\text{diastole}}$<br>*NOTE*                                                          | n/a                                         |
| A/B Ratio     | Velocities Ratio             | two Doppler blood flow peak velocities                         | $A/B = V_1/V_2$                                                                                                   | n/a                                         |
| PI            | Pulsatility Index            | two Doppler blood flow peak velocities and TAMAX               | $PI = (V_{\text{max}} - V_{\text{diastole}})/TAMAX$<br>*NOTE*                                                     | n/a                                         |
| RI            | Resistivity Index            | two Doppler blood flow peak velocities                         | $RI = (V_{\text{max}} - V_{\text{diastole}})/V_{\text{max}}$<br>*NOTE*                                            | n/a                                         |
| EFW #1        | Estimated Fetal Weight       | Average of BPD and Average of AC                               | $EFW = 10^{(-1.7492 + 0.166 * BPD + (0.046 * AC) - ((2.646 * AC * BPD)/1000))}$                                   | Shepard, AJOG, 142:47, 1982                 |
| EFW #2        | Estimated Fetal Weight       | Average of FL and Average of AC                                | $EFW = 10^{(1.304 + (0.05281 * AC) + (0.1938 * FL) - (0.004 * AC * FL))}$                                         | Hadlock-Radiology 150:535, 1984             |
| EFW #3        | Estimated Fetal Weight       | Average of BPD, Average of AC and Average of FL                | $EFW = 10^{(1.335 - 0.0034 * AC * FL) + (0.0316 * BPD) + 0.0457 * AC + (0.158 * FL)}$                             | Hadlock, AJOG, 151:333, 1985                |
| EFW #4        | Estimated Fetal Weight       | Average of FL, Average of AC and Average of HC                 | $EFW = 10^{(1.326 - (0.00326 * AC * FL)) + (0.0170 * HC) + (0.0438 * AC) + (0.158 * FL)}$                         | Hadlock, AJOG, 151:333, 1985                |
| EFW #5        | Estimated Fetal Weight       | Average of FL, Average of AC, Average of HC and Average of BPD | $EFW = 10^{(1.3596 - (0.00386 * AC * FL)) + (0.0064 * HC) + (0.00061 * BPD * AC) + (0.0424 * AC) + (0.174 * FL)}$ | Hadlock, AJOG, 151:333, 1985                |
| EFW Graph     | Estimated Fetal Weight Graph | EFW Calculation                                                | $LOG\ EFW\ Mean\ (g) = 0.578 + (0.332 * MA) - (0.00354 * MA^2)$<br><br>$1SD = 0.127 * EFW\ Mean$                  | Hadlock, Radiology, 1991; 181               |

$EFW[g] = \exp(0.578 + 0.332 * MA - 0.00354 * MA^2)$  where MA is menstrual age [weeks]

\*NOTE\*:  $V_{\text{diastole}} = V_{\text{min}}$  or  $V_{\text{end-diastole}}$  (depends on preset selection)

Table 9–3. OB Calculation Formulas Part 1—USA Version Factory Defaults (cont'd)

USA Method (cont'd)

| Calc Mnemonic | Calc Name                | Formula                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | Author Reference                                                                                                                                                         |
|---------------|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| CUA**         | Composite Ultrasound Age | 1. $CUA(BPD)=9.54+1.482*(BPD)+0.1676*(BPD)^2$<br>2. $CUA(HC)=8.96+0.540*(HC)+0.0003*(HC)^3$<br>3. $CUA(AC)=8.14+0.753*(AC)+0.0036*(AC)^2$<br>4. $CUA(FL)=10.35+2.460*(FL)+0.170*(FL)^2$<br>5. $CUA(BPD,AC)=9.57+0.524*(AC)+0.1220*(BPD)^2$<br>6. $CUA(BPD,HC)=10.32+0.009*(HC)^2+1.3200*(BPD)+0.00012*(HC)^3$<br>7. $CUA(BPD,FL)=10.50+0.197*(BPD)*(FL)+0.9500*(FL)+0.7300*(BPD)$<br>8. $CUA(HC,AC)=10.31+0.012*(HC)^2+0.3850*(AC)$<br>9. $CUA(HC,FL)=11.19+0.070*(HC)*(FL)+0.2630*(HC)$<br>10. $CUA(AC,FL)=10.47+0.442*(AC)+0.3140*(FL)^2-0.0121*(FL)^3$<br>11. $CUA(BPD,AC,FL)=10.61+0.175*(BPD)*(FL)+0.2970*(AC)+0.7100*(FL)$<br>12. $CUA(HC,BPD,FL)=11.38+0.070*(HC)*(FL)+0.9800*(BPD)$<br>13. $CUA(HC,AC,FL)=10.33+0.031*(HC)*(FL)+0.3610*(HC)+0.0298*(AC)*(FL)$<br>14. $CUA(HC,AC,BPD)=10.58+0.005*(HC)^2+0.3635*(AC)+0.02864*(BPD)*(AC)$<br>15. $CUA(BPD,HC,AC,FL)=10.85+0.060*(HC)*(FL)+0.6700*(BPD)+0.1680*(AC)$<br>Unit: CUA (week), BPD (cm), HC (cm), AC (cm), FL (cm) | Hadlock, Radiology, 1984 152:497–501                                                                                                                                     |
| FL/AC         | FL/AC ratio              | FL/AC                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | Hadlock, AJR 141:979, 1983                                                                                                                                               |
| FL/HC         | FL/HC ratio              | FL/HC                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | Hadlock, JUM 3:439, 1984                                                                                                                                                 |
| FL/BPD        | FL/BPD ratio             | FL/BPD                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | Hadlock, AJOG 141:759, 1987                                                                                                                                              |
| CI            | Cephalic Index           | BPD/OFD short distance HC/long distance HC                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | Hadlock, AJR 137:83, 1981                                                                                                                                                |
| HC/AC         | HC/AC ratio              | HC/AC                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | Campbell, BRJ. OG 84:165, 1977                                                                                                                                           |
| AFI*          | Amniotic Fluid Index     | AFI=AFI1 (distance) + AFI2 (distance) + AFI3 (distance) + AFI4 (distance)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | Dr. Rutherford/Dr. Phelan, Obstetrics & Gynecology Volume 70, No. 3, Part 1, p. 353–6, Sept. 1987. Dr. C. C. Smith, The Female Patient, Volume 15, p. 85-97, March 1990. |

\* The skip function using the Ellipse rocker switch is not available for AFI.

\*\* Formulas are used only if Hadlock HC, FL, AC and BPD are used and the “CUA/AUA for Hadlock” preset, found on System Parameters page 3, is set to “CUA”. If other authors are used, CUA automatically changes to AUA and an average value is displayed.

Table 9–4. OB Calculation Formulas Part 2—USA Version Factory Defaults



## OB Measurements and Formulas

### European Method

|               |               |              |                |
|---------------|---------------|--------------|----------------|
|               |               |              |                |
| GS<br>TOKYO   | BPD<br>JEANTY | AC<br>JEANTY | TAD<br>ERIKSEN |
| CRL<br>JEANTY | FL<br>JEANTY  | HC<br>JEANTY | OB<br>Report   |

Page 1

|                |             |                   |              |
|----------------|-------------|-------------------|--------------|
|                |             |                   |              |
| BD<br>JEANTY   | Ft<br>PARIS |                   | OB<br>Graph  |
| OFD<br>HANSMAN | USER1       | ANATOMY<br>SURVEY | OB<br>Report |

Page 2

|                  |              |              |                 |
|------------------|--------------|--------------|-----------------|
|                  |              |              |                 |
| Select<br>Locatn | D/S<br>Ratio | A/B<br>Ratio | TRACE<br>AUTO   |
| ▲<br>▼<br>1      | RI           | PI           | Transf<br>CALCS |

Page 3

|       |       |               |              |
|-------|-------|---------------|--------------|
|       |       |               |              |
| USER2 | USER4 | HEART<br>RATE | OB<br>Graph  |
| USER3 | USER5 |               | OB<br>Report |

Page 4

Transf CALCS is only available with Realtime Doppler Calculation option installed

Figure 9–4. European Sub-Menus

| Calc Mnemonic | Calc Name                 | Input Measurements                                       | Formula                                                  | Author Reference                                  |
|---------------|---------------------------|----------------------------------------------------------|----------------------------------------------------------|---------------------------------------------------|
| HR            | Heart Rate (beats/minute) | one 2 beat time (Measured manually or automatically)     | $HR[BPM] = \frac{120[sec]}{2beat\ time[sec]}$            | n/a                                               |
| GS            | Gestational Sac           | one distance                                             | Curve data is available in the Advanced Reference Manual | Tokyo University Method 1986, 6 by Univ. of Tokyo |
| CRL           | Crown Rump Length         | one distance                                             |                                                          | JEANTY : Radiology, 143 : 513, 1982               |
| BPD           | Biparietal Diameter       | one distance                                             |                                                          |                                                   |
| HC            | Head Circumference        | circumference by trace, ellipse, circle or two distances |                                                          |                                                   |
| AC            | Abdominal Circumference   | circumference by trace, ellipse, circle or two distances |                                                          |                                                   |

Table 9–5. OB Calculation Formulas—European Version Factory Defaults

European Method (cont'd)

| Calc Mnemonic | Calc Name                     | Input Measurements                                | Formula                                                                  | Author Reference                                                      |
|---------------|-------------------------------|---------------------------------------------------|--------------------------------------------------------------------------|-----------------------------------------------------------------------|
| FL            | Femur Length                  | one distance                                      | Curve data is available in the Advanced Reference Manual                 | JEANTY : Radiology, 143 : 513, 1982                                   |
| BD            | Binocular Distance            | one distance                                      |                                                                          | ERIKSEN                                                               |
| TAD           | Transverse Abdominal Diameter | one distance                                      |                                                                          |                                                                       |
| OFD           | Occipito Frontal Diameter     | one distance                                      | Known LMP Case:                                                          | Hansmann: Ultraschalldiag-nostik in Begurtshilfe und Gynakologie 1985 |
|               |                               |                                                   | Unknown LMP Case:                                                        |                                                                       |
| Ft            | Foot Distance                 | one distance                                      | Curve data is available in the Advanced Reference Manual                 | PARIS                                                                 |
| EFW           | Estimated Fetal Weight        | Average of BPD and Average of AC [cm]             | $EFW[g] = 10^{(1.7288 + 0.09184*BPD + (0.02581*AC) + (0.00011*BPD*AC))}$ | Shepard/ Warsoff                                                      |
| FL/AC         | FL/AC ratio                   | Average of FL and Average of AC                   | FL/AC                                                                    | Hadlock, AJR, 141 : 979, 1983                                         |
| FL/BPD        | FL/BPD ratio                  | Average of FL and Average of BPD                  | FL/BPD                                                                   | Hadlock, AJOG, 141 : 759, 1987                                        |
| CI            | Cephalic Index                | Average of BPD and Average of OFD                 | BPD/OFD short distance HC/long distance HC                               | Hadlock, AJR, 137 : 83, 1981                                          |
| HC/AC         | HC/AC ratio                   | Average of HC and Average of AC                   | HC/AC                                                                    |                                                                       |
| FL/Ft         | FL/Ft ratio                   | Average of FL and Average of Ft                   | FL/Ft                                                                    |                                                                       |
| BD/BPD        | BD/BPD ratio                  | Average of BPD and Average of BD                  | BD/BPD                                                                   |                                                                       |
| D/S Ratio     | Diastolic/Systolic Ratio      | two Doppler blood flow peak velocities            | $D/S = V_{diastole} / V_{systole}$ *NOTE*                                | n/a                                                                   |
| A/B Ratio     | Velocities Ratio              | two Doppler blood flow peak velocities            | $A/B = V_1 / V_2$                                                        | n/a                                                                   |
| PI            | Pulsatility Index             | two Doppler blood flow peak velocities and TA_MAX | $PI = (V_{max} - V_{diastole}) / TA\_MAX$ *NOTE*                         | n/a                                                                   |
| RI            | Resistivity Index             | two Doppler blood flow peak velocities            | $RI = (V_{max} - V_{diastole}) / V_{max}$ *NOTE*                         | n/a                                                                   |

\*NOTE\*:  $V_{diastole} = V_{min}$  or  $V_{end-diastole}$  (depends on preset selection)

Table 9–5. OB Calculation Formulas—European Version Factory Defaults (cont'd)



## OB Measurements and Formulas

### Other OB Calculation Formulas

| Calc Mnemonic | Calc Name                     | Input Measurements                                       | Formula                                                  | Author Reference                                                            |
|---------------|-------------------------------|----------------------------------------------------------|----------------------------------------------------------|-----------------------------------------------------------------------------|
| BPD           | Biparietal Diameter           | one distance                                             | Curve data is available in the Advanced Reference Manual | CAMPBELL : King's College Hosp. London (Am.J. obst gynecol). Oct. 1, 1982   |
| FL            | Femur Length                  | one distance                                             |                                                          |                                                                             |
| BD            | Binocular Distance            | one distance                                             |                                                          |                                                                             |
| GS            | Gestational Sac               | one distance                                             |                                                          | HANSMANN: M and AI : Geburtsch, u, Frauenheilk 39 : 656, 1979               |
| CRL           | Crown Rump Length             | one distance                                             |                                                          | ROBINSON: Robins on and AI: Br J Gynecol, 82 : 702, 1975                    |
| CRL           | Crown Rump Length             | one distance                                             |                                                          | PARIS                                                                       |
| BPD           | Biparietal Diameter           | one distance                                             |                                                          |                                                                             |
| TAD           | Transverse Abdominal Diameter | one distance                                             |                                                          |                                                                             |
| FL            | Femur Length                  | one distance                                             |                                                          |                                                                             |
| BPD           | Biparietal Diameter           | one distance                                             |                                                          |                                                                             |
| HC            | Head Circumference            | circumference by trace, ellipse, circle or two distances |                                                          | SOSTOA: Hospital de la Santa Cruz y San Pablo, servieio de obst. y gynecol. |
| AC            | Abdominal Circumference       | circumference by trace, ellipse, circle or two distances |                                                          |                                                                             |
| FL            | Femur Length                  | one distance                                             |                                                          |                                                                             |
| BD            | Binocular Distance            | one distance                                             |                                                          |                                                                             |
| OFD           | Occipito Frontal Diameter     | one distance                                             |                                                          |                                                                             |
| FL            | Femur Length                  | one distance                                             |                                                          | Hadlock, AJR, 138 : 875, 1982                                               |

Table 9–6. OB Calculation Formulas—Other Selections

Other OB Calculation Formulas (cont'd)

| Calc Mnemonic | Calc Name               | Input Measurements                                       | Formula                                                                                                              | Author Reference                    |
|---------------|-------------------------|----------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|-------------------------------------|
| BPD           | Biparietal Diameter     | one distance                                             | Curve data is available in the Advanced Reference Manual                                                             | Hadlock, JUM 1:97, 1982             |
| AC            | Abdominal Circumference | circumference by trace, ellipse, circle or two distances |                                                                                                                      | Hadlock, AJR, 139 : 367, 1982       |
| HC            | Head Circumference      | circumference by trace, ellipse, circle or two distances |                                                                                                                      | Hadlock, AJR, 138 : 649, 1982       |
| CRL           | Crown Rump Length       | one distance                                             |                                                                                                                      | NELSON                              |
| BPD           | Biparietal Diameter     | one distance                                             |                                                                                                                      | KURTZ                               |
| BD            | Binocular Distance      | one distance                                             |                                                                                                                      | BERKOWITZ                           |
| EFW           | Estimated Fetal Weight  | Average of BPD and Average of TAD [mm]                   | $EFW[kg] = 0.515263 - (0.105775 * BPD) + (0.000930707 * (BPD)^2) + (0.0649145 * TAD) - (0.00020562 * (TAD[mm])^2)$   | German                              |
| EFW           | Estimated Fetal Weight  | Average of BPD and Average of AC [cm]                    | $EFW [g] = 10^{(3 - 1.7492 + 0.046 * AC + 0.166 * BPD - 0.002646 * AC * BPD)}$                                       | Shepard : Richards/ Berkowitz       |
| EFW           | Estimated Fetal Weight  | Average of FL, Average of AC and Average of HC [cm]      | $EFW [g] = 10^{(1.5662 - (0.0108 * HC) + (0.0468 * AC) + (0.171 * FL) + (0.00034 * (HC)^2) - (0.003685 * AC * FL))}$ | Hadlock, Radiology, 150 : 535, 1984 |
| FL/HC         | FL/HC ratio             | Average of FL and Average of HC                          | FL/HC                                                                                                                | Hadlock, JUM, 3 : 439, 1984         |

Table 9–6. OB Calculation Formulas—Other Selections (cont'd)



## OB Measurements and Formulas

### Other OB Calculation Formulas (cont'd)

| Calc Mnemonic    | Calc Name                                                        | Input Measurements               | Formula                                                          | Author Reference                                                      |
|------------------|------------------------------------------------------------------|----------------------------------|------------------------------------------------------------------|-----------------------------------------------------------------------|
| GS               | Gestational Sac                                                  | one distance                     | Known LMP Case:<br>Unknown LMP Case:                             | Rempen                                                                |
| CRL              | Crown Rump Length                                                | one distance                     | Known LMP Case:<br>Unknown LMP Case:                             |                                                                       |
| BPD              | Biparietal Diameter                                              | one distance                     | Known LMP Case:<br>Unknown LMP Case:                             |                                                                       |
| BPD              | Biparietal Diameter                                              | one distance                     | Known LMP Case:<br>Unknown LMP Case:                             |                                                                       |
| CRL              | Crown Rump Length                                                | one distance                     | Known LMP Case:<br>Unknown LMP Case:                             | Hansmann: Ultraschalldiag–nostik in Begurtshilfe und Gynakologie 1985 |
| HC               | Head Circumference                                               | Average of BIP<br>Average of FRO | HC [mm] = $2.325 * (BPD [mm]^2 + OFD [mm]^2)^{0.5}$              |                                                                       |
| ThD              | Transverse Thorax Diameter                                       | one distance                     | Known LMP Case:<br>Unknown LMP Case:                             |                                                                       |
| FL               | Femur Length                                                     | one distance                     | Known LMP Case:<br>Unknown LMP Case:                             |                                                                       |
| AC               | Abdominal Circumference                                          | one distance                     | Known LMP Case:<br>Unknown LMP Case:                             |                                                                       |
| EFW              | Estimated Fetal Weight                                           | Average of BPD<br>Average of AC  | EFW[g] = $-3200.40479 + 157.07186 * AC[cm] + 15.90391 * BPD[cm]$ | Merz, et al., Womens Hospital, University of Mainz                    |
| BPD              | Biparietal Diameter                                              | one distance                     | Curve data is available in the Advanced Reference Manual         | Tokyo Shinozuka Method, 1996                                          |
| APTD x TTD (AxT) | Antero Posterior Trunk Diameter<br><br>Transverse Trunk Diameter | two distance                     |                                                                  |                                                                       |

Table 9–6. OB Calculation Formulas—Other Selections (cont'd)

Other OB Calculation Formulas (cont'd)

| Calc Mnemonic | Calc Name               | Input Measurements                                      | Formula                                                                                                                                                | Author Reference                |
|---------------|-------------------------|---------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|
| AC            | Abdominal Circumference | circumference by trace, ellipse, circle or two distance | Curve data is available in the Advanced Reference Manual                                                                                               | Tokyo Shinozuka Method, 1996    |
| FL            | Femur Length            | one distance                                            |                                                                                                                                                        |                                 |
| CRL           | Crown Rump Length       | one distance                                            |                                                                                                                                                        |                                 |
| EFW1          | Estimated Fetal Weight  | four distance                                           | EFW1 [gram] = $1.07 * \text{BPD}^3 + 3.42 * \text{APTD} * \text{TTD} * \text{FL}$<br>G = [day], BPD = [cm], FL = [cm], APTD * TTD = [cm <sup>2</sup> ] |                                 |
| EFW2          | Estimated Fetal Weight  | 2 distance and 1 circumference                          | EFW2 [gram] = $1.07 * \text{BPD}^3 + 0.30 * \text{AC}^2 * \text{FL}$<br>G = [day], BPD = [cm], AC = [cm], FL = [cm]                                    |                                 |
| EFW3          | Estimated Fetal Weight  | four distance                                           | EFW3 [gram] = $1.07 * \text{BPD}^3 + 2.91 * \text{APTD} * \text{TTD} * \text{LV}$<br>G = [day], BPD = [cm], LV = [cm], APTD * TTD = [cm <sup>2</sup> ] |                                 |
| BPD           | Biparietal Diameter     | one distance                                            | Curve data is available in the Advanced Reference Manual                                                                                               | ASUM, D. Robinson et al, 1990.8 |
| CRL           | Crown Rump Length       | one distance                                            |                                                                                                                                                        | ASUM, Silva et al, 1991.6       |
| AC            | Abdominal Circumference | circumference by trace, ellipse, circle or two distance | AC (mm) = $-56.582 + 11.402 * (\text{GA} - 2)$<br>GA = [week]                                                                                          | ASUM, Deler                     |

Table 9–6. OB Calculation Formulas—Other Selections (cont'd)

GS Measurement

GS Measurement for the Tokyo University and European methods requires only one distance measurement.

GS Measurement for the USA method requires 3 distance measurements obtained from longitudinal and transverse images. The 3 measurements are anterior–posterior, transverse and longitudinal. It is advantageous to use the dual-image (split screen) feature for this GS measurement.





### Amniotic Fluid Index (AFI)

The Amniotic Fluid Index is determined by taking distance measurements of the 4 quadrants of the uterine cavity. These 4 measurements are added together to calculate the Amniotic Fluid Index.

Press the **Measurement** key once to display a “×” cursor on the screen, and to display distance in cm on the bottom part of screen. Select **AFI** from the OB Measurement Soft Menu. The distance measurement changes to AFI1 (2, 3 or 4).

Perform a standard distance measurement as described in 7–7.

When the measurement of the first quadrant is completed, unfreeze and move to the second quadrant. After the image is obtained, press **Freeze** and then **Measurement**. The user is then prompted to continue with the AFI measurements. Perform a standard distance measurement for the second quadrant.

Repeat these steps for quadrants 3 and 4.

When all 4 quadrants have been measured, the AFI will be calculated and displayed in the Calculation Result Area of the display.



### Hints

The skip function using the Ellipse rocker switch is not available for AFI.

Unfreezing the image after doing an AFI measurement does not delete the previous measurements. Unfreeze and change scan planes as necessary. Make the remaining measurements to arrive at the AFI value by following the basic distance measurement procedure.

If any generic distance measurements are unassigned before invoking the Amniotic Fluid Index these measurements will be used. Only the remaining distance measurements need to be made to arrive at an AFI.

---

## OB Measurements and Formulas

---

### Amniotic Fluid Index (cont'd)


The normal values are considered to be:

36–40 weeks  
0–5 cm = very low  
5.1–8.0 cm = low  
8.1 – 18.0 cm = normal  
>18.0 = high


Dr. Rutherford/Dr. Phelan, *Obstetrics and Gynecology*, Volume 70, No. 3, Part 1, p. 353–6, Sept. 1987.

28–40 weeks  
15.0 cm = average  
>20.0 – 24.0 = hydramnios  
<5.0–6.0 = Oligohydramnios

Dr. C.C. Smith, *The Female Patient*, Volume 15, p. 85–97, March 1990.



Sagittal



Transverse

Select **AFI** from the OB Measurement Soft Menu.



## OB Measurements and Formulas

---

### User-Programmed Calculations (Tables)

User-programmed calculation tables are displayed in the measurement menus only after they are entered by the user.

The User Table Editor is found in the Set Up/Preset Program Menu page 7. Refer to 14–77 on how to define these user selections and program the user tables.

### Select Locatn

This choice is used to select 1 of 3 locations at which S/D (D/S) Ratio, A/B Ratio, RI and PI can be taken and recorded on the report page.

The locations are numbered 1, 2 and 3 as a default from the factory. However, the user can customize each location designation with a 6 character name. This customization edit is done in any report page that displays S/D (D/S) Ratio, A/B Ratio, RI or PI. Changes will be permanently saved even after system power down.

Up to three measurements can be taken for S/D (D/S) Ratio, A/B Ratio, RI or PI. The latest value or measurement average can be reported. Averaging depends upon the average activity preset selection in the Set Up/Preset Program page 3.

Press the **Select Locatn** rocker switch to change the location number as necessary.

|                                                                                       |              |              |                 |
|---------------------------------------------------------------------------------------|--------------|--------------|-----------------|
|                                                                                       |              |              |                 |
| Select<br>Locatn                                                                      | D/S<br>Ratio | A/B<br>Ratio | TRACE<br>AUTO   |
|  1 | RI           | PI           | Transf<br>CALCS |

Transf CALCS is only available with Realtime Doppler Calculation option installed

Figure 9–5. European OB Calculation Sub-Menu



---

## OB Measurements and Formulas

---

### S/D (D/S) Ratio, A/B Ratio, RI or PI

Refer to the procedure found on 10–29.

### Trace Auto

Use the Trace Auto Sub-Menu selection to turn on/off the Doppler Auto Trace function. See 10–28 for more details.

### Transf Calcs

This selection is only available if the Realtime Doppler Calculation option is installed and Realtime Doppler Calculations are displayed in the upper left corner of the screen. See 7–19 for details on RT Doppler operation.

### Fetal Heart Rate Measurement

The fetal heart rate is typically taken in M-Mode, though it also can be performed in Doppler Mode. Refer to the procedure found on 10–30.

## Helpful Hints



### Hints

- Any measurement can be repeated by selecting that measurement again from the Soft-Menu.
- The system displays several measurements, but the summary report retains only the last 3 measurements of each type.
- The 3 report page measurements can be averaged and the average used in other calculations. This averaging selection is made in the Set Up/Preset Program page 3. Average activity on/off and average number (number of measurements averaged) can be set to the user's preference. Other presets on this page that affect OB measurements for the European method are Add 1 week to EDD and Display EDD with GA value.
- If results are displayed as question marks, the measurement is not within the table range.



# OB Summary Reports

## Starting an Exam

Accurate and complete Summary Report presentation starts at the beginning of the patient exam. Always begin an exam by entering as much new patient information as possible.



*NOTE: Refer to the Beginning an Exam section on 4–3 for instructions.*

After selecting the OB exam type, fill in the following information after PT Name, PT ID, Notes, and Oper ID:

1. **AGE.** Type the woman's age.
2. **PREGNANCY ORIGIN DATA SELECTION** Enter 1 for LMP, 2 for BBT (not in USA version), 3 for EDC (EDD in Europe) and 4 for GA.  
Type in the required information for category selected.
3. **GRAVIDA.** Type the number of pregnancies.
4. **PARA.** Type the number of live births.
5. **AB.** Type the number of abortions, miscarriages, etc.
6. **ECTOPIC.** Type the number of pregnancies outside the uterus.
7. **REF. MD.** Type the name of the woman's physician.
8. **COMMENT** Type any pertinent medical history with regard to this pregnancy, e.g., diabetes, bleeding.
9. **EXIT** Press **Return** after having completed filling in this information.



*NOTE: If patient information needs to be edited, use the **ID/Name** key to display the patient entry menu for editing information.*

## OB Report Page Layout

The basic OB Report format consists of 2 pages for each region available (USA, Europe, Tokyo University and Osaka University).

The first page consists of 4 basic areas. The first, or top portion, is generally patient data that was entered on the patient entry menu at the beginning of the exam. This area cannot be edited.

The second is the measurements area. It varies by region format. Generally this area contains measurement titles, the last 3 measurements made, the measurement average (if averaging is used), estimated gestational age and the ability to choose if the measurement is used in the calculation of the composite age. This area can be edited by deleting measurements, selecting yes/no for composite age or rearranging the vertical location of a measurement using the Report Format Editor in the Set Up/Preset Program page 6.

The third area is calculations. It varies by region format. These items are results of measurements or patient data entered on the Patient Entry Menu. This area cannot be edited.

The fourth area is comments. The same amount of area is available for user comments for each region. The user can edit this area freely. Comments made on page 1 are not reflected on page 2. Therefore, comments on each page can be edited individually.



## OB Summary Reports

### OB Report Page 1

|               |               |            |    |             |          |
|---------------|---------------|------------|----|-------------|----------|
| HOSPITAL NAME |               | Oper ID:   |    | 98/05/09    |          |
| ID:           | NAME:         |            |    | AGE: 0      |          |
| LMP (OPE)     | GA (LMP)      | EDD (LMP)  |    | G0 P0 A0 E0 |          |
| Ref MD:       | NOTE:         |            |    |             |          |
| POS:          | PLAC:         |            |    | ← 1/2       |          |
| MEASUREMENTS  | CUA           | LAST       | 1  | 2           | 3        |
| BPD (HADLOCK) | y             | (          |    |             | )        |
| HC (HADLOCK)  | y             | (          |    |             | )        |
| OFD (HC)      | y             | (          |    |             | )        |
| AC (HADLOCK)  | y             | (          |    |             | )        |
| FL (HADLOCK)  | y             | (          |    |             | )        |
| CRL (HADLOCK) | y             | (          |    |             | )        |
| GS (HELLMAN)  | y             | (          |    |             | )        |
| CALCULATIONS  |               |            |    |             |          |
| CI            | (70-86)       | EFW        | g± | g           | ( lb oz) |
| FL/BPD        | (71-87)       | Based On:  | (  |             | )        |
| FL/AC         | (20-24)       | AF1 (cm)   |    |             | HR (BPM) |
| FL/HC         | (??-??-??)    | LMP: (OPE) |    |             |          |
| HC/AC         | (?. ??-?. ??) | AGE: LMP   |    |             | CUA      |
|               |               | EDD: LMP   |    |             | CUA      |
| COMMENTS:     |               |            |    |             |          |

Figure 9–6. USA Basic OB Summary Report Page 1



*NOTE: When the head circumference (HC) is measured by the ellipse or 2 distance method, the longer axis will automatically be assigned as the OFD measurement.*

<AVG> or <LAST> indicates that the Average or Latest value will be displayed in this column. This depends on the preset, *Average Activity*, in the Set Up/Preset Program menu page 3.

|                   |        |             |   |          |       |
|-------------------|--------|-------------|---|----------|-------|
| HOSPITAL NAME     |        | Oper ID:    |   | 98/05/09 |       |
| ID:               | NAME:  |             |   | AGE: 0   |       |
| LMP:              | BBT:   | G0 P0 A0 E0 |   |          |       |
| Ref MD:           | NOTE:  |             |   |          |       |
| MEASUREMENTS      | AGE    | CGA         | 1 | 2        | ← 1/2 |
| GS (TOKYO)        |        | y           | ( |          | )     |
| CRL (JEANTY)      |        | y           | ( |          | )     |
| BPD (JEANTY)      |        | y           | ( |          | )     |
| HC (JEANTY)       |        | y           | ( |          | )     |
| FL (JEANTY)       |        | y           | ( |          | )     |
| TAD (ERIKSEN)     |        | y           | ( |          | )     |
| AC (JEANTY)       |        | y           | ( |          | )     |
| BD (JEANTY)       |        | y           | ( |          | )     |
| OFD (HANSMANN)    |        | y           | ( |          | )     |
| Ft (PARIS)        |        | y           | ( |          | )     |
| CALCULATIONS      |        |             |   |          |       |
| GA (LMP)          | CI     | HR          |   |          |       |
| EDD (LMP)         | FL/BPD | 1           | 2 | 3        |       |
| LMP (GA)          | FL/AC  | P1          |   |          |       |
| CGA               | HC/AC  | R1          |   |          |       |
| EDD (CGA)         | FL/Ft  | D/S         |   |          |       |
| EFW (Sheep/Wares) | BD/BPD |             |   |          |       |
| COMMENTS:         |        |             |   |          |       |

Figure 9–7. European Basic OB Report Page 1

OB Report Page 1 (cont'd)

|              |       |             |                 |
|--------------|-------|-------------|-----------------|
| ID:          | NAME: | Oper ID:    | 99/03/18        |
| LMP:         | BBT:  | EDC(LMP):   | AGE: 0          |
| Ref MD:      | NOTE: | G0 P0 A0 E0 |                 |
| MEASUREMENTS | GA    | EDC         | CGA 1 2 + 1/2 3 |
| GS(TOKYO)    |       |             | y ( )           |
| CRL(TOKYO)   |       |             | y ( )           |
| BPD(TOKYO)   |       |             | y ( )           |
| FL(TOKYO)    |       |             | y ( )           |
| LV(TOKYO)    |       |             | y ( )           |
| APTD(TOKYO)  |       |             | ( )             |
| TTD(TOKYO)   |       |             | ( )             |
| UT01(YA)     |       |             | y ( )           |
| UT10()       |       |             | y ( )           |
| CALCULATIONS |       |             |                 |
| GA(LMP)      | CTAR  | HR          | AFI 3           |
| EDC(LMP)     | EDD   | 1 2         |                 |
| CGA          | EF    | PI          |                 |
| EDC(CGA)     | FS    | RI          |                 |
| EFBW(TOKYO)  | Vmax  | S/D         |                 |
| COMMENTS:    | PLI   |             |                 |

Figure 9–8. Tokyo University Basic OB Report Page 1

|              |                |             |                 |
|--------------|----------------|-------------|-----------------|
| ID:          | NAME:          | Oper ID:    | 99/03/18        |
| LMP:         | BBT:           | EDC(LMP):   | AGE: 0          |
| Ref MD:      | NOTE:          | G0 P0 A0 E0 |                 |
| MEASUREMENTS | GA             | EDC         | CGA 1 2 + 1/2 3 |
| CRL(OSAKA)   |                |             | y ( )           |
| BPD(OSAKA)   |                |             | y ( )           |
| FTR(OSAKA)   |                |             | y ( )           |
| HL(OSAKA)    |                |             | y ( )           |
| FL(OSAKA)    |                |             | y ( )           |
| CALCULATIONS |                |             |                 |
| GA(LMP)      | CTAR           | HR          | AFI 3           |
| EDC(LMP)     | EDD            | 1 2         |                 |
| CGA          | EF             | PI          |                 |
| EDC(CGA)     | FS             | RI          |                 |
| GA(EFBW)     | Vmax           | S/D         |                 |
| GA(IUGR)     | PLI            |             |                 |
| EFBW(OSAKA)  | IU-EFBW(OSAKA) |             |                 |
| COMMENTS:    |                |             |                 |

Figure 9–9. Osaka University Basic OB Report Page 1





## OB Summary Reports

---

### Composite Age Determination

The composite fetal age is displayed as CUA (Composite Ultrasound Age) on the USA Report and CGA (Composite Gestational Age) on the European, Tokyo University and Osaka University Reports.



*NOTE: There is a preset selection to choose how gestational age is derived for the USA Hadlock formulas. The “CUA/AUA for Hadlock” preset, found on System Parameters page 3, allows gestational age to be derived by linear regression (CUA) or averaging (AUA).*

The CUA or CGA field can be edited on the report page by the designation ‘Y’ for yes or ‘N’ for no.

‘Y’ will include the last measurement or measurement average in the calculation composite age. Changing the CUA/CGA designation to ‘N’ will exclude the measurement from the composite age calculation. The default value is ‘Y’.



*NOTE: There is a preset selection on Set Up/System Parameters page one that allows for the resetting of the CUA fields at New Patient selection (Reset on or off).*

### Fetal Weight Estimates

The estimated fetal birth weight calculation varies by region. EFW or EFBW is calculated for:

USA Version—The system automatically selects 1 of the 5 available formulas based on the measurements taken.

Tokyo or European Version—It is based on the formula selected in the Set Up/System Parameter menu page 2.

Osaka University Version—It is based on a single formula listed on 9–10. Standard Deviation is also displayed with EFBW in the Osaka University OB Report.

OB Report Page 2

<AVG> or <LAST> indicates that the Average or Latest value will be displayed in this column. This depends on the preset, *Average Activity*, in the Set Up/Preset Program menu page 3.

SITE represents a Doppler site that the user can edit with a maximum of 6 characters. Factory default labels are 1, 2 and 3.

|               |           |               |   |             |   |
|---------------|-----------|---------------|---|-------------|---|
| HOSPITAL NAME |           | Oper ID:      |   | 98/05/09    |   |
| ID:           | NAME:     | AGE: 0        |   |             |   |
| LMP(OPE)      | GA(LMP)   | EDD(LMP)      |   | G0 P0 A0 E0 |   |
| Ref MD:       | POS:      | NOTE:         |   | + 2/2       |   |
|               | FLUID:    | PLAC:         |   |             |   |
| BIOPHYSICAL:  | MOVEMENT: | PREVIA?:      |   | GRADE:      |   |
|               | FLUID:    | TONE:         |   | BREATHING:  |   |
|               |           | REACTIVE NST: |   | SCORE: /10  |   |
| AFI(cm)       | SUM       | =             | 1 | 2           | 3 |
|               |           |               | + | +           | + |
|               |           |               |   |             |   |
| HR(BPM)       | SITE      | LAST          | 1 | 2           | 3 |
| PI            | 1         |               | ( |             | ) |
|               | 2         |               | ( |             | ) |
|               | 3         |               | ( |             | ) |
| RI            | 1         |               | ( |             | ) |
|               | 2         |               | ( |             | ) |
|               | 3         |               | ( |             | ) |
| S/D           | 1         |               | ( |             | ) |
|               | 2         |               | ( |             | ) |
|               | 3         |               | ( |             | ) |
| A/B           | 1         |               | ( |             | ) |
|               | 2         |               | ( |             | ) |
|               | 3         |               | ( |             | ) |
| COMMENTS:     |           |               |   |             |   |

Figure 9–10. USA Basic OB Summary Report Page 2

|               |           |               |   |            |   |
|---------------|-----------|---------------|---|------------|---|
| HOSPITAL NAME |           | Oper ID:      |   | 98/05/09   |   |
| ID:           | NAME:     | AGE: 0        |   |            |   |
| LMP:          | BBT:      | G0 P0 A0 E0   |   |            |   |
| Ref MD:       | NOTE:     | + 2/2         |   |            |   |
| BIOPHYSICAL:  | MOVEMENT: | TONE:         |   | BREATHING: |   |
|               | FLUID:    | REACTIVE NST: |   | SCORE: /10 |   |
| HR(BPM)       | SITE      | LAST          | 1 | 2          | 3 |
| PI            | 1         |               | ( |            | ) |
|               | 2         |               | ( |            | ) |
|               | 3         |               | ( |            | ) |
| RI            | 1         |               | ( |            | ) |
|               | 2         |               | ( |            | ) |
|               | 3         |               | ( |            | ) |
| D/S           | 1         |               | ( |            | ) |
|               | 2         |               | ( |            | ) |
|               | 3         |               | ( |            | ) |
| A/B           | 1         |               | ( |            | ) |
|               | 2         |               | ( |            | ) |
|               | 3         |               | ( |            | ) |
| COMMENTS:     |           |               |   |            |   |

Figure 9–11. Basic OB Report Page 2 for Europe, Tokyo University and Osaka University



## OB Summary Reports

---

### Biophysical Profile

A biophysical profile source is available to evaluate the well being of the fetus. The user can rate the fetus with a score of 0, 1 or 2 on the following items:

- Movement                      Rate fetal physical activity
- Tone                              Rate fetal heart tone
- Breathing                      Rate fetal breathing movement
- Fluid                              Rate fetal amniotic fluid
- Reactive NST                  Rate a Reactive Non-Stress Test

A sum of the above 5 items is then displayed as an overall biophysical profile score.

If any value is entered other than 0, 1 or 2 for each of the 5 items, an error message "Invalid, Input a proper data (0–2)" is displayed.

### HR, PI, RI, S/D (D/S) Ratio and A/B Ratio

The largest portion of page 2 displays the last 3 measurements and latest/average value for HR, PI, RI, S/D (D/S) Ratio and A/B Ratio. The last 4 measurements can be taken at 3 different sites on the fetus.

The factory default for these site fields is 1, 2 and 3. However, the user can freely edit these site fields on the Report Page. For some regions, these site fields may be common between pages 1 and 2. If this is the case, edits on either page are reflected on the other. These edits are remembered after power down. The maximum length of the site field is 6 characters.

The comments section can be edited freely using the alphanumeric keyboard or annotation library. Simply move the highlight cursor into the comment area and press the **Comment** key. The system is now in the annotation mode. When the highlight cursor leaves the comment area, the system exits the annotation mode. These comments are separate from those on page 1.



## OB Summary Reports

---

### USA Only

Two items are different on the OB Summary Report page 2 for the USA version.

1. The third line of the top portion will display all of the LMP, GA and EDD data based on the operator's input. For example, if LMP was entered:

LMP (OPE) mm/dd/yy GA (LMP) ##W#D EDD (LMP) mm/dd/yy  
will be displayed.

2. Amniotic Fluid Index (AFI) displays the sum of the 4 measurements and the 4 individual measurement values.

### Editing the Report

#### Cursor Movement

Movement of the edit cursor on the report page can be accomplished by using the Trackball or the up/down/left/right arrow keys on the keyboard.

#### Edit Actions

When positioned on a measurement field the **Back Space** or **Space** keys will delete the current value in that field. After pressing **Back Space** or **Space**, the only options remaining are to press **Set** or **Esc**.

**Set** deletes the old value.

**Esc** returns the old value to the report page.



*NOTE: Any calculations relating to the deleted measurement field are recomputed.*

#### Logic Field Edits

The CUA/CGA field requires a "Y" or "N" to be entered in either upper or lower case.

When the cursor is positioning in the desired field, press "Y" or "N". Once yes or no has been entered, press **Return** to record that entry and proceed to the next edit field.



*NOTE: Calculations affected are recomputed.*



## OB Summary Reports

---

### Comment Edits

The comment field consists of 2 lines with a total of 120 characters. Use normal keyboard functions to type necessary comments.

Comments made on each report page are not reflected on the other report pages. Therefore, comments on each page can be edited individually.

When finished, press **Return** to complete the comment edit function. Since this is the last edit field on the report page, the cursor moves to the first edit field.

### Measurement Edits

Measurements on the OB Report Page can only be edited by deletion. After a measurement is deleted, the average or other related calculations will be recomputed and the report page displayed again.

### Operator Entered Values

In the calculations section of the report page, if the value is a result of operator entered data (generally from the patient entry menu) and not measurements, the designation "OPE" is seen after the calculation name.

An example would be EDD (OPE). This means the estimate date of delivery was a value entered by the operator on the patient entry menu.

### Recording Summary Reports

The summary report can be saved like any ultrasound image. Once it is displayed on the screen, it can be recorded on the VCR, printed on the B/W or color page printer, photographed by the multi-image camera, stored on a MOD with the Image Archive option or placed on regular paper with a line printer.

### Line Printer Operation

#### USA Only

Refer to 13–21 for details on line printer operation.



# Anatomical Survey

## Overview

The Anatomical Survey page provides a checklist that promotes a thorough routine reporting of OB exams.

The header of the page is patient data similar to the corresponding OB Report page. This information can only be edited by using the **New Patient** or **ID Name** keys.

The third line of the header will be different depending on the type of OB Report Format selected.

There are 13 pre-programmed and a maximum of 6 user programmable anatomical features.

|                   |                   |                   |             |
|-------------------|-------------------|-------------------|-------------|
| HOSPITAL NAME     |                   | Oper ID: 98/05/09 |             |
| ID:               | NAME:             | AGE: 0            |             |
| LMP(OPE)          | GA(LMP)           | EDD(LMP)          | G0 P0 A0 E0 |
| Ref MD:           | NOTE:             |                   |             |
|                   | ANATOMICAL SURVEY |                   |             |
|                   | IMAGED?           |                   |             |
| HEAD              | YES               | NO                | APPEARANCE: |
| LATERAL VENT      | YES               | NO                |             |
| CEREBELLUM        | YES               | NO                |             |
| FACE              | YES               | NO                |             |
| HEART             | YES               | NO                |             |
| SPINE             | YES               | NO                |             |
| STOMACH           | YES               | NO                |             |
| KIDNEYS           | YES               | NO                |             |
| BLADDER           | YES               | NO                |             |
| CORD              | YES               | NO                |             |
| CORD INSERTION    | YES               | NO                |             |
| UPPER EXTREMITIES | YES               | NO                |             |
| LOWER EXTREMITIES | YES               | NO                |             |
|                   | YES               | NO                |             |
|                   | YES               | NO                |             |
|                   | YES               | NO                |             |
|                   | YES               | NO                |             |
|                   | YES               | NO                |             |
| COMMENTS:         |                   |                   |             |

Figure 9–12. Anatomical Survey Page



## Anatomical Survey

---

### Editing

When the Anatomical Survey Page is enabled from the OB calculation menu, the highlight cursor appears at the YES/NO field of the first feature on the checklist. Both Yes and No are displayed for all items on the list.

The highlight cursor is moved by the **Trackball** or the up/down/left/right arrow keys.

YES/NO will appear in the "IMAGED?" column even if there is no entry in the user programmed fields. Select yes by pressing 'Y' on the keyboard. When yes is selected, no will be erased. The cursor automatically moves to the comments field.

In the comments field, comments can be entered by typing or selecting annotations from the Soft-Menu. One line of comments is available for each anatomical feature. Press **Return** when completed.



*NOTE: Any meaningless entries will display the message: "Please select 'y' or 'n'." or "Please select 'n' or 'a'."*

### User Programmed Features

The user can program 6 anatomical feature fields. The maximum number of characters is 19.

To enter a user feature:

**Trackball** the highlight cursor to a User Field.

Type in the desired feature name and press **Return**. The highlight cursor moves to the Yes/No field.

Repeat the steps as necessary for any additional feature fields.

The entries will be saved and retained even after the system is powered down.



### Hints

While in a user programmed field, rolling the Trackball up/down or pressing the up/down arrow keys have no effect on cursor position. Only motion of the Trackball to the left/right or pressing the left/right arrow keys will work.



# OB Graphs

## Overview

The LOGIQ™ 500 can display a fetal growth curve graph (OB Graph) from data in each measurement table for fetal age.

The graphic displays a horizontal scale for gestational weeks and a vertical scale for the measurement value. The curve graphic shows the standard values with 2 additional curves representing the standard deviation.

If information was entered on the patient data entry menu that allows for the estimation of a gestational age, an asterisk (\*) is displayed on the graph to show the current status of the fetus relative to the measurement table.

The OB graph page also has patient data fields, measurement/calculation data fields and a comment field.

The system does not allow any field on the OB graph to be edited.



*NOTE: In many cases, the LMP (or any form of estimated gestational age) is required in order to properly compute the OB Graph. Ensure that the LMP is properly entered in the Patient Entry Menu.*



## OB Graph Selection

OB Graph is a factory default selection in the measurement menus. Select **OB Graph** from the Measurement Menu.




Figure 9–13. Basic OB Graph Display (Osaka University)



*NOTE: The label in the upper left corner of the graph indicates the standard deviation used to plot the graph.*

## OB Graph Labeling

Table 9–7 is a summary of the labeling found in the OB Graph for deviation.

## OB Graph Labeling (cont'd)

| Author              | Measurement table                    | Deviation  | Deviation Label |
|---------------------|--------------------------------------|------------|-----------------|
| Hadlock             | BPD, HC, AC, FL                      | Vertical   | +2SD/-2SD       |
|                     | CRL                                  | —          | No Label        |
|                     | EFW                                  | Vertical   | +2SD/-2SD       |
| Hellman             | GS                                   | —          | No Label        |
| ASUM                | BPD, AC                              | Horizontal | +2SD/-2SD(WEEK) |
|                     | CRL                                  | —          | No Label        |
| Tokyo               | GS, CRL, BPD, FL, LV                 | Horizontal | Deviation(WEEK) |
| Tokyo-S (Shinozuka) | CRL, BPD, FL, AC, AxT(APTDxTTD), EFW | Vertical   | +1.64SD/-1.64SD |
| Osaka               | CRL, BPD, FL, HL, FTA, EFBW          | Vertical   | +1.5SD/-1.5SD   |
| Hansmann            | BPD, CRL, OFD, HC, TThd, FL          | Vertical   | 95%/5%          |
|                     | GS, TAD                              | Vertical   | 90%/10%         |
|                     | AC                                   | —          | No Label        |
| Rempen              | GS, CRL, BPD                         | Vertical   | 95%/5%          |
| Campbell            | BPD, CRL, FL                         | Vertical   | 90%/10%         |
|                     | BD                                   | —          | No Label        |
| Berkowitz           | BD                                   | —          | No Label        |
| Eriksen             | TAD                                  | —          | No Label        |
| Kurtz               | BPD                                  | Vertical   | 90%/10%         |
| Nelson              | CRL                                  | Vertical   | 90%/10%         |
| Robinson            | CRL                                  | Vertical   | 90%/10%         |
| Jeanty              | CRL, BPD, FL, AC, HC                 | Vertical   | 90%/10%         |
|                     | BD                                   | —          | No Label        |
| Paris               | CRL, BPD, FL, Ft                     | Vertical   | 90%/10%         |
|                     | TAD                                  | —          | No Label        |
| Sostoa              | BPD, FL, AC, HC                      | Vertical   | 90%/10%         |
|                     | BD, OFD                              | —          | No Label        |
| User Table USA      |                                      | Horizontal | Deviation(WEEK) |
| User Table TOKYO    |                                      | Horizontal | Deviation(WEEK) |
| User Table OSAKA    |                                      | Vertical   | +1.5SD/-1.5SD   |
| User Table Europe   |                                      | Vertical   | 90%/10%         |

Table 9-7. OB Graph Deviation Labeling



## OB Graphs

---

### Changing OB Graph Selection

After the OB graph has been displayed, the user can change or select which measurement parameter is to be used to display a graph.

Highlight **CHANGE GRAPH** and press **Return**. The graphic is deleted and a measurement selection menu appears.

Select the desired measurement value and press **Return**. The new measurements, calculations and graphic are displayed.

### Exit OB Graph

Select **OB Graph** in the OB Calculation Sub-Menu.



#### Hints

The user is not allowed to edit any field on the OB Graph display except the Comments section. Comments on the OB Graph page are different from any other OB Report Page.

If the user does not enter an estimated gestational age on the Patient Entry Menu screen, the mark (\*) indicating the current status of the fetus will not be displayed on the graph.

# Fetal Trend Management

## (software option)

### Overview

Fetal Trend Management is an option to the LOGIQ™ 500 OB Calculation package that enhances the user's ability to monitor the development of the fetus.

If patient data, measurements and calculations are saved during the initial examination, this information can be compared to results of follow up examinations. The OB Graph function can be used to display the current data or combine the current data with past data to show a fetal growth trend.

When previous measurement data is saved, all other information, which is input at the Patient Entry Menu, is retained. This information is automatically displayed in the Patient Entry Menu on subsequent exams when the patient ID is entered in the Patient Entry Menu.

If the previous measurement data was saved on MOD, the MOD must be inserted before inputting the patient ID in the Patient Entry Menu.

### Storing Patient Information

After an OB examination, the user can save the resultant patient data to two different types of media.

The user can save the patient data to the system hard disk. While this offers the convenience of not maintaining removable media, the system hard disk's storage capacity is very limited.

A Magnetic Optical Disk (MOD) would be the best storage media. Although this removable disk must be maintained by the user, it offers far greater storage capacity. Approximately one year's worth of information can be stored on one disk.



### **Data Storage Estimations**

If the following standard assumptions are made:

- Work days/week = 5
- Patients/day = 30
- Diagnosis/patient = 5
- Measurements/patient = 5

then approximately one year's worth of studies can be stored on the hard drive.

When hard drive capacity is reached, the message "Data is full. Delete needless data" is displayed.

However, the average MOD can store more than 30 years worth of data using the same assumptions.

### **Media Selection Preset Parameter**

The user must choose the type of media that will store the patient OB information. This is done in the Set Up/Preset Program sub-menu page 3. The parameter is:

Media Selection for Fetal Trend : HD MO

Select HD for the system hard drive (limited storage space) or MOD for a removable magnetic optical disk (much greater storage space).

## Saving Data

Before the diagnosis is complete, ensure that all patient information such as Name, ID, Ref MD and EDC has been entered. If it has not, use the ID/Name key to enter this necessary information. Select the OB Graph function from the OB Calculation soft menu.

Five additional commands are displayed on the OB Graph that relate to the Fetal Trend Management option as shown in Figure 9–14. These commands are:

**LIST-ID**—Displays a list by patient ID number.

**LIST DATA**—Displays data for a specified patient.

**TREND-BOTH**—Shows data points on the OB graph based on current and past data.

**TREND-PRESENT**—Shows data points on the OB graph based on current data only.

**SAVE**—Saves Patient information to specified media.

**CHANGE GRAPH**—Changes the measurement value graphed (performs the same function as in the basic OB calculation package).




Figure 9–14. OB Graph Display



## Fetal Trend Management

---

### Type of Data Saved

The type of data that is recorded during the SAVE function is:

1. Date and Time
2. Patient Name and ID
3. Calculated EDC or EDD
4. Measurement Author's Name
5. Measured or Calculated Data

### Save Command



*NOTE: To avoid any trouble in searching later, it is better to have all patient information entered before saving.*

Use the **Trackball** to highlight the SAVE command and press **Set**.

If patient information has been entered, the patient data, measurements and calculations will be saved to the designated media (HD or MOD).

If there are similar files in the data list (i.e. same ID but Name and EDC are different), the system displays the Patient List Menu (LIST-ID). The user will have to select the file to which the data will be added.

If the user does not input the ID, the system displays the List-Data Menu. The user will then have to input the ID.

If the user does not input other patient data, the system displays the List-Data menu. The user will then have to input the patient data.

If more than one candidate is available on the list, the system first displays the Patient List Menu. The user then has to pick the proper file in which the current data will be added.



*NOTE: If measurement averaging is turned on, only the average value is saved. If measurement averaging is turned off, only the last measured or calculated value is saved.*

## Save Function Messages

If patient information has not been entered, the message displayed is:

**“Input patient’s information.”**

If a patient ID number is not entered, the message displayed is:

**“Input ID:”**



*NOTE: The ID number is necessary or the SAVE function CANNOT be accomplished. It is possible to save data without a patient name or EDC but the information displayed with the LIST-DATA function is limited.*

If all of the data is the same as a file on the data list, the message displayed is:

**“Overwrite existing data? ‘y’ ‘n’.”**

where ‘y’ replaces the archived file with the new file.  
‘n’ causes the system to do nothing. This cancels the SAVE function request.

If the storage media is full, especially in the case of the limited capacity of the hard drive, the message displayed is:

**“Data is full. Delete needless data.”**

Eliminating needless spaces in the data entered will help to conserve storage space.





## Fetal Trend Management

### Growth Trending

The default command after the OB Graph is displayed is TREND-PRESENT. This will display an OB Graph for the current author/measurement selected.

However, after additional examinations, the user can recall past data and display it with the current data to show a trend in fetal growth.

When the user selects the TREND-BOTH command from the OB Graph display, the system automatically searches the storage media and gets the data and displays the OB Trend Graph as shown in Figure 9-15.




Figure 9-15. OB Trend Graph

The system will search for like data displayed on the OB Graph. If Hadlock BPD is the data displayed, then the current patient's Hadlock BPD data is used with the past Hadlock BPD data to display the trend on the graph.

If the user wants to trend different author/measurement data, that graph must be displayed first before the TREND-BOTH command is selected.

An asterisk (\*) is used to show the present data point on the graph while the pound sign (+) is used for all past data points.

## List ID Management

The system allows the user to manage the data that was previously stored using the Patient ID List.

Choose LIST-ID from the OB Graph display. The system displays all data found on the storage media. The oldest data will have the first order number. The newest data is displayed on the last page as shown in Figure 9-16.

```

[PATIENT LIST MENU]

PT NAME      :
PT ID        :
EDD(YY/MM/DD) :  /  /

NUMBER : ██████████

No      ID      NAME      PAGE : 0 / 0
      EDD

Ctrl+N:Next Page      Ctrl+L:Data List      Ctrl+M:Modify
Ctrl+P:Prvs Page      Ctrl+D:Delete      Ctrl+K:Save
Ctrl+A:Arrange         Ctrl+S:Search      Ctrl+C:Cancel
CLEAR:Clear entered characters

```

Figure 9–16. LIST-ID Data Display

## List ID Commands

The commands at the bottom of the List-ID page allow the user to perform the following functions:

- CTRL, N—Displays the next page  
CTRL, P—Displays the previous page  
CTRL, A—Allows user to rearrange the order of the display list  
CTRL, L—Displays the selected patients data list  
CTRL, D—Deletes saved patient data  
CTRL, S—Searches the ID List for desired patient information  
CTRL, M—Allows user to modify patient data of a selected file  
CTRL, K—Saves the current patient data  
CTRL, C—Cancels the current process  
CLEAR—Clears entered characters



## Fetal Trend Management

### Control, A

Used to control the order in which the LIST-ID data is displayed.

Press **Ctrl, A** simultaneously. The cursor moves to the number window as shown in Figure 9–17.

[PATIENT LIST MENU]

PT NAME :  
PT ID :  
EDD (MM / DD / YY) : 95/03/03

NUMBER :  Arrange : 1-ID 2-Name 3-EDC +R-REV

| NO    | ID             | NAME                  | PAGE :#### | EDC |
|-------|----------------|-----------------------|------------|-----|
| 00001 | A1234567890123 | AAAAAAAAAAAAAAAAAAAAA | 95/09/15   |     |
| 00002 | B9876543210987 | BBBBBBBBBBBBBBBBBBBBB | 95/09/23   |     |
| 00003 | C6754321897653 | CCCCCCCCCCCCCCCCCCCCC | 95/09/28   |     |
|       |                |                       |            |     |
| 00015 | Z3456188273653 | ZZZZZZZZZZZZZZZZZZZZZ | 95/09/13   |     |

Ctrl+N:Next Page      Ctrl+L:Data List      Ctrl+M:Modify  
Ctrl+P:Prvs Page      Ctrl+D:Delete      Ctrl+K:Save  
Ctrl +A:Arrange      Ctrl +S:Search      Ctrl +C:Cancel

CLEAR:Clear entered characters

Operator Message Area

Figure 9–17. List Arrangement Selections

The following message is displayed:

**“Arrange : 1-ID 2-Name 3-EDC +R-REV”**

- 1 = Rearrange by ID number smallest to largest
- 1R = Rearrange by ID number largest to smallest
- 2 = Rearrange by name A to Z.
- 2R = Rearrange by name Z to A.
- 3 = Rearrange by EDD newest to oldest
- 3R = Rearrange by EDD oldest to newest.

“R” added to the number selection means reverse order.

Type in the number or number/letter for the list display arrangement desired and press **Return**.

The temporary file sort for display begins. The time required depends on the amount of information in the data base.



---

## Fetal Trend Management

---

### Control, D

The Control, D function allows the user to delete specified files from the LIST-ID data base (HD or MOD).

Press **Ctrl, D** simultaneously. The cursor moves to the number window.

Type in the number from the list displayed on the screen and press **Return**.

Before deleting any data, the following message is displayed:

**“Are you sure to delete the data ? (y/n)”**

where ‘y’ means to proceed with deleting the data.

‘n’ means do nothing.

The data selected is then deleted from the data base.

If the user is cleaning up the media, especially in the case of the Hard Drive, the following can be entered instead of a number:

|          |                                                 |
|----------|-------------------------------------------------|
| # Select | Delete selected file numbers                    |
| A-All    | Delete ALL files from the media data base.      |
| P-Page   | Delete ALL files on the current page displayed. |

### Control, K

Operates like the SAVE command on the OB Graph display. Pressing these two keys simultaneously from the PATIENT LIST MENU display causes the system to save the current patient data. See Storing Patient Information on 9-39 for details.



## Fetal Trend Management

### Control, L

The Control, L function allows the user to display the data list for a selected patient file.

Press **Ctrl, L** simultaneously. The cursor moves to the number window.

Type in the number from the list displayed on the screen and press **Return**.

The data list for the patient selected is displayed on the screen as shown in Figure 9–18.

| [DATA LIST MENU]      |       |                               |       |                       |       |       |       |                 |       |  |  |  |  |  |  |
|-----------------------|-------|-------------------------------|-------|-----------------------|-------|-------|-------|-----------------|-------|--|--|--|--|--|--|
| PT NAME               | :     | 12345678901234567890123456789 |       |                       |       |       |       |                 |       |  |  |  |  |  |  |
| PT ID                 | :     | 12345678901234                |       |                       |       |       |       |                 |       |  |  |  |  |  |  |
| EDD (MM / DD / YY)    | :     | 12 / 31 / 93                  |       |                       |       |       |       |                 |       |  |  |  |  |  |  |
| NUMBER :              |       |                               |       |                       |       |       |       |                 |       |  |  |  |  |  |  |
| List-Page :####       |       | Meas-Page :####               |       |                       |       |       |       |                 |       |  |  |  |  |  |  |
| No                    | GA    | MEAS1                         | MEAS2 | MEAS3                 | MEAS4 | MEAS5 | MEAS6 | MEAS7           | MEAS8 |  |  |  |  |  |  |
|                       | UNIT1 | UNIT2                         | UNIT3 | UNIT4                 | UNIT5 | UNIT6 | UNIT7 | UNIT8           |       |  |  |  |  |  |  |
|                       | AUTH1 | AUTH2                         | AUTH3 | AUTH4                 | AUTH5 | AUTH6 | AUTH7 | AUTH8           |       |  |  |  |  |  |  |
| 123                   | ##W#D | ###.#                         | ###.# | ###.#                 | ###.# | ###.# | ###.# | ###.#           | ###.# |  |  |  |  |  |  |
| 123                   | ##W#D | ###.#                         | ###.# | ###.#                 | ###.# | ###.# | ###.# | ###.#           | ###.# |  |  |  |  |  |  |
| 123                   | ##W#D | ###.#                         | ###.# | ###.#                 | ###.# | ###.# | ###.# | ###.#           | ###.# |  |  |  |  |  |  |
|                       |       |                               |       |                       |       |       |       |                 |       |  |  |  |  |  |  |
| 123                   | ##W#D | ###.#                         | ###.# | ###.#                 | ###.# | ###.# | ###.# | ###.#           | ###.# |  |  |  |  |  |  |
| 123                   | ##W#D | ###.#                         | ###.# | ###.#                 | ###.# | ###.# | ###.# | ###.#           | ###.# |  |  |  |  |  |  |
|                       |       |                               |       |                       |       |       |       |                 |       |  |  |  |  |  |  |
| Ctrl+N: Next Page     |       |                               |       | Ctrl+G: Graph         |       |       |       | Ctrl+K: Save    |       |  |  |  |  |  |  |
| Ctrl+P: Prvs Page     |       |                               |       | Ctrl+D: Delete        |       |       |       | Ctrl +C: Cancel |       |  |  |  |  |  |  |
| Ctrl+F: Forward Page  |       |                               |       | Ctrl+B: Backward Page |       |       |       |                 |       |  |  |  |  |  |  |
| Operator Message Area |       |                               |       |                       |       |       |       |                 |       |  |  |  |  |  |  |

Figure 9–18. Patient List-Data Menu

Page two of the Data List Menu shows the Number of studies, Gestational Age (by LMP), Tech ID and Reference MD information.



## Fetal Trend Management

---

### Control, M

The Control, M function allows the user to modify the patient data (Name, ID & EDD) on any file in the data base. This will be helpful if the patient information was not entered EXACTLY the same for each exam.

Press **Ctrl, M** simultaneously. The cursor moves to the number window as shown in Figure 9–19.

| [PATIENT LIST MENU] |                      |
|---------------------|----------------------|
| PT NAME             | :                    |
| PT ID               | :                    |
| EDD(Y Y/M M/D D)    | : / /                |
| NUMBER :            | <input type="text"/> |

Figure 9–19. Modified Data Entry

Type in the number from the list displayed on the screen and press **Return**.

The selected patient information is displayed and the cursor moves to the PT NAME window.

Use the **Trackball** to highlight the desired window and change the necessary information with the alphanumeric keys. Once the **Return** key is pressed while in the EDD window or the **Trackball** is used to move the cursor outside of the patient information area, the following message appears:

**“ Are you sure to modify the previous data? (y/n) ”**

where ‘y’ means modify the old data and save it. The cursor moves to the PT Name window.

‘n’ means do nothing. The cursor moves to the PT Name window.

If the desired modifications will create a file exactly like a current existing file, the following message appears:

**“Same data exists.”**



## Fetal Trend Management

---

### Control, M (cont'd)

At this point the user may:

#C – Combine a selected number data file with the one being modified. The file being modified will be added to the file number selected.

Q – This will quit the sub mode operation and return to the previous modification mode.



*NOTE: Caution should be used when modifying files so as not to change/delete existing files or create false ones.*

### Control, S

The Control, S function allows the user to quickly search the PATIENT LIST MENU data base for a specific patient file.

Press **Ctrl, S** simultaneously. The cursor moves to the PT Name window as shown in Figure 9–20.

| [PATIENT LIST MENU] |       |
|---------------------|-------|
| PT NAME             | :     |
| PT ID               | :     |
| EDD (YY/MM/DD)      | : / / |
| NUMBER              | :     |

Figure 9–20. Search Data Entry

Use the **Trackball** to highlight the desired window and type the necessary information with the alphanumeric keys. Once the **Return** key is pressed while in the EDD window or the **Trackball** is used to move the cursor outside of the patient information area, the search functions begin.

When the search has ended, the list is displayed on the screen.

If the search ends in failure, the system may beep or the cursor remains at the PT Name window and the following message is displayed:

**“No data exists. Input other information.”**

At the end of the search, the patient information is placed in the windows (Name, ID & EDD) of the file that most closely matches the search criteria.

## Data List Management

The LIST-DATA function allows the user to manage the measurement/calculation data stored in a specific patient file.

Choose LIST-DATA from the OB Graph or the Control, L function from the LIST-ID display. The system displays past gestational age data for the same patient criteria displayed on the OB Graph or for the information that was entered in the List-ID display. Refer to Figure 9–21.

The measurement/calculation list order is based on gestational age and cannot be rearranged.

| [PATIENT DATA MENU]            |       |                               |       |       |       |                |       |       |       |  |  |  |
|--------------------------------|-------|-------------------------------|-------|-------|-------|----------------|-------|-------|-------|--|--|--|
| PT NAME                        | :     | 12345678901234567890123456789 |       |       |       |                |       |       |       |  |  |  |
| PT ID                          | :     | 12345678901234                |       |       |       |                |       |       |       |  |  |  |
| EDD (MM / DD / YY)             | :     | 12 / 31 / 93                  |       |       |       |                |       |       |       |  |  |  |
| NUMBER :                       |       | <input type="text"/>          |       |       |       |                |       |       |       |  |  |  |
| List-Page : #####              |       | Meas-Page : #####             |       |       |       |                |       |       |       |  |  |  |
| No                             | GA    | 1                             | 2     | 3     | 4     | 5              | 6     | 7     | 8     |  |  |  |
|                                |       | MEAS1                         | MEAS2 | MEAS3 | MEAS4 | MEAS5          | MEAS6 | MEAS7 | MEAS8 |  |  |  |
|                                |       | UNIT1                         | UNIT2 | UNIT3 | UNIT4 | UNIT5          | UNIT6 | UNIT7 | UNIT8 |  |  |  |
|                                |       | AUTH1                         | AUTH2 | AUTH3 | AUTH4 | AUTH5          | AUTH6 | AUTH7 | AUTH8 |  |  |  |
| 123                            | ##W#D | ###.#                         | ###.# | ###.# | ###.# | ###.#          | ###.# | ###.# | ###.# |  |  |  |
| 123                            | ##W#D | ###.#                         | ###.# | ###.# | ###.# | ###.#          | ###.# | ###.# | ###.# |  |  |  |
| 123                            | ##W#D | ###.#                         | ###.# | ###.# | ###.# | ###.#          | ###.# | ###.# | ###.# |  |  |  |
|                                |       |                               |       |       |       |                |       |       |       |  |  |  |
| 123                            | ##W#D | ###.#                         | ###.# | ###.# | ###.# | ###.#          | ###.# | ###.# | ###.# |  |  |  |
| 123                            | ##W#D | ###.#                         | ###.# | ###.# | ###.# | ###.#          | ###.# | ###.# | ###.# |  |  |  |
|                                |       |                               |       |       |       |                |       |       |       |  |  |  |
| Ctrl+N:Next Page               |       | Ctrl+L:Data List              |       |       |       | Ctrl+M:Modify  |       |       |       |  |  |  |
| Ctrl+P:Prvs Page               |       | Ctrl+D:Delete                 |       |       |       | Ctrl+K:Save    |       |       |       |  |  |  |
| Ctrl +A:Arrange                |       | Ctrl +S:Search                |       |       |       | Ctrl +C:Cancel |       |       |       |  |  |  |
| CLEAR:Clear entered characters |       |                               |       |       |       |                |       |       |       |  |  |  |
| Operator Message Area          |       |                               |       |       |       |                |       |       |       |  |  |  |

Figure 9–21. List-Data Menu

## Data List Commands

Some functions at the bottom of the data list screen operate the same as those on the LIST-ID screen:

CTRL, N—Displays the next page  
 CTRL, P—Displays the previous page  
 CTRL, D—Deletes saved patient data.  
 CTRL, K—Saves the current patient data  
 CTRL, C—Cancels the current process  
 CLEAR—Clears entered characters

Functions specific to the LIST-ID display are:

CTRL, F—Moves to the right for the displayed page  
 CTRL, B—Moves to the left for the displayed page  
 CTRL, G—Displays the OB Graph





## *Fetal Trend Management*

---

### **Control, F and Control, B**

Since only eight measurements/calculations can be displayed on the screen at one time, Control F and Control B allow the user to move to the right and left, respectively, to display additional measurements/calculations.

### **Control, G**

Control, G enables the OB Graph function for the number of files selected.

Press **Ctrl, G** simultaneously. The cursor moves to the number window.

Type in the number from the list displayed on the screen and press **Return**.

The selected OB Graph will be displayed and the TREND–BOTH selection is highlighted.

# OB–Multigestational

## (software option)

### Overview

The LOGIQ™ 500 offers an optional calculation package that allows the user to measure and report multiple fetus development. The system is capable of reporting a maximum of four fetuses.

### Patient Entry Menu

If the Multigestational Option has been purchased, an extra entry appears on the Patient Entry Menu to the right of “RefMD:”. This entry is called “FETUS NUMBER:”. The factory default number is one.

### Entering Fetus Number

It is generally during the first exam that the user recognizes a multiple gestation. If during the first exam more than one fetus is imaged, the user can use the ID/Name key to return to the Patient Entry Menu and enter the necessary number of fetuses.

For subsequent exams the user enters the correct number of fetuses when the New Patient key is used.



*NOTE: Data Management Center transfer function will not work if the fetus number is greater than one. An error message will appear:*

*“Multigestation data transfer is not supported”*

*DMC data transfer will function normally if the fetus number is set to one.*



## Distinguishing Each Fetus

For measurements/calculations and report page displays, the fetuses are labeled A, B, C and D.

Each fetus is distinguished on the report pages by its letter and total number of fetuses. For example, a report page could be noted as FETUS:A/3. This is fetus A from a total of 3.

Following this fetal number designation, the user can use twelve characters to describe the position and twelve characters to describe placenta location of the fetus. This is displayed as:

POS: ##### PLAC: #####

## Measurements/Calculations

During measurements/calculations, all measurements are performed on fetus A first. The user then switches to fetus B and performs the same measurements, etc.

The user can switch between fetuses by pressing **Blue Shift** and then **N**.



*NOTE: A fetus change can be accomplished at anytime during the exam. However, the Blue Shift, N function is NOT available while a measurement cursor is active.*

After changing to the next fetus, any measurements made will be recorded and reported to the new designated fetus. Any active measurement or calculation not completed will be cancelled when the fetus number is changed.



## **Change the Number of Fetuses**

In the ID/Name function the user can increase or decrease the fetus number or cancel the Multigestational function. This would need to be accomplished due to the demise of a fetus or error in the definition of the fetus number.

### **Number Increase**

If the fetus number is increased, the message:

“Are you sure to change the number? (y/n)”

is displayed.

‘y’ will create the report pages for the new fetus.

‘n’ returns to the previous number.

If a fetus number increase is not available, the message: “Data range over. Input proper data.” is displayed.

### **Number Decrease**

If the fetus number is decreased because some measurement has already been completed or the fetus number was fixed in a previous exam, the message

“Are you sure to change the number? (y/n)”

is displayed.

‘y’ prompts the user to select which fetus to delete (A, B, C or D).

‘n’ returns to the previous number.

The data of deleted fetuses is deleted for the measurement result window and report pages. If the fetus number is reduced to one, the multigestation function is cancelled.

After the decrease in fetus number, there is no need to adjust the IDs (A, B, C or D) of the remaining fetuses. The same ID for the remaining fetuses can be maintained and used.



## Report Page Layout

The Multigestational Option adds additional pages to the OB Report displays.

Pages one and two of the Multigestational Option appear the same as the Basic OB Reporting package; except, pages one and two will be generated for each fetus. They will have the fetus identifying information between the heading and measurements. The fetal information displayed is:

“Fetus : \$/% POS: ##### PLAC: #####”

\$ = Fetal ID (A, B, C or D)

% = Total number of fetuses

The second to last page, depending on the number of fetuses, is a summary of the measurements and calculations found on page one. However, all fetuses (maximum of 4) are displayed on this page. Since page one is different by region, page three is also different by region selected (i.e. USA, European, Tokyo University or Osaka University).

The last page, depending on the number of fetuses, is a summary that compares all fetuses. The heading is different by region. It contains the following:

|                 |                                                                              |
|-----------------|------------------------------------------------------------------------------|
| Fetus Pos:      | Fetus position. This 12 character string is common between pages 1, 2 and 4. |
| Fluid:          | An estimate of the amount of amniotic fluid. Common between pages 2 and 4.   |
| Biophys:        | Total biophysical profile score from page 2.                                 |
| HR (BPM):       | Heart rate in beats per minute from page 2.                                  |
| AFI (USA Only): | Amniotic Fluid Index summary from page 2.                                    |
| Placenta:       | Shows location fetus, grade and previa (yes, no or partial)                  |
| Membrane:       | Shows location between fetus to fetus.                                       |

## OB Graph

With the Multigestational Option, the OB Graph selection includes pages to display graphs for individual fetuses as well as a graph to plot all fetuses simultaneously.

For individual fetus reports, the fetus is designated by fetus ID and total number of fetuses (i.e. Fetus: A/3 or Fetus: B/2)

Change the fetus pages by moving the highlight cursor to ← or → and press **Set**.




Figure 9–22. OB Graph Multigestational Option (Osaka University Version)



## OB–Multigestational

### OB Graph (cont'd)

If all fetuses are displayed simultaneously, different symbols are used to mark each fetus. The symbols are:

| <u>Fetus</u> | <u>Present</u> |
|--------------|----------------|
| A            | *              |
| B            | ⊗              |
| C            | ⊠              |
| D            | ⋄              |




Figure 9–23 OB Graph Multigestational Option (Osaka University Version)

Only gestational age, based on operator's input, is displayed for simultaneous plotting.

## Fetal Trend Management (Multigestational Option)

Fetal Trend Management with the Multigestational Option operates much like the basic fetal trend package described earlier. Figure 9–24 shows an example of the Multigestational Fetal Trend Graph.




Figure 9–24. Fetal Trend Graph Multigestational Option (Hadlock-USA Version)

## Save Function

The **SAVE** command on the trend graph display or **Ctrl + K** in the Data List display will cause the system to compare the current fetus combination between present and past data. If patient data and fetus number are compatible, the data is saved.

If there is any mismatch in data, the message:

“Incompatible data. Do you save? (y/n)”

is displayed.

‘y’ = Saves data as it is

‘n’ = Cancel save function





## Data List Menu

The Data List Menu of fetal trend information is slightly different with the Multigestational Option. Each fetus has it's own display page as shown in Figure 9–25.

[DATA LIST MENU]

PT NAME :  
PT ID :  
EDD(MM/DD/YY) : / /

NUMBER :

List-Page:  
No    GA

Meas-Page:

Ctrl+N: Next Page    Ctrl+G: Graph    Ctrl+K: Save  
Ctrl+P: Prvs Page    Ctrl+D: Delete    Ctrl+C: Cancel  
Ctrl+F: Forward Page    Ctrl+B: Backward Page

Figure 9–25. Multigestational Data List Menu

# Data Management Center (DMC)

## Overview

The LOGIQ™ 500 is capable of interfacing with a personal computer (PC). This interface is used to send OB, GYN and Vascular measurement and calculation data from the LOGIQ™ 500 to the PC for display, processing and evaluation.

Additionally, the LOGIQ™ 500 can accept patient data from a personal computer (PC) and transfer it to the patient entry menu.



*NOTE: Data Management Center transfer function will not work if the fetus number is greater than one with the Multigestation Option. An error message appears:  
"Multigestation data transfer is not supported"*

*DMC data transfer will function normally if the fetus number is set to one.*

## Operational Setup

The PC must be connected to Port A or B of the LOGIQ™ 500 by an RS-232C interface cable. "Computer" must be selected for Port A or B on the Set Up/System Parameters page 5.

The **Record 1** or **Record 2** keys must be assigned to the computer selection. The Record 1 and Record 2 presets are found on Set Up/Systems Parameters page 5.

## Operational Setup (cont'd)




Figure 9–26. System Parameters Page 5

## Transferring OB Data

Data transfer from the LOGIQ™ 500 to the DMC computer (PC) will only occur while operating in the OB, GYN or Vascular Exam categories.

The user should input all necessary information into the Patient Entry Menu. Scan the patient. Measure and calculate all the desired parameters.

Press **Record 1** or **Record 2**, depending on the key assignment. Communication with the PC will start by displaying the message: "Communicating with computer now".



**NOTE:** Normal system operation is suspended during the transfer of data.

When all data has been transferred with no errors, the message "Communication was completed" is displayed.

Refer to the DMC User Manual for processing exam data on the PC.

## Error Messages

If an anomaly occurs in the communication or transfer of data between the LOGIQ™ 500 and the PC, an error message will be displayed on the LOGIQ™ 500 screen. Possible error messages are:

| Message                                         | Possible Cause                                                                                       |
|-------------------------------------------------|------------------------------------------------------------------------------------------------------|
| Check Computer. No power.                       | PC power is off. Interface cable broken or disconnected.                                             |
| Check Computer. Bad condition.                  | PC is busy. LOGIQ™ 500 senses the PC cannot accept data.                                             |
| Check Computer. Abnormal response.              | LOGIQ™ 500 senses an abnormal response from the PC.                                                  |
| Check Computer. Checksum Error occurred.        | LOGIQ™ 500 senses an abnormal acknowledgement from the PC after all data has been sent.              |
| Check Computer. No response.                    | LOGIQ™ 500 senses no acknowledgement from the computer after all data has been sent.                 |
| Can't send data to Computer in this category.   | User is trying to send data to the PC in an exam category other than OB, GYN or Vascular.            |
| System doesn't have option.                     | The OB calculation option for the LOGIQ™ 500 has been turned off.                                    |
| SYSTEM ERROR: CANNOT COMMUNICATE WITH COMPUTER. | A LOGIQ™ 500 system error has occurred that will not allow communication with the computer to begin. |
| This function is not available.                 | This function is not available in some situations.                                                   |

Table 9–8. Data Management Center Error Messages



## Error Messages (cont'd)



### Hints

- Data transfer can only be accomplished in the OB, GYN and Vascular Exam categories.
- Date is always transferred as Year–Month–Day. Time is transferred as 24 hour time (i.e. 00:00 to 23:59).
- Gestational Age Data is always transferred in weeks and days.
- The LOGIQ™ 500 can store three values of the same measurement. Only the average of selected measurements will be transferred.
- The system keyboard is locked out during the transfer of data.
- Only the Record 1 or Record 2 keys can start the data transfer.
- The following DMC items are not sent by the LOGIQ™ 500:
  - Serial No.
  - DateRep
  - TimeRep
  - RepPtName
  - RefMDID
  - History
  - FetusPos
  - FetusSex.



*NOTE: Currently, the transfer of multi-fetus information is not compatible with the DMC. If it is desired to use the DMC with multiple fetuses, each fetus should be treated individually by using a unique Patient ID number for each.*



---

## *Data Management Center (DMC)*

---

### **Patient Data Input**

The PC (personal computer) can be used to read patient data cards and transfer that data to the Patient Entry menu via the DMC serial port connection.

### **Preparation**

Connect the PC and LOGIQ™ 500 using an isolated RS232C cable and serial ports equivalent to the DMC cabling.

Assign “Computer” to the proper LOGIQ™ 500 serial port in the Setup/System Parameter Menu page 5.

Display the New Patient Entry Menu on the LOGIQ™ 500.

### **Send Data**

On the PC, read the hospital ID card with patient data and transfer that data to the LOGIQ™ 500. The transfer starts by displaying the following message on the LOGIQ™ 500 monitor:

“Communicating with Computer now.”

When the reception is complete with no errors, the the following message on the LOGIQ™ 500 monitor:

“Communication was Completed.”

The patient data received is then displayed on the Patient Entry Menu.



## Data Management Center (DMC)

---

### Data Transferred

The data items transferred and accepted by the LOGIQ™ 500 are:

- Patient Name
- Patient ID Number
- Patient Age
- Patient Sex
- Patient Birthday
- Patient Height
- Patient Weight
- Patient Body Surface Area

### Error Messages

Refer to Table 9–8 on 9–63 for error messages in patient data transfer.

If the New Patient Entry Menu is not displayed when the PC tries to send data, the error message:

“Can’t receive data. Return to New Patient Menu.”

is displayed.



### Hints

If data is omitted, only the items sent by the PC will be displayed in the Patient Entry Menu.

If both age and birthday are transferred to the LOGIQ™ 500, the system will use the birth date to calculate the age.

All LOGIQ™ 500 keys cannot be used while the data transfer is in progress.



# GYN Measurements

## B-Mode

### Ovarian Length, Height, and Width

The length, height and width of the left and right ovaries can be measured and recorded on the GYN summary report.

Each measurement is a typical distance measurement made on a B-Mode image in the appropriate scan plane.

Typically, the height and width are measured on the axial plane while the length is measured on the sagittal plane.

To measure ovarian length:

Scan the patient's right or left ovary in the sagittal plane.

|         |         |         |            |
|---------|---------|---------|------------|
|         |         |         |            |
| Lt Ov-L | Lt Ov-H | Lt Ov-W |            |
| Rt Ov-L | Rt Ov-H | Rt Ov-W | GYN Report |

Figure 9–27. GYN Calculation Sub-Menu

Select **Rt Ov-L** or **Lt Ov-L** from the GYN Calculation Sub-Menu. Perform a standard distance measurement as described in 7–7.

To measure ovarian height and width:

Scan the patient's right or left ovary in the axial plane.

Select **Rt Ov-H**, **Rt Ov-W**, **Lt Ov-H** or **Lt Ov-W** from the GYN Calculation Sub-Menu. Perform the necessary distance measurements as described in 7–7.





## GYN Measurements

---

### Uterine Length, Height, and Width

The GYN Calculation Sub-Menu page 2 allows for measurements to be recorded as uterine length (Ut-L), height (Ut-H) or width (Ut-W).

Length is typically measured in the sagittal plane while height and width are measured in the axial plane.

Scan the patient in the appropriate scan plane.

|      |      |      |            |
|------|------|------|------------|
|      |      |      |            |
| Ut-L | Ut-H | Ut-W |            |
| Endo |      |      | GYN Report |

Figure 9–28. GYN Calculation Sub-Menu

Select the desired measurement from the GYN Calculation Sub-Menu (Ut-L, Ut-H or Ut-W).

Measurements are recorded on the GYN Summary Report. The last 3 measurements of each type can be averaged by the Summary Report Page.

### Endometrium Thickness

An endometrium thickness measurement is also available on the GYN Calculation Sub-Menu.

Scan the patient.

Select **ENDO** from the GYN Calculation Sub-Menu. Perform a standard distance measurement as described in 7–7.

### Doppler Mode

Three GYN related measurements can be taken in Doppler Mode and recorded on the GYN Summary Report. These are resistive index measurements for a left ovarian vessel, uterine vessel and right ovarian vessel.



Resistive index

RI measurements are velocity values measured on the Doppler Spectrum for vessels in the ovaries and uterus.


To make GYN Resistive Index measurements:

Scan the desired vessel in Doppler mode.

Select the desired **RI** measurement from the GYN Calculation Sub-Menu.

|              |           |             |               |
|--------------|-----------|-------------|---------------|
|              |           |             |               |
| Lt Ov-<br>RI | Ut-<br>RI | R Ov-<br>RI | TRACE<br>AUTO |
|              |           |             | GYN<br>Report |

Figure 9–29. GYN Calculation Sub-Menu



Use the **Trackball** to move the cursor to the peak point of the Doppler waveform where it is desired to start the RI measurement.

Press **Set** to fix the start-point cursor. An end-point cursor appears.

Use the **Trackball** to move the end-point cursor to the diastolic point of the Doppler waveform.

Press **Set** to complete the measurement. A calculated value is displayed and stored in the GYN Summary Report.



# GYN Summary Report

## GYN Report Layout

|                                                                          |       |          |    |          |    |
|--------------------------------------------------------------------------|-------|----------|----|----------|----|
| HOSPITAL NAME                                                            |       | Oper ID: |    | 05/09/98 |    |
| ID:                                                                      | NAME: |          |    | AGE:     |    |
| LMP:                                                                     |       | G0       | P0 | A0       | E0 |
| Ref MD:                                                                  | NOTE: |          |    |          |    |
| MEASUREMENTS                                                             | AVE   | 1        | 2  | 3        |    |
| Uterus Length                                                            | (     |          |    |          | )  |
| Uterus Width                                                             | (     |          |    |          | )  |
| Uterus Height                                                            | (     |          |    |          | )  |
| Right Ovary Length                                                       | (     |          |    |          | )  |
| Right Ovary Width                                                        | (     |          |    |          | )  |
| Right Ovary Height                                                       | (     |          |    |          | )  |
| Left Ovary Length                                                        | (     |          |    |          | )  |
| Left Ovary Width                                                         | (     |          |    |          | )  |
| Left Ovary Height                                                        | (     |          |    |          | )  |
| Endometrium Thickness                                                    | (     |          |    |          | )  |
| CALCULATIONS                                                             | AVE   | 1        | 2  | 3        |    |
| Uterus RI                                                                |       |          |    |          |    |
| Right Ovary RI                                                           |       |          |    |          |    |
| Left Ovary RI                                                            |       |          |    |          |    |
| COMMENTS:                                                                |       |          |    |          |    |
| <div style="background-color: black; height: 1.2em; width: 100%;"></div> |       |          |    |          |    |

Figure 9–30. GYN Report Page

The top portion of the report contains patient information entered into the Patient Entry Menu.

The middle portion of the report (measurements) shows the last 3 measurement values made. <AVG> or <LAST> indicates that the Average or Latest value will be displayed in this column. This depends on the preset, *Average Activity*, in the Setup/Preset Program menu page 3.

The calculation section shows Doppler measurements and their average, while the comments area provides 2 lines of comment data.

## GYN Calculation Formulas

| Calc Mnemonic | Calc Name                            | Input Measurements                     | Formula                                                        |
|---------------|--------------------------------------|----------------------------------------|----------------------------------------------------------------|
| UT-L          | Uterine Length                       | one distance                           | Ut-L[cm or mm]=d1                                              |
| UT-H          | Uterine Height                       | one distance                           | Ut-H[cm or mm]=d1                                              |
| UT-W          | Uterine Width                        | one distance                           | Ut-W[cm or mm]=d1                                              |
| Endo          | Endometrium Thick-ness               | one distance                           | Endo[cm or mm]=d1                                              |
| Lt. Ov-L      | Left Ovarian Length                  | one distance                           | Lt. Ov-L[cm or mm]=d1                                          |
| Lt. Ov-H      | Left Ovarian Height                  | one distance                           | Lt. Ov-H[cm or mm]=d1                                          |
| Lt. Ov-W      | Left Ovarian Width                   | one distance                           | Lt. Ov-W[cm or mm]=d1                                          |
| Rt. Ov-L      | Right Ovarian Length                 | one distance                           | Rt. Ov-L[cm or mm]=d1                                          |
| Rt. Ov-H      | Right Ovarian Height                 | one distance                           | Rt. Ov-H[cm or mm]=d1                                          |
| Rt. Ov-W      | Right Ovarian Width                  | one distance                           | Rt. Ov-W[cm or mm]=d1                                          |
| Lt. Ov-RI     | Left Ovarian Vessel Resistive Index  | two Doppler blood flow peak velocities | Lt. Ov-RI= $(V_{\max}-V_{\text{diastole}})/V_{\max}$<br>*NOTE* |
| Ut-RI         | Uterine Vessel Resistive Index       | two Doppler blood flow peak velocities | Ut-RI= $(V_{\max}-V_{\text{diastole}})/V_{\max}$<br>*NOTE*     |
| Rt. Ov-RI     | Right Ovarian Vessel Resistive Index | two Doppler blood flow peak velocities | Rt. Ov-RI= $(V_{\max}-V_{\text{diastole}})/V_{\max}$<br>*NOTE* |

\*NOTE\*:  $V_{\text{diastole}} = V_{\min}$  or  $V_{\text{end-diastole}}$  (depends on preset selection)

Table 9–9. GYN Calculation Formulas



## *GYN Summary Report*

---

This page left blank intentionally.

# Introduction

## Overview

The basic cardiac measurement package offers a choice of LV Measurement methods. Analysis of the left ventricle can be performed by 6 different methods.

- Cubed method
- Teichholz method
- Bullet method
- Modified Simpson's Rule method
- Single-plane method
- Biplane ellipsoid method

Each selection requires a series of related measurements specific to that method.

Selection of the LV measurement method is made from the initial Cardiac Calculation Sub-Menu.

|             |               |              |              |
|-------------|---------------|--------------|--------------|
|             |               |              |              |
| LV<br>CUBED | LV<br>BULLET  | LV<br>SP-ELP |              |
| LV<br>TEICH | LV<br>SIMPSON | LV<br>BP-ELP | LV<br>Report |

Figure 10–1. Cardiac Calculation Sub-Menu



## Introduction

---

### Overview (cont'd)

Once the LV measurement method is chosen, a Sub-Menu with the available measurement selections appears.

In general, measurement techniques are the same as those described in the *General Measurements and Calculations* chapter. After selecting a measurement, be sure to follow the prompts on the image screen.

Measurements are organized in sequence to streamline system usage and give easy access to calculated values on the report page. Some dependent calculations require several measurements to be taken before the result can be calculated.

## Report Pages

<AVG> or <LAST> indicates that the Average or Latest value will be displayed in this column. This depends on the preset, *Average Activity*, in the Set Up/Preset Program menu page 3.

The 3 report page measurements can be averaged and the average used in other calculations.

After displaying an LV Summary Report, the calculation method can be changed by moving the highlighted cursor to “Change Method” and pressing **Set**.

The current report will be cleared and a method selection menu appears. Use the **Trackball** to select a new LV method and press **Set**.

After pressing **Set** the new LV method report is displayed.

|                                                  |         |                                                                                                                                                                              |      |          |
|--------------------------------------------------|---------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------|----------|
| HOSPITAL NAME                                    |         | Oper ID:                                                                                                                                                                     |      | 05/09/98 |
| ID:                                              | NAME:   |                                                                                                                                                                              |      | AGE:     |
| SEX:                                             | HEIGHT: | WEIGHT:                                                                                                                                                                      | BSA: |          |
| Ref MD:                                          | NOTE:   |                                                                                                                                                                              |      |          |
| METHOD: Biplane Ellipsoid                        |         | SELECT                                                                                                                                                                       |      |          |
| MEASUREMENTS                                     |         | <div style="border: 1px solid black; padding: 5px;"> Cubed<br/> Teichholz<br/> Bullet<br/> Modified Simpson's Rule<br/> Single Plane Ellipsoid<br/> Biplane Ellipsoid </div> |      |          |
| LVMLd<br>LVMLs<br>LVAd<br>LVAs<br>LVAMd<br>LVAMs |         |                                                                                                                                                                              |      |          |
| CALCULATIONS                                     |         |                                                                                                                                                                              |      |          |
| HR<br>EdV<br>EsV<br>SV<br>CO<br>EF               |         |                                                                                                                                                                              |      |          |
| CHANGE METHOD                                    |         |                                                                                                                                                                              |      |          |
| COMMENTS:                                        |         |                                                                                                                                                                              |      |          |

Figure 10–2. Report Page Change Method





## Introduction

---

### Recording Summary Reports

The summary report can be saved like any ultrasound image. Once it is displayed on the screen, it can be recorded on the VCR, printed on the B/W or color page printer, photographed by the multi-image camera, stored on a MOD with the Image Archive option or placed on regular paper with a line printer.

#### USA Only

Line Printer Operation—Refer to 13–21 for details on line printer operation.

### BSA Calculation Methods

The LOGIQ™ 500 can calculate BSA (Body Surface Area) by 2 different formulas, depending on the method selected in the Set Up/System Parameter Menu page 2.

#### Oriental Formula

$$\text{BSA [m}^2\text{]} = (\text{Height [cm]})^{.03} \times (\text{Weight [g]})^A \times 3.207 \times 10^{-4}$$
$$A = 0.7285 - 0.0188 \times \log (\text{Weight [g]})$$

#### Occidental Formula

$$\text{BSA [m}^2\text{]} = (0.425 \times \log (\text{Weight [g]}) + 0.725 \times \log (\text{Height [cm]}) + 1.8654) \times 10^{-4}$$



NOTE: 1 inch = 2.54 cm  
1 lb = 453.59243 g  
g = grams  
cm = centimeters

#### Heart Rate (HR)

Refer to the procedure found on 10–30.

# LV Measurement Methods

## Cubed Method

|       |      |       |            |
|-------|------|-------|------------|
|       |      |       |            |
| LVIDd | IVSd | LVPWd | HEART RATE |
| LVIDs | IVSs | LVPWs | Return     |

Figure 10–3. Cardiac Calculation Sub-Menu—Cubed Method

|                      |         |          |      |          |   |
|----------------------|---------|----------|------|----------|---|
| HOSPITAL NAME        |         | Oper ID: |      | 05/09/98 |   |
| ID:                  | NAME:   |          |      | AGE:     |   |
| SEX:                 | HEIGHT: | WEIGHT:  | BSA: |          |   |
| Ref MD:              | NOTE:   |          |      |          |   |
| METHOD: Cubed        |         |          |      |          |   |
| MEASUREMENTS         |         | AVE      | 1    | 2        | 3 |
| LVIDd                |         |          | {    |          | } |
| LVIDs                |         |          | {    |          | } |
| IVSd                 |         |          | {    |          | } |
| IVSs                 |         |          | {    |          | } |
| LVPWd                |         |          | {    |          | } |
| LVPWs                |         |          | {    |          | } |
| CALCULATIONS         |         |          |      |          |   |
| HR                   |         |          |      |          |   |
| Edv                  |         |          |      |          |   |
| EsV                  |         |          |      |          |   |
| SV                   |         |          |      |          |   |
| CO                   |         |          |      |          |   |
| EF                   |         |          |      |          |   |
| FS                   |         |          |      |          |   |
| <b>CHANGE METHOD</b> |         |          |      |          |   |
| COMMENTS:            |         |          |      |          |   |

Figure 10–4. Cubed Report Page Display

## LV Measurement Methods

### Cubed Method (cont'd)




Figure 10–5. Cubed Method Measurements

| Calc Mnemonic | Calc Name                                         | Input Measurements                         | Formula                                                                 |
|---------------|---------------------------------------------------|--------------------------------------------|-------------------------------------------------------------------------|
| LVIDd         | Left Ventricular Internal Diameter, Diastole      | one distance                               | $LVIDd = d1[\text{cm or mm}]$                                           |
| LVIDs         | Left Ventricular Internal Diameter, Systole       | one distance                               | $LVIDs = d1[\text{cm or mm}]$                                           |
| IVSd          | Interventricular Septal Thickness, Diastole       | one distance                               | $IVSd = d1[\text{cm or mm}]$                                            |
| IVSs          | Interventricular Septal Thickness, Systole        | one distance                               | $IVSs = d1[\text{cm or mm}]$                                            |
| LVPWd         | Left Ventricle Posterior Wall Thickness, Diastole | one distance                               | $LVPWd = d1[\text{cm or mm}]$                                           |
| LVPWs         | Left Ventricle Posterior Wall Thickness, Systole  | one distance                               | $LVPWs = d1[\text{cm or mm}]$                                           |
| EdV           | End Diastole Volume                               | one distance                               | $EdV[\text{ml}] = LVIDd^3$                                              |
| EsV           | End Systole Volume                                | one distance                               | $EsV[\text{ml}] = LVIDs^3$                                              |
| SV            | Stroke Volume                                     | two distance                               | $SV[\text{ml}] = EdV - EsV$                                             |
| CO            | Cardiac Output                                    | two distances and one 2 beat time interval | $CO[1/\text{min}] = SV \times HR / 1000$                                |
| EF            | Ejection Fraction                                 | two distances                              | $EF = SV / EdV \times 100$                                              |
| FS            | Fractional Shortening                             | two distances                              | $FS = (1 - LVIDs / LVIDd) \times 100$                                   |
| HR            | Heart Rate (beats/minute)                         | one 2 beat time interval                   | $HR[\text{BPM}] = 120 [\text{sec}] / 2 \text{ beat time } [\text{sec}]$ |
| PHT           | Pressure Half Time                                | one time interval                          | $PHT = t1[\text{ms or sec}]$                                            |
| MVA           | Mitral Valve Area                                 | one pressure half time interval            | $MVA[\text{cm}^2] = 220 / PHT$                                          |
| ET            | Ejection Time                                     | one time interval                          | $ET = t1[\text{ms or sec}]$                                             |

Table 10–1. Cubed Method Formulas

## Teichholz Method

|       |      |       |            |
|-------|------|-------|------------|
|       |      |       |            |
| LVIDd | IVSd | LVPWd | HEART RATE |
| LVIDs | IVSs | LVPWs | Return     |

Figure 10–6. Cardiac Calculation Sub-Menu—Teichholz Method

|                      |         |          |      |          |
|----------------------|---------|----------|------|----------|
| HOSPITAL NAME        |         | Oper ID: |      | 05/09/98 |
| ID:                  | NAME:   |          |      | AGE:     |
| SEX:                 | HEIGHT: | WEIGHT:  | BSA: |          |
| Ref MD:              | NOTE:   |          |      |          |
| METHOD: Teichholz    |         |          |      |          |
| MEASUREMENTS         | AVE     | 1        | 2    | 3        |
| LVIDd                |         | {        |      | }        |
| LVIDs                |         | {        |      | }        |
| IVSd                 |         | {        |      | }        |
| IVSs                 |         | {        |      | }        |
| LVPWd                |         | {        |      | }        |
| LVPWs                |         | {        |      | }        |
| CALCULATIONS         |         |          |      |          |
| HR                   |         |          |      |          |
| Edv                  |         |          |      |          |
| EsV                  |         |          |      |          |
| SV                   |         |          |      |          |
| CO                   |         |          |      |          |
| EF                   |         |          |      |          |
| FS                   |         |          |      |          |
| <b>CHANGE METHOD</b> |         |          |      |          |
| COMMENTS:            |         |          |      |          |

Figure 10–7. Teichholz Report Page Display




Figure 10–8. Teichholz Methods Measurements



## LV Measurement Methods

---

### Teichholz Method (cont'd)

| Calc Mnemonic | Calc Name                                         | Input Measurements                         | Formula                                                  |
|---------------|---------------------------------------------------|--------------------------------------------|----------------------------------------------------------|
| LVIDd         | Left Ventricular Internal Diameter, Diastole      | one distance                               | $LVIDd = d1[\text{cm or mm}]$                            |
| LVIDs         | Left Ventricular Internal Diameter, Systole       | one distance                               | $LVIDs = d1[\text{cm or mm}]$                            |
| IVSd          | Interventricular Septal Thickness, Diastole       | one distance                               | $IVSd = d1[\text{cm or mm}]$                             |
| IVSs          | Interventricular Septal Thickness, Systole        | one distance                               | $IVSs = d1[\text{cm or mm}]$                             |
| LVPWd         | Left Ventricle Posterior Wall Thickness, Diastole | one distance                               | $LVPWd = d1[\text{cm or mm}]$                            |
| LVPWs         | Left Ventricle Posterior Wall Thickness, Systole  | one distance                               | $LVPWs = d1[\text{cm or mm}]$                            |
| EdV           | End Diastole Volume                               | one distance                               | $EdV[\text{ml}] = LVIDd^3 \times 7 / (2.4 + LVIDd)$      |
| EsV           | End Systole Volume                                | one distance                               | $EsV[\text{ml}] = LVIDs^3 \times 7 / (2.4 + LVIDs)$      |
| SV            | Stroke Volume                                     | two distances                              | $SV[\text{ml}] = EdV - EsV$                              |
| CO            | Cardiac Output                                    | two distances and one 2 beat time interval | $CO[1/\text{min}] = SV \times HR / 1000$                 |
| EF            | Ejection Fraction                                 | two distances                              | $EF = SV / EdV \times 100$                               |
| FS            | Fractional Shortening                             | two distances                              | $FS = (1 - LVIDs / LVIDd) \times 100$                    |
| HR            | Heart Rate (beats/minute)                         | one 2 beat time interval                   | $HR = 120[\text{sec}] / 2 \text{ beat time}[\text{sec}]$ |
| PHT           | Pressure Half Time                                | one time interval                          | $PHT = t1[\text{ms or sec}]$                             |
| MVA           | Mitral Valve Area                                 | one pressure half time interval            | $MVA[\text{cm}^2] = 220 / PHT$                           |
| ET            | Ejection Time                                     | one time interval                          | $ET = t1[\text{ms or sec}]$                              |

Table 10–2. Teichholz Method Formulas

## Bullet Method

|      |       |  |               |
|------|-------|--|---------------|
|      |       |  |               |
| LVLd | LVAMd |  | HEART<br>RATE |
| LVLs | LVAMs |  | Return        |

Figure 10–9. Cardiac Calculation Sub-Menu—Bullet Method

|                      |         |          |   |          |   |
|----------------------|---------|----------|---|----------|---|
| HOSPITAL NAME        |         | Oper ID: |   | 05/09/98 |   |
| ID:                  | NAME:   | BSA:     |   | AGE:     |   |
| SEX:                 | HEIGHT: | WEIGHT:  |   |          |   |
| Ref MD:              | NOTE:   |          |   |          |   |
| METHOD: Bullet       |         |          |   |          |   |
| MEASUREMENTS         |         | AVE      | 1 | 2        | 3 |
| LVLd                 |         |          | ( |          | ) |
| LVLs                 |         |          | ( |          | ) |
| LVAMd                |         |          | ( |          | ) |
| LVAMs                |         |          | ( |          | ) |
| CALCULATIONS         |         |          |   |          |   |
| HR                   |         |          |   |          |   |
| EdV                  |         |          |   |          |   |
| EsV                  |         |          |   |          |   |
| SV                   |         |          |   |          |   |
| CO                   |         |          |   |          |   |
| EF                   |         |          |   |          |   |
| <b>CHANGE METHOD</b> |         |          |   |          |   |
| COMMENTS:            |         |          |   |          |   |

Figure 10–10. Bullet Report Page Display




Figure 10–11. Bullet Method Measurements



## LV Measurement Methods

---

### Bullet Method (cont'd)

| Calc Mnemonic | Calc Name                                     | Input Measurements                                                                     | Formula                                                                 |
|---------------|-----------------------------------------------|----------------------------------------------------------------------------------------|-------------------------------------------------------------------------|
| LVLd          | Left Ventricular Length, Diastole             | one distance                                                                           | $LVLd = d1[\text{cm or mm}]$                                            |
| LVLs          | Left Ventricular Length, Systole              | one distance                                                                           | $LVLs = d1[\text{cm or mm}]$                                            |
| LVAMd         | Left Ventricular Area, Mitral Valve, Diastole | one area (by ellipse, trace or circle)                                                 | $LVAMd = a1[\text{cm}^2]$                                               |
| LVAMs         | Left Ventricular Area, Mitral Valve, Systole  | one area (by ellipse, trace or circle)                                                 | $LVAMs = a1[\text{cm}^2]$                                               |
| EdV           | End Diastole Volume                           | one distance and one area (by ellipse, trace or circle)                                | $EdV[\text{ml}] = 5 \times LVLd \times LVAMd / 6$                       |
| EsV           | End Systole Volume                            | one distance and one area (by ellipse, trace or circle)                                | $EsV[\text{ml}] = 5 \times LVLs \times LVAMs / 6$                       |
| SV            | Stroke Volume                                 | two distances and two areas (by ellipse, trace or circle)                              | $SV[\text{ml}] = EdV - EsV$                                             |
| CO            | Cardiac Output                                | two distances and two areas (by ellipse, trace or circle) and one 2 beat time interval | $CO[1/\text{min}] = SV \times HR / 1000$                                |
| EF            | Ejection Fraction                             | two distances and two areas (by ellipse, trace or circle)                              | $EF = SV / EdV \times 100$                                              |
| HR            | Heart Rate (beats/minute)                     | one 2 beat time interval                                                               | $HR[\text{BPM}] = 120 [\text{sec}] / 2 \text{ beat time } [\text{sec}]$ |
| PHT           | Pressure Half Time                            | one time interval                                                                      | $PHT = t1[\text{ms or sec}]$                                            |
| MVA           | Mitral Valve Area                             | one pressure half time interval                                                        | $MVA[\text{cm}^2] = 220 / PHT$                                          |
| ET            | Ejection Time                                 | one time interval                                                                      | $ET = t1[\text{ms or sec}]$                                             |

Table 10–3. Bullet Method Formulas

## Modified Simpson's Rule Method

|      |       |       |            |
|------|-------|-------|------------|
|      |       |       |            |
| LVLd | LVAMd | LVAPd | HEART RATE |
| LVLs | LVAMs | LVAPs | Return     |

Figure 10–12. Cardiac Calculation Sub-Menu—Modified Simpson's Rule Method

|                                 |         |          |      |          |   |
|---------------------------------|---------|----------|------|----------|---|
| HOSPITAL NAME                   |         | Oper ID: |      | 05/09/98 |   |
| ID:                             | NAME:   |          |      | AGE:     |   |
| SEX:                            | HEIGHT: | WEIGHT:  | BSA: |          |   |
| Ref MD:                         | NOTE:   |          |      |          |   |
| METHOD: Modified Simpson's Rule |         |          |      |          |   |
| MEASUREMENTS                    |         | AVE      | 1    | 2        | 3 |
| LVLd                            |         |          | (    |          | ) |
| LVLs                            |         |          | (    |          | ) |
| LVAMd                           |         |          | (    |          | ) |
| LVAMs                           |         |          | (    |          | ) |
| LVAPd                           |         |          | (    |          | ) |
| LVAPs                           |         |          | (    |          | ) |
| CALCULATIONS                    |         |          |      |          |   |
| HR                              |         |          |      |          |   |
| EdV                             |         |          |      |          |   |
| EsV                             |         |          |      |          |   |
| SV                              |         |          |      |          |   |
| CO                              |         |          |      |          |   |
| EF                              |         |          |      |          |   |
| <b>CHANGE METHOD</b>            |         |          |      |          |   |
| COMMENTS:                       |         |          |      |          |   |

Figure 10–13. Modified Simpson's Rule Report Page Display




Figure 10–14. Modified Simpson's Rule Method Measurements





## LV Measurement Methods

### Modified Simpson's Rule Method (cont'd)

| Calc Mnemonic | Calc Name                                          | Input Measurements                                                                      | Formula                                                                |
|---------------|----------------------------------------------------|-----------------------------------------------------------------------------------------|------------------------------------------------------------------------|
| LVLd          | Left Ventricular Length, Diastole                  | one distance                                                                            | $LVLd = d1[\text{cm or mm}]$                                           |
| LVLs          | Left Ventricular Length, Systole                   | one distance                                                                            | $LVLs = d1[\text{cm or mm}]$                                           |
| LVAMd         | Left Ventricular Area, Mitral Valve, Diastole      | one area (by ellipse, trace or circle)                                                  | $LVAMd = a1[\text{cm}^2]$                                              |
| LVAMs         | Left Ventricular Area, Mitral Valve, Systole       | one area (by ellipse, trace or circle)                                                  | $LVAMs = a1[\text{cm}^2]$                                              |
| LVAPd         | Left Ventricular Area, Papillary Muscles, Diastole | one area (by ellipse, trace or circle)                                                  | $LVAPd = a1[\text{cm}^2]$                                              |
| LVAPs         | Left Ventricular Area, Papillary Muscles, Systole  | one area (by ellipse, trace or circle)                                                  | $LVAPs = a1[\text{cm}^2]$                                              |
| EdV           | End Diastole Volume                                | one distance and two areas (by ellipse, trace or circle)                                | $EdV[\text{ml}] = LVLd/3 \times (LVAMd + (LVAMd + LVAPd)/2 + LVAPd/3)$ |
| EsV           | End Systole Volume                                 | one distance and two areas (by ellipse, trace or circle)                                | $EsV[\text{ml}] = LVLs/3 \times (LVAMs + (LVAMs + LVAPs)/2 + LVAPs/3)$ |
| SV            | Stroke Volume                                      | two distances and four areas (by ellipse, trace or circle)                              | $SV[\text{ml}] = EdV - EsV$                                            |
| CO            | Cardiac Output                                     | two distances and four areas (by ellipse, trace or circle) and one 2 beat time interval | $CO[1/\text{min}] = SV \times HR/1000$                                 |
| EF            | Ejection Fraction                                  | two distances and four areas (by ellipse or trace or circle)                            | $EF = SV/EdV \times 100$                                               |
| HR            | Heart Rate (beats/minute)                          | one 2 beat time interval                                                                | $HR = 120[\text{sec}] / 2 \text{ beat time}[\text{sec}]$               |
| PHT           | Pressure Half Time                                 | one time interval                                                                       | $PHT = t1[\text{ms or sec}]$                                           |
| MVA           | Mitral Valve Area                                  | one pressure half time interval                                                         | $MVA[\text{cm}^2] = 220/PHT$                                           |
| ET            | Ejection Time                                      | one time interval                                                                       | $ET = t1[\text{ms or sec}]$                                            |

Table 10–4. Modified Simpson's Rule Method Formulas



Single Plane Ellipsoid Method

|      |      |  |               |
|------|------|--|---------------|
|      |      |  |               |
| LVLd | LVAd |  | HEART<br>RATE |
| LVLs | LVAs |  | Return        |

Figure 10–15. Cardiac Calculation Sub-Menu—Single Plane Ellipsoid Method

|                                |         |          |   |          |   |
|--------------------------------|---------|----------|---|----------|---|
| HOSPITAL NAME                  |         | Oper ID: |   | 05/09/98 |   |
| ID:                            | NAME:   | BSA:     |   | AGE:     |   |
| SEX:                           | HEIGHT: | WEIGHT:  |   |          |   |
| Ref MD:                        | NOTE:   |          |   |          |   |
| METHOD: Single Plane Ellipsoid |         |          |   |          |   |
| MEASUREMENTS                   |         | AVE      | 1 | 2        | 3 |
| LVLd                           |         |          | { |          | } |
| LVLs                           |         |          | { |          | } |
| LVAd                           |         |          | { |          | } |
| LVAs                           |         |          | { |          | } |
| CALCULATIONS                   |         |          |   |          |   |
| HR                             |         |          |   |          |   |
| EdV                            |         |          |   |          |   |
| EsV                            |         |          |   |          |   |
| SV                             |         |          |   |          |   |
| CO                             |         |          |   |          |   |
| EF                             |         |          |   |          |   |
| <b>CHANGE METHOD</b>           |         |          |   |          |   |
| COMMENTS:                      |         |          |   |          |   |

Figure 10–16. Single Plane Ellipsoid Report Page Display




Figure 10–17. Single Plane Ellipsoid Method Measurements



## LV Measurement Methods

---

### Single Plane Ellipsoid Method (cont'd)

| Calc Mnemonic | Calc Name                         | Input Measurements                                                                      | Formula                                                                |
|---------------|-----------------------------------|-----------------------------------------------------------------------------------------|------------------------------------------------------------------------|
| LVLd          | Left Ventricular Length, Diastole | one distance                                                                            | $LVLd = d1[\text{cm or mm}]$                                           |
| LVLs          | Left Ventricular Length, Systole  | one distance                                                                            | $LVLs = d1[\text{cm or mm}]$                                           |
| LVAd          | Left Ventricular Area, Diastole   | one area (by ellipse, trace or circle)                                                  | $LVAd = a1[\text{cm}^2]$                                               |
| LVAs          | Left Ventricular Area, Systole    | one area (by ellipse, trace or circle)                                                  | $LVAs = a1[\text{cm}^2]$                                               |
| EdV           | End Diastole Volume               | one distance and two areas (by ellipse, trace or circle)                                | $EdV[\text{ml}] = 8 / (3 \pi) \times (LVAd)^2 / LVLd$                  |
| EsV           | End Systole Volume                | one distance and two areas (by ellipse, trace or circle)                                | $EsV[\text{ml}] = 8 / (3 \pi) \times (LVAs)^2 / LVLs$                  |
| SV            | Stroke Volume                     | two distances and four areas (by ellipse, trace or circle)                              | $SV[\text{ml}] = EdV - EsV$                                            |
| CO            | Cardiac Output                    | two distances and four areas (by ellipse, trace or circle) and one 2 beat time interval | $CO[1/\text{min}] = SV \times HR / 1000$                               |
| EF            | Ejection Fraction                 | two distances and four areas (by ellipse, trace or circle)                              | $EF = SV / EdV \times 100$                                             |
| HR            | Heart Rate (beats/minute)         | one 2 beat time interval                                                                | $HR[\text{BPM}] = 120 [\text{sec}] / 2 \text{ beat time} [\text{sec}]$ |
| PHT           | Pressure Half Time                | one time interval                                                                       | $PHT = t1[\text{ms or sec}]$                                           |
| MVA           | Mitral Valve Area                 | one pressure half time interval                                                         | $MVA[\text{cm}^2] = 220 / PHT$                                         |
| ET            | Ejection Time                     | one time interval                                                                       | $ET = t1[\text{ms or sec}]$                                            |

Table 10–5. Single Plane Ellipsoid Method Formulas



Bi Plane Ellipsoid Methods

|       |       |      |            |
|-------|-------|------|------------|
|       |       |      |            |
| LVMLd | LVAMd | LVAd | HEART RATE |
| LVMLs | LVAMs | LVAs | Return     |

Figure 10–18. Cardiac Calculation Sub-Menu—Bi Plane Ellipsoid Method

HOSPITAL NAME  
ID:  
SEX: HEIGHT: NAME: WEIGHT: Oper ID: 05/14/98  
Ref MD: NOTE: BSA: AGE:  
METHOD: Biplane Ellipsoid

MEASUREMENTS  
LVMLd  
LVMLs  
LVAd  
LVAs  
LVAMd  
LVAMs

CALCULATIONS  
HR  
EdV  
EsV  
SV  
CO  
EF

CHANGE METHOD

COMMENTS:

| AVE | 1 | 2 | 3 |
|-----|---|---|---|
| (   |   |   | ) |
| (   |   |   | ) |
| (   |   |   | ) |

Figure 10–19. Bi Plane Ellipsoid Report Page Display




Figure 10–20. Bi Plane Ellipsoid Methods Measurements



## LV Measurement Methods

### Bi Plane Ellipsoid Method (cont'd)

| Calc Mnemonic | Calc Name                                         | Input Measurements                                                                      | Formula                                                                |
|---------------|---------------------------------------------------|-----------------------------------------------------------------------------------------|------------------------------------------------------------------------|
| LVMLd         | Left Ventricle Medial-Lateral Dimension, Diastole | one distance                                                                            | $LVMLd = d1[\text{cm or mm}]$                                          |
| LVMLs         | Left Ventricle Medial-Lateral Dimension, Systole  | one distance                                                                            | $LVMLs = d1[\text{cm or mm}]$                                          |
| LVAMd         | Left Ventricular Area, Mitral Valve, Diastole     | one area (by ellipse, trace or circle)                                                  | $LVAMd = a1[\text{cm}^2]$                                              |
| LVAMs         | Left Ventricular Area, Mitral Valve, Systole      | one area (by ellipse, trace or circle)                                                  | $LVAMs = a1[\text{cm}^2]$                                              |
| LVAd          | Left Ventricle Area, Diastole                     | one area (by ellipse, trace or circle)                                                  | $LVAd = a1[\text{cm}^2]$                                               |
| LVAs          | Left Ventricle Area, Systole                      | one area (by ellipse, trace or circle)                                                  | $LVAs = a1[\text{cm}^2]$                                               |
| EdV           | End Diastole Volume                               | one distance and two areas (by ellipse, trace or circle)                                | $EdV[\text{ml}] = 8/(3\pi) \times (LVAd \times LVAMd) / LVMLd$         |
| EsV           | End Systole Volume                                | one distance and two areas (by ellipse, trace or circle)                                | $EsV[\text{ml}] = 8/(3\pi) \times (LVAs \times LVAMs) / LVMLs$         |
| SV            | Stroke Volume                                     | two distances and four areas (by ellipse, trace or circle)                              | $SV[\text{ml}] = EdV - EsV$                                            |
| CO            | Cardiac Output                                    | two distances and four areas (by ellipse, trace or circle) and one 2 beat time interval | $CO[1/\text{min}] = SV \times HR / 1000$                               |
| EF            | Ejection Fraction                                 | two distances and four areas (by ellipse, trace or circle)                              | $EF = SV / EdV \times 100$                                             |
| HR            | Heart Rate (beats/minute)                         | one 2 beat time interval                                                                | $HR[\text{BPM}] = 120 [\text{sec}] / 2 \text{ beat time} [\text{sec}]$ |
| PHT           | Pressure Half Time                                | one time interval                                                                       | $PHT = t1[\text{ms or sec}]$                                           |
| MVA           | Mitral Valve Area                                 | one pressure half time interval                                                         | $MVA[\text{cm}^2] = 220 / PHT$                                         |
| ET            | Ejection Time                                     | one time interval                                                                       | $ET = t1[\text{ms or sec}]$                                            |

Table 10–6. Bi Plane Ellipsoid Method Formulas

# Additional Cardiology Calculations

## Volume

The volume calculation can be made from 1, 2 or 3 distances, from 1 distance and an ellipse, from a single ellipse or from 2 ellipse measurements.

*General Measurements and Calculations* details how to make distance and ellipse measurements.



***IMPORTANT NOTE:*** When a volume calculation is desired, the necessary measurements or combination of measurements must be made **BEFORE** selecting **VOLUME** in the Sub-menu.

|        |            |  |              |
|--------|------------|--|--------------|
|        |            |  |              |
| VOLUME | %<br>Steno |  |              |
| ANGLE  |            |  | LV<br>Report |

Figure 10–21. Cardiac Calculation Sub-Menu



## Additional Cardiac Calculations

---

### Examples

When a volume calculation is desired, do one of the following:

- Make 1 distance measurement.
- Make 2 distance measurements.
- Make 1 ellipse measurement.
- Make 3 distance measurements.



*NOTE: Three distances should be done in the dual format mode (side by side images). One measurement is usually made in the sagittal plane and 2 measurements in the axial plane.*

- Make 1 distance and 1 ellipse measurement.
- Make 2 ellipse measurements.

Formulas used are found in the table on 10–21.



### Hints

- Volumes are most accurate when measurements are taken in the sagittal and axial scan planes.
- Use the side by side dual format option to display sagittal and axial plane images simultaneously.
- Make all necessary measurements before selecting volume.
- See *General Measurements and Calculations* for details on the steps to take when doing basic distance and area measurements.
- An overflow error message (OOR) will occur during a 2 ellipse, or 1 ellipse/1 distance volume calculation if: the minor (small) axis of the large ellipse is equal to/less than the major (large) axis of the small ellipse (d1), or equal to/less than the single distance measurement (d2).




Figure 10–22. Volume Calculation Error Possibilities



## Additional Cardiac Calculations

### Volume Calculation Formulas

| Calc Name                   | Input Measurements                                                                                                                           | Formula                                                                                                                    |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------|
| Volume (spherical)          | one distance                                                                                                                                 | $\text{Vol}[\text{ml}] = (\pi/6) \times d^3$                                                                               |
| Volume (prolate spheroidal) | two distances, $d_1 > d_2$                                                                                                                   | $\text{Vol}[\text{ml}] = (\pi/6) \times d_1 \times d_2^2$                                                                  |
| Volume (prolate spheroidal) | one ellipse, $d_1$ major axis, $d_2$ minor axis                                                                                              | $\text{Vol}[\text{ml}] = (\pi/6) \times d_1 \times d_2^2$                                                                  |
| Volume (spheroidal)         | three distances                                                                                                                              | $\text{Vol}[\text{ml}] = (\pi/6) \times d_1 \times d_2 \times d_3$                                                         |
| Volume (spheroidal)         | one distance $d_1$ , one ellipse, $d_2$ major axis, $d_3$ minor axis                                                                         | $\text{Vol}[\text{ml}] = (\pi/6) \times d_1 \times d_2 \times d_3$                                                         |
| Volume (spheroidal)         | two ellipse, ellipse 1 with axes $d_1$ and $d_2$ , ellipse 2 with axes $d_3$ and $d_4$ , with $ d_2 - d_3  \leq  d_1 - d_4 $ and $d_2 > d_3$ | $\text{Vol}[\text{ml}] = (\pi/6) \times d_1 \times d_2 \times d_4$<br>( $d_3$ is not used, assuming it is close to $d_2$ ) |

Table 10–7. Volume Calculation Formulas




Figure 10–23. Volume Calculation Examples



## Additional Cardiac Calculations

### Angle


This function is intended to measure the angle between 2 intersecting planes:

|        |         |  |           |
|--------|---------|--|-----------|
|        |         |  |           |
| VOLUME | % Steno |  |           |
| ANGLE  |         |  | LV Report |

Figure 10–24. Cardiac Calculations Sub-Menu

Select **ANGLE** from the Calculations Sub-Menu. A “×” cursor with a vertical dashed line appears. The displayed angle is zero degrees.

Use the **Trackball** to position the line cursor and the **Zoom Size/Rotation** control to adjust the angle of the line cursor.




Press **Set** to fix the position of the first cursor.

Use the **Trackball** to position the second line cursor and the **Zoom Size/Rotation** control to adjust the angle of the second line cursor. The angle between the 2 cursors is constantly updated on the display.



## Additional Cardiac Calculations

### Angle (cont'd)



If a fine adjustment of the 2 cursors is needed, press **Measurement** to toggle which line cursor is active. Adjust the active cursor with the **Trackball** and **Zoom Size/Rotation** control.

Press **Set** to complete the angle measurement.

### % stenosis (stenosis ratio)

To calculate percent stenosis in a B-Mode image:

The area of the blood vessel can be measured by the ellipse, trace, circle, 2 distance or 1 distance methods. See *General Measurements and Calculations* for details on the steps to take when doing basic distance and area measurements. The Set Up/Preset Program page 3 defines the default method.

Select **% Steno** from the Calculations Sub-Menu.

|        |         |  |           |
|--------|---------|--|-----------|
|        |         |  |           |
| VOLUME | % Steno |  |           |
| ANGLE  |         |  | LV Report |

Figure 10–25. Cardiac Calculations Sub-Menu

The first cursor is the residual area of the blood vessel and the second cursor is the lumen area of the blood vessel.

### CAUTION



When using a distance measurement to calculate area for % stenosis, the measurement should always be taken from a cross sectional view of the vessel.

Do NOT take a distance measurement from a longitudinal view for an area calculation. This may lead to an inaccurate assessment of % stenosis.



## Additional Cardiac Calculations

---

### PHT (Pressure Half Time)

A PHT measurement is used to calculate Mitral Valve Area (MVA). The measurement is taken on the Doppler spectral display of the mitral valve.




Figure 10–26. PHT Measurement on Mitral Valve Waveform

Select **PHT** from the Calculation Sub-Menu. Perform a time interval measurement as described 7–18. Set the first cursor on the peak of the mitral valve waveform and the second on to a point halfway down the mitral valve waveform slope.

|     |    |         |               |
|-----|----|---------|---------------|
|     |    |         |               |
| PHT | ET | Max PG  | TRACE<br>AUTO |
| MVA |    | Mean PG | LV<br>Report  |

Figure 10–27. Cardiac Calculation Sub-Menu

### MVA (Mitral Valve Area)

A Mitral Valve Area is calculated from the pressure half time measurement. If PHT was selected, MVA is also calculated. If MVA is selected, PHT must be measured in order to calculate MVA.

## Additional Cardiac Calculations

### ET (Ejection Time)

Ejection time for the left ventricle is the time during which the aortic valve is open. It is measured on a Doppler Spectrum display. The time between the two cursors is ejection time (ET).




Figure 10–28. Ejection Time for Aortic Valve

Select **ET** from the Calculation Sub-Menu. Perform a time interval measurement as described on 7–18. Set the first cursor at the opening of the aortic valve and the second at the closing of the aortic valve.

|     |    |         |               |
|-----|----|---------|---------------|
|     |    |         |               |
| PHT | ET | Max PG  | TRACE<br>AUTO |
| MVA |    | Mean PG | LV<br>Report  |

Figure 10–29. Cardiac Calculation Sub-Menu



## Additional Cardiac Calculations

---

### Max PG

To measure MAX PG (Maximum Pressure Gradient):

|     |    |         |               |
|-----|----|---------|---------------|
|     |    |         |               |
| PHT | ET | Max PG  | TRACE<br>AUTO |
| MVA |    | Mean PG | LV<br>Report  |

Figure 10–30. Cardiac Calculations Sub-Menu

Select **MAX PG** from the Calculations Sub-Menu. Perform peak velocity measurements as described on 7–14. Set the first cursor on  $V_{MAX}$  and the second on Vd.

### Automatic Calculation of Max PG

If velocity was previously measured and calculated in Doppler Mode, MAX PG is calculated from the  $V_{MAX}$  and Vd measurements taken for velocity when MAX PG is first selected in the Calculation Sub-Menu.



---

## Additional Cardiac Calculations

---

### Mean PG

To measure Mean PG (Mean Pressure Gradient):

|     |    |         |               |
|-----|----|---------|---------------|
|     |    |         |               |
| PHT | ET | Max PG  | TRACE<br>AUTO |
| MVA |    | Mean PG | LV<br>Report  |

Figure 10–31. Cardiac Calculation Sub-Menu

Select **MEAN PG** from the Calculation Sub-Menu. Perform the measurement described on 7–15. Set the first cursor on the waveform trace starting point.

If Trace Auto is selected, the waveform is automatically traced after the second cursor is fixed.

If Trace Auto is not selected, manually trace the desired portion of the waveform.

Press **Set** to complete the the Mean PG measurement.

### Automatic Calculation of Mean PG

If velocity was previously measured and calculated in Doppler Mode, Mean PG is calculated from the velocity value when Mean PG is first selected in the Calculation Sub-Menu.

The waveform is automatically traced. Velocity and Mean PG are also displayed on the screen.



## Additional Cardiac Calculations

---

### Trace Auto

To turn Doppler Auto Trace on or off:

|     |    |         |               |
|-----|----|---------|---------------|
|     |    |         |               |
| PHT | ET | Max PG  | TRACE<br>AUTO |
| MVA |    | Mean PG | LV<br>Report  |

Figure 10–32. Cardiac Calculation Sub-Menu

Select **TRACE AUTO** from the Calculation Sub-Menu.

Use the Trace Auto Sub-Menu selection to turn on/off the Doppler Auto Trace function.

If Trace Auto is highlighted, the Doppler envelope will automatically be traced as determined by parameters set in Custom Display page 13 and Preset Program page 7.

If Trace Auto is not highlighted, the Auto Trace function is off and traces can be performed manually.

### Trace Method

The system will trace the peak value, floor and peak, mean and peak or mode and peak, depending on the method selected by the user.



## Additional Cardiac Calculations

### S/D (D/S) Ratio, RI, A/B Ratio or PI

To measure the systolic/diastolic (diastolic/systolic) ratio, resistance index, velocity ratio and pulsatility index:

|           |           |            |              |
|-----------|-----------|------------|--------------|
|           |           |            |              |
| S/D Ratio | A/B Ratio | HEART RATE | TRACE AUTO   |
| PI        | RI        |            | Transf CALCS |

Transf CALCS is only available with Realtime Doppler Calculation option installed

Figure 10–33. Cardiac Calculation Sub-Menu

Select **S/D (D/S) RATIO**, **RI**, **A/B RATIO** or **PI** from the Calculations Sub-Menu. Perform peak velocity measurements as described on 7–14.

The first cursor is the start point on the Doppler waveform. This would be systole for S/D ratio (diastole for D/S ratio), peak velocity for RI, “A” velocity for A/B ratio and  $V_{MAX}$  for PI.

The second cursor is the end-point cursor to the end point of the Doppler waveform. This would be diastole for S/D ratio (systole for D/S ratio), minimum velocity for RI, “B” velocity for A/B ratio and Vd for PI.



*NOTE: For the PI calculation, if Trace Auto is not selected, manually trace the waveform between  $V_{MAX}$  and Vd.*

*If 2 velocity measurements are made prior to selecting S/D (D/S) ratio, A/B ratio or RI, the system automatically calculates S/D (D/S) ratio, A/B ratio or RI from the last two velocity measurements recorded.*



*NOTE: For the PI calculation, if Trace Auto is on, the system automatically traces the waveform when **Set** is pressed to fix Vd.*



### Hints

With S/D (D/S), PI and RI, if “Vd” is selected as the diastole velocity in Preset Program page 4, the end point is defined as Vd if Trace Auto is OFF.





## Additional Cardiac Calculations

---

### Heart Rate (HR)

To calculate the Heart Rate:

|              |              |               |                 |
|--------------|--------------|---------------|-----------------|
|              |              |               |                 |
| S/D<br>Ratio | A/B<br>Ratio | HEART<br>RATE | TRACE<br>AUTO   |
| PI           | RI           |               | Transf<br>CALCs |

Transf CALCS is only available with Realtime Doppler Calculation option installed

Figure 10–34. Cardiac Calculations Sub-Menu

Select **HEART RATE** from the Calculations Sub-Menu.  
Measure the time between two heartbeats.




Figure 10–35. Two Heart Beat Reference

### Transf Calcs

This selection is only available if the Realtime Doppler Calculation option is installed and Realtime Doppler Calculations are displayed in the upper left corner of the screen.

Use the Trans Calcs selection to turn on/off the transfer of Automatic Doppler Calculations displayed on the screen to the report page. See 7–19 for details on Realtime Doppler operation.

# ECG Option

## Overview

A physiological input panel is available for the LOGIQ™ 500. This panel has inputs for ECG, physiological and auxiliary signals.

The ECG Soft-Menu controls the signals connected to the Physiological panel.

Approved accessory cables provide the proper signals to the Physiological Panel.




Figure 10–36. Optional Physiological Input Panel



**NOTE:** The LOGIQ™ 500 can calculate heart rate in BPM from the ECG waveform. See Set Up/Custom Display page 17 to turn ECG Heart Rate Display on or off.

## WARNING



Do not use with defibrillator.



## *ECG Option*

---

### **Physio Sweep Speed**

The sweep speed of the physio signal on the B-Mode image can be set independent of the timeline (M-Mode and Doppler) sweep speed.

This is accomplished using the preset parameter “Physio Sweep Speed on B”. It is found on Set Up/Custom Display page 3.

### **ECG Sub-Menu**

The 3 pages of the ECG Sub-Menus provide for control of the physiological input signals. In each case, select the proper sub-menu page by pressing the **Sub-Menu Select** rocker switch.

## ECG Lead Placement

The 3 patient ECG leads are color coded white, black and green.

White is connected to the patient's right arm, black to the left arm, and green (ground) to the right foot (often placed on the right side of the abdomen). The right arm connection may change if the patient is in the decubitus position. Refer to Figure 10–37.




Figure 10–37. Common ECG Lead Placement

Once the leads are connected to the patient and the LOGIQ™ 500, the ECG waveform amplifier needs about 10 seconds to stabilize the waveform on the CRT before adjustments are made.

## ECG Sync Mark Display

A flashing heart can be displayed as a R-wave sync mark when a valid ECG trace is present. It can be turned on or off in Set Up/Custom Display page 17. This sync mark is added after system memory. Therefore, it will not be able to be recalled during cine review.



## ECG Sub-Menu Page 1

| B      | Preset           | Set Up      | ECG         |
|--------|------------------|-------------|-------------|
| Single | Sync.<br>Selectn | R<br>Delay  | PCG<br>Wave |
| Dual   | Ref.<br>Scan     | ECG<br>Wave | Aux<br>Wave |

Figure 10–38. ECG Sub-Menu Page 1

|                                  | Single                                                                                                                                                                                                                                    | Dual                                                                                                                                                                                                                                                                                                                                                  |
|----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | Updates the B-Mode image at the R1 time. Each time the R1 point is reached on the ECG, PCG, or AUX Waveform, the B-Mode image changes. In between R1 trigger points, the B-Mode image is frozen. The display shows the delay value as S1. | Updates the B-Mode image at the R1 and R2 times. Each time the R1 and R2 points are reached on the ECG, PCG, or Aux waveform, the B-Mode image changes. In between each trigger time, the B-Mode image is frozen.<br><br><i>NOTE: R<sub>2</sub> on the ECG will be marked with a vertical dotted line. The display shows the delays as S1 and S2.</i> |
| <b>Benefits</b>                  | Allows for the acquisition of an image at a specific point in the Cardiac Cycle.                                                                                                                                                          | Allows for the acquisition of an image at two points in the Cardiac Cycle.                                                                                                                                                                                                                                                                            |
| <b>Values</b>                    | If Single is highlighted, the one-trigger function is on. If Single is not highlighted, the function is off.<br><br><i>NOTE: The Trackball can be used to change the value of the active delay.</i>                                       | If Dual is highlighted, the two-trigger function is on. If Dual is not highlighted, the two-trigger function is off.<br><br><i>NOTE: The Trackball can be used to change the value of the active delay.</i>                                                                                                                                           |
| <b>Affects on Other Controls</b> | If Dual was highlighted when single is pressed, the Dual function is turned off and the Single function turned on.                                                                                                                        | If Single is highlighted when Dual is pressed, the single function is turned off and the Dual function turned on.                                                                                                                                                                                                                                     |

## ECG Sub-Menu Page 1

|                                  | Sync Selectn                                                                                                                                                                                                                                    | Ref Scan (Reference Scan)                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>               | In the Single trigger mode, Sync Selection allows for the choice of triggering on R1 or R2.<br><br>In the Dual trigger mode, Sync Selection toggles trackball control of the S1 or S2 delays displayed in the upper left corner of the monitor. | Reference Scan provides the ability to display both an active and trigger updated image simultaneously in the Dual B-Mode display format. Reference Scan functions in the Single trigger Mode only. A real-time image can not be obtained in Dual Trigger Mode.<br><br>In dual image mode, press <b>Left</b> to display a trigger updated image on the left and active image on the right. Press <b>Right</b> to display a trigger updated image on the right and active image on the left. |
| <b>Benefits</b>                  | Provides trackball control of delay and selection of update trigger point.                                                                                                                                                                      | Provides the ability to view the real-time image while simultaneously acquiring images at a trigger point on the ECG, PCG or Aux waveform.<br><br>These images are best recorded on video tape.                                                                                                                                                                                                                                                                                             |
| <b>Values</b>                    | R1 or R2 trigger select or variable delay adjustment.                                                                                                                                                                                           | The Reference Scan function is on when the Sub-Menu selection is highlighted.<br><br>The function is off when the Sub-Menu selection is not highlighted.<br><br><i>NOTE: The real-time image will only update at each trigger point. The image is paused between the trigger points.</i>                                                                                                                                                                                                    |
| <b>Affects on Other Controls</b> | Changes trigger point in Single Trigger Mode. Changes Trackball control of delays in Dual Mode.                                                                                                                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |

|                    | R Delay                                                                                                                                                                                                                                                                          | ECG Wave, PCG Wave, Aux Wave                                                                     |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|
| <b>Description</b> | R Delay assigns Trackball Control to the adjustment of the trigger delays for R1(S1) or R2 (S2). The graphic display of the delay is found in the upper left corner of the monitor.                                                                                              | These 3 menu selections provide the ability to turn each waveform on for display on the monitor. |
| <b>Benefits</b>    | Allows for the easy adjustment of the R1 (S1) or R2 (S2) delays while scanning.                                                                                                                                                                                                  | Allows for the choice to display a waveform.                                                     |
| <b>Values</b>      | When R Delay is highlighted, the Trackball controls the adjustment of the delays.<br><br>If R Delay is not highlighted, the delays cannot be adjusted with the Trackball.<br><br><i>NOTE: R Delay is automatically activated when SINGLE or DUAL trigger modes are selected.</i> | On or Off.                                                                                       |



## ECG Option

### ECG Gain Pages 2 and 3

| B                                                                                    | Preset                                                                            | Set Up                                                                               | ECG                                                                               |
|--------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|
| ECG Gain                                                                             | ECG Positn                                                                        | PCG Gain                                                                             | PCG Positn                                                                        |
|  10 |  |  10 |  |

Page 2

| B                                                                                    | Preset                                                                            | Set Up | ECG |
|--------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|--------|-----|
| Aux Gain                                                                             | Aux Positn                                                                        |        |     |
|  10 |  |        |     |

Page 3

Figure 10–39. ECG Sub-Menus Pages 2 and 3

|                    | ECG Gain, PCG Gain, Aux Gain                                                                                                                                                                                                                              | ECG Positn, PCG Positn, Aux Positn (Position)                                                                                                         |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Allows for the amplitude control of the ECG, PCG or Aux waveform.                                                                                                                                                                                         | Allows for the vertical positioning of the ECG, PCG or Aux waveform on the image display.                                                             |
| <b>Benefits</b>    | Allows for amplitude adjustment to compensate for different levels of ECG, PCG or Aux output.                                                                                                                                                             | Allows for positioning of the ECG, PCG or Aux waveform to minimize impact on the scan image. Adjustable to user requirements.                         |
| <b>Values</b>      | The PCG & AUX gain can be adjusted from –10 to +10 in 2-digit increments. ECG gain can be adjusted from –20 to +10 in 2-digit increments.<br><i>NOTE: ECG, PCG and Aux default gain setting can be changed in the Set Up Custom Display menu page 17.</i> | Changes display position of the waveform.<br><i>NOTE: ECG, PCG and Aux default position can be changed in the Set Up Custom Display menu page 17.</i> |

### ECG/Cine Gauge/Image Tracking

When an ECG waveform is displayed with the cine gauge, an arrow pointer will appear above the ECG waveform.

As the cine gauge marker is moved with the Cine Scroll control, the arrow above the ECG waveform will move to indicate where on the ECG cycle the displayed image was taken.



# Advanced Cardiac Calculations

## (AMCAL option)

### Overview

This option provides additional measurement, calculation and report capabilities to the Left-Ventricular calculations not found in the basic cardiac option package.

Added to the basic package is the Gibson, Single Plane–DISC and Bi-Plane–DISC methods of LV calculations. A page of measurements is also added to the cardiac calculation menu. This page varies with scan mode. Additional scan mode measurement categories are:

|        |                                                                                                                                                         |
|--------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| B-Mode | Parasternal Long Axis (PLAX), Parasternal Short Axis (PSAX) Aortic Valve, PSAX Mitral Valve, PSAX Papillary Muscles, Apical 4 chamber, Apical 2 chamber |
| M-Mode | Left/Right Ventricles, Mitral Valve, Aortic Valve, Pulmonic Valve, Tricuspid Valve                                                                      |
| D-Mode | Mitral Valve, Aortic Valve, Pulmonic Valve, Tricuspid Valve                                                                                             |

The number and type of measurements vary for each calculation. The formulas used are shown in the calculation specification tables in this chapter.





## Advanced Cardiac Calculations

---

### Measurement Sequences

Each of the titles found in the first layer of the cardiac menu consists of a sequence of measurements. The second layer Sub-Menu consists of individual measurements/calculations that can be performed in sequence or individually.

If the measurements are performed in sequence, the display automatically returns to the first layer menu when all are complete.

If the second layer Sub-Menu measurements are performed individually, only those items measured are recorded on the report page.

### Sequence Philosophy

The primary philosophy of the cardiac calculation package is for the user to select an item from the first menu layer. The system then prompts the user to perform a series of measurements in sequence. This measurement sequence can be modified by the user if the factory sequence is not satisfactory.

This operation should not be confused with the Auto Sequence top menu selection. Items programmed into an auto sequence can consist of first layer sequences as well as individual measurements. Together these form a sequence of measurements and calculations that may be performed automatically.



### Hints

If a measurement in a sequence has already been made, it may be skipped. The system will invoke the next measurement in the sequence.

Pressing the top of the **Ellipse** rocker switch skips a measurement in the sequence.

Pressing the bottom of the **Ellipse** rocker switch moves back to measurements previously skipped.

Press the **Clear** key to quit a sequence at any time. The final calculation result can not be obtained.



### Re-measurement

#### Overview

Values obtained during Advanced Cardiac calculations or after a calculation sequence has been completed can be re-measured as long as those measurement results remain displayed on the screen.

However, once the measured values are deleted from the display by the CLEAR function or by scrolling off the top of the measurement display area, re-measurement is not possible.

#### Operation Method

1. Complete the calculation or calculation sequence.
2. Measurements and resultant calculations will be displayed in the designation area. Only those measurements displayed are available for re-measurement.
3. Press the **Measurement** key and the top of the **Ellipse** rocker switch. The first measurement available to be re-measured will be highlighted and the name is duplicated at the bottom of the measurement display area.
4. Use the **Ellipse** rocker switch to select the desired value to re-measure.

Calculations and generic measurements are not available to be re-measured and will be skipped in the selection process.



*NOTE: The re-measurement function is not available with the Single Plane DISC and Bi-Plane DISC methods.*



## Advanced Cardiac Calculations

---

### Operation Method (cont'd)

5. Perform the necessary measurement steps for the selected value using the Trackball and Set key.

When the **Set** key is pressed to complete the re-measurement:

- The resultant value of the selected item and all related calculations will be recalculated and displayed on the screen.
- The new result will be entered on the report page.
- The old measurement is deleted from the screen and report page as it is replaced by the new value.

If additional re-measurements need to be performed, repeat the previous steps.



### Hints

Each value of an M-Mode continuous sequence can be remeasured separately.

If a measurement is selected from the soft-menu during the re-measure process, the re-measurement process is cancelled and the old value is restored.

After a re-measurement is complete (fixed), the old graphic is erased and the same cursor symbols are used for the remeasured value.

If two or more identical measurement items are displayed in the calculation result area, only the last one supports re-measurement.

The re-measurement of Heart Rate is always manual, even though Auto Heart Rate is selected as a preset.

There is no automatic compare function for Systole/Diastole.

It is not possible to do re-measurements on VCR playback values in Image Memory.

Re-measurement is not possible during a calculation sequence. Only after the sequence is compiled.

Once cleared or erased, measurements cannot be recalled for re-measurement.



## Automatic Determination of Systole and Diastole

The system can be programmed to automatically determine the cardiac phase as systole or diastole. The parameter in the Setup/Preset Program Sub-Menu page 3 used to make this selection is:

Diastole/Systole Determination : Manual    Auto

The auto determination of systole/diastole can be fine tuned using the following parameters also found in Set Up/Preset Program page 3:

R Delay Time of End Systole Prior Edge : 0–2000msec  
End Systole Period from Prior Edge : 0–2000msec

### Automatic

If Auto is selected, the determination is made by comparing the displayed B-Mode or M/D Mode image frames with the R-point on the ECG wave. Therefore, 's' for systole and 'd' for diastole is added automatically to the measurement or calculation name.

If phase cannot automatically be determined, the system displays the following message:

**“Diastole ('d') or Systole ('s') ?”**

Entering 'd' or 's' are the only two acceptable inputs.

|     |                                                                                                                     |
|-----|---------------------------------------------------------------------------------------------------------------------|
| 's' | Erases the message, recognizes the phase as systole and adds a 's' to the end of the measurements or calculations.  |
| 'd' | Erases the message, recognizes the phase as diastole and adds a 'd' to the end of the measurements or calculations. |

The user will not be able to perform measurements of the opposite phase until the next image change is made.



## Advanced Cardiac Calculations

---

### Automatic (cont'd)

When the measurement phase requested by the system is different than the phase of the cursor position, the following message is displayed:

**“Display ##### image”**

##### represents the phase opposite of the one just completed.

Auto determination of systole/diastole is only available on a frozen image. It cannot be done on a recalled image or one from an external video input.

### Manual

If Manual is selected, each measurement has a ‘s’ or ‘d’ at the end of its name. The system will not determine the phase or ask the operator to make that determination.

### Auto Trace Measurements

Doppler measurement menus (D-AV, D-MV, D-PV and D-TV) each have a selection called Auto Trace. When the user selects Auto Trace, the start point and end point of a measurement is set. After the end point is set, the system traces the waveform and automatically calculates the values as follows:

Ao Auto Trace: (Set start and end points for FVI-AV)  
Calculates: FVI-AV, MaxPG, MeanPG, Ejection Time, AccT, DecT, PFV and Mean-V

MV Auto Trace: (Set start and end points for FVI)  
Calculates: FVI, MaxPG, MeanPG, Pressure Half Time, MV Area, AccT, DecT, PFV, Mean-V and Ejection Time

PV Auto Trace: (Set start and end points for FVI-PV)  
Calculates: FVI-PV, MaxPG, MeanPG, AccT, DecT, PFV, Mean-V and Ejection Time

TV Auto Trace: (Set start and end points for FVI)  
Calculates: FVI, MaxPG, MeanPG, Pressure Half Time, AccT, DecT, PFV, Mean-V and Ejection Time

## Continuous M-Mode Measurements

When M-Mode measurements can be taken on the same timeline, the system will support a continuous measurement of all values along that line.

This can be done because the end point of one measurement is the start point of the next measurement. After the end point of the last measurement is set, keystrokes are reduced by not having to set the next start point.




Figure 10-40. Continuous M-Mode Measurements

This type of M-Mode measurements is supported by M-LV/RV, Teichholz, Cubed and Gibson menu selections.

After selecting the measurement sequence, use the **Trackball** to position the first cursor and press **Set**. The name of the first value appears; use the **Trackball** to position the cursor at the end of the first measurement and press **Set**. This is the start point of the next measurement. The name of the next measurement appears; use the **Trackball** to position the end point of the second measurement and press **Set**.

This process continues until all measurements are taken.



## Advanced Cardiac Calculations

---

### Continuous M-Mode Measurements (cont'd)



#### Hints

The **Ellipse** rocker switch is used to skip measurements or move back to measurements that were skipped.

Before the end point is set, uncompleted measurements can be removed by pressing **Clear** once. That measurement can then be started again.

Measurements can be taken in Zoom mode as long as the magnification is not changed during the measurement process.

Measurements can be programmed not to appear in the continuous measurement sequence in the Set Up/Preset Program menu pages 5 and 6.

If the measurement is set to OFF in the Preset Program menus, it will not appear in the continuous sequence. Additional Set points will have to be made if measurements are skipped in the continuous sequence.

### Advanced Cardiac Calculations Measurement Menus

The illustrations that follow show the Advanced Cardiac Calculation measurement Sub-Menus that are available with the AMCAL option.

Each selection in the first layer of Sub-Menus will yield a second layer of Sub-Menus. This second layer contains measurements and calculations necessary to complete the category or goal selected in the first layer.

The illustrations shown are those that result if the determination of systole/diastole is set to Automatic and then Manual.



---

## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (First Layer)

|             |               |              |                |
|-------------|---------------|--------------|----------------|
|             |               |              |                |
| LV<br>CUBED | LV<br>SP-DISC | LV<br>SP-ELP | LV<br>GIBSON   |
| LV<br>TEICH | LV<br>BP-DISC | LV<br>BP-ELP | CARD<br>Report |

Page One

|               |  |  |  |
|---------------|--|--|--|
|               |  |  |  |
| LV<br>BULLET  |  |  |  |
| LV<br>SIMPSON |  |  |  |

Page Two

|             |              |        |                |
|-------------|--------------|--------|----------------|
|             |              |        |                |
| PLAX        | PSAX<br>-MV  | AP-4CH | TRACE<br>AUTO  |
| PSAX<br>-AV | PSAX<br>-PAP | AP-2CH | CARD<br>Report |

Page Three

|      |      |  |                |
|------|------|--|----------------|
|      |      |  |                |
| D-MV | D-PV |  | TRACE<br>AUTO  |
| D-AV | D-TV |  | CARD<br>Report |

Page Four

|         |      |      |                |
|---------|------|------|----------------|
|         |      |      |                |
| M-LV/RV | M-AV | M-TV |                |
| M-MV    | M-PV |      | CARD<br>Report |

Page Five

Illustration 1. First Layer Sub-Menus





## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Automatic Determination of Systole/Diastole

|      |      |      |             |
|------|------|------|-------------|
|      |      |      |             |
| MEAS | IVS  | LVPW | Return      |
| RVID | LVID | ET   | CARD Report |

Page One

|     |            |     |             |
|-----|------------|-----|-------------|
|     |            |     |             |
|     | FS         | EdV | Return      |
| LVM | HEART RATE | EsV | CARD Report |

Page Two

|    |      |    |             |
|----|------|----|-------------|
|    |      |    |             |
| SV | EF   | SI | Return      |
| CO | MVCF | CI | CARD Report |

Page Three

|  |        |        |             |
|--|--------|--------|-------------|
|  |        |        |             |
|  | IVS/PW | %STIVS | Return      |
|  |        | %STPW  | CARD Report |

Page Four

Illustration 2. Cubed, Teichholz and Gibson Sub-Menus

|      |  |  |             |
|------|--|--|-------------|
|      |  |  |             |
| LVL  |  |  | Return      |
| LVAM |  |  | CARD Report |

Page One

|  |            |     |             |
|--|------------|-----|-------------|
|  |            |     |             |
|  |            | EdV | Return      |
|  | HEART RATE | EsV | CARD Report |

Page Two

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO |    | CI | CARD Report |

Page Three

Illustration 3. LV Bullet Sub-Menus



---

## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Automatic Determination of Systole/Diastole (cont'd)

|     |  |  |             |
|-----|--|--|-------------|
|     |  |  |             |
| LVL |  |  | Return      |
| LVA |  |  | CARD Report |

Page One

|  |            |     |             |
|--|------------|-----|-------------|
|  |            |     |             |
|  |            | EdV | Return      |
|  | HEART RATE | EsV | CARD Report |

Page Two

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO |    | CI | CARD Report |

Page Three

Illustration 4. LV SP-DISC Sub-Menus

|      |      |  |             |
|------|------|--|-------------|
|      |      |  |             |
| LVL2 | LVA2 |  | Return      |
| LVL4 | LVA4 |  | CARD Report |

Page One

|  |            |     |             |
|--|------------|-----|-------------|
|  |            |     |             |
|  |            | EdV | Return      |
|  | HEART RATE | EsV | CARD Report |

Page Two

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO |    | CI | CARD Report |

Page Three

Illustration 5. LV BP-DISC Sub-Menus



## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Automatic Determination of Systole/Diastole (cont'd)

|      |      |  |             |
|------|------|--|-------------|
|      |      |  |             |
| LVL  | LVAP |  | Return      |
| LVAM |      |  | CARD Report |

Page One

|  |            |     |             |
|--|------------|-----|-------------|
|  |            |     |             |
|  |            | EdV | Return      |
|  | HEART RATE | EsV | CARD Report |

Page Two

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO |    | CI | CARD Report |

Page Three

Illustration 6. Modified Simpson's Rule Sub-Menus

|     |  |  |             |
|-----|--|--|-------------|
|     |  |  |             |
| LVL |  |  | Return      |
| LVA |  |  | CARD Report |

Page One

|  |            |     |             |
|--|------------|-----|-------------|
|  |            |     |             |
|  |            | EdV | Return      |
|  | HEART RATE | EsV | CARD Report |

Page Two

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO |    | CI | CARD Report |

Page Three

Illustration 7. Single Plane Ellipsoid Sub-Menus



---

## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Automatic Determination of Systole/Diastole (cont'd)

|      |     |  |             |
|------|-----|--|-------------|
|      |     |  |             |
| LVML | LVA |  | Return      |
| LVAM |     |  | CARD Report |

Page One

|  |            |     |             |
|--|------------|-----|-------------|
|  |            |     |             |
|  |            | EdV | Return      |
|  | HEART RATE | EsV | CARD Report |

Page Two

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO |    | CI | CARD Report |

Page Three

Illustration 8. Biplane Ellipsoid Sub-Menus

|      |       |      |             |
|------|-------|------|-------------|
|      |       |      |             |
| AO   | LVID  | ALSs | Return      |
| IVSd | LVPWd | LADs | CARD Report |

Page One

|       |       |            |             |
|-------|-------|------------|-------------|
|       |       |            |             |
| LVOTs | AAs   |            | Return      |
|       | LVOTA | HEART RATE | CARD Report |

Page Two

Illustration 9. PLAX Sub-Menus



## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Automatic Determination of Systole/Diastole (cont'd)

|      |       |  |             |
|------|-------|--|-------------|
|      |       |  |             |
| PAD  | AO    |  | Return      |
| RVOT | LADML |  | CARD Report |

Page One

|      |       |            |             |
|------|-------|------------|-------------|
|      |       |            |             |
| PAAs | RVOTA |            | Return      |
| AAs  |       | HEART RATE | CARD Report |

Page Two

Illustration 10. PSAX-AV Sub-Menus

|      |  |            |             |
|------|--|------------|-------------|
|      |  |            |             |
| MVOA |  |            | Return      |
| LVAM |  | HEART RATE | CARD Report |

Page One

Illustration 11. PSAX-MV Sub-Menus

|      |  |            |             |
|------|--|------------|-------------|
|      |  |            |             |
| LVAP |  |            | Return      |
|      |  | HEART RATE | CARD Report |

Page One

Illustration 12. PSAX-PAP Sub-Menus



---

## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Automatic Determination of Systole/Diastole (cont'd)

|     |      |      |             |
|-----|------|------|-------------|
|     |      |      |             |
| LVA | LVML | RVL  | Return      |
| LVL | RVA  | RVML | CARD Report |

Page One

|      |       |            |             |
|------|-------|------------|-------------|
|      |       |            |             |
| TAML | LADML | RAA        | Return      |
| LAAs | LADSI | HEART RATE | CARD Report |

Page Two

|       |  |  |             |
|-------|--|--|-------------|
|       |  |  |             |
| RADML |  |  | Return      |
| RADSI |  |  | CARD Report |

Page Three

Illustration 13. AP-4CH Sub-Menus

|     |      |            |             |
|-----|------|------------|-------------|
|     |      |            |             |
| LVA | LVML |            | Return      |
| LVL |      | HEART RATE | CARD Report |

Page One

Illustration 14. AP-2CH Sub-Menus



## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Automatic Determination of Systole/Diastole (cont'd)

|      |      |      |             |
|------|------|------|-------------|
|      |      |      |             |
| MEAS | RVDd | LVID | Return      |
|      | IVS  | LVPW | CARD Report |

Page One

|        |        |            |             |
|--------|--------|------------|-------------|
|        |        |            |             |
| IVS/PW | %STIVS | LVM        | Return      |
| FS     | %STPW  | HEART RATE | CARD Report |

Page Two

Illustration 15. M-LV/RV Sub-Menus

|      |      |      |             |
|------|------|------|-------------|
|      |      |      |             |
| EPSS | V_DE | D_CE | Return      |
| V_EF | D_DE | T_AC | CARD Report |

Page One

|     |       |            |             |
|-----|-------|------------|-------------|
|     |       |            |             |
| P-R | A/E   |            | Return      |
|     | PR-AC | HEART RATE | CARD Report |

Page Two

Illustration 16. M-MV Sub-Menus

|     |       |       |             |
|-----|-------|-------|-------------|
|     |       |       |             |
| AOd | LADs  | LVPEP | Return      |
| ALS | RVOTs | LVET  | CARD Report |

Page One

|        |  |            |             |
|--------|--|------------|-------------|
|        |  |            |             |
| LA/AO  |  |            | Return      |
| PEP/ET |  | HEART RATE | CARD Report |

Page Two

Illustration 17. M-AV Sub-Menus



---

## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Automatic Determination of Systole/Diastole (cont'd)

|       |       |        |             |
|-------|-------|--------|-------------|
|       |       |        |             |
| PADs  | RVET  | PAA's  | Return      |
| RVPEP | aWave | PEP/ET | CARD Report |

Page One

|  |  |            |             |
|--|--|------------|-------------|
|  |  |            |             |
|  |  |            | Return      |
|  |  | HEART RATE | CARD Report |

Page Two

Illustration 18. M-PV Sub-Menus

|       |       |        |             |
|-------|-------|--------|-------------|
|       |       |        |             |
| TV_DE | TV_AC | TV-A/E | Return      |
| TV_EF | P-R   | TPR-AC | CARD Report |

Page One

|  |  |            |             |
|--|--|------------|-------------|
|  |  |            |             |
|  |  |            | Return      |
|  |  | HEART RATE | CARD Report |

Page Two

Illustration 19. M-TV Sub-Menus

|      |       |        |             |
|------|-------|--------|-------------|
|      |       |        |             |
| MVOA | PFVLA | FVI-MV | Return      |
|      | PFVMV |        | CARD Report |

Page One

|        |     |               |             |
|--------|-----|---------------|-------------|
|        |     |               |             |
| ET-MV  | MVA | MV Auto Trace | Return      |
| PHT-MV |     | HEART RATE    | CARD Report |

Page Two

|           |         |            |             |
|-----------|---------|------------|-------------|
|           |         |            |             |
| MeanV -MV | AccT-MV | MaxPG -MV  | Return      |
| PkFVI -MV | DecT-MV | MeanPG -MV | CARD Report |

Page Three

|  |       |       |             |
|--|-------|-------|-------------|
|  |       |       |             |
|  | SV-MV | SI-MV | Return      |
|  | CO-MV | CI-MV | CARD Report |

Page Four

Illustration 20. D-MV Sub-Menus





## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Automatic Determination of Systole/Diastole (cont'd)

|     |       |        |             |
|-----|-------|--------|-------------|
|     |       |        |             |
| AOs | PFVAV | FVI-AV | Return      |
|     | PFVLA | FVI-LA | CARD Report |

Page One

|        |       |               |             |
|--------|-------|---------------|-------------|
|        |       |               |             |
| ET-AV  | AV_A  | Ao Auto Trace | Return      |
| PHT-AV | LVOTA | HEART RATE    | CARD Report |

Page Two

|          |         |           |             |
|----------|---------|-----------|-------------|
|          |         |           |             |
| MeanV-AV | AccT-AV | MaxPG-AV  | Return      |
| PkFVI-AV | DecT-AV | MeanPG-AV | CARD Report |

Page Three

|       |       |       |             |
|-------|-------|-------|-------------|
|       |       |       |             |
| Qp:Qs | SV-AV | SI-AV | Return      |
|       | CO-AV | CI-AV | CARD Report |

Page Four

Illustration 21. D-AV Sub-Menus

|      |       |        |             |
|------|-------|--------|-------------|
|      |       |        |             |
| TAML | PFVRA | FVI-TV | Return      |
|      | PFVTV |        | CARD Report |

Page One

|        |  |               |             |
|--------|--|---------------|-------------|
|        |  |               |             |
| ET-TV  |  | TV Auto Trace | Return      |
| PHT-TV |  | HEART RATE    | CARD Report |

Page Two

|          |         |           |             |
|----------|---------|-----------|-------------|
|          |         |           |             |
| MeanV-TV | AccT-TV | MaxPG-TV  | Return      |
| PkFVI-TV | DecT-TV | MeanPG-TV | CARD Report |

Page Three

|  |       |       |             |
|--|-------|-------|-------------|
|  |       |       |             |
|  | SV-TV | SI-TV | Return      |
|  | CO-TV | CI-TV | CARD Report |

Page Four

Illustration 22. D-TV Sub-Menus



---

## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Automatic Determination of Systole/Diastole (cont'd)

|        |       |        |             |
|--------|-------|--------|-------------|
|        |       |        |             |
| PADs   | PFVPA | FVI-PV | Return      |
| FVI-RV | PFVRV | FVI-PA | CARD Report |

Page One

|        |      |               |             |
|--------|------|---------------|-------------|
|        |      |               |             |
| ET-PV  | PV_A | PV Auto Trace | Return      |
| PHT-PV |      | HEART RATE    | CARD Report |

Page Two

|          |         |           |             |
|----------|---------|-----------|-------------|
|          |         |           |             |
| MeanV-PV | AccT-PV | MaxPG-PV  | Return      |
| PkFVI-PV | DecT-PV | MeanPG-PV | CARD Report |

Page Three

|       |       |       |             |
|-------|-------|-------|-------------|
|       |       |       |             |
| Qp:Qs | SV-PV | SI-PV | Return      |
|       | CO-PV | CI-PV | CARD Report |

Page Four

Illustration 23. D-PV Sub-Menus



## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Manual Determination of Systole/Diastole

|       |       |      |             |
|-------|-------|------|-------------|
|       |       |      |             |
| MEASd | LVIDd | IVSd | Return      |
| MEASs | LVIDs | IVSs | CARD Report |

Page One

|       |            |     |             |
|-------|------------|-----|-------------|
|       |            |     |             |
| LVPWd | ET         | EdV | Return      |
| LVPWs | HEART RATE | EsV | CARD Report |

Page Two

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO | FS | CI | CARD Report |

Page Three

|      |        |        |             |
|------|--------|--------|-------------|
|      |        |        |             |
| MVCF | IVS/PW | %STIVS | Return      |
| LVM  | RVIDd  | %STPW  | CARD Report |

Page Four

Illustration 24. Cubed, Teichholz and Gibson Sub-Menus

|      |       |  |             |
|------|-------|--|-------------|
|      |       |  |             |
| LVLd | LVAMd |  | Return      |
| LVLs | LVAMs |  | CARD Report |

Page One

|  |            |     |             |
|--|------------|-----|-------------|
|  |            |     |             |
|  |            | EdV | Return      |
|  | HEART RATE | EsV | CARD Report |

Page Two

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO |    | CI | CARD Report |

Page Three

Illustration 25. Bullet Sub-Menus



---

## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Manual Determination of Systole/Diastole (cont'd)

|      |       |  |             |
|------|-------|--|-------------|
|      |       |  |             |
| LVLd | LVA d |  | Return      |
| LVLs | LVA s |  | CARD Report |

Page One

|  |            |     |             |
|--|------------|-----|-------------|
|  |            |     |             |
|  |            | EdV | Return      |
|  | HEART RATE | EsV | CARD Report |

Page Two

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO |    | CI | CARD Report |

Page Three

Illustration 26. LV SP-DISC Sub-Menus

|       |       |  |             |
|-------|-------|--|-------------|
|       |       |  |             |
| LVL2d | LVA2d |  | Return      |
| LVL2s | LVA2s |  | CARD Report |

Page One

|       |       |  |             |
|-------|-------|--|-------------|
|       |       |  |             |
| LVL4d | LVA4d |  | Return      |
| LVL4s | LVA4s |  | CARD Report |

Page Two

|  |            |     |             |
|--|------------|-----|-------------|
|  |            |     |             |
|  |            | EdV | Return      |
|  | HEART RATE | EsV | CARD Report |

Page Three

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO |    | CI | CARD Report |

Page Four

Illustration 27. LV BP-DISC Sub-Menus



## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Manual Determination of Systole/Diastole (cont'd)

|      |       |       |             |
|------|-------|-------|-------------|
|      |       |       |             |
| LVLd | LVAMd | LVAPd | Return      |
| LVLs | LVAMs | LVAPs | CARD Report |

Page One

|  |            |     |             |
|--|------------|-----|-------------|
|  |            |     |             |
|  |            | EdV | Return      |
|  | HEART RATE | EsV | CARD Report |

Page Two

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO |    | CI | CARD Report |

Page Three

Illustration 28. Modified Simpson's Rule Sub-Menus

|      |       |  |             |
|------|-------|--|-------------|
|      |       |  |             |
| LVLd | LVAAd |  | Return      |
| LVLs | LVAAs |  | CARD Report |

Page One

|  |            |     |             |
|--|------------|-----|-------------|
|  |            |     |             |
|  |            | EdV | Return      |
|  | HEART RATE | EsV | CARD Report |

Page Two

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO |    | CI | CARD Report |

Page Three

Illustration 29. Single Plane Ellipsoid Sub-Menus



---

## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Manual Determination of Systole/Diastole (cont'd)

|       |       |      |             |
|-------|-------|------|-------------|
|       |       |      |             |
| LVMLd | LVAMd | LVAd | Return      |
| LVMLs | LVAMs | LVAs | CARD Report |

Page One

|  |            |     |             |
|--|------------|-----|-------------|
|  |            |     |             |
|  |            | EdV | Return      |
|  | HEART RATE | EsV | CARD Report |

Page Two

|    |    |    |             |
|----|----|----|-------------|
|    |    |    |             |
| SV | EF | SI | Return      |
| CO |    | CI | CARD Report |

Page Three

Illustration 30. Biplane Ellipsoid Sub-Menus

|      |       |      |             |
|------|-------|------|-------------|
|      |       |      |             |
| AOd  | LVIDd | AOs  | Return      |
| IVSd | LVPWd | ALSs | CARD Report |

Page One

|       |       |            |             |
|-------|-------|------------|-------------|
|       |       |            |             |
| LVIDs | LVOTs | AAs        | Return      |
| LADs  | LVOTA | HEART RATE | CARD Report |

Page Two

Illustration 31. PLAX Sub-Menus Manual Determination



## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Manual Determination of Systole/Diastole (cont'd)

|      |       |      |             |
|------|-------|------|-------------|
|      |       |      |             |
| AOd  | RVOTd | AOs  | Return      |
| PADd |       | PADs | CARD Report |

Page One

|       |       |            |             |
|-------|-------|------------|-------------|
|       |       |            |             |
| LADML | RVOTA | AAs        | Return      |
| RVOTs | PAAs  | HEART RATE | CARD Report |

Page Two

Illustration 32. PSAX-AV Sub-Menus Manual Determination

|       |  |  |             |
|-------|--|--|-------------|
|       |  |  |             |
| MVOA  |  |  | Return      |
| LVAMd |  |  | CARD Report |

Page One

|       |  |            |             |
|-------|--|------------|-------------|
|       |  |            |             |
|       |  |            | Return      |
| LVAMs |  | HEART RATE | CARD Report |

Page Two

Illustration 33. PSAX-MV Sub-Menus Manual Determination

|       |  |  |             |
|-------|--|--|-------------|
|       |  |  |             |
| LVAPd |  |  | Return      |
|       |  |  | CARD Report |

Page One

|       |  |            |             |
|-------|--|------------|-------------|
|       |  |            |             |
| LVAPs |  |            | Return      |
|       |  | HEART RATE | CARD Report |

Page Two

Illustration 34. PSAX-PAP Sub-Menus Manual Determination



Advanced Cardiac Calculations

AMCAL Sub-Menus (Second Layer)—Manual Determination of Systole/Diastole (cont'd)

|                  |                   |                   |             |
|------------------|-------------------|-------------------|-------------|
|                  |                   |                   |             |
| LVA <sub>d</sub> | LVML <sub>d</sub> | RVL <sub>d</sub>  | Return      |
| LVL <sub>d</sub> | RVA <sub>d</sub>  | RVML <sub>d</sub> | CARD Report |

Page One

|      |                  |            |             |
|------|------------------|------------|-------------|
|      |                  |            |             |
| LVAs | LVMLs            | RVLs       | Return      |
| LVLs | RVA <sub>s</sub> | HEART RATE | CARD Report |

Page Two

|       |       |       |             |
|-------|-------|-------|-------------|
|       |       |       |             |
| RVMLs | LAAs  | LADSI | Return      |
| TAML  | LADML | RAA   | CARD Report |

Page Three

|       |  |  |             |
|-------|--|--|-------------|
|       |  |  |             |
| RADML |  |  | Return      |
| RADSI |  |  | CARD Report |

Page Four

Illustration 35. AP-4CH Sub-Menus Manual Determination

|                  |                   |  |             |
|------------------|-------------------|--|-------------|
|                  |                   |  |             |
| LVA <sub>d</sub> | LVML <sub>d</sub> |  | Return      |
| LVL <sub>d</sub> |                   |  | CARD Report |

Page One

|      |       |            |             |
|------|-------|------------|-------------|
|      |       |            |             |
| LVAs | LVMLs |            | Return      |
| LVLs |       | HEART RATE | CARD Report |

Page Two

Illustration 36. AP-2CH Sub-Menus Manual Determination





## Advanced Cardiac Calculations

---

### AMCAL Sub-Menus (Second Layer)—Manual Determination of Systole/Diastole (cont'd)

|       |      |       |             |
|-------|------|-------|-------------|
|       |      |       |             |
| MEASd | RVDd | LVIDd | Return      |
| MEASs | IVSd | LVPWd | CARD Report |

Page One

|      |       |            |             |
|------|-------|------------|-------------|
|      |       |            |             |
| IVSs | LVIDs | LVM        | Return      |
|      | LVPWs | HEART RATE | CARD Report |

Page Two

|        |        |  |             |
|--------|--------|--|-------------|
|        |        |  |             |
| IVS/PW | %STIVS |  | Return      |
| FS     | %STPW  |  | CARD Report |

Page Three

Illustration 37. M-LV/RV Sub-Menus Manual Determination

## Cardiac Measurements

Measurements taken in the Advanced Cardiac Calculation option are the same type of measurements taken in the basic package.

The user should review and become familiar with how to make measurements as shown in chapters titled *Basic Measurements/Calculations* and *Cardiology* of this manual set.

The system will prompt the user for any measurements necessary to complete a selected sequence.

## Customizing Measurement Sequences

The sequence in which the second layer of measurements/ Cardiology calculations is presented may be customized by the user. This is done by using the Set Up/Preset Program Sub-Menu page 7. The selection is called:

Cardiac Calculation : Submenu

Selecting this preset displays the calculation sequence Set Up menu.

|                                                                         |           |                           |  |
|-------------------------------------------------------------------------|-----------|---------------------------|--|
| HOSPITAL NAME                                                           |           | 05/09/98 06:01:21<br>???? |  |
| <Calculation Sequence Set up Menu>                                      |           |                           |  |
| [1] Select a calculation                                                |           |                           |  |
| CALCULATIONS                                                            | [B Mode   | ][Parasternal Long Axis ] |  |
| AAs<br>LVOTa<br>*PLAX                                                   |           |                           |  |
| [2] Set up the sequence order                                           |           | ←01/01→                   |  |
| MEASUREMENTS                                                            |           |                           |  |
| <01>A0d                                                                 | <07>LADs  |                           |  |
| <02>IVSd                                                                | <08>LVOTs |                           |  |
| <03>LVPWd                                                               | <09>LVIDs |                           |  |
| <04>LVIDd                                                               | <10>HR    |                           |  |
| <05>A0s                                                                 | <11>AAs   |                           |  |
| <06>ALSs                                                                | <12>LVOTa |                           |  |
| [ORDER] [RESET] [SAVE] [EXIT] ]<br>'TRACKBALL' and 'SET' are available. |           |                           |  |

Figure 10–41. Calculation Sequence Set Up Menu



## *Advanced Cardiac Calculations*

---

### **Pop up Menus**

The main menu (B-Mode, M-Mode, etc.) and the Sub-Menu (Aortic Valve, Mitral Valve, etc.) are pop-up menus that provide the ability to quickly display the measurements and their sequence.

To display a pop up menu, use the **Trackball** to move the highlight cursor to the main or Sub-Menu area. Press **Set**. The pop up menu is displayed.

While the pop up menu is displayed, use the up/down **Trackball** motion or the **Up/Down** arrow keys to select the desired menu item.

Use the **Set** key to activate the menu choice.

### **Page Numbers**

The page number area not only shows the Present Page/Total Pages, but this area is used to navigate through the display of available pages.

Use the **Trackball** to move the arrow cursor to highlight the navigation arrows on either side of the page numbers.

Press **Set** to display the previous or next page.



---

## Advanced Cardiac Calculations

---

### New Measurement Sequence

#### ORDER Command

After selecting the desired Main Menu and Sub-Menu, the calculations are displayed in the top half of the menu.

Use the **Trackball** to highlight the desired calculation and press **Set**. Measurement and calculation items are displayed on the bottom half of the menu. The order in which they are executed is shown to their left.

To change the order or omit undesired items from the sequence, use the ORDER command. **Trackball** the highlight cursor to ORDER and press **Set**.

All measurements and calculations on the bottom half of the menu no longer have a number to their left.

Use the **Trackball** to move the highlight cursor to the first measurement to be made in the customized sequence. Press **Set**. The number one is placed to the left of that item. Use the **Trackball** to move the highlighted cursor to the second item in the new sequence and press **Set**.

Continue the process until all desired items are designated.

*NOTE: Items not designated will be deleted from the sequence.*

Use the SAVE command to store the new custom sequence.



#### Hints

- The customized sequence is saved until the ORDER command is used to generate a new sequence.
- Items not designated will not automatically appear as a measurement prompt.

#### RESET Command

The RESET command resets the measurement sequence to the factory defaults. Prior to resetting, the following message is displayed:

**“Reset? Press ‘ESC’ to cancel. ‘SET’ to confirm”**

**Set** resets the order to factory defaults. **Esc** causes nothing to happen.



## Advanced Cardiac Calculations

---

### SAVE Command

The SAVE command will save the changes as the new measurement sequence. Prior to saving, the following message is displayed:

**“Overwrite? Press ‘ESC’ to cancel. ‘SET’ to confirm.”**

**Set** saves the order changes. **Esc** causes nothing to happen.

If the rearrangement process was not completed, the message displayed is:

**“Complete the re-arrangement first.”**

### EXIT Command

If the exit command is selected before choosing, completing or saving the order, the following message is displayed:

**“Quit? Press ‘ESC’ to cancel. ‘SET’ to confirm.”**

**Set** exits from this screen. **Esc** causes nothing to happen.

If a new sequence programming was completed but not saved, the message displayed is:

**“Overwrite? Press ‘ESC’ to cancel. ‘SET’ to confirm.”**

**Set** saves the order changes and exits. **Esc** causes nothing to happen.

### Auto Sequence Programming



#### Hints

If there is a series of measurements/calculations that are routinely taken in a specific sequence, the auto sequence feature of the system can be used to program these measurements/calculations to automatically appear on the display in the desired order. Refer to 14–81 for details on programming an auto sequence.

## Advanced Cardiac Specification Tables

The following tables show each first layer menu selection with the measurements and calculations found in it's second layer menu.

The tables show the mnemonic, a description of that mnemonic (calc name), input measurements required and formula used.

### LV Calculation Formulas (Cubed Method)




Figure 10-42. Cubed Method Measurements

Refer to 10-7 for the following Cubed Method Formulas:

LVIDd, LVIDs, IVSd, IVSs, LVPWd, LVPWs, HR, ET, EdV, EsV, FS, SV, EF and CO.

| Calc Mnemonic | Calc Name                   | Input Measurements                | Formula                                           |
|---------------|-----------------------------|-----------------------------------|---------------------------------------------------|
| SI            | Stroke Volume Index         | body surface area and SV          | $SI[ml/m^2]=SV/BSA$                               |
| CI            | Cardiac Index               | body surface area and CO          | $CI[ml/min/m^2]=CO/BSA$                           |
| MVCF          | Mean Vcf                    | two distances & one time interval | $MVCF[circ/s]=(LVIDd-LVIDs)/(LVIDd \times ET)$    |
| LVM           | Left Ventricle Cardiac Mass | three distances                   | $LVM[g]=1.04*[(IVSd+LVPWd+LVIDd)^3-LVIDd^3]-13.6$ |

Table 10-8. Cubed Method Calculation Formulas

## LV Calculation Formulas (Teichholz Method)




Figure 10-43. Teichholz Method Measurements

Refer to 10-9 for the following Teichholz Method Formulas:

LVIDd, LVIDs, IVSd, IVSs, LVPWd, LVPWs, HR, ET, EdV, EsV, FS, SV, EF and CO.

| Calc Mnemonic | Calc Name                   | Input Measurements                | Formula                                                                                                 |
|---------------|-----------------------------|-----------------------------------|---------------------------------------------------------------------------------------------------------|
| SI            | Stroke Volume Index         | body surface area and SV          | $SI[\text{ml}/\text{m}^2] = \text{SV}/\text{BSA}$                                                       |
| CI            | Cardiac Index               | body surface area and CO          | $CI[\text{ml}/\text{min}/\text{m}^2] = \text{CO}/\text{BSA}$                                            |
| MVCF          | Mean Vcf                    | two distances & one time interval | $\text{MVCF}[\text{circ}/\text{s}] = (\text{LVIDd} - \text{LVIDs}) / (\text{LVIDd} \times \text{ET})$   |
| LVM           | Left Ventricle Cardiac Mass | three distances                   | $\text{LVM}[\text{g}] = 1.04 * [(\text{IVSd} + \text{LVPWd} + \text{LVIDd})^3 - \text{LVIDd}^3] - 13.6$ |

Table 10-9. Teichholz Method Calculation Formulas

## LV Calculation Formulas (Bullet Method)




Figure 10-44. Bullet Method Measurements

Refer to 10-11 for the following Bullet Method Formulas:

LVLd, LVLs, LVAMd, LVAMs, HR, EdV, EsV, SV, EF and CO.

| Calc Mnemonic | Calc Name           | Input Measurements       | Formula                 |
|---------------|---------------------|--------------------------|-------------------------|
| SI            | Stroke Volume Index | body surface area and SV | $SI[ml/m^2]=SV/BSA$     |
| CI            | Cardiac Index       | body surface area and CO | $CI[ml/min/m^2]=CO/BSA$ |

Table 10-10. Bullet Method Calculation Formulas

## LV Calculation Formulas (LV SP-DISC Method)




Figure 10-45. LV SP-DISC Method Measurements

| Calc Mnemonic | Calc Name                         | Input Measurements | Formula                      |
|---------------|-----------------------------------|--------------------|------------------------------|
| LVLd          | Left Ventricular Length, Diastole | one distance       | $LVLd=d1[cm \text{ or } mm]$ |
| LVLs          | Left Ventricular Length, Systole  | one distance       | $LVLs=d1[cm \text{ or } mm]$ |

Table 10-11. LV SP-DISC Method Formulas





## Advanced Cardiac Calculations

### LV Calculation Formulas (LV SP-DISC Method) (cont'd)

| Calc Mnemonic | Calc Name                       | Input Measurements                                                                     | Formula                                                                |
|---------------|---------------------------------|----------------------------------------------------------------------------------------|------------------------------------------------------------------------|
| LVAd          | Left Ventricular Area, Diastole | one area (by trace only) is divided into 20 elliptical disk segments                   | $LVAd = a1[\text{cm}^2]$                                               |
| LVAs          | Left Ventricular Area, Systole  | one area (by trace only) is divided into 20 elliptical disk segments                   | $LVAs = a1[\text{cm}^2]$                                               |
| HR            | Heart Rate (beats/minute)       | one 2 beat time interval                                                               | $HR[\text{BPM}] = 120 [\text{sec}] / 2 \text{ beat time} [\text{sec}]$ |
| EdV           | End Diastole Volume             | one distance and two areas (by ellipse, trace or circle)                               | $EdV = \sum_{i=1}^{20} (\pi/4) * ((LVLd/20) a_{di}^2) \quad (*1)$      |
| EsV           | End Systole Volume              | one distance and two areas (by ellipse, trace or circle)                               | $EsV = \sum_{i=1}^{20} (\pi/4) * ((LVLs/20) a_{si}^2) \quad (*2)$      |
| SV            | Stroke Volume                   | two distances and two areas (by ellipse, trace or circle)                              | $SV[\text{ml}] = EdV - EsV$                                            |
| EF            | Ejection Fraction               | two distances and two areas (by ellipse, trace or circle)                              | $EF = SV / EdV \times 100$                                             |
| CO            | Cardiac Output                  | two distances and two areas (by ellipse, trace or circle) and one 2 beat time interval | $CO[1/\text{min}] = SV \times HR / 1000$                               |
| SI            | Stroke Volume Index             | body surface area and SV                                                               | $SI[\text{ml}/\text{m}^2] = SV / BSA$                                  |
| CI            | Cardiac Index                   | body surface area and CO                                                               | $CI[\text{ml}/\text{min}/\text{m}^2] = CO / BSA$                       |

(\*1) where  $a_{di}$  is the  $i^{\text{th}}$  disc diameter of LVAd.

(\*2) where  $a_{si}$  is the  $i^{\text{th}}$  disc diameter of LVAs.

Reference: Journal of American Society of Echocardiography, Vol. 2 Nov 5, 1989.

Table 10–11. LV SP-DISC Method Formulas (cont'd)

## LV Calculation Formulas (LV BP-DISC Method)




Figure 10-46. LV BP-DISC Method Measurements

| Calc Mnemonic | Calc Name                                           | Input Measurements                                                   | Formula                                            |
|---------------|-----------------------------------------------------|----------------------------------------------------------------------|----------------------------------------------------|
| LVA2d         | Left Ventricular Area, Apical 2 Chamber, Diastole   | one area (by trace only) is divided into 20 elliptical disk segments | $LVA2d = a1[cm^2]$                                 |
| LVA2s         | Left Ventricular Area, Apical 2 Chamber, Systole    | one area (by trace only) is divided into 20 elliptical disk segments | $LVA2s = a1[cm^2]$                                 |
| LVA4d         | Left Ventricular Area, Apical 4 Chamber, Diastole   | one area (by trace only) is divided into 20 elliptical disk segments | $LVA4d = a1[cm^2]$                                 |
| LVA4s         | Left Ventricular Area, Apical 4 Chamber, Systole    | one area (by trace only) is divided into 20 elliptical disk segments | $LVA4s = a1[cm^2]$                                 |
| LVL2d         | Left Ventricular Length, Apical 2 Chamber, Diastole | one distance                                                         | $LVL2d = d1[cm \text{ or } mm]$                    |
| LVL2s         | Left Ventricular Length, Apical 2 Chamber, Systole  | one distance                                                         | $LVL2s = d1[cm \text{ or } mm]$                    |
| LVL4d         | Left Ventricular Length, Apical 4 Chamber, Diastole | one distance                                                         | $LVL4d = d1[cm \text{ or } mm]$                    |
| LVL4s         | Left Ventricular Length, Apical 4 Chamber, Systole  | one distance                                                         | $LVL4s = d1[cm \text{ or } mm]$                    |
| HR            | Heart Rate (beats/minute)                           | one 2 beat time interval                                             | $HR[BPM] = 120 [sec] / 2 \text{ beat time } [sec]$ |

Table 10-12. LV BP-DISC Method Formulas



## Advanced Cardiac Calculations

### LV Calculation Formulas (LV BP-DISC Method) (cont'd)

| Calc Mnemonic | Calc Name           | Input Measurements                                                                     | Formula                                                                                            |
|---------------|---------------------|----------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|
| EdV           | End Diastole Volume | one distance and two areas (by ellipse, trace or circle)                               | $\text{EdV[ml]} = \frac{(\pi/4) * ((\text{LVLd}/20) \sum_{i=1}^{20} a_{di} b_{di})}{1} \quad (*1)$ |
| EsV           | End Systole Volume  | one distance and two areas (by ellipse, trace or circle)                               | $\text{EsV[ml]} = \frac{(\pi/4) * ((\text{LVLS}/20) \sum_{i=1}^{20} a_{si} b_{si})}{1} \quad (*2)$ |
| SV            | Stroke Volume       | two distances and two areas (by ellipse, trace or circle)                              | $\text{SV[ml]} = \text{EdV} - \text{EsV}$                                                          |
| EF            | Ejection Fraction   | two distances and two areas (by ellipse, trace or circle)                              | $\text{EF} = \text{SV} / \text{EdV} \times 100$                                                    |
| CO            | Cardiac Output      | two distances and two areas (by ellipse, trace or circle) and one 2 beat time interval | $\text{CO[1/min]} = \text{SV} \times \text{HR} / 1000$                                             |
| SI            | Stroke Volume Index | body surface area and SV                                                               | $\text{SI[ml/m}^2] = \text{SV} / \text{BSA}$                                                       |
| CI            | Cardiac Index       | body surface area and CO                                                               | $\text{CI[ml/min/m}^2] = \text{CO} / \text{BSA}$                                                   |

(\*1) where  $a_{di}$  is the  $i^{\text{th}}$  disc diameter of LVA2d,  $b_{di}$  is the  $i^{\text{th}}$  disc diameter of LVA4d, LVLd is the max (LVL2d, LVL4d).

(\*2) where  $a_{si}$  is the  $i^{\text{th}}$  disc diameter of LVA2s,  $b_{si}$  is the  $i^{\text{th}}$  disc diameter of LVA4s, LVLS is the max (LVL2s, LVL4s).

Reference: Journal of American Society of Echocardiography, Vol. 2 Nov 5, 1989.

Table 10–12. LV BP-DISC Method Formulas (cont'd)

## LV Calculation Formulas (Modified Simpson's Rule Method)




Figure 10–47. Modified Simpson's Rule Method Measurements

Refer to 10–13 for the following Modified Simpson's Rule Method Formulas:

LVLd, LVLs, LVAMd, LVAMs, LVAPd, LVAPs, HR, EdV, EsV, SV, EF and CO.

| Calc Mnemonic | Calc Name           | Input Measurements       | Formula                 |
|---------------|---------------------|--------------------------|-------------------------|
| SI            | Stroke Volume Index | body surface area and SV | $SI[ml/m^2]=SV/BSA$     |
| CI            | Cardiac Index       | body surface area and CO | $CI[ml/min/m^2]=CO/BSA$ |

Table 10–13. Modified Simpson's Rule Method Calculation Formulas

## LV Calculation Formulas (Single Plane Ellipsoid Method)




Figure 10–48. Single Ellipsoid Method Measurements

Refer to 10–15 for the following Single Plane Ellipsoid Method Formulas:

LVLd, LVLs, LVAd, LVAs, HR, EdV, EsV, SV, EF and CO.

| Calc Mnemonic | Calc Name           | Input Measurements       | Formula                 |
|---------------|---------------------|--------------------------|-------------------------|
| SI            | Stroke Volume Index | body surface area and SV | $SI[ml/m^2]=SV/BSA$     |
| CI            | Cardiac Index       | body surface area and CO | $CI[ml/min/m^2]=CO/BSA$ |

Table 10–14. Single Plane Ellipsoid Method Calculation Formulas

### LV Calculation Formulas (Bi Plane Ellipsoid Method)




Figure 10–49. Bi Plane Ellipsoid Method Measurements

Refer to 10–17 for the following Bi Plane Ellipsoid Method Formulas:

LVMLd, LVMLs, LVAMd, LVAMs, LVAd, LVAs, HR, EdV, EsV, SV, EF and CO.

| Calc Mnemonic | Calc Name           | Input Measurements       | Formula                                                      |
|---------------|---------------------|--------------------------|--------------------------------------------------------------|
| SI            | Stroke Volume Index | body surface area and SV | $SI[\text{ml}/\text{m}^2] = \text{SV}/\text{BSA}$            |
| CI            | Cardiac Index       | body surface area and CO | $CI[\text{ml}/\text{min}/\text{m}^2] = \text{CO}/\text{BSA}$ |

Table 10–15. Bi Plane Ellipsoid Method Calculation Formulas

## LV Calculation Formulas (Gibson Method)




Figure 10-50. Gibson Method Measurements

| Calc Mnemonic | Calc Name                                         | Input Measurements       | Formula                                                                 |
|---------------|---------------------------------------------------|--------------------------|-------------------------------------------------------------------------|
| LVIDd         | Left Ventricular Internal Diameter, Diastole      | one distance             | $LVIDd = d1[\text{cm or mm}]$                                           |
| LVIDs         | Left Ventricular Internal Diameter, Systole       | one distance             | $LVIDs = d1[\text{cm or mm}]$                                           |
| IVSd          | Interventricular Septal Thickness, Diastole       | one distance             | $IVSd = d1[\text{cm or mm}]$                                            |
| IVSs          | Interventricular Septal Thickness, Systole        | one distance             | $IVSs = d1[\text{cm or mm}]$                                            |
| LVPWd         | Left Ventricle Posterior Wall Thickness, Diastole | one distance             | $LVPWd = d1[\text{cm or mm}]$                                           |
| LVPWs         | Left Ventricle Posterior Wall Thickness, Systole  | one distance             | $LVPWs = d1[\text{cm or mm}]$                                           |
| HR            | Heart Rate (beats/minute)                         | one 2 beat time interval | $HR[\text{BPM}] = 120 [\text{sec}] / 2 \text{ beat time } [\text{sec}]$ |
| ET            | Ejection Time                                     | one time interval        | $ET = t1[\text{ms or sec}]$                                             |

Table 10-16. Gibson Method Measurement Formulas



## Advanced Cardiac Calculations

---

### LV Calculation Formulas (Gibson Method) (cont'd)

| Calc Mnemonic | Calc Name                   | Input Measurements                         | Formula                                                              |
|---------------|-----------------------------|--------------------------------------------|----------------------------------------------------------------------|
| EdV           | End Diastole Volume         | three distances                            | $EdV[ml] = \pi/6 \times (LVIDd)^2 \times (0.98 \times LVIDd + 0.59)$ |
| EsV           | End Systole Volume          | three distances                            | $EsV[ml] = \pi/6 \times (LVIDs)^2 \times (1.14 \times LVIDs + 4.18)$ |
| FS            | Fractional Shortening       | two distances                              | $FS = (1 - LVIDs/LVIDd) \times 100$                                  |
| SV            | Stroke Volume               | two distances                              | $SV[ml] = EdV - EsV$                                                 |
| EF            | Ejection Fraction           | two distances                              | $EF = SV/EdV \times 100$                                             |
| CO            | Cardiac Output              | two distances and one 2 beat time interval | $CO[1/min] = SV \times HR/1000$                                      |
| SI            | Stroke Volume Index         | body surface area and SV                   | $SI[ml/m^2] = SV/BSA$                                                |
| CI            | Cardiac Index               | body surface area and CO                   | $CI[ml/min/m^2] = CO/BSA$                                            |
| MVCF          | Mean Vcf                    | two distances & one time interval          | $MVCF[circ/s] = (LVIDd - LVIDs) / (LVIDd \times ET)$                 |
| LVM           | Left Ventricle Cardiac Mass | three distances                            | $LVM[g] = 1.04 \times [(IVSd + LVPWd + LVIDd)^3 - LVIDd^3] - 13.6$   |

Table 10–17. Gibson Method Calculation Formulas

## B-Mode Analysis – Parasternal Long Axis



Figure 10–51. PLAX Measurements

| Calc Mnemonic | Calc Name                                             | Input Measurements       | Formula                                                            |
|---------------|-------------------------------------------------------|--------------------------|--------------------------------------------------------------------|
| AOd           | Aortic Root Dimension, Diastole                       | one distance             | $AOd=d1[\text{cm or mm}]$                                          |
| IVSd          | Interventricular Septal Dimension at Diastole         | one distance             | $IVSd=d1[\text{cm or mm}]$                                         |
| LVPWd         | Left Ventricular Posterior Wall Thickness at Diastole | one distance             | $LVPWd=d1[\text{cm or mm}]$                                        |
| LVIDd         | Left Ventricular Interior Dimension at Diastole       | one distance             | $LVIDd=d1[\text{cm or mm}]$                                        |
| AOs           | Aortic Root Dimension at Systole                      | one distance             | $AOs=d1[\text{cm or mm}]$                                          |
| LVIDs         | Left Ventricular Interior Dimension at Systole        | one distance             | $LVIDs=d1[\text{cm or mm}]$                                        |
| ALSS          | Aortic Valve Leaflet Separation, Systole              | one distance             | $ALSS=d1[\text{cm or mm}]$                                         |
| LADs          | Left Atrial Dimension at Systole                      | one distance             | $LADs=d1[\text{cm or mm}]$                                         |
| LVOTs         | LV Outflow Tract Diameter at Systole                  | one distance             | $LVOTs=d1[\text{cm or mm}]$                                        |
| AAs           | Aortic Area at Systole                                | one distance (AOs)       | $AAs[\text{cm}^2]=(\pi/4) \times (AOs)^2$                          |
| LVOTA         | LV Outflow Tract Area at Systole                      | one distance             | $LVOTA[\text{cm}^2]=(\pi/4) \times (LVOTs)^2$                      |
| HR            | Heart Rate (beats/minute)                             | one 2 beat time interval | $HR[\text{BPM}]=120[\text{sec}] / 2 \text{ beat time}[\text{sec}]$ |

Table 1. Parasternal Long Axis (PLAX)





## B-Mode Analysis – Parasternal Short Axis (PSAX-AV)



Figure 10-52. PSAX-AV Measurements

| Calc Mnemonic | Calc Name                                         | Input Measurements       | Formula                                                               |
|---------------|---------------------------------------------------|--------------------------|-----------------------------------------------------------------------|
| AOd           | Aortic Root Dimension, Diastole                   | one distance             | $AOd=d1[\text{cm or mm}]$                                             |
| PADd          | Pulmonary Artery Diameter, Diastole               | one distance             | $PADd=d1[\text{cm or mm}]$                                            |
| RVOTd         | Right Ventricle Outflow Tract Diameter, Diastole  | one distance             | $RVOTd=d1[\text{cm or mm}]$                                           |
| AOs           | Aortic Root Dimension at Systole                  | one distance             | $AOs=d1[\text{cm or mm}]$                                             |
| PADs          | Pulmonary Artery Diameter at Systole              | one distance             | $PADs=d1[\text{cm or mm}]$                                            |
| LADML         | Left Atrium Medial-Lateral Diameter at Systole    | one distance             | $LADML=d1[\text{cm or mm}]$                                           |
| RVOTs         | Right Ventricle Outflow Tract Diameter at Systole | one distance             | $RVOTs=d1 [\text{cm or mm}]$                                          |
| RVOTA         | Right Ventricle Outflow Tract Area at Systole     | one distance             | $RVOTA [\text{cm}^2]= (\pi/4) \times (RVOTs)^2$                       |
| PAAs          | Pulmonary Artery Area at Systole                  | one distance (PADs)      | $PAAs[\text{cm}^2]= (\pi/4) \times (PADs)^2$                          |
| AAs           | Aortic Area at Systole                            | one distance (AOs)       | $AAs[\text{cm}^2]= (\pi/4) \times (AOs)^2$                            |
| HR            | Heart Rate (beats/minute)                         | one 2 beat time interval | $HR[\text{BPM}]=120 [\text{sec}] / 2 \text{ beat time } [\text{sec}]$ |

Table 2. Parasternal Short Axis - Aortic Valve (PSAX-AV)

**B-Mode Analysis – Parasternal Short Axis (PSAX-MV)**




Figure 10-53.PSAX–MV Measurements

| Calc Mnemonic | Calc Name                                     | Input Measurements                     | Formula                                        |
|---------------|-----------------------------------------------|----------------------------------------|------------------------------------------------|
| LVAMd         | Left Ventricular Area, Mitral Valve, Diastole | one area (by ellipse, trace or circle) | $LVAMd=a1[cm^2]$                               |
| LVAMs         | Left Ventricular Area, Mitral Valve, Systole  | one area (by ellipse, trace or circle) | $LVAMs=a1[cm^2]$                               |
| MVOA          | Mitral Valve Maximum Orifice Area at Diastole | one area (by ellipse, trace or circle) | $MVOA=a1[cm^2]$                                |
| HR            | Heart Rate (beats/minute)                     | one 2 beat time interval               | $HR[BPM]=120 [sec]/ 2 \text{ beat time [sec]}$ |

Table 3. Parasternal Short Axis - Mitral Valve (PSAX-MV)



## B-Mode Analysis – Parasternal Short Axis (PSAX-PAP)



Figure 10-54. PSAX-PAP Measurements

| Calc Mnemonic | Calc Name                                            | Input Measurements                     | Formula                                            |
|---------------|------------------------------------------------------|----------------------------------------|----------------------------------------------------|
| LVAPd         | Left Ventricular Area, Papillary Muscles at Diastole | one area (by ellipse, trace or circle) | $LVAPd = a1[cm^2]$                                 |
| LVAPs         | Left Ventricular Area, Papillary Muscles at Systole  | one area (by ellipse, trace or circle) | $LVAPs = a1[cm^2]$                                 |
| HR            | Heart Rate (beats/minute)                            | one 2 beat time interval               | $HR[BPM] = 120 [sec] / 2 \text{ beat time } [sec]$ |

Table 4. Parasternal Short Axis - Papillary Muscles (PSAX-PAP)

## B-Mode Analysis – Apical 4 Chamber (AP-4CH)




Figure 10–55. AP–4CH Measurements (Diastole)

| Calc Mnemonic | Calc Name                                            | Input Measurements                     | Formula                       |
|---------------|------------------------------------------------------|----------------------------------------|-------------------------------|
| LVAd          | Left Ventricle Area, Diastole                        | one area (by ellipse, trace or circle) | $LVAd=a1[cm^2]$               |
| LVLd          | Left Ventricular Length Dimension, Diastole          | one distance                           | $LVLd=d1[cm \text{ or } mm]$  |
| LVMLd         | Left Ventricular Medial-Lateral Dimension, Diastole  | one distance                           | $LVMLd=d1[cm \text{ or } mm]$ |
| RVAd          | Right Ventricle Area, Diastole                       | one area (by ellipse, trace or circle) | $RVAd=a1[cm^2]$               |
| RVLd          | Right Ventricular Length Dimension, Diastole         | one distance                           | $RVLd=d1[cm \text{ or } mm]$  |
| RVMLd         | Right Ventricular Medial-Lateral Dimension, Diastole | one distance                           | $RVMLd=d1[cm \text{ or } mm]$ |

Table 5. Apical 4 Chamber (AP-4CH) Diastole



## B-Mode Analysis – Apical 4 Chamber (AP-4CH) (cont'd)




Figure 10-56. AP-4CH Measurements (Systole)

| Calc Mnemonic | Calc Name                                           | Input Measurements                     | Formula                                         |
|---------------|-----------------------------------------------------|----------------------------------------|-------------------------------------------------|
| LVAs          | Left Ventricle Area, Systole                        | one area (by ellipse, trace or circle) | $LVAs=a1[cm^2]$                                 |
| LVLS          | Left Ventricular Length Dimension, Systole          | one distance                           | $LVLS=d1[cm \text{ or } mm]$                    |
| LVMLS         | Left Ventricular Medial-Lateral Dimension, Systole  | one distance                           | $LVMLS=d1[cm \text{ or } mm]$                   |
| RVAs          | Right Ventricle Area, Systole                       | one area (by ellipse, trace or circle) | $RVAs=a1[cm^2]$                                 |
| RVLS          | Right Ventricular Length Dimension at Systole       | one distance                           | $RVLS=d1[cm \text{ or } mm]$                    |
| RVMLS         | Right Ventricular Medial-Lateral Dimension, Systole | one distance                           | $RVMLS=d1[cm \text{ or } mm]$                   |
| TAML          | Tricuspid Annulus Medial-Lateral Dimension, Systole | one distance                           | $TAML=d1[cm \text{ or } mm]$                    |
| LAA           | Left Atrial Area, Systole                           | one area (by ellipse, trace or circle) | $LAA=a1[cm^2]$                                  |
| LADML         | Left Atrial Medial-Lateral Dimension, Systole       | one distance                           | $LADML=d1[cm \text{ or } mm]$                   |
| LADSI         | Left Atrial Superior-Inferior Dimension, Systole    | one distance                           | $LADSI=d1[cm \text{ or } mm]$                   |
| RAA           | Right Atrial Area, Systole                          | one area (by ellipse, trace or circle) | $RAAs=a1[cm^2]$                                 |
| RADML         | Right Atrial Medial-Lateral Dimension, Systole      | one distance                           | $RADML=d1[cm \text{ or } mm]$                   |
| RADSI         | Right Atrial Superior-Inferior Dimension, Systole   | one distance                           | $RADSI=d1[cm \text{ or } mm]$                   |
| HR            | Heart Rate (beats/minute)                           | one 2 beat time interval               | $HR[BPM]=120 [sec]/ 2 \text{ beat time } [sec]$ |

Table 6. Apical 4 Chamber (AP-4CH) Systole

## B-Mode Analysis – Apical 2 Chamber (AP-2CH)




Figure 10-57. AP-2CH Measurements

| Calc Mnemonic | Calc Name                                           | Input Measurements                     | Formula                                            |
|---------------|-----------------------------------------------------|----------------------------------------|----------------------------------------------------|
| LVAd          | Left Ventricle Area, Diastole                       | one area (by ellipse, trace or circle) | $LVAd = a1[cm^2]$                                  |
| LVLd          | Left Ventricular Length Dimension, Diastole         | one distance                           | $LVLd = d1[cm \text{ or } mm]$                     |
| LVMLd         | Left Ventricular Medial-Lateral Dimension, Diastole | one distance                           | $LVMLd = d1[cm \text{ or } mm]$                    |
| LVAs          | Left Ventricle Area, Systole                        | one area (by ellipse, trace or circle) | $LVAs = a1[cm^2]$                                  |
| LVLs          | Left Ventricular Length Dimension at Systole        | one distance                           | $LVLs = d1[cm \text{ or } mm]$                     |
| LVMLs         | Left Ventricular Medial-Lateral Dimension, Systole  | one distance                           | $LVMLs = d1[cm \text{ or } mm]$                    |
| HR            | Heart Rate (beats/minute)                           | one 2 beat time interval               | $HR[BPM] = 120 [sec] / 2 \text{ beat time } [sec]$ |

Table 7. Apical 2 Chamber (AP-2CH)



## M-Mode Analysis – Left/Right Ventricle (M-LV/RV)




Figure 10-58. M-LV/RV Measurements

| Calc Mnemonic | Calc Name                                                          | Input Measurements                    | Formula                                                                 |
|---------------|--------------------------------------------------------------------|---------------------------------------|-------------------------------------------------------------------------|
| RVDd          | Right Ventricular Dimension, Diastole                              | one distance                          | $RVDd = d1[\text{cm or mm}]$                                            |
| IVSd          | Interventricular Septal Thickness, Diastole                        | one distance                          | $IVSd = d1[\text{cm or mm}]$                                            |
| LVIDd         | Left Ventricular Interior Dimension, Diastole                      | one distance                          | $LVIDd = d1[\text{cm or mm}]$                                           |
| LVPWd         | Left Ventricular Posterior Wall Thickness at Diastole              | one distance                          | $LVPWd = d1[\text{cm or mm}]$                                           |
| IVSs          | Interventricular Septal Thickness, Systole                         | one distance                          | $IVSs = d1[\text{cm or mm}]$                                            |
| LVPWs         | Left Ventricular Posterior Wall Thickness at Systole               | one distance                          | $LVPWs = d1[\text{cm or mm}]$                                           |
| LVIDs         | Left Ventricular Interior Dimension at Systole                     | one distance                          | $LVIDs = d1[\text{cm or mm}]$                                           |
| HR            | Heart Rate (beats/minute)                                          | one 2 beat time interval              | $HR[\text{BPM}] = 120 [\text{sec}] / 2 \text{ beat time } [\text{sec}]$ |
| IVS/PW        | Interventricular Septum/Left Ventricle Posterior Ratio at Diastole | two distances (IVSd & LVPWd)          | $IVS/PW[\%] = IVSd / LVPWd \times 100$                                  |
| FS            | Left Ventricle Internal Dimension Fractional Shortening            | two distances (LVIDd & LVIDs)         | $FS[\%] = [(LVIDd - LVIDs) / LVIDd] \times 100$                         |
| %STIVS        | Interventricular Shortening                                        | two distances (IVSd & IVSs)           | $\%STIVS[\%] = [(IVSs - IVSd) / IVSd] \times 100$                       |
| %STPW         | Left Ventricle Posterior Wall Shortening                           | two distances (LVPWd & LVPWs)         | $\%STPW[\%] = [(LVPWs - LVPWd) / LVPWd] \times 100$                     |
| LVM           | Left Ventricle Cardiac Mass                                        | three distances (LVPWd, LVIDd & IVSd) | $LVM[g] = 1.04 \times [(IVSd + LVPWd + LVIDd)^3 - LVIDd^3] - 13.6$      |

Table 8. Left/Right Ventricle (M-LV/RV)

## M-Mode Analysis – Mitral Valve (M-MV)




Figure 10-59. M-MV Measurements

| Calc Mnemonic | Calc Name                   | Input Measurements        | Formula                                                                 |
|---------------|-----------------------------|---------------------------|-------------------------------------------------------------------------|
| EPSS          | E point Septal Separation   | one distance              | $EPSS = d1[\text{cm or mm}]$                                            |
| V_EF          | Mitral Valve E-F Velocity   | one slope                 | $V\_EF = s1[\text{cm/s}]$                                               |
| V_DE          | Mitral Valve D-E Velocity   | one slope                 | $V\_DE = s1[\text{cm/s}]$                                               |
| D_DE          | Mitral Valve D-E Separation | one distance              | $V\_DE = d1[\text{cm or mm}]$                                           |
| D_CE          | Mitral Valve C-E Separation | one distance              | $D\_CE = d1[\text{cm or mm}]$                                           |
| T_AC          | Mitral Valve A-C Interval   | one time interval         | $T\_AC = t1[\text{msec or sec}]$                                        |
| P-R           | P-R Interval                | one time interval         | $P-R = t1[\text{msec or sec}]$                                          |
| HR            | Heart Rate (beats/minute)   | one 2 beat time interval  | $HR[\text{BPM}] = 120 [\text{sec}] / 2 \text{ beat time } [\text{sec}]$ |
| A/E           | Mitral Valve A/E Ratio      | two distances (A-C & D-E) | $A/E[\%] = (D\_AC / D\_DE) \times 100$                                  |
| PR-AC         | Mitral Valve PR-AC Interval | two time intervals        | $PR-AC[\text{msec or sec}] = P-R - T\_AC$                               |

Table 9. Mitral Valve (M-MV)



*V\_DE measurement calculates D\_DE simultaneously.  
D\_DE measurement calculates V\_DE simultaneously.  
T\_AC measurement calculates D\_AC simultaneously.*

*V\_EF, V\_DE, D\_DE, D\_CE and T\_AC measurement results should be performed through the slope measurement.*





## M-Mode Analysis – Aortic Valve (M-AV)




Figure 10–60. M-AV Measurements

| Calc Mnemonic | Calc Name                                           | Input Measurements              | Formula                                                               |
|---------------|-----------------------------------------------------|---------------------------------|-----------------------------------------------------------------------|
| AOd           | Aortic Root Dimension, Diastole                     | one distance<br>A-B             | $AOd=d1[\text{cm or mm}]$                                             |
| ALSS          | Aortic Valve Leaflet Separation, Systole            | one distance<br>E-F             | $ALSS=d1[\text{cm or mm}]$                                            |
| LADs          | Left Atrium Dimension, Systole                      | one distance<br>C-D             | $LADs=d1[\text{cm or mm}]$                                            |
| RVOTs         | Right Ventricular Outflow Tract Diameter at Systole | one distance<br>G-H             | $RVOTs=d1[\text{cm or mm}]$                                           |
| LVPEP         | Left Ventricle Pre-Ejection Period                  | one time interval<br>Q-I        | $LVPEP=t1[\text{msec or sec}]$                                        |
| LVET          | Left Ventricle Ejection Time                        | one time interval<br>I-K        | $LVET=t1[\text{msec or sec}]$                                         |
| HR            | Heart Rate (beats/minute)                           | one 2 beat time interval        | $HR[\text{BPM}]=120 [\text{sec}] / 2 \text{ beat time } [\text{sec}]$ |
| LA/AO         | LADs/AOd Ratio                                      | two distances<br>(LADs & AOd)   | $LA/AO[\%]=LADs/AOd \times 100$                                       |
| PEP/ET        | Left Ventricle Systole Time Interval Ratio          | two distances<br>(LVPEP & LVET) | $PEP/ET[\%]= (LVPEP/LVET) \times 100$                                 |

Table 10. Aortic Valve (M-AV)

## M-Mode Analysis – Pulmonic Valve (M-PV)




Figure 10–61. M–PV Measurements

| Calc Mnemonic | Calc Name                                   | Input Measurements           | Formula                                                                 |
|---------------|---------------------------------------------|------------------------------|-------------------------------------------------------------------------|
| PADs          | Pulmonic Artery Diameter at Systole         | one distance                 | $PADs = d1[\text{cm or mm}]$                                            |
| RVPEP         | Right Ventricle Pre-Ejection Period         | one time interval            | $RVPEP = t1[\text{msec or sec}]$                                        |
| RVET          | Right Ventricle Ejection Time               | one time interval            | $RVET = t1[\text{msec or sec}]$                                         |
| aWAVE         | Pulmonic Valve a-Wave Amplitude             | one distance                 | $aWAVE = d1[\text{cm or mm}]$                                           |
| HR            | Heart Rate (beats/minute)                   | one 2 beat time interval     | $HR[\text{BPM}] = 120 [\text{sec}] / 2 \text{ beat time } [\text{sec}]$ |
| PAAs          | Pulmonary Artery Area at Systole            | one distance (PADs)          | $PAAs[\text{cm}^2] = (\pi/4) \times (PADs)^2$                           |
| PEP/ET        | Right Ventricle Systole Time Interval Ratio | two distances (RVPEP & RVET) | $PEP/ET[\%] = (RVPEP/RVET) \times 100$                                  |

Table 11. Pulmonic Valve (M-PV)



## M-Mode Analysis – Tricuspid Valve (M-TV)




Figure 10-62. M-TV Measurements

| Calc Mnemonic | Calc Name                      | Input Measurements       | Formula                                                                |
|---------------|--------------------------------|--------------------------|------------------------------------------------------------------------|
| TV_EF         | Tricuspid Valve E-F Velocity   | one slope                | $TV\_EF = v1[\text{cm/s or m/s}]$                                      |
| TV_DE         | Tricuspid Valve D-E Amplitude  | one distance             | $TV\_DE = d1[\text{cm or mm}]$                                         |
| TV_CE         | Tricuspid Valve C-E Amplitude  | one distance             | $TV\_CE = d1[\text{cm or mm}]$                                         |
| TV_AC         | Tricuspid Valve A-C Interval   | one time interval        | $TV\_AC = t1[\text{msec or sec}]$                                      |
| P-R           | P-R Interval                   | one time interval        | $P-R = t1[\text{msec or sec}]$                                         |
| HR            | Heart Rate (beats/minute)      | one 2 beat time interval | $HR[\text{BPM}] = 120 [\text{sec}] / 2 \text{ beat time} [\text{sec}]$ |
| TV-A/E        | Tricuspid Valve A/E Ratio      | one distance (TV<D-E>)   | $TV-A/E[\%] = (D < TV\_AC > / TV\_DE) \times 100$                      |
| TPR-AC        | Tricuspid Valve PR-AC Interval | two time intervals       | $TPR-AC[\text{ms or sec}] = P-R - TV\_AC$                              |

Table 12. Tricuspid Valve (M-TV)



*TV\_AC measurement calculates D\_AC simultaneously.*

*TV\_DE, TV\_EF, TV\_AC and TV\_CE measurement results should be performed through the slope measurement.*




Figure 10–63. D-MV Measurements

| Calc Mnemonic | Calc Name                              | Input Measurements                     | Formula                                            |
|---------------|----------------------------------------|----------------------------------------|----------------------------------------------------|
| MVOA          | Mitral Valve Orifice Area at Diastole  | one area (by ellipse, trace or circle) | $MVOA = a[cm]^2$                                   |
| PFVMV         | Peak Flow Velocity at the Mitral Valve | one velocity                           | $PFV - MV = v[cm/s \text{ or } m/s]$               |
| PFVLA         | Peak Flow Velocity at Left Atrium      | one velocity                           | $PFV - LA = v[cm/s \text{ or } m/s]$               |
| ET            | Ejection Time                          | one time interval                      | $ET = t[ms \text{ or } sec]$                       |
| FVI           | Flow Velocity Integral                 | flow velocities                        | $FVI[mm \text{ or } cm] = \int V dt$               |
| PHT           | Pressure Half Time                     | one time interval                      | $PHT = t[ms \text{ or } sec]$                      |
| HR            | Heart Rate (beats/minute)              | one 2 beat time interval               | $HR[BPM] = 120 [sec] / 2 \text{ beat time } [sec]$ |

Table 13. Mitral Valve (D-MV)



## Doppler Analysis – Mitral Valve (D-MV) (cont'd)

| Calc Mnemonic | Calc Name                   | Input Measurements                                                                  | Formula                                                                                  |
|---------------|-----------------------------|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|
| AccT          | Flow Acceleration Time      | one time interval                                                                   | $\text{AccT}[\text{msec or sec}] = t_{\text{PFV-MV}} - t_{\text{ET}}$                    |
| DecT          | Flow Deceleration Time      | one time interval                                                                   | $\text{DecT}[\text{msec or sec}] = t_{\text{ET}} - t_{\text{PFV-MV}}$                    |
| PkFVI         | Peak Flow Velocity Intergal | one velocity,<br>one time interval                                                  | $\text{PkFVI}[\text{cm/s or m/s}] = 1.14 \times \text{PFV-MV} \times \text{ET}/20 + 0.3$ |
| MeanFV        | Mean Flow Velocity          | flow velocities, one area<br>(by ellipse, trace or circle)                          | $\text{MeanFV}[\text{cm/s or m/s}] = \text{FVI}/t_{\text{FVI}}$                          |
| SV            | Stroke Volume               | flow velocities and one area<br>(by ellipse, trace or circle)                       | $\text{SV}[\text{ml}] = \text{FVI} \times \text{MVOA}$                                   |
| SI            | Stroke Volume Index         | body surface area and SV                                                            | $\text{SI}[\text{ml}/\text{m}^2] = \text{SV}/\text{BSA}$                                 |
| CO            | Cardiac Output              | flow velocities, one area<br>(by ellipse, trace or circle)<br>and one time interval | $\text{CO}[\text{l/min}] = \text{SV} \times \text{HR}/1000$                              |
| CI            | Cardiac Index               | body surface area and CO                                                            | $\text{CI}[\text{ml/min}/\text{m}^2] = \text{CO}/\text{BSA}$                             |
| MaxPG         | Maximum Pressure Gradient   | two velocities                                                                      | $\text{MaxPG}[\text{mmHg}] = 4 \times (\text{PFV-MV}^2)$                                 |
| MeanPG        | Mean Pressure Gradient      | flow velocities                                                                     | $\text{MeanPG}[\text{mmHg}] = \frac{n}{4 \times \sum_{i=1}^n (V^2/n)}$                   |
| MVA           | Mitral Valve Area           | one time interval                                                                   | $\text{MVA}[\text{cm}^2] = 220/\text{PHT}$                                               |

Table 13. Mitral Valve (D-MV) (cont'd)

## Doppler Analysis – Aortic Valve (D-AV)




Figure 10-64. D-AV Measurements

| Calc Mnemonic | Calc Name                                                | Input Measurements       | Formula                                                              |
|---------------|----------------------------------------------------------|--------------------------|----------------------------------------------------------------------|
| AOS           | Aortic Root Dimension, Systole                           | one distance             | $AOS=d1[\text{cm or mm}]$                                            |
| FVI-AV        | Flow Velocity Intergal at Aortic Valve                   | flow velocities          | $FVI-AV[\text{mm or cm}]=\int Vdt$                                   |
| FVI-LV        | Flow Velocity Intergal at Left Ventricular Outflow Tract | flow velocities          | $FVI-LV[\text{mm or cm}]=\int Vdt$                                   |
| PFVAV         | Peak Flow Velocity at Aortic Valve                       | one velocity             | $PFVAV=v1[\text{cm/s or m/s}]$                                       |
| PFVLV         | Peak Flow Velocity at Left Ventricular Outflow Tract     | one velocity             | $PFVLV=v1[\text{cm/s or m/s}]$                                       |
| ET            | Ejection Time                                            | one time interval S-E    | $ET=t1[\text{ms or sec}]$                                            |
| PHT           | Pressure Half Time                                       | one time interval        | $PHT=t1[\text{ms or sec}]$                                           |
| HR            | Heart Rate (beats/minute)                                | one 2 beat time interval | $HR[\text{BPM}]=120 [\text{sec}]/ 2 \text{ beat time } [\text{sec}]$ |

Table 14. Aortic Valve (D-AV)

## Doppler Analysis – Aortic Valve (D-AV) (cont'd)

| Calc Mnemonic | Calc Name                                      | Input Measurements                                                                  | Formula                                                                                                         |
|---------------|------------------------------------------------|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|
| AV-A          | Aortic Valve Area by cont'Eq                   | one distance<br>two velocities                                                      | $AV-A[mm \text{ or } cm] = LVOT-A \times PFV-LV/PFV-AV$                                                         |
| LVOTA         | Left Ventricular Outflow Tract Area at Systole | LVOTs of PLAX                                                                       | $LVOTA[cm^2] = \pi \times (LVOTs)^2 / 4$                                                                        |
| AccT          | Flow Acceleration Time                         | two time intervals                                                                  | $AccT[msec \text{ or } sec] = t<PFV-AV> - t<ET>$                                                                |
| DecT          | Flow Deceleration Time                         | two time intervals                                                                  | $DecT[msec \text{ or } sec] = t<ET> - t<PFV-AV>$                                                                |
| PkFVI         | Peak Flow Velocity Integral                    | one velocity,<br>one time interval                                                  | $PkFVI[mm \text{ or } cm] = 1.14 \times PFV-AV \times ET / 20 \pm .03$                                          |
| MeanFV        | Mean Flow Velocity                             | flow velocities,<br>one time interval                                               | $MeanFV[cm/s \text{ or } m/s] = FVI / t<FVI>$                                                                   |
| QP:QS         | QP:QS Ratio                                    | two CO                                                                              | $QP:QS = CO(PV) / CO(AV)$                                                                                       |
| MaxPG         | Maximum Pressure Gradient                      | two velocities                                                                      | $MaxPG[mmHg] = 4 \times (PFV-AV^2)$                                                                             |
| MeanPG        | Mean Pressure Gradient                         | flow velocities                                                                     | $MeanPG[mmHg] = \frac{1}{4} \times \sum_{i=1}^n (V_i^2)$                                                        |
| SV            | Stroke Volume                                  | flow velocities, one area<br>(by ellipse or trace or circle)                        | AV-A is calculated.<br>$SV[ml] = FVI-AV \times AV-A$<br>else<br>$SV[ml] = FVI-AV \times (\pi/4) \times (AOs)^2$ |
| SI            | Stroke Volume Index                            | body surface area and SV                                                            | $SI[ml/m^2] = SV/BSA$                                                                                           |
| CO            | Cardiac Output                                 | flow velocities, one area<br>(by ellipse, trace or circle)<br>and one time interval | $CO[ml] = SV \times HR / 1000$                                                                                  |
| CI            | Cardiac Index                                  | body surface area and CO                                                            | $CI[ml/min/m^2] = CO/BSA$                                                                                       |

Table 14. Aortic Valve (D-AV) (cont'd)

## Doppler Analysis – Pulmonic Valve (D-PV)




Figure 10-65. D-PV Measurements

| Calc Mnemonic | Calc Name                                               | Input Measurements             | Formula                                                              |
|---------------|---------------------------------------------------------|--------------------------------|----------------------------------------------------------------------|
| PADs          | Pulmonary Artery Diameter at Systole                    | one distance<br>B-Mode         | $PADs=d1[\text{cm or mm}]$                                           |
| FVI-PV        | Flow Velocity Integral at Pulmonic Valve                | one FVI                        | $FVI-PV=\int Vdt$                                                    |
| FVI-PA        | Flow Velocity Integral at Pulmonic Aorta                | one FVI                        | $FVI-PA=\int Vdt$                                                    |
| FVI-RV        | Flow Velocity Integral at Right Ventricle Outflow Tract | one FVI                        | $FVI-RV=\int Vdt$                                                    |
| PFVPA         | Peak Flow Velocity at Pulmonic Aorta                    | one velocity                   | $PFVPA=v1[\text{cm/s or m/s}]$                                       |
| PFVRV         | Peak Flow Velocity at Right Ventricle Outflow Tract     | one velocity                   | $PFVRV=v1[\text{cm/s or m/s}]$                                       |
| ET            | Ejection Time                                           | one time interval S–E          | $ET=t1[\text{ms or sec}]$                                            |
| PHT           | Pressure Half Time                                      | one time interval              | $PHT=t1[\text{ms or sec}]$                                           |
| HR            | Heart Rate (beats/minute)                               | one 2 beat time interval       | $HR[\text{BPM}]=120 [\text{sec}]/ 2 \text{ beat time } [\text{sec}]$ |
| PV-A          | Pulmonic Valve Area by cont'Eq                          | one distance<br>two velocities | $PV-A[\text{cm}^2]=(\pi/4) \times (PADs)^2 \times PFV-PA/PFV-RV$     |
| AccT          | Flow Acceleration Time                                  | two time intervals             | $AccT[\text{msec or sec}]=t<PFV-PA>-ts<ET>$                          |
| DecT          | Flow Deceleration Time                                  | two time intervals             | $DecT[\text{msec or sec}]=te<ET>-t<PFV-PA>$                          |

Table 15. Pulmonic Valve (D-PV)





## Doppler Analysis – Pulmonic Valve (D-PV) (cont'd)

| Calc Mnemonic | Calc Name                   | Input Measurements                                                            | Formula                                                        |
|---------------|-----------------------------|-------------------------------------------------------------------------------|----------------------------------------------------------------|
| PkFVI         | Peak Flow Velocity Integral | one velocity and one time interval                                            | $PkFVI[\text{mm or cm}] = 1.14 \times PFV-PA \times ET/20+0.3$ |
| MeanFV        | Mean Flow Velocity          | flow velocities and one time interval                                         | $MeanFV[\text{mm or cm}] = FVI/t<FVI>$                         |
| QP:QS         | QP:QS Ratio                 | two CO                                                                        | $QP:QS = CO(PV)/CO(AV)$                                        |
| CO            | Cardiac Output              | flow velocities, one area (by ellipse, trace or circle) and one time interval | $CO[\text{ml}] = SV \times HR/1000$                            |
| CI            | Cardiac Index               | body surface area and CO                                                      | $CI[\text{ml/min/m}^2] = CO/BSA$                               |
| MaxPG         | Maximum Pressure Gradient   | two velocities                                                                | $MaxPG[\text{mmHg}] = 4 \times (PFV-PA^2)$                     |
| MeanPG        | Mean Pressure Gradient      | flow velocities                                                               | $MeanPG[\text{mmHg}] = 4 \times \sum_{i=1}^n (V^2/n)$          |
| SV            | Stroke Volume               | flow velocities<br>RVOTA of PLAX–AV                                           | $SV[\text{ml}] = FVI-PV \times RVOTA$                          |
| SI            | Stroke Volume Index         | body surface area and SV                                                      | $SI[\text{ml/m}^2] = SV/BSA$                                   |

Table 15. Pulmonic Valve (D-PV) (cont'd)

## Advanced Cardiac Calculations

### Doppler Analysis – Tricuspid Valve (D-TV)




Figure 10-66. D-TV Measurements

| Calc Mnemonic | Calc Name                                           | Input Measurements       | Formula                                                              |
|---------------|-----------------------------------------------------|--------------------------|----------------------------------------------------------------------|
| TAML          | Tricuspid Annulus Medial-Lateral Dimension, Systole | one distance             | $TAML=d1[\text{cm or mm}]$                                           |
| FVI           | Flow Velocity Integral                              | flow velocities          | $FVI[\text{mm or cm}]=\int Vdt$                                      |
| PFVTV         | Peak Flow Velocity at Tricuspid Valve               | one velocity             | $PFVTV=v1[\text{cm/s or m/s}]$                                       |
| PFVRA         | Peak Flow Velocity at Right Atrial                  | one velocity             | $PFVRA=v1[\text{cm/s or m/s}]$                                       |
| ET            | Ejection Time                                       | one time interval        | $ET=t1[\text{ms or sec}]$                                            |
| PHT           | Pressure Half Time                                  | one time interval        | $PHT=t1[\text{ms or sec}]$                                           |
| HR            | Heart Rate (beats/minute)                           | one 2 beat time interval | $HR[\text{BPM}]=120 [\text{sec}]/ 2 \text{ beat time } [\text{sec}]$ |
| AccT          | Flow Acceleration Time                              | two time intervals       | $AccT[\text{msec or sec}]=t<PFVTV>-ts<ET>$                           |
| DecT          | Flow Deceleration Time                              | two time intervals       | $DecT[\text{msec or sec}]=te<ET>-t<PFVTV>$                           |

Table 16. Tricuspid Valve (D-TV)



## Doppler Analysis – Tricuspid Valve (D-TV) (cont'd)

| Calc Mnemonic | Calc Name                   | Input Measurements                                  | Formula                                                                 |
|---------------|-----------------------------|-----------------------------------------------------|-------------------------------------------------------------------------|
| PkFVI         | Peak Flow Velocity Integral | one velocity, one time interval                     | $PkFVI[cm/s \text{ or } m/s] = 1.14 \times PFV_{TV} \times ET/20 + 0.3$ |
| MeanFV        | Mean Flow Velocity          | one velocity, one time interval                     | $MeanFV[cm/s \text{ or } m/s] = FVI / t < FVI >$                        |
| MaxPG         | Maximum Pressure Gradient   | two velocities                                      | $MaxPG[mmHg] = 4 \times (PFV - TV)^2$                                   |
| MeanPG        | Mean Pressure Gradient      | flow velocities                                     | $MeanPG[mmHg] = \frac{n}{4 \times \sum_{i=1}^n (V^2/n)}$                |
| SV            | Stroke Volume               | flow velocities, one distance                       | $SV[m] = (\pi/4) \times FVI \times (TAML)^2$                            |
| SI            | Stroke Volume Index         | body surface area and SV                            | $SI[m^3/m^2] = SV/BSA$                                                  |
| CO            | Cardiac Output              | flow velocities, one distance and one time interval | $CO[1/min] = SV \times HR/1000$                                         |
| CI            | Cardiac Index               | body surface area and CO                            | $CI[l/min/m^2] = CO/BSA$                                                |

Table 17. Tricuspid Valve (D-TV) (cont'd)

## Advanced Cardiac Reports

### Overview

The averaging function for AMCAL measurements can be set to a maximum of five measurements. The "Average Number" preset, located on Preset Program page 3, sets the default for the AMCAL measurements.

Two types of report formats are supported by the LOGIQ™ 500 Advanced Cardiac Calculation option.

The first format is the standard display that shows a separate page for each of the first layer menu choices (i.e. Cubed Method, Bullet Method, PSAX-AV, D-MV, etc.). Along with the patient information are the measured values, their average and the calculated values.

Two examples are shown in Figure 10–67 and Figure 10–68.

|               |         |          |        |          |            |
|---------------|---------|----------|--------|----------|------------|
| HOSPITAL NAME |         | Oper ID: |        | 05/14/98 |            |
| ID:           | NAME:   | BSA:     |        | AGE:     |            |
| SEX:          | HEIGHT: | WEIGHT:  |        |          |            |
| Ref MD:       | NOTE:   |          |        |          |            |
| [LV           | ][Cubed | ]        | ←??/?→ |          |            |
| MEASUREMENTS  |         | AVE      | 1      | 2        | 3          |
| RVIDd         |         |          | (      |          | )          |
| IVSd          |         |          | (      |          | )          |
| LVIDd         |         |          | (      |          | )          |
| LVPWd         |         |          | (      |          | )          |
| IVSs          |         |          | (      |          | )          |
| LVIDs         |         |          | (      |          | )          |
| LVPWs         |         |          | (      |          | )          |
| ET            |         |          | (      |          | )          |
| HR            |         |          | (      |          | )          |
| CALCULATIONS  |         |          |        |          |            |
| Edv           |         |          | SI     |          |            |
| EsV           |         |          | CI     |          |            |
| FS            |         |          | MVCF   |          |            |
| SV            |         |          | LVM    |          |            |
| EF            |         |          | IVS/PW |          |            |
| CO            |         |          | %STIVS |          |            |
|               |         |          | %STPW  |          |            |
| COMMENTS:     |         | [LIST    | ][EDIT | ][S-RP   | ][AUTO-P ] |

Figure 10–67. Standard Cubed Method Report Page

<AVG> or <LAST> indicates that the Average or Latest value will be displayed in this column. This depends on the preset, *Average Activity*, in the Setup/Preset Program menu page 3.

## Overview (cont'd)

|                                         |                      |          |        |          |   |
|-----------------------------------------|----------------------|----------|--------|----------|---|
| HOSPITAL NAME                           |                      | Oper ID: |        | 05/09/98 |   |
| ID:                                     | NAME:                | AGE:     |        |          |   |
| SEX:                                    | HEIGHT:              | WEIGHT:  | BSA:   |          |   |
| Ref MD:                                 | NOTE:                |          | ~??/?~ |          |   |
| [B Mode                                 | ] [PSAX Aortic Valve |          | ]      |          |   |
| MEASUREMENTS                            |                      | AVE      | 1      | 2        | 3 |
| PADd                                    |                      |          | {      |          | } |
| RVOTd                                   |                      |          | {      |          | } |
| AAd                                     |                      |          | {      |          | } |
| PADs                                    |                      |          | {      |          | } |
| LADML                                   |                      |          | {      |          | } |
| AAs                                     |                      |          | {      |          | } |
| RVOTs                                   |                      |          | {      |          | } |
| HR                                      |                      |          | {      |          | } |
| CALCULATIONS                            |                      |          |        |          |   |
| PARs                                    |                      |          |        |          |   |
| AAs                                     |                      |          |        |          |   |
| RVOTA                                   |                      |          |        |          |   |
| COMMENTS: [LIST] [EDIT] [S-RP] [AUTO-P] |                      |          |        |          |   |

Figure 10–68. Standard PSAX-AV Report Page

The second report format is the list of measured and calculated values. When a measurement is made, it is entered into each report page that displays that measurement. If the standard report format is used, many unnecessary report pages could be generated. The List type of reporting can be customized to the user's needs, eliminating unnecessary printing requirements.

A List type report page is shown in Figure 10–69.

|                                          |                     |                   |                          |              |  |
|------------------------------------------|---------------------|-------------------|--------------------------|--------------|--|
| Hospital Name                            |                     | Tech. ID: 1 2 3 4 |                          | 10 / 18 / 93 |  |
| ID: Patient ID                           | NAME: Patient Name  | AGE: ###          |                          |              |  |
| SEX: 1                                   | HEIGHT: 123cm       | WEIGHT: 123.45kg  | BSA: 12.34m <sup>2</sup> |              |  |
| Ref MD:                                  | NOTE:               |                   |                          |              |  |
| [ M-Mode                                 | ] [ Aortic Valve    |                   | ] <-- ## / ## -->        |              |  |
| AAd                                      | ###.###mm           | (###.##           | ###.##                   | ###.##       |  |
| ALSS                                     | ###.###mm           | (###.##           | ###.##                   | ###.##       |  |
| LADs                                     | ###.###mm           | (###.##           | ###.##                   | ###.##       |  |
| RVOTs                                    | ###.###mm           | (###.##           | ###.##                   | ###.##       |  |
| LA/AO                                    | ###.###%            | (###.##           | ###.##                   | ###.##       |  |
| [ D-Mode                                 | ] [ Tricuspid Valve |                   | ]                        |              |  |
| TAML                                     | ###.###mm           | (###.##           | ###.##                   | ###.##       |  |
| FVI                                      | ###.###cm           | (###.##           | ###.##                   | ###.##       |  |
| ET                                       | ###.###ms           | (###.##           | ###.##                   | ###.##       |  |
| PHT                                      | ###.###ms           | (###.##           | ###.##                   | ###.##       |  |
| MeanV                                    | ###.###cm/s         | (###.##           | ###.##                   | ###.##       |  |
| MeanPG                                   | ###.###mmHg         | (###.##           | ###.##                   | ###.##       |  |
| SV                                       | ###.###ml           | (###.##           | ###.##                   | ###.##       |  |
| CO                                       | ###.###l/min        | (###.##           | ###.##                   | ###.##       |  |
| COMMENTS: [ RETURN ] [ S-RP ] [ AUTO-P ] |                     |                   |                          |              |  |
| Comments Area                            |                     |                   |                          |              |  |
| Operator Message Area                    |                     |                   |                          |              |  |

Figure 10–69. List Type Report Page



---

## Advanced Cardiac Calculations

---

### Standard Report

The fifth line of the standard report page layout has three important areas of information.

|             |                                                       |
|-------------|-------------------------------------------------------|
| Main Menu   | Shows the main or first layer category (D-Mode)       |
| Sub-Menu    | Shows the sub or second layer category (Aortic Valve) |
| Page Number | Shows present page number and total pages             |

### Pop up Menus

The Main and Sub-Menus are pop-up menus that provide the ability to quickly display reports.

To display a pop up menu, use the **Trackball** to move the highlight cursor to the main or Sub-Menu area. Press **Set**. The pop up menu is displayed.

While the pop up menu is displayed, use the up/down **Trackball** motion or the **Up/Down** arrow keys to select the desired menu item.

Use the **Set** key to activate the menu choice.

### Page Numbers

The page number area not only shows the Present Page/Total Pages, but this area is used to navigate through the display of available pages.

Use the **Trackball** to move the arrow cursor to highlight the navigation arrows on either side of the page numbers.

Press **Set** to display the previous or next page.



## Advanced Cardiac Calculations

---

### Report Page Commands

Four commands at the bottom of the report page allow the user to:

|        |                                                                                           |
|--------|-------------------------------------------------------------------------------------------|
| List   | Change the report format to the LIST type.                                                |
| Edit   | Allows for deleting measurement data on the report page                                   |
| S–RP   | Used to select the titles and order of report page printing.                              |
| Auto–P | Automatically prints out all reports, selected by the S–RP function, in sequential order. |

Use the **Trackball** to highlight the desired command and press **Set**.

### List Type Reports

The List type report displays only those measurements or calculations actually performed. The section order of the list is the same that was programmed by the S–RP command. The item order is the same as the standard report.

All commands are the same as the standard report. The difference is the Return command. Selecting **Return** changes the display to the standard format.



### Hints

When the amount of section information will not fit on one page, a new page is generated so that all will fit on the same page.

### Select Reports (S–RP)

The customization of the cardiac calculation reporting is accomplished with this function.

Select S–RP from any standard report display. The screen for selecting the report order is displayed as shown in Figure 10–70.

Select Reports (S–RP)

|                                   |                           |          |                            |          |
|-----------------------------------|---------------------------|----------|----------------------------|----------|
| HOSPITAL NAME                     |                           | Oper ID: |                            | 05/09/98 |
| ID:                               | NAME:                     |          |                            | AGE:     |
| SEX:                              | HEIGHT:                   | WEIGHT:  | BSA:                       |          |
| Ref MD:                           | NOTE:                     |          |                            |          |
| B Mode Analysis                   |                           |          |                            |          |
| 1                                 | Parasternal Long Axis     |          | D Mode Analysis            |          |
| 2                                 | PSAX Aortic Valve         |          | 12 Mitral Valve            |          |
| 3                                 | PSAX Mitral Valve         |          | 13 Aortic Valve            |          |
| 4                                 | PSAX Papillary Muscles    |          | 14 Pulmonic Valve          |          |
| 5                                 | Apical 4 Chamber Diastole |          | 15 Tricuspid Valve         |          |
| 6                                 | Apical 4 Chamber Systole  |          |                            |          |
| 7                                 | Apical 2 Chamber          |          | LV Analysis                |          |
| M Mode Analysis                   |                           |          |                            |          |
| 8                                 | Left/Right Ventricle      |          | 19 Single Plane Ellipsoid  |          |
| 9                                 | Mitral Valve              |          | 20 Biplane Ellipsoid       |          |
| 10                                | Aortic Valve              |          | 21 Modified Simpson's Rule |          |
| 11                                | Pulmonic Valve            |          | 16 Cubed                   |          |
|                                   | Tricuspid Valve           |          | 22 Gibson                  |          |
|                                   |                           |          | 17 Teichholz               |          |
|                                   |                           |          | 18 Bullett                 |          |
|                                   |                           |          | 23 Single Plane Disc       |          |
|                                   |                           |          | 24 Biplane Disc            |          |
| [ORDER ] [RESET ] [SAVE ] [EXIT ] |                           |          |                            |          |

Figure 10–70. Select Report Display Order

Commands at the bottom of this Select Report Display allow the user to:

- |       |                                                                                 |
|-------|---------------------------------------------------------------------------------|
| ORDER | Erases all previous print order numbers and starts the order selection process. |
| RESET | Resets the print order to the factory default.                                  |
| SAVE  | Saves the chosen combination.                                                   |
| EXIT  | Exits the S–RP Process and changes to the previous screen.                      |





## Advanced Cardiac Calculations

---

### Operation and Messages

#### ORDER

The ORDER command is used to erase all previous designations and commence the order selection process. Select the ORDER command from the bottom menu. The message displayed is:

**“Reset? Press ‘ESC’ to cancel. ‘SET’ to confirm”**

**Set** erases the old order and starts the order selection process. **Esc** causes nothing to happen.

To designate a new order, use the **Trackball** to move the cursor to the desired selection.

Press **Set**. The order number is placed to the left of the measurement group name. The highlighted cursor stays in its position until it is moved by the **Trackball** to the next selection.

Continue this process until the desired selections have been made in the preferred order.



*NOTE: If Apical 4 Chamber–Diastole is selected, the same order number will be assigned to Apical 4 Chamber–Systole and vice-versa.*

*Measurement groups not selected will not be printed.*

#### RESET

The RESET command resets the displayed combination to the factory defaults. Prior to resetting, the following message is displayed:

**“Reset? Press ‘ESC’ to cancel. ‘SET’ to confirm”**

**Set** resets the order to factory defaults. **Esc** causes nothing to happen.



## *Advanced Cardiac Calculations*

---

### **Operation and Messages (cont'd)**

#### **SAVE**

The SAVE command will save the changes as the new report printing configuration. Prior to saving, the following message is displayed:

**“Overwrite? Press ‘ESC’ to cancel. ‘SET’ to confirm”**

**Set** saves the order changes. **Esc** causes nothing to happen.

If the Order process was not completed, the following message is displayed:

**“Complete the present set up first.”**

#### **EXIT**

If the exit command is selected before choosing, completing or saving a new sequence program, the following message is displayed:

**“Quit? Press ‘ESC’ to cancel. ‘SET’ to confirm”**

**Set** exits from this screen. **Esc** causes nothing to happen.

If a new sequence program was completed but not saved, the following message is displayed:

**“Overwrite? Press ‘ESC’ to cancel. ‘SET’ to confirm”**

**Set** saves the order changes and exits. **Esc** causes nothing to happen.



## *Advanced Cardiac Calculations*

---

This page left blank intentionally.



# Exam Preparation

## Introduction

Measurements and calculations derived from ultrasound images are intended to supplement other clinical procedures available to the attending physician. The accuracy of measurements is not only determined by the system accuracy, but also by use of proper medical protocols by the user. When appropriate, be sure to note any protocols associated with a particular measurement or calculation. Formulas and databases used within the system software that are associated with specific investigators are so noted. Be sure to refer to the original article describing the investigator's recommended clinical procedures.

## General Guidelines

Any measurement can be repeated by selecting that measurement again from the Sub-Menu.

The system will retain as many as 8 measurements, but the summary report retains only the last 3 measurements of each type.

The 3 report page measurements can be averaged and the average used in other calculations.



# Measurements

## Carotid Artery Measurements

Page 1 of the Vascular Sub-Menu contains measurements for the patient's right side. Page 2 of the Sub-Menu contains identical left side measurements.

### ECA, CCA, Bifurc and ICA

To measure the velocity of the blood through the right/left external carotid artery, common carotid artery, carotid artery bifurcation or internal carotid artery in the Doppler spectrum:

Image the appropriate artery.

|           |              |               |                |
|-----------|--------------|---------------|----------------|
|           |              |               |                |
| Rt<br>ECA | Rt<br>BIFURC | Rt<br>ICA/CCA | %<br>Steno     |
| Rt<br>CCA | Rt<br>ICA    |               | VASC<br>Report |

Page 1

|           |              |               |                |
|-----------|--------------|---------------|----------------|
|           |              |               |                |
| Lt<br>ECA | Lt<br>BIFURC | Lt<br>ICA/CCA | %<br>Steno     |
| Lt<br>CCA | Lt<br>ICA    |               | VASC<br>Report |

Page 2

Figure 11–1. Vascular Calculation Sub-Menu

Select the desired measurement from the Vascular Calculation Sub-Menu. The calculation mode is set and a “×” cursor with a vertical dotted line appears.

Move the cursor to the desired position on the Doppler spectrum and press **Set**.

The velocity measurement is completed and recorded on the Vascular Summary Report.



## Measurements

---

### RT (LT) ICA/CCA

To calculate a velocity ratio between the right (left) Internal carotid artery and common carotid artery:

|           |              |               |                |
|-----------|--------------|---------------|----------------|
|           |              |               |                |
| Rt<br>ECA | Rt<br>BIFURC | Rt<br>ICA/CCA | %<br>Steno     |
| Rt<br>CCA | Rt<br>ICA    |               | VASC<br>Report |

Page 1

|           |              |               |                |
|-----------|--------------|---------------|----------------|
|           |              |               |                |
| Lt<br>ECA | Lt<br>BIFURC | Lt<br>ICA/CCA | %<br>Steno     |
| Lt<br>CCA | Lt<br>ICA    |               | VASC<br>Report |

Page 2

Figure 11–2. Vascular Calculation Sub-Menus

Image the left or right internal carotid artery first. Select **RT ICA/CCA** from page 1 or **LT ICA/CCA** from page 2 of the Vascular Calculation Sub-Menu.

Move the cursor to the peak point of the Internal carotid artery Doppler spectrum and press **Set**.

Image the common carotid artery and press **Freeze**.

Press the **Measurement** key.

Move the cursor to the peak point (velocity) on the Doppler spectrum of the common carotid artery and press **Set**. The ICA/CCA calculation is complete and a calculated value is displayed.



*NOTE: If individual ICA and CCA measurements were already made, the ICA/CCA will be automatically calculated from previous measurements.*



## S/D (D/S) Ratio, RI, A/B Ratio or PI

Refer to the procedure found on 10–29.

|           |           |            |              |
|-----------|-----------|------------|--------------|
|           |           |            |              |
| S/D Ratio | A/B Ratio | % Steno    | TRACE AUTO   |
| PI        | RI        | HEART RATE | Transf CALCS |

Transf CALCS is only available with Realtime Doppler Calculation option installed

Figure 11–3. Vascular Calculation Sub-Menu

## Stenosis Ratio (% stenosis)

Refer to the procedure found on 10–23.

## Heart Rate (HR)

Refer to the procedure found on 10–30.

## Trace Auto

Use the Trace Auto Sub-Menu selection to turn on/off the Doppler Auto Trace function. See 10–28 for more details.

## Transf Calcs

This selection is only available if the Realtime Doppler Calculation option is installed and Realtime Doppler Calculations are displayed in the upper left corner of the screen.

Use the Trans Calcs selection to turn on/off the transfer of Automatic Doppler Calculations displayed on the screen to the report page. See 7–19 for details on Realtime Doppler operation.

## Helpful Hints



### Hints

If there is a series of measurements/calculations that are routinely taken in a specific sequence, the auto sequence feature of the system can be used to program these measurements/calculations to automatically appear on the display in the desired order. Refer to 14–81 for details on programming an auto sequence.



# Vascular Summary Report

## Introduction

The vascular summary report is structured to automatically display vascular measurements made at specific anatomical sites. An average of the 3 measurements can also be displayed. Calculated ratios are automatically summarized and displayed.

## Displaying the Summary Report

The summary report can be displayed at any time during the exam by selecting **VASC REPORT** from the Vascular Calculation Sub-Menu.

Data entered as new patient data is automatically summarized at the top of the summary report. Measurements are summarized below the patient data. Calculations and ratios appear at the bottom of the summary report.

To exit the summary report display and resume scanning, select **VASC REPORT** from the Sub-Menu to deactivate it.

|                   |         |         |          |
|-------------------|---------|---------|----------|
| HOSPITAL NAME     | Oper    | ID:     | 05/09/98 |
| ID:               | NAME:   |         | AGE:     |
| SEX:              | HEIGHT: | WEIGHT: | BSA:     |
| Ref MD:           | NOTE:   |         |          |
| MEASUREMENTS      | AVE     | 1       | 2 3      |
| Left ICA          | (       |         | )        |
| Left ECA          | (       |         | )        |
| Left BIFURC       | (       |         | )        |
| Left CCA          | (       |         | )        |
| Left ICA/CCA      | (       |         | )        |
| Right ICA         | (       |         | )        |
| Right ECA         | (       |         | )        |
| Right BIFURC      | (       |         | )        |
| Right CCA         | (       |         | )        |
| Right ICA/CCA     | (       |         | )        |
| CALCULATIONS      |         |         |          |
| S/D Ratio         | (       |         | )        |
| A/B Ratio         | (       |         | )        |
| % Stenosis        | (       |         | )        |
| Pulsatility Index | (       |         | )        |
| Resistivity Index | (       |         | )        |
| COMMENTS:         |         |         |          |

Figure 11–4. Vascular Report Page

<AVG> or <LAST> indicates that the Average or Latest value will be displayed in this column. This depends on the preset, *Average Activity*, in the Set Up/Preset Program menu page 3.





## *Vascular Summary Report*

---

### **Editing the Summary Report**

New patient data cannot be edited. Measurements can be edited by the operator.

To edit a measurement:

Select **VASC REPORT** from any page of the Vascular Calculation Sub-Menu.

Use the **Trackball** or **Return** key to move the cursor to the measurement that is to be edited. The measurement can be deleted by pressing the **BackSpace** or **Space** bar and then **Set**.

The average measurements, calculations and ratios are automatically updated to reflect the edited values.

When editing is complete, select **Exit** to leave the summary report and resume imaging.

### **Recording Summary Reports**

The summary report can be saved like any ultrasound image. Once it is displayed on the screen, it can be recorded on the VCR, printed on the B/W or color page printer, photographed by the multi-image camera, stored on a MOD with the Image Archive option or placed on regular paper with a line printer.

### **Line Printer Operation**

#### **USA Only**

Refer to 13–21 for details on line printer operation.



## Vascular Summary Report

### Vascular Calculation Formulas

| Calc Mnemonic                                 | Calc Name                                                                   | Input Measurements                                            | Formula                                               |
|-----------------------------------------------|-----------------------------------------------------------------------------|---------------------------------------------------------------|-------------------------------------------------------|
| RT ECA                                        | Right External Carotid Artery Velocity                                      | one Doppler blood flow peak velocity                          | RT ECA=v1<br>[cm/s or m/s]                            |
| RT CCA                                        | Right Common Carotid Artery Velocity                                        | one Doppler blood flow peak velocity                          | RT CCA=v1<br>[cm/s or m/s]                            |
| RT BIFURC                                     | Right Carotid Bifurcation Velocity                                          | one Doppler blood flow peak velocity                          | RT BIFURC=v1<br>[cm/s or m/s]                         |
| RT ICA                                        | Right Internal Carotid Artery Velocity                                      | one Doppler blood flow peak velocity                          | RT ICA=v1<br>[cm/s or m/s]                            |
| RT ICA/CCA                                    | Right Internal Carotid Artery Velocity/Common Carotid Artery Velocity Ratio | two Doppler blood flow peak velocities                        | RT ICA/CCA= $V_{ICA}/V_{CCA}$                         |
| LT ECA, LT CCA, LT BIFURC, LT ICA, LT ICA/CCA | Same as above, for Left Carotid Artery                                      | Same as above                                                 | Same as above                                         |
| A/B Ratio                                     | Velocities Ratio                                                            | two Doppler blood flow peak velocities                        | $A/B = V_1/V_2$                                       |
| % Stenosis                                    | Stenosis Ratio                                                              | two areas (by ellipse, trace, circle or distance)             | % Stenosis= $[1-(A_{residual}/A_{lumen})] \times 100$ |
| S/D Ratio                                     | Systolic Velocity/Diastolic Velocities Ratio                                | two Doppler blood flow peak velocities                        | $S/D = V_{systole}/V_{diastole}$<br>*NOTE*            |
| PI                                            | Pulsatility Index                                                           | two Doppler blood flow peak velocities and TAMAX              | $PI = (V_{max} - V_{diastole}) / TAMAX$<br>*NOTE*     |
| RI                                            | Resistivity Index                                                           | two Doppler blood flow peak velocities                        | $RI = (V_{max} - V_{diastole}) / V_{max}$<br>*NOTE*   |
| HR                                            | Heart Rate (beats/minute)                                                   | one 2 beat time interval (measured manually or automatically) | HR[BPM]= $120[sec] / 2beat\ time[sec]$                |

\*NOTE\*:  $V_{diastole} = V_{min}$  or  $V_{end-diastole}$  (depends on preset selection)

Table 11–1. Vascular Calculation Formulas

# Advanced Vascular

(software option)

## Overview

The Advanced Vascular software option provides increased measurement and reporting capabilities for vascular and venous studies.

## Menu Selections

The Advanced Vascular menu selections offer more flexibility than the Basic Software package.

In the Advanced Menus, the user can select

- left/right side on menu page one.
- three Doppler sites on page three.
- three stenosis sites on page four.

Each Doppler and Stenosis site can be edited (5 characters maximum) on the Vascular Report pages.



### % Stenosis

*The % Stenosis measurement can be performed by either a distance or area measurement in Advanced Vascular. The “% Stenosis Calculation” preset, found on Set Up/Preset Program page 4, erases all previous % Stenosis measurement results when this preset is changed. Depending on the measurement method used, the Vascular Report shows either an A (for area) or a D (for distance) after %Sten:.*

The menus displayed will depend on the preset for determination of systole/diastole (auto/manual) in the Set Up/Preset Program page 3 as shown in Figure 11–5 and Figure 11–6.



*NOTE: If the Realtime Doppler Calculation option is installed and activated, the user can choose to disable the automatic transfer of the realtime results to the result window and report page. Do this by selecting the manual calc option in page one of the Advanced Vascular Calculation Sub-Menu.*



## Advanced Vascular

### Menu Selections (cont'd)

|                 |          |          |                |
|-----------------|----------|----------|----------------|
|                 |          |          |                |
| Select Side     | Rt DICAd | Rt MICAd | Venous Report  |
| ▲<br>▼<br>RIGHT | Rt DICAs | Rt MICAs | Carotid Report |

Page One

|          |         |  |                |
|----------|---------|--|----------------|
|          |         |  |                |
| Rt PICAd | Rt ECAd |  | Venous Report  |
| Rt PICAs | Rt ECAs |  | Carotid Report |

Page Two

|          |          |          |                |
|----------|----------|----------|----------------|
|          |          |          |                |
| Rt DCCAd | Rt MCCAd | Rt PCCAd | Venous Report  |
| Rt DCCAs | Rt MCCAs | Rt PCCAs | Carotid Report |

Page Three

|         |         |              |                |
|---------|---------|--------------|----------------|
|         |         |              |                |
|         |         | Rt ICA /CCAd | Venous Report  |
| Rt VERT | Rt SUBC | Rt ICA /CCAs | Carotid Report |

Page Four

|                      |              |              |               |
|----------------------|--------------|--------------|---------------|
|                      |              |              |               |
| Select Location      | Rt S/D Ratio | Rt A/B Ratio | TRACE AUTO    |
| ▲<br>▼<br>\$\$\$\$\$ | Rt RI        | Rt PI        | Transf. CALCs |

Page Five

|                 |          |            |                |
|-----------------|----------|------------|----------------|
|                 |          |            |                |
| Select Location | Rt %Sten | HEART RATE | Venous Report  |
| ▲<br>▼<br>##### |          |            | Carotid Report |

Page Six

\$\$\$\$\$ = Doppler Sample Site      ##### = Stenosis Site

Transf CALCs is only available with Realtime Doppler Calculation option installed

Figure 11–5. Vascular Calculation Sub-Menu Manual Determination Systole/Diasatole

Menu Selections (cont'd)

|              |         |         |             |
|--------------|---------|---------|-------------|
|              |         |         |             |
| Select Side  | Rt DICA | Rt PICA | Rt DCCA     |
| ▲<br>▼ RIGHT | Rt MICA | Rt ECA  | Manual Calc |

Page One

|         |             |             |                |
|---------|-------------|-------------|----------------|
|         |             |             |                |
| Rt MCCA | Rt VERT     | Rt SUBC     | Venous Report  |
| Rt PCCA | Rt ICA/CCAs | Rt ICA/CCAd | Carotid Report |

Page Two

|                 |              |              |               |
|-----------------|--------------|--------------|---------------|
|                 |              |              |               |
| Select Location | Rt S/D Ratio | Rt A/B Ratio | TRACE AUTO    |
| ▲<br>▼ \$\$\$\$ | Rt RI        | Rt PI        | Transf. CALCs |

Page Three

|                 |          |            |                |
|-----------------|----------|------------|----------------|
|                 |          |            |                |
| Select Location | Rt %Sten | HEART RATE | Venous Report  |
| ▲<br>▼ #####    |          |            | Carotid Report |

Page Four

\$\$\$\$ = Doppler Sample Site      ##### = Stenosis Site

Transf CALCs is only available with Realtime Doppler Calculation option installed

Figure 11–6. Vascular Calculation Sub-Menu Auto Determination Systole/Diastole



## Advanced Vascular

### Menu Selections (cont'd)

|                                  |              |              |               |                       |             |             |                |
|----------------------------------|--------------|--------------|---------------|-----------------------|-------------|-------------|----------------|
|                                  |              |              |               |                       |             |             |                |
| Select Side                      | Rt DCIA      | Rt PICA      | Rt DCCA       | Rt MCCA               | Rt VERT     | Rt SUBC     | Venous Report  |
| ▲<br>▼<br>RIGHT                  | Rt MICA      | Rt ECA       | Manual Calc   | Rt PCCA               | Rt ICA/CCAs | Rt ICA/CCAd | Carotid Report |
| Page One                         |              |              |               | Page Two              |             |             |                |
| Select Location                  | Rt S/D Ratio | Rt A/B Ratio | TRACE AUTO    | Select Location       | Rt %Sten    | HEART RATE  | Venous Report  |
| ▲<br>▼<br>\$\$\$\$\$             | Rt RI        | Rt PI        | Transf. CALCS | ▲<br>▼<br>#####       |             |             | Carotid Report |
| Page Three                       |              |              |               | Page Four             |             |             |                |
| \$\$\$\$\$ = Doppler Sample Site |              |              |               | ##### = Stenosis Site |             |             |                |

Figure 11–7. Advanced Vascular Calculation Sub-Menu with Realtime Doppler ON



*NOTE: The Advanced Vascular Sub-Menu changes when the Realtime Doppler option is installed and turned ON. Shaded areas denote the differences.*



### Hints

With the Realtime Doppler Calculation Option, the DICA, MICA, PICA, ECA, DCCA, MCCA, PCCA, VERT and SUBC measurements will use  $V_{\max}$  and  $V_d$  derived from the Realtime Doppler Calculations.

If it is desired not to use the Realtime Doppler Calculations, select **MANUAL CALC**. When Manual Calc is selected, the user must select the measurement. The Trackball and Set key are then used to measure  $V_{\max}$  and  $V_d$ .

All results, automatic or manual, are displayed on the monitor and transferred to the report page.

## Report Page Layout

Unlike the Basic Vascular Option, Advanced Vascular offers three report pages.

### Advanced Vascular Report Page 1

Page one is a summary of the left and right carotid measurements plus Doppler measurements and % Stenosis at three different sites.

|               |       |            |         |                                                                                   |     |           |           |
|---------------|-------|------------|---------|-----------------------------------------------------------------------------------|-----|-----------|-----------|
| HOSPITAL NAME |       | NAME:      |         | Oper ID:                                                                          |     | 05/09/98  |           |
| ID:           | SEX:  | HEIGHT:    | WEIGHT: | BSA:                                                                              | HR: | AGE:      | BPM       |
| Ref MD:       | NOTE: |            |         |                                                                                   |     |           |           |
|               | SYST  | RIGHT DIAS |         |                                                                                   |     | LEFT SYST | 1/3- DIAS |
| DICA(m/s)     |       |            |         |  |     |           |           |
| MICA(m/s)     |       |            |         |                                                                                   |     |           |           |
| PICA(m/s)     |       |            |         |                                                                                   |     |           |           |
| ECA(m/s)      |       |            |         |                                                                                   |     |           |           |
| DCCA(m/s)     |       |            |         |                                                                                   |     |           |           |
| MCCA(m/s)     |       |            |         |                                                                                   |     |           |           |
| PCCA(m/s)     |       |            |         |                                                                                   |     |           |           |
| ICA/CCA       |       |            |         |                                                                                   |     |           |           |
| VERT(m/s)     |       | A          | R       | B                                                                                 |     | A         | R         |
| SUBC(m/s)     |       | T          | B       | M                                                                                 |     | T         | B         |
| %Stenosis(%)  | 1     | 2          | 3       |                                                                                   | 4   | 5         | 6         |
|               | 1     | 2          | 3       |                                                                                   | 4   | 5         | 6         |
| PI            |       |            |         |                                                                                   |     |           |           |
| RI            |       |            |         |                                                                                   |     |           |           |
| A/B           |       |            |         |                                                                                   |     |           |           |
| S/D           |       |            |         |                                                                                   |     |           |           |
| COMMENTS:     |       |            |         |                                                                                   |     |           |           |

Figure 11–8. Vascular Report Page 1

Provisions are made to annotate the vertebral artery in diastole as: A = Antegrade, R = Retrograde or B = Biphasic. The Subclavian Artery can be labelled as: T= Triphasic, B = Biphasic or M = Monophasic.


Comments can be added in the Vascular diagram area.



## Advanced Vascular

### Advanced Vascular Report Page 1 (cont'd)

To note one of the three choices for Vert and Subc:



Use the **Trackball** or **Arrow** keys to place the cursor on top of the desired letter.

Press **Set**.

The desired letter will be placed in brackets (i.e. A<R>B or T B<M>).

The user can freely edit the stenosis location names (%A1%%) and Doppler site names (%B1%%). A maximum of five characters can be used for each site name.

### Advanced Vascular Report Page 2

Page two is a list for the left and right sides of the last three calculations with maximum average or latest values for % stenosis and the three designated Doppler sites.

Figure 11–9 shows an example of the layout for page two.

|               |         |          |       |            |      |
|---------------|---------|----------|-------|------------|------|
| HOSPITAL NAME |         | Oper ID: |       | 05/09/98   |      |
| ID:           | NAME:   | BSA:     |       | AGE:       |      |
| SEX:          | HEIGHT: | WEIGHT:  | NOTE: |            |      |
| Ref MD:       |         |          |       |            |      |
| Stn(%)MAX     | RIGHT   | 2        | 3     | Stn(%) MAX | LEFT |
| 1             | 1       |          |       | 4          | 1    |
| 2             | 2       |          |       | 5          | 2    |
| 3             | 3       |          |       | 6          | 3    |
| PI            |         |          |       | PI         |      |
| 1             |         |          |       | 4          |      |
| 2             |         |          |       | 5          |      |
| 3             |         |          |       | 6          |      |
| RI            |         |          |       | RI         |      |
| 1             |         |          |       | 4          |      |
| 2             |         |          |       | 5          |      |
| 3             |         |          |       | 6          |      |
| A/B           |         |          |       | A/B        |      |
| 1             |         |          |       | 4          |      |
| 2             |         |          |       | 5          |      |
| 3             |         |          |       | 6          |      |
| S/D           |         |          |       | S/D        |      |
| 1             |         |          |       | 4          |      |
| 2             |         |          |       | 5          |      |
| 3             |         |          |       | 6          |      |
| COMMENTS:     |         |          |       |            |      |

Figure 11–9. Vascular Report Page 2



### Advanced Vascular Report Page 3

Page three is a list of the last three measurements, with the maximum, average, or latest, for the carotid measurements taken in systole/diastole on the left and right sides. Figure 11–10 shows an example of the layout for page three.

|               |         |          |   |          |     |          |   |   |
|---------------|---------|----------|---|----------|-----|----------|---|---|
| HOSPITAL NAME |         |          |   | Oper ID: |     | 05/09/98 |   |   |
| ID:           | NAME:   |          |   | WEIGHT:  |     | AGE:     |   |   |
| SEX:          | HEIGHT: | WEIGHT:  |   | BSA:     |     |          |   |   |
| Ref MD:       | NOTE:   |          |   |          |     |          |   |   |
| SYSTOLE       |         | DIASTOLE |   | +3/3     |     |          |   |   |
| RIGHT         | MAX     | 1        | 2 | 3        | MAX | 1        | 2 | 3 |
| DICA(m/s)     |         |          |   |          |     |          |   |   |
| MICA(m/s)     |         |          |   |          |     |          |   |   |
| PICA(m/s)     |         |          |   |          |     |          |   |   |
| ECA(m/s)      |         |          |   |          |     |          |   |   |
| DCCA(m/s)     |         |          |   |          |     |          |   |   |
| MCCA(m/s)     |         |          |   |          |     |          |   |   |
| PCCA(m/s)     |         |          |   |          |     |          |   |   |
| ICA/CCA       |         |          |   |          |     |          |   |   |
| VERT(m/s)     |         |          |   |          |     |          |   |   |
| SUBC(m/s)     |         |          |   |          |     |          |   |   |
|               |         |          |   |          |     |          |   |   |
| LEFT          | MAX     | 1        | 2 | 3        | MAX | 1        | 2 | 3 |
| DICA(m/s)     |         |          |   |          |     |          |   |   |
| MICA(m/s)     |         |          |   |          |     |          |   |   |
| PICA(m/s)     |         |          |   |          |     |          |   |   |
| ECA(m/s)      |         |          |   |          |     |          |   |   |
| DCCA(m/s)     |         |          |   |          |     |          |   |   |
| MCCA(m/s)     |         |          |   |          |     |          |   |   |
| PCCA(m/s)     |         |          |   |          |     |          |   |   |
| ICA/CCA       |         |          |   |          |     |          |   |   |
| VERT(m/s)     |         |          |   |          |     |          |   |   |
| SUBC(m/s)     |         |          |   |          |     |          |   |   |

Figure 11–10. Vascular Report Page 3

On pages two and three, the user can choose to display the maximum, average or latest value. The preset “Average Activity: OFF/ON” in the Set Up/Preset Program menu page 3 controls AVG/LAST.

ON = AVG value  
OFF = LAST value

The preset “MAX Velocity (Advanced Vascular) : OFF/ON” overrides Average Activity.

ON = MAX value  
OFF = Average Activity selection



## Venous Comments

The Advanced Vascular Option also provides a Venous Comments page to report observations made during an exam. Figure 11–11 shows the layout of the Venous Comments page.

|                                                        |                         |          |         |                         |       |
|--------------------------------------------------------|-------------------------|----------|---------|-------------------------|-------|
| HOSPITAL NAME                                          |                         | Oper ID: |         | 05/09/98                |       |
| ID:                                                    | NAME:                   | BSA:     |         | AGE:                    |       |
| SEX:                                                   | HEIGHT:                 | WEIGHT:  |         |                         |       |
| Ref MD:                                                | NOTE:                   |          |         |                         |       |
| RIGHT                                                  |                         |          |         |                         |       |
| Comp?                                                  | Doppler                 | LEFT     |         |                         |       |
| CFV                                                    | Y N P Sp Ag Ph Rf Cn PI | Comp?    | Doppler | Y N P Sp Ag Ph Rf Cn PI |       |
| PSEV                                                   | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| MSFV                                                   | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| DSFV                                                   | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| PROFV                                                  | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| PPopV                                                  | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| MPopV                                                  | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| DPopV                                                  | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| PPTV                                                   | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| MPTV                                                   | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| DPTV                                                   | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| PPerV                                                  | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| MPerV                                                  | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| DPerV                                                  | Y N P Sp Ag Ph Rf Cn PI | Y N P    | Sp Ag   | Ph Rf                   | Cn PI |
| PGSAPH Y N P Sp Ag Ph Rf Cn PI Y N P Sp Ag Ph Rf Cn PI |                         |          |         |                         |       |
| MGSAPH Y N P Sp Ag Ph Rf Cn PI Y N P Sp Ag Ph Rf Cn PI |                         |          |         |                         |       |
| DGSAPH Y N P Sp Ag Ph Rf Cn PI Y N P Sp Ag Ph Rf Cn PI |                         |          |         |                         |       |
| COMMENTS:                                              |                         |          |         |                         |       |

Figure 11–11. Venous Comment Report Page

## Abbreviations

Abbreviations found on the Venous Comments page designate the following:

|       |   |                          |
|-------|---|--------------------------|
| CFV   | = | Common Femoral Vein      |
| SFV   | = | Superficial Femoral Vein |
| PROFV | = | Profunda Femoris Vein    |
| PopV  | = | Popliteal Vein           |
| PTV   | = | Posterior Tibial Vein    |
| PerV  | = | Peroneal Vein            |
| GSAPH | = | Greater Saphenous Vein   |

Prefixes found on some abbreviations designate:

|   |   |          |
|---|---|----------|
| P | = | Proximal |
| M | = | Mid      |
| D | = | Distal   |



## Observations

Observations can be noted for each vein for both the left and right sides. The first section of observations denotes compressability (Comp?) of the vein. The user can note if compression could be accomplished on the vein by designating one of the following:

Y = Yes  
N = No  
P = Partial

In the second part, the user can comment on observations made during the Doppler study. The user can select one or more of the following designations:

|    |   |             |                                                                                              |
|----|---|-------------|----------------------------------------------------------------------------------------------|
| Sp | = | Spontaneous | Indicates the vein is easily imaged without aids.                                            |
| Ag | = | Augmented   | Velocity signal achieved by distal limb compression or release of proximal limb compression. |
| Ph | = | Phasicity   | Indicates the waxing and waning of the velocity signal with respiration.                     |
| Rf | = | Reflux      | Indicates reverse or backward venous flow was observed.                                      |
| Cn | = | Continuous  | Indicates no change in velocity signal after compression.                                    |
| Pl | = | Pulsatility | Indicates the vein responds to the cardiac cycle and not to the respiratory cycle.           |

## Selecting Abbreviations

In order to select abbreviations on this comment page:

- Move the highlight cursor to the desired abbreviation and press **Set**.
- The abbreviation will be enclosed by brackets "< >".
- To cancel the brackets, move the highlight cursor to the bracketed abbreviation and press **Set**.
- The brackets are deleted.

## Comments

General Comments can be written at the bottom of each Vascular Report page. The comments are separate for each page. Individual notes can be entered on each one.



## Advanced Vascular

### Vascular Calculation Formulas

| Calc Mnemonic | Calc Name                                                 | Input Measurements                   | Formula                   |
|---------------|-----------------------------------------------------------|--------------------------------------|---------------------------|
| Rt DICAd      | Right Distal Internal Carotid Artery Velocity, diastole   | one Doppler blood flow peak velocity | RT DICAd=v1 [cm/s or m/s] |
| Rt DICAs      | Right Distal Internal Carotid Artery Velocity, systole    | one Doppler blood flow peak velocity | RT DICAs=v1 [cm/s or m/s] |
| Rt MICAd      | Right Mid-Internal Carotid Artery Velocity, diastole      | one Doppler blood flow peak velocity | RT MICAd=v1 [cm/s or m/s] |
| Rt MICAs      | Right Mid-Internal Carotid Artery Velocity, systole       | one Doppler blood flow peak velocity | RT MICAs=v1 [cm/s or m/s] |
| Rt PICAd      | Right Proximal Internal Carotid Artery Velocity, diastole | one Doppler blood flow peak velocity | RT PICAd=v1 [cm/s or m/s] |
| Rt PICAs      | Right Proximal Internal Carotid Artery Velocity, systole  | one Doppler blood flow peak velocity | RT PICAs=v1 [cm/s or m/s] |
| Rt ECAd       | Right External Carotid Artery Velocity, diastole          | one Doppler blood flow peak velocity | RT ECAd=v1 [cm/s or m/s]  |
| Rt ECAs       | Right External Carotid Artery Velocity, systole           | one Doppler blood flow peak velocity | RT ECAs=v1 [cm/s or m/s]  |
| Rt DCCAd      | Right Distal Common Carotid Artery Velocity, diastole     | one Doppler blood flow peak velocity | RT DCCAd=v1 [cm/s or m/s] |
| Rt DCCAs      | Right Distal Common Carotid Artery Velocity, systole      | one Doppler blood flow peak velocity | RT DCCAs=v1 [cm/s or m/s] |
| Rt MCCAd      | Right Mid-Common Carotid Artery Velocity, diastole        | one Doppler blood flow peak velocity | RT MCCAd=v1 [cm/s or m/s] |
| Rt MCCAs      | Right Mid-Common Carotid Artery Velocity, systole         | one Doppler blood flow peak velocity | RT MCCAs=v1 [cm/s or m/s] |
| Rt PCCAd      | Right Proximal Common Carotid Artery Velocity, diastole   | one Doppler blood flow peak velocity | RT PCCAd=v1 [cm/s or m/s] |
| Rt PCCAs      | Right Proximal Common Carotid Artery Velocity, systole    | one Doppler blood flow peak velocity | RT PCCAs=v1 [cm/s or m/s] |

Table 11–2. Advanced Vascular Calculation Formulas

Vascular Calculation Formulas (cont'd)

| Calc Mnemonic                              | Calc Name                                                                   | Input Measurements                               | Formula                                                 |
|--------------------------------------------|-----------------------------------------------------------------------------|--------------------------------------------------|---------------------------------------------------------|
| Rt ICA/CCAd<br>Rt ICA/CCAs                 | Right Internal Carotid Artery Velocity/Common Carotid Artery Velocity Ratio | two Doppler blood flow peak velocities           | $RT\ XICA/XCCA = V_{XICA}/V_{XCCA}$                     |
| Rt VERT                                    | Right Vertebra Artery Velocity, systole                                     | one Doppler blood flow peak velocity             | $RT\ VERT = v1$<br>[cm/s or m/s]                        |
| Rt SUBC                                    | Right Subclavian Artery Velocity, systole                                   | one Doppler blood flow peak velocity             | $RT\ SUBC = v1$<br>[cm/s or m/s]                        |
| Rt PI                                      | Right Pulsatility Index                                                     | two Doppler blood flow peak velocities and TAMAX | $RT\ PI = (V_{max} - V_{diastole}) / TAMAX$<br>*NOTE*   |
| Rt RI                                      | Right Resistivity Index                                                     | two Doppler blood flow peak velocities           | $RT\ RI = (V_{max} - V_{diastole}) / V_{max}$<br>*NOTE* |
| Rt S/D Ratio                               | Right Systolic Velocity/Diastolic Velocities Ratio                          | two Doppler blood flow peak velocities           | $RT\ S/D = V_{systole} / V_{diastole}$<br>*NOTE*        |
| Rt A/B Ratio                               | Right side Velocities Ratio                                                 | two Doppler blood flow peak velocities           | $RT\ A/B = V_1 / V_2$                                   |
| Same as above, but for Left Carotid Artery | Same as above, but for Left Carotid Artery                                  | Same as above                                    | Same as above                                           |

\*NOTE\*:  $V_{diastole} = V_{min}$  or  $V_{end-diastole}$  (depends on preset selection)

Table 11–2. Advanced Vascular Calculation Formulas (cont'd)



## *Advanced Vascular*

---

This page left blank intentionally.



# Urology Basic Calculations

## Overview

The Urology Calculations factory default Sub-Menus consists of 3 pages. The calculations available are the same as the general menus for Radiology/Abdomen and Small Parts. They are:

|        |            |  |                   |
|--------|------------|--|-------------------|
|        |            |  |                   |
| VOLUME | %<br>STENO |  |                   |
| ANGLE  |            |  | Urology<br>Report |

Urology Report is only available with Urology Calculation option installed

|              |              |               |                 |
|--------------|--------------|---------------|-----------------|
|              |              |               |                 |
| S/D<br>Ratio | A/B<br>Ratio | HEART<br>RATE | TRACE<br>AUTO   |
| PI           | RI           |               | Transf<br>CALCS |

Transf CALCS is only available with Realtime Doppler Calculation option installed

|         |    |               |               |
|---------|----|---------------|---------------|
|         |    |               |               |
| Max PG  | CO | HEART<br>RATE | TRACE<br>AUTO |
| Mean PG | SV | FV            | FVO           |

Figure 12–1. Urology Calculation Menus

Refer to the chapters *Abdomen and Small Parts (8)* and *Cardiology (10)* for steps to perform these measurements and calculations.

## Presumed Circle Area Ratio (PCAR)

Presumed Circle Area Ratio (PCAR) Measurement is available with the Urology Calculation Option.

This measurement must first be selected in the Measurement VFD Softmenu Setup Submenu, which is accessed by selecting the "Measurement Menu B" preset on Preset Program page 7. It will then be accessible from the B-Mode Measurement Sub-Menu.

|      |  |  |  |
|------|--|--|--|
|      |  |  |  |
| PCAR |  |  |  |
|      |  |  |  |

Figure 12–2. Urology Calculation Sub-Menu

Select **PCAR** from the Calculation Sub-Menu. The trace cursor appears in B-Mode.

Move the cursor to the measurement start point and press **Set**. Trace the measurement area and press **Set** to fix the PCAR value displayed on the bottom of the screen and calculate the area.

## Stepper Volume (STVOL)

The ERB7 probe is designed to be used with the mechanical stepping device option and needle placement guide. The needle placement guide matches the electronic needle placement grid display on the monitor as the biopsy guidezone is displayed.

This measurement must first be selected in the Measurement VFD Softmenu Setup Submenu, which is accessed by selecting the "Measurement Menu B" preset on Preset Program page 7. It will then be accessible from the B-Mode Measurement Sub-Menu.

Refer to 12–7 for further details.





# Urology Calculation

(software option)

## Urology Summary Report

### Urology Report Page

The Urology Report Page is provided to improve productivity and facilitate consistency in Urology procedures. It is available by selecting Urology Report from the measurement Sub-Menu.

|              |          |          |          |
|--------------|----------|----------|----------|
| ID:          | NAME:    | Oper ID: | 25/02/99 |
|              |          | AGE:     |          |
| Ref MD:      | NOTE:    |          |          |
| MEASUREMENTS |          |          |          |
| VOLUME #1    |          |          |          |
| VOLUME #2    |          |          |          |
| VOLUME #3    |          |          |          |
| LESION SITE  |          | BIOPSY   |          |
| 1            |          | YES      | NO       |
| 2            |          | YES      | NO       |
| 3            |          | YES      | NO       |
| 4            |          | YES      | NO       |
| PAP :        |          |          |          |
| PSA :        | ng/ml    |          |          |
| PSAD:        | BASED ON | VOLUME   | #N       |
| PSAD:        | BASED ON | VOLUME   | #N       |
| DRE :        |          |          |          |
| :            |          |          |          |
| COMMENTS:    |          |          |          |

Figure 3–45. Urology Report Page



---

## *Urology Calculations (Option)*

---

### **Urology Report Page (cont'd)**

The first, or top portion, is generally patient data that was entered on the patient entry menu at the beginning of the exam. This area cannot be edited.

The second is the measurements area. There are 3 fields for volume measurements. The first, second and third volume measurements taken are recorded as Volume #1, Volume #2 and Volume #3, respectively. If the volume measurements were performed by STVOL, "STVOL" is displayed.

The system stores the 3 most recent volume measurements. If a fourth measurement is made, Volume #1 is erased. Volume #2 becomes #1. Volume #3 becomes #2 and the fourth measurement is Volume #3. The very first measurement is lost.

Comments, pertinent for that specific volume measurement, may be added to the right of the volume measurement.

The third area is for Lesion Sites. A description can be input (maximum 20 characters each) for each lesion site. Whether the lesion site was biopsied needs to be selected ("Y" for yes or "N" for no).

The fourth area is for calculations. The PAP (Prostatic Acid Phosphatase concentration) and PSA (Prostatic Specific Antigen concentration) need to be input.

The PSAD is automatically calculated once the single digit field following the words "BASED ON VOLUME" is entered (1, 2, or 3). This represents which volume measurement (corresponding to the entire gland) should be used to compute PSAD. The calculation is dependent on the default selected in the "PSA Calculation" preset found on Set Up/System Parameters page 3.

The DRE (Digital Rectal Exam) result are input below the PPSA calculations. Any other exam which is desired to be reported can be input under DRE.

The fifth area is for comments and findings.



## *Urology Calculations (Option)*

---

### **PSA Measurement**

The “*PSA Calculation*” preset, located on System Parameters page 3, allows PSA to be calculated two ways. The two selections are:

PSAD: Prostatic Specific Antigen (PSA) Density – defined as:

$$\text{PSAD} = \text{PSA} / \text{Volume}$$

PPSA: Predicted Prostate Specific Antigen – defined as:

$$\text{PPSA} = \text{Volume} \times \text{PPSA Coefficient}$$

PPSA Coefficient is entered in the “*PPSA Coefficient (1)*” and “*PPSA Coefficient (2)*” presets, located on Preset Program page 4.

PSAD and PPSA are calculated on the Urology Report Page.

## Stepper Volume Calculation

### Overview

Using the ERB7 probe with the stepping device option, a new volume measurement, named STVOL, is available with the Urology Calculation Option.

A stepper volume calculation (STVOL) is now available with the LOGIQ™ 500 when using the ERB7 probe with the stepping device option.

### ERB7 Probe Preparation

The optional 7.0 MHz Transaxial (ERB7) Probe is a micro-convex and linear bi-plane probe for transaxial prostate scanning.





Figure 12–3. ERB7 Probe and cable

### ERB7 Probe Preparation (cont'd)

1. Remove the probe from the box and carefully examine it for any damage.
2. Clean and then disinfect/sterilize the ERB7 Probe.
3. Inspect a sterile/sanitary sheath. Place a small amount of ultrasound gel inside the sheath tip (the gel is between the sheath inner surface and the probe aperture).
4. Place the sheath tip over the probe aperture and then pull the sheath end toward the probe handle.
5. Place a rubberband/twist lock over the sheath at the end of the probe shaft. Ensure the rubberband/twist lock is tight around the sheath. Rub your finger over the tip of the probe to ensure all air bubbles have been eliminated.
6. Place a small amount of ultrasound gel on the gel-filled sheath tip outer surface.



**NOTE:** Remember to rinse all probe sheaths of powder before placing on the probe. Powder can degrade the scan image displayed.



- 1 Twist Lock
- 2 Sterile/Sanitary Outer Sheath
- 3 Optional Finger Cot for Water Path
- 4 Optional Water Path Inlet
- 5 Optional Finger Cot Rubberband for Water Path

Figure 12–4. Sterile/Sanitary Sheath Placement

## Performing STVOL Measurement

Stepper Volume (STVOL) is the method used to calculate the volume of an organ using the LOGIQ™ 500 Urology software, Transaxial Probe and a mechanical stepping device that moves the probe in fixed increments.

### CAUTION



The stepper volume measurement works only in B-Mode with Convex Scan (NOT available Linear scan).

### Stepper Volume Formula

Stepper Volume (STVOL) is the method used to calculate the volume of an organ using the LOGIQ™ 500 Urology software, Transaxial Probe and a mechanical stepping device that moves the probe in fixed increments.

The calculation is based on the fact that each step or area measurement is taken at equal stepper increments. The area measured at each slice is then used to compute the total volume of the organ, according to the following model:

For **N** slices, there are **N-1** volumes between the slices plus a small volume at each end.

It is assumed that the small volumes of the end caps are cones with a base equal to the measured area of the end slice and height equal to the slice spacing. The cones will have volumes  $V=A_1 d/3$  and  $A_N d/3$  for slices **1** and **N** respectively, where **d** is the spacing between the slices.

The volumes between the slices are assumed to be segments of a paraboloid of revolution, where the volume between any slice **n** and **n+1** is  $(A_n+A_{n+1})d/2$ .

The sum of the volumes between slices is:

$$V=d [(A_1+A_2/2)+(A_2+A_3/2)+....+(A_{n-1}+A_n/2)]$$

Therefore, the total volume of the organ is calculated as:

$$V=d [(5/6)A_1+A_2+A_3+....+(5/6)A_N]$$

Press **Freeze** to have an Image of the organ to be measured. The probe should be inserted to its most cephalic position (i.e base of the prostate).

Press **Measurement** to fix the STVOL measurement.



## *Urology Calculations (Option)*

---

### **Mechanical Stepper/Needle Placement Guide Option**

This mechanical stepping device is intended to be used with the LOGIQ™ 500 and the 7MHz Transaxial (ERB7) Probe.

The variable position mechanical stepping device is used to hold the 7MHz Transaxial (ERB7) Probe, lock the probe in place and move the probe in small increments throughout an organ in cross-sectional slice. Used with the LOGIQ™ 500 Stepper Volume Calculation, a complete three-dimensional volume calculation can be made of an organ.

The needle placement guide template matches the LOGIQ™ 500 electronic needle placement grid for the 7MHz Transaxial (ERB7) Probe to aid in needle placement.

The universal table mounting system holds the stepper and is compatible with most radiology and urologic tables.

#### **CAUTION**



Standard stepper increment is 5.0mm. Assembly for 2.5mm increment movement is available. Therefore it is necessary to match the stepper movement to the ultrasound system software.

Refer to the Stepping Device Instructions for details.



# Recording Images

## Image Memory

### Storage

Without the optional Cine Memory extension, the following number of images or proper combination of formats can be saved in memory:

- Single Format 8 images
- Dual Format 4 images
- Timeline Format 4 (NTSC) or 3 (PAL)

With the optional Cine Memory extension, the following number of images or proper combinations of formats can be saved in memory:

- Single Format 8 images
- Dual Format 8 images
- Timeline Format 4 (NTSC) or 3 (PAL)

The LOGIQ™ 500 has storage space for images in system memory.

This storage is temporary. The images are erased when the **New Patient** key is pressed or power is turned off.






## Storage (cont'd)

To save images in system memory:

Press **Freeze** to stop image acquisition.



Use the **Cine Scroll** control, if necessary, to display the best image.

Press the **Image Memory** key. The displayed image is saved to system memory and the number of images currently stored is displayed at the bottom of the screen.

The system gives a Memory Full warning beep each time an image is about to be stored in a full memory.

Pressing **Image Memory** after the warning deletes the oldest image in system memory and saves the current displayed image.



## Recording Images

---

### Recall

#### Image Recall



To recall images stored in system memory for review or archival:

- Press the **Image Recall** key. The soft-menu displays the Image Recall Sub-Menu with the display mode and time the images were stored.
- The last image to be placed into Image Memory is displayed.

Use the appropriate Sub-Menu rocker switch to select the desired image to recall.

After the recalled image is displayed, no image post processing can be performed. Only Comment, Measurement, Body Pattern, Clear, Ext Video and Mic functions are available.

Press **Freeze** to continue scanning or press **Image Recall** to return to the previous system status before Image Recall was initially pressed.

### Helpful hints



#### Hints

- Images stored in system memory will be lost if the **New Patient** key is pressed or power is turned off.
- Images are stored in system memory on a “first in, first out” basis (maximum 8 images).
- Use the highest quality videotape possible and record in the slowest speed possible.



## Peripheral Devices

The LOGIQ™ 500 can save scan images to a variety of optional devices. Available archiving options are:

- Black/white video page printer
- Color video page printer
- Video cassette recorder (VCR or VTR)
- Multi-image Camera
- Line Printer – USA Only
- Laser camera

| Device                           | Manufacturer          | Model       | Video Signal |
|----------------------------------|-----------------------|-------------|--------------|
| Black & White Video Page Printer | SONY                  | UP-890MD    | NTSC         |
|                                  | SONY                  | UP-890CE    | PAL          |
| Color Video Page Printer         | SONY                  | UP-1800     | NTSC         |
|                                  | SONY                  | UP-2800     | PAL          |
|                                  | SONY                  | UP-2850     | PAL          |
|                                  | SONY                  | UP-2900     | NTSC         |
|                                  | SONY                  | UP-2950     | NTSC         |
| VCR                              | SONY                  | SVO-9500MD  | NTSC         |
|                                  | SONY                  | SVO-9500MDP | PAL          |
| Multi-image Camera               | International Imaging | IIE-460     | NTSC         |

Table 8–2. Approved Optional Recording Devices

### CAUTION



**DO NOT** connect any probes or accessories without approval by GE.

All peripherals are considered options to the console. Those listed in this section have been tested and verified to be compatible with the LOGIQ™ 500 system. Please refer to the local Sales or Field Service representative for the latest list of approved peripherals.

### WARNING



When using any peripheral device or accessory, observe all danger messages, warnings, and cautions given in peripheral operator manuals.



## Recording Images

---

### Black/White Video Page Printer


#### Operations

Remote control of the B/W printer is limited to the print function only. Adjustments to print quality are done on the page printer. No status or error messages are available to be displayed on the LOGIQ™ 500.

There is a six second delay built in to the **Record** key recognition. This eliminates improper operation if the **Record** key is pressed more than once in rapid succession.

To print an image:

Press **Freeze** to stop image acquisition.



Freeze



Record 1

*NOTE: Remember that **Cine Scroll** may be used to look at previous image frames to obtain the best image.*

Press **Record 1** to activate the print function on the standard B/W Video Page Printer.



*NOTE: The **Record 1** and **Record 2** keys can be preset to activate almost any peripheral device. This is done by assigning the key control to the output port to which the peripheral is attached in the Set Up/System Parameters menu page 5.*


For details on the page printer operation, consult the Sony Operator Manual provided with the printer.



## Black/White Video Page Printer (cont'd)

### Configurations

The Sony UP-890 Black and White Video Page Printer is mounted on the left side of the LOGIQ™ 500 under the keyboard. The Sony connections are Video, AC power, and Remote Shutter (Expose).



1 Black/White Page Printer

Figure 13–1. B/W Page Printer Location



## Recording Images

### Black/White Video Page Printer (cont'd)

Front Panel of LOGIQ 500




Figure 13-2. B/W Printer Connections



## Color Video Page Printer

### Operations

The color video page printer allows for recording 1 or 4 images on a single sheet of paper in color or Black/White.

When the LOGIQ™ 500 is powered on, the Preset/Set Up selections set the memory mode (1 or 4 images) and the input signal select (RGB, Video, S-Video) on the color printer. All other functions need to be selected at the printer control panels.

The **Record 1** or **Record 2** keys can be programmed in the Preset/Set Up function to activate the print function on the Sony color page printer.

|   |
|---|
| 1 |
|---|

One image format: After the **Record** key is pressed, the image is first stored in printer memory and then printed.

|   |   |
|---|---|
| 1 | 2 |
| 3 | 4 |

Four Image Format: The first, second and third time the **Record** key is pressed, the image is stored in printer memory. The fourth time the **Record** key is pressed, the image is stored in printer memory and all four are printed on one sheet.

The color video page printer provides some feedback if there is an error in the print process. The following error messages could be displayed due to color printer problems:

- Paper is jamming. Check paper.
- Check ribbon cassette setting.
- Check paper cassette setting.
- No paper. Set paper.
- Place paper print side up.

For details on the page printer operation, consult the Sony Operator Manual provided with the printer.

## Recording Images

---

### Color Video Page Printer (cont'd)

#### Configurations

The Sony UP-1800 Color Video Page Printer is mounted behind the keyboard under the monitor. It is mounted facing the front of the system. The Sony UP-1800 connections are Red, Green, Blue, Sync, and RS-232 Remote Control.




Figure 13–3. Color Page Printer Location, Front Mount



Color Video Page Printer (cont'd)




Figure 13-4. Connections to LOGIQ™ 500




Figure 13-5. Connections to Color Page Printer



## Recording Images

---

### Video Cassette Recorder (VCR)

#### Operations

An optional video cassette recorder (VCR or VTR) is available for the LOGIQ™ 500. The optional VCR is S-VHS compatible for improved recording quality.

The LOGIQ™ 500 supports remote control of the Sony SVO-9500MD/MDP S-VHS Video Cassette Recorder only. Other VCRs can record and playback video with the LOGIQ™ 500, however remote control, advanced search functions and auto calibration will not be supported. Control of approved VCRs other than the Sony SVO-9500MD/MDP must be done by the VCR front panel.



*NOTE: Measurements from VCR Playback video requires the optional VCR measurement playback board for proper operation.*

Explanation of VCR operations assumes that the Sony SVO-9500 is being used. It also assumes that the "Heading VCR Playback" option is installed.

#### CAUTION



The system can keep track of patient and tape information. The system can search a tape for patient images. However, the system CANNOT:

- Stop or indicate when it comes to the end of the last study on the tape.
- Prevent recording over previous studies if record is pressed while positioned in the middle of a study.

Without the Heading VCR Playback option:

- The hard drive does not store patient or tape id information.
- Auto calibration for playback measurements is not available.
- Image, Patient and Tape Search functions are not available.

The Record keys and foot switch can be programmed in the Preset/Set Up menu to activate the Record/Pause function for the VCR.



## Video Cassette Recorder (VCR) (cont'd)

### Configurations

The Sony SVO-9500MD is mounted behind the keyboard and under the monitor. This recorder takes advantage of all the advanced tracking features of the LOGIQ™ 500. The Sony SVO-9500MD connections are S-VHS Video In/Out, Audio In/Out, and RS-232 Remote Control.




Figure 13–6. VCR Location



Recording Images

Video Cassette Recorder (VCR) (cont'd)




Figure 13-7. VCR Wiring Diagram

## Video Cassette Recorder (VCR) (cont'd)




Figure 13–8. LOGIQ™ 500 VCR Connections




Figure 13–9. VCR Connections

## Recording Images

---

### VCR and Color Page Printer

Both the Sony UP-1800 Color video Page Printer and Sony SVO-9500MD can be mounted behind the keyboard, under the monitor.



*NOTE: Wiring for the combination of peripherals is the same as for individual peripherals shown in previous sections.*




Figure 13–10. VCR and Color Page Printer Locations



## **Multi-Image Camera (MIC) (IIE Model 460)**

### **Operations**

The multi-image camera provides for the storage of still images on an 8½x11 sheet of x-ray film.

The Record keys can be programmed in the Set Up/System Parameters menu page 5 to control the print (expose) function of the camera. All other camera functions are adjusted by the controls mounted on the camera.

The LOGIQ™ 500 beeps after each exposure. The camera controls the placement of the images.

Some error messages could be displayed if a problem occurs in camera operation. Possible error messages are:

- Maximum number of prints have been made.
- No video signal. Check connections.
- Check Multi-Image camera.

For details on camera operation, consult the IIE Operator Manual provided with the camera.

## Recording Images

---

### Multi-Image Camera (MIC) (IIE Model 460) (cont'd)

#### Configurations

An International Imaging Electronics Multi-Image Camera can be attached to the LOGIQ™ 500. Contact a Sales or Service representative for details.

The Multi-Image Camera Connections are Composite Video and RS-232 Remote Control.




Figure 13–11. Multi-Image Camera Location



Multi-Image Camera (MIC) (IIE Model 460) (cont'd)




Figure 13–12. Multi-Image Camera Wiring Diagram



Recording Images

Multi-Image Camera (MIC) (IIE Model 460) (cont'd)




Figure 13-13. LOGIQ™ 500 Connections for Multi-Image Camera




Figure 13-14. Multi-Image Camera Connections



## Line Printer (USA ONLY)

### Operations

Display the desired report page:

- OB
- GYN
- OB Graph
- Cardiology
- Vascular
- Anatomical Survey

While the report page is displayed, press **Ctrl + R** to start the printing process. This Ctrl + R function is not available unless the report page is displayed on the screen.

The message “Preparing for print now” is displayed. After two seconds the message changes to “Newpt/Power Off will delete print data”. This message is displayed until the print operation is complete.



### Hints

If more than one report page is to be printed, the desired report page must be displayed and Ctrl + R pressed.

If printing was not completed due to some sort of error, the failed report must be displayed and Ctrl + R pressed again.

Ensure all reports have printed prior to deleting report data by pressing **New Patient** or power off.

### Configurations

The line printer will be connected to a serial port on the system. This can be accomplished at Port MIC, Port A or Port B on the rear panel.

The appropriate parameter in Set Up/System Parameters page 5 must be set to “Line Printer”.



## Recording Images

---

### Laser Camera

The LOGIQ™ 500 can print images to a Laser Camera for archival.

Contact the local GE Service Representative for details about hook-up and operation.

### Video Signal Specifications

Table 8–3 gives the video specifications of the Ultrasound system. This data may be needed when interfacing the system with OEM cameras or recording devices.

|                                    | NTSC           | PAL            |
|------------------------------------|----------------|----------------|
| Total Lines per Frame [line]       | 525            | 625            |
| Vertical Field Frequency [Hz]      | 60             | 50             |
| Horizontal Scanning Frequency [Hz] | 15.733         | 15.625         |
| Displayed Image Pixels [mm]        | 207.0 by 157.3 | 207.0 by 157.3 |
| Total Horizontal Line Time [ms]    | 63.56          | 64.00          |
| Horizontal Display [ms]            | 49.54          | 48.81          |
| Front Porch Width [ms]             | 2.76           | 3.09           |
| Sync Pulse Width [ms]              | 4.73           | 4.68           |
| Back Porch Width [ms]              | 6.53           | 7.42           |
| Total Horizontal Blanking [ms]     | 14.02          | 15.19          |
| Vertical Blanking Interval [H]     | 31.50          | 38.50          |
| Vertical Front Porch Width [H]     | 6.5            | 9.0            |
| Vertical Sync Width [H]            | 3              | 2.5            |
| Vertical Back Porch Width [H]      | 22             | 27             |

Table 8–3. Video Signal Specifications

### Maintenance

Consult the peripheral device operator manual for necessary routine maintenance.



## **MOD Image Archive (option)**

### **Overview**

Image Archive is an option to store and recall scan images using a 3.5 inch Magneto Optical Disk (128 or 230 Megabytes). The Magneto Optical Disk (MOD) allows for much faster and greater storage capacity than a Floppy Disk Drive (FDD).

Along with the increased speed and storage capacity, the user can perform measurements and calculations on images recalled from a MOD. Images recorded as hard copy (film) do not allow for additional measurements at a later date. Images recalled from a MOD have much better resolution than VCR playback images, providing better detail for additional measurements.

In addition to the image archive function, the MOD allows for a System Backup disk to be made in the unlikely event of a hard disk failure. System software updates and service diagnostics may also be accomplished much faster.

### **Archive Functions**

The Image Archive option allows the user to perform the following functions:

- Store Image to MOD
- Recall Image from MOD
- Patient File search of MOD
- MOD Media file search
- Delete a file from MOD
- Format the MOD disk
- Eject the MOD disk from the drive



## Recording Images

---

### Related Preset Parameters

Two selections must be made in the Set Up/Preset Program Sub-Menu page 2 that affect image archiving.

Image Archive Compression:

- **No** is no DEFF format image compression.
- **Yes** is for DEFF format image compression. Storage time for a compressed image is more than double.



*NOTE: If image compression is used, images saved by the LOGIQ™ 500 cannot be read on other TIFF readers, DEFF devices or the LOGIQ™ 700. However, they can be read on the LOGIQ™ 400.*

B/W Image Archive & Color Image Archive:

- **B/W** is for storing in black and white format (less storage space required).
- **Color** is for storing in color format (more storage space required).



*NOTE: An External Video Image is recognized as a color image and will use the Color Image Archive preset.*

### System ID Entry/Display

#### Overview

The LOGIQ™ 500 system ID number is assigned at the factory. A System ID should be unique for a system. This number is used by the system when initializing a MOD (image archive option) for image storage.

If there is more than one system at a facility, the ID numbers will be different. A scanner will write to a MOD if the media was initialized on a different scanner (different ID number). However, images can only be deleted on the system that initialized the MOD.

## Enter/Display System ID Number

To display or enter a System ID number:

- Select the **Set Up** Top Menu.
- Select **Utility** from the Set Up Sub Menu.




Figure 13–15. Utility Menu

- Enter number **10** (System ID Entry/Display) and press **Return**.
- Enter number **1** (ID Entry) or **2** (ID Display) and press **Return**.
- Enter an ID number in the range of **00000–16383** and press **Return**. If the ID was displayed, make a mental note or record the number for future reference.
- Press **Ctrl, R** simultaneously to exit the Utility menu and return to the previous scan mode.



*NOTE: Images which are not stored on the GE LOGIQ™ series systems will not be listed by the patient search function or displayed on a LOGIQ™ series machine.*



*NOTE: If there is more than one LOGIQ™ 400 or LOGIQ™ 500 in a facility or group of hospitals, each system should have it's own unique System ID number, unless it is necessary to delete MOD images on any system.*



## Recording Images

---

### Media Format (DEFF)

Display the Archive sub-menu.

| Archive        | DICOM       | Auto Seq    | CINE     |
|----------------|-------------|-------------|----------|
| Patient Search | Store Image |             | MO Eject |
| Media Search   |             | DEFF Format |          |

Figure 13–16. Archive Sub-Menu

Select **DEFF Format**. The following messages are displayed during the disk initialization and formatting process:

**“Do you continue ? (Y/N)”**

Insert the disk and press ‘Y’ to continue or ‘N’ to quit. If ‘Y’ is pressed,

**“Initializing now”**

appears. The disk format function is in progress.

The MO (disk) ID and System ID will be written on the disk media. One system can register a maximum of 10,000 disks.

If the number of registered disks is greater than 10,000, previously registered disk information will be overwritten in order to register the new disk information.

To avoid confusion, it is essential that care is taken not to attempt image storage to a previous (old) disk who’s MO ID and System ID is not available (IDs have been deleted from the media file).



*NOTE: The media list on the Media Search Menu contains both previous (old) and current (new) disk information.*





## **Disk Verification**

In order to protect images from being deleted by other LOGIQ™ 500 systems, the disk is always verified by the system before it can be used.

If the MO disk is not inserted into the drive or an unformatted MO disk is inserted, the following message is displayed:

“MEDIA IS EITHER UNMOUNTED OR UNFORMATTED.”

When a disk is inserted into a drive, the system first reads the MO (disk ID) and System ID from the media. These two IDs are displayed at the bottom of the monitor while the system verifies that they are registered in the verification list on the system hard drive.

## **Storing Images**

Display the Archive sub-menu. Select **Store Image** in the sub-menu. The following message is displayed during the storage process:

“In Progress. Please wait.”

System operation is interrupted for five seconds before scan operation can continue. The image storage operation is then performed in the “background”.

## **Image Comments**

With the comment cursor in the home position, four characters can be entered as a comment that is stored with the patient name, id and date of the scan. Refer to Figure 13–17 on 13–28 for examples.

## **Optional Method**

In addition to the Store Image function in the Soft Menu, the Record 1 & 2 keys can be programmed in the Set Up/System Parameter Sub-Menu page 5.

By selecting Image Archive for the Record1 B&W/Color or Record2 B&W/Color selection, these keys automatically store the B/W or Color image to MOD when pressed.



## Recording Images

### Patient Search

Select **Patient Search** from the Archive Sub-Menu. The Image Archive Patient Search Menu is displayed on the monitor as shown in Figure 13–17.

```
[IMAGE ARCHIVE PATIENT SEARCH MENU]

PT NAME      :
PT ID        :
DATE(YY/MM/DD) : / /

PAGE : 1 / 2

FILE  ID      NAME      DATE
----  -
AABH 98765    Johnson, Mary Jane  072897
AABG 98765    Johnson, Mary Jane  072897abcd
AABF 98765    Johnson, Mary Jane  072897
AABE 98765    Johnson, Mary Jane  072897test
AABD 98765    Johnson, Mary Jane  072897
AABC 23456    Doe, Jane B.         072897
AABB 23456    Doe, Jane B.         072897
AABA 23456    Doe, Jane B.         072897
AAAZ 23456    Doe, Jane B.         072897
AAAY 23456    Doe, Jane B.         072897
AAXX 12345    Smith, Mary A.       072897abcd
AAAW 12345    Smith, Mary A.       072897abcd
AAAV 12345    Smith, Mary A.       072897test
AAAU 12345    Smith, Mary A.       072897
AAPT 12345    Smith, Mary A.       072897

Search  Lock   List All  Delete
Tag All  Unlock Exit
Recall

MO : IMAGE
FREE : 112 MB

Ctrl+C:Cancel  Ellipse:Page Scroll
```

Figure 13–17. Patient Search Menu

|          |                                                                                                                                                            |
|----------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Comments | If four characters were typed at the comment home position prior to storing the image, the characters (comment) are displayed to the right after the date. |
| MO       | Shows the type of disk in the drive as an image disk.                                                                                                      |
| FREE     | Shows the amount of memory space available on the MOD in MegaBytes.                                                                                        |



## Recording Images

**Search Criteria** A specific MOD can be searched for patient files that meet certain criteria. The hard disk data base and backup MODs can also be searched for a list of MODs and files that meet the criteria listed below.

|              |                                               |
|--------------|-----------------------------------------------|
| Patient Name | Maximum 29 characters                         |
| Patient ID   | Maximum 14 characters                         |
| Date         | Format specified in System Parameters page 1. |

**Criteria Input Region** Both the Patient Search Menu and Media Search Menu have a criteria input region as shown in Figure 13–18.

|                |   |     |
|----------------|---|-----|
| PT NAME        | : |     |
| PT ID          | : |     |
| DATE(YY/MM/DD) | : | / / |

Figure 13–18. Criteria Input Region

Type the desired name in the highlighted region. Press **Return** or **Trackball** to the PT ID region.

Input the Patient ID. Press **Return** or **Trackball** to the Date region.

Input the desired date. Press **Return** or **Trackball** to move the cursor out of the criteria input region. The search function starts immediately.



### Hints

Spaces can be used as a wild card in the date entry. For example:

|         |                                                       |
|---------|-------------------------------------------------------|
| MM/ /   | lists everything for a specific month.                |
| / /YY   | lists all in a specific year.                         |
| MM/ /YY | lists all within a specific month of a specific year. |

If nothing is entered in one or more of the search criteria, a list of files is displayed that match the remaining criteria entered.



## Recording Images

---

**Searching the MOD** The search function is initiated when the **Trackball** is used to move the highlighted cursor out of the criteria area or the **Return** key is pressed while in the Date entry area.

The first 15 files are displayed on the monitor. If the search is not complete, the message:

**“In Progress, Please wait”**

appears at the bottom of the monitor. This message remains until the search process is totally complete.



*NOTE: Images which are not stored on the GE LOGIQ™ series systems will not be listed by the patient search function or displayed on a LOGIQ™ series machine.*



### Hints

In order to perform any other function, the user must wait until the search is complete or press **Ctrl, C** simultaneously to cancel the search function.

If no search criteria was entered by the user, a complete file listing of the entire MOD is executed. Only 15 files are displayed at one time on the monitor.

Use the top or bottom of the Ellipse rocker switch to scroll to the next or previous patient list page.



## Recording Images

---

### Menu Commands

The commands listed at the bottom of the Patient Search Menu can be used to manipulate the file information displayed.

|          |                                                                                                                                                       |
|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| Search   | Changes the current mode to the search criteria input mode.                                                                                           |
| Tag All  | Selects all images on the displayed page.                                                                                                             |
| Recall   | Recalls an image from the MOD for display on the monitor.                                                                                             |
| Lock     | Locks an image file.<br>A locked file has an asterisk (*) after its file number and cannot be deleted while locked.                                   |
| Unlock   | Unlocks an image file.                                                                                                                                |
| List All | Displays all images on the MOD in the Patient Search Menu.                                                                                            |
| Exit     | Exits the Patient Search Menu.                                                                                                                        |
| Delete   | Erases an unlocked image file from the list.<br>A file must be unlocked before being erased.<br><i>The image is not deleted from the affected MO.</i> |



**NOTE:** An image file will still be listed in the Media Search Menu even though the file has been deleted from the Patient Search Menu.

|         |                                                      |
|---------|------------------------------------------------------|
| Ellipse | Scrolls between the available pages to be displayed. |
| CTRL+C  | Cancels the search process.                          |



## Recording Images

---

### Image Selection for Recall

To select an image for recall from the Patient Search Menu:

- **Trackball** to highlight the desired image file and press **Set**. The image file data will be outlined with the box.
- **Trackball** to highlight the next image file and press **Set**.
- Repeat the procedure to “TAG” all of the desired image files to be recalled.



*NOTE: If all the files are required for recall, Trackball to TAG ALL and press **Set**.*

Once all the desired image files are tagged, **Trackball** to RECALL and press **Set**.

The first image is displayed and the Ellipse rocker switch lights up. The **Ellipse** rocker switch is used to display the next or previous image.

### Image Recall Process

During the Image Recall process, from the Patient Search Menu, the message:

**“In progress. Please wait.”**

is displayed; changing modes and inputting scan parameters are not allowed.

To quit the recall process, press the **Freeze** key. The system returns to the previous menu displayed.



### Hints

Storage of a recalled image can be accomplished as previously described. A new image is stored on the MOD and the original recalled image remains on the MOD unchanged.

## Media Search

The Image Archive option has the ability to store a file of all registered (Own–System) disks on the hard drive. This file consists of MO (disk) ID, Patient ID, Patient Name and Date with note information.

## Hard Disk Capacity

The capacity of the hard drive to store all media information is limited.

When the hard disk has reached its capacity to store media information, the user is asked to make a Backup MOD. The following message is displayed:

**“Too much data. Store failed.”**

The user should perform the Backup function at this time.

The Media Search function can search this data base or back-up MODs and indicate which Disk ID number contains information that meets the search criteria.

The media search menu is shown in Figure 13–19.

[IMAGE ARCHIVE MEDIA SEARCH MENU]

PT NAME :  
PT ID :  
DATE(Y Y/M M/D D) : / /

Search : PAGE : 1 / 3

| MO   | ID    | NAME               | DATE       |
|------|-------|--------------------|------------|
| 0000 | 98765 | Johnson, Mary Jane | 072897     |
| 0000 | 98765 | Johnson, Mary Jane | 072897     |
| 0000 | 98765 | Johnson, Mary Jane | 072897abcd |
| 0000 | 98765 | Johnson, Mary Jane | 072897     |
| 0000 | 98765 | Johnson, Mary Jane | 072897test |
| 0000 | 98765 | Johnson, Mary Jane | 072897     |
| 0001 | 56789 | White, July C.     | 072897abcd |
| 0001 | 56789 | White, July C.     | 072897test |
| 0001 | 56789 | White, July C.     | 072897rstv |
| 0001 | 56789 | White, July C.     | 072897     |
| 0001 | 34567 | Black, June K.     | 072897wxyz |
| 0001 | 34567 | Black, June K.     | 072897wxyz |
| 0001 | 34567 | Black, June K.     | 072897wxyz |
| 0001 | 34567 | Black, June K.     | 072897wxyz |
| 0001 | 34567 | Black, June K.     | 072897     |

Search(Ctrl+S) BackupSearch(Ctrl+B) Exit

Ctrl+C:Cancel Ellipse:Page Scroll

Figure 13–19. Media Search Menu



## Recording Images

---

### Menu Commands

The commands listed at the bottom of the Media Search Menu can be used to manipulate the file information displayed.

|                        |                                                                       |
|------------------------|-----------------------------------------------------------------------|
| Ellipse                | Displays the next page or previous page of files                      |
| CTRL+C                 | Cancels the search process                                            |
| Search (CTRL+S)        | Changes the current mode to the search criteria input mode            |
| Backup Search (CTRL+B) | Executes search of media backup MODs                                  |
| CLEAR                  | Clears (delete) entered characters in highlighted criteria input area |
| Exit                   | Exits the Media Search menu                                           |

### Messages

Messages displayed during the Media Search are similar to those displayed during Store Image, Image Recall and Patient Search. These messages are:

**“In progress. Please wait.”**

**“Search failed.”**

**“MEDIA IS EITHER UNMOUNTED OR UNFORMATTED.”**

The following message is displayed when the last page of the media list is displayed on the monitor:

**“There may be file in Backup Media fits this criteria.”**

The user should then search any Media backup disks for desired files.



*NOTE: When both NTSC and PAL images are stored on the same MO (disk), the following message is displayed:  
“Incompatible image are not listed up.”*

### Media File Recall

It is important to note that image files CANNOT be recalled from the media search menu. This function simply gives the user a list of all MO (disk) ID numbers, Patient Names, Patient IDs and Dates that meet the search criteria.





## MO Eject

To eject a disk from the drive, the user must select MO Eject from the Archive sub-menu.

**“Eject MO from the drive now.”**

is displayed on the monitor. The inserted MO is then ejected from the drive automatically.



*NOTE: After a MO is inserted into the drive and the search process is executed from the Patient Search Menu or the Media Search Menu, selecting MO EJECT from the soft menu is the only method to eject the MO. The MO cannot be ejected by the eject button on the front of the drive.*



## *Recording Images*

---

This page left blank intentionally.



# Time Adjustment

## Overview

The system clock for the LOGIQ™ 500 is set by service personnel to the local time during the installation process.

If any other time adjustment is needed during the year, the clock can be changed in the Set Up/Utility Menu.

## Time Adjustment

To adjust the system clock:

- Select the **Set Up** Top Menu.
- From the Set Up Sub-Menus, select Utility. The Utility Menu appears as shown in Figure 14–1.




Figure 14–1. Utility Menu



## Time Adjustment

---

### Time Adjustment (cont'd)

- Enter menu number “1” or “01” (Time Adjustment) in the Utility Menu and press **Return**. The Time Adjustment Menu appears as shown in Figure 14–2.

|                                                       |                   |
|-------------------------------------------------------|-------------------|
| HOSPITAL NAME                                         | 98/05/09 01:00:24 |
| Time Adjustment(1/1)                                  | S317              |
| [Caution] You can change only time within +,-12 hours |                   |
| Input the difference of time (HH:MM) ?                |                   |

Figure 14–2. Time Adjustment Menu



*NOTE: Time is changed by adjusting the current time a maximum of plus or minus 12 hours.*

- To Advance (add) the time, type in the number of hours and minutes to add to the currently displayed time. HH:MM is the format to use.
- To Retard (subtract) the time, type a minus sign (–) then the hours and minutes to be subtracted from the currently displayed time. –HH:MM is the format to use.
- Press **Return**. This completes the time adjustment and the Utility Menu reappears with the new time displayed.
- Press the Utility rocker switch to return to the previous scan mode.

The 12 hour display format has no am/pm indicator. The date will change at noon if the set time is 12 hours off.

# Preset Parameters

## Overview

LOGIQ™ 500 presets provide the user with a powerful tool to customize initial system operation for a particular exam type.

The organization structure of the presets in the system is as follows:

- Within each of the 7 exam categories to choose from in the patient entry menu, there are 8 user programmable and up to 8 factory applications presets.
- Parameters used in these presets are divided into 3 groups: Preset Program, Custom Display, and System Parameters.
- Changes to System Parameters affect all applications presets in all exam categories. These are the only parameters that affect factory presets. Nothing else can be changed in factory presets by the user.
- User programmable application presets can also be affected by parameters found in Preset Program and Custom Display. These parameters can customize the preset within an exam category or every user preset within an exam category.
- For additional flexibility, the first 9 pages of the Custom Display parameters are specific to the probe in use.

Decreased setup time and increased productivity can be accomplished by establishing an application preset for types of studies that are frequently performed by customizing all parameters to values that are normally used at the start of an exam.



## *Preset Parameters*

---

### **Overview** (cont'd)

Each parameter may be affected by a change in exam category, user preset, or probe. Each parameter description in this section is coded to show how it is affected. The codes are as follows:

|    |                                                   |
|----|---------------------------------------------------|
| EC | Exam Category dependent                           |
| UP | User Preset dependent                             |
| P  | Probe dependent                                   |
| R  | Reset to preset value with new patient selection. |



# Custom Display

## Overview

Custom Display in the Set Up Top Menu is 18 pages of preset parameters that are exam category/preset name dependent. The first 9 pages are also probe dependent within the exam category and preset name. The last 9 pages detail mode-specific items for the exam category/preset name.

| B              | Preset         | Set Up      | ECG         |
|----------------|----------------|-------------|-------------|
| Custom Display | Preset Program | Save Values | User Define |
| System Paramtr |                | Utility     | Diag.       |

Figure 14–3. Setup Sub-Menus

A brief description of each preset parameter and the choices available for that preset follows.



*NOTE: If it is necessary to quickly return to the Custom Display Sub-Menu page 1, press **Ctrl** and **1** simultaneously.*

*Refer to 14–6 for the key to the codes used to describe each parameter's affect.*



## Parameter Menu Command Lines

An example of a Custom Display parameter menu is shown in Figure 14–4.

```
03/28/99 14:37:38
GE Medical Systems C551
USER PARAM SYSTEM PARAM PRESET PROG CUSTOM DISP
CATEGORY : RAD/ABDOMEN PRESET NAME : ABDOMEN
PAGE : 17/18 PRIOR NEXT CONTENTS
COMMAND : SAVE RESET DELETE EXIT

***** Imaging Parameter 17 *****
ECG Auto Display Off On
ECG Gain -10
ECG Position -40
ECG Heart Rate Display Off On
ECG R Delay 1 [ms] 0
ECG R Delay 1 to 2 [ms] 300
ECG Audio Beat Sound Off On
ECG Audio Beat Sound Tone Tone-1 Tone-2 Tone-3
ECG Sync Mark Display Off On
PCG Auto Display Off On
PCG Gain -6
PCG Position -25
AUX Auto Display Off On
AUX Gain -6
AUX Position -20

'TRACKBALL' and 'SET' are available.
```

Figure 14–4. Custom Display Menu

The top protected lines show the date, time, hospital name and active probe. The next four lines provide information about the menu displayed plus commands to manipulate the menu pages and parameters.

### Command line 1

This line shows the parameter categories available with the current one highlighted. The categories listed are:

- User Param is not active at this time.
- System Param has 7 pages of programmed parameters.
- Preset Program has 10 pages of programmable parameters.
- Custom Display has 18 pages of programmable parameters.





### **Command line 2**

This line shows the exam category and application preset name.

### **Command line 3**

This line shows the menu page displayed (1/18 is page one of eighteen).

PRIOR will display the previous menu page. NEXT will display the following menu page. CONTENTS will show a summary of all Custom Display pages.

- Use the **Trackball** to place the arrow cursor on PRIOR or NEXT.
- Press **Set** to change menu pages.

### **Command line 4**


This line lists 4 commands that affect the parameters in the category displayed. They are:

|        |                                                                                                                           |
|--------|---------------------------------------------------------------------------------------------------------------------------|
| SAVE   | is used to save any preset parameter change(s).                                                                           |
| RESET  | if parameters were changed after the SAVE function was used, RESET will return those parameters to the saved values.      |
| DELETE | is used to erase the currently displayed user preset. Factory presets cannot be erased.                                   |
| EXIT   | returns the system to the last scan mode. Parameter changes will only be temporary if they were not saved before exiting. |



## Changing a Parameter Value

To change a parameter value:




**Trackball** the arrow cursor to the desired parameter.

Press **Set** to highlight the choice.

If the parameter lists several values, **Trackball** the arrow cursor to the desired value and press **Set**.

If the parameter needs the value to change, use the **Ellipse** rocker switch to increase or decrease the value.



When the desired value is selected or changed, **Trackball** the arrow cursor to SAVE.

Press **Set**.

A prompt at the bottom of the screen asks:

“Overwrite Existing Data? ‘y’ or ‘n’”.

Press ‘y’ to save the new values for the preset displayed.

Press ‘n’ and the prompt “Input User Preset Number (1–8)” appears. Enter the desired number.






*NOTE: If no is selected to overwrite current data, the system is asking to make a new user application preset. The eight available spaces are displayed by selecting the preset top menu.*



## Changing a Parameter Value (cont'd)

To customize the name of the new preset (1–8):

**Trackball** to the preset name in command line 2.



Press **Set**. Type in the personalized preset name.

**Trackball** to Save.

Press **Set**. All parameter data is saved as a new preset.



### Hints

- Parameters on all of the menu pages can be customized before selecting Save.
- Parameter changes will be lost if Save is not selected.
- A maximum of 8 user exam application presets can be saved.
- Only a change in System Parameters will affect factory presets.

## Custom Display Contents




Figure 14–5. Custom Display Contents Page

The custom display contents page shows a summary of the 18 menu pages.

Note that the first 9 pages of presets are probe dependent as well as category/preset dependent like the last 9 pages.

The 3D Mode menus (listed on pages 9 and 18) are available with the option installed.

**Trackball** the cursor to the desired page and press **Set**.



**Page 1 of 18 (Imaging Parameter 1 – Probe Dependent 1)**

|                                                                                                                                                                                                                                             |                     |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|
| <b>Enter Depth</b>                                                                                                                                                                                                                          | <b>EC, UP, P, R</b> |
| Choose the desired default display depth for the probe indicated.<br>4, 5, 6, 7, 8, 10, 12, 14, 16, 18, 20, 22 or 24 cm.                                                                                                                    |                     |
| <b>Focus Depth [mm]</b>                                                                                                                                                                                                                     | <b>EC, UP, P, R</b> |
| Choose the desired default focus depth for the probe indicated.<br>Depth from 1 to 240 mm in 1 mm increments.                                                                                                                               |                     |
| <b>Image Rotation</b>                                                                                                                                                                                                                       | <b>EC, UP, P, R</b> |
| Choose the desired default image orientation.<br><br>0 (N)      Normal top to bottom orientation.<br>90 (W)    90° left to right orientation.<br>180 (S)   Inverted bottom to top orientation.<br>270 (E)   270° right to left orientation. |                     |
| <b>Image Reverse</b>                                                                                                                                                                                                                        | <b>EC, UP, P, R</b> |
| Choose the initial default setting of image reverse (left-right orientation).<br><br>On or Off.                                                                                                                                             |                     |
| <b>Zoom Auto Start</b>                                                                                                                                                                                                                      | <b>EC, UP, P</b>    |
| Choose to automatically place the image into the zoom mode after a probe selection change.<br><br>Off      No auto zoom at probe change.<br>On      Auto zoom mode at probe change.                                                         |                     |
| <b>Zoom Reference</b>                                                                                                                                                                                                                       | <b>EC, UP, P</b>    |
| Choose the zoom reference image default setting.<br><br>On or Off.                                                                                                                                                                          |                     |
| <b>Zoom Factor</b>                                                                                                                                                                                                                          | <b>EC, UP, P, R</b> |
| Choose the default zoom factor for the read (freeze) zoom.<br><br>1.0, 1.2, 1.5, 2.0, 2.5, 3.0 or 4.0.                                                                                                                                      |                     |



**Page 1 of 18 (Imaging Parameter 1 – Probe Dependent 1) (cont'd)**

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                                  |                                      |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------|--------------------------------------|
| <b>Standoff Setting EC, UP, P</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |                                  |                                      |
| Choose a setting for the type of probe standoff, if used, for the 739L or L764 (LH) probes. When any other probe is selected the preset is "Not Available".                                                                                                                                                                                                                                                                                                                                                                                  |                                  |                                      |
| Off (No standoff), WP (Water path) or CPL (Couplant)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                                  |                                      |
| <b>Needle Guide Type EC, UP, P</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                  |                                      |
| Choose the type of biopsy guidezone to be displayed for the guide attachment used. For most probes, choose either Single Fixed Angle (SGL) or Multi Angle (MBX). For the E721/MTZ, the probe guidezone choices are:                                                                                                                                                                                                                                                                                                                          |                                  |                                      |
| TV0° Reusable metal guide with a 0 degree offset angle.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |                                  |                                      |
| TR5° Civco disposable guide with a 5 degree offset angle.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                  |                                      |
| <b>Guide Attachment</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | <b>Needle Guide Type Setting</b> | <b>Biopsy Zone Soft Menu Choices</b> |
| Multi Angle Attachment                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | MBX                              | MBX1, MBX2, MBX3, OFF                |
| Fixed Angle Attachment                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | SGL                              | SGL, OFF                             |
| E721/MTZ Metal Guide                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | TV0°                             | TV0°, OFF                            |
| E721/MTZ Disposable Guide                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | TR5°                             | TR5°, OFF                            |
| ERB7 Disposable Guide                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | BX                               | 1~8, OFF                             |
| ERB7 Needle Placement Guide                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | GBX                              | GBX, OFF                             |
| <div><div><b>DANGER</b></div><div>Failure to match the guidezone displayed to the guide may cause the needle to track a path outside the displayed zone.</div><div>It is extremely important that when using the adjustable angle biopsy guides, the angle displayed on the screen matches the angle set on the guide, otherwise the needle will not follow the displayed guidezone which could result in repeated biopsies or patient injury.</div></div> |                                  |                                      |
| <b>Acoustic Power B/M [%] EC, UP, P, R</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |                                  |                                      |
| Set the initial percentage of acoustic power to be used when in B- or B/M-Modes.                                                                                                                                                                                                                                                                                                                                                                                                                                                             |                                  |                                      |
| Percentage from 0 to 100% in 10% increments.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                                  |                                      |
| <b>Acoustic Power D/CFM/PDI [%] EC, UP, P, R</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |                                  |                                      |
| Set the initial percentage of acoustic power to be used when in Doppler, Color Flow or Power Doppler imaging modes.                                                                                                                                                                                                                                                                                                                                                                                                                          |                                  |                                      |
| Percentage from 0 to 100% in 10% increments.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                                  |                                      |



Page 2 of 18 (Imaging Parameter 2 – Probe Dependent 2)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                     |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|
| <b>B Gain</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | <b>EC, UP, P, R</b> |
| Choose the default gain value for B-Mode.<br><br>Value from 0 to 98 in 2 digit increments. The maximum value is probe dependent and may not reach 98.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                     |
| <b>B Dynamic Range</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | <b>EC, UP, P, R</b> |
| Choose the default dynamic range value for B-Mode.<br><br>Value from 30 to 90 in 6 digit increments.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |                     |
| <b>B Edge Enhance</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | <b>EC, UP, P, R</b> |
| Choose the default edge enhancement value for B-Mode.<br><br>Off, Low, Mid or High.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                     |
| <b>B Image Softener</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | <b>EC, UP, P, R</b> |
| Choose the B-Mode image softener function default.<br><br>On or Off.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |                     |
| <b>B Softener Level</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | <b>EC, UP, P, R</b> |
| Choose the B-Mode image softener default level.<br><br>D1, D2, D3, Low, Mid or High.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |                     |
| <b>B Frame Average</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | <b>EC, UP, P</b>    |
| Choose the B-Mode frame averaging default value.<br><br>Off (no frame averaging), 1, 2, 3, 4, 5 or 6.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                     |
| <b>B Image Frequency</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | <b>EC, UP, P, R</b> |
| Choose the B-Mode image frequency optimization setting.<br><br>Resolution                      Process image for best resolution.<br>Normal                          Compromise between resolution and penetration.<br>Penetration                      Process image for penetration.<br>Further Penetration <b>C358 Probe Only</b> —Process image for more penetration<br>Harmonic1 <b>S222, S317 and C358 Probes Only</b> —Process image with<br>tissue harmonic frequency.<br>Harmonic2 <b>S222, S317 and C358 Probes Only</b> —Process image with<br>tissue harmonic frequency.<br><br>The soft menu display shows this parameter as an imaging frequency depending on the active probe. |                     |



Page 2 of 18 (Imaging Parameter 2 – Probe Dependent 2) (cont'd)

B Gray Scale Map    EC, UP, P, R

Choose the B-Mode gray scale default map.  
1, 2, 3, 4, 5, 6, 7, 8, 9,10, 11, 12, 13, 14, 15 or 16.




Figure 14–6. Gray Scale Map Graphs





Page 2 of 18 (Imaging Parameter 2 – Probe Dependent 2) (cont'd)

|                                                                                                                                                                                                                                                                                                                                                                      |                               |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------|
| <b>B Rejection</b>                                                                                                                                                                                                                                                                                                                                                   | <b>EC, UP, P, R</b>           |
| Choose the B-Mode rejection default level.<br>Level from 0 to 40 in 1 digit increments.                                                                                                                                                                                                                                                                              |                               |
| <b>B Focus Combi Number</b>                                                                                                                                                                                                                                                                                                                                          | <b>EC, UP, P, R</b>           |
| Choose the default number of focus points in B-Mode.<br>1, 2, 3 or 4.<br>If Range is selected as B Combi Focus Width, the maximum number is 3.                                                                                                                                                                                                                       |                               |
| <b>B Combi Focus Width</b>                                                                                                                                                                                                                                                                                                                                           | <b>EC, UP, P</b>              |
| Choose the default width between focus points in B-Mode.<br><br>Narrow    Focal points close together.<br>Normal    Focal point at nominal distance.<br>Wide       Focal points far apart.<br>Range     Improves near field resolution as well as that of the transmit focal zone set by 'Focus Position'. Maximum number of displayed, moveable focus markers is 3. |                               |
| <b>B Angle/Width</b>                                                                                                                                                                                                                                                                                                                                                 | <b>[deg, mm] EC, UP, P, R</b> |
| Choose the amount of the scan angle to display as a default for B-Mode.<br>The amount or degrees available depends on the probe.                                                                                                                                                                                                                                     |                               |
| <b>B Target Frame Rate</b>                                                                                                                                                                                                                                                                                                                                           | <b>EC, UP, P</b>              |
| Choose the default value for B-Mode frame rate.<br><br>Low        Low frame rate, better resolution.<br>Mid        Mid range frame rate and resolution.<br>High       Higher frame rate, less resolution.                                                                                                                                                            |                               |



## Page 3 of 18 (Imaging Parameter 3 – Probe Dependent 3)

|                                                                                           |                           |
|-------------------------------------------------------------------------------------------|---------------------------|
| <b>Timeline Sweep Speed</b>                                                               | <b>EC, UP, P, R</b>       |
| Choose the default value for the timeline (M-Mode and Doppler spectrum) sweep speed.      |                           |
| Slow                                                                                      | 16 microseconds per pixel |
| Mid                                                                                       | 8 microseconds per pixel  |
| Fast                                                                                      | 4 microseconds per pixel  |
| Very-Fast*                                                                                | 2 microseconds per pixel  |
| * Works with Doppler only, not M-Mode or M/D-Modes.                                       |                           |
| <b>Physio Sweep Speed on B</b>                                                            | <b>EC, UP, P, R</b>       |
| Choose the default value for the physio signal sweep speed displayed on the B-Mode image. |                           |
| Slow                                                                                      | 16 microseconds per pixel |
| Mid                                                                                       | 8 microseconds per pixel  |
| Fast                                                                                      | 4 microseconds per pixel  |
| <b>M Gain (Delta from B)</b>                                                              | <b>EC, UP, P, R</b>       |
| Choose the default value for M-Mode gain difference from the B-Mode Gain.                 |                           |
| Value from –98 to 98 in 2 digit increments.                                               |                           |
| <b>M Dynamic Range</b>                                                                    | <b>EC, UP, P, R</b>       |
| Choose the M-Mode dynamic range default value.                                            |                           |
| Value from 30 to 90 in 6 digit increments.                                                |                           |
| <b>M Edge Enhance</b>                                                                     | <b>EC, UP, P, R</b>       |
| Choose the M-Mode edge enhancement default level.                                         |                           |
| Off, Low, Mid or High.                                                                    |                           |
| <b>M Image Softener</b>                                                                   | <b>EC, UP, P, R</b>       |
| Choose the M-Mode image softener function default.                                        |                           |
| On or Off.                                                                                |                           |
| <b>M Softener Level</b>                                                                   | <b>EC, UP, P, R</b>       |
| Choose the M-Mode image softener default level.                                           |                           |
| Low, Mid or High.                                                                         |                           |



Page 3 of 18 (Imaging Parameter 3 – Probe Dependent 3)

|                                                         |                                                                                       |
|---------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>M Image Frequency</b> EC, UP, P, R                   |                                                                                       |
| Choose the M-Mode image frequency optimization setting. |                                                                                       |
| Resolution                                              | Process image for best resolution.                                                    |
| Normal                                                  | Compromise between resolution and penetration.                                        |
| Penetration                                             | Process image for penetration.                                                        |
| Further Penetration                                     | <b>C358 Probe Only</b> —Process image for more penetration.                           |
| Harmonic1                                               | <b>S222, S317 and C358 Probes Only</b> —Process image with tissue harmonic frequency. |
| Harmonic2                                               | <b>S222, S317 and C358 Probes Only</b> —Process image with tissue harmonic frequency. |
| This parameter is not selectable from the soft menu.    |                                                                                       |
| <b>M Gray Scale Map</b> EC, UP, P, R                    |                                                                                       |
| Choose a M-Mode gray scale default map.                 |                                                                                       |
| 1, 2, 3, 4, 5, 6, 7, 8, 9 or 10.                        |                                                                                       |
| <b>M Rejection</b> EC, UP, P, R                         |                                                                                       |
| Choose the M-Mode rejection default level.              |                                                                                       |
| Level from 0 to 40 in 2 digit increments.               |                                                                                       |
| <b>M Peak Hold</b> EC, UP, P, R                         |                                                                                       |
| Choose the default M-Mode peak hold function status.    |                                                                                       |
| On or Off.                                              |                                                                                       |



## Page 4 of 18 (Imaging Parameter 4 – Probe Dependent 4)

|                                                                                                                                                                |                                              |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------|
| <b>PWD/CFM/PDI Penetration</b>                                                                                                                                 | <b>EC, UP, P, R</b>                          |
| Choose to have the default for penetration, using the lower probe frequency, for pulsed Doppler, color flow mapping and power Doppler imaging modes on or off. |                                              |
| Off                                                                                                                                                            | Use higher probe frequency.                  |
| On                                                                                                                                                             | Use lower probe frequency.                   |
| <b>PWD/CFM/PDI Slant Auto Start</b>                                                                                                                            | <b>EC, UP, P, R</b>                          |
| Choose the default position of the slant scan window.                                                                                                          |                                              |
| Off                                                                                                                                                            | Vertical (no slant)                          |
| Other-Side                                                                                                                                                     | Slanted away from the probe orientation mark |
| Marker-Side                                                                                                                                                    | Slanted toward the probe orientation mark    |
| <b>PWD Gain</b>                                                                                                                                                | <b>EC, UP, P, R</b>                          |
| Choose the default gain value for pulsed Doppler Mode.                                                                                                         |                                              |
| Value from 0 to 32 in 2 digit increments.                                                                                                                      |                                              |
| <b>PWD Dynamic Range</b>                                                                                                                                       | <b>EC, UP, P, R</b>                          |
| Choose the default dynamic range value for pulsed Doppler Mode.                                                                                                |                                              |
| Value from 18 to 48 in 6 digit increments.                                                                                                                     |                                              |
| <b>PWD Sample Volume [mm]</b>                                                                                                                                  | <b>EC, UP, P, R</b>                          |
| Choose the default pulsed Doppler sample volume cursor size (in millimeters).                                                                                  |                                              |
| 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14 or 16.                                                                                                                   |                                              |
| <b>PWD Gray Scale Map</b>                                                                                                                                      | <b>EC, UP, P, R</b>                          |
| Choose the pulsed Doppler Mode gray scale default map.                                                                                                         |                                              |
| 1, 2, 3, 4, 5, 6, 7, 8, 9 or 10.                                                                                                                               |                                              |
| <b>PWD Rejection</b>                                                                                                                                           | <b>EC, UP, P, R</b>                          |
| Choose the pulsed Doppler Mode rejection default level.                                                                                                        |                                              |
| Level from 0 to 40 in 2 digit increments.                                                                                                                      |                                              |
| <b>PWD Velocity [cm/s]</b>                                                                                                                                     | <b>EC, UP, P, R</b>                          |
| Choose the pulsed Doppler Mode velocity scale default value.                                                                                                   |                                              |
| Value from 10 to 1000 in 10 digit increments.                                                                                                                  |                                              |
| <b>PWD Wall Filter [cm/s]</b>                                                                                                                                  | <b>EC, UP, P, R</b>                          |
| Choose the pulsed Doppler Mode wall filter default value.                                                                                                      |                                              |
| 0.3, 0.5, 0.7, 1.0, 1.5, 2.0, 3.0, 5.0, 7.0, 10, 15, 20, 30 or 50.                                                                                             |                                              |

Page 4 of 18 (Imaging Parameter 4 – Probe Dependent 4) (cont'd)

|                                                                                                                                          |                                                                                              |
|------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|
| <b>PWD Alternative Scan</b> <b>EC, UP, P</b>                                                                                             |                                                                                              |
| Choose the action of the B-Pause key while in PWD Mode.                                                                                  |                                                                                              |
| On                                                                                                                                       | Pressing B-Pause toggles between B-Mode frozen/PWD real-time or B-Mode real-time/PWD frozen. |
| Off                                                                                                                                      | Pressing B-Pause freezes or unfreezes the B-Mode image only.                                 |
| <b>PWD/CWD Flow Angle [deg]</b> <b>EC, UP, P, R</b>                                                                                      |                                                                                              |
| Choose the default for the Doppler flow $\theta$ angle.                                                                                  |                                                                                              |
| Value from -80 to +80 in 5° increments.                                                                                                  |                                                                                              |
| <b>CFM/PWD Ratio in Triplex</b> <b>EC, UP, P, R</b>                                                                                      |                                                                                              |
| Choose the velocity ratio between the CFM display and Pulsed Doppler.                                                                    |                                                                                              |
| 1/1                                                                                                                                      | Same velocity display.                                                                       |
| 1/2                                                                                                                                      | Pulsed Doppler velocity 2 times CFM.                                                         |
| 1/4                                                                                                                                      | Pulsed Doppler velocity 4 times CFM.                                                         |
| <b>PDI/PWD Ratio in Triplex</b> <b>EC, UP, P, R</b>                                                                                      |                                                                                              |
| Choose the velocity ratio between the PDI display and Pulsed Doppler.                                                                    |                                                                                              |
| 1/1                                                                                                                                      | Same velocity display.                                                                       |
| 1/2                                                                                                                                      | Pulsed Doppler velocity 2 times PDI.                                                         |
| 1/4                                                                                                                                      | Pulsed Doppler velocity 4 times PDI.                                                         |
| <b>CFM/PDI Shrink in Triplex</b> <b>EC, UP, P, R</b>                                                                                     |                                                                                              |
| Choose to turn CFM/PDI shrink function on or off. Allows for shrinking of the CFM/PDI window to a specified value to improve frame rate. |                                                                                              |
| On or Off.                                                                                                                               |                                                                                              |
| <b>CFM/PDI Shrink Angle/Width [deg, mm]</b> <b>EC, UP, P, R</b>                                                                          |                                                                                              |
| Choose the shrink value for the CFM/PDI window in degrees for convex and sector probes or millimeters for linear probes.                 |                                                                                              |
| Value range depends on probe type.                                                                                                       |                                                                                              |



## Page 5 of 18 (Imaging Parameter 5 – Probe Dependent 5)

|                                                                                                                                          |                            |
|------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|
| <b>CWD Gain</b>                                                                                                                          | <b>EC, UP, P, R</b>        |
| Choose the default gain value for continuous wave Doppler Mode.<br>Value from 0 to 32 in 2 digit increments.                             |                            |
| <b>CWD Dynamic Range</b>                                                                                                                 | <b>EC, UP, P, R</b>        |
| Choose the continuous wave Doppler Mode dynamic range default value.<br>Value from 18 to 48 in 6 digit increments.                       |                            |
| <b>CWD Gray Scale Map</b>                                                                                                                | <b>EC, UP, P, R</b>        |
| Choose the continuous wave Doppler Mode gray scale default map.<br>1, 2, 3, 4, 5, 6, 7, 8, 9 or 10.                                      |                            |
| <b>CWD Rejection</b>                                                                                                                     | <b>EC, UP, P, R</b>        |
| Choose the continuous wave Doppler Mode rejection default level.<br>Level from 0 to 40 in 2 digit increments.                            |                            |
| <b>CWD Velocity</b>                                                                                                                      | <b>[cm/s] EC, UP, P, R</b> |
| Choose the continuous wave Doppler Mode velocity scale default value.<br>Value from 10 to 1000 in 10 digit increments.                   |                            |
| <b>CWD Wall Filter</b>                                                                                                                   | <b>[cm/s] EC, UP, P, R</b> |
| Choose the continuous wave Doppler Mode wall filter default value.<br>0.3, 0.5, 0.7, 1.0, 1.5, 2.0, 3.0, 5.0, 7.0, 10, 15, 20, 30 or 50. |                            |

## Page 6 of 18 (Imaging Parameter 6 – Probe Dependent 6)

|                                                                                                                                                                                                           |                     |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|
| <b>CFM Color Threshold</b>                                                                                                                                                                                | <b>EC, UP, P, R</b> |
| Choose a priority level for B-Mode display during CFM mode. The higher the number, the greater the emphasis on CFM Mode (color) versus B-Mode (gray scale).<br>Level from 0 to 100 in 1 digit increments. |                     |
| <b>CFM Threshold Velocity</b>                                                                                                                                                                             | <b>EC, UP, P, R</b> |
| Choose the level above which color will no longer be used to display velocity.<br>Level from 0 to 25 in 1 digit increments.                                                                               |                     |
| <b>CFM Threshold Turbulence</b>                                                                                                                                                                           | <b>EC, UP, P, R</b> |
| Choose the level above which color will no longer be used to display turbulence.<br>Level from 0 to 100 in 1 digit increments.                                                                            |                     |
| <b>CFM Wall Echo Cancellor</b>                                                                                                                                                                            | <b>EC, UP, P, R</b> |
| Choose the initial default setting of the wall echo cancel function.<br>On or Off.                                                                                                                        |                     |

Page 6 of 18 (Imaging Parameter 6 – Probe Dependent 6) (cont'd)

|                                                                                                                                                      |                                          |
|------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------|
| <b>CFM/PDI Spatial Filter</b> <b>EC, UP, P</b>                                                                                                       |                                          |
| Choose the desired amount of averaging between color pixels.                                                                                         |                                          |
| Off                                                                                                                                                  | No spatial filtering.                    |
| Low                                                                                                                                                  | Low spatial filtering.                   |
| Mid                                                                                                                                                  | Medium spatial filtering.                |
| High                                                                                                                                                 | Maximum spatial filtering.               |
| <b>CFM Target Frame Rate</b> <b>EC, UP, P</b>                                                                                                        |                                          |
| Choose the desired frame rate range while in CFM mode.                                                                                               |                                          |
| Density                                                                                                                                              | Lowest frame rate for high color density |
| Low                                                                                                                                                  | Low frame rate for CFM.                  |
| Mid                                                                                                                                                  | Medium frame rate for CFM.               |
| High                                                                                                                                                 | High frame rate for low color density.   |
| <b>CFM Target Frame Rate Dtl</b> <b>EC, UP, P</b>                                                                                                    |                                          |
| Choose the desired frame rate range while in the Detail (Dtl) Mode of CFM presentation.                                                              |                                          |
| Density                                                                                                                                              | Lowest frame rate for high color density |
| Low                                                                                                                                                  | Low frame rate for CFM Detail Mode.      |
| Mid                                                                                                                                                  | Medium frame rate for CFM Detail Mode.   |
| High                                                                                                                                                 | High frame rate for CFM Detail Mode.     |
| <b>CFM Packet Size Survey, Map, SrvyDtl and MpDtl</b> <b>EC, UP, P</b>                                                                               |                                          |
| Choose the default packet size for each CFM diagnostic mode display. This will be for Survey Mode, Map Mode, Survey-Detail Mode and Map-Detail Mode. |                                          |
| Small                                                                                                                                                | Small packet size.                       |
| Mid                                                                                                                                                  | Medium packet size.                      |
| Large                                                                                                                                                | Large packet size.                       |
| <i>NOTE: The higher the packet size, the better the color presentation but with a slower frame rate.</i>                                             |                                          |
| <b>CFM Velocity [cm/s]</b> <b>EC, UP, P, R</b>                                                                                                       |                                          |
| Choose the default CFM velocity scale in centimeters per second.                                                                                     |                                          |
| Value from 0 to 1000 in 10 digit increments.                                                                                                         |                                          |
| <b>CFM B Gain</b> <b>EC, UP, P, R</b>                                                                                                                |                                          |
| Choose the default B-Mode gain setting while in CFM Mode.                                                                                            |                                          |
| Value from 0 to 31 in 1 digit increments.                                                                                                            |                                          |



**Page 6 of 18 (Imaging Parameter 6 – Probe Dependent 6) (cont'd)**

|                                                                                               |                            |
|-----------------------------------------------------------------------------------------------|----------------------------|
| <b>CFM B High Resolution    EC, UP, P, R</b>                                                  |                            |
| Choose the initial default setting of the high resolution function in CFM Mode.<br>On or Off. |                            |
| <b>CFM B Frame Average    EC, UP, P, R</b>                                                    |                            |
| Choose the level of frame averaging for CFM while in combination B-Mode with CFM.             |                            |
| Off                                                                                           | No B-Mode frame averaging. |
| Low                                                                                           | Low frame averaging.       |
| Mid                                                                                           | Medium frame averaging.    |
| High                                                                                          | High frame averaging.      |
| <b>CFM B Wall Filter    EC, UP, P, R</b>                                                      |                            |
| Choose the default level of wall filtering for B-Mode while in CFM.                           |                            |
| Low, M1, M2, M3, M4 and High.                                                                 |                            |
| <b>CFM B Angle/Width    [deg, mm]    EC, UP, P, R</b>                                         |                            |
| Choose the default B-Mode angle or width of the display while in CFM.                         |                            |
| Value from 1 to 124 in one digit increments. Maximum is determined by the probe type.         |                            |
| <b>CFM/PDI Start Depth    [mm]    EC, UP, P, R</b>                                            |                            |
| Set the start depth of the CFM/PDI window area.                                               |                            |
| Depth in millimeters ranging from 0 to 230.                                                   |                            |
| <b>CFM/PDI Vertical Size    [mm]    EC, UP, P, R</b>                                          |                            |
| Set the vertical size of the CFM/PDI window.                                                  |                            |
| Size in millimeters ranging from 10 to 240.                                                   |                            |
| <b>CFM/PDI Focus Depth    [%]    EC, UP, P, R</b>                                             |                            |
| Set the desired depth of the focal zone in CFM/PDI Mode window.                               |                            |
| Percentage ranging from 0 to 100.                                                             |                            |



Page 7 of 18 (Imaging Parameter 7 – Probe Dependent 7)

|                                                                                                                                                                                                          |                     |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|
| <b>CFM M Gain (Delta from B)</b>                                                                                                                                                                         | <b>EC, UP, P, R</b> |
| Choose the default level for M-Mode gain while in CFM.<br>Level from –31 to 31 in 1 digit increments.                                                                                                    |                     |
| <b>CFM M Wall Filter</b>                                                                                                                                                                                 | <b>EC, UP, P, R</b> |
| Choose the default level for the M-Mode Wall Filter while in CFM Mode.<br>Low, M1, M2, M3, M4, or High.                                                                                                  |                     |
| <b>CFM ACE</b>                                                                                                                                                                                           | <b>EC, UP, P, R</b> |
| Choose the default setting for Adaptive Color Enhancement, which reduces clutter noise or color flash.<br>On or Off.                                                                                     |                     |
| <b>CFM Noise Blanker</b>                                                                                                                                                                                 | <b>EC, UP, P, R</b> |
| Choose the default setting for Noise Blanker, which reduces random color noise.<br>OFF, LOW, MID or HIGH.                                                                                                |                     |
| <b>CFM Persistence</b>                                                                                                                                                                                   | <b>EC, UP, P, R</b> |
| Choose the default setting for Persistence, which retains highest color pixel value for a specified period of time or until a higher value is detected.<br>OFF, SHORT, MID or HIGH.                      |                     |
| <b>CFM MR-Flow</b>                                                                                                                                                                                       | <b>EC, UP, P, R</b> |
| Choose the default setting of the MR-Flow option (combination focus in CFM area).<br>Limited to LA39, C551 (CAE), C358, E721, 739L and 546L probes only. Must have ACE-2 option installed.<br>On or Off. |                     |
| <b>CFM MR-Flow Combi Number</b>                                                                                                                                                                          | <b>EC, UP, P</b>    |
| Set the number of focal zones to be used for the combination focus for MR-Flow option.<br>2 or 3.                                                                                                        |                     |



## Page 8 of 18 (Imaging Parameter 8 – Probe Dependent 8)

|                                                                                                                                                                                                                                                                                                                                                                                                           |  |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| <b>PDI Color Threshold</b> <b>EC, UP, P, R</b>                                                                                                                                                                                                                                                                                                                                                            |  |
| Choose a priority level for B-Mode display during PDI mode. The higher the number, the greater the emphasis on B-Mode versus PDI Mode.<br>Choose 0 to 100 in 1 digit increments.                                                                                                                                                                                                                          |  |
| <b>PDI Threshold Power</b> <b>EC, UP, P, R</b>                                                                                                                                                                                                                                                                                                                                                            |  |
| Set the threshold value (%) at which power information will start to be displayed.<br>Set from 0 to 50 percent.<br>Example: Setting = 20%; Power values less than 20% will not be displayed. They will be set to the background color. Only power values over 20% will be displayed as color.                                                                                                             |  |
| <b>PDI Target Frame Rate</b> <b>EC, UP, P</b>                                                                                                                                                                                                                                                                                                                                                             |  |
| Choose the desired frame rate range while in PDI mode.<br><br>Density              Lowest frame rate for high color density<br>Low                  Low frame rate for PDI.<br>Mid                  Medium frame rate for PDI.<br>High                 High frame rate for PDI.                                                                                                                           |  |
| <b>PDI Target Frame Rate Dtl</b> <b>EC, UP, P</b>                                                                                                                                                                                                                                                                                                                                                         |  |
| Choose the desired frame rate range while in the Detail (Dtl) Mode of PDI presentation.<br><br>Density              Lowest frame rate for high color density<br>Low                  Low frame rate for PDI Detail Mode.<br>Mid                  Medium frame rate for PDI Detail Mode.<br>High                 High frame rate for PDI Detail Mode.                                                      |  |
| <b>PDI Packet Size Survey, Map, SrvyDtl and MpDtl</b> <b>EC, UP, P</b>                                                                                                                                                                                                                                                                                                                                    |  |
| Choose the default packet size for each PDI diagnostic mode display. This will be for Survey Mode, Map Mode, Survey-Detail Mode and Map-Detail Mode.<br><br>Small                Small packet size.<br>Mid                  Medium packet size.<br>Large                Large packet size.<br><br><i>NOTE: The higher the packet size the better the color presentation but with a slower frame rate.</i> |  |
| <b>PDI Velocity [cm/s]</b> <b>EC, UP, P, R</b>                                                                                                                                                                                                                                                                                                                                                            |  |
| Choose the default PDI velocity scale in centimeters per second.<br>Value from 0 to 1000 in 10 digit increments.                                                                                                                                                                                                                                                                                          |  |

Page 8 of 18 (Imaging Parameter 8 – Probe Dependent 8) (cont'd)

|                                                                                                                      |                            |
|----------------------------------------------------------------------------------------------------------------------|----------------------------|
| <b>PDI B Gain</b>                                                                                                    | <b>EC, UP, P, R</b>        |
| Choose the default B-Mode gain setting while in PDI Mode.<br>Value from 0 to 31 in 1 digit increments.               |                            |
| <b>PDI B High Resolution</b>                                                                                         | <b>EC, UP, P, R</b>        |
| Choose the initial default setting of the high resolution function in PDI Mode.<br>On or Off.                        |                            |
| <b>PDI B Frame Average</b>                                                                                           | <b>EC, UP, P, R</b>        |
| Choose the level of frame averaging for B-Mode while in PDI.                                                         |                            |
| Off                                                                                                                  | No B-Mode frame averaging. |
| Low                                                                                                                  | Low frame averaging.       |
| Mid                                                                                                                  | Medium frame averaging.    |
| High                                                                                                                 | High frame averaging.      |
| <b>PDI B Wall Filter</b>                                                                                             | <b>EC, UP, P, R</b>        |
| Choose the default level of wall filtering for B-Mode while in PDI.<br>Low, M1, M2, M3, M4 and High.                 |                            |
| <b>PDI M Gain (Delta from B)</b>                                                                                     | <b>EC, UP, P, R</b>        |
| Choose the default level for M-Mode gain while in PDI.<br>Level from –31 to 31 in 1 digit increments.                |                            |
| <b>PDI M Wall Filter</b>                                                                                             | <b>EC, UP, P, R</b>        |
| Choose the default level for the M-Mode Wall Filter while in PDI Mode.<br>Low, M1, M2, M3, M4, or High.              |                            |
| <b>PDI ACE</b>                                                                                                       | <b>EC, UP, P, R</b>        |
| Choose the default setting for Adaptive Color Enhancement, which reduces clutter noise or color flash.<br>ON or OFF. |                            |
| <b>PDI Noise Blanker</b>                                                                                             | <b>EC, UP, P, R</b>        |
| Choose the default setting for Noise Blanker, which reduces random color noise.<br>OFF, LOW, MID or HIGH.            |                            |



**Page 8 of 18 (Imaging Parameter 8 – Probe Dependent 8) (cont'd)**

|                                                                                                                                                                                                                |                     |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|
| <b>PDI Persistence</b>                                                                                                                                                                                         | <b>EC, UP, P, R</b> |
| Choose the default setting for Persistence, which retains highest color pixel value for a specified period of time or until a higher value is detected.<br>OFF, SHORT, MID or HIGH.                            |                     |
| <b>PDI MR-Flow</b>                                                                                                                                                                                             | <b>EC, UP, P, R</b> |
| Choose the default setting of the MR-Flow option (combination focus in PDI color area).<br>Limited to LA39, C551 (CAE), C358, E721, 739L and 546L probes only. Must have ACE-2 option installed.<br>On or Off. |                     |
| <b>PDI MR-Flow Combi Number</b>                                                                                                                                                                                | <b>EC, UP, P</b>    |
| Set the number of focal zones to be used for the combination focus for MR-Flow option.<br>2 or 3.                                                                                                              |                     |



Page 9 of 18 (Imaging Parameter 9 – Probe Dependent 9)

|                                                                                                                                                   |                     |
|---------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|
| <b>3D Aspect Ratio</b>                                                                                                                            | <b>EC, UP, P, R</b> |
| Choose the ratio between depth and width in 3D-Mode. Depth is the first number.<br>Width is the second number.<br>2/1, 1/1, 1/2, 1/3, 1/4 or 1/5. |                     |
| <b>3D Swing Angle [deg]</b>                                                                                                                       | <b>EC, UP, P</b>    |
| Choose the swinging angle of the MIP image display.<br>Angle from 5 to 60 degrees in 5 degree increments.                                         |                     |
| <b>3D Swing Center [%]</b>                                                                                                                        | <b>EC, UP, P</b>    |
| Choose the center position to swing the MIP image display.<br>Percentage from 0 to 100% in 10% increments.                                        |                     |
| <b>3D View Direction</b>                                                                                                                          | <b>EC, UP, P, R</b> |
| Choose the viewpoint for the construction of the MIP images when in 3D-Mode.<br>Normal or Reverse.                                                |                     |




Figure 14–7. MIP Image Principles

**CAUTION**



The marker that indicates the viewpoint for the MIP image is displayed on the left side of the monitor.

Normal is indicated by the marker “↗ ↑ ↖”.

Reverse is indicated by the marker “↘ ↓ ↙”.

Ensure that this marker is consulted for MIP image viewpoint.



## Page 10 of 18 (Imaging Parameter 10)

|                                                                                                                                                                                                                                         |                                                                                                     |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|
| <b>Image Mode</b>                                                                                                                                                                                                                       | <b>EC, R</b>                                                                                        |
| Choose the default imaging mode for this preset.                                                                                                                                                                                        |                                                                                                     |
| B                                                                                                                                                                                                                                       | B-Mode only                                                                                         |
| B/B                                                                                                                                                                                                                                     | Dual B-Mode                                                                                         |
| B/M                                                                                                                                                                                                                                     | B- and M-Modes                                                                                      |
| B/D                                                                                                                                                                                                                                     | B- and Doppler Modes                                                                                |
| CFM-B                                                                                                                                                                                                                                   | B-Mode with Color Flow                                                                              |
| CFM-B/B                                                                                                                                                                                                                                 | B-Mode with Color Flow and B-Mode                                                                   |
| CFM-B/M                                                                                                                                                                                                                                 | B- and M-Modes with Color Flow                                                                      |
| CFM-B/D                                                                                                                                                                                                                                 | B- and Doppler Modes with Color Flow                                                                |
| <b>Acoustic Power Knob Effect</b>                                                                                                                                                                                                       | <b>EC</b>                                                                                           |
| Set how the Acoustic Output knob will affect the scan image.                                                                                                                                                                            |                                                                                                     |
| <b>Both B and D</b> —Knob changes power for B/M/D/CFM/PDI all the time regardless of scan mode.                                                                                                                                         |                                                                                                     |
| <b>Only B</b> —Knob changes power for B or B/M mode only. D/CFM/PDI modes cannot be adjusted by Acoustic power knob. If B or M-Modes are in combination with D/CFM/PDI they also cannot be adjusted. It is fixed at setup preset value. |                                                                                                     |
| <b>Independent B and D</b> —Knob changes power for B/M-Mode while in B/M-Mode. Knob changes power for D/CFM/PDI mode while in D/CFM/PDI mode.                                                                                           |                                                                                                     |
| <b>No Effect</b> —Acoustic output knob has no effect on changing power. Values fixed at set up preset numbers.                                                                                                                          |                                                                                                     |
| <b>Displayed TI</b>                                                                                                                                                                                                                     | <b>EC</b>                                                                                           |
| Choose the type of thermal index display.                                                                                                                                                                                               |                                                                                                     |
| TIS                                                                                                                                                                                                                                     | Thermal Index – Soft Tissue                                                                         |
| TIC                                                                                                                                                                                                                                     | Thermal Index – Cranial Bone                                                                        |
| TIB                                                                                                                                                                                                                                     | Thermal Index – Bone                                                                                |
| <b>Image Gamma [%]</b>                                                                                                                                                                                                                  | <b>EC, UP</b>                                                                                       |
| Choose the level of gamma correction applied to the image during photography.                                                                                                                                                           |                                                                                                     |
| Choose from 50 to 200 in 5 digit increments.                                                                                                                                                                                            |                                                                                                     |
| <b>Auto Focus Control</b>                                                                                                                                                                                                               | <b>EC</b>                                                                                           |
| Choose to have the focus control follow depth changes as a percent of depth or have the focus depth be a fixed value.                                                                                                                   |                                                                                                     |
| On                                                                                                                                                                                                                                      | Focus is a percent of the current depth (Focus Pos % Depth preset) until changed by focus position. |
| Off                                                                                                                                                                                                                                     | Focus is at a fixed depth until changed by focus position control (Focus Depth [mm] Preset).        |



Page 10 of 18 (Imaging Parameter 10) (cont'd)

|                                                                                                            |                                                           |
|------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|
| <b>Zoom Reference Update    EC, UP</b>                                                                     |                                                           |
| Choose to update timing of reference image in Zoom.                                                        |                                                           |
| Off<br>0.5s, 1.0s, or 2.0s                                                                                 | No reference update,<br>Update at time interval selected. |
| <b>TGC Depth Remap    EC</b>                                                                               |                                                           |
| Choose the default setting for the TGC graphic to be recomputed and displayed when parameters are changed. |                                                           |
| On or Off.                                                                                                 |                                                           |
| <b>Cine Gauge Auto Display    EC</b>                                                                       |                                                           |
| Choose the initial default setting for the Cine gauge to be displayed when Cine is activated.              |                                                           |
| On or Off.                                                                                                 |                                                           |
| <b>Cine Speed    EC, R</b>                                                                                 |                                                           |
| Choose the default Cine loop playback speed.                                                               |                                                           |
| 1, 1/2, 1/4, 1/8 or 1/16.                                                                                  |                                                           |
| <b>Image Back Color    EC, UP</b>                                                                          |                                                           |
| Choose the image area background color.                                                                    |                                                           |
| BLK                                                                                                        | Black                                                     |
| GR1                                                                                                        | 10% Gray                                                  |
| GR2                                                                                                        | 20% Gray                                                  |
| GR3                                                                                                        | 30% Gray                                                  |
| DRD                                                                                                        | Dark Red                                                  |
| BDX                                                                                                        | Bordeaux Red                                              |
| GGR                                                                                                        | Green Gray                                                |
| DGR                                                                                                        | Dark Green                                                |
| <b>Frame Rate Display    EC</b>                                                                            |                                                           |
| Choose the default setting of the frame rate display on the screen.                                        |                                                           |
| On or Off.                                                                                                 |                                                           |
| <b>Probe Direction Mark Display    EC</b>                                                                  |                                                           |
| Choose the default setting of the probe direction marker display.                                          |                                                           |
| On or Off.                                                                                                 |                                                           |



## Page 10 of 18 (Imaging Parameter 10) (cont'd)

| Audio Tone Equalizer    EC                                                                 |                                             |            |               |  |
|--------------------------------------------------------------------------------------------|---------------------------------------------|------------|---------------|--|
| Choose the desired type of audio output response                                           |                                             |            |               |  |
| Flat                                                                                       | Normal setting                              |            |               |  |
| 1                                                                                          | Reduced high tones, increased mid tones     |            |               |  |
| 2                                                                                          | Reduced high tones, increased mid/low tones |            |               |  |
| 3                                                                                          | Reduced high tones, increased mid tones     |            |               |  |
|                                                                                            | <u>Bass</u>                                 | <u>Mid</u> | <u>Treble</u> |  |
| Flat                                                                                       | 0db                                         | 0db        | 0db           |  |
| 1                                                                                          | 0db                                         | 2db        | -4db          |  |
| 2                                                                                          | 2db                                         | 2db        | -2db          |  |
| 3                                                                                          | 0db                                         | 4db        | -2db          |  |
| Display TGC Curve    EC                                                                    |                                             |            |               |  |
| Choose the default setting of the TGC Curve display for both large and small image format. |                                             |            |               |  |
| On or Off.                                                                                 |                                             |            |               |  |

## Page 11 of 18 (Imaging Parameter 11)

| B Focus Dual Marker    EC, UP                                                                            |                                                                                 |
|----------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------|
| Choose to display focus markers on both dual B-Mode images or just to the right side of the right image. |                                                                                 |
| Off                                                                                                      | Only one focus marker is displayed.                                             |
| On                                                                                                       | Focus markers will be displayed to the right of both the right and left images. |
| B Color Auto Start    EC, R                                                                              |                                                                                 |
| Choose the initial default setting for B color (colorized gray scale).                                   |                                                                                 |
| On or Off.                                                                                               |                                                                                 |
| B Color Map    EC, UP                                                                                    |                                                                                 |
| Choose a B color default map.                                                                            |                                                                                 |
| 1, 2, 3, 4, 5 or 6.                                                                                      |                                                                                 |
| B Tag Auto Start    EC, R                                                                                |                                                                                 |
| Choose the initial default setting for B color tag.                                                      |                                                                                 |
| On or Off.                                                                                               |                                                                                 |





Page 11 of 18 (Imaging Parameter 11) (cont'd)

|                                                                                                                                                                                               |              |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|
| <b>B Tag Center [%]</b>                                                                                                                                                                       | <b>EC, R</b> |
| Choose the center point of the B color tag window as a default.<br>Percentage from 0 to 100%. A percentage of the entire displayed B color window (colorized gray scale).                     |              |
| <b>B Tag Width [%]</b>                                                                                                                                                                        | <b>EC, R</b> |
| Choose a default for the desired B color tag width.<br>Percentage from 0 to 100% for the B color window size.                                                                                 |              |
| <b>B Format B Single</b>                                                                                                                                                                      | <b>EC</b>    |
| Choose a default display size for a single B-Mode image.<br>Large or small.                                                                                                                   |              |
| <b>B Format B Dual</b>                                                                                                                                                                        | <b>EC</b>    |
| Choose a default display size for a dual B-Mode image.<br>Large or small.                                                                                                                     |              |
| <b>B Scale Mark</b>                                                                                                                                                                           | <b>EC</b>    |
| Choose as a default what the scale markers will represent in B-Mode.<br><br>Depth/Width    Depth and width marks.<br>Depth            Depth marks only.<br>Comb             Combination mode. |              |
| <b>B Video Inverse</b>                                                                                                                                                                        | <b>EC</b>    |
| Choose as a default how to display B-Mode video.<br><br>Off                White on black.<br>On                 Black on white.                                                              |              |



## Page 12 of 18 (Imaging Parameter 12)

|                                                                                                                                                                         |                                                                        |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------|
| <b>Timeline Sweep Method</b>                                                                                                                                            | <b>EC</b>                                                              |
| Choose as a default between two types of timeline presentations.                                                                                                        |                                                                        |
| Scroll                                                                                                                                                                  | Continuous update at edge of the screen.                               |
| Moving-Bar                                                                                                                                                              | Timeline is updated by a bar moving across the screen.                 |
| <b>Timeline Format</b>                                                                                                                                                  | <b>EC</b>                                                              |
| Choose as a default how the timeline (M-Mode or Doppler spectrum) will be displayed. Only one type of display format is available at one time.                          |                                                                        |
| Side/Side                                                                                                                                                               | B-Mode on left half of the display and the timeline on the right half. |
| Top/Bottom                                                                                                                                                              | B-Mode on the top of the display and the timeline across the bottom.   |
| <b>Auto B Melt at Unfreeze</b>                                                                                                                                          | <b>EC, UP</b>                                                          |
| When the B Pause function is used to freeze B-Mode while in B/M or B/D, the preset sets the default to have the B-Mode image remain paused when the system is unfrozen. |                                                                        |
| Off                                                                                                                                                                     | B-Mode remains paused at unfreeze.                                     |
| On                                                                                                                                                                      | B-Mode exits pause condition at unfreeze.                              |
| <b>M Mode Simul PD Activity</b>                                                                                                                                         | <b>EC</b>                                                              |
| Choose as a default to have the ability to display simultaneous M-Mode and pulsed Doppler spectrum.                                                                     |                                                                        |
| Off                                                                                                                                                                     | No simultaneous M and PD display.                                      |
| On                                                                                                                                                                      | Simultaneous display capable.                                          |
| <b>M Color Auto Start</b>                                                                                                                                               | <b>EC, R</b>                                                           |
| Choose the initial default setting for M color (colorized gray scale).                                                                                                  |                                                                        |
| On or Off.                                                                                                                                                              |                                                                        |
| <b>M Color Map</b>                                                                                                                                                      | <b>EC, UP</b>                                                          |
| Choose a M color default map.                                                                                                                                           |                                                                        |
| 1, 2, 3, 4, 5 or 6.                                                                                                                                                     |                                                                        |
| <b>M Tag Auto Start</b>                                                                                                                                                 | <b>EC, R</b>                                                           |
| Choose the initial default setting for M color tag.                                                                                                                     |                                                                        |
| On or Off.                                                                                                                                                              |                                                                        |
| <b>M Tag Center [%]</b>                                                                                                                                                 | <b>EC, R</b>                                                           |
| Choose the center point of the M color tag window as a default.                                                                                                         |                                                                        |
| Percentage from 0 to 100%. A percentage of the entire displayed M color window (colorized gray scale).                                                                  |                                                                        |



Page 12 of 18 (Imaging Parameter 12) (cont'd)

|                                                                                                                                                                                                                                              |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>M Tag Width [%]    EC, R</b>                                                                                                                                                                                                              |
| Choose a default for the desired M color tag width.<br>Percentage from 0 to 100% for the M color window size.                                                                                                                                |
| <b>M Format (Side/Side, Full)    EC</b>                                                                                                                                                                                                      |
| Choose the default size for the M-Mode display format when it is full screen or side by side. Side/Side must be selected for "Timeline Format" on the Set Up/Custom Display Menu page 12 in order to display this format.<br>Large or small. |
| <b>M Mode FULL B + Cursor    EC</b>                                                                                                                                                                                                          |
| Choose the default setting for the full B plus M-Mode cursor display format.<br>On or Off.                                                                                                                                                   |
| <b>M Mode Side/Side B Large    EC</b>                                                                                                                                                                                                        |
| Choose the default setting for the large B-Mode (left side) and M-Mode (right side) format. Side/Side must be selected for "Timeline Format" on the Set Up/Custom Display Menu page 12 in order to display this format.<br>On or Off.        |
| <b>M Mode Side/Side B Small    EC</b>                                                                                                                                                                                                        |
| Choose the default setting for the small B-Mode (left side) and M-Mode (right side) format. Side/Side must be selected for "Timeline Format" on the Set Up/Custom Display Menu page 12 in order to display this format.<br>On or Off.        |
| <b>M Mode Top/Bottom B Large    EC</b>                                                                                                                                                                                                       |
| Choose the default setting for the large B-Mode (top) and M-Mode (bottom) format. Top/Bottom must be selected for "Timeline Format" on the Set Up/Custom Display Menu page 12 in order to display this format.<br>On or Off.                 |
| <b>M Mode Top/Bottom B Mid    EC</b>                                                                                                                                                                                                         |
| Choose the default setting for the medium size B-Mode (top) and M-Mode (bottom) format. Top/Bottom must be selected for "Timeline Format" on the Set Up/Custom Display Menu page 12 in order to display this format.<br>On or Off.           |



## Page 12 of 18 (Imaging Parameter 12) (cont'd)

|                                                                                                                                                                                                                     |                                     |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------|
| <b>M Mode Top/Bottom B Small</b>                                                                                                                                                                                    | <b>EC</b>                           |
| Choose the default setting for the small size B-Mode (top) and M-Mode (bottom) format. Top/Bottom must be selected for "Timeline Format" on the Set Up/Custom Display Menu page 12 in order to display this format. |                                     |
| On or Off.                                                                                                                                                                                                          |                                     |
| <b>M Mode Full M</b>                                                                                                                                                                                                | <b>EC</b>                           |
| Choose the default setting for the full screen M-Mode format.                                                                                                                                                       |                                     |
| On or Off.                                                                                                                                                                                                          |                                     |
| <b>M Scale Mark</b>                                                                                                                                                                                                 | <b>EC</b>                           |
| Choose the method to display scale markers in M-Mode.                                                                                                                                                               |                                     |
| Time/Depth                                                                                                                                                                                                          | Time horizontal and depth vertical. |
| Time                                                                                                                                                                                                                | Time horizontal only.               |
| <b>M Video Inverse</b>                                                                                                                                                                                              | <b>EC</b>                           |
| Choose the display video type for the M-Mode portion of the display.                                                                                                                                                |                                     |
| Off                                                                                                                                                                                                                 | White on black (inverse off).       |
| On                                                                                                                                                                                                                  | Black on white (inverse on).        |

## Page 13 of 18 (Imaging Parameter 13)

|                                                                                          |                                                            |
|------------------------------------------------------------------------------------------|------------------------------------------------------------|
| <b>D/CFM Scale Mark</b>                                                                  | <b>EC</b>                                                  |
| Choose the default value for the scale markers in Doppler with color flow.               |                                                            |
| Velo                                                                                     | Velocity                                                   |
| Frq                                                                                      | Frequency                                                  |
| Velo/Frq                                                                                 | Velocity/Frequency                                         |
| Frq/Velo                                                                                 | Frequency/Velocity                                         |
| <b>D Type for Wall Filter</b>                                                            | <b>EC</b>                                                  |
| Choose the type of wall filter display while in Doppler.                                 |                                                            |
| Velo                                                                                     | Velocity                                                   |
| Freq                                                                                     | Frequency                                                  |
| <b>D Realtime Auto Trace</b>                                                             | <b>EC, R</b>                                               |
| Choose the initial default setting for the auto trace function during real-time Doppler. |                                                            |
| Off                                                                                      | No Real-time trace function.                               |
| On                                                                                       | Real-time trace function enabled.                          |
| Calc                                                                                     | Real-time trace function with calculation display enabled. |

Page 13 of 18 (Imaging Parameter 13) (cont'd)

|                                                                                                                             |                                           |
|-----------------------------------------------------------------------------------------------------------------------------|-------------------------------------------|
| <b>D Realtime Trace Method</b>                                                                                              | <b>EC, R</b>                              |
| Choose the value that will be traced during real-time Doppler.                                                              |                                           |
| Peak                                                                                                                        | Peak or highest values (TAMAX).           |
| Floor                                                                                                                       | Bottom or lowest values and peak (TAMIN). |
| Mean                                                                                                                        | Average values and peak (TAMEAN).         |
| Mode                                                                                                                        | Brightest echo trace and peak (TAMODE).   |
| <b>D Realtime Calc Dir.</b>                                                                                                 | <b>EC, R</b>                              |
| Choose the default for the portion of the Doppler signal to be auto traced.                                                 |                                           |
| Compo                                                                                                                       | Forward and reverse flow.                 |
| Forwd                                                                                                                       | Only forward flow.                        |
| Revrs                                                                                                                       | Only reverse flow.                        |
| <b>D Color Auto Start</b>                                                                                                   | <b>EC, R</b>                              |
| Choose the initial default setting for D color (colorized gray scale).                                                      |                                           |
| On or Off.                                                                                                                  |                                           |
| <b>D Color Map</b>                                                                                                          | <b>EC, UP</b>                             |
| Choose a D color default map.                                                                                               |                                           |
| 1, 2, 3, 4, 5 or 6.                                                                                                         |                                           |
| <b>D Tag Auto Start</b>                                                                                                     | <b>EC, R</b>                              |
| Choose the initial default setting for D color tag.                                                                         |                                           |
| On or Off.                                                                                                                  |                                           |
| <b>D Tag Center [%]</b>                                                                                                     | <b>EC, R</b>                              |
| Choose the center point of the D color tag window as a default.                                                             |                                           |
| Percentage from 0 to 100%. A percentage of the entire displayed D color window (colorized gray scale).                      |                                           |
| <b>D Tag Width [%]</b>                                                                                                      | <b>EC, R</b>                              |
| Choose a default for the desired D color tag width.                                                                         |                                           |
| Percentage from 0 to 100% for the D color window size.                                                                      |                                           |
| <b>D Format (Side/Side, Full)</b>                                                                                           | <b>EC</b>                                 |
| Choose the default Doppler display size for the side by side or full Doppler formats.                                       |                                           |
| Side/Side must be selected for "Timeline Format" on the Set Up/Custom Display Menu page 12 in order to display this format. |                                           |
| Large, small or very-small.                                                                                                 |                                           |



**Page 13 of 18 (Imaging Parameter 13) (cont'd)**

|                                                                                                                                                                                                                                         |           |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| <b>D Mode Full B + Cursor</b>                                                                                                                                                                                                           | <b>EC</b> |
| Choose the default setting for the full B-Mode plus Doppler cursor format display.<br>On or Off.                                                                                                                                        |           |
| <b>D Mode Side/Side B Large</b>                                                                                                                                                                                                         | <b>EC</b> |
| Choose the default setting for the large B-Mode (left side) and Doppler (right side) format. Side/Side must be selected for "Timeline Format" on the Set Up/Custom Display Menu page 12 in order to display this format.<br>On or Off.  |           |
| <b>D Mode Side/Side B Small</b>                                                                                                                                                                                                         | <b>EC</b> |
| Choose the default setting for the small B-Mode (left side) and Doppler (right side) format. Side/Side must be selected for "Timeline Format" on the Set Up/Custom Display Menu page 12 in order to display this format.<br>On or Off.  |           |
| <b>D Mode Top/Bottom B Large</b>                                                                                                                                                                                                        | <b>EC</b> |
| Choose the default setting for the large B-Mode (top) and Doppler spectrum (bottom) format. Top/Bottom must be selected for "Timeline Format" on the Set Up/Custom Display Menu page 12 in order to display this format.<br>On or Off.  |           |
| <b>D Mode Top/Bottom B Mid</b>                                                                                                                                                                                                          | <b>EC</b> |
| Choose the default setting for the medium B-Mode (top) and Doppler spectrum (bottom) format. Top/Bottom must be selected for "Timeline Format" on the Set Up/Custom Display Menu page 12 in order to display this format.<br>On or Off. |           |
| <b>D Mode Top/Bottom B Small</b>                                                                                                                                                                                                        | <b>EC</b> |
| Choose the default setting for the small B-Mode (top) and Doppler spectrum (bottom) format. Top/Bottom must be selected for "Timeline Format" on the Set Up/Custom Display Menu page 12 in order to display this format.<br>On or Off.  |           |
| <b>D Mode Full D</b>                                                                                                                                                                                                                    | <b>EC</b> |
| Choose the default setting for the full screen Doppler spectrum format.<br>On or Off.                                                                                                                                                   |           |

Page 14 of 18 (Imaging Parameter 14)

|                                                                                                            |                                                                                     |
|------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <b>D Image Processing    EC</b>                                                                            |                                                                                     |
| Choose the desired type of Doppler spectrum image processing.                                              |                                                                                     |
| Black Holes Filling                                                                                        | Reduce the noise in the Doppler spectrum by filling in vertically.                  |
| Normal                                                                                                     | Standard Processing. No additional filtering to reduce noise.                       |
| Smoothing                                                                                                  | Reduce the noise in the Doppler spectrum by filling in vertically and horizontally. |
| Leveling                                                                                                   | Reduce the noise for Doppler spectrum direction.                                    |
| <b>D Time Resolution    EC, UP</b>                                                                         |                                                                                     |
| Choose the amount of emphasis to be placed on Doppler time resolution as compared to frequency resolution. |                                                                                     |
| 1                                                                                                          | Time resolution much greater than frequency resolution.                             |
| 2                                                                                                          | Time resolution slightly greater than frequency resolution.                         |
| 3                                                                                                          | Frequency resolution slightly greater than time resolution.                         |
| 4                                                                                                          | Frequency resolution much greater than time resolution.                             |
| 5                                                                                                          | Frequency resolution the greatest compared to time resolution.                      |
| <b>D Audio to Speaker    EC</b>                                                                            |                                                                                     |
| Choose the default type of Doppler audio to be sent to the speakers.                                       |                                                                                     |
| Stereo                                                                                                     | Forward & Reverse signals separate.                                                 |
| Mono-F                                                                                                     | Forward signals only.                                                               |
| Mono-R                                                                                                     | Reverse signals only.                                                               |
| <b>D Audio Volume [%]    EC, R</b>                                                                         |                                                                                     |
| Choose the default percentage value for the Doppler volume to the speakers.                                |                                                                                     |
| Percentage from 0 to 100% in 4 digit steps.                                                                |                                                                                     |
| <b>D Video Inverse    EC</b>                                                                               |                                                                                     |
| Choose the default type of video used to display the Doppler spectrum.                                     |                                                                                     |
| Off                                                                                                        | White on black.                                                                     |
| On                                                                                                         | Black on white.                                                                     |
| <b>Display Dop Frequency    EC</b>                                                                         |                                                                                     |
| Choose the default setting for the Doppler Frequency display.                                              |                                                                                     |
| On or Off.                                                                                                 |                                                                                     |
| <b>Display Dop PRF    EC</b>                                                                               |                                                                                     |
| Choose the default setting for the Doppler PRF display on the left side of the monitor.                    |                                                                                     |
| On or Off.                                                                                                 |                                                                                     |
| <b>Display Spectrm Ref. Value    EC</b>                                                                    |                                                                                     |
| Choose the default setting for the display of the Doppler velocity increments on the spectrum.             |                                                                                     |
| On                                                                                                         | Only Doppler velocity reference scale is displayed.                                 |
| Off                                                                                                        | Only Doppler velocity maximum value is displayed.                                   |
| <b>Spec. Inv in Triplex    EC</b>                                                                          |                                                                                     |
| Choose the default setting for the color scale bar when Doppler spectrum is inverted.                      |                                                                                     |
| Independent                                                                                                | The color scale bar remains the same when the Doppler spectrum is inverted.         |
| Simultaneous                                                                                               | The color scale bar inverts when the Doppler spectrum is inverted.                  |



## Page 15 of 18 (Imaging Parameter 15)

|                                                                                                                                                                                                                                                                                                                                                                     |               |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|
| <b>PWD HPRF</b>                                                                                                                                                                                                                                                                                                                                                     | <b>EC, UP</b> |
| Choose the default setting for the HPRF selection in PWD Mode.<br>On or Off.                                                                                                                                                                                                                                                                                        |               |
| <b>PWD Base Line [%]</b>                                                                                                                                                                                                                                                                                                                                            | <b>EC, R</b>  |
| Choose the default percentage value for forward pulsed Doppler spectrum display.<br>Percentage from –75% to +75% in 25% increments (0 is 50% forward and 50% reverse).                                                                                                                                                                                              |               |
| <b>CWD Base Line [%]</b>                                                                                                                                                                                                                                                                                                                                            | <b>EC, R</b>  |
| Choose the default percentage value for forward continuous wave Doppler spectrum display.<br>Percentage from –75% to +75% in 25% increments (0 is 50% forward and 50% reverse).                                                                                                                                                                                     |               |
| <b>CWD Non-aliasing mode</b>                                                                                                                                                                                                                                                                                                                                        | <b>EC</b>     |
| Choose how to process the CWD signal for aliasing.<br>Off                      Non-aliasing mode off.<br>On                        Increase signal to noise ratio to reduce aliasing.                                                                                                                                                                               |               |
| <b>CFM Capture Interval [s]</b>                                                                                                                                                                                                                                                                                                                                     | <b>EC, R</b>  |
| Choose the default value (time) for the CFM capture interval.<br>Off                      No time interval update.<br>1/2                      Update at .5 second intervals.<br>1                         Update at 1 second intervals.<br>2                         Update at 2 second intervals.<br>ECG                     Update synchronized with ECG signal. |               |





Page 15 of 18 (Imaging Parameter 15) (cont'd)

|                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>CFM Velo Dynamic Range</b> <b>EC, UP</b>                                                                                                                                                                                                                                                                                                                        |
| Choose the default value for the CFM velocity dynamic range display.<br>Value from 10 to 100 in 10 digit increments.                                                                                                                                                                                                                                               |
| <b>CFM/PDI Sequence</b> <b>EC, UP</b>                                                                                                                                                                                                                                                                                                                              |
| Use the Ellipse rocker switch to choose the default diagnostic mode sequence for Color Flow Mapping or Power Doppler Imaging.<br><br>Srvy>Mp>MpDtl            Survey, map and map/detail.<br>Srvy>SrvyDtl            Survey and survey/detail.<br>SrvyDtl>MpDtl            Survey/detail and map/detail.<br>Srvy>SrvyDtl>Map        Survey, survey/detail and map. |
| <b>CFM/PDI Initial Sequence</b> <b>EC</b>                                                                                                                                                                                                                                                                                                                          |
| Choose the initial default diagnostic mode for Color Flow Mapping or Power Doppler Imaging.<br><br>Survey, map, survey/detail or map/detail.                                                                                                                                                                                                                       |
| <b>CFM Color Map 1, 2, 3, 4, 5 and 6</b> <b>EC, UP</b>                                                                                                                                                                                                                                                                                                             |
| Choose the default color map type for each of the four selection positions in the CFM mode.<br><br>Map V 1, 2, 3, 4, 5, 6, 7 or 8. (velocity)<br>Map VT 1, 2, 3, 4, 5, 6, 7 or 8. (velocity/turbulence)<br>Map T 1 or 2. (turbulence)<br>Map A 1, 2, 3 or 4. (advanced velocity maps)*<br><br>* ACE-2 option needed for advanced velocity maps.                    |
| <b>CFM Tag Auto Start</b> <b>EC, R</b>                                                                                                                                                                                                                                                                                                                             |
| Choose the initial default setting for CFM tag.<br><br>On or Off.                                                                                                                                                                                                                                                                                                  |
| <b>CFM Tag Center [cm/s]</b> <b>EC, R</b>                                                                                                                                                                                                                                                                                                                          |
| Choose the default center point in cm/sec for the CFM color tag window.<br><br>Value from 0 to 500.                                                                                                                                                                                                                                                                |
| <b>CFM Tag Width [cm/s]</b> <b>EC, R</b>                                                                                                                                                                                                                                                                                                                           |
| Choose the default width in cm/sec for the CFM color tag window.<br><br>Width from 0 to 1000.                                                                                                                                                                                                                                                                      |
| <b>CFM Base Line [%]</b> <b>EC, R</b>                                                                                                                                                                                                                                                                                                                              |
| Choose a default baseline for the CFM color map.<br><br>-75, -50, -25, 0, 25, 50 or 75.                                                                                                                                                                                                                                                                            |
|  75<br>0<br>-75                                                                                                                                                                                                                                                                 |



## Page 16 of 18 (Imaging Parameter 16)

|                                                                         |                                      |
|-------------------------------------------------------------------------|--------------------------------------|
| <b>CFM/PDI Initial Mode</b>                                             | <b>EC, R</b>                         |
| Choose the default mode for color CFM or PDI.                           |                                      |
| CFM                                                                     | Regular CFM mode.                    |
| PDI                                                                     | Power Doppler Imaging mode.          |
| <b>PDI B/W Display</b>                                                  | <b>EC</b>                            |
| Set the default parameter for the PDI transparent mode Maps 1 and 2.    |                                      |
| On or Off.                                                              |                                      |
| <b>PDI Default Map</b>                                                  | <b>EC</b>                            |
| Choose the default color map for the PDI mode.                          |                                      |
| P1, P2, P3, P4, P5, P6, P7 or P8                                        |                                      |
| <b>PDI Dynamic Range</b>                                                | <b>EC, UP</b>                        |
| Choose the default value for the Power Doppler dynamic range display.   |                                      |
| Value from 10 to 100 in 10 digit increments.                            |                                      |
| <b>PDI Tag Auto Start</b>                                               | <b>EC, R</b>                         |
| Choose the initial default setting for PDI tag.                         |                                      |
| On or Off.                                                              |                                      |
| <b>PDI Tag Center [cm/s]</b>                                            | <b>EC, R</b>                         |
| Choose the default center point in cm/sec for the PDI color tag window. |                                      |
| Value from 0 to 500.                                                    |                                      |
| <b>PDI Tag Width [cm/s]</b>                                             | <b>EC, R</b>                         |
| Choose the default width in cm/sec for the PDI color tag window.        |                                      |
| Width from 0 to 1000.                                                   |                                      |
| <b>PDI Capture Interval [s]</b>                                         | <b>EC, R</b>                         |
| Choose the default value (time) for the PDI capture interval.           |                                      |
| Off                                                                     | No time interval update.             |
| .5                                                                      | Update at .5 second intervals.       |
| 1                                                                       | Update at 1 second intervals.        |
| 2                                                                       | Update at 2 second intervals.        |
| ECG                                                                     | Update synchronized with ECG signal. |

Page 17 of 18 (Imaging Parameter 17)

|                                                                                                                                           |                  |
|-------------------------------------------------------------------------------------------------------------------------------------------|------------------|
| <b>ECG Auto Display</b>                                                                                                                   | <b>EC</b>        |
| Choose the initial default setting for the ECG waveform to be displayed on the monitor.<br>On or Off.                                     |                  |
| <b>ECG Gain</b>                                                                                                                           | <b>EC, R</b>     |
| Choose the default value for ECG gain.<br>Value from –20 to +10 in 2 digit increments.                                                    |                  |
| <b>ECG Position</b>                                                                                                                       | <b>EC, R</b>     |
| Choose the default value for the ECG display position (horizontal).<br>Value from –50 to +50 in 5 digit increments.                       |                  |
| <b>ECG Heart Rate Display</b>                                                                                                             | <b>EC</b>        |
| Choose the initial default setting for the heart rate, calculated from the ECG input, to be displayed on the monitor.<br>On or Off.       |                  |
| <b>ECG R Delay 1 [ms]</b>                                                                                                                 | <b>EC, UP, R</b> |
| Choose the default value in milliseconds for the R wave Delay 1.<br>Value from 0 to 2000 in 5 digit increments.                           |                  |
| <b>ECG R Delay 1 to 2 [ms]</b>                                                                                                            | <b>EC, UP, R</b> |
| Choose the default value in milliseconds for the R wave delay between Delay 1 and Delay 2.<br>Value from 0 to 2000 in 5 digit increments. |                  |
| <b>ECG Audio Beat Sound</b>                                                                                                               | <b>EC, UP</b>    |
| Choose the default setting to represent the ECG input as an audio tone.<br>On or Off.                                                     |                  |
| <b>ECG Audio Beat Sound Tone</b>                                                                                                          | <b>EC</b>        |
| Choose as a default the type of tone to be used to represent the ECG input.<br>Tone-1, tone-2 or tone-3.                                  |                  |



**Page 17 of 18 (Imaging Parameter 17) (cont'd)**

|                                                                                                                                                                           |               |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|
| <b>ECG Sync Mark Display</b>                                                                                                                                              | <b>EC, UP</b> |
| Choose the default setting for the flashing heart mark display. The heart mark is displayed in sync with the R-wave only when a valid ECG trace is present.<br>On or Off. |               |
| <b>PCG Auto Display</b>                                                                                                                                                   | <b>EC</b>     |
| Choose the initial default setting for the PCG input displayed on the monitor.<br>On or Off.                                                                              |               |
| <b>PCG Gain</b>                                                                                                                                                           | <b>EC, R</b>  |
| Choose the default value for PCG gain.<br>Value from –10 to +10 in 2 digit increments.                                                                                    |               |
| <b>PCG Position</b>                                                                                                                                                       | <b>EC, R</b>  |
| Choose as a default the horizontal position of the PCG waveform on the display.<br>Value from –50 to +50 in 5 digit increments.                                           |               |
| <b>AUX Auto Display</b>                                                                                                                                                   | <b>EC</b>     |
| Choose the initial default setting for the auxiliary signal input to be displayed on the monitor.<br>On or Off.                                                           |               |
| <b>AUX Gain</b>                                                                                                                                                           | <b>EC, R</b>  |
| Choose the default gain value for the auxiliary input signal.<br>Value from –10 to +10 in 2 digit increments.                                                             |               |
| <b>AUX Position</b>                                                                                                                                                       | <b>EC, R</b>  |
| Choose as a default the horizontal position of the auxiliary waveform on the display.<br>Value from –50 to +50 in 5 digit increments.                                     |               |



Page 18 of 18 (Imaging Parameter 18)

|                                                                                                                                                                                                                                                             |                  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
| <b>3D Swing Speed</b>                                                                                                                                                                                                                                       | <b>EC, UP, R</b> |
| Choose the swing speed of the displayed MIP image.<br>1/1, 1/2, 1/4, 1/8 or 1/16 with 1/1 being the fastest swing speed.                                                                                                                                    |                  |
| <b>3D Scan Line Resolution</b>                                                                                                                                                                                                                              | <b>EC, UP</b>    |
| Choose the default setting for 3D-Mode Scan Line Resolution.<br>Normal or High.                                                                                                                                                                             |                  |
| <b>3D Number of Frames in Build</b>                                                                                                                                                                                                                         | <b>EC, UP</b>    |
| Choose the number of frames to be constructed for the final MIP image. The higher the number the better the 3D resolution but with longer rendering time.<br>Number from 5 to 21 in 2 frame steps.                                                          |                  |
| <b>3D B Percent of Source Frames</b>                                                                                                                                                                                                                        | <b>EC, UP</b>    |
| Choose the number of source frames to be used for constructing the B/W MIP image.<br>Number from 1/1, 1/2, 1/3, 1/4 or 1/5.                                                                                                                                 |                  |
| <b>3D B Projection Technique</b>                                                                                                                                                                                                                            | <b>EC, UP, R</b> |
| Choose the method to be used in constructing the MIP image.<br><br>Min Minimum Intensity Projection.<br>Max Maximum Intensity Projection.<br>Grad Backward images are faded.<br>Inv Same as Min but the gray levels of the constructed images are inverted. |                  |
| <b>3D B Gradient Threshold [%]</b>                                                                                                                                                                                                                          | <b>EC, UP</b>    |
| Choose the gray level threshold percentage to be used in the B/W Gradient Mode.<br>Percentage from 0 to 100% in 2% increments.                                                                                                                              |                  |
| <b>3D B Background Level [%]</b>                                                                                                                                                                                                                            | <b>EC, UP</b>    |
| Choose the percentage of the gray level for the background of the B/W MIP images.<br>Percentage from 0 to 100% in 2% increments.                                                                                                                            |                  |
| <b>3D B Frame Line Level [%]</b>                                                                                                                                                                                                                            | <b>EC, UP</b>    |
| Choose the percentage of the gray level to be used for the frame line of the B/W MIP image.<br>Percentage from 0 to 100% in 2% increments.                                                                                                                  |                  |



## Page 18 of 18 (Imaging Parameter 18) (cont'd)

|                                                                                                                                                                                                                        |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>3D CFM/PDI Percent of Source Frames [%] EC, UP *</b>                                                                                                                                                                |
| Choose the number of source frames to be used for constructing the CFM MIP image.<br>1/1, 1/2, 1/3, 1/4 or 1/5.                                                                                                        |
| <b>3D CFM/PDI Projection Technique EC, UP, R *</b>                                                                                                                                                                     |
| Choose the method to be used in constructing the MIP images.<br><br>Norm            Normal Maximum Intensity Construction.<br>Grad            Backward images are faded.<br>Shaded        Object's edges are enhanced. |
| <b>3D CFM/PDI Gradient Threshold [%] EC, UP *</b>                                                                                                                                                                      |
| Choose the color threshold percentage to be used in the Gradient Mode.<br><br>Percentage from 0 to 100% in 2% increments.                                                                                              |
| <b>3D CFM/PDI Shade Object Threshold [%] EC, UP *</b>                                                                                                                                                                  |
| Choose the color threshold percentage to be used for the construction in Shaded Mode.<br><br>Percentage from 0 to 100% in 2% increments.                                                                               |
| <b>3D CFM/PDI Shade Boundary Threshold [%] EC, UP *</b>                                                                                                                                                                |
| Choose the color threshold percentage to be used for the boundary in Shaded Mode.<br><br>Percentage from 0 to 100% in 2% increments.                                                                                   |
| <b>3D CFM/PDI Frame Line Level [%] EC, UP *</b>                                                                                                                                                                        |
| Choose the percentage of the gray level to be used for the frame line of the B/W MIP images.<br><br>Percentage from 0 to 100% in 2% increments.                                                                        |



*\* NOTE: 3D images cannot be rendered in PDI-Mode.*

# System Parameters

## Overview

System Paramtr (parameters) in the Set Up Top Menu is 7 pages of preset parameters that when set, are the same for all exam categories and preset names. Changing system parameters also affects these exam application presets found in factory default menu.

Generally, these items are common to modes, probes, formats and applications.

| B              | Preset         | Set Up      | ECG         |
|----------------|----------------|-------------|-------------|
| Custom Display | Preset Program | Save Values | User Define |
| System Paramtr |                | Utility     | Diag.       |

Figure 14–8. Setup Sub-Menu (System Parameters)

A brief description of each preset parameter and the choices available for that preset follows.



*NOTE: If it is necessary to quickly return to page one of the System Parameters Sub-Menus, press **Ctrl** and **1** simultaneously.*



## System Parameters

---

### System Parameters Contents




Figure 14–9. System Parameters Contents Page

The System Parameters Contents page shows a summary of the 7 menu pages.

**Trackball** the cursor to the desired page and press **Set**.



Page 1 of 7 (System Setup)

|                                                                                                               |                                                                                                                        |
|---------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|
| <b>Hospital Name</b>                                                                                          |                                                                                                                        |
| Enter the desired hospital or facility name. This appears at the top of all display screens and report pages. |                                                                                                                        |
| Maximum of 29 alphanumeric characters.                                                                        |                                                                                                                        |
| <b>Date Format</b>                                                                                            |                                                                                                                        |
| Choose the desired presentation format for all dates displayed.                                               |                                                                                                                        |
| Year/Month/Day, Month/Day/Year, Day/Month/Year                                                                |                                                                                                                        |
| <b>Time Format</b>                                                                                            |                                                                                                                        |
| Choose the desired presentation format for all time graphics displayed.                                       |                                                                                                                        |
| 24 H                                                                                                          | 24 hour clock                                                                                                          |
| 12 H                                                                                                          | 12 hour clock                                                                                                          |
| <b>Language</b>                                                                                               |                                                                                                                        |
| Choose the desired language for all displayed graphics by using the <b>Ellipse</b> rocker switch.             |                                                                                                                        |
| English, French, German, Spanish and Italian.                                                                 |                                                                                                                        |
| <b>Diag Category</b>                                                                                          |                                                                                                                        |
| Choose the desired default diagnostic category.                                                               |                                                                                                                        |
| Rad/Abd                                                                                                       | Radiology/Abdomen                                                                                                      |
| Obst                                                                                                          | Obstetrical                                                                                                            |
| Gyne                                                                                                          | Gynecology                                                                                                             |
| Card                                                                                                          | Cardiology                                                                                                             |
| Vasc                                                                                                          | Vascular                                                                                                               |
| Urol                                                                                                          | Urology                                                                                                                |
| SmlPts                                                                                                        | Small Parts                                                                                                            |
| <b>Change Probe and Category</b>                                                                              |                                                                                                                        |
| Choose to manually select the desired category with the probe.                                                |                                                                                                                        |
| Independent                                                                                                   | Manual selection of exam category and manual selection of probe.                                                       |
| Simultaneous                                                                                                  | When the probe is manually selected, the system automatically defaults to the assumed exam application for that probe. |
| <b>Report Cursor Blink</b>                                                                                    |                                                                                                                        |
| Choose to have the report cursor blinking or solid.                                                           |                                                                                                                        |
| Off                                                                                                           | Solid                                                                                                                  |
| On                                                                                                            | Blinking                                                                                                               |
| <b>Report Cursor Type</b>                                                                                     |                                                                                                                        |
| Choose the style of the report cursor desired.                                                                |                                                                                                                        |
| Underscore                                                                                                    |                                     |
| Block                                                                                                         |                                     |



## System Parameters

---

### Page 1 of 7 (System Setup) (cont'd)

|                                                                                                                    |                                  |
|--------------------------------------------------------------------------------------------------------------------|----------------------------------|
| <b>CUA Selection Reset on New Patient</b>                                                                          |                                  |
| Choose to reset all Composite Ultrasound Age flags on the OB Report pages to yes when New Patient is selected.     |                                  |
| On or Off.                                                                                                         |                                  |
| <b>Insite Access Enable</b>                                                                                        |                                  |
| Choose the default setting to access remote diagnostics via telephone line.                                        |                                  |
| On or Off.                                                                                                         |                                  |
| <b>Keyboard LED Bright</b>                                                                                         |                                  |
| Choose the amount of illumination desired for the keyboard keys and switches.                                      |                                  |
| Dark                                                                                                               | Dim for dark rooms.              |
| Bright                                                                                                             | Brighter for well lighted rooms. |
| <b>Keyboard Menu Bright</b>                                                                                        |                                  |
| Choose the amount of illumination desired for the Soft Menu Vacuum Florescent Display.                             |                                  |
| Dark                                                                                                               | Dim setting.                     |
| Mid1                                                                                                               | Low medium setting.              |
| Mid2                                                                                                               | High medium setting.             |
| Bright                                                                                                             | Brightest setting.               |
| <b>Keyboard Tab</b>                                                                                                |                                  |
| Choose the comment cursor movement desired when using the Tab key.                                                 |                                  |
| Normal                                                                                                             | Eight character spaces.          |
| Word                                                                                                               | Next word (after a space).       |
| <b>Operation Error Beep</b>                                                                                        |                                  |
| Choose to disable the audio error beep or set the volume of the tone that is heard when a keystroke error is made. |                                  |
| Off                                                                                                                | No error beep.                   |
| Low                                                                                                                | Low volume beep on error.        |
| Loud                                                                                                               | High volume beep on error.       |
| <b>Standby Time [min]</b>                                                                                          |                                  |
| The system can be set to go into stand-by mode if no keyboard activity is sensed for a specific period of time.    |                                  |
| 5, 10, 15, 30, 45, 60 minutes or INF (infinity-never)                                                              |                                  |

Page 1 of 7 (System Setup) (cont'd)

|                                                                                                                                   |                           |
|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------|
| <b>Power on Memory Test</b>                                                                                                       |                           |
| Choose to execute the system memory test at power up.                                                                             |                           |
| On or Off.                                                                                                                        |                           |
| <b>Power on Status</b>                                                                                                            |                           |
| Choose to reset exam/system parameters at power up or set to the parameters saved when power was turned off.                      |                           |
| Keep-Latest                                                                                                                       | Set to latest values.     |
| Reset                                                                                                                             | Reset to system defaults. |
| <b>System Error Erase Interval</b>                                                                                                |                           |
| Set the time that a system error message is displayed on the monitor before it is erased.                                         |                           |
| 0.1, 1, 10 minutes or infinity.                                                                                                   |                           |
| Infinity is until the next error occurs.                                                                                          |                           |
| <b>Display Probe Name</b>                                                                                                         |                           |
| Choose to display the names of the probes attached to the console in the soft menu when an available Probe Select key is pressed. |                           |
| On or Off.                                                                                                                        |                           |

Page 2 of 7 (System Setup)

|                                                                                                |                               |
|------------------------------------------------------------------------------------------------|-------------------------------|
| <b>Setup Blue Back Start</b>                                                                   |                               |
| Choose to display a color for the background of the preset menus. Color is set by Col BackGnd. |                               |
| Off                                                                                            | Background is transparent.    |
| On                                                                                             | Display a colored background. |
| <b>Mask Image Characters</b>                                                                   |                               |
| Choose to mask (disable) the graphic characters on the image display.                          |                               |
| On or Off.                                                                                     |                               |
| <b>PCG Filter</b>                                                                              |                               |
| Choose the desired level of PCG (Physiological Signal) filtering.                              |                               |
| Low                                                                                            | Low level filtering.          |
| Mid1                                                                                           | Medium level filtering.       |
| <b>Delay Trigger Source</b>                                                                    |                               |
| Choose the source for the cardiac trigger delay.                                               |                               |
| ECG                                                                                            | Electrocardiogram input.      |
| AUX                                                                                            | Auxiliary signal input.       |



## System Parameters

---

### Page 2 of 7 (System Setup) (cont'd)

|                                                                                                                                                                               |                                               |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|
| <b>ID/Name Prohibition After Measurement</b>                                                                                                                                  |                                               |
| Choose to be able to use the ID/Name Key to change patient information after a measurement has been taken or lockout the use of the ID/Name key after a measurement is taken. |                                               |
| Off                                                                                                                                                                           | ID/Name Key functional after measurement.     |
| On                                                                                                                                                                            | ID/Name Key not functional after measurement. |
| <b>BSA Type</b>                                                                                                                                                               |                                               |
| Choose the formula used to calculate Body Surface Area.                                                                                                                       |                                               |
| Oriental                                                                                                                                                                      | Asian formula.                                |
| Occidental                                                                                                                                                                    | USA formula.                                  |
| <b>OB Report Format</b>                                                                                                                                                       |                                               |
| Choose the format for calculation and display of OB measurements.                                                                                                             |                                               |
| Tokyo                                                                                                                                                                         | Tokyo University Method.                      |
| Osaka                                                                                                                                                                         | Osaka University Method.                      |
| USA                                                                                                                                                                           | United States Method.                         |
| Europe                                                                                                                                                                        | European Method.                              |
| <b>EFBW Form (Europe)</b>                                                                                                                                                     |                                               |
| Choose the method for calculating estimate fetal birth weight for the European method.                                                                                        |                                               |
| Shep/Wars                                                                                                                                                                     | Shepherd/Warshoff Method.                     |
| Rich/Berk                                                                                                                                                                     | Richards/Berkowitz Method.                    |
| Hadlock                                                                                                                                                                       | Hadlock Method.                               |
| German                                                                                                                                                                        | German Accepted Method.                       |
| Merz                                                                                                                                                                          | German Method                                 |
| <b>EFBW Form (Tokyo)</b>                                                                                                                                                      |                                               |
| Choose the formula to be used for calculating the estimate fetal birth weight for the Tokyo method.                                                                           |                                               |
| Tokyo                                                                                                                                                                         | Tokyo University                              |
| Tokyo-S1                                                                                                                                                                      | Tokyo Sinozuka EFW1                           |
| Tokyo-S2                                                                                                                                                                      | Tokyo Sinozuka EFW2                           |
| Tokyo-S3                                                                                                                                                                      | Tokyo Sinozuka EFW3                           |
| <b>Col BackGnd</b>                                                                                                                                                            |                                               |
| Choose the color for the background of Preset and Patient Entry menus displayed on the monitor.                                                                               |                                               |
| Color options are:                                                                                                                                                            |                                               |
| BLK                                                                                                                                                                           | Black                                         |
| GRY                                                                                                                                                                           | Gray                                          |
| ULB                                                                                                                                                                           | Ultramarine Blue                              |
| BLU                                                                                                                                                                           | Blue                                          |
| PWB                                                                                                                                                                           | Periwinkle Blue                               |
| MGR                                                                                                                                                                           | Moss Green                                    |
| MNK                                                                                                                                                                           | Mink                                          |
| MRN                                                                                                                                                                           | Maroon                                        |

Page 2 of 7 (System Setup) (cont'd)

|                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|-------|-----|--------|-----|-------|-----|-------------|-----|-----------|-----|-----------|-----|-----------|-----|---------------|
| <b>Col Sys Inf</b>  | <p>Choose the color for all alphanumeric system information displayed on the monitor.</p> <p>Color options are:</p> <table> <tr><td>WHT</td><td>White</td></tr> <tr><td>YLW</td><td>Yellow</td></tr> <tr><td>GRN</td><td>Green</td></tr> <tr><td>LGR</td><td>Light Green</td></tr> <tr><td>PBL</td><td>Pale Blue</td></tr> <tr><td>PNK</td><td>Baby Pink</td></tr> <tr><td>RED</td><td>Vermilion</td></tr> <tr><td>YOR</td><td>Yellow Orange</td></tr> </table> | WHT | White | YLW | Yellow | GRN | Green | LGR | Light Green | PBL | Pale Blue | PNK | Baby Pink | RED | Vermilion | YOR | Yellow Orange |
| WHT                 | White                                                                                                                                                                                                                                                                                                                                                                                                                                                           |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| YLW                 | Yellow                                                                                                                                                                                                                                                                                                                                                                                                                                                          |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| GRN                 | Green                                                                                                                                                                                                                                                                                                                                                                                                                                                           |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| LGR                 | Light Green                                                                                                                                                                                                                                                                                                                                                                                                                                                     |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| PBL                 | Pale Blue                                                                                                                                                                                                                                                                                                                                                                                                                                                       |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| PNK                 | Baby Pink                                                                                                                                                                                                                                                                                                                                                                                                                                                       |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| RED                 | Vermilion                                                                                                                                                                                                                                                                                                                                                                                                                                                       |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| YOR                 | Yellow Orange                                                                                                                                                                                                                                                                                                                                                                                                                                                   |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| <b>Col Imag Crs</b> | <p>Choose the color for the image characters on the display (Probe Orientation, Zoom Depth, and Image Recall graphics). See color options under <i>Col Sys Inf</i>.</p>                                                                                                                                                                                                                                                                                         |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| <b>Col Scale</b>    | <p>Choose the color for the scale markers on the display. See color options under <i>Col Sys Inf</i>.</p>                                                                                                                                                                                                                                                                                                                                                       |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| <b>Col Comment</b>  | <p>Choose the color for comments and annotations entered on the display. See color options under <i>Col Sys Inf</i>.</p>                                                                                                                                                                                                                                                                                                                                        |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| <b>Col Mes Crs</b>  | <p>Choose the color for the fixed measurement cursor on the display. See color options under <i>Col Sys Inf</i>.</p>                                                                                                                                                                                                                                                                                                                                            |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| <b>Col Applcat</b>  | <p>Choose the color for completed measurement information displayed on the monitor. See color options under <i>Col Sys Inf</i>.</p>                                                                                                                                                                                                                                                                                                                             |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| <b>Col Moving</b>   | <p>Choose a color for the active movable cursor on the display. See color options under <i>Col Sys Inf</i>.</p>                                                                                                                                                                                                                                                                                                                                                 |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| <b>Col Aux</b>      | <p>Choose a color for all auxiliary signal information displayed on the monitor (ECG, AUX, and PCG). See color options under <i>Col Sys Inf</i>.</p>                                                                                                                                                                                                                                                                                                            |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |
| <b>Col Dop Trac</b> | <p>Choose a color for the Doppler Spectrum Trace function. See color options under <i>Col Sys Inf</i>.</p>                                                                                                                                                                                                                                                                                                                                                      |     |       |     |        |     |       |     |             |     |           |     |           |     |           |     |               |



## System Parameters

---

### Page 3 of 7 (System Setup)

|                                                                                                                                                                                                                                                                                                                                 |                                                                                                                                           |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Use A/N Keys for User Define</b>                                                                                                                                                                                                                                                                                             |                                                                                                                                           |
| Choose to use the alphanumeric keys A~Z and 0~9 with the user define function.<br>On or Off.                                                                                                                                                                                                                                    |                                                                                                                                           |
| <b>User Define Key Names &amp; Locks</b>                                                                                                                                                                                                                                                                                        |                                                                                                                                           |
| Use this sub-menu to review, name or lock/unlock all available user define keys. The 11 character name of each user defined key is displayed.<br><br>An asterisk to the left of the name means the key is locked. The user define function cannot be deleted or modified when the key is locked.<br><br>See 14–101 for details. |                                                                                                                                           |
| <b>CUA/AUA for HADLOCK</b>                                                                                                                                                                                                                                                                                                      |                                                                                                                                           |
| Choose how gestational age is derived for the Hadlock formulas.                                                                                                                                                                                                                                                                 |                                                                                                                                           |
| CUA                                                                                                                                                                                                                                                                                                                             | Composite Ultrasound Age—gestational age derived by linear regression.                                                                    |
| AUA                                                                                                                                                                                                                                                                                                                             | Average Ultrasound Age—gestational age derived by averaging.                                                                              |
| <b>PSA Calculation</b>                                                                                                                                                                                                                                                                                                          |                                                                                                                                           |
| Choose how the PSA Calculation is derived.                                                                                                                                                                                                                                                                                      |                                                                                                                                           |
| PSAD                                                                                                                                                                                                                                                                                                                            | Prostate Specific Antigen Density—defined as PSA divided by volume.                                                                       |
| PPSA                                                                                                                                                                                                                                                                                                                            | Predicted Prostate Specific Antigen—a number which gives the normal level of PSA that would be expected for a prostate of a given volume. |
| <b>Active Measurement Cursor</b>                                                                                                                                                                                                                                                                                                |                                                                                                                                           |
| Choose the Measurement active cursor type.<br><br>X–mark or Plus–mark.                                                                                                                                                                                                                                                          |                                                                                                                                           |



Page 4 of 7 (System Setup – Body Pattern)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Body Pattern Package 1–8</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <p>Choose up to 16 body patterns (2 sub-menu pages) for each of the 8 packages. Eight custom body pattern packages are then available for selection to be displayed in the Set Up/Preset Program Menu page one.</p> <p>An illustration of each body pattern available is shown on 6–24. To select the body pattern package to be used for each user preset within an exam category, refer to 14–66.</p> <p>See the following pages for details on Customizing or deleting a body pattern package .</p> |
| <b>Body Pattern Probe Change Erase</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <p>Choose to erase the displayed body pattern when a probe is changed.</p> <p>On or Off.</p>                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Body Pattern Unfreeze Erase</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <p>Choose to erase the body pattern when the system exits the freeze function.</p> <p>On or Off.</p>                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Body Pattern Copy To Active Side</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <p>Choose to have the selected body pattern copied to the active side during dual format displays.</p> <p>On or Off.</p>                                                                                                                                                                                                                                                                                                                                                                               |




## System Parameters

---

### Body pattern package (1–8) customization

Choose the **Set Up** Top Menu, then choose **System Paramtr** (system parameters) from the Set Up Sub-Menu and go to page 4.

Use the **Trackball** to move the arrow cursor to the word “Submenu” of the desired body pattern package number to be customized.





Press **Set**. The body pattern selection Sub-Menu appears.



Figure 14–10. Initial Body Pattern Sub-Menu

The rectangle in the display represents space for two rows of body patterns (2 Sub-Menu display pages). There is a space for 8 body patterns in each row.

Use the **Trackball** to move the arrow cursor to a desired body pattern location.



Press **Set**. A box appears around a Body Pattern or blank space if no pattern is programmed.

(continued)



## Body pattern package (1–8) customization (cont'd)



Figure 14–11. Adding Body Patterns



Use the **Ellipse** rocker switch to cycle through the body pattern choices.

When the desired body pattern is displayed, use the **Trackball** and **Zoom Size/Rotation** knob to position the probe orientation marker to the desired default position on the body pattern.

Press **Set** to complete the selection for that menu position.



**NOTE:** As many as 16 selections (2 rows of 8) can be made for each body pattern package.




## System Parameters

---

### Body Pattern Deletion




Use the **Trackball** to position the arrow cursor over the pattern to be deleted.



Press **Set**. A box appears around the Body Pattern.



Use the **Ellipse** rocker switch to cycle to the blank body pattern selection.



Press **Set** to record a blank space for that body pattern location.

Page 5 of 7 (System Setup – Recording)

|                                                                                                                                                                                                                                                                                                                          |                                       |                                  |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------|----------------------------------|
| <b>Report Video Inverse to Printer</b>                                                                                                                                                                                                                                                                                   |                                       |                                  |
| Choose to send calculation summary report video to the print device in inverse (reverse) video.<br>On or Off.                                                                                                                                                                                                            |                                       |                                  |
| <b>Record 1 B&amp;W or Color</b>                                                                                                                                                                                                                                                                                         |                                       |                                  |
| Choose the type of device to be controlled by the Record 1 key for black and white or color images. Use the Ellipse rocker switch to cycle through the selections.<br><br>Black and White Polaroid Camera, Color Polaroid Camera, Black and White Printer, Color Printer, Multi Image Camera, Image Archive or Computer. |                                       |                                  |
| <b>Record 2 B&amp;W or Color</b>                                                                                                                                                                                                                                                                                         |                                       |                                  |
| Choose the type of device to be controlled by the Record 2 key for black and white or color images. Use the Ellipse rocker switch to cycle through the selections.<br><br>Black and White Polaroid Camera, Color Polaroid Camera, Black and White Printer, Color Printer, Multi Image Camera, Image Archive or Computer. |                                       |                                  |
| <b>Report</b>                                                                                                                                                                                                                                                                                                            |                                       |                                  |
| Choose the type of device that will be used to print the calculation summary report pages.<br><br>Black and White Polaroid Camera, Color Polaroid Camera, Black and White Printer, Color Printer, Multi Image Camera or Image Archive.                                                                                   |                                       |                                  |
| <b>External Video Signal</b>                                                                                                                                                                                                                                                                                             |                                       |                                  |
| Choose the type of video signal, from an external source, that is being put into the LOGIQ™ 500.                                                                                                                                                                                                                         |                                       |                                  |
| S-Video                                                                                                                                                                                                                                                                                                                  | Two component video.                  |                                  |
| RGB                                                                                                                                                                                                                                                                                                                      | Red, green, blue and sync.            |                                  |
| Col-Comp                                                                                                                                                                                                                                                                                                                 | Composite color signal.               |                                  |
| <b>Color Printer Type</b>                                                                                                                                                                                                                                                                                                |                                       |                                  |
| Choose the type of color printer (option) that is attached to the LOGIQ™ 500.                                                                                                                                                                                                                                            |                                       |                                  |
| UP3030                                                                                                                                                                                                                                                                                                                   | Sony UP-3030MD.                       |                                  |
| UP1850                                                                                                                                                                                                                                                                                                                   | Sony UP-1800/UP-1850/UP-2950/UP-2850. |                                  |
| UP1850P                                                                                                                                                                                                                                                                                                                  | Sony UP-1850P.                        |                                  |
| <b>Color Printer Memory</b>                                                                                                                                                                                                                                                                                              |                                       |                                  |
| Choose how printer memory will be divided and thus how the color print will be formatted.                                                                                                                                                                                                                                |                                       |                                  |
| Single                                                                                                                                                                                                                                                                                                                   | One picture per sheet.                |                                  |
| Quad-Frm                                                                                                                                                                                                                                                                                                                 | Four pictures per sheet.              |                                  |
| Quad-FrmIs                                                                                                                                                                                                                                                                                                               | Not available.                        |                                  |
| UP-1800 Series/UP-2800 Series/UP-2900 Series                                                                                                                                                                                                                                                                             | <b>Single</b><br>Available            | <b>Quad-FRM</b><br>Not Available |
| UP-1850 Series/UP-2850 Series/UP-2950 Series                                                                                                                                                                                                                                                                             | Available                             | Available                        |



## System Parameters

---

### Page 5 of 7 (System Setup – Recording) (cont'd)

|                                                                                                 |                                     |
|-------------------------------------------------------------------------------------------------|-------------------------------------|
| <b>Color Printer Signal</b>                                                                     |                                     |
| Choose the type of video signal being sent to the color printer.                                |                                     |
| RGB                                                                                             | Red, green, blue and sync.          |
| Video                                                                                           | Composite Video.                    |
| S-Video                                                                                         | Two component video.                |
| <b>Plug 1–BW</b>                                                                                |                                     |
| Choose which device has its record/expose cable connected to plug one.                          |                                     |
| B&W–Prnt                                                                                        | Black and white video page printer. |
| B&W–Pola                                                                                        | Black and white Polaroid camera.    |
| Col–Pola                                                                                        | Color Polaroid camera.              |
| <b>Plug 2–CLR</b>                                                                               |                                     |
| Choose which device has its record/expose cable connected to plug two.                          |                                     |
| B&W–Prnt                                                                                        | Black and white video page printer. |
| B&W–Pola                                                                                        | Black and white Polaroid camera.    |
| Col–Pola                                                                                        | Color Polaroid camera.              |
| <b>Port MIC</b>                                                                                 |                                     |
| Choose the type of device that has an interface cable connected to the multi-image camera port. |                                     |
| Off                                                                                             | No connection.                      |
| Multi Image Camera                                                                              | Multi Image Camera.                 |
| VCR                                                                                             | Video Cassette Recorder.            |
| Color Printer                                                                                   | Color Printer.                      |
| Computer                                                                                        | Data Management Center Computer.    |
| Line Printer–USA Only                                                                           | Paper Line Printer.                 |
| <b>Port A/B</b>                                                                                 |                                     |
| Choose the type of device that has an interface cable connected to VCR/printer ports A and B.   |                                     |
| Off                                                                                             | No connection.                      |
| Multi Image Camera                                                                              | Multi Image Camera.                 |
| VCR                                                                                             | Video Cassette Recorder.            |
| Color Printer                                                                                   | Color Printer.                      |
| Computer                                                                                        | Data Management Center Computer.    |
| Line Printer–USA Only                                                                           | Paper Line Printer.                 |
| <b>Maskline Record</b>                                                                          |                                     |
| Choose to display the operator message area during the record function.                         |                                     |
| On or Off.                                                                                      |                                     |

Page 5 of 7 (System Setup – Recording) (cont'd)

|                                                                                                                                 |
|---------------------------------------------------------------------------------------------------------------------------------|
| <b>B&amp;W Printer Exposure Pulse Length [msec]</b>                                                                             |
| Set the exposure pulse length needed to trigger the B&W printer print function.<br>Length from 50 to 500 ms in 50 ms intervals. |
| <b>B&amp;W Printer Exposure Min Interval [sec]</b>                                                                              |
| Set the minimum interval in seconds between B&W printer exposure capabilities.<br>Interval from 1 to 10 in 0.5 sec intervals.   |
| <b>B&amp;W Image Archive</b>                                                                                                    |
| Choose how a Black & White image will be archived<br>B&W or Color.                                                              |
| <b>Color Image Archive</b>                                                                                                      |
| Choose how a Color image will be archived.<br>B&W or Color.                                                                     |

Page 6 of 7 (System Setup – User ID and Password)

|                                                                                                                                                                                                                                                                                                 |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Password Ask</b>                                                                                                                                                                                                                                                                             |
| Choose to use the password security option at power up by asking for the user ID and password.<br>On or Off.<br><i>NOTE: Refer to 3–18 for usage of the password function.</i>                                                                                                                  |
| <b>User ID 1–8</b>                                                                                                                                                                                                                                                                              |
| As many as eight user ID codes can be entered into the system parameters. These can be used for security purposes. When used with an accompanying password, user ID establishes limited access to the system at power up.<br>Maximum of 8 alphanumeric characters for each user ID.             |
| <b>Password 1–8</b>                                                                                                                                                                                                                                                                             |
| As many as 8 passwords can be entered to correspond with the 8 user ID codes. These must be used when password ask is on.<br>Maximum of 8 alphanumeric characters for each password.<br><i>Initially, User ID Codes are blank. A User ID must be entered before a password can be recorded.</i> |



## System Parameters

---

### Selecting an ID

- Select the desired User ID number on page 6 of the Set Up/System Parameter Menu, by moving the arrow cursor over the User ID number and pressing **Set**.
  - The User ID entry area is highlighted.
  - Type in a maximum of 8 characters for the User ID and press **Return**.
- The User ID is now recorded and displayed on the System Parameter Menu.

### Entering a Password

Once a User ID is established or changed, a password can be entered.

- Move the arrow cursor to the desired password and press **Set**.
- Type in the new password. Press **Return**.
- Type in the new password a second time for re-confirmation. Press **Return**.



*NOTE: The password function is upper and lower case sensitive.*

*The Password function is not active unless "Password Ask" in the Set Up/System Parameter menu page 6 is set to On.*

### Changing/Deleting a Password

To change or delete a password:

- Move the arrow cursor to the desired password and press **Set**.
- In the highlighted password area, type the old password. Press **Return**.
- To change a password, type in the new characters. Press **Return**. To delete a password, type in 8 spaces. Press **Return**.
- Type the new password or 8 spaces a second time for confirmation. Press **Return**.

The password is now changed or deleted.

## Changing/Deleting a User ID

To change or delete a User ID:

- Move the arrow cursor to the desired User ID and press **Set**.
- Input the old password. Press **Return**.
- Type in the new User ID to change or type in spaces to delete the User ID. Press **Return**.

The User ID is now changed or deleted.



**NOTE:** *If the User ID is deleted (all spaces), the corresponding password is also deleted.*

*If no password is selected with the User ID, the ID can be changed or deleted directly. There is no requirement to enter the old password.*

## Page 7 of 7 (DICOM)

| Menu Selection at New Patient                                                            |                                                                        |
|------------------------------------------------------------------------------------------|------------------------------------------------------------------------|
| Choose what screen will be displayed when the New Patient key is pressed.                |                                                                        |
| Patient_Entry                                                                            | Displays standard patient entry screen for the default exam category.  |
| Schedule                                                                                 | Displays the DICOM patient schedule screen.                            |
| Sort Criteria for Schedule List                                                          |                                                                        |
| Choose how to sort the entries on the DICOM patient schedule screen.                     |                                                                        |
| Date&Time                                                                                | Arrange by date and time of patient session.                           |
| Name                                                                                     | Arrange by Patient Name.                                               |
| Patient Name Format                                                                      |                                                                        |
| Choose how the Patient Name entry will be displayed for the Patient Entry Menu.          |                                                                        |
| Full_Name                                                                                | One block for entering patient name.                                   |
| Last &First                                                                              | Separate blocks for entering patient last name and patient first name. |
| Auto Deletion of Transferred Queue                                                       |                                                                        |
| Choose to automatically delete the image from the queue after successful DICOM transfer. |                                                                        |
| Yes or No.                                                                               |                                                                        |



# Preset Program

## Overview

Preset Program in the Set Up top menu is 10 pages of preset parameters that are generally considered application parameters. They are exam category and preset name dependent.

Special features in the Preset Program Sub-Menus are:

- Annotation Library script designation
- Summary report display format editing
- User programmable OB table editing
- Calculation sequence programming
- Mode measurement menu customizing

| B              | Preset         | Set Up      | ECG         |
|----------------|----------------|-------------|-------------|
| Custom Display | Preset Program | Save Values | User Define |
| System Paramtr |                | Utility     | Diag.       |

Figure 14–12. Set Up Top Sub-Menus (Preset Program)

A brief description of each preset parameter and the choices available for that preset follows.



*NOTE: If it is necessary to quickly return to the Preset Program Sub-Menu page 1, press **Ctrl** and **1** simultaneously.*

*Refer to 14–6 for the key to the codes used to describe each parameter's affect.*



## Preset Program Contents




Figure 14–13. Preset Program Contents Page

The Preset Program Contents pages shows a summary of the 10 menu pages.

**Trackball** the cursor to the desired page and press **Set**.



## Preset Program

---

### Page 1 of 10 (Application)

|                                                                                                                                                                                                                                                                                                                              |  |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| <b>Body Pattern Probe Mark Preset    EC, UP</b>                                                                                                                                                                                                                                                                              |  |
| Choose to activate the preset position of the probe marker graphic.<br>On or Off.                                                                                                                                                                                                                                            |  |
| <b>Body Pattern Package    EC, UP</b>                                                                                                                                                                                                                                                                                        |  |
| Choose one of the 8 customized body pattern packages to be displayed with the current application preset.<br>Body pattern package 1, 2, 3, 4, 5, 6, 7 or 8.<br>Refer to 14–55 for details on programming the contents of each body pattern package.                                                                          |  |
| <b>Comment Erase at Display Change    EC</b>                                                                                                                                                                                                                                                                                 |  |
| Choose to erase all comments when a change is made in scan mode.<br>On or Off.                                                                                                                                                                                                                                               |  |
| <b>Comment Home Position Horizontal    EC, R</b>                                                                                                                                                                                                                                                                             |  |
| Choose a point on the horizontal axis for the comment cursor home position.<br>Enter 1 to 56 in 1 digit increments.                                                                                                                                                                                                          |  |
| <b>Comment Home Position Vertical    EC, R</b>                                                                                                                                                                                                                                                                               |  |
| Choose a point on the vertical axis for the comment cursor home position.<br>Enter 1 to 26 in 1 digit increments.                                                                                                                                                                                                            |  |
| <b>Comment Clear Key Function    EC</b>                                                                                                                                                                                                                                                                                      |  |
| Choose to erase all comments with the Clear key or just the line of comments that the cursor is on.<br><br>All        Erase all comments on the screen.<br>Line      Erase the comments on the cursor line. The second time Clear is pressed, all comments will be erased.                                                   |  |
| <b>Trackball Function at Freeze    EC</b>                                                                                                                                                                                                                                                                                    |  |
| Choose the function assigned to the Trackball when the image is frozen.<br><br>No                No function assigned.<br>Measure        Trackball assigned to the measurement function.<br>Comment        Trackball assigned to the comment function.<br>Body              Trackball assigned to the body pattern function. |  |



Page 2 of 10 (Application)

|                                                                                                                                                                                              |                                                         |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------|
| <b>VCR Rec Audio Doppler</b>                                                                                                                                                                 | <b>EC</b>                                               |
| Choose the format for audio recording in Doppler mode on the VCR.                                                                                                                            |                                                         |
| Stereo                                                                                                                                                                                       | Right/left channel separation.                          |
| Dop-L/Voice-R                                                                                                                                                                                | Doppler on left channel and voice on the right channel. |
| <b>VCR Play Audio</b>                                                                                                                                                                        | <b>EC</b>                                               |
| Choose how to playback VCR recorded audio.                                                                                                                                                   |                                                         |
| Stereo                                                                                                                                                                                       | Left and right separation.                              |
| Mono-L                                                                                                                                                                                       | Left speaker only.                                      |
| Mono-R                                                                                                                                                                                       | Right speaker only.                                     |
| <b>VCR Play Audio Volume [%]</b>                                                                                                                                                             | <b>EC, R</b>                                            |
| Choose the default percentage of volume level for VCR playback audio.                                                                                                                        |                                                         |
| Percentage from 0 to 100 in 4 digit increments.                                                                                                                                              |                                                         |
| <b>VCR Color Intensity</b>                                                                                                                                                                   | <b>EC</b>                                               |
| Choose the desired Color Intensity processing for VCR playback.                                                                                                                              |                                                         |
| 1                                                                                                                                                                                            | Optimized Chrominance Gain.                             |
| 2                                                                                                                                                                                            | Normal Processing.                                      |
| <b>Image Archive Compression</b>                                                                                                                                                             | <b>EC</b>                                               |
| Choose to compress the image archive data for the DEFF format.                                                                                                                               |                                                         |
| Yes or No.                                                                                                                                                                                   |                                                         |
| <i>NOTE: If image compression is used, images saved by the LOGIQ™ 500 cannot be read on other TIFF Readers, DEFF devices or the LOGIQ™ 700. However, they can be read on the LOGIQ™ 400.</i> |                                                         |



## Page 3 of 10 (Application – Measurement)

|                                                                    |                                                                |
|--------------------------------------------------------------------|----------------------------------------------------------------|
| <b>Time Unit</b>                                                   | <b>EC</b>                                                      |
| Choose the unit of display for all time measurements.              |                                                                |
| ms                                                                 | milliseconds.                                                  |
| sec                                                                | seconds.                                                       |
| <b>Length Unit</b>                                                 | <b>EC</b>                                                      |
| Choose the unit of display for all length (distance) measurements. |                                                                |
| mm                                                                 | millimeters.                                                   |
| cm                                                                 | centimeters.                                                   |
| <b>Velocity Unit</b>                                               | <b>EC, UP</b>                                                  |
| Choose the unit of display for all velocity measurements.          |                                                                |
| cm/s                                                               | centimeters per second.                                        |
| m/s                                                                | meters/second.                                                 |
| auto                                                               | displays in either cm/s or m/s depending on velocity measured. |
| <b>Circ/Area Method</b>                                            | <b>EC</b>                                                      |
| Choose the method to be used for circumference/area measurements.  |                                                                |
| Trace                                                              | Manually trace with the Trackball.                             |
| Ellipse                                                            | Ellipse method.                                                |
| Circle                                                             | Circle method.                                                 |
| 2 Dist                                                             | Two diameter (distance) method.                                |
| Diameter                                                           | Single distance method.                                        |
| <b>Echo Level Method</b>                                           | <b>EC</b>                                                      |
| Choose the size of the echo level measurement box cursor.          |                                                                |
| Box 3mm                                                            | 3mm x 3mm box.                                                 |
| Box 5mm                                                            | 5mm x 5mm box.                                                 |
| Box 7mm                                                            | 7mm x 7mm box.                                                 |
| <b>Heart Rate Method</b>                                           | <b>EC</b>                                                      |
| Choose the method for measuring heart rate.                        |                                                                |
| Auto                                                               | Automatic Calibration.                                         |
| Manual                                                             | Manual Calibration.                                            |
| <b>VCRPB Calib Method</b>                                          | <b>EC</b>                                                      |
| Choose the default method for VCR playback calibration.            |                                                                |
| Auto                                                               | Automatic Calibration.                                         |
| Manual                                                             | Manual Calibration.                                            |

Page 3 of 10 (Application – Measurement) (cont'd)

|                                                                                                                                                                                                                                                  |               |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|
| <b>Average Activity</b>                                                                                                                                                                                                                          | <b>EC</b>     |
| Choose to enable the calculation summary report page measurement averaging function.<br>On or Off.                                                                                                                                               |               |
| <b>Max Velocity (Advanced Vascular)</b>                                                                                                                                                                                                          | <b>EC</b>     |
| This selection is used with the Advanced Vascular option only. When turned on it overrides the Average Activity preset to cause the Advanced Vascular report page to display maximum velocity rather than latest or average.<br>On or Off.       |               |
| <b>Average Number</b>                                                                                                                                                                                                                            | <b>EC</b>     |
| Choose how many measurements will be used in the averaging function. The Cardiology exam category with AMCAL has a choice of 5 measurements whereas the other exam categories have a choice of 3 measurements.<br>1, 2, 3, 4 or 5.               |               |
| <b>Mask Image</b>                                                                                                                                                                                                                                | <b>EC, UP</b> |
| Choose to put the measurement/calculation results in a colored box over the image or just the graphic characters over the image.<br>Off Graphics are placed on top of the image.<br>On Graphics are placed in a colored box on top of the image. |               |
| <b>Trace Auto</b>                                                                                                                                                                                                                                | <b>EC</b>     |
| Choose to enable the automatic measurement of TAMAX, TAMIN, TAMEAN or TAMODE. See "D Realtime Trace Method" in Set Up/Custom Display page 13 for trace method selection on 14–37.<br>On or Off.                                                  |               |
| <b>Add 1 Week to EDD</b>                                                                                                                                                                                                                         | <b>EC</b>     |
| Choose to add an additional week to the displayed estimated date of delivery.<br>On or Off.                                                                                                                                                      |               |
| <b>Display EDD with GA Value</b>                                                                                                                                                                                                                 | <b>EC</b>     |
| Choose to display the estimated EDD with GA.<br>On or Off.                                                                                                                                                                                       |               |



**Page 3 of 10 (Application – Measurement) (cont'd)**

|                                                                                                                                                                                       |                                                 |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------|
| <b>Diastole/Systole Determination    EC</b>                                                                                                                                           |                                                 |
| Used with Advanced Cardiac Calculation option. Determine if measurement is systole or diastole automatically from the measurements relationship to the ECG wave.                      |                                                 |
| Manual                                                                                                                                                                                | Automatic determination disabled.               |
| Auto                                                                                                                                                                                  | Automatic determination enabled.                |
| <b>R Delay Time of End Systole Prior Edge    EC, UP</b>                                                                                                                               |                                                 |
| Used with automatic determination of systole/diastole in the Advanced Cardiac Calculation option. Sets the delay time on the ECG wave to ensure systole will be properly determined.  |                                                 |
| Use the Ellipse rocker switch to enter a number between 0 and 2000 msec.                                                                                                              |                                                 |
| <b>End Systole Period from Prior Edge    EC, UP</b>                                                                                                                                   |                                                 |
| Used with Automatic determination of Systole/Diastole in the Advanced Cardiac Calculation option. Sets the delay time on the ECG wave to ensure diastole will be properly determined. |                                                 |
| Use the Ellipse rocker switch to enter a number between 0 and 2000 msec.                                                                                                              |                                                 |
| <b>Media Selection for Fetal Trend    EC</b>                                                                                                                                          |                                                 |
| Used with the Fetal Growth Trend option. Sets the media to be used to store patient data. The hard drive has a much limited capacity.                                                 |                                                 |
| HD                                                                                                                                                                                    | System Hard Drive.                              |
| MO                                                                                                                                                                                    | Magnetic Optical Disk (removeable).             |
| <b>Measurement Clear Operation    EC</b>                                                                                                                                              |                                                 |
| When Clear is pressed, decide if it will erase measurements only or measurements and comments.                                                                                        |                                                 |
| Meas.-only                                                                                                                                                                            | Erase only measurements not comments.           |
| with-Comment                                                                                                                                                                          | Erase measurements and comments simultaneously. |

Page 4 of 10 (Application – Measurement)

|                                                                                                                                      |                                                                                               |
|--------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------|
| <b>Diastole Velocity for PI EC, UP</b>                                                                                               |                                                                                               |
| Choose which velocity value will be used to calculate end diastole.                                                                  |                                                                                               |
| Vmin                                                                                                                                 | Minimum value of the whole trace.                                                             |
| Vd                                                                                                                                   | Minimum value detected before Vmax (end diastole). Minimum value part of a heart rate period. |
| If Vd is not detected, use Vmin.                                                                                                     |                                                                                               |
| <b>Diastole Velocity for RI EC, UP</b>                                                                                               |                                                                                               |
| Choose which velocity value will be used to calculate end diastole.                                                                  |                                                                                               |
| Vmin                                                                                                                                 | Minimum value of the whole trace.                                                             |
| Vd                                                                                                                                   | Minimum value detected before Vmax (end diastole). Minimum value part of a heart rate period. |
| If Vd is not detected, use Vmin.                                                                                                     |                                                                                               |
| <b>Diastole Velocity for S/D (D/S) EC, UP</b>                                                                                        |                                                                                               |
| Choose which velocity value will be used to calculate end diastole.                                                                  |                                                                                               |
| Vmin                                                                                                                                 | Minimum value of the whole trace.                                                             |
| Vd                                                                                                                                   | Minimum value detected before Vmax (end diastole). Minimum value part of a heart rate period. |
| If Vd is not detected, use Vmin.                                                                                                     |                                                                                               |
| <b>Averaging Number for Doppler Realtime Calc EC, UP</b>                                                                             |                                                                                               |
| Choose how many heart rate periods will be averaged prior to being displayed in the realtime calculation window.                     |                                                                                               |
| 1, 2, 3, 4 or 5.                                                                                                                     |                                                                                               |
| <b>RI Calculation to bidirectional Flow EC, UP</b>                                                                                   |                                                                                               |
| When the flow direction for systolic velocity and diastolic velocity are different, set this parameter to on for the RI calculation. |                                                                                               |
| Off                                                                                                                                  | Bidirectional flow not considered in RI calculation.                                          |
| On                                                                                                                                   | Bidirectional flow considered in RI calculation.                                              |
| <b>Hip Orientation EC, UP</b>                                                                                                        |                                                                                               |
| Choose the type of Hip Orientation method to be used for scanning the infant hip.                                                    |                                                                                               |
| Cranial-left                                                                                                                         | Cranial structures are on the left; caudal structures are on the right.                       |
| Caudal-left                                                                                                                          | Caudal structures are on the left; cranial structures are on the right.                       |
| <b>Stepper Increment [mm] EC, UP</b>                                                                                                 |                                                                                               |
| Choose the stepper increment for area measurements in millimeters.                                                                   |                                                                                               |
| 5 or 2.5.                                                                                                                            |                                                                                               |



## Preset Program

### Page 4 of 10 (Application – Measurement)

|                                                                                                                                                        |                                                                                 |
|--------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------|
| <b>PPSA Coefficient (1) EC, UP</b>                                                                                                                     |                                                                                 |
| Set the desired default PPSA Coefficient for the first PPSA calculation on the Urology Report page.                                                    |                                                                                 |
| Coefficient from 0.01 to 0.99.                                                                                                                         |                                                                                 |
| <b>PPSA Coefficient (2) EC, UP</b>                                                                                                                     |                                                                                 |
| Set the desired default PPSA Coefficient for the second PPSA calculation on the Urology Report page.                                                   |                                                                                 |
| Coefficient from 0.01 to 0.99.                                                                                                                         |                                                                                 |
| <b>Display MaxPG with Velocity EC, UP</b>                                                                                                              |                                                                                 |
| Choose the default for displaying MaxPG with peak velocity.                                                                                            |                                                                                 |
| Off or On.                                                                                                                                             |                                                                                 |
| <b>Arrow Mark Display EC, UP</b>                                                                                                                       |                                                                                 |
| Choose the default for displaying the arrow mark when the Measurement key is pressed a certain number of times. See Table 7–1 on 7–5 for more details. |                                                                                 |
| Off or On.                                                                                                                                             |                                                                                 |
| <b>% Stenosis Calculation EC, UP</b>                                                                                                                   |                                                                                 |
| In Advanced Vascular, choose the default measurement type for calculating % Stenosis.                                                                  |                                                                                 |
| Area or Distance.                                                                                                                                      |                                                                                 |
| <b>CAUTION:</b> If this preset is changed, the previous % Stenosis measurement results will be erased.                                                 |                                                                                 |
| <b>B Auto Caliper Control EC, UP</b>                                                                                                                   |                                                                                 |
| Choose the location of the caliper for subsequent measurements while in B-Mode.                                                                        |                                                                                 |
| Off                                                                                                                                                    | The caliper is not displayed after a measurement is performed.                  |
| On                                                                                                                                                     | The caliper is in the home position.                                            |
| 1'st                                                                                                                                                   | The caliper is positioned at the previous measurement's first cursor position.  |
| 2'nd                                                                                                                                                   | The caliper is positioned at the previous measurement's second cursor position. |
| <b>M Auto Caliper Control EC, UP</b>                                                                                                                   |                                                                                 |
| Choose the location of the caliper for subsequent measurements while in M-Mode.                                                                        |                                                                                 |
| Off                                                                                                                                                    | The caliper is not displayed after a measurement is performed.                  |
| On                                                                                                                                                     | The caliper is in the home position.                                            |
| 1'st                                                                                                                                                   | The caliper is positioned at the previous measurement's first cursor position.  |
| 2'nd                                                                                                                                                   | The caliper is positioned at the previous measurement's second cursor position. |
| <b>D Auto Caliper Control EC, UP</b>                                                                                                                   |                                                                                 |
| Choose the location of the caliper for subsequent measurements while in Doppler Mode.                                                                  |                                                                                 |
| Off                                                                                                                                                    | The caliper is not displayed after a measurement is performed.                  |
| On                                                                                                                                                     | The caliper is in the home position.                                            |
| 1'st                                                                                                                                                   | The caliper is positioned at the previous measurement's first cursor position.  |
| 2'nd                                                                                                                                                   | The caliper is positioned at the previous measurement's second cursor position. |





Page 5 & 6 of 10 (Application – Measurement Sub-Menu)

| Cardiac Measurement Choices    EC, UP                                                                                                  |                                                 |
|----------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------|
| Used with Advanced Cardiac Calculation option. Determine if the measurement will be used in the automatic M-Mode calculation sequence. |                                                 |
| Off                                                                                                                                    | Automatic measurement callout will NOT be made. |
| On                                                                                                                                     | Automatic measurement callout WILL be made.     |



## Page 7 of 10 (Application – Measurement Submenu)

| Report Format                                                                                                                                                                                                                               | EC |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----|
| If the current exam category has an associated summary report page, select this menu option to change the vertical (y) position of the measurement elements on the report page. The elements that can be edited depend on each report page. |    |

### Advanced Report Formatting

The LOGIQ™ 500 allows for user control of each exam category summary report measurement display format.

Except for the protected 5 lines of each report page (that has hospital and patient information) and the bottom area for calculations and comments, the vertical position (“y”) of the measurements can be customized by the user.

Factory default items can be deleted or user desired items added.



Example report format change

Exam Category: OB

OB measurements are limited to changes in the vertical position (“y”) on the report page.

To change the position of a measurement:



From the Set Up/Preset Program Menu page 7, use the **Trackball** to place the arrow cursor over Report Format and press **Set**. The report format editor screen appears.

GE Medical Systems98/05/12 14:12:24  
C364

REPORT FORMAT EDITOR FOR OB USA

| MEASUREMENTS  |      |               |      |
|---------------|------|---------------|------|
| NAME          | Ypos | NAME          | Ypos |
| BPD (HADLOCK) | 6    | HC (HADLOCK)  | 7    |
| OFD (HC)      | 8    | AC (HADLOCK)  | 9    |
| FL (HADLOCK)  | 10   | CRL (HADLOCK) | 11   |
| GS (HELLMAN)  | 12   |               |      |

| CALCULATIONS |      |      |           |      |      |
|--------------|------|------|-----------|------|------|
| NAME         | Xpos | Ypos | NAME      | Xpos | Ypos |
| CI           | 1    | 20   | EFW       | 29   | 20   |
| FL/BPD       | 1    | 21   | SD        | 39   | 20   |
| FL/AC        | 1    | 22   | SD        | 29   | 20   |
| FL/HC        | 1    | 23   | LMP (GA)  | 26   | 23   |
| HC/AC        | 1    | 24   | GA (LMP)  | 32   | 24   |
| CUA          | 46   | 24   | EDD (LMP) | 32   | 25   |
| EDD (CUA)    | 46   | 25   |           |      |      |

SAVE

EXIT

'TRACKBALL' and 'SET' are available.

Figure 14–14. Report Format Editor



Use the **Trackball** or **Return** to position the highlighted area over the “Ypos” of the measurement to be changed.

Type the new position number.



Use the **Trackball** to highlight any additional positions to be changed.


(continued)




## Preset Program

---

### Example report format change (cont'd)



When all changes are complete, highlight Save and press **Set**.



To exit the Report Format Editor and return to the Set Up/  
Preset Program Menu, highlight Exit and press **Set**.

To Add or Delete a measurement:

Enter the Report Format Editor in the same manner as “To  
change the position”.

Use the **Trackball** to highlight the desired measurement name  
(to delete or change) or highlight a blank area (to add).



To delete a name, type in all spaces.


To change a name, type over the current name.

To add a name, type in the new name.




*NOTE: When adding or changing a measurement name, it must  
contain 2 important elements:*

- *Measurement type. For example, GS, BPD, CRL, HC.*
- *Author's name. The author's name must appear in  
parentheses, for example (Hadlock), (Hellman).*



Highlight Save and press **Set** to record all changes.



Highlight Exit and press **Set** to return to the Set Up/Preset  
Program menu.



*NOTE: If 2 measurement names are placed in the same  
position, an error message appears.*

*“More than 1 Measurement Field Name at the SAME Ypos”.*



Page 7 of 10 (Application – Measurement Submenu) (cont'd)

User Table Editor

If the factory preset OB tables are not suitable, the user can input data for five OB tables. The titles given to these tables appear as selections in the OB Calculations Soft-Menu when displayed.


To program a user table

Select the **Set Up** Top Menu.

| B              | Preset         | Set Up      | ECG         |
|----------------|----------------|-------------|-------------|
| Custom Display | Preset Program | Save Values | User Define |
| System Paramtr |                | Utility     | Diag.       |

Figure 14–15. Set Up Top Menu

Select **Preset Program** from the Set Up Sub-Menu. The preset program selections are displayed over the image on the monitor.



Use the **Trackball** to place the arrow cursor over User Table Editor.

Press **Set** to display the User Table Programming Main Menu.

The User Table Programming main menu appears as shown in Figure 14–16.



## Preset Program

### To program a user table (cont'd)

99/03/17 16:29:42  
C358

USER TABLE PROGRAMMING

Table Select : **UT01**

Dimension : **Dist (mm)**  
Circ (mm)  
Area (mm<sup>2</sup>)


Standard Dev : SD\_Age

Title : UT01  
Author Name :  
Note :  
Min. Data : 1  
Interval : 1

TABLE\_EDIT                      EXIT

'TRACKBALL' and 'SET' are available.

Figure 14–16. User Table Programming



Use the **Ellipse** rocker switch to select the desired table number from UT01 to UT10.

Press **Set** to select the desired table.

Use the **Trackball** to place the arrow cursor on the desired table dimension. The choices are Dist (mm) for distance in millimeters, Circ (mm) for circumference in millimeters and Area for area in square millimeters. Press **Set** when the cursor is over the correct dimension.

Use the **Trackball** to move the arrow cursor to the title block. Type in the desired table title (maximum 4 characters). Press **Return**.

Type in the author's name (maximum 8 characters). Press **Return**.



## Preset Program

---

### To program a user table (cont'd)

Type in any desired notes (maximum 30 characters). Press **Return**.

Type in the minimum or beginning data value for the table (maximum 4 digits). Press **Return**.

Type in the interval value between each table data point (maximum 4 digits). Press **Return**.

The selection "Table-Edit" should now be highlighted. Press **Set**.



*NOTE: The **Return** key or **Trackball** can be used to cycle through the menu selections. To make a choice, press **Set**.*



Preset Program

To enter table values

The Table-Edit Menu shown in Figure 14–17 allows for adding data, editing data or clearing (deleting) a User Table.

|                    |        |                        |                   |         |      |
|--------------------|--------|------------------------|-------------------|---------|------|
| GE Medical Systems |        |                        | 98/05/12 14:12:58 |         |      |
|                    |        |                        | C364              |         |      |
| TABLE : TBL5       |        | UNIT : mm <sup>2</sup> | PAGE : 1/2        |         |      |
| 1                  | 0.0000 | 21                     | W D±              | 41      | W D± |
| 2                  | W D±   | 22                     | W D±              | 42      | W D± |
| 3                  | W D±   | 23                     | W D±              | 43      | W D± |
| 4                  | W D±   | 24                     | W D±              | 44      | W D± |
| 5                  | W D±   | 25                     | W D±              | 45      | W D± |
| 6                  | W D±   | 26                     | W D±              | 46      | W D± |
| 7                  | W D±   | 27                     | W D±              | 47      | W D± |
| 8                  | W D±   | 28                     | W D±              | 48      | W D± |
| 9                  | W D±   | 29                     | W D±              | 49      | W D± |
| 10                 | W D±   | 30                     | W D±              | 50      | W D± |
| 11                 | W D±   | 31                     | W D±              | 51      | W D± |
| 12                 | W D±   | 32                     | W D±              | 52      | W D± |
| 13                 | W D±   | 33                     | W D±              | 53      | W D± |
| 14                 | W D±   | 34                     | W D±              | 54      | W D± |
| 15                 | W D±   | 35                     | W D±              | 55      | W D± |
| 16                 | W D±   | 36                     | W D±              | 56      | W D± |
| 17                 | W D±   | 37                     | W D±              | 57      | W D± |
| 18                 | W D±   | 38                     | W D±              | 58      | W D± |
| 19                 | W D±   | 39                     | W D±              | 59      | W D± |
| 20                 | W D±   | 40                     | W D±              | 60      | W D± |
| SAVE               |        | CLEAR                  | SETUP             | NEXT-PG | EXIT |

Figure 14–17. Table Edit Menu

For each data point, the weeks (2 digits), days (1 digit) and deviation (2 digits) must be entered.

The menu choices at the bottom allow:

- SAVE
- Saves the table in its current form.
- CLEAR
- Clears or deletes the entire table (and title).
- SETUP
- Returns to the User Table Programming Setup menu.
- NEXT-PG
- Goes to the next page to continue data entry.
- EXIT
- Exits the User Programmable Table function.

**Trackball** to one of the above choices and press **Set**. When exit is selected without doing a “Save”, all additions or edits are lost. The table reverts to its previous state.



*NOTE: If a new table was added, the measurement can be added to the summary report. Follow the instructions under Report Format in the Set Up/Preset program.*



Page 7 of 10 (Application – Measurement Submenu) (cont'd)


**Sequence 1–8 EC, UP (Auto Sequence Programming)**

The user can program a maximum of 8 calculation sequences with a maximum of eight measurements per sequence.

These are labeled Sequence 1–8 in the Set Up/Preset Program Menu as USERSEQUENCE1–8 in the Measurement Sub-Menu or in the Auto Sequence Top Menu.

When Sequence 1–8 is chosen from the Auto Sequences Top Menu or the Measurement Sub-Menu, the system automatically prompts the user through the calculations previously programmed.

**Program a User Sequence**



Use the **Trackball** to move the arrow cursor to the desired sequence number (1–8) on page 7 of the Set Up/Preset Program Menus.

Press **Set**. The first page of the 8 Measurement Sequence Sub-Menus is displayed.

```

07/30/97 08:40:04
S317
USER PARAM      SYSTEM PARAM  PRESET PROG  CUSTOM DISP
CATEGORY : OBSTETRICS  PRESET NAME : General  OB
PAGE : 7/10      PRIOR  NEXT  CONTENTS
COMMAND : $AVE  RESET  DELETE  EXIT

***** Measurement Sequence Submenu *****
USERSEQUENCE1  1  UNDEFINED  NEXT CALC
UNDEFINED      OB CTAR  OB EF      OB PLI
OB Vmax        OB S/D  OB D/S     OB RI
OB PI          AFI     USA  OB USER1 USA  OB USER2
USA  OB USER3 USA  OB USER4 USA  OB USER5 EUROP OB USER1
EUROP OB USER2 EUROP OB USER3 EUROP OB USER4 EUROP OB USER5
TOKYO OB USER1 TOKYO OB USER2 TOKYO OB USER3 TOKYO OB USER4
TOKYO OB USER5 OSAKA OB USER1 OSAKA OB USER2 OSAKA OB USER3
OSAKA OB USER4 OSAKA OB USER5 GS HELLMAN  GS HANSMANN
GS TOKYO      GS REMPEN  CRL HADLOCK  CRL JEANTY
CRL CAMPBEL   CRL HANSMANN CRL ROBINSN  CRL PARIS
CRL NELSON    CRL TOKYO   CRL TOKYO-S  CRL OSAKA
CRL REMPEN    CRL ASUM   BPD HADLOCK  BPD JEANTY
BPD CAMPBEL   BPD HANSMANN BPD PARIS    BPD SOSTOA
BPD KURTZ     BPD TOKYO   BPD TOKYO-S  BPD OSAKA
BPD REMPEN    BPD ASUM   FL HADLOCK   FL JEANTY
FL CAMPBEL    FL HANSMANN FL PARIS     FL SOSTOA
FL TOKYO      FL TOKYO-S  FL OSAKA     AC HADLOCK
AC JEANTY     AC HANSMANN AC SOSTOA    AC TOKYO-S
AC ASUM       HC HADLOCK  HC JEANTY    HC HANSMANN
'TRACKBALL' and 'SET' are available.
    
```

Figure 14–18. Measurement Sequence Sub-Menu



**NOTE:** The first row under the sub-menu title line shows 3 important items:



## Preset Program

### Program a User Sequence (cont'd)




Figure 14-19. Sequence Sub-Menu Status Line

**8** is the order in which the measurement or calculation will be made in USERSEQUENCE 2. 1 is first, 2 is second.....8 is last.

UNDEFINED is the current measurement in that sequence location. Undefined means it is blank.

NEXTCALC allows the user to cycle through the 8 pages of measurements and calculations available to be programmed.

Use the **Ellipse** rocker switch to cycle to the desired usersequence measurement number (1-8).



Use the **Trackball** to move the arrow cursor to the desired measurement/calculation.

Move the arrow cursor to Next Calc and press **Set** if a sub-menu page change is necessary.

Press **Set** when the arrow cursor is over the desired measurement/calculation.

The newly selected measurement/calculation appears as the current measurement for the highlighted USERSEQUENCE number replacing undefined.

Repeat the previous steps as necessary to program additional sequence number locations for any of the 8 user sequences.



## Program a User Sequence (cont'd)



### Hints

USERSEQUENCE1–8 are choices that can be programmed in the Measurement Menus for each mode (B, M, D, and CFM) per exam category.



Auto Sequence Top Menu choices are exam category dependent. Sequences programmed in Vascular are not available when the Cardiac exam category is selected.

UNDEFINED means no measurement/calculation is in that location.

Use UNDEFINED to delete a desired measurement.

To customize the Usersequence/Auto Sequence name:

Move the **Trackball** to the USERSEQUENCE# at the top of one of the 8 Calculation menus.



Press **Set**.

Type in the name of the Usersequence. The space available is 2 lines of 7 characters each.

Move the **Trackball** to SAVE.

Press **Set**.

## Preset Program

---

### Performing an Auto Sequence (USERSEQUENCE)

Select an Auto Sequence, or a user programmable mode measurement menu selection, the system automatically prompts the user to do the first measurement in the sequence.

After the first measurement/calculation is completed, the system automatically prompts the user to do the second measurement in the sequence.

This continues until all measurements are completed.

### Sequence Modification

During the Auto Sequence function:



Press the top of the **Ellipse** rocker switch to skip a desired measurement and advance to the next measurement in the sequence.

Press the bottom of the **Ellipse** rocker switch to skip back to a previous measurement in the sequence.

Measurements skipped will appear again at the end of the sequence. This gives the user a second chance to make the measurements prior to exiting the Auto Sequence function.

Page 7 of 10 (Application – Measurement Submenu) (cont'd)

| Measurement Menu B/M/D/CFM EC, UP                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Each exam category has a set of factory default measurements/calculations. The 4 measurement menu modes allow the user to add a maximum of 4 new pages to the measurement sub-menus per exam category per mode within a user application preset. With 8 items per page, a total of 32 measurements/calculations can be added for each exam category in B, M, D, which is distinct and separate from those that can be programmed into CFM-Mode.</p> <p><i>NOTE: Measurements and Calculations setup for Measurement Menu B are not transferred to Menu M, D, or CFM. Each one is separate.</i></p> |

Organization

There are 4 pages of space available for the user to add measurements and calculations to an exam category measurement menu.

These pages do not appear in the Soft Menu display unless items are placed in the number location of that page.

|   |   |   |   |
|---|---|---|---|
|   |   |   |   |
| 1 | 3 | 5 | 7 |
| 2 | 4 | 6 | 8 |

Sub Menu page 1 of 4

|    |    |    |    |
|----|----|----|----|
|    |    |    |    |
| 9  | 11 | 13 | 15 |
| 10 | 12 | 14 | 16 |

Sub Menu page 2 of 4

|    |    |    |    |
|----|----|----|----|
|    |    |    |    |
| 17 | 19 | 21 | 23 |
| 18 | 20 | 22 | 24 |

Sub Menu page 3 of 4

|    |    |    |    |
|----|----|----|----|
|    |    |    |    |
| 25 | 27 | 29 | 31 |
| 26 | 28 | 30 | 32 |


Sub Menu page 4 of 4

Figure 14–20. Measurement Menu Key Numbers

The spaces available for calculations are called Key Numbers. Key Numbers 1 through 8 are on page 1, 9 through 16 on page 2 etc. (See Figure 14–20).

## Preset Program

### Programming Measurements



Use the **Trackball** to move the arrow cursor to the desired Measurement Menu (B, M, D, or CFM) on page 7 of the Set Up/Preset Program menus.

Press **Set**. The last page viewed of the 4 measurement pages is displayed.

```
97/09/19 14:10:03
C364
USER PARAM      SYSTEM PARAM  PRESET PROG  CUSTOM DISP
CATEGORY : RAD/ABDOMEN  PRESET NAME : 1
PAGE : 7/10      PRIOR      NEXT      CONTENTS
COMMAND : SAVE    RESET    DELETE    EXIT

***** Measurement VFD Softmenu Setup Submenu *****
KEY NUMBER      1      UNDEFINED      NEXT CALC
UNDEFINED      USERSEQUENCE1  USERSEQUENCE2  USERSEQUENCE3
USERSEQUENCE4  USERSEQUENCE5  USERSEQUENCE6  USERSEQUENCE7
USERSEQUENCE8  TRACE AUTO      ANATOMY SURVEY  OB GRAPH
OB REPORT      GYN REPORT      VASC REPORT     LV REPORT
CARD REPORT    Venous Report   Carotid Report  OB CTAR
OB EF          OB PLI          OB Vmax          OB S/D
OB D/S         OB RI          OB PI            AFI
USA OB USER1  USA OB USER2  USA OB USER3  USA OB USER4
USA OB USER5  EUROP OB USER1  EUROP OB USER2  EUROP OB USER3
EUROP OB USER4  EUROP OB USER5  TOKYO OB USER1  TOKYO OB USER2
TOKYO OB USER3  TOKYO OB USER4  TOKYO OB USER5  OSAKA OB USER1
OSAKA OB USER2  OSAKA OB USER3  OSAKA OB USER4  OSAKA OB USER5
GS HELLMAN      GS HANSMANN      GS TOKYO        GS REMPEN
CRL HADLOCK     CRL JEANTY       CRL CAMPBEL     CRL HANSMANN
CRL ROBINSN     CRL PARIS        CRL NELSON      CRL TOKYO
CRL TOKYO-S     CRL OSAKA        CRL REMPEN      CRL ASUM
BPD HADLOCK     BPD JEANTY       BPD CAMPBEL     BPD HANSMANN
BPD PARIS       BPD SOSTOA       BPD KURTZ       BPD TOKYO
BPD TOKYO-S     BPD OSAKA        BPD REMPEN      BPD ASUM
'TRACKBALL' and 'SET' are available.
```

Figure 14–21. Measurement Setup Sub-Menu



**NOTE:** The first row under the sub-menu title line shows 3 important items.



## Programming Measurements (cont'd)

```
***** Measurement VFD Softmenu Setup Submenu *****  
KEY NUMBER    1 USERSEQUENCE1 NEXT CALC
```

Figure 14–22. Sequence Sub-Menu Status Line

**1** is the position in the soft menu where the measurement for key number one will be displayed.

USERSEQUENCE1, or the customized name, is the current measurement sequence in that key number location.

NEXTCALC allows the user to cycle through the four pages of measurements and calculations available.

Remembering how the key numbers are organized (see Figure 14–20), use the **Ellipse** rocker switch to change to the desired key number.

Move the arrow cursor to NEXTCALC and press **Set** if a sub-menu page change is necessary.

Press **Set** when the arrow cursor is over the desired measurement/calculation.

The newly selected measurement/calculation appears in the highlighted key number location.



*NOTE: If key number is desired to be empty (blank, deleted) choose UNDEFINED.*

*Factory default measurements/calculations are protected and cannot be changed.*

*Adding authors to a measurement menu does NOT automatically add the measurement to the report page.*

Repeat the previous steps to program any key number in any mode menu. Be sure to use the Save function prior to exiting the Calculation menus and the Preset Program menus.



Page 7 of 10 (Application – Measurement Submenu) (cont'd)

|                                                                                                                                                                                                                |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Cardiac Calculation    EC, UP</b>                                                                                                                                                                           |
| This sub menu can be used to rearrange the order of the automatic cardiac measurement prompts for the Advanced Cardiac Calculation option.<br>Refer to 10–63 for details on measurement sequence modification. |
| <b>D Realtime Calc.    EC, UP</b>                                                                                                                                                                              |
| This sub menu is used to select and arrange the parameters that are displayed in the D Realtime Calc. window in the upper left of the monitor.                                                                 |

|                                                  |              |                       |                |
|--------------------------------------------------|--------------|-----------------------|----------------|
| 07/28/97 11:28:11                                |              |                       |                |
| C551                                             |              |                       |                |
| USER PARAM                                       | SYSTEM PARAM | PRESET PROG           | CUSTOM DISP    |
| CATEGORY : OBSTETRICS                            |              | PRESET NAME : General | OB             |
| PAGE : 7/10                                      |              | PRIOR                 | NEXT           |
| COMMAND : SAVE    RESET                          |              | DELETE                | EXIT           |
| ***** Realtime Doppler Calculation Submenu ***** |              |                       |                |
| KEY NUMBER                                       | 1            | Vmax(PEAK)            |                |
| UNDEFINED                                        | Vmax(PEAK)   | Vmin(PEAK)            | Vd(PEAK)       |
| TAMAX                                            | TAMEAN       | TAMOD                 | TAMIN          |
| HR (FFT)                                         | GEN PI       | GEN RI                | GEN S/D        |
| OB PI                                            | OB RI        | OB S/D                | OB D/S         |
| VASC PI                                          | VASC RI      | VASC S/D              | AVERAGE NUMBER |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| UNDEFINED                                        | UNDEFINED    | UNDEFINED             | UNDEFINED      |
| 'TRACKBALL' and 'SET' are available.             |              |                       |                |

Figure 14–23. Realtime Doppler Calculation Sub-Menu

A maximum of 8 lines (parameters) can be displayed.

|                                                  |                                    |
|--------------------------------------------------|------------------------------------|
| ***** Realtime Doppler Calculation Submenu ***** |                                    |
| KEY NUMBER                                       | 1                      Vmax (PEAK) |

Figure 14–24. Realtime Doppler Calculation Sub-Menu Status Line





Page 7 of 10 (Application – Measurement Submenu) (cont'd)

1 is the position number or line of the Doppler Realtime display.

Use the **Ellipse** rocker switch to change to the desired line number (1–8).

Move the arrow cursor to the desired measurement/calculation to be displayed on that line.

Press **Set** when the arrow cursor is over the desired measurement/calculation.

The newly selected measurement/calculation appears to the right of the position number.



*NOTE: Remember that this setup is Exam Category and User Preset dependent. If changes are desired to the factory defaults, a new User Preset will be created.*



## Preset Program

---

### Page 8 & 9 of 10 (Application – Annotation Library)

|                                                                                                                                                                                                                                                                                                                                                                                                                                                  |               |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|
| <b>Annotation Library 1–24</b>                                                                                                                                                                                                                                                                                                                                                                                                                   | <b>EC, UP</b> |
| <p>Type in the desired (most commonly used) display annotations. There is space for a maximum of 20 characters (numbers, letters and spaces) for each library script. A maximum of 24 annotation library scripts is available per exam category.</p> <p>Simply place the arrow cursor across from the desired annotation library number and press <b>Set</b>.</p> <p>Type the desired script.</p> <p>Press <b>Set</b> to complete the entry.</p> |               |

### Page 10 of 10 (Application – Patient Information)

|                                                                                                                                                                                                                                                                                                                                                                                                                      |                        |      |                       |       |                        |           |                       |        |                      |        |             |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------|------|-----------------------|-------|------------------------|-----------|-----------------------|--------|----------------------|--------|-------------|
| <b>New Patient Y/N Ask</b>                                                                                                                                                                                                                                                                                                                                                                                           | <b>EC</b>              |      |                       |       |                        |           |                       |        |                      |        |             |
| <p>In the new patient entry menu, choose to ask to erase all current patient data.</p> <p>On or Off.</p>                                                                                                                                                                                                                                                                                                             |                        |      |                       |       |                        |           |                       |        |                      |        |             |
| <b>Display Unit Age 1–4</b>                                                                                                                                                                                                                                                                                                                                                                                          | <b>EC</b>              |      |                       |       |                        |           |                       |        |                      |        |             |
| <p>Choose the priority for the units to be used to display patient age. Age 1 is the highest priority. Age 4 is the lowest priority.</p> <table><tr><td>year</td><td>Display age in years.</td></tr><tr><td>month</td><td>Display age in months.</td></tr><tr><td>week</td><td>Display age in weeks.</td></tr><tr><td>day</td><td>Display age in days.</td></tr><tr><td>nodisp</td><td>No display.</td></tr></table> |                        | year | Display age in years. | month | Display age in months. | week      | Display age in weeks. | day    | Display age in days. | nodisp | No display. |
| year                                                                                                                                                                                                                                                                                                                                                                                                                 | Display age in years.  |      |                       |       |                        |           |                       |        |                      |        |             |
| month                                                                                                                                                                                                                                                                                                                                                                                                                | Display age in months. |      |                       |       |                        |           |                       |        |                      |        |             |
| week                                                                                                                                                                                                                                                                                                                                                                                                                 | Display age in weeks.  |      |                       |       |                        |           |                       |        |                      |        |             |
| day                                                                                                                                                                                                                                                                                                                                                                                                                  | Display age in days.   |      |                       |       |                        |           |                       |        |                      |        |             |
| nodisp                                                                                                                                                                                                                                                                                                                                                                                                               | No display.            |      |                       |       |                        |           |                       |        |                      |        |             |
| <b>Display Unit Height 1–3</b>                                                                                                                                                                                                                                                                                                                                                                                       | <b>EC</b>              |      |                       |       |                        |           |                       |        |                      |        |             |
| <p>Choose the priority to be used for the units to display the patient's height.</p> <table><tr><td>cm</td><td>Centimeters</td></tr><tr><td>inch</td><td>Inches</td></tr><tr><td>feet-inch</td><td>Feet and inches</td></tr><tr><td>nodisp</td><td>No display</td></tr></table>                                                                                                                                      |                        | cm   | Centimeters           | inch  | Inches                 | feet-inch | Feet and inches       | nodisp | No display           |        |             |
| cm                                                                                                                                                                                                                                                                                                                                                                                                                   | Centimeters            |      |                       |       |                        |           |                       |        |                      |        |             |
| inch                                                                                                                                                                                                                                                                                                                                                                                                                 | Inches                 |      |                       |       |                        |           |                       |        |                      |        |             |
| feet-inch                                                                                                                                                                                                                                                                                                                                                                                                            | Feet and inches        |      |                       |       |                        |           |                       |        |                      |        |             |
| nodisp                                                                                                                                                                                                                                                                                                                                                                                                               | No display             |      |                       |       |                        |           |                       |        |                      |        |             |
| <b>Display Unit Weight 1–2</b>                                                                                                                                                                                                                                                                                                                                                                                       | <b>EC</b>              |      |                       |       |                        |           |                       |        |                      |        |             |
| <p>Choose the priority to be used for the units to display the patient's weight.</p> <table><tr><td>kg</td><td>Kilograms</td></tr><tr><td>lbs</td><td>Pounds</td></tr><tr><td>none</td><td>None</td></tr></table>                                                                                                                                                                                                    |                        | kg   | Kilograms             | lbs   | Pounds                 | none      | None                  |        |                      |        |             |
| kg                                                                                                                                                                                                                                                                                                                                                                                                                   | Kilograms              |      |                       |       |                        |           |                       |        |                      |        |             |
| lbs                                                                                                                                                                                                                                                                                                                                                                                                                  | Pounds                 |      |                       |       |                        |           |                       |        |                      |        |             |
| none                                                                                                                                                                                                                                                                                                                                                                                                                 | None                   |      |                       |       |                        |           |                       |        |                      |        |             |



# Save Values

## Overview

Save Values is a selection in the Set Up Menu that allows for the quick retention of a current system settings and preset parameters.

It works like the Save command in the Custom Display, System Parameter or Preset Program sub-menus.

| B              | Preset         | Set Up      | ECG         |
|----------------|----------------|-------------|-------------|
| Custom Display | Preset Program | Save Values | User Define |
| System Paramtr |                | Utility     | Diag.       |

Figure 14–25. Set Up Sub-Menu

## Saving Scan Values

If at any time during a scan it is felt that the system settings are optimized and should be saved as an Exam Applications Preset, simply select **Save Values**.

The system will start the process of saving all of the current settings as a User Exam Applications Preset.

If a factory default applications preset is active when Save Values is selected, the message:

“Input User Preset Number (1-8)”

appears at the bottom of the screen. Enter the preset position number 1 through 8. The data is saved and available for recall at the designated preset location.



## *Save Values*

---

### **Saving Scan Values (cont'd)**

If a User Preset is active when Save Values is selected, the message:

“Overwrite existing data ‘Y’ or ‘N’”

appears at the bottom of the screen.

Choose ‘Y’ to overwrite data in the current preset location or ‘N’ to designate a new preset location 1 through 8.



# Exam Applications Presets

## Overview

The top menu Preset selection will display a maximum of eight factory application preset choices for each exam category.

The user may define a maximum of eight application presets for display on the Preset sub-menu designated as page one.

Preset parameters that are exam category dependent and available for customization are found in the Preset Program and Custom Display menus. Some parameters are also probe dependent within the application preset.

Items found under System Parameters affect each exam category, application preset or probe.

## Defining a User Preset

User presets appear on page one of the Preset Sub-Menu. They will appear on the Soft Menu Display, as shown in Figure 14–26.

| B             | Preset        | Set Up        | ECG           |
|---------------|---------------|---------------|---------------|
| User Preset-1 | User Preset-3 | User Preset-5 | User Preset-7 |
| User Preset-2 | User Preset-4 | User Preset-6 | User Preset-8 |

Figure 14–26. User Preset Configuration



## Defining a User Preset (cont'd)

To define a User Preset (1–8), first ensure that:

- the system is set to the proper exam category by using **New Patient** or **ID/Name** keys
- the desired probe is connected and activated.
- all keyboard controls (mode, depth, gain, etc.) are set to the desired values.


| B              | Preset         | Set Up      | ECG         |
|----------------|----------------|-------------|-------------|
| Custom Display | Preset Program | Save Values | User Define |
| System Paramtr |                | Utility     | Diag.       |

Figure 14–27. Set Up Top Menu Selection

Select the **Set Up** Top Menu. Select Custom Display or Preset Program.

Make all desired changes to the parameters in the Custom Display and Preset Parameter Sub Menus.

Use the **Trackball** to move the arrow cursor to Save.



(continued)

Press **Set**. The message “Input User Preset Number (1–8)” appears at the bottom of the monitor.

Enter the desired preset number.

The Preset Name at the top of the menu now displays the number entered (1–8).

## Defining a User Preset (cont'd)



*NOTE: If a user application preset was used at the start of the process, the system prompt states:*


*"Overwrite Existing Data? 'y' or 'n'"*

*Enter 'y' to write over the existing user preset.*

*Enter 'n' to get the message: "Input User Preset Number (1-8)" and create a new user preset.*



Use the **Trackball** to move the arrow cursor over the Preset Name area.




Press **Set**.



*NOTE: If new data is not saved, all changes will be lost when another preset or New Patient is selected.*

## Naming a User Preset



To give the user preset a name instead of a number:

While in a Custom Display or Preset Program menu page, use the **Trackball** to move the arrow cursor over the Preset Name area. Press **Set**.

Type in the desired name for the User Preset in the designated exam category.

Use the **Trackball** to move the arrow cursor to SAVE and press **Set**.

Answer Yes to the question “ Overwrite existing data? ‘Y’ or ‘N’ ” to save the new User Preset name.

**IMPORTANT NOTE:** To ensure that all parameters and the user preset name are saved by the LOGIQ™ 500 before exiting the Set Up function:

- select System Parameters, any menu pages 1/7 through 7/7.
- **Trackball** to the SAVE function in the System Parameters menu and press **Set**.
- **IMPORTANT:** Answer ‘Y’ to overwrite existing data again.

*These steps are important to ensure that the new user preset name is saved.*





## **Deleting User Presets and Names**

Display any Set Up/Custom Display or Preset Program sub menu page.

- Use the **Trackball** to move the arrow cursor to the preset name and press **Set**.
- Type in the User Preset position number (1–8) and use the **Space Bar** to delete the rest of the name.
- Use the **Trackball** to move the arrow cursor to Delete menu selection and press **Set**.
- Answer Yes to the question. “ Delete data? ‘Y’ or ‘N’ ”.
- Select System Parameters, any menu page 1–7.
- **Trackball** to the SAVE function in the System Parameters menu and press **Set**.
- Answer Yes to the question “Overwrite existing data? ‘Y’ or ‘N’” to save the new User Preset name.

The last three steps are important to ensure that the User Preset name is cleared.

## **Recall Preset**

Recall Preset is a selection found on page 3 of the Preset Sub-Menu. Recall Preset can be used during an examination to return the majority of the parameters to the latest saved preset values.


If front panel and soft menu values were changed while scanning, Recall Preset quickly resets these values.

Recall Preset does not recover preset data that was changed with the Save function or was deleted.



# User Define Function

## Overview



Six keys in the middle of the right side of the console keyboard as well as alphanumeric (A/N) keys A~Z 1~0 may be programmed to duplicate a single keystroke or a sequence of 64 keystrokes maximum.

A keystroke is considered to be a press of a key, one movement of a knob or movement of the Trackball.

## Programming the User Define Function

- Access the User Define function by pressing the **Set Up** Top Menu Select key.
- Select User Define from the Set Up Sub-Menu.
- Perform the keystroke or sequence of keystrokes to be programmed.
- When complete, press the desired User Define key number or A/N key to record the programming.
- The system will ask: "Save data? 'Y' or 'N'." Press 'Y' on the keyboard to save the user define sequence. Press 'N' to cancel the user define programming function.
- If a sequence is already assigned to the selected key, the message will be displayed:

"Overwrite Existing data? 'y' or 'n'"

User Define keys that have been programmed will be half lighted (back lit).



## Programming the User Define Function (cont'd)

- When 'y' is the replay for Save Data or Overwrite Existing data,

"Input Name: " is displayed

A user define function name must be entered with a maximum of 11 characters.

- Press the Mode Top Menu Select function to return to the previous scan mode, if necessary.

The next time a programmed User Define or A/N key is pressed, the recorded keystroke or keystroke sequence will be performed by a single key operation.



### Hints

#### Ctrl, 1

If selecting Set Up Top Menu with Custom Display, System Parameter or Preset Program as part of the User Define Sequence, press **Ctrl** and **1** simultaneously to set the sub-menu page to one. Press **Set** the appropriate number of times to ensure proper movement to the desired sub-menu page.

#### Blue Shift, S

A pause can be valuable in a User Define or A/N key sequence, especially if it is desired to program a key to toggle functions on/off, between two different settings or between selections.

Enable the **Blue Shift** key and press **S** to add a pause in the User Define or A/N key programming sequence.

With a pause in the user define sequence, the first time the User Define or A/N key is pressed, all steps will be executed up to the pause. The next time the key is pressed, all steps will be executed from the pause to the end (or next pause). This is especially useful in programming "toggle" functions.



### User Define Key Program Example

The following is a brief example of programming a user define sequence showing the use of Ctrl, 1 and Blue Shift, S.

As an example, assume that it is desired to turn the measurement “Average Activity” on and off. This parameter is in the Set Up/Preset Program menu page 3. Perform the following steps to program a User Define function to toggle Average Activity:

- Press the **Set Up** top menu key.
- Select and highlight **User Define** in the Set Up Sub-menu.
- Select and highlight **Presets Program** in the Set Up Sub-menu.
- Press **Ctrl, 1** simultaneously (this ensures that the menu page displayed is number one).
- Press **Set** twice.
- **Trackball** to “Average Activity OFF” and press **Set**.
- **Trackball** to “Exit” and press **Set**.
- Enable the **Blue Shift** key (key is lighted) and press **S**. This adds the pause to the sequence.
- Disable the **Blue Shift** key (key is not lighted).
- Select and highlight **Presets Program** in the Set Up Sub-menu.
- Press **Ctrl, 1** simultaneously (this ensures that the menu page displayed is number one).
- Press **Set** twice.
- **Trackball** to “Average Activity ON” and press **Set**.
- **Trackball** to “Exit” and press **Set**.
- Press the desired **User Define** or **A/N** key. This ends the programming operation.



User Define Key Program Example (cont'd)

The first time the User Define or A/N key is pressed, Average Activity will be set to OFF. The second time the User Define or A/N key is pressed, Average Activity will be set to ON. The third time the User Define or A/N key is pressed, Average Activity will be set to OFF. Each press of the **User Define** or **A/N** key toggles the Average Activity setting.

User Define Names & Lock/Unlock

The User Define Key Names & Locks Sub-Menu in the Set Up/System Parameters menu page 3 will allow the user to lock a key, keeping it from being deleted or reprogrammed. This Sub-Menu also allows the user to name or change the name of any user defined key.




Figure 14–28. User Define Key Names & Locks Sub-Menu



## User Define Function

---

### Names

To enter or edit the name of a User Define key:

- Place the arrow cursor to the right of the desired key in the sub-menu.
- Press **Set**.
- In the highlighted area that appears, type in the desired name. Eleven (11) characters maximum.

### Lock/Unlock

An asterisk next to the User Define key indicates the key is locked. The key cannot be programmed or the function defined for that key cannot be edited.

To lock or unlock a key:

- Place the arrow cursor on the desired key in the sub-menu.
- Press **Set**.
- An asterisk appears to indicate it is locked. No asterisk indicates it is unlocked.



### Hints

Press Shift + User Define key to display the key name on the monitor.



## **Deleting User Define Functions**

To delete a User Define or A/N key sequence:

- Ensure that the User Define or A/N key to be deleted is unlocked.
- Select **Set Up** Top Menu.
- Select User Define function.
- Press the User Define or A/N key number to be deleted.
- Press 'Y' (for Yes) to overwrite current data.

A User Define or A/N key can also be re-programmed by writing over (replacing) the existing sequence.

## **Saving User Define Functions**

All user define functions can be saved to MOD by using Setup/Utilities/User Data Back Up/Save. The data can be loaded back on to the system in the event of a hard drive failure.



# User Data Back-up

## Overview

The Set Up/Utility Menu, option #04 allows the user to back-up (save) all user preset data to MOD. In the event of a system upgrade or catastrophic hard drive failure, the user application presets can be easily loaded back to the system if they have been previously saved to MOD.

Once presets have been backed up, the user data back-up should be performed periodically to update any changes that have been made to the user preset parameters.

## Saving Presets

Before saving presets, the MOD must first be initialized (formatted) prior to saving user application preset data.

- Select Set Up/Utility Menu.
- Select #09 Media Initialize and press **Return**.
- Insert a MOD into the drive and press **Return**.

The system then takes the time to format the MOD (approx. 12 minutes).



*NOTE: If presets have already been saved on the MOD, the MOD should be initialized again to ensure that old preset data has been deleted from the MOD.*

To save presets:

- Enter #04 User Data Backup from the Utility Menu and press **Return**.
- Enter # 1 (Save) and press **Return**.





## **Saving Presets (cont'd)**

- Set media to drive (ensure the MOD is in the drive) and press **Return**.
- Select the exam category name to be saved (01–08) or “All System & Application Data” (09), enter its number and press **Return**.

The system takes a few minutes to save the selected item to the MOD. A succeed or fail message will be displayed.

Eject the MOD, label it and store it in a secure place for future reference.

## **Loading Presets**

To load presets back onto the system:

- Type #04 User Data Backup from the Utility Menu and press **Return**.
- Type #02 Load and press **Return**.
- Set the media to the drive (ensure the MOD is in the drive) and press **Return**.
- Select the application destination and press **Return**.
- Eject the MOD and store it in a secure place.

The system will take a few minutes to load the presets from the MOD to the system hard drive.



## *User Data Backup*

---

This page left blank intentionally.



# Probe Overview

## Ergonomics

Probes have been ergonomically designed to:

- Handle and manipulate with ease
- Connect to the system with one hand
- Be lightweight and balanced
- Have rounded edges and smooth surfaces.

Cables have been designed to:

- Connect to system with appropriate cable length
- Stand up to typical wear by cleaning and disinfectant agents, contact with approved gel, etc.


## Cable handling

Take the following precautions with probe cables:

- Keep free from wheels
- Do not bend the cable acutely
- Avoid crossing cables between probes.

## Probe orientation

Each probe is provided with an orientation marking (refer to Figure 15–1). This mark is used to identify the end of the probe corresponding to the side of the image having the orientation mark on the display.



1 Orientation Marking

Figure 15–1. Orientation Marking on Probe

## Labeling

Each probe is labeled with the following information:

- Seller's name and manufacturer
- Operating frequency
- GE part number
- Probe serial number
- Month and year of manufacture
- Probe designation—provided on the probe grip and the top of the connector housing, so it is easily read when mounted on the system and is also automatically displayed on the screen when the probe is selected.



Labeling (cont'd)




Figure 15-2. Probe Adapter Label




Figure 15-3. Probe Handle Labels

## Labeling (cont'd)




Figure 15-4. Probe Connector Labels




Figure 15-5. Displayed Probe Information



## Probe Overview

### Applications

Below is a list of probes and their intended applications.

| Probe Application                                                                   | B510 | C358 | C364 | C386 | C551 | C721 | E721 | ERB7 | I739 | LA39 | 546L | 739L |
|-------------------------------------------------------------------------------------|------|------|------|------|------|------|------|------|------|------|------|------|
| Abdomen                                                                             |      | ●    | ●    | ●    | ●    |      |      |      |      |      | ●    |      |
| Small Parts                                                                         |      |      |      |      |      |      |      |      | ○    | ●    | ○    | ●    |
| Periph. Vasc.                                                                       |      |      |      |      |      |      |      |      | ○    | ●    | ●    | ●    |
| Obstetrics                                                                          |      | ●    | ●    | ●    | ●    |      | ●    |      |      |      | ○    |      |
| Gynecology                                                                          |      | ●    | ○    | ○    | ●    |      | ●    |      |      |      | ○    |      |
| Pediatrics                                                                          |      | ○    |      |      | ○    | ○    |      |      |      | ○    |      | ○    |
| Neonatal                                                                            |      |      |      |      |      | ●    |      |      |      | ○    |      |      |
| Urology                                                                             |      | ○    | ○    |      | ○    |      | ●    | ●    |      |      |      | ○    |
| Surgery                                                                             |      |      |      |      |      |      |      |      | ●    | ○    |      |      |
| Cardiac                                                                             | ●    |      |      |      |      |      |      |      |      |      |      |      |
| Endocavity                                                                          |      |      |      |      |      |      | ●    | ●    |      |      |      |      |
| Biopsy                                                                              |      | ✓    | ✓    | ✓    | ✓    | ✓    | ✓    | ✓    |      | ✓    | ✓    | ✓    |
| ● = Main Application      ○ = Alternative Application      ✓ = Option Kit available |      |      |      |      |      |      |      |      |      |      |      |      |

| Probe Application                                                                   | L764 | LD | P509 | S220 | S222 | S316 | S317 | S611 | T739 | CWD2 | CWD5 |
|-------------------------------------------------------------------------------------|------|----|------|------|------|------|------|------|------|------|------|
| Abdomen                                                                             |      |    |      | ○    | ○    | ●    | ●    | ●    |      |      |      |
| Small Parts                                                                         | ●    |    |      |      |      |      |      |      | ○    |      |      |
| Periph. Vasc.                                                                       | ○    |    |      |      |      |      |      |      | ○    |      | ●    |
| Obstetrics                                                                          |      |    |      |      |      |      |      |      |      |      |      |
| Gynecology                                                                          |      |    |      |      |      |      |      |      |      |      |      |
| Pediatrics                                                                          |      |    |      |      |      | ○    | ○    | ●    |      |      |      |
| Neonatal                                                                            |      |    |      |      |      |      |      |      |      |      |      |
| Urology                                                                             |      |    |      | ○    | ○    | ○    | ○    |      |      |      |      |
| Surgery                                                                             |      |    |      |      |      |      |      |      | ●    |      |      |
| Cardiac                                                                             |      |    | ●    | ●    | ○    | ●    | ●    | ●    |      | ●    |      |
| Transcranial                                                                        |      |    |      |      | ●    |      |      |      |      |      |      |
| Biopsy                                                                              | ✓    | ●  |      |      |      | ✓    | ✓    | ✓    | ✓    |      |      |
| ● = Main Application      ○ = Alternative Application      ✓ = Option Kit available |      |    |      |      |      |      |      |      |      |      |      |

Table 15–1. Probe Indications for Use

## Specifications

| Probe Designation                                                                                                                                                                                                                                                                                            | Principal Clinical Uses             | Image Frequency (MHz) | Doppler Frequency (MHz) |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------|-----------------------|-------------------------|
| B510                                                                                                                                                                                                                                                                                                         | Cardiac                             | 5.0                   | 4.0                     |
| C358                                                                                                                                                                                                                                                                                                         | Abdom./OB/GYN                       | 3.8                   | 2.5                     |
| C364                                                                                                                                                                                                                                                                                                         | Radiology/Obstetrics                | 3.3                   | 2.5                     |
| C386                                                                                                                                                                                                                                                                                                         | OB/GYN/Abdom.                       | 3.5                   | 2.5                     |
| C551                                                                                                                                                                                                                                                                                                         | Radiology/Obstetrics                | 5.0                   | 4.0                     |
| C721                                                                                                                                                                                                                                                                                                         | Neonatal                            | 6.6                   | 5.0                     |
| E721                                                                                                                                                                                                                                                                                                         | Transvaginal<br>Transrectal         | 6.6                   | 5.0                     |
| ERB7                                                                                                                                                                                                                                                                                                         | Urology                             | 7.0                   | 5.0                     |
| I739                                                                                                                                                                                                                                                                                                         | Intra-operative                     | 6.7                   | 5.0                     |
| 546L                                                                                                                                                                                                                                                                                                         | Deep Vascular/<br>Small Parts       | 5.2                   | 4.0                     |
| 739L                                                                                                                                                                                                                                                                                                         | Peripheral Vascular/<br>Small Parts | 6.7                   | 5.0                     |
| L764                                                                                                                                                                                                                                                                                                         | Mammo, PV & Small<br>Parts          | 6.6                   | 5.0                     |
| LA39                                                                                                                                                                                                                                                                                                         | Small Parts                         | 8.7                   | 5.0                     |
| LD                                                                                                                                                                                                                                                                                                           | Biopsy                              | 3.5                   | n/a                     |
| P509                                                                                                                                                                                                                                                                                                         | Cardiac                             | 5.0                   | 3.5                     |
| S220                                                                                                                                                                                                                                                                                                         | Cardiac                             | 2.5                   | 2.2                     |
| S222                                                                                                                                                                                                                                                                                                         | Transcranial,<br>Cardiac            | 2.9                   | 2.0                     |
| S316                                                                                                                                                                                                                                                                                                         | Cardiac                             | 3.3                   | 2.5                     |
| S317                                                                                                                                                                                                                                                                                                         | Cardiac/Abdom.                      | 3.3                   | 2.5                     |
| S611                                                                                                                                                                                                                                                                                                         | Cardiac                             | 5.7                   | 4.0                     |
| T739                                                                                                                                                                                                                                                                                                         | Intra-operative                     | 6.7                   | 5.0                     |
| <div> <div> <b>B</b> = Bipolar Sector<br/> <b>E</b> = Endocavitary<br/> <b>L</b> = Linear Array<br/> <b>S</b> = Phased Array Sector </div> <div> <b>C</b> = Curved Array (Convex)<br/> <b>I</b> = Intra-operative Linear<br/> <b>P</b> = Multiplane TEE<br/> <b>T</b> = Intra-operative Linear </div> </div> |                                     |                       |                         |

Table 15–2. System Probe Definitions





## Probe Overview

### Specifications (cont'd)

| Adapter Designation | Adaptable Probe | System Operates Probe As: |
|---------------------|-----------------|---------------------------|
| PA 51               | CBF             | C364                      |
|                     | CAE             | C551                      |
|                     | LH              | L764                      |
|                     | MTZ             | E721                      |
| 5S                  | UC              | S316                      |
|                     | W               | S220                      |

Table 15–3. Probe Adapter Usage

| Probe Name | Material of Headshell | Use                | Type               | Catalog Number | Mfg. by | Probe Family (Headshell) | Part No.  |
|------------|-----------------------|--------------------|--------------------|----------------|---------|--------------------------|-----------|
| B510       | PU                    | Intercavity        | Biplane Sector     | H45202BT       | GEYMS   |                          | 2123593   |
| C358       | NORYL                 | Abdom./ OB/GYN     | Convex             | H45202CD       | GEYMS   |                          | 2193617   |
| C364       | PES                   | Abdom.             | Convex             | H45202CF       | GEYMS   | CBF                      | P9607AB   |
| C386       | NORYL                 | Abdom./ OB/GYN     | Convex             | H45202CC       | GEYMS   |                          | 2147187–2 |
| C551       | PES                   | Abdom.             | Convex             | H45202CE       | GEYMS   | CAE                      | P9607AD   |
| C721       | NORYL                 | Neonatal           | Convex             | H45202MN       | GEYMS   | MTZ                      | 2121267-2 |
| E721       | PES                   | Intercavity        | Convex             | H45202MT       | GEYMS   | MTZ                      | P9607AF   |
| ERB7       | NORYL                 | Urology            | Convex/ Linear     | H40392L        | GEYMS   |                          | 2239590   |
| I739       | NORYL                 | Intra–op.          | Linear             | H45202JG       | GEYMS   |                          | 2147189–2 |
| 546L       | NORYL                 | Abdom. Superficial | Linear             | H45202LE       | GEYMS   |                          | 2144266–2 |
| 739L       | NORYL                 | Superficial        | Linear             | H45202AG       | GEYMS   |                          | 2107460–2 |
| L764       | PES                   | Sm. Parts          | Linear             | H45202HP       | GEYMS   | LH                       | 2121377-2 |
| LA39       | NORYL                 | Superficial        | Linear             | H45202LA       | GEYMS   |                          | 2155078–2 |
| LD         | PES                   | Biopsy             | Linear             | H45202LD       | GEYMS   |                          | 2124318   |
| P509       | PU                    | Intercavity        | Multi plane Sector | H45202RT       | GEYMS   |                          | 2169773   |
| S220       | PES                   | Cardiac            | Sector             | H45202WG       | GEYMS   | W                        | 2121793–2 |
| S222       | NORYL                 | Trans-cranial      | Sector             | H45202TC       | GEYMS   |                          | 2159263   |

Table 15–4. LOGIQ™ 500 Probe List

Specifications (cont'd)

| Probe Name | Material of Headshell | Use                                | Type          | Catalog Number | Mfg. by | Probe Family (Headshell) | Part No.  |
|------------|-----------------------|------------------------------------|---------------|----------------|---------|--------------------------|-----------|
| S316       | PES                   | Cardiac/<br>Abdom.                 | Sector        | H45202SC       | GEYMS   | UC                       | P9606AB   |
| S317       | NORYL                 | Cardiac/<br>Abdom.                 | Sector        | H45202SD       | GEYMS   |                          | 2144268-2 |
| S611       | NORYL                 | Cardiac/<br>Neonatal/<br>Pediatric | Sector        | H45202SF       | GEYMS   |                          | 2144267-2 |
| T739       | NORYL                 | Intra-op.                          | Linear        | H45202TG       | GEYMS   |                          | 2147188-2 |
| CWD2       | NORYL                 | Cardiac                            | Single<br>CWD | H45202DB       | GEYMS   |                          | 2123594   |
| CWD5       | NORYL                 | Peripheral<br>Vascular             | Single<br>CWD | H45202DE       | GEYMS   |                          | 2123595   |
| CBF        | PES                   | Abdom.                             | Convex        | H46022CB       | GEYMS   |                          | P9603AD   |
| CAE        | PES                   | Abdom.                             | Convex        | H46022CA       | GEYMS   |                          | P9603AE   |
| MTZ        | PES                   | Intercavity                        | Convex        | H46022MT       | GEYMS   |                          | P9603AU   |
| UC         | PES                   | Cardiac/<br>Abdom.                 | Sector        | H4163B         | GEYMS   |                          | P9600BE   |
| LH         | PES                   | Sm. Parts                          | Linear        | H46022LH       | GEYMS   |                          | P9601AS   |
| W          | PES                   | Cardiac                            | Sector        | H4162C         | GEYMS   |                          | P9600BH   |

Table 15-4. LOGIQ™ 500 Probe List (cont'd)



## Probe Overview

---

### Probe Usage

For details on connecting, activating, deactivating, disconnecting, transporting and storing the probes, refer to 3–27.

### Care and Maintenance

#### Inspecting probes



Inspect the probe's lens, cable, and casing after each use. Look for any damage that would allow liquid to enter the probe. If any damage is found, do not use the probe until it has been inspected and repaired/replaced by a GE Service Representative.



Keep a log of all probe maintenance, along with a picture of any probe malfunction.

#### Environmental Requirements

Probes should be operated, stored, or transported within the parameters outlined below.

|             | Operational                 | Storage                      | Transport                     |
|-------------|-----------------------------|------------------------------|-------------------------------|
| Temperature | 10° - 40° C<br>50° - 104° F | -10° - 60° C<br>14° - 140° F | -40° - 60° C<br>-40° - 140° F |
| Humidity    | 30-85%<br>non-condensing    | 30-90%<br>non-condensing     | 30-90%<br>non-condensing      |
| Pressure    | 700-1060hPa                 | 700-1060hPa                  | 700-1060hPa                   |

Table 15–5. Probe Environmental Requirements

## Probe Safety

### Handling precautions

#### WARNING



Ultrasound probes are highly sensitive medical instruments that can easily be damaged by improper handling. Use care when handling and protect from damage when not in use. **DO NOT** use a damaged or defective probe. Failure to follow these precautions can result in serious injury and equipment damage.

### Electrical shock hazard



#### Electrical Hazard

The probe is driven with electrical energy that can injure the patient or user if live internal parts are contacted by conductive solution:

- **DO NOT** immerse the probe into any liquid beyond the level indicated by the immersion level diagram. Refer to Figure 15–6 on 15–15. Never immerse the probe connector or probe adaptors into any liquid.
- **DO NOT** drop the probes or subject them to other types of mechanical shock or impact. Degraded performance or damage such as cracks or chips in the housing may result.
- Inspect the probe before and after each use for damage or degradation to the housing, strain relief, lens, and seal. A thorough inspection should be conducted during the cleaning process.
- **DO NOT** kink, tightly coil, or apply excessive force on the probe cable. Insulation failure may result.
- Electrical leakage checks should be performed on a routine basis by GE Service or qualified hospital personnel. Refer to the service manual for leakage check procedures.

### Mechanical hazards

#### CAUTION



A defective probe or excessive force can cause patient injury or probe damage:

- Observe depth markings and do not apply excessive force when inserting or manipulating intercavitary probes.
- Inspect probes for sharp edges or rough surfaces that could injure sensitive tissue.

**Special handling instructions****Using protective sheaths****CAUTION****Biological  
Hazard**

Protective barriers may be required to minimize disease transmission. Probe sheaths are available for use with all clinical situations where infection is a concern. Use of legally marketed, sterile probe sheaths is strongly recommended for intra-cavitary and intra-operative procedures. Use of legally marketed, sterile, pyrogen free probe sheaths is **REQUIRED** for neurological intra-operative procedures.

**Instructions.** Custom made sheaths are available for each probe. Each probe sheath kit consists of a flexible sheath used to cover the probe and cable and elastic bands used to secure the sheath.

Sterile probe sheaths are supplied as part of biopsy kits for those probes intended for use in biopsy procedures. In addition to the sheath and elastic bands, there are associated accessories for performing a biopsy procedure which are included in the kit. Refer to the biopsy instructions for the specific probes in the *Discussion* section of this chapter for further information.

**Reordering.** To reorder sheaths, refer to 16–42.

**CAUTION**

Devices containing latex may cause severe allergic reaction in latex sensitive individuals. Refer to FDA's March 29, 1991 Medical Alert on latex products.

**CAUTION**

Do not use pre-lubricated condoms as a sheath. In some cases, they may damage the probe. Lubricants in these condoms may not be compatible with probe construction.



## Probe handling and infection control

This information is intended to increase user awareness of the risks of disease transmission associated with using this equipment and provide guidance in making decisions directly affecting the safety of the patient as well as the equipment user.

Diagnostic ultrasound systems utilize ultrasound energy that must be coupled to the patient by direct physical contact. Depending on the type of examination, this contact occurs with a variety of tissues ranging from intact skin in a routine exam to recirculating blood in a surgical procedure. The level of risk of infection varies greatly with the type of contact.

One of the most effective ways to prevent transmission between patients is with single use or disposable devices. However, ultrasound transducers are complex and expensive devices that must be reused between patients. It is very important, therefore, to minimize the risk of disease transmission by using barriers and through proper processing between patients.

### CAUTION



Adequate cleaning and disinfection are necessary to prevent disease transmission. It is the responsibility of the equipment user to verify and maintain the effectiveness of the infection control procedures in use. Always use sterile, legally marketed probe sheaths for intra-cavitary and intra-operative procedures.

For neurological intra-operative procedures, use of a legally marketed, sterile, pyrogen free probe sheath is **REQUIRED**. Probes for neuro surgical use must not be sterilized with liquid chemical sterilants because of the possibility of neuro toxic residues remaining on the probe.



## Probe Overview

---

### Probe Cleaning Process



To clean the probe:

1. After each use, disconnect the probe from the ultrasound console and remove all coupling gel from the probe by wiping with a soft cloth and rinsing with flowing water.
2. Wash the probe with mild soap in lukewarm water. Scrub the probe as needed using a soft sponge, gauze, or cloth to remove all visible residue from the probe surface. Prolonged soaking or scrubbing with a soft bristle brush (such as a toothbrush) may be necessary if material has dried onto the probe surface.
3. Rinse the probe with enough clean potable water to remove all visible soap residue.
4. Air dry or dry with a soft cloth.



Special Cleaning Instructions for the E721(MTZ): When cleaning the E721(MTZ) probe, it is important to be sure that all surfaces are thoroughly cleaned. This probe has an adjustable two-part handle that must be disassembled to gain access to all surfaces. To disassemble the handle, completely remove the handle adjustment screw located mid-way between the cable entry and probe tip. The two handle halves and adjustment screw must be thoroughly cleaned along with the main probe shaft as described earlier in step 2. After rinsing and drying is completed, the probe handle can be loosely reassembled for the disinfection process.

Probe Cleaning Process (cont'd)




Figure 15–6. Probe Immersion Levels





## Probe Overview

---

### Disinfecting probes



After Each Use

Ultrasound probes can be disinfected using liquid chemical germicides. The level of disinfection is directly related to the duration of contact with the germicide. Increased contact time produces a higher level of disinfection.

2% Glutaraldehyde-based solutions have been shown to be very effective for this purpose. Cidex is the only germicide that has been evaluated for compatibility with the material used to construct the probes.

### CAUTION



In order for liquid chemical germicides to be effective, all visible residue must be removed during the cleaning process. Thoroughly clean the probe, as described earlier before attempting disinfection.

1. Prepare the germicide solution according to the manufacturer's instructions. Be sure to follow all precautions for storage, use and disposal.
2. Place the cleaned and dried probe in contact with the germicide for the time specified by the germicide manufacturer. High-level disinfection is recommended for surface probes and is required for endocavitary and intraoperative probes (follow the germicide manufacturer's recommended time).

Probes for neuro surgical intra-operative use must NOT be sterilized with liquid chemical sterilants because of the possibility of neuro toxic residues remaining on the probe. Neurological procedures must be done with the use of legally marketed, sterile, pyrogen free probe sheaths.

3. After removing from the germicide, rinse the probe following the germicide manufacturer's rinsing instructions. Flush all visible germicide residue from the probe and allow to air dry.

**Disinfecting probes (cont'd)**

Special Disinfecting Instructions for the E721(MTZ): To properly disinfect the E721(MTZ) probe, the probe handle can be reassembled loosely so that the entire probe with handle can be immersed in the germicide solution. The adjustment screw must be kept loose so that germicide can penetrate to all surfaces. After immersing, rotate and shake the probe while it is below the surface of the germicide to eliminate air pockets. Allow the germicide to remain in contact with the fully immersed probe, for high level disinfection, according to the germicide manufacturer's recommended time. To remove all germicide residue, final rinsing should be done following the germicide manufacturer's instructions. Remove excess water by shaking and allow to air dry.

**CAUTION****CREUTZFELD–JACOB DISEASE**

Neurological use on patients with this disease must be avoided. If a probe becomes contaminated, there is no adequate disinfecting means.

**Biological  
Hazard**

**Disinfecting probes (cont'd)****WARNING**

Ultrasound transducers can easily be damaged by improper handling and by contact with certain chemicals. Failure to follow these precautions can result in serious injury and equipment damage.

- Do not immerse the probe into any liquid beyond the level specified for that probe. Never immerse the transducer connector or probe adapters into any liquid.
- Avoid mechanical shock or impact to the transducer and do not apply excessive bending or pulling force to the cable.
- Transducer damage can result from contact with inappropriate coupling or cleaning agents:
  - Do not soak or saturate transducers with solutions containing alcohol, bleach, ammonium chloride compounds or hydrogen peroxide
  - Avoid contact with solutions or coupling gels containing mineral oil or lanolin
  - Avoid temperatures above 60° C.
- Inspect the probe prior to use for damage or degeneration to the housing, strain relief, lens and seal. Do not use a damaged or defective probe.

## Coupling gels

### CAUTION



Do not use unrecommended gels (lubricants). They may damage the probe and void the warranty.

### Applying

In order to assure optimal transmission of energy between the patient and probe, a conductive gel or couplant must be applied liberally to the patient where scanning will be performed.

### Precautions

Coupling gels should not contain the following ingredients as they are known to cause probe damage:

- Methanol, ethanol, isopropanol, or any other alcohol-based product
- Mineral oil
- Iodine
- Lotions
- Lanolin
- Aloe Vera
- Olive Oil
- Methyl or Ethyl Parabens (para hydroxybenzoic acid)
- Dimethylsilicone

## Planned Maintenance

### CAUTION



The following maintenance schedule is suggested for the system and probes to ensure optimum operation and safety.

| Do the Following   | Daily | After Each Use | As Necessary |
|--------------------|-------|----------------|--------------|
| Inspect the Probes | X     |                |              |
| Clean the Probes   |       | X              |              |
| Disinfect Probes   |       | X              |              |

Table 15–6. Planned Maintenance Program



# Probe Discussion

## Introduction

The LOGIQ™ 500 supports four types of probes:

- **Curved Array (Convex).** Curved Array (Convex) probes, including 'micro' convex, are usually designated by the prefix/suffix "C"; the endocavity probe is designated by the prefix/suffix "E".
- **Linear Array.** Linear Array probes are designated by the prefix/suffix "L"; the linear intra-operative probes are designated by the prefix/suffix "I" or "T".
- **Phased Array Sector.** Phased Array Sector probes are designated by the prefix/suffix "S"; the biplane TEE probe is designated by the prefix/suffix "B"; the multi-plane TEE probe is designated by the prefix/suffix "P".
- **Continuous Wave Doppler.** Continuous Wave Doppler probes are designated by the prefix/suffix "CWD".

## Probe naming conventions

|                                                                                                                                                                                                                                                             | TYPE                                                                                                                                                                                                                                                                                                                                          |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <div><div><div>C364L</div><div>— *Type</div><div>— Array Aperture (mm)</div><div>— Center Frequency</div><div>— Type</div></div><div><div>CWD2</div><div>— Center Frequency</div><div>— Type</div></div><div>* Type can be a prefix or a suffix</div></div> | <div><div>B</div><div>Biplane TEE</div><div>C</div><div>Convex</div><div>CWD</div><div>Pencil CWD</div><div>E</div><div>Endo (TR/TV)</div><div>I</div><div>I-shaped intra-operative</div><div>L</div><div>Linear</div><div>P</div><div>Multiplane TEE</div><div>S</div><div>Sector</div><div>T</div><div>T-shaped intra-operative</div></div> |

Table 15–7. Probe Naming Conventions

## Curved Array (Convex) Probes

The convex probes for the system are the C358, C364 (CBF), C386, C551 (CAE), C721, E721 (MTZ) and ERB7 probes. The ERB7 probe is a bi-plane (convex and linear) probe.

**ERB7 probe not currently available in the USA.**

## Biopsy Guidezone

The type of guide displayed for the convex probes is illustrated below.





Figure 15–7. Typical Convex Probe Biopsy Guidezone

| Probe             | Description                                                                                   | Intended Uses                                                                                        | Capabilities and Features                                                                                                                                                                                                                                            | Illustration |
|-------------------|-----------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|
| <b>C358</b>       | The C358 probe is a general purpose probe for obtaining B-Mode, Doppler, and Color Flow data. | <ul style="list-style-type: none"> <li>General Purpose</li> </ul>                                    | <ul style="list-style-type: none"> <li>Wide field of view</li> <li>Penetration</li> <li>Good image uniformity</li> <li>Good B-Mode Resolution</li> <li>CFM/Doppler detectability</li> <li>Ergonomics for scanning and cleaning</li> <li>Biopsy capability</li> </ul> |              |
| <b>C364 (CBF)</b> | The C364 probe is a general purpose probe for obtaining B-Mode, Doppler, and Color Flow data. | <ul style="list-style-type: none"> <li>General Purpose</li> <li>Radiology</li> <li>OB/GYN</li> </ul> | <ul style="list-style-type: none"> <li>Wide field of view</li> <li>Penetration</li> <li>Good B-Mode Resolution</li> <li>CFM/Doppler detectability</li> <li>Ergonomics for scanning and cleaning</li> <li>Biopsy capability</li> </ul>                                |              |



## Probe Discussion

### Curved Array (Convex) Probes (cont'd)

| Probe             | Description                                                                                         | Intended Uses                                                                                                                 | Capabilities and Features                                                                                                                                                                                                                                                                             | Illustration                                                                          |
|-------------------|-----------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>C386</b>       | The C386 probe is a general purpose convex probe for obtaining B-Mode, Doppler and Color Flow data. | <ul style="list-style-type: none"> <li>• OB/GYN</li> <li>• General Abdominal</li> </ul>                                       | <ul style="list-style-type: none"> <li>• Wideband for B-Mode resolution &amp; homogeneity</li> <li>• CFM/Doppler sensitivity</li> <li>• Wide field of view</li> <li>• Large radius for better surface contact</li> <li>• Ergonomics for scanning and cleaning</li> <li>• Biopsy capability</li> </ul> |    |
| <b>C551 (CAE)</b> | The C551 probe is a general purpose probe for obtaining B-Mode, Doppler, and Color Flow data.       | <ul style="list-style-type: none"> <li>• General Purpose</li> <li>• Radiology</li> <li>• Obstetrics and Gynecology</li> </ul> | <ul style="list-style-type: none"> <li>• Wide field of view</li> <li>• Good B-Mode Resolution</li> <li>• CFM/Doppler detectability</li> <li>• Ergonomics for scanning and cleaning</li> <li>• Biopsy capability</li> </ul>                                                                            |   |
| <b>C721</b>       | The C721 probe is a general purpose probe for obtaining B-Mode, Doppler, and Color Flow data.       | <ul style="list-style-type: none"> <li>• Neonatal</li> <li>• Pediatrics</li> </ul>                                            | <ul style="list-style-type: none"> <li>• Small Footprint</li> <li>• Wide field of view</li> <li>• Good B-Mode Resolution</li> <li>• CFM/Doppler detectability</li> <li>• Ergonomics for scanning and cleaning</li> <li>• Biopsy capability</li> </ul>                                                 |  |



## Probe Discussion

### Curved Array (Convex) Probes (cont'd)

| Probe             | Description                                                                                                                                                                                                                                                   | Intended Uses                                                                           | Capabilities and Features                                                                                                                                                                                                                                                                              | Illustration                                                                         |
|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| <b>E721 (MTZ)</b> | The E721 probe is an endocavitary probe for obtaining B-Mode, Doppler, and Color Flow data.                                                                                                                                                                   | <ul style="list-style-type: none"> <li>• Transvaginal</li> <li>• Transrectal</li> </ul> | <ul style="list-style-type: none"> <li>• Wide field of view</li> <li>• Small headshell and probe shaft</li> <li>• Adjustable handle angle</li> <li>• Good B-Mode Resolution</li> <li>• CFM/Doppler sensitivity</li> <li>• Ergonomics for scanning and cleaning</li> <li>• Biopsy capability</li> </ul> |   |
| <b>ERB7</b>       | The ERB7 probe is a bi-plane (convex and linear) probe for endocavity scanning and stepper volume measurements. The ERB7 can also be utilized for obtaining B-Mode, Doppler, and Color Flow images.<br><br><b>The ERB7 probe is not available in the USA.</b> | <ul style="list-style-type: none"> <li>• Transrectal</li> </ul>                         | <ul style="list-style-type: none"> <li>• Wide field of view</li> <li>• Small headshell and probe shaft</li> <li>• Adjustable handle angle</li> <li>• Good B-Mode Resolution</li> <li>• CFM/Doppler sensitivity</li> <li>• Ergonomics for scanning and cleaning</li> <li>• Biopsy capability</li> </ul> |  |





## Probe Discussion

### Linear Array Probes

The linear probes for the system are the I739, 546L, 739L, L764 (LH), T739, LA39, LD and ERB7 probes. The ERB7 probe is a bi-plane (convex and linear) probe, see 15–23 for details.

**The LD probe is available only in Europe.**

### Biopsy Guidezone

The type of guide displayed for the linear probe is illustrated below.







Figure 15–8. Typical Linear Probe Biopsy Guidezone

| Probe       | Description                                                                                         | Intended Uses                                                                                                                                | Capabilities and Features                                                                                                                                                                                                                                                                                                             | Illustration |
|-------------|-----------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|
| <b>I739</b> | The I739 probe is a linear probe intended for intra-operative applications.                         | <ul style="list-style-type: none"> <li>• Intra-operative imaging</li> <li>• Superficial imaging for organs with space constraints</li> </ul> | <ul style="list-style-type: none"> <li>• Low height “I” shaped micro-case</li> <li>• Long and light probe cable</li> <li>• Wideband for B-Mode resolution &amp; homogeneity</li> <li>• Slant Scan</li> <li>• CFM Doppler sensitivity</li> <li>• Ergonomics for scanning and cleaning</li> </ul>                                       |              |
| <b>546L</b> | The 546L probe is a general purpose linear probe for obtaining B-Mode, Doppler and Color flow data. | <ul style="list-style-type: none"> <li>• Small Parts</li> <li>• Deep vessel vascular</li> <li>• General Abdominal and OB/GYN</li> </ul>      | <ul style="list-style-type: none"> <li>• Wide field of view</li> <li>• Slant Scan</li> <li>• Good penetration and resolution for deep vascular</li> <li>• Wideband for B-Mode resolution &amp; homogeneity</li> <li>• CFM Doppler sensitivity</li> <li>• Ergonomics for scanning and cleaning</li> <li>• Biopsy capability</li> </ul> |              |

## Linear Array Probes (cont'd)

| Probe            | Description                                                                                         | Intended Uses                                                                                                                                | Capabilities and Features                                                                                                                                                                                                                                                                                                    | Illustration                                                                          |
|------------------|-----------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>739L</b>      | The 739L probe is a general purpose probe for obtaining B-Mode, Doppler, and Color Flow data.       | <ul style="list-style-type: none"> <li>• Small Parts</li> <li>• Peripheral Vascular</li> </ul>                                               | <ul style="list-style-type: none"> <li>• Wide field of view</li> <li>• Slant scan</li> <li>• Wideband for B-Mode resolution &amp; homogeneity</li> <li>• CFM Doppler sensitivity</li> <li>• Ergonomics for scanning and cleaning</li> <li>• Biopsy capability</li> </ul>                                                     |    |
| <b>L764 (LH)</b> | The L764 probe is a general purpose probe for obtaining B-Mode, Doppler, and Color Flow data.       | <ul style="list-style-type: none"> <li>• Mammography</li> <li>• Small Parts</li> <li>• Peripheral Vascular</li> </ul>                        | <ul style="list-style-type: none"> <li>• Wide field of view</li> <li>• Good B-Mode Resolution</li> <li>• CFM Doppler detectability</li> <li>• Ergonomics for scanning and cleaning</li> <li>• Biopsy capability</li> </ul>                                                                                                   |   |
| <b>T739</b>      | The T739 probe is a linear probe intended for intra-operative applications.                         | <ul style="list-style-type: none"> <li>• Intra-operative imaging</li> <li>• Superficial imaging for organs with space constraints</li> </ul> | <ul style="list-style-type: none"> <li>• Low height "T" shaped micro-case</li> <li>• Long and light probe cable</li> <li>• Slant scan</li> <li>• Wideband for B-Mode resolution &amp; homogeneity</li> <li>• CFM Doppler sensitivity</li> <li>• Ergonomics for scanning and cleaning</li> <li>• Biopsy capability</li> </ul> |  |
| <b>LA39</b>      | The LA39 probe is a general purpose linear probe for obtaining B-Mode, Doppler and Color Flow data. | <ul style="list-style-type: none"> <li>• Intra-operative imaging</li> <li>• Superficial imaging for organs with space constraints</li> </ul> | <ul style="list-style-type: none"> <li>• Slant scan</li> <li>• Wideband for B-Mode resolution &amp; homogeneity</li> <li>• CFM Doppler sensitivity</li> <li>• Ergonomics for scanning and cleaning</li> <li>• Biopsy capability</li> </ul>                                                                                   |  |
| <b>LD</b>        | The LD probe is a linear probe for biopsies<br><br><b>The LD probe is available only in Europe.</b> | <ul style="list-style-type: none"> <li>• Biopsy</li> </ul>                                                                                   | <ul style="list-style-type: none"> <li>• Ergonomics for scanning and cleaning</li> <li>• High sensitivity of needle echo</li> <li>• Wide field of view</li> </ul>                                                                                                                                                            |   |



## Probe Discussion

### Sector Probes

The sector probes for the system are the S220 (W), S222, S316 (UC), S317, S611, B510 and P509 probes.

**The P509 probe is not available in the USA.**

### Biopsy Guidezone

The type of guide displayed for the sector probes is illustrated below.




Figure 15–9. Typical Sector Probe Biopsy Guidezone

| Probe           | Description                                                                                   | Intended Uses                                                                           | Capabilities and Features                                                                                                                                                                                                                                    | Illustration                                                                          |
|-----------------|-----------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>S220 (W)</b> | The S220 probe is a general purpose probe for obtaining B-Mode, Doppler, and Color Flow data. | <ul style="list-style-type: none"> <li>Cardiology</li> <li>General Abdominal</li> </ul> | <ul style="list-style-type: none"> <li>Small footprint</li> <li>Penetration</li> <li>Steerable Doppler</li> <li>Good B-Mode Resolution</li> <li>CFM Doppler detectability</li> <li>Ergonomics for scanning and cleaning</li> </ul>                           |  |
| <b>S222</b>     | The S222 probe is a general purpose probe for obtaining B-Mode, Doppler, and Color Flow data. | <ul style="list-style-type: none"> <li>Transcranial</li> <li>Cardiology</li> </ul>      | <ul style="list-style-type: none"> <li>Small footprint</li> <li>Dedicated handle design for Transcranial</li> <li>Wideband for B-Mode resolution &amp; homogeneity</li> <li>CFM Doppler sensitivity</li> <li>Ergonomics for scanning and cleaning</li> </ul> |  |

## Sector Probes (cont'd)

| Probe            | Description                                                                                                                                                                                                                          | Intended Uses                                                                                                   | Capabilities and Features                                                                                                                                                                                                                                                                           | Illustration                                                                          |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>S316 (UC)</b> | The S316 probe is a general purpose probe for obtaining B-Mode, Doppler, and Color Flow data.                                                                                                                                        | <ul style="list-style-type: none"> <li>• Cardiology</li> <li>• General Abdominal</li> </ul>                     | <ul style="list-style-type: none"> <li>• Small footprint</li> <li>• Steerable Doppler</li> <li>• Good B-Mode Resolution</li> <li>• CFM Doppler detectability</li> <li>• Ergonomics for scanning and cleaning</li> <li>• Biopsy capability</li> </ul>                                                |    |
| <b>S317</b>      | The S317 probe is a general purpose sector probe for obtaining B-Mode, Doppler and Color Flow data.                                                                                                                                  | <ul style="list-style-type: none"> <li>• Cardiology</li> <li>• General Abdomen</li> </ul>                       | <ul style="list-style-type: none"> <li>• Small footprint</li> <li>• Steerable Doppler</li> <li>• Wideband for B-Mode resolution &amp; homogeneity</li> <li>• CFM Doppler sensitivity</li> <li>• Ergonomically designed micro-case for scanning and cleaning</li> <li>• Biopsy capability</li> </ul> |   |
| <b>S611</b>      | The S611 probe is a general purpose sector probe for obtaining B-Mode, Doppler and Color Flow data.                                                                                                                                  | <ul style="list-style-type: none"> <li>• Cardiology, Pediatric and Neonatal</li> <li>• Neonatal head</li> </ul> | <ul style="list-style-type: none"> <li>• Small footprint</li> <li>• Steerable Doppler</li> <li>• Wideband for B-Mode resolution &amp; homogeneity</li> <li>• CFM Doppler sensitivity</li> <li>• Ergonomically designed micro-case for scanning and cleaning</li> <li>• Biopsy capability</li> </ul> |  |
| <b>B510</b>      | <p>The B510 probe is a Biplane Transesophageal probe for obtaining B-Mode, Doppler, and Color Flow data.</p> <p><i>NOTE: Refer to the User Manual that comes with the probe for details on handling, operation and cleaning.</i></p> | <ul style="list-style-type: none"> <li>• Cardiology from the Esophagus</li> </ul>                               | <ul style="list-style-type: none"> <li>• Small footprint</li> <li>• Biplane imaging</li> <li>• Steerable Doppler</li> <li>• Good B-Mode Resolution</li> <li>• CFM/Doppler sensitivity</li> <li>• Ergonomics for one hand manipulation of the handle</li> </ul>                                      |  |



## Probe Discussion

---

### Sector Probes (cont'd)



| Probe       | Description                                                                                                                                                                                                                                                                                        | Intended Uses                                                                   | Capabilities and Features                                                                                                                                                                                                                                  | Illustration                                                                        |
|-------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <b>P509</b> | <p>The P509 probe is a Multi-plane Transesophageal probe for obtaining B-Mode, Doppler, and Color Flow data.</p> <p><b>The P509 probe is not available in the USA.</b></p> <p><i>NOTE: Refer to the User Manual that comes with the probe for details on handling, operation and cleaning.</i></p> | <ul style="list-style-type: none"><li>• Cardiology from the Esophagus</li></ul> | <ul style="list-style-type: none"><li>• Small footprint</li><li>• Multiplane imaging</li><li>• Steerable Doppler</li><li>• Good B-Mode Resolution</li><li>• CFM/Doppler sensitivity</li><li>• Ergonomics for one hand manipulation of the handle</li></ul> |  |



## Probe Discussion

### CWD Probes

The CWD probes for the system are the CWD2 and CWD5 probes.

| Probe       | Description                                                                      | Intended Uses                                                                           | Capabilities and Features                                                                                                                             | Illustration                                                                         |
|-------------|----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| <b>CWD2</b> | The CWD2 probe is a pencil CWD probe for obtaining Continuous Wave Doppler data. | <ul style="list-style-type: none"><li>• CW Doppler for Cardiology and Abdomen</li></ul> | <ul style="list-style-type: none"><li>• High sensitivity</li><li>• T-bar support for holding</li><li>• Ergonomics for scanning and cleaning</li></ul> |    |
| <b>CWD5</b> | The CWD5 probe is a pencil CWD probe for obtaining Continuous Wave Doppler data. | <ul style="list-style-type: none"><li>• CW Doppler for Peripheral Vascular</li></ul>    | <ul style="list-style-type: none"><li>• High sensitivity</li><li>• T-bar support for holding</li><li>• Ergonomics for scanning and cleaning</li></ul> |  |



## Biopsy Special Concerns

### Precautions Concerning the Use of Biopsy Procedures

#### WARNING



Do not freeze the image during a biopsy procedure. The image must be live to avoid a positioning error.

Biopsy guidezones are intended to assist the user in determining optimal probe placement and approximate the needle path. However, actual needle movement is likely to deviate from the guideline. Always monitor the relative positions of the biopsy needle and the subject mass during the procedure.

#### CAUTION



The use of biopsy devices and accessories that have not been evaluated for use with this equipment may not be compatible and could result in injury.

#### CAUTION



The invasive nature of biopsy procedures requires proper preparation and technique to control infection and disease transmission. Equipment must be cleaned as appropriate for the procedure prior to use.

- Follow the probe cleaning and disinfection procedures and precautions to properly prepare the probe.
- Follow the manufacturer's instructions for the cleaning of biopsy devices and accessories.
- Use protective barriers such as gloves and probe sheaths.
- After use, follow proper procedures for decontamination, cleaning, and waste disposal.

#### CAUTION



Improper cleaning methods and the use of certain cleaning and disinfecting agents can cause damage to the plastic components that will degrade imaging performance or increase the risk of electric shock. Refer to probe safety and handling precautions on 15–11.

# Preparing for a Biopsy

## Displaying the Guidezone

Activate the Biopsy Zone by selecting it from the B-Mode Sub-Menu.

| B           | Preset       | Set Up  | ECG          |
|-------------|--------------|---------|--------------|
| Biopsy Zone | Image Rotatn | Rejectn | Edge Enhance |
| ▲▼ SGL      | ▲▼ 0DEG      | ▲▼ 40   | ▲▼ MID       |

Figure 15–10. B-Mode Sub-Menu (Biopsy Zone)

If the single fixed angle (SGL) was selected as the needle guide type in Set Up/Custom Display page 1, the fixed zone angle is displayed.

If the needle guide type selected is multiple angle (MULTI), the Biopsy Zone rocker switch can be used to cycle through the angle selections MBX1, MBX2, MBX3 and OFF.

The Biopsy Guidezone adjusts along with image adjustments, such as image inversion/rotations, zoom and depth changes.




Figure 15–11. Biopsy Guidezone Example






## Preparing for a Biopsy

### Determining Needle Length



**NOTE:** Press the **Measurement** key once to display the integrated biopsy depth cursor and center line while the guidezone is present. Use the **Trackball** to measure the needle length. Needle length is from the top of needle barrel to the target.



- 1 Needle Length
- 2 Target

Figure 15–12. Needle Length Measurement

### CAUTION



Pressing measurement while the biopsy guidezone is displayed will add a line down the center of the guidezone with a depth marker that can be adjusted with the Trackball.

The needle may vary from the center line or guidezone for various reasons:

- Needle barrel to needle clearance or strength.
- Bracket manufacturing tolerance.
- Needle deflection due to tissue resistance.
- Needle size chosen. Thinner needles may deflect more.
- Wrong needle guide type selection in Set Up/Custom Display page 1.

The display should be carefully monitored during a biopsy for any needle deviation from the center line or guidezone.

## Needle Guide Type Preset Selection

A preset is provided in Set Up/Custom Display page 1 to choose single (SGL) or Multiple (MBX) guidezone angle capabilities. Available selections for respective probes are shown in the chart below.

|          | Depth in cm at Center Channel |                   |      |      |
|----------|-------------------------------|-------------------|------|------|
|          | Fixed-Angle                   | Multi-Angle (MBX) |      |      |
| Probe    | SGL                           | MBX1              | MBX2 | MBX3 |
| C358     | n/a                           | 4.0               | 6.0  | 8.0  |
| C364/CBF | 8.0                           | 4.0               | 6.0  | 8.0  |
| C386     | n/a                           | 4.0               | 6.0  | 8.0  |
| C551/CAE | 7.0                           | 4.0               | 5.5  | 7.0  |
| C721     | 3.0                           | 2.0               | 3.0  | 4.5  |
| 546L     | n/a                           | 4.0               | 5.5  | 7.0  |
| L764/LH  | 2.0                           | 2.0               | 4.0  | 7.0  |
| 739L     | 2.0                           | 2.0               | 4.0  | 6.0  |
| LA39     | 1.5                           | n/a               | n/a  | n/a  |
| S316/UC  | 8.0                           | n/a               | n/a  | n/a  |
| S317     | 8.0                           | 4.0               | 6.0  | 8.0  |
| S611     | n/a                           | 4.0               | 5.5  | 7.0  |
| T739     | 2.0                           | 2.0               | 4.0  | 6.0  |

Table 15–8. Biopsy Guide Attachment Selection

### DANGER



Failure to match the guidezone displayed to the guide may cause the needle to track a path outside the zone.

It is extremely important that when using the adjustable angle biopsy guides, the angle displayed on the screen matches the angle set on the guide, otherwise the needle will not follow the displayed guidezone which could result in repeated biopsies or patient injury.



**NOTE:** Although the multi-angle guides are compatible with the Civco Ultrapro and Ultrapro II, it is recommended the multi-angle guides only be used with the Ultrapro II.



## Needle Guide Type Preset (cont'd)




Figure 15–13. Biopsy Guide Depth Selection

### DANGER



When the biopsy guidezone is displayed, the message:

“Confirm BX type of Bracket”


is displayed at the bottom of the screen with the angle selected.

Ensure that the MBX type selected for each probe is the same as the angle selected on the actual biopsy guide.


## E721 Type Selection

When the E721 probe is attached and active, the needle guide type selection choices in Set Up/Custom Display page 1 are:

- TV0° Reusable metal guide with a 0° offset angle.
- TR5° Civco disposable guide with a 5° offset angle.



TV0°



TR5°

Figure 15–14. Biopsy Guides

## ERB7 Type Selection

When the ERB7 probe is attached and active, the needle guide:

- BX Disposable guide with 5 mm increments for vertical orientation (1~8 selections are available).
- GBX Needle Placement guide with 5 mm increments for seed implant template grid alignment.

## LD Type Selection

| LD Probe Angle (degree) |              |      |              |
|-------------------------|--------------|------|--------------|
| MBX1                    | MBX2         | MBX3 | MBX4         |
| -30 (30)                | -16.1 (16.1) | 0    | 16.1 (-16.1) |

Table 15–9. LD Probe Angle




Figure 15–15. LD Probe Alignment




Figure 15–16. LD Probe with Needle Guide



## Preparing for a Biopsy

---

### Preparing the Biopsy Guide Attachment

Convex, Sector and Linear probes have optional biopsy guide attachments for each probe. The guide consists of a non-disposable bracket to attach to the probe, disposable needle clip to attach to the bracket, sheath, gel (sterile gel if necessary) and disposable needle barrels.

The disposable needle barrels are available for a variety of needle sizes.

#### CAUTION



Please refer to the manufacturer's instructions included in the biopsy kit.

### Fixed Needle Guide Assembly

Identify the appropriate biopsy guide bracket by matching the label on the bracket with the probe to be used.

Orient the bracket so that the needle clip attachment will be on the same side as the probe orientation mark (ridge).




Figure 15–17. Probe/Bracket Alignment

Attach the biopsy bracket to the probe by sliding the bracket over the end of the probe until it clicks or locks in place.

### Fixed Needle Guide Assembly (cont'd)

Place an adequate amount of coupling gel on the face of the probe.

Place the proper sanitary sheath over the probe and biopsy bracket. Use the rubber bands supplied to hold the sheath in place.





Figure 15–18. Applying Sanitary Sheath

Snap the fixed or adjustable needle clip onto the biopsy guide bracket.



1 Sheath

Figure 15–19. Fixed Needle Clip Attachment

Choose the desired gauge (size) needle barrel. Twist it back and forth to remove it from the plastic tree.



## Preparing for a Biopsy

---

### Fixed Needle Guide Assembly (cont'd)




Figure 15-20. Needle Barrels

Place the needle barrel into the needle clip with the desired gauge facing the needle clip and snap into place.




Figure 15-21. Needle Barrel Installation

### CAUTION



Ensure that all guide parts are seated properly prior to performing a biopsy.

## Multi-Angle Biopsy Guide Assembly

### WARNING



*DO NOT attempt to use the multi-angle biopsy bracket and needle guide until the manufacturer's instructions, provided with the biopsy bracket and needle guide in the kit, have been read and thoroughly understood.*

Scan the patient and identify the target for biopsy. Move the probe to locate the target to the center of the image. Enable the system biopsy guidezone and try guidezone angles MBX1 to MBX3 to decide the best angle setting for needle path.



*NOTE: The biopsy guidezone type can be selected in the Set Up/Custom Display menu page 1.*

Identify the appropriate biopsy guide bracket by matching the label on the bracket with the probe to be used.




Figure 15–22. S317 Multi-Angle Biopsy Guide Bracket

Orient the bracket so that the needle clip attachment will be on the same side as the probe orientation mark (ridge).




Figure 15–23. Probe/Bracket Alignment

Attach the biopsy bracket to the probe by sliding the bracket over the end of the probe until it clicks or locks in place.





## Preparing for a Biopsy

---

### Multi-Angle Biopsy Guide Assembly (cont'd)

Pull up on the knob to freely move the needle guide attachment. Align the knob with the selected position of the needle guide attachment from MBX1, MBX2 and MBX3, to match the guidezone display on the ultrasound system.




Figure 15–24. Select the angle position

Push the knob down into the desired slot to secure the angle position of the needle guide attachment.




Figure 15–25. Fix the angle position

### CAUTION



Hold the bracket in place on the probe when pushing the knob to secure the angle position of the needle guide attachment. Excessive force may cause the bracket to release from the probe.

### Multi-Angle Biopsy Guide Assembly (cont'd)

Place an adequate amount of coupling gel on the face of the probe.

Place the proper sanitary sheath tightly over the probe and biopsy bracket. Use the rubber bands supplied to hold the sheath in place.





Figure 15–26. Applying Sanitary Sheath

Snap the needle clip onto the biopsy guide bracket.



1 Sheath

Figure 15–27. Fixing the Needle Clip Attachment



## Preparing for a Biopsy

### Multi-Angle Biopsy Guide Assembly (cont'd)

Push the locking mechanism towards the bracket to secure the lock. Make sure the needle guide is firmly attached to the bracket.




Figure 15-28. Locking the Needle Clip

Choose the desired gauge (size) needle barrel. Twist it back and forth to remove it from the plastic tree.




Figure 15-29. Needle Barrel Selection

Place the needle barrel into the needle clip with the desired gauge facing the needle clip and snap into place.




Figure 15-30. Needle Barrel Installation

### CAUTION



Ensure that all guide parts are seated properly prior to performing a biopsy.

## Biopsy Need Path Verification

To verify that the path of the needle is accurately indicated within the guidezone on the system monitor, perform the following:

- Properly install the bracket and biopsy guide.
- Scan in a container filled with water.
- Display the biopsy guidezone on the monitor.
- Ensure that the needle echo falls within the guidezone markers.



**NOTE:** The guidezone center line can be displayed by pressing the **Measurement** key.

## The Biopsy Procedure

Place coupling gel on the scanning surface of the probe/sheath/biopsy guide assembly.

Activate the biopsy guidezone on the system through the B-Mode Sub-Menu. When using multi-angle guides, ensure that the proper guidezone angle is displayed.

| B           | Preset       | Set Up  | ECG          |
|-------------|--------------|---------|--------------|
| Biopsy Zone | Image Rotatn | Rejectn | Edge Enhance |
| ▲▼ SGL      | ▲▼ 0DEG      | ▲▼ 40   | ▲▼ MID       |

Figure 15–31. B-Mode Sub-Menu (Biopsy Zone)

Scan to locate the target. Center the target in the electronic guidezone path.



**NOTE:** Enabling color flow would allow for visualization of the vascular structure around the area to be biopsied.

Place the needle in the guide between the needle barrel and needle clip. Direct it into the area of interest for specimen retrieval.



## Preparing for a Biopsy

---

### Post Biopsy

When the biopsy is complete, remove the needle barrel, needle clip and probe sheath. Properly dispose of these items in accordance with current facility guidelines.


The biopsy bracket can be sterilized in a recommended disinfecting agent and reused.

### E721 Probe Biopsy Guide

#### Preparation

To prepare the E721 for use:

1. Remove the probe from the box and carefully examine it for any damage.
2. If the biopsy guide is to be attached, use the filling removal tool (Figure 15–32) to clean out the attachment area on the probe head.



1 Probe Head

2 Attachment

3 Filling Removal Tool

Figure 15–32. Attachment Filling Removal

3. Clean, then disinfect/sterilize the probe.



**NOTE:** Ensure that protective gloves are worn.

## Preparation (cont'd)

Install the sheath:

1. Remove the sheath from its package. Do not unroll the sheath.




*NOTE: Remember to rinse all sanitary probe sheaths of powder before placing on the probe. Powder can degrade the displayed image.*

2. Place an adequate amount of ultrasound gel inside the sheath tip (the gel is between the sheath inner surface and the probe aperture).



*NOTE: Ensure that only acoustic coupling gel is used for this purpose.*

3. Place the sheath tip over the probe aperture and then pull the sheath end toward the probe handle.
4. Inspect the sheath for nicks, cuts or tears.



1 Probe Handle

2 Probe Body

3 Sanitary Sheath

Figure 15–33. Probe with Sheath




## Preparing for a Biopsy

---

### Preparation (cont'd)


5. Rub a finger over the tip of the probe to ensure all air bubbles have been removed.



1 Biopsy Guide

Figure 15–34. Reusable Metal Biopsy Guide 0° Angle

6. If a biopsy is to be performed, snap the metal or plastic biopsy guide on to the probe over the sheath (Figure 15–34).
7. Place an adequate amount of ultrasound gel on the gel-filled sheath tip **outer** surface.



1 Sanitary Sheath

Figure 15–35. Civco Disposable Biopsy Guide 5° Angle

8. Ensure the guide is properly seated and secure by pushing forward on the needle insertion end of the guide until the attachment node is firmly in place in it's hole.

## Scanning

1. Scan the patient. The probe handle orientation mark indicates the image scan plane. Be sure that the **Image Reverse** function is **Off**.




Figure 15–36. Probe and Guidezone Alignment

2. Rotate, retract, or advance the probe, as necessary, to see all pertinent anatomy.
3. If a biopsy is being performed, activate the biopsy guidezone.

## CAUTION



Scan the patient to determine the correct puncture depth and site **before** inserting the needle.

Ensure that the proper guidezone (TV0 or TR5) is displayed.

## Post Biopsy

If the exam is over:

1. Remove the biopsy guide and twist lock/clamp. Remove and properly dispose of the sheath.
2. Thoroughly clean the probe and equipment. Refer to your institution's infection control guidelines for disinfection/sterilization protocols.
3. After sterilization, return the probe to its carrying case.





## Preparing for a Biopsy

### Biopsy Probes

| Biopsy Guide Attachment |                                  |       |                            |                 |
|-------------------------|----------------------------------|-------|----------------------------|-----------------|
| Probe                   | Application                      | Mfg.  | Needle Sizes               | Repl. Kit       |
| C364/CBF                | Abdominal, OB/GYN                | Civco | 14, 16, 18, 20, 22, 25 AWG | Civco Ultra-Pro |
| C551/CAE                | OB/GYN, Abdominal                | Civco | 14, 16, 18, 20, 22, 25 AWG | Civco Ultra-Pro |
| C721                    | Neonatal                         | Civco | 14, 16, 18, 20, 22, 25 AWG | Civco Ultra-Pro |
| E721                    | Endocavitary                     | YMS   | Requires 25 cm needles     | H46222AD        |
|                         |                                  | Civco | 18 AWG                     | H4550BG         |
| ERB7                    | Endocavitary                     | YMS   |                            |                 |
|                         |                                  | Civco |                            |                 |
| 739L                    | Small Parts, Peripheral Vascular | Civco | 14, 16, 18, 20, 22, 25 AWG | Civco Ultra-Pro |
| L764                    | Breast, Small Parts              | Civco | 14, 16, 18, 20, 22, 25 AWG | Civco Ultra-Pro |
| LA39                    | Small Parts, Mammography         | Civco | 14, 16, 18, 20, 22, 25 AWG | Civco Ultra-Pro |
| S316/UC                 | Abdominal                        | Civco | 14, 16, 18, 20, 22, 25 AWG | Civco Ultra-Pro |
| S317                    | Abdominal                        | Civco | 14, 16, 18, 20, 22, 25 AWG | Civco Ultra-Pro |
| T739                    | Intra-operative                  | Civco | 14, 16, 18, 20, 22, 25 AWG | Civco Ultra-Pro |

Table 15–10. Probe Biopsy Chart—Fixed Angle

| Multi Angle Biopsy Guide Attachment |                                  |       |                                |                    |
|-------------------------------------|----------------------------------|-------|--------------------------------|--------------------|
| Probe                               | Application                      | Mfg.  | Needle Sizes<br>(not included) | Repl. Kit          |
| C358                                | Abdominal, OB/GYN                | Civco | 14 through 23 AWG              | Civco Ultra-Pro II |
| C364                                | Abdominal, OB/GYN                | Civco | 14 through 23 AWG              | Civco Ultra-Pro II |
| C386                                | OB/GYN, Abdominal                | Civco | 14 through 23 AWG              | Civco Ultra-Pro II |
| C551                                | OB/GYN, Abdominal                | Civco | 14 through 23 AWG              | Civco Ultra-Pro II |
| C721                                | Small Parts                      | Civco | 14 through 23 AWG              | Civco Ultra-Pro II |
| 546L                                | Small Parts, Peripheral Vascular | Civco | 14 through 23 AWG              | Civco Ultra-Pro II |
| 739L                                | Small Parts, Peripheral Vascular | Civco | 14 through 23 AWG              | Civco Ultra-Pro II |
| L764                                | Breast, Small Parts              | Civco | 14 through 23 AWG              | Civco Ultra-Pro II |
| S317                                | Abdominal, Cardiac               | Civco | 14 through 23 AWG              | Civco Ultra-Pro II |
| S611                                | Abdominal, Pediatric             | Civco | 14 through 23 AWG              | Civco Ultra-Pro II |
| T739                                | Intra-operative                  | Civco | 14 through 23 AWG              | Civco Ultra-Pro II |

Table 15–11. Probe Biopsy Chart—Multi Angle

# System Data

## Features/Specifications

### System Dimensions

- Height: 47–59 in (120–150cm)
- Width: 21 in (54 cm)
- Depth: 37 in (94 cm)
- Weight: approx. 397 lbs (180 kg)

### Electrical Power

- Voltage: 100–120 V<sub>AC</sub> (220–240 V<sub>AC</sub> available)
- Frequency: 60 Hz (50 Hz available)
- Power 1350 VA Max

### System Operating/Display Modes

- B-Mode
- M-Mode
- PW Doppler
- HPRF Doppler
- Color flow imaging
- Color flow M-Mode
- Power Doppler imaging
- Steerable CW Doppler
- Non-imaging CW Doppler
- Real Time Triplex simultaneous B/D/CFI and B/D/PDI
- *SmartTrac* cursor control
- Colorized B, M and PW Doppler Display
- User selectable/programmable presets

### Transducer Types

- Phased array sector
- Convex array
- Linear array
- Micro-convex array
- Endocavitary (TV/TR)
- Non-imaging CW Doppler

### Display Monitor

- Single 12" color display
- NTSC color formats (PAL available)
- Full swivel/tilt/height adjust monitor

### Operator Interface

- Full alphanumeric keyboard
- Context sensitive backlit keys
- Three active probe ports
- Dynamic soft menu display
- User definable keys
- Body patterns

### Peripheral Interface

- Multi-image cameras
- Color thermal printers
- B/W thermal printer
- S-VHS video recorder
- Laser cameras

### Imaging Formats

- Sector
  - Variable angle
  - Steerable, selectable FOV
- Linear
  - Steerable Doppler and CF
- Convex
  - Variable angle
  - Steerable, selectable FOV
- Micro-convex
  - Variable angle
  - Steerable, selectable FOV
- Image reverse and rotate
- Split screen for B, M, D, CFI & PDI

### Digital Image Processing

- Micron imaging
- 256 shades of gray
- Adjustable dynamic range
- *SmartZoom* acoustic and display zooms
- Dynamic focusing
- Multi-frequency
- User selectable gray maps
- Parallel processing



## System Data

---

### Features/Specifications (cont'd)

#### Digital Doppler

- Steerable sample volume
- Variable velocity scale
- Variable sample size
- Variable wall filter
- Scrolling memory
- Doppler angle-correction with calibrated velocity scaling
- Dual Doppler display
- Spectrum invert
- Stereo audio output
- HPRF

#### Digital Color Flow

- Adaptive color enhancement (ACE)
- MTI motion suppression
- Variable velocity scale
- Color velocity tag
- Frame averaging
- Color baseline shift
- Color capture
- Selectable color maps
- Variable wall filter
- Variable FOV size and position

#### Digital Power Doppler

- Adaptive color enhancement (ACE)
- Variable FOV size and position
- Selectable color maps
- Transparent maps
- Topographic maps

#### InSite™ Remote Diagnostic Capability

- Live application support
- Remote execution of system diagnostics
- *Virtual Operator* capability
- Image transfer
- Error log analysis
- Scan data analysis
- Keystroke analysis

#### Measurements

- Distance
- Circumference (Ellipse/Trace)
- Area (Ellipse/Trace)
- Volume (Ellipsoid)
- Angle Between Two Lines (B-Mode)
- Time & Slope (M-Mode)
- Single Velocity (Doppler)
- Time Average Maximum Velocity

#### Calculations

- OB Calculations package
- Vascular Calculations package
- Cardiac Calculations package
- Advanced OB package
- Advanced Vascular Calculations package
- Advanced Cardiac Calculations package

#### Other Features

- Measurements/Calculations and annotations on cine and VCR
- Dual image cine display
- Variable speed cine display
- VCR playback auto calibration
- Acoustic output displays: MI/TIS/TIB/TIC/AP
- Cine memory 36–160 frames
- ECG/Physio module (optional)



## LOGIQ™ 500 Clinical Measurement Accuracy

### Basic Measurements

The following information is intended to provide guidance to the user in determining the amount of variation or measurement error that should be considered when performing clinical measurements with this equipment. Error can be contributed by equipment limitations and improper user technique. Be sure to follow all measurement instructions and develop uniform measurement techniques among all users to minimize the potential operator error. Also, in order to detect possible equipment malfunctions that could affect measurement accuracy, a quality assurance (QA) plan should be established for the equipment that includes routine accuracy checks with tissue mimicking phantoms.

Please be advised that all distance and Doppler related measurements through tissue are dependent upon the propagation velocity of sound within the tissue. The propagation velocity usually varies with the type of tissue, but an average velocity for soft tissue is assumed. This equipment is designed for, and the accuracy statements listed on 16–6 are based on, an assumed average velocity of 1540 m/s. The percent accuracy when stated applies to the measurement obtained (not the full scale range). Where the accuracy is stated as a percent with a fixed value, the expected inaccuracy is the greater of the two.



## System Data

### Basic Measurements (cont'd)

| Measurement                                         | Units                              | Useful Range                                             | Accuracy                                                                         | Limitations or Conditions                       |
|-----------------------------------------------------|------------------------------------|----------------------------------------------------------|----------------------------------------------------------------------------------|-------------------------------------------------|
| Depth                                               | mm                                 | Full Screen                                              | $\pm 5\%$ or 1 mm                                                                |                                                 |
| Distance:<br>Axial<br>Lateral<br>Lateral<br>Lateral | mm<br>mm<br>mm<br>mm               | Full Screen<br>Full Screen<br>Full Screen<br>Full Screen | $\pm 5\%$ or 1 mm<br>$\pm 5\%$ or 2 mm<br>$\pm 5\%$ or 4 mm<br>$\pm 5\%$ or 4 mm | Linear Probes<br>Convex Probes<br>Sector Probes |
| Circumference:<br>Trace<br>Ellipse                  | mm<br>mm                           | Full Screen<br>Full Screen                               | $\pm 10\%$ or 1 mm<br>$\pm 5\%$ or 1 mm                                          |                                                 |
| Area:<br>Trace<br>Ellipse                           | mm <sup>2</sup><br>mm <sup>2</sup> | Full Screen<br>Full Screen                               | $\pm 5\%$ or 1 mm <sup>2</sup><br>$\pm 5\%$ or 1 mm <sup>2</sup>                 |                                                 |
| Time                                                | s                                  | Timeline Display                                         | $\pm 5\%$ or 10 ms                                                               | M or Doppler                                    |
| Slope                                               | mm/s                               | Timeline Display                                         | $\pm 5\%$ or 1 mm/s                                                              | M-Mode Only                                     |
| Doppler SV Position                                 | mm                                 | Full Screen                                              | $\pm 2$ mm                                                                       | Any Direction                                   |
| Velocity                                            | cm/s                               | From 0 to 100 cm/s<br>From 100 to 130 cm/s               | $\pm 10\%$ or 1 cm/s<br>$\pm 5\%$ or 1 cm/s<br>$\pm 50\%$                        | PW Doppler<br>Color Flow                        |
| Doppler Angle<br>Correction                         | cm/s<br>cm/s                       | From 0° to 60°<br>From 60° to 80°                        | $\pm 5\%$<br>$\pm 12\%$                                                          |                                                 |

Table 16–1. System Measurements and Accuracies

### LOGIQ™ 500 Clinical Calculation Accuracy

Estimate the overall inaccuracy of a combined measurement and calculation by including the stated inaccuracy from the basic measurement accuracy statements.

#### CAUTION



Diagnostic errors may result from the inappropriate use of clinical calculations. Review the referenced source of the stated formula or method to become familiar with the intended uses and possible limitations of the calculation.



Calculation formulas and databases are provided as a tool to assist the user, but should not be considered an undisputed database, in making a clinical diagnosis. The user is encouraged to research the literature and judge the equipment capabilities on an ongoing basis in order to assess its utility as a clinical tool.



# Warranties

## Scope and Duration of Warranties

### Product warranties

GE warrants that the ultrasound products are:

1. Free from defects in material, workmanship, and title.
2. Conform to our published product specifications in effect on the date of shipment of the products. The product specifications are available on request.

### Patent and copyright warranty

GE warrants that when delivered, the products will not be subject to any valid patent or copyright infringement claim.

The warranty period for all warranties, except the warranty of title and the patent and copyright warranty, is limited in time as shown below:

| Ultrasound Part                                       | Time Limit |
|-------------------------------------------------------|------------|
| Ultrasound systems, components, modules, and upgrades | 12 Months  |
| Ultrasound probes and transducers                     | 12 Months  |
| Ultrasound water path attachment kit                  | 3 Months   |

Table 16–2. Warranty Time Limits



*NOTE: If a maintenance agreement is desired to extend coverage beyond the manufacturer's warranty, please contact the local service or sales representative.*



## Warranties

---

### Patent and copyright warranty (cont'd)

The warranty period begins on the date the products are delivered. But, if GE assembles the products, the warranty period begins on the earlier of:

1. Five (5) days after the date GE notifies you that we have completed assembly and the products are operating in accordance with our published product specifications

OR

2. The date you first use the products for patient use.

If assembly is delayed for 30 days or more after the date of delivery for a reason beyond our reasonable control, the warranty period will begin on the thirtieth day after the date of delivery.

The warranty period for any product or part furnished to correct a warranty failure will be the unexpired term of the warranty applicable to the repaired or replaced products.

## Warranty Exclusions

These warranties are exclusive and in lieu of all other warranties, whether written, oral, expressed, implied, or statutory.



*No warranty of merchantability or fitness for a particular purpose applies.*

The warranties **DO NOT** cover:

1. Any defect or deficiency (including failure to conform to product specifications) which results, in whole or in part, from:
  - a. Any alteration, improper storage, handling, use or maintenance, or any extraordinary use of the products, by anyone other than GE.
  - b. Failure to follow any of GE's written recommendations or instructions.
  - c. Combining the products with products of others or with incompatible GE products.
  - d. Any of your designs, specifications or instructions.

OR

  - e. Any cause external to the products as furnished by GE or beyond our reasonable control.
2. Ultrasound supplies and accessories identified by catalog numbers which start with the letter "E" (which are covered by a separate printed warranty).
3. Products which are not listed in our price pages at the time of sale. Non-listed products are provided with the manufacturer's warranties, if any, GE is permitted to pass on to you. Otherwise, non-listed products are provided AS IS.
4. The payment and reimbursement of any facility costs arising from repair or replacement of the products or parts.





## Warranties

---

### Exclusive Warranty Remedies

#### Product warranties

If you promptly notify us of your warranty claim and make the product available for service, GE will at our option, repair, adjust, or replace (with new or exchange replacement parts) the non-conforming product or parts of the product. Warranty service will be performed without charge from 8:00 am to 5:00 pm, Monday through Friday, excluding our holidays, and outside those hours at our prevailing service rates and subject to the availability of personnel.

#### Patent and copyright warranty

GE will defend or settle any suit against you to the extent it is based on an infringement claim which would be a breach of the patent and copyright warranty. If the infringement claim is valid, GE will pay all damages and costs awarded against you due to the breach. In addition, GE will obtain a license for you to continue using the infringing product, provide a non-infringing replacement, alter the product so that it is non-infringing, or remove the infringing product and refund the price (less reasonable depreciation) and any return transportation costs paid by you.

The above describes your exclusive remedies and our sole liability for any warranty claims. You agree that GE and our representatives have no liability to you for:

1. Any penal, incidental, or consequential damages such as lost profit or revenue.
2. Any assistance not required under our quotation.
3. Anything occurring after the warranty period ends.



# System Care and Maintenance

## Overview

Refer to Section 7 of the LOGIQ™ 500 Service Manual (P9030TA) for any additional maintenance guidance.

Contact the local Service Representative for parts or planned maintenance inspections. It is recommended that planned maintenance be performed on the system every six months.

## Inspecting the System



Examine the following on a monthly basis:

- Connectors on cables for any mechanical defects.
- Entire length of electrical and power cables for cuts or abrasions.
- Equipment for loose or missing hardware.
- Control panel and keyboard for defects.
- Casters for proper locking operation.

### DANGER



To avoid electrical shock hazard, do not remove panels or covers from console. This servicing must be performed by qualified service personnel. Failure to do so could cause serious injury.



### Electrical Hazard

If any defects are observed or malfunctions occur, do not operate the equipment but inform a qualified service person. Contact a Service Representative for information.



## System Care and Maintenance

---

### Weekly Maintenance



The LOGIQ™ 500 system requires weekly care and maintenance to function safely and properly. Clean the following:

- System cabinet
- Monitor
- Operator control panel
- Foot switch
- Video Cassette Recorder (VCR)
- Multi Imaging Camera (MIC)
- Video Page Printer

Failure to perform required maintenance may result in unnecessary service calls.

### Cleaning the system

Prior to cleaning any part of the system:

1. Turn off the system power.

#### System cabinet



To clean the system cabinet:

1. Moisten a soft, non-abrasive folded cloth with a mild, general purpose, non-abrasive soap and water solution.
2. Wipe down the top, front, back, and both sides of the system cabinet.



**NOTE:** Do not spray any liquid directly into the unit.

## Cleaning the system (cont'd)

### Monitor

To clean the monitor face and filter:

Remove the monitor filter as shown in Figure 16-1.

Daily  
Weekly  
Monthly  
Other

☐ ☒ ☐ ☐




Figure 16-1. Monitor Filter Removal

1. Slide the filter clamps towards the monitor sides.
2. Pull the filter clamps out.



## System Care and Maintenance

---

### Cleaning the system (cont'd)

#### Monitor (cont'd)

Use a soft, folded cloth and a glass cleaner solution. Apply the glass cleaner to the cloth. Gently wipe the monitor face and filter.

#### CAUTION



*Do NOT use a glass cleaner that has a hydrocarbon base (such as Benzene, Methyl Alcohol or Methyl Ethyl Ketone) on monitors with the filter (anti-glare shield). Prolonged use of such cleaners will damage the filter (anti-glare shield). Hard rubbing will also damage the filter.*



*NOTE: When cleaning the monitor, make sure not to scratch the monitor.*

After cleaning the monitor face and filter, re-install the filter on the monitor by reversing the removal procedure in Figure 16–1. Insert the filter clamps and slide them inwards securely.



*NOTE: Make sure that the filter is securely fixed by the filter clamps. These are important to prevent a hazard from the filter falling. A “click” sound will be heard when the filter clamp is secured by the lock mechanism.*

#### Operator Controls

To clean the operator control panel:



1. Moisten a soft, non-abrasive folded cloth with a mild, general purpose, non-abrasive soap and water solution.
2. Wipe down operator control panel.
3. Use a cotton swab to clean around keys or controls. Use a toothpick to remove solids from between keys and controls.



*NOTE: When cleaning the operator control panel, make sure not to spill or spray any liquid on the controls, into the system cabinet, or in the probe connection receptacle.*



## System Care and Maintenance

---

### Cleaning the system (cont'd)

#### Foot Switch



To clean the foot switch:

1. Moisten a soft, non-abrasive folded cloth with a mild, general purpose, non-abrasive soap and water solution.
2. Wipe the external surfaces of the unit then dry with a soft, clean, cloth.

#### VCR



To clean the VCR:

1. Turn off the VCR power. If possible, disconnect the power cord.
2. Wipe the external surfaces of the unit with a soft, clean, dry cloth.



**NOTE:** *Do not use a wet cloth or any cleaning fluid because it may enter and damage the unit.*

3. Clean the record and playback heads with a soft, non-abrasive cleaning system, according to the manufacturer's instructions.

#### For more information

See the VCR's Operator Manual.

#### Multi Image Camera



To clean the MIC:

1. Moisten a soft, non-abrasive folded cloth with a mild, general purpose, non-abrasive soap and water solution.
2. Wipe down the top, front, back, and both sides of the unit.

#### For more information

Review the Multi Image Camera's Operator Manual for details.



## System Care and Maintenance

---

### Cleaning the system (cont'd)

#### Video Page Printer

To clean the external surface of the video page printer:



1. Turn off the power. If possible, disconnect the power cord.
2. Wipe the external surfaces of the unit with a soft, clean, dry cloth.
3. Remove stubborn stains with a cloth lightly dampened with a mild detergent solution.



*NOTE: Never use strong solvents, such as thinner or benzine, or abrasive cleansers because they will damage the cabinet.*

No further maintenance, such as lubrication, is required.

To clean the surface of the print head:

1. Run the cleaning sheet (provided with the printer) through the printer.

**For more information**

Review the Video Page Printer's Operator Manual for details.

### Other Maintenance

#### Replacing illuminated key caps/lamps



Contact a local Service Representative when a key cap or lamp needs to be replaced.

As Necessary

## Cleaning the air filters

Daily  
Weekly  
Monthly  
Other

Quarterly

Clean the system's air filters to ensure that a clogged filter does not cause the system to overheat and reduce system performance and reliability. It is recommended the filters be cleaned *quarterly* (once every three months).

## Locating

The 3 air filters are located in the front and back of the system:




Figure 16–2. Location of Air Filters





## System Care and Maintenance

---

### Removing

To remove the console air filters:

1. Lift up on the Velcro tab securing the air filter.
2. Pull out the filter.

To remove the power supply air filter:

1. Pull out and up on the filter cover.
2. Pull the filter screen out.




Figure 16-3. Single Filter Screen Removal

### Cleaning

To clean the filter:

1. Be sure to shake the filter in an area away from the system.
2. Wash the filter in a mild soapy solution, rinse and air dry or dry with a cloth.



*NOTE: Allow wet filter to dry thoroughly before installing.*

3. Slide the filter back into the system.
4. Secure the filter with the Velcro tab or cover.



# Troubleshooting

## Introduction

Listed in this section are problem or system messages that may be encountered, possible causes for the problem or message, and the appropriate action to take to correct the situation. If additional information or assistance is needed contact your local Applications, Sales or Service Representative.

## Trouble images

The LOGIQ™ 500 has the capacity to store as many as eight trouble images. The trouble image files may be saved to MOD for archival or loaded back into the system hard drive for display and review. The trouble images on file can also be retrieved for analysis through Insite™, the exclusive remote diagnostic capabilities of the LOGIQ™ family of products.



### To Store a Trouble Image to the Hard Drive

Should the LOGIQ™ 500 ever have a problem imaging, the system has the capability of storing up to 8 trouble images on the hard disk drive. To store a trouble image:

1. Freeze the trouble image.
2. Press the **Ctrl, W** keys simultaneously.

The message *"In progress. Please wait."* is displayed while the image is being saved to the hard drive.

After the image has been saved, the message "Input comment? 'y' or 'n' " is displayed. Press '**n**' for no comments. Press '**y**' to add a maximum of 42 characters as comments to the trouble image file. Press **Return** after commenting is complete.

After the eighth image is saved, there is a warning to indicate the disk space is full. The ninth attempt to save a trouble image will require writing over the first image, etc. (first in, first out).

The message *"Overwrite the old image (DIAGIG00) ? 'y' or 'n'"* is displayed. Press '**y**' to replace the old trouble image with the new one. Press '**n**' to cancel the trouble image store action.



*NOTE: Trouble images do not display Patient Name, ID or Hospital Name information. That portion of the image is blanked out and will not be shown when the file is recalled for display.*

## To Save Trouble Images to MOD

To save trouble images to MOD:

- Select the Set Up Top menu.
- From the Set Up Sub-Menu, select Utility. The Utility Menu appears as shown in Figure 16–4.



*NOTE: MOD media used to save trouble images must first be initialized (formatted) in the utility menu selection “09 Media Initialization”.*




Figure 16–4. Utility Menu

- From the Utility menu, select “03 Trouble Image Save/Load/Display”. Enter menu number 03 and press **Return**.



## Troubleshooting

### To Save Trouble Images to MOD (cont'd)

- From the Trouble Image Save/Load/Display menu select "1: Save". Enter menu number 1 and press **Return**.




Figure 16–5. Trouble Image Save/Load/Display Menu

- A message will be displayed to place an MOD disk in the drive and press any key to continue.

The system will read all trouble images from the hard drive and display the image menu in 2 pages of 4 files each. The line available under the file number, date/time and mode information is the 42 character comment line.




Figure 16–6. Saved Trouble Image Menu



To Save Trouble Images to MOD (cont'd)

Trouble image files must be saved one at a time.

- Enter the number of the file (i.e. "01") and press **Return**.

The message "Saving now" is displayed.

The system will take a few moments to save the image file.

After the save operation is complete, the message "Press 'Y' to continue, or 'Ctrl+R' to return." is displayed. Pressing 'Y' will prompt for another file menu number to be saved. Pressing Ctrl+R will return to the previous menu screen.

```
HOSPITAL NAME                               98/05/09 02:24:48
Trouble Image Save/Load/Display              S317
      1:Save   2:Load   3:Display

Enter Menu Number ?1
Set a media to drive and then press any key.

No  Data Name      Date      Mode      (1/1)
01  DIAGIG00.VR   1998/05/09 01:08      B
    TESTING
02  DIAGIG01.VR   1998/05/09 01:09      BD
    TEST 2
03  DIAGIG02.VR   1998/05/09 01:10      BD
    TEST 3
04  DIAGIG03.VR   1998/05/09 01:10      BM
    TEST 4

Enter Menu Number ?1
Saving now
Succeeded
Press 'Y' to continue, or 'Ctrl+R' to return._
```

Figure 16–7. Trouble Image Save Complete



### To Load Trouble Images from MOD

Trouble images can be loaded from the MOD to the hard drive.  
To load trouble images from the MOD:

- From the Set Up/Utility menu select *"03 Trouble Image Save/Load/Display"*. Enter menu number 03 and press **Return**.
- From the Trouble Image Save/Load/Display menu select *"2: Load"*. Enter menu number 2 and press **Return**.
- A message will be displayed to place a MOD disk in the drive and press any key to continue.

The system will read all trouble images from the MOD and display the image menu in 2 pages of 4 files each. The line available under the file number, date/time and mode information is the 42 character comment line.

Trouble image files must be loaded one at a time.

- Enter the number of the file to be loaded (i.e. "01") and press **Return**.

The message *"Loading now"* is displayed.

The system will take a few moments to load the image file.

After the load operation is complete, the message *"Press 'Y' to continue, or 'Ctrl+R' to return."* is displayed. Pressing 'Y' will prompt for another file menu number to be loaded. Pressing Ctrl+R will return to the previous menu screen.

### To Display Trouble Images

Trouble images can only be displayed from the files stored on the hard drive.

- From the Set Up/Utility menu select *"03 Trouble Image Save/Load/Display"*. Enter menu number 03 and press **Return**.
- From the Trouble Image Save/Load/Display menu select *"3: Display"*. Enter menu number 3 and press **Return**.

A menu of the trouble Image files on the hard drive is displayed in 2 pages of 4 files each if necessary.



## To Display Trouble Images (cont'd)

```
HOSPITAL NAME          98/05/09 02:25:22
Trouble Image Save/Load/Display      S317
    1: Save    2: Load    3: Display

Enter Menu Number ?3

No  Data Name      Date      Mode      (1/1)
01  DIAGIG00.VR  1998/05/09 01:08    B
    TESTING
02  DIAGIG01.VR  1998/05/09 01:09    BD
    TEST 2
03  DIAGIG02.VR  1998/05/09 01:10    BD
    TEST 3
04  DIAGIG03.VR  1998/05/09 01:10    BM
    TEST 4
Enter Menu Number ?_
```



Figure 16–8. Trouble Image Display Menu

- Enter the number of the file to be displayed (i.e. “01”) and press **Return**.  
The image will be displayed on the monitor for analysis.

## Loose cables

If ECG cables are loose:

1. Check the ECG cable lead connection  
OR
2. Check the connection on the back panel.

If peripheral cables are loose:

1. Check the connections on the back of the peripheral  
OR
2. Check the connection on the back panel.





### Display Messages

The LOGIQ™ 500 provides a variety of messages concerning the status of the system's operation.

- System Error Messages are displayed when a hardware problem is detected.
- Operation Error Messages are displayed when the user selects a control that is not appropriate to the current scan mode or function.
- Operation Guide Messages are displayed to inform the operator of a completed task or required action.
- Warning Messages are displayed for patient safety or system failure reasons.
- System Preparation Message is displayed when the system powers up.

These messages are displayed in the single line message area at the bottom of the display. The message display time varies depending on the type of message. Some messages may be accompanied by an audio tone.

For Operation Error Messages, the user can choose to turn on/off the audio tone that accompanies the display in the Set Up/System Parameters page 1.

The User can turn on/off the recording of the error message line in the Set Up/System Parameters page 5 (Maskline Record).

The following tables outline the error messages, their possible cause and cure.

## System Error Message Description

| Category | Problem/<br>Message                                                 | Possible<br>Cause                              | Possible<br>Corrective Action |
|----------|---------------------------------------------------------------------|------------------------------------------------|-------------------------------|
| General  | SYS ERROR: BUS<br>ERROR.<br>PC: #####<br>AD: #####                  | Bus error detected.                            | *                             |
|          | SYS ERROR: FILE<br>NOT FOUND.<br>FILE: #####                        | Requested file not found.                      | *                             |
|          | SYS ERROR: BOARD<br>NON EXISTENT.<br>BOARD: #####                   | Hardware board is non<br>existent.             | *                             |
|          | SYS ERROR:<br>UNEXPECTED CPU<br>EXCEPTION.                          | Unexpected CPU<br>exception is<br>encountered. | *                             |
|          | SYS ERROR: ILLEGAL<br>FILE FORMAT.<br>FILE: #####                   | File format is illegal.                        | *                             |
|          | SYS ERROR: SCSI<br>ERROR.<br>DEVICE: ###<br>CODE: #####             | SCSI device error is<br>encountered.           | *                             |
|          | SYS ERROR: LOCAL<br>CPU ERROR.<br>BOARD: #####                      | Local CPU error is<br>encountered.             | *                             |
|          | SYS ERROR:<br>INTERRUPT TIME OUT.<br><br>BOARD: #####<br>VECTOR: ## | Time out of waiting<br>interrupt.              | *                             |

Table 16–3. System Error Message Description

- \* If System Error Message occurs, record the message and contact a Service Representative. Be aware that the system may be operating at reduced capabilities, some functions inoperative or not operating at all.



## Operation Error Message Description

| Category                                    | Problem/<br>Message                   | Possible<br>Cause                                                             | Possible<br>Corrective Action              |
|---------------------------------------------|---------------------------------------|-------------------------------------------------------------------------------|--------------------------------------------|
| General                                     | Register is filled.                   | Cannot register any more because registration function has reached its limit. | *                                          |
|                                             | Register is empty.                    | Cannot delete any more because registration function is empty.                | *                                          |
|                                             | Save data? 'y' or 'n'.                | System requires operator reaction to save some data.                          | Y to save data,<br>N to not save data.     |
|                                             | Overwrite existing data? 'y' or 'n'.  | There is pre-existing customized data when saving.                            | Y to save changes,<br>N to ignore changes. |
|                                             | Please Select 'y' or 'n'.             | A key other than those specified was pressed.                                 | Enter Y or N.                              |
|                                             | Improper Data Selected.               | The data input is not correct.                                                | Enter correct data.                        |
| Recording<br>Devices<br>Color Video Printer | Check Printer. Paper is jammed.       | Paper in the printer is jammed.                                               | Locate and clear the paper jam.            |
|                                             | Check Printer. Ribbon Cassette error. | Something is wrong with the ribbon cassette or its setting.                   | Remove ribbon and re-install.              |
|                                             | Check Printer. Paper Cassette error.  | The paper cassette setting is wrong.                                          | Remove cassette, set and re-install.       |
|                                             | Check Printer, Paper Supply empty.    | The printer is out of paper.                                                  | Add paper.                                 |
|                                             | Shut the top hatch of printer.        | The top hatch of the printer is open.                                         | Open/shut top hatch.                       |
|                                             | Shut the bottom hatch of printer.     | The bottom hatch of the printer is open.                                      | Open/shut bottom hatch.                    |
|                                             | Check Printer, Paper is incorrect.    | The placing of paper is incorrect.                                            | Remove cassette, adjust and reinstall.     |
|                                             | Check Printer.                        | There is no response from the printer.                                        | Turn power off and on again.               |

Table 16–4. Operation Error Message Description

- \* If System Error Message occurs, record the message and contact a Service Representative. Be aware that the system may be operating at reduced capabilities, some functions inoperative or not operating at all.

## Operation Error Message Description (cont'd)

| Category                                        | Problem/Message                              | Possible Cause                                                                            | Possible Corrective Action                                                                                 |
|-------------------------------------------------|----------------------------------------------|-------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Recording Devices<br/>Multi Image Camera</b> | Check M.I.C., No Exposures Left.             | The maximum number of prints have been made.                                              | Reverse cassette or install a new cassette.                                                                |
|                                                 | Check M.I.C., No Video Input.                | There is no video signal being input to the M.I.C.                                        | Check all video connections.                                                                               |
|                                                 | Check M.I.C.                                 | There is no response from the M.I.C.                                                      | Turn power off, turn power on.                                                                             |
| <b>Recording Devices<br/>VCR</b>                | Check VCR.                                   | There is no response from the VCR.                                                        | Check all VCR controls.                                                                                    |
|                                                 | Check VCR. No Cassette.                      | The VCR has no tape.                                                                      | Insert tape.                                                                                               |
|                                                 | This function is not available.              | The function selected is not available in some situations.                                | Perform from VCR front panel.                                                                              |
|                                                 | Protected. Check Tape.                       | This tape has been write protected.                                                       | Use non protected tape.                                                                                    |
| <b>Line Printer</b>                             | Check Line Printer. The port isn't set.      | The line printer has not been assigned to Port MIC, A or B.                               | Assign the line printer to Port MIC, A or B.                                                               |
|                                                 | Check Line Printer. No Power.                | The line printer is powered off or the serial cable is not connected to port MIC, A or B. | Power on the line printer or connect the serial cable to port MIC, A or B.                                 |
|                                                 | Check Line Printer. Printing Error.          | An error has occurred in sending data to the line printer or the printer is out of paper. | Resend the report data or add paper to the printer.                                                        |
|                                                 | Check Line Printer. Printing isn't finished. | Report Print Buffer is full.                                                              | Ctrl + R cannot be used until space is available. The system can store 10 reports for background printing. |

Table 16–4. Operation Error Message Description (cont'd)



## Operation Error Message Description (cont'd)

| Category                         | Problem/<br>Message                    | Possible<br>Cause                                                                                   | Possible<br>Corrective Action       |
|----------------------------------|----------------------------------------|-----------------------------------------------------------------------------------------------------|-------------------------------------|
| Measurement &<br>Calculation     | Measurement(s)<br>Required.            | The system requires the<br>operator to input<br>measurement(s).                                     | Input necessary<br>measurements.    |
|                                  | Press 'FREEZE' to<br>complete.         | The system requires the<br>operator to press<br><b>FREEZE</b> before<br>completing<br>measurements. | Freeze before all<br>measurements.  |
|                                  | No correct data left.                  | The system no longer<br>has correct data.                                                           | Enter new<br>measurements.          |
|                                  | Manual Calibration<br>Required.        | The system cannot<br>calibrate automatically<br>because there is no<br>barcode data on tape.        | Reply Yes to manual<br>calibration. |
|                                  | Barcode sum check<br>error.            | There is something<br>wrong with the barcode<br>data.                                               | Retry, use manual<br>calibration.   |
| Scanning<br>Function             | Transducer Not<br>Recognized.          | An invalid probe has<br>been connected.                                                             | Remove invalid probe.               |
|                                  | Transducer Calibration in<br>Process.  | Probe change is taking<br>place.                                                                    | Wait for process to<br>complete.    |
| Option Setting                   | Option Not Available.                  | The option to perform<br>the requested function is<br>not available.                                | Install desired option.             |
| New Patient,<br>ID/Name Function | Invalid, Input a proper<br>data.       | Improper data has been<br>input.                                                                    | Input valid data.                   |
|                                  | Press 'RETURN' to exit.                | An invalid key has been<br>pressed.                                                                 | Press <b>Return</b> .               |
|                                  | Data range over, Input<br>proper data. | The input data is too<br>large.                                                                     | Numbers too large.                  |
|                                  | Please answer 'y' or 'n'.              | A key other than those<br>specified was pressed.                                                    | Y for yes, N for no.                |
|                                  | Improper Data Selected.                | Improper data has been<br>input.                                                                    | Input proper data.                  |

Table 16–4. Operation Error Message Description (cont'd)

## Operation Error Message Description (cont'd)

| Category         | Problem/Message                             | Possible Cause                                                                                                                   | Possible Corrective Action                                                                          |
|------------------|---------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|
| Storage Function | Mismatched System ID. Insert correct media. | The disk is not a registered disk. Images cannot be deleted.                                                                     | Take disk to system with proper ID.                                                                 |
|                  | The MO is full, change MO.                  | There is not enough space to store additional images on the disk.                                                                | Use a different disk.                                                                               |
|                  | Media is locked.                            | The disk has been write protected.                                                                                               | Disable the write protection or use another disk.                                                   |
|                  | Media is either unmounted or unformatted.   | The disk is either not formatted or inserted in the drive.                                                                       | Either format the disk or insert it into the drive.                                                 |
|                  | Too much data. Store failed.                | The hard drive capacity of the hard drive to store additional MOD media search data is full.                                     | Backup MOD media search data from the hard drive to a separate MOD before attempting image storage. |
|                  | This function is not available.             | During the execution of the background storage process, this message appears if Ctrl + W, Record 1 or Record 2 keys are pressed. |                                                                                                     |
| Miscellaneous    | No function is defined.                     | No function is defined for the User Define key that was pressed.                                                                 | Use setup menu to define.                                                                           |
|                  | (Y/N)                                       | Operator Decision Requested                                                                                                      | Type in Y for Yes or N for No                                                                       |

Table 16–4. Operation Error Message Description (cont'd)



## Operation Guide Message Description

| Category         | Problem/<br>Message                                  | Possible<br>Cause                                                      | Possible<br>Corrective Action               |
|------------------|------------------------------------------------------|------------------------------------------------------------------------|---------------------------------------------|
| Scan Guide       | Invalid CFM image was cleared.                       | The PRF is changed, and invalid CFM image is cleared.                  | NA                                          |
| Miscellaneous    | In progress. Please wait.                            | The system needs several seconds to do something.                      | NA                                          |
| VCR Search Guide | Searching.                                           | The VCR is searching the tape.                                         | NA                                          |
|                  | Search was completed.                                | The VCR search is completed.                                           | NA                                          |
|                  | Requested information doesn't exist on the tape.     | The patient image doesn't exist on the tape.                           | Perform a system search.                    |
|                  | Requested information doesn't exist in System.       | The information doesn't exist on the hard drive.                       | Check all patient information.              |
|                  | Can not read Tape ID. Register the new tape? (y/n)   | The tape ID cannot be read.                                            | "Yes" to register new tape, "No" to cancel. |
|                  | Retry to read Tape ID? (y/n)                         | The tape ID cannot be read.                                            | "Yes" to retry, "No" to cancel.             |
|                  | Please play VCR image and freeze at the other frame. | The operator has retried to read the tape ID.                          | Press <b>Play</b> and <b>Freeze</b> .       |
|                  | Save patient information on VCR Tape? (y/n)          | Operator pressed key for saving patient information on the VCR tape.   | "Yes" to continue, "No" to cancel.          |
|                  | Save patient information to System? (y/n)            | Operator pressed key for saving patient information on the hard drive. | "Yes" to save, "No" to cancel.              |
|                  | System ID = #####<br>Tape ID = #####                 | Tape was mounted or the new tape was registered.                       | NA                                          |
|                  | System ID = #####<br>New Tape ID = #####             | The tape was mounted.                                                  | NA                                          |

Table 16–5. Operation Guide Message Description

## Operation Guide Message Description (cont'd)

| Category                | Problem/Message                      | Possible Cause                                                                    | Possible Corrective Action            |
|-------------------------|--------------------------------------|-----------------------------------------------------------------------------------|---------------------------------------|
| Preset & Setup Function | Data is saved.                       | The system completed saving the data that was requested.                          | NA                                    |
|                         | Input value by Rocker Button.        | Operator can use rocker switch to input the value.                                | Use the <b>Ellipse</b> rocker switch. |
|                         | Input character strings.             | Operator needs to input some character strings.                                   | Type in characters.                   |
|                         | Select an initial data.              | Operator needs to select initial data.                                            | Input proper data.                    |
|                         | 'TRACKBALL' and 'SET' are available. | Operator can use the <b>Trackball</b> and <b>Set</b> key to operate the function. | NA                                    |
|                         | Recommended value are copied.        | The system completed copying the data that was requested.                         | NA                                    |

Table 16–5. Operation Guide Message Description (cont'd)

## Warning Message Description

| Category                                                                                                                 | Problem/Message                                  | Possible Cause                                                                       | Possible Corrective Action             |
|--------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------|--------------------------------------------------------------------------------------|----------------------------------------|
|  <b>WARNING</b><br>Temperature Safety | WARNING: Temperature error.                      | System internal temperature is at a maximum, but not caused by fan failure.          | Reduce room temperature.               |
|                                                                                                                          | WARNING: Cooling Fan Failure.                    | System fan has stopped working properly.                                             | Call service.                          |
|                                                                                                                          | WARNING: Cooling Fan Failure. Temperature error. | System internal temperature is at a maximum, caused by the fan not working properly. | Reduce room temperature; call service. |
|                                                                                                                          | WARNING: HV Abnormal.                            | The High Voltage is abnormal.                                                        | Call service.                          |
|                                                                                                                          | WARNING: Probe Temperature error.                | The probe temperature is at a maximum.                                               | Call service.                          |
|                                                                                                                          | WARNING: NOW START THE POWER OFF PROCESS.        | The system shutdown process has begun.                                               | Call service.                          |

Table 16–6. Warning Message Description





# Operator Diagnostics

## Introduction

Refer to Section 4–3 of the LOGIQ™ 500 Service Manual (P9030TA).



The user is able to run diagnostic tests such as:

System Test 1 (reduced version)

Black & White Test Pattern

Color Test Pattern

Graphics Test Pattern

The test patterns can be used to adjust the display monitor or peripheral displays. The system test provides a reduced (pass/fail) version of the diagnostic used to check the LOGIQ™ 500 software.

## Probe Selection

Connect and select the C364 probe for diagnostic test purposes.

## Accessing Diagnostics

To access the available diagnostics, select the Set Up Top Menu.

| B              | Preset         | Set Up      | ECG         |
|----------------|----------------|-------------|-------------|
| Custom Display | Preset Program | Save Values | User Define |
| System Paramtr |                | Utility     | Diag.       |

Figure 16–9. Set Up Sub-Menu

Select **DIAG** from the Set Up Sub-Menu. The diagnostic test menu appears on the screen.

```

HOSPITAL NAME                               98/05/09 02:26:14
Diagnosis Menu(1/1)                         S317

No Time      Menu Name
01 30M       System Test(Reduced)
02 0M        Test Pattern Black & White
03 0M        Test Pattern Color
04 0M        Test Pattern Graphics

Enter Menu Number ?
  
```

Figure 16–10. Diagnostic Test-Menu

Type in the desired test number and press **Return**. The test begins running or the test pattern is displayed.

EXE TIME is the test execution time. (The time needed to run the test.)



## System Test 1 (reduced)

Select 01 from the diagnostic test menu.

System Test 1 is a collection of 24 individual tests that are run in sequence. These are the full or reduced versions of the hardware tests.

Three possible messages will be displayed to the right of the test name.

RUNNING—Test is in progress.

FAILED—Test was not completed successfully for some reason.

PASSED—Test was run and completed successfully.

Figure 16–11 shows the screen display during system Test 1 operation.

|                                  |                     |
|----------------------------------|---------------------|
| GE MEDICAL SYSTEMS               | 03/21/94 01: 31: 50 |
| System Test                      |                     |
| 03 Beam Former Test 1            | Running             |
| 04 Beam Former Test 2            |                     |
| 05 Log Test                      |                     |
| 06 Frequency Characteristic Test |                     |
| 08 Doppler Test 1                |                     |
| 09 Doppler Test 2                |                     |
| 11 Color Flow Test 1             |                     |
| 12 Color Flow Test 2             |                     |
| 13 TLMP Test (Reduced)           |                     |
| 14 DDSC Test (Reduced)           |                     |
| 15 CINE Test (Reduced)           |                     |
| 16 VIDO Test (Reduced)           |                     |
| 17 VIDO & VPBM Test (Reduced)    |                     |
| 18 EXBO Test (Reduced)           |                     |
| 19 RCTL Test (Reduced)           |                     |
| 20 PIOP Test (Reduced)           |                     |
| 21 Keyboard RAM Test (Reduced)   |                     |
| 27 Power Supply & Thermal Test   |                     |
| 29 B Mode Noise Floor Test       |                     |
| 30 D Mode Noise Floor Test       |                     |
| 31 CFM Mode Noise Floor Test     |                     |

Figure 16–11. System Test 1 Screen Display

Should any test fail, consider that the system may be operating at reduced capabilities. Record the System Test 1 results and contact a local service representative.

Total run time for System Test 1 is approximately 15 minutes.



### **Test Pattern Black & White**

Select 02 from the Diagnostic Test Menu.

A black/white test pattern is displayed.

This pattern can be used to adjust the LOGIQ™ 500 monitor or peripheral displays.

### **Test Pattern Color**

Select 03 from the Diagnostic Test Menu.

A color video test pattern is displayed.

This pattern can be used to adjust the LOGIQ™ 500 or peripheral output devices.

### **Test Pattern Graphics**

Select 04 from the Diagnostic Test Menu.

A graphic video test pattern is displayed.

This checks the graphics generator of the LOGIQ™ 500 and can be used for monitors or peripheral devices.



# Assistance

## Clinical Questions

For information, call your local Applications, Sales or Service Representative.



**For USA  
Only**

GE Clinical Answer Center at 1-800-682-5327 or  
414-524-5698.

## Service Questions

For service, call your local Service Representative.



**For USA  
Only**

GE CARES at 1-800-437-1171

## Literature

To request the latest GE Accessories catalog or equipment brochures, call your local Applications, Sales or Service Representative.



**For USA  
Only**

Response Center at 1-800-643-6439

## Accessories

To place an order, call your local Applications, Sales or Service Representative.



**For USA  
Only**

GE Access Center at 1-800-472-3666.

## Supplies/Accessories

### CAUTION



**DO NOT** connect any probes or accessories without approval by GE.

The following supplies/accessories have been verified to be compatible with the LOGIQ™ 500 system:

### CAUTION



The LOGIQ™ 500 ultrasound system complies with regulatory requirements of **CE**<sub>0459</sub>, except Gel, Disinfectant, Biopsy Kits, Sheath Sets and B510 probe.

## Peripherals

| Accessory                               | Units                 | Catalog Number |
|-----------------------------------------|-----------------------|----------------|
| S-VHS VCR Tapes                         | 60 Minute S-VHS tape  | E7010GG        |
|                                         | 126 Minute S-VHS tape | E7010GF        |
| Sony VCR<br>Model SVO-9500MD            | Each                  | H4550MA        |
| Sony B&W Printer<br>Model UP-890        | Each                  | E8310KA        |
| Sony Color Printer<br>Model UP-1800     | Each                  | H4550MB        |
| IIE Multi Image Camera<br>Model IIE 460 | Each                  | H4550KF        |
| Sony UP-1800 Paper                      | Color                 | E8310JH        |
|                                         | B/W                   | E8310JJ        |
| Sony UP-890 Paper                       |                       | E8310FA        |

Table 16–7. Peripherals and Accessories

## Console

| Accessory            | Units | Catalog Number |
|----------------------|-------|----------------|
| Foot Switch          | Each  | H40582L        |
| Vinyl Cover          | Each  | H40642L        |
| Extended Cine Memory | Kit   | H40592L        |

Table 16–8. Console Accessories



## Assistance

---

### Probes

| Accessory            | Units | Catalog Number |
|----------------------|-------|----------------|
| B510                 | Each  | H45202BT       |
| C358                 | Each  | H45202CD       |
| C364                 | Each  | H45202CF       |
| C386                 | Each  | H45202CC       |
| C551                 | Each  | H45202CE       |
| C721                 | Each  | H45202MN       |
| E721                 | Each  | H45202MT       |
| ERB7                 | Each  | H40392LC       |
| I739                 | Each  | H45202JG       |
| LA39                 | Each  | H45202LA       |
| 546L                 | Each  | H45202LE       |
| 739L                 | Each  | H45202AG       |
| L764                 | Each  | H45202HP       |
| P509                 | Each  | H45202RT       |
| S220                 | Each  | H45202WG       |
| S222                 | Each  | H45202TC       |
| S316                 | Each  | H45202SC       |
| S317                 | Each  | H45202SD       |
| S611                 | Each  | H45202SF       |
| T739                 | Each  | H45202TG       |
| CWD2                 | Each  | H45202DB       |
| CWD5                 | Each  | H45202DE       |
| PA51 Adapter         | Each  | H45202PA       |
| S5 Adapter           | Each  | H46322AC       |
| Probe Cable Arm      | Each  | H40622L        |
| Left Side Holder     | Each  | H40602L        |
| Special TV/TR Holder | Each  | H40612L        |

Table 16–9. Probes and Accessories

## Gel

| Accessory              | Units                                         | Catalog Numbers |
|------------------------|-----------------------------------------------|-----------------|
| Thermasonic Gel Warmer | Holds three plastic bottles (250 ml or 8 oz). | E8365BH         |
| Aquasonic 100 Scan Gel | 5 liter jug                                   | E8365AF         |
|                        | 250 ml plastic bottles (12/case)              | E8365BA         |
| Scan Ultraound Gel     | 8 oz plastic bottles (12/case)                | E8365BC         |
|                        | 1 gallon plastic jug                          | E8365BD         |
|                        | Four 1-gallon plastic jugs                    | E8365BK         |

Table 16–10. Gel

## Disinfectant

| Accessory                  | Units             | Catalog Number |
|----------------------------|-------------------|----------------|
| Cidex Activated Dialdehyde | 16/1 qt. bottles  | E8386EB        |
|                            | 4/1 gal. bottles  | E8386EC        |
|                            | 2/2.5 gal bottles | E8386ED        |

Table 16–11. Disinfectant

## Civco Biopsy Starter Kits (includes bracket)

| Accessory             | Units       | Catalog Number |
|-----------------------|-------------|----------------|
| C358 Probe Biopsy Kit | Kit         | E8385RK        |
| C364 Probe Biopsy Kit | Kit         | E8385MD        |
| C551 Probe Biopsy Kit | Kit         | E8385MF        |
| S316 Probe Biopsy Kit | Kit         | E8385MK        |
| 739L Probe Biopsy Kit | Kit         | E8385MC        |
| E721 Probe Biopsy Kit | Kit (GEYMS) | H46222AD       |
|                       | Kit (Civco) | H4550BG        |
| ERB7 Probe Biopsy Kit | Kit         |                |
|                       | Kit (Civco) |                |
| C721 Probe Biopsy Kit | Kit         | H40252L        |
| L764 Probe Biopsy Kit | Kit         | E8385LA        |
| T739 Probe Biopsy Kit | Kit         | H45262L        |
| LA39 Probe Biopsy Kit | Kit         | E8385MM        |

Table 16–12. Probe Biopsy Brackets





## Assistance

---

### Multi-Angle Brackets

| Accessory  | Units | Catalog Number |
|------------|-------|----------------|
| S317 Probe | Each  | E8385PD        |
| C364 Probe | Each  | E8385PA        |
| C386 Probe | Each  | E8385PB        |
| 546L Probe | Each  | E8385PC        |
| T739 Probe | Each  | E8385PF        |
| 739L Probe | Each  | E8385PG        |

Table 16–13. Multi-Angle Brackets

### Biopsy Replacement Kits

| Accessory                                  | Units             | Catalog Number |
|--------------------------------------------|-------------------|----------------|
| Civco UltraPro                             | Kit               | E8385LC        |
| Civco UltraPro II                          | Kit               | E8385PE        |
| Civco Cleaning Brush for E721 Biopsy Guide | 5 Brushes Per Set | E8323GE        |

Table 16–14. Biopsy Kits

### Ultrasound Probe and Cord Sheath Sets

| Accessory                                                                         | Units      | Catalog Number |
|-----------------------------------------------------------------------------------|------------|----------------|
| Sterile Ultrasound Probe Sheath Set                                               | 20 Per Set | E8385CA        |
| Sterile Ultrasound Cord Sheath Set                                                | 20 Per Set | E8385CB        |
| Sanitary Rectal/Vaginal Probe Cover                                               | 20 Per Set | E8385CC        |
| Sterile Combination Probe and Cord Cover Set                                      | 12 Per Set | E8385CE        |
| Sterile Ultrasound Probe Sheath Set for Wide (2.5 and 3.5) Aperture Sector Probes | 20 Per Set | E8385CF        |

Table 16–15. Biopsy Sheath Kits

## Physio Accessories

| Accessory          | Units | Catalog Number |
|--------------------|-------|----------------|
| Physio Input Panel | Kit   | H40552L        |
| ECG Cables         | Set   | H40562L        |
| PCG Sensor         | Each  | H40572L        |

Table 16–16. Patient ECG Lead Wires

## Patient Electrodes

| Accessory            | Units    | Catalog Number |
|----------------------|----------|----------------|
| Adult                | Box/300  | E8811EE        |
| Pediatric            | Box/300  | E8811EN        |
| Adult/off center     | Box/500  | E8365DA        |
| Pediatric/off center | Box/1000 | E8365DB        |

Table 16–17. Patient ECG Electrodes

Contact the distributor, GE affiliate or sales representative for approved peripherals.



# Quality Assurance

## Introduction

A good Quality Assurance Evaluation program consists of planned systematic actions that provide the user with adequate confidence that their diagnostic ultrasound system will produce consistently high quality images and quantitative information.

Therefore, it is in the best interests of every ultrasound user to routinely monitor equipment performance.

The frequency of Quality Assurance evaluations should be based on user's specific needs and clinical practice.

Periodic monitoring is essential in order to detect the performance changes that occur through normal aging of system components. Routine equipment evaluations may also reduce the duration of exams, number of repeat exams, and maintenance time required.



Refer to the *System Care and Maintenance* section of this chapter for system and peripheral routine preventive maintenance instructions.



## Typical Tests to Perform

Quality assurance measurements provide results relating to system performance. Typically these are:

- Axial Measurement Accuracy
- Lateral Measurement Accuracy
- Axial and Lateral Resolution
- Penetration
- Functional & Contrast Resolution
- Gray Scale Photography.

With these tests, a performance baseline can be set at installation with the phantom in your department. Future test results can be compared to the baseline in order to maintain a record of system performance trends.

## Frequency of tests



Every 3 Months  
Or  
Every 400 Patients  
Or  
If system  
performance  
in question

Quality assurance tests are used to determine whether a scanner is providing the same level of performance from day to day.

The frequency of testing varies with the amount of system usage and modes to be tested. It is recommended that the user perform quality assurance tests at least every three months or every 400 patient studies. Tests should also be performed when a question about system performance exists.

*A mobile system may require more frequent tests.*

Image quality should also be tested immediately after the following events:

- Service calls
- System upgrades/modifications
- Dropped probe, power surge, etc.



## Quality Assurance

---

### Phantoms

Quality Assurance Evaluations should be done with phantoms and test objects that are applicable to the parameters being evaluated or to the user's clinical practice.


Typical phantoms are composed of material that acoustically mimic human tissue. Pins, anechoic and echogenic targets are physically positioned to provide information for a variety of tests.

Doppler phantoms are currently expensive and complicated to deal with on the user level. If a problem with any Doppler parameters or measurement is suspected, contact a local service representative for evaluation.



*The RMI 403GS phantom is still available. Due to the superior penetration and resolution capabilities of GE LOGIQ systems, the RMI 405GSX is recommended. It is the most current one available to our field service personnel and will provide the targets and extended life necessary for consistent system testing.*

Phantoms (cont'd)



1. Penetration
2. Axial Distance Measurement
3. Functional Resolution
4. Lateral Resolution
5. Lateral Distance Measurement
6. Axial Resolution
7. Contrast Resolution and Gray Scale Photography
8. Gray Scale Plane Targets

Figure 16–12. Phantoms



### Baselines

An absolute necessity for a quality assurance program is establishing baselines for each test or check. Baselines are established after the system has been verified to be working properly at installation or after a repair. If a probe or major assembly is replaced, new baselines should be generated.

Baselines can be made by adjusting system parameters to prescribed levels or to the best possible image. The key factor to remember is reproducibility. The same conditions must be reproduced for each periodic check.

All system parameters not displayed on the monitor should be recorded for the permanent record.

### Periodic Checks

Periodic checks are performed on a regular basis as previously recommended. For the data to be valid, periodic checks should mimic the baseline setup parameters.

The resulting image, when scanning the phantom exactly as before, should be recorded and compared to the baseline. When a matching image is obtained, it can be assumed that the system performance has not degraded from the baseline.

If a significant difference between the baseline and periodic check is noted, double check the system setup and repeat the test. If the difference between the baseline and periodic check persists, contact a local Service Representative.

*Failing to reproduce the control settings as in the baselines will introduce errors in the data and potentially invalidate the results.*



## **Results**

Lack of standardization among test instruments, the wide range of acceptance criteria, and incomplete knowledge regarding the significance of certain performance parameters prohibit the establishment of absolute performance criteria for these tests.

Quality Assurance Evaluation results should be compared to previously-recorded results.

Performance trends can then be detected. Unacceptable performance or diminishing trends should be identified for maintenance or repair before a malfunction or inappropriate diagnosis occurs.

The user should determine the best method for recording and archiving the baseline and periodic checks. In most cases the choice is hard copy.

It is important to maintain good consistent records for inspections that may arise, as well as to detect system performance trends.





### System Setup

The user should tailor the tests to their particular needs. It is certainly not necessary to make all checks with all probes. A representative example, with the probes used most often by the customer, should be adequate in judging system performance trends.

Use a gray scale phantom as the scan object for the tests. Commercial phantoms are supplied with its own operator manual. Be familiar with proper phantom operating procedures prior to use for quality assurance evaluations.

1. Adjust image monitor. Brightness and contrast should be set to the normal viewing of a good gray scale image.
2. Check all recording devices for proper duplication of image monitor. Ensure that what is seen is what is recorded. Check the B/W or Color page printer, VCR, MIC, or Laser camera.
3. Annotate non-displayed image processing controls.
4. Set TGC slide pots to center (detent) position.
5. Place focal zone marker(s) in area of interest for an optimum image.

### Test Procedures

The following are recommended Quality Assurance tests. A brief description of the test, the benefit it provides and steps to accomplish the test are supplied.

The importance of recording scan parameters and consistent record keeping cannot be stressed enough. Reproducibility to monitor system trends is the key to quality assurance evaluations.

Using the system's dual image display format is often very convenient and saves recording media.



## Axial distance measurements

|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Axial measurements are the distance measurements obtained along the sound beam. Refer to Figure 16–12 on 16–47 for details.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Benefit</b>     | The accurate measurement of the size, depth and volume of a structure is a critical factor in determining a proper diagnosis. Most imaging systems use depth markers and/or electronic calipers for this purpose.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Method</b>      | Axial distance should be measured in the near, mid and far fields as well as in zoom. If necessary, different depths or fields of view can be tested.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Procedure</b>   | <p>To measure axial distance:</p> <ol style="list-style-type: none"><li>1. Scan a test phantom with precisely-spaced vertical pin targets. Adjust all scan controls, as necessary, for the best image of the pin targets to typical depths for the probe being used.</li><li>2. Press <b>Freeze</b> to stop image acquisition and perform a standard distance measurement between the pins at different points in the image. Record all images for archiving.</li><li>3. Scan the vertical pins in zoom or at different depth/scale factors.</li><li>4. Press <b>Freeze</b> to stop image acquisition; repeat the distance measurements between pins and record the images for archiving.</li><li>5. Document the measurements for reference and future comparison.</li></ol> |



Contact a Service Engineer if vertical measurements differ by more than 1.5% of the actual distance.



## Lateral distance measurements

|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Lateral measurements are distance measurements obtained perpendicular to the axis of the sound beam. Refer to Figure 16–12 on 16–47 for details.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Benefit</b>     | The purpose is the same as vertical measurements. Precisely-spaced horizontal pin targets are scanned and results compared to the known distance in the phantom.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Method</b>      | Lateral distance should be measured in the near, mid and far fields as well as in zoom. If necessary, different depths of fields of view can be tested.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Procedure</b>   | <p>To measure lateral distance:</p> <ol style="list-style-type: none"><li>1. Scan a test phantom with precisely-spaced horizontal pin targets. Adjust all scan controls, as necessary, for the best image of the pin targets from side to side.</li><li>2. Press <b>Freeze</b> to stop image acquisition and perform a standard distance measurement between the pins at different points in the image. Record all images for archiving.</li><li>3. Scan the horizontal pins in zoom or at different depth/scale factors.</li><li>4. Press <b>Freeze</b> to stop image acquisition; repeat the distance measurements between pins and record the images for archiving.</li><li>5. Document the measurements for reference and future comparison.</li></ol> |



Contact a Service Engineer if horizontal measurements differ by more than 3mm or 3% of that depth, whichever is greater.



## Axial resolution

### Description

Axial resolution is the minimum reflector separation between two closely-spaced objects to produce discrete reflections along the axis of the sound beam. It can also be monitored by checking the vertical size of known pin targets. Refer to Figure 16–12 on 16–47 for details.

Axial resolution is affected by the transmitting section of the system and the probe.

### Benefit

In clinical imaging, poor axial resolution displays small structures lying close together as a single dot. This may lead to improper interpretation of the ultrasound image.

### Procedure

To measure Axial resolution:

1. Scan a test phantom with precisely-spaced vertical pin targets.
2. Adjust all scan controls, as necessary, for the best image of the pin targets to typical depths for the probe being used.
3. Press **Freeze** to stop image acquisition.
4. Perform a standard distance measurement of the pin vertical thickness at different points in the image. Record all images for archiving.
5. Scan the vertical pins in zoom or at different depth/scale factors.
6. Press **Freeze** to stop image acquisition; repeat the vertical thickness measurements of the pins and record the images for archiving.
7. Document the measurements for reference and future comparison.



Axial resolution should remain stable over time. Contact a Service Engineer if any changes are observed.



## Lateral resolution

### Description

Lateral resolution is the minimum reflector separation between two closely spaced objects to produce discrete reflections perpendicular to the axis of the sound beam. It can also be monitored by checking the horizontal size of known pin targets. Refer to Figure 16–12 on 16–47 for details.

Lateral resolution is dependent upon the beam width produced by the probe. The narrower the beam, the better the lateral resolution.

The beam width is affected by the frequency, degree of focusing, and distance of the object from the face of the probe.

### Benefit

Clinically, poor lateral resolution will display small structures lying close together as a single dot. This may lead to improper interpretation of the ultrasound image.

### Procedure

To measure lateral resolution:

1. Scan a test phantom with precisely-spaced horizontal pin targets.
2. Adjust all scan controls, as necessary, for the best image of the pin targets from side to side.
3. Press **Freeze** to stop image acquisition and perform a standard distance measurement of the horizontal thickness of a pin at different points in the image. Record all images for archiving.
4. Scan the horizontal pins in zoom or at different depth/scale factors.
5. Press **Freeze** to stop image acquisition; repeat the horizontal thickness measurements of the pins and record the images for archiving.
6. Document the measurements for reference and future comparison.



Pin width should remain relatively constant over time ( $\pm 1\text{mm}$ ). Dramatic changes in pin width may indicate beamforming problems. Contact a Service Engineer if beam width changes consistently over 2 to 3 periodic tests.



## Penetration

|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | <p>Penetration is the ability of an imaging system to detect and display weak echoes from small objects at large depths. Refer to Figure 16–12 on 16–47 for details.</p> <p>Penetration can be affected by the system's:</p> <ul style="list-style-type: none"><li>• Transmitter/receiver</li><li>• Degree of probe focusing</li><li>• Attenuation of the medium</li><li>• Depth and shape of reflecting object</li><li>• Electromagnetic interference from local surroundings.</li></ul>                                                                                                                                                        |
| <b>Benefit</b>     | <p>Weak reflecting echoes are commonly produced from the internal structure of organs. Definition of this tissue texture is important in the interpretation of the ultrasound findings.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Method</b>      | <p>Scan a phantom to see how echoes begin to fade as depth is increased. The maximum depth of penetration is the point at which homogeneous material in the phantom begins to lose brightness.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Procedure</b>   | <p>To measure penetration:</p> <ol style="list-style-type: none"><li>1. Set the front panel TGC slide pots to their center (detent) position.</li><li>2. Gain and acoustic output can be adjusted, as necessary, since these values are displayed on the monitor.</li><li>3. Scan a test phantom along the vertical pin targets to typical depths for the probe being used.</li><li>4. Perform a standard distance measurement from the top of the image displayed to the point at which homogeneous material in the phantom begins to lose brightness.</li><li>5. Document the depth measurement for reference and future comparison.</li></ol> |



Contact a Service Engineer if the depth of penetration shifts more than one centimeter (1cm) when using the same probe and same system settings.



### Functional resolution

**Description** Functional resolution is an imaging system's ability to detect and display the size, shape, and depth of an anechoic structure, as opposed to a pin target. Refer to Figure 16–12 on 16–47 for details.

The very best possible image is somewhat less important than reproducibility and stability over time. Routine tests at the same settings should produce the same results.

**Benefit** The data obtained will give a relative indication of the smallest structure the system is capable of resolving at a given depth.

**Procedure** To measure functional resolution:

1. Set the front panel TGC slide pots to their center (detent) position.
2. Gain and acoustic output can be adjusted as necessary, since these values are displayed on the monitor.
3. Scan a test phantom with a vertical row of anechoic cyst targets to typical depths for the probe being used.
4. Evaluate the cysts at various depths for a good (round) shape, well-defined borders and no fill in. Remember, TGC slide pots are centered and should remain fixed. This may **NOT** provide optimal cystic clearing.
5. Document all results for future reference and comparison.



Contact a Service Engineer if a greatly distorted image is obtained.



## Contrast resolution

**Description** Contrast resolution is the ability of an imaging system to detect and display the shape and echogenic characteristics of a structure. Refer to Figure 16–12 on 16–47 for details.

Specific values measured are less important than stability over time. Routine tests at the same settings should produce the same results.

**Benefit** A correct diagnosis is dependent upon an imaging system's ability to differentiate between a cystic or solid structure versus echo patterns from normal surrounding tissue.

**Method** A phantom with echogenic targets of different sizes and depths should be used.

**Procedure** To measure contrast resolution:

1. Set the front panel TGC slide pots to their center (detent) position. Set dynamic range to 54 db.
2. Gain and acoustic output can be adjusted, as necessary, since these values are displayed on the monitor.
3. Scan a test phantom with echogenic targets at the depths available.
4. Evaluate the echogenic targets for contrast between each other and between the surrounding phantom material. Remember, TGC slide pots are centered and should remain fixed. This may **NOT** provide an optimal scan image.
5. Document all results for future reference and comparison.



Contact a Service Engineer if the echogenic characteristics or shapes of the targets appear distorted.





## Gray Scale photography

|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Poor photography will cause loss of low level echoes and the lack of contrast between large amplitude echoes. Refer to Figure 16–12 on 16–47 for details.                                                                                                                                                                                                                                                                                           |
| <b>Benefit</b>     | When photographic controls and film processors are properly adjusted, weak echoes, as well as strong echoes, are accurately recorded on film.                                                                                                                                                                                                                                                                                                       |
| <b>Procedure</b>   | <ol style="list-style-type: none"><li>1. Adjust the camera according to the manufacturer's instructions until the hard copy and video display are equal.</li><li>2. Scan the phantom and it's echogenic contrast targets.</li><li>3. Make a hard copy photograph of the display and compare it to the image on the video monitor for contrast and weak echo display.</li><li>4. Document all results for future reference and comparison.</li></ol> |



Contact a Service Engineer if camera cannot duplicate what is on the image monitor.

*NOTE: Optimization of brightness/contrast controls on the display monitor is imperative in order to make sure that the hardcopy and monitor look alike.*

*The display monitor is adjusted first. The hardcopy camera or printer is adjusted to match the display monitor.*



## Setting up a Record Keeping System

### Preparation

The following is needed:

- Quality Assurance binder
- Magnetic photo pages for hard copy images, or Floppy disk holder for archived images
- Quality Assurance Checklists.

Display the following information while testing quality assurance:

- Acoustic Output
- Gain
- Depth
- Probe
- Dynamic Range
- Set up new patient to be the name of the test.

Annotate the following:

- Any control where its value is **NOT** displayed.
- Significant phantom information.

### Record Keeping

Complete the following:

1. Fill out the Ultrasound Quality Assurance Checklist for each probe, as scheduled.
2. Make a hard copy or archive the image.
3. Compare images to baseline images and acceptable values.
4. Evaluate trends over previous test periods.
5. File hard copy/archive floppy and checklist in Quality Assurance binder.



## Ultrasound Quality Assurance Checklist

Performed By \_\_\_\_\_ Date \_\_\_\_\_

System \_\_\_\_\_ Serial Number \_\_\_\_\_

Probe Type \_\_\_\_\_ Probe Model \_\_\_\_\_ Serial Number \_\_\_\_\_

Phantom Model \_\_\_\_\_ Serial Number \_\_\_\_\_ Room Temperature \_\_\_\_\_

Acoustic Output \_\_\_\_\_ Gain \_\_\_\_\_ Focal Zone \_\_\_\_\_

Gray Map \_\_\_\_\_ TGC \_\_\_\_\_ Depth \_\_\_\_\_

Monitor Settings \_\_\_\_\_

Peripheral Settings \_\_\_\_\_

Other Image Processing Control Settings \_\_\_\_\_

| Test                            | Base-line Value Range | Tested Value | Image Hard Copy/ Archived | Acceptable? Yes/No | Service Called (Date) | Date Resolved |
|---------------------------------|-----------------------|--------------|---------------------------|--------------------|-----------------------|---------------|
| Vertical Measurement Accuracy   |                       |              |                           |                    |                       |               |
| Horizontal Measurement Accuracy |                       |              |                           |                    |                       |               |
| Axial Resolution                |                       |              |                           |                    |                       |               |
| Lateral Resolution              |                       |              |                           |                    |                       |               |
| Penetration                     |                       |              |                           |                    |                       |               |
| Functional Resolution           |                       |              |                           |                    |                       |               |
| Contrast Resolution             |                       |              |                           |                    |                       |               |
| Gray Scale Photography          |                       |              |                           |                    |                       |               |



# Bioeffects

## Concerns Surrounding the Use of Diagnostic Ultrasound

During a diagnostic ultrasound examination, high frequency sound penetrates and interacts with tissue in and around the area of anatomy to be imaged. Only a small portion of this sound energy is reflected back to the probe for use in constructing the image while the remainder is dissipated within the tissue. The interaction of sound energy with tissue at sufficiently high levels can produce biological effects (aka bioeffects) of either a mechanical or thermal nature. Although the generation of bioeffect is intentional with therapeutic ultrasound, it is generally undesired in diagnostic applications and may be harmful in some conditions.



*NOTE: The American Institute of Ultrasound in Medicine has published a document entitled "Medical Ultrasound Safety". This three part document covers Bioeffects and Biophysics, Prudent Use and Implementing ALARA.*

*Ultrasound users should order this document from the AIUM to become more familiar with Ultrasound safety.*

In the USA, contact the AIUM by telephone at 1-800-638-5352. To write them concerning their publications, use the following address:

AIUM  
14750 Sweitzer Lane  
Suite 100  
Laurel, MD, USA 20707-5906



## **Thermal Bioeffect**

As with most forms of energy, ultrasound is attenuated as it passes through tissue and is converted to heat, which, if produced at sufficient rates, will increase tissue temperature to a point where tissue damage may result. Major factors contributing to thermal bioeffect can be categorized as tissue characteristics or control parameters:

- Physical tissue characteristics like acoustic impedance, attenuation, absorption, and perfusion determine the rates of heat production and heat transfer. The susceptibility of some tissues to injury from heat, such as developing fetal tissue, further complicate the concern for long-term effects.
- The time-average density of available ultrasound energy is mainly determined by acoustic parameters like output frequency, pulse amplitude, pulse duration, duty cycle, beam shape, and beam motion. These parameters are controlled by the operator through equipment selections such as probe type, operating mode, focal depth, sample volume location, and output control settings. The operator also has significant influence by controlling probe motion and dwell time. These “Control Parameters” form the means through which the operator can minimize thermal bioeffect.



### Mechanical Bioeffect

In a similar manner, the interaction of ultrasound energy with tissue can produce a number of non-thermal or mechanical effects. The most significant is cavitation which results from the action of the oscillating ultrasound pressure on tiny gas bubbles within the tissues. Cavitation has caused mechanical damage on a cellular level such as microscopic tears and hemorrhage in laboratory tests with small animals. The major contributing factors can again be categorized as either tissue characteristics or control parameters:

- The physical characteristics of tissue such as the presence and size of microscopic gas bubbles and the sensitivity of the tissue to the effects of cavitation will influence the potential for and magnitude of cavitation.
- Acoustic field parameters like output frequency, peak pulse amplitudes, and perhaps pulse length are the primary parameters affecting the onset of cavitation. These are controllable by the operator through appropriate equipment selections.

Although it is generally accepted that no harmful biological effects have been demonstrated at the frequency, intensity, and exposure times used in diagnostic examinations, research into the potential for harmful effects continues. The operator is encouraged to survey the literature for future developments on bioeffects and to become familiar with the references at the end of this section.



## **Operator Awareness and Actions to Minimize Bioeffect**

The operator must be aware of the particular conditions that exist during the examination to recognize the potential for bioeffect and then take appropriate action to reduce the risk. The recognition of potential harm comes from an understanding of tissue characteristics and a real-time knowledge of acoustic output. Taking appropriate action requires familiarity with equipment operation and examination skills like implementing alternative techniques for obtaining the same diagnostic information.

### **Tissue characteristics**

Tissue characteristics vary considerably throughout the body. They influence the acoustic field and determine the heating/cooling rates and cavitation potential. Ultrasound energy dissipates as it passes through the tissue causing the deeper tissue to encounter much lower levels. Some tissues like bone readily convert ultrasound energy to heat, while others like blood and amniotic fluid pass the energy on to adjacent tissue relatively unattenuated.

A particular situation that represents a tissue combination requiring extra precaution is a third trimester transabdominal fetal examination where there is a very thin abdominal wall and a long fluid path. The relative lack of attenuating tissue along the acoustic path will significantly increase the available energy in the fetal tissue. Additionally, fetal tissues are more susceptible to long term injury due to nature of developing tissue. Focusing the ultrasound beam on or near fetal bone further increases the risk.

Other than fetal tissue, there is increased susceptibility for heating in any tissue that cannot easily conduct or distribute heat due to low blood perfusion. As the examination progresses, the operator must be aware of changing tissue conditions.



### Acoustic output

Awareness of the acoustic output level can be a difficult task for the operator, especially when the objective is to obtain a quality image. Older ultrasound equipment had limited means, if any, for indicating the acoustic output level. In most cases, the operator had to be familiar with the output intensities as described in the operator manual. To improve operator awareness of acoustic output, this system incorporates an output display that directly indicates the potential for mechanical and thermal bioeffects as equipment controls are adjusted.

The output display consists of four numeric index values that indicate the potential for producing bioeffects (three indices are for heating effect and one for cavitation). As the user changes equipment settings that alter the acoustic output, the output display indices are immediately updated to reflect the change in potential for producing bioeffect. The indices are based on mathematical models and each is normalized so that the potential for bioeffect becomes more significant as the indices reach a value 1.0 or larger.

A mechanical index (MI) provides an indication of the potential for the possible onset of transient cavitation within tissue while the three thermal indices provide an indication of the potential for heat generation within tissue. The different thermal indices may be used depending on the type of tissue being examined:

- Soft Tissue Thermal Index (TIS) is used as an indicator of the potential to generate heat within soft tissues. This is the most used thermal indicator.
- Bone Thermal Index (TIB) is used as an indicator of the potential to generate heat at the beam focus when focusing on or near bone that is adjacent to very sensitive tissue. This index is intended as a thermal indicator for second and third trimester fetal examination or transfontanelle neonatal cephalic exams.
- Cranial Bone Thermal Index (TIC) is used as an indicator of the potential to generate heat in the near-field when the beam passes through bone at the surface as with adult or pediatric cranial applications.





## Acoustic output (cont'd)

Tissue heating is more of a concern when the acoustic beam is stationary, so the thermal index is likely to increase when Doppler or M-Modes are selected. The influence of specific operator controls on acoustic output is described along with the functional purpose of the control throughout the user manual and a summary is provided in *Safety*.

The operator now has easy access to the status of acoustic output and, when combined with the knowledge of tissue characteristics and beam location, the risk of potential bioeffect can be readily assessed. This display conforms to the AIUM/NEMA Output Display Standard [1] for ultrasound imaging equipment.

## Operator intervention

When conditions indicate a potential for harmful bioeffect, the operator should take action promptly to reduce the risk by changing equipment settings or altering procedural techniques:

- Output display index values much greater than 1.0 represent an increased risk for tissues in particular beam locations. The potential for heating will normally only occur near the surface or at the focus, while the potential for cavitation is reduced away from the focus. Selecting non-scan operating modes such as PW or CW Doppler and M-Mode will significantly increase the thermal index because the beam is stationary.
- Optimize gain and other image enhancement features before increasing the acoustic output control or other equipment controls that significantly affect the output level. Become thoroughly familiar with all controls that affect output and observe the output display for results. Controls affecting output are described throughout the user manual.



### Operator intervention (cont'd)

- Develop and practice skills to localize anatomy and optimize image quality rapidly, then freeze the image as soon as the necessary diagnostic information is obtained. It takes time for tissue temperatures to increase, so reducing exposure time can significantly reduce the potential for injury.
- Avoid susceptible tissues, if possible, by changing probe position, entrance angles or probe type. Higher frequency probes will not penetrate as deep while linear probes have lower near-field energy density. Avoid focusing on bone or poorly perfused tissue. Do not allow the acoustic beam to penetrate or focus on or near the eye.

Although choices like probe selection, mode of operation and other control adjustments have a significant affect on output levels, the ability to change these selections is often restricted by the type of examination or clinical objectives. Therefore, some examinations may require relatively high output levels to achieve success.

The decision to raise acoustic output to potentially harmful levels must include an assessment of the risk/benefit potential. Such decisions are routine with imaging modalities incorporating ionizing radiation such as Nuclear Medicine, X-ray and CT. The principle of ALARA is widely used in these modalities for minimizing the exposure risk and is now a recommended practice with high-level diagnostic ultrasound.

### CAUTION



During each ultrasound examination, the clinical user is expected to weigh the medical benefit of the diagnostic information obtained against the risk of harmful effects. Once an optimal image is achieved the need for increasing acoustic output or prolonging the exposure can not be justified. It is important, therefore, for the user to be familiar with system controls that affect image quality as well as acoustic output. Complete descriptions of image optimization and acoustic output controls are provided in the user instructions.

## Implementing ALARA Methods

The primary objective for any ultrasound examination is to obtain diagnostic information of sufficient quality to benefit the patient. Image quality can usually be improved by increasing the acoustic output or taking more time to refine the image. These same actions, however, will also increase the risk of harmful bioeffects when imaging sensitive tissues or when high output levels are used. The operator is therefore encouraged to use the lowest acoustic output setting necessary to produce clinically acceptable data.

The principle of ALARA, which stands for **As Low As Reasonably Achievable**, is to keep the radiation exposure at the minimum level necessary to obtain the diagnostic information. This principle is widely practiced in medical x-ray protection where exposure at any level is potentially harmful. Historically, ALARA was initiated as a cautious approach for dealing with uncertain hazards but has since become the principle method for reducing the risk of injury from hazards that do not have safe minimum threshold.

While no minimum thresholds for harmful bioeffects have been established with the use of diagnostic ultrasound, the principle of ALARA can be readily implemented on equipment incorporating an output display. As the operator adjusts the equipment to optimize the image quality, the display interactively updates to indicate the effect on output.

Controls that have no noticeable impact on image quality should be set to minimize the output while controls that improve the image quality and also increase acoustic output should be set no higher than needed to achieve a diagnostic quality image. If the output display indicates values much greater than 1.0, the operator should reduce the exposure time and freeze the image as soon as possible.

At very low levels ( $< 0.4$ ), the display is inactive and the potential for harmful bioeffect is negligible. More detailed information concerning the use of ALARA in medical practice can be found in NCRP Report No. 107<sup>2</sup>.



## Clinical instructions for fetal use

Following are illustrative examples of clinical instructions for fetal Doppler use which were prepared by Harold Schulman, M.D. and John Hobbins, M.D.<sup>3</sup>

### Umbilical artery

To obtain a signal:

For continuous wave Doppler instrument, place the pencil probe with appropriate conducting gel on the maternal abdomen. The volume should be at a comfortable level and the power output and gain at a mid-setting. Slowly move the pencil probe by changing the angle or location on the maternal abdomen until the characteristic sound of the fetal umbilical flow can be heard. This is a swishing sound with a rate usually between 120 and 160 BPM. It should be lacking in discernible clicks or sounds similar to valve movement on cardiac auscultation. When the image is obtained, adjust the controls in order to optimize the image by reducing the power output to the lowest setting at which a good quality signal is obtained. The gain setting should also be reduced to make the signal appear crisp with minimal background noise. If diastolic flow signals are not obtained, the angle of the beam incident to the umbilical cord may be too high and another area of cord should be examined to determine whether diastolic flow is present. When a good quality signal is obtained all the way across the display screen, the image may be frozen.

Prior to taking measurements, the signals should be examined to ensure that the variations in waveform size that may be caused by fetal breathing movements are not present. Fetal breathing movements will invalidate any measurements due to the variations they cause in umbilical flow.



### **Umbilical artery (cont'd)**

If a duplex Doppler system is being used, then the umbilical cord should be visualized and attempted to be seen in as much of its length as is feasible, considering its usual coiling. The sample volume of the pulsed Doppler or the sample line if a duplex continuous wave instrument is being used should be placed over one of the two smaller vessels in the umbilical cord. The sample volume size ("gate") should be adjusted to encompass the entire vessel. If this vessel is being seen only in cross section, and diastolic flow is not seen, this may be an artifact caused by the angle between the beam and the vessel and another sampling site should be examined. When this cross sectional view is obtained, rotating the transducer 90 degrees may bring a greater length of the cord into view and permit examination of an area where the artery can be seen at a more advantageous angle. Indices of pulsatility, if calculated, should be determined for each of several heartbeats and averaged.

### **Uterine artery**

For continuous wave Doppler systems, place the pencil probe in the lower lateral portion of the mother's abdomen, generally just above the groin and directed toward the cervix. Adjust the angle or position until the characteristic waveform of the uterine/arcuate artery is obtained. Without visualization, it is important to search carefully for the waveform with the greatest amount of diastolic flow as proximal vessels in the uterine system (internal iliac and hypogastric arteries) may have diminished diastolic flow compared to the systolic flow in a way that might appear abnormal if obtained from the uterine artery itself. The power and gain settings should be adjusted to the lowest levels at which adequate signals are obtained with a minimal amount of noise relative to the signal. The image may be frozen when the desired number of waveforms are present on the screen.



## Bioeffects

---

### Uterine artery (cont'd)

For pulsed Doppler systems, place the transducer in the lower outer portion of the maternal abdomen. The orientation should be towards the parametria with some parts of the amniotic sac visualized. The parametria area can then be examined by maneuvering the sample volume until the characteristic waveforms are obtained. Examination should be confined to the area that is composed of the myometrium and may be most fruitful close to the placenta. Indices of pulsatility, if calculated, should be determined for each of several heartbeats and averaged.

### Fetal heart

Fetal cardiac Doppler studies can be performed only with duplex ultrasound systems. The fetal heart should be carefully examined in all standard planes to determine whether or not the anatomy is normal. These planes should include: four-chamber view, long axis left ventricles, short axis of the ventricles, short axis of the great vessels, aortic arch, and pulmonary artery-ductus views. The area for sampling should be visualized as clearly as possible with the anticipated flow direction at as low an angle as possible to the ultrasound beam. The sample volume should be placed in this area with an appropriate size selected to sample as desired. Power output control should be at the lowest setting compatible with obtaining an adequate image and the gain setting should be maintained to keep an adequate signal without excessive noise. Excluding the actual valves from the sample volume will help to minimize unnecessary noise from the signals obtained.

## Efficacy of Fetal Doppler

The following are clinical obstetrical conditions where there are experimental data that demonstrate the efficacy of Doppler.  
Provided by John C. Hobbins, M.D. and Peter Burns, Ph.D.<sup>4</sup>

### IUGR

Many studies have shown a good correlation between abnormal waveform (or decreased blood flow) and increased perinatal mortality<sup>5</sup>, fetal hypoxia<sup>6,7,8</sup>, and neonatal events such as necrotizing enterocolitis (NEC), and intraventricular hemorrhage (IVH).<sup>8</sup> It must be pointed out that the Doppler studies initially designed to identify altered fetal growth had a predictably low sensitivity<sup>9</sup> because the category of small for dates frequently includes genetically small but completely normal fetuses and neonates. Thus far, all studies have shown that the rare Doppler pattern of reverse diastolic flow in the umbilical artery has been very highly correlated with adverse outcome<sup>10,11</sup> and may warrant immediate intervention. Conversely, a normal waveform in the umbilical artery has been rarely associated with stillbirth in a high risk pregnancy.

### Cardiac Anomalies

Doppler has become an integral part of fetal echocardiographic studies. Its use has been directly derived from well-documented pediatric and adult Doppler research. The etiology and seriousness of fetal arrhythmias is relatively easy to determine with pulsed Doppler interrogation of ventricular diastolic filling patterns.<sup>12–15</sup> Benign patterns such as premature atrial extrasystoles can be more easily diagnosed with Doppler than with the older approach of M-Mode echocardiography. Doppler has proven useful as well in the understanding and interpretation of structural heart disease in the fetus. A close correlation has been found between atrioventricular valve regurgitation in the fetus with structural heart disease and the appearance of non-immune hydrops<sup>16</sup>. Improved diagnostic and prognostic ability has permitted more accurate counseling of parents whose fetuses have structural cardiac abnormalities.



## Bioeffects

---

### Summary

It is important to realize that current investigation DOES NOT support the concept that Doppler waveform analysis in the above obstetrical conditions provides information that is meant to replace other conventional tests, or the biophysical profile, also, it should not replace classical non-Doppler ultrasound scanning methods. Rather, it appears that Doppler is of value when used in conjunction with these other tests. Lastly, there is little data to suggest that an isolated Doppler examination of the fetus can be used as a screening tool in a low-risk population to identify the compromised fetus.

As with the evolution of any new diagnostic technique, the indications for use may change when new experimental data become available.

### Suggested Guidelines

Following are suggested guidelines prepared by Harold Schulman, M.D.<sup>17</sup>

### Overview

There are two types of Doppler instruments, continuous wave (CW) and pulsed. Modern instruments use directional Dopplers, that is they portray forward and reverse flow on a split screen. The CW Doppler transverses the entire vessel diameter and summarizes a variety of reflections, including those from the vessel wall and neighboring vessels. The pulsed Doppler is focused and may sample the red cell flow in different loci within the radius, but has the advantage of rejecting extraneous signals from other vessels. Its disadvantage is in that some current instrumentation utilize output energies which exceed guidelines for safety on obstetric ultrasound.





## **Methodology**

Doppler studies on the umbilical arteries are carried out after manual palpation of the uterus and fetus. In pregnancies before 28 weeks, the fetus cannot be easily felt hence the probe is placed on the upper third of the uterus at the midline. The probe is gently rotated until the umbilical signal is heard and clearly identified on the screen. There should be 3-4 clear dense waveforms of equal amplitude. If there is fetal breathing, the waveforms will not be equal in intensity or frequency; hence, measurements must be postponed until the fetus is quiet. If there is tachycardia or bradycardia, the measurements have no clear meaning because now the waveform represents a complex mixture of resistance, cardiac output, and time needed to empty a vascular region. The image is frozen and the measurements are taken. We have found it useful to continuously record the examination on audiotape because the signal is immediately captured and not lost when the fetus or cord moves.

This process is repeated at least twice more and the results from the three studies are averaged. Although the S/D ratio is, in theory, an angle-independent measurement, there may be differences in the ratio when the artery is studied near the abdomen or near the insertion in the placenta. When three different angle measurements are taken, we found an average experimental error of 6% and a maximum error of 16%. The uterine artery signals may be more difficult to obtain and require more training to develop expertise. The probe is directed in the lower quadrants at the paracervical area where the uterine artery enters the uterus. Three potential signals may be encountered, hence the examiner should be persistent until the correct waveform is identified. The mature waveform is generally present by 20 weeks, but should be achieved no later than 26 weeks. The pulse wave frequency should coincide with the maternal pulse and 4 equal and dense images should be seen on each side of the uterus. The results from both sides are averaged to give a single number. The S/D ratio appears to provide a simple and reproducible measurement, but a number of other calculations have also been proposed. Interobserver error with this technique is 4% with a maximum error of 10%.



## Variance studies for fetal Doppler measurements

The following guidance<sup>17</sup> is provided to help users assess the potential variation of results when performing fetal Doppler measurements among the same or between different users. In order to complete this assessment, each clinic intending to use this system for fetal Doppler should conduct a clinical trial to establish the intra- and inter-operator variances.

The recommended protocol for such a clinical trial is:

1. Select at least 20 high-risk patients.
2. For each patient, select at least one Doppler waveform from the same fetal artery.
3. For each waveform, have each of at least three trained sonographers independently complete at least three repeat measurements for each of the Doppler calculations (e.g., S/D Ratio, PI, RI, etc.).
4. Perform an appropriate Analysis of Variance (ANOVA) over all 20 patient waveforms and also for each patient waveform separately.
5. Use the results of the ANOVA to estimate both intra- and inter-operator error variance components. (The SAS procedure VARCOMP may be used.)



## **Training and User Assistance**

Maintaining awareness of potentially harmful bioeffects and being able to recognize contributing conditions is essential to minimize the risk. Gaining experience with the system and becoming familiar with controls affecting output by observing the output display will improve the user's confidence to determine the presence of risk and how to reduce it. The user manual instructions and applications training are the best methods for learning these basic skills. Get fully acquainted with the user manual and review it frequently. Contact a representative at any time to request additional training or assistance.

As indicated by their titles, the following references are intended to provide additional detailed information concerning bioeffects.

1. Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment, AIUM/NEMA, 1992<sup>1</sup>
2. Implementation of the Principle of As Low As Reasonably Achievable (ALARA) for Medical and Dental Personnel, National Council on Radiation Protection and Measurements (NCRP), Report No.107, December 31, 1990.<sup>2</sup>
3. Biological Effects of Ultrasound: Mechanisms and Clinical Implications, NCRP Report No. 74, December 30, 1983.<sup>18</sup>
4. Exposure Criteria for Medical Diagnostic Ultrasound: I. Criteria Based on Thermal Mechanisms, NCRP Report No. 113, June 1, 1992.<sup>19</sup>
5. Bioeffects Considerations for the Safety of Diagnostic Ultrasound, Journal of Ultrasound in Medicine, AIUM, September 1988.<sup>20</sup>
6. Geneva Report on Safety and Standardization in Medical Ultrasound, WFUMB, May 1990.<sup>21</sup>



## FDA Acoustic Output Tables

### Maximum output summary

The following tables list the typical maximum acoustic output levels achievable with the LOGIQ™ 500 for all probes and operational modes. It is intended that this information be useful in making ALARA decisions and selecting the most appropriate probe for the application. In accordance with US FDA Guidelines, the overall maximum acoustic SPTA intensity for LOGIQ™ 500 is limited to 720 mW/cm<sup>2</sup> and MI is limited to 1.9. Modes for which TI does not exceed 0.1 are indicated by <0.1.

#### B510

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|
| B-Mode                                       | 0.42 | 0.25         | –            | –            | 0.27 | 8.90                |
| M-Mode                                       | 0.42 | –            | 0.02         | 0.04         | 0.02 | 9.60                |
| B/M-Mode                                     | 0.42 | 0.17         | 0.02         | 0.04         | 0.20 | 15.53               |
| Pulsed Wave Doppler                          | 0.50 | –            | 0.21         | 0.32         | 0.24 | 101.80              |
| B/Pulsed Wave Doppler                        | 0.50 | 0.17         | 0.21         | 0.32         | 0.41 | 107.73              |
| B/Color Flow Doppler                         | 0.42 | 0.28         | –            | –            | 0.30 | 46.93               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.50 | 0.22         | 0.11         | 0.16         | 0.36 | 77.33               |
| M-Mode/Color M-Mode                          | 0.35 | –            | 0.11         | 0.11         | 0.13 | 19.20               |
| Color Flow with M Doppler                    | 0.35 | –            | 0.11         | 0.11         | 0.13 | 19.20               |
| CW Doppler                                   | 0.02 | –            | 0.36         | 0.67         | 0.40 | 147.26              |

Table 17–1. B510 Maximum Acoustic Output Levels

#### CWD2

| Mode                      | MI    | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> |
|---------------------------|-------|--------------|--------------|--------------|------|---------------------|
| Color Flow with M Doppler | 0.098 | –            | 0.49         | 3.296        | 1.21 | 280.805             |

Table 17–2. CWD2 Maximum Acoustic Output Levels

## CWD5

| Mode                      | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC   | I <sub>SPTA.3</sub> |
|---------------------------|------|--------------|--------------|--------------|-------|---------------------|
| Color Flow with M Doppler | 0.06 | –            | 1.919        | 3.549        | 3.557 | 485.1999            |

Table 17–3. CWD5 Maximum Acoustic Output Levels

## C364

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> | I <sub>SPPA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|---------------------|
| B-Mode                                       | 0.68 | 0.27         | –            | –            | 0.34 | 22.8                | 142.2               |
| M-Mode                                       | 0.68 | –            | 0.02         | 0.20         | 0.09 | 63.0                | 142.2               |
| B/M-Mode                                     | 0.68 | 0.18         | 0.02         | 0.20         | 0.31 | 78.3                | 142.2               |
| Pulsed Wave Doppler                          | 0.95 | –            | 1.06         | 3.23         | 2.08 | 616.1               | 158.2               |
| B/Pulsed Wave Doppler                        | 0.95 | 0.18         | 1.06         | 3.23         | 2.30 | 631.4               | 158.2               |
| B/Color Flow Doppler                         | 0.78 | 1.01         | –            | –            | 1.85 | 156.5               | 142.2               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.95 | 0.60         | 0.53         | 1.61         | 2.08 | 393.9               | 158.2               |
| M-Mode/Color M-Mode                          | 0.68 | –            | 0.33         | 1.92         | 0.67 | 700.1               | 142.2               |
| Color Flow with M Doppler                    | 0.68 | –            | 1.06         | 3.23         | 2.08 | 616.1               | 142.2               |

Table 17–4. C364 Maximum Acoustic Output Levels

## C358

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> | I <sub>SPPA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|---------------------|
| B-Mode                                       | 1.23 | 1.12         | –            | –            | 2.13 | 124.8               | 343.5               |
| M-Mode                                       | 1.23 | –            | 0.29         | 1.08         | 0.53 | 268.9               | 343.5               |
| B/M-Mode                                     | 1.23 | 0.75         | 0.29         | 1.08         | 1.95 | 352.1               | 343.5               |
| Pulsed Wave Doppler                          | 0.91 | –            | 0.32         | 2.05         | 1.36 | 508.7               | 310.3               |
| B/Pulsed Wave Doppler                        | 1.23 | 0.75         | 0.32         | 2.05         | 2.78 | 591.9               | 343.5               |
| B/Color Flow Doppler                         | 1.23 | 1.94         | –            | –            | 2.46 | 107.9               | 343.5               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.23 | 1.34         | 0.16         | 1.03         | 2.62 | 349.9               | 343.5               |
| M-Mode/Color M-Mode                          | 1.23 | –            | 0.80         | 3.38         | 2.78 | 717.3               | 343.5               |
| Color Flow with M Doppler                    | 1.23 | –            | 0.80         | 3.38         | 2.78 | 717.3               | 343.5               |

Table 17–5. C358 Maximum Acoustic Output Levels



## FDA Acoustic Output Data

---

### C386

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|
| B-Mode                                       | 0.62 | 0.25         | –            | –            | 0.52 | 4.5                 |
| M-Mode                                       | 0.62 | –            | <0.1         | 0.17         | 0.10 | 51.7                |
| B/M-Mode                                     | 0.62 | 0.17         | <0.1         | 0.17         | 0.45 | 54.7                |
| Pulsed Wave Doppler                          | 0.61 | –            | 0.15         | 0.93         | 0.60 | 376.1               |
| B/Pulsed Wave Doppler                        | 0.62 | 0.17         | 0.15         | 0.93         | 0.95 | 379.1               |
| B/Color Flow Doppler                         | 0.85 | 1.39         | –            | –            | 1.93 | 140.3               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.85 | 0.78         | <0.1         | 0.46         | 1.44 | 259.7               |
| M-Mode/Color M-Mode                          | 0.85 | –            | 0.45         | 0.97         | 1.58 | 117.2               |
| Color Flow with M Doppler                    | 0.85 | –            | 0.45         | 0.97         | 1.58 | 117.2               |

Table 17–6. C386 Maximum Acoustic Output Levels

### C551

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|
| B-Mode                                       | 0.76 | 0.64         | –            | –            | 0.92 | 51.0                |
| M-Mode                                       | 0.76 | –            | 0.26         | 0.30         | 0.32 | 59.6                |
| B/M-Mode                                     | 0.76 | 0.43         | 0.26         | 0.30         | 0.93 | 93.6                |
| Pulsed Wave Doppler                          | 0.98 | –            | 0.54         | 2.38         | 1.83 | 656.9               |
| B/Pulsed Wave Doppler                        | 0.98 | 0.43         | 0.54         | 2.38         | 2.45 | 690.9               |
| B/Color Flow Doppler                         | 1.18 | 1.25         | –            | –            | 1.38 | 111.9               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.18 | 0.84         | 0.27         | 1.19         | 1.91 | 401.4               |
| M-Mode/Color M-Mode                          | 0.76 | –            | 0.37         | 1.40         | 0.68 | 641.7               |
| Color Flow with M Doppler                    | 0.76 | –            | 0.37         | 1.40         | 0.68 | 641.7               |

Table 17–7. C551 Maximum Acoustic Output Levels

## C721

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC   | I <sub>SPTA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|-------|---------------------|
| B-Mode                                       | 0.46 | 0.35         | –            | –            | 0.34  | 10.3                |
| M-Mode                                       | 0.35 | –            | < 0.1        | < 0.1        | < 0.1 | 15.6                |
| B/M-Mode                                     | 0.46 | 0.23         | < 0.1        | < 0.1        | 0.26  | 22.4                |
| Pulsed Wave Doppler                          | 0.42 | –            | 0.31         | 0.24         | 0.28  | 91.8                |
| B/Pulsed Wave Doppler                        | 0.46 | 0.23         | 0.31         | 0.24         | 0.50  | 98.7                |
| B/Color Flow Doppler                         | 0.46 | 0.85         | –            | –            | 0.96  | 43.8                |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.46 | 0.54         | 0.15         | 0.12         | 0.73  | 71.2                |
| M-Mode/Color M-Mode                          | 0.30 | –            | 0.62         | 0.52         | 0.73  | 87.5                |
| Color Flow with M Doppler                    | 0.30 | –            | 0.62         | 0.52         | 0.73  | 87.5                |

Table 17–8. C721 Maximum Acoustic Output Levels

## E721

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> | I <sub>SPPA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|---------------------|
| B-Mode                                       | 0.81 | 0.24         | –            | –            | 0.28 | 39.0                | 283.4               |
| M-Mode                                       | 0.81 | –            | 0.10         | 0.12         | 0.09 | 34.8                | 283.4               |
| B/M-Mode                                     | 0.81 | 0.16         | 0.10         | 0.12         | 0.28 | 60.7                | 283.4               |
| Pulsed Wave Doppler                          | 0.74 | –            | 0.67         | 1.42         | 0.70 | 621.8               | 253.9               |
| B/Pulsed Wave Doppler                        | 0.81 | 0.16         | 0.67         | 1.42         | 0.88 | 647.8               | 283.4               |
| B/Color Flow Doppler                         | 0.81 | 0.39         | –            | –            | 0.42 | 36.5                | 283.4               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.81 | 0.27         | 0.33         | 0.71         | 0.65 | 342.1               | 283.4               |
| M-Mode/Color M-Mode                          | 0.81 | –            | 0.73         | 0.68         | 0.76 | 150.2               | 283.4               |
| Color Flow with M Doppler                    | 0.81 | –            | 0.73         | 0.68         | 0.76 | 150.2               | 283.4               |

Table 17–9. E721 Maximum Acoustic Output Levels



## FDA Acoustic Output Data

---

### ERB7 (Convex)

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|
| B-Mode                                       | 1.01 | 0.22         | –            | –            | 0.33 | 41.1                |
| M-Mode                                       | 1.01 | –            | 0.13         | 0.18         | 0.16 | 54.3                |
| B/M-Mode                                     | 1.01 | 0.14         | 0.13         | 0.18         | 0.38 | 81.7                |
| Pulsed Wave Doppler                          | 1.03 | –            | 0.17         | 0.54         | 0.24 | 254.5               |
| B/Pulsed Wave Doppler                        | 1.03 | 0.14         | 0.17         | 0.54         | 0.46 | 281.9               |
| B/Color Flow Doppler                         | 1.01 | 0.22         | –            | –            | 0.32 | 32.2                |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.03 | 0.18         | < 0.1        | 0.27         | 0.39 | 157.1               |
| M-Mode/Color M-Mode                          | 1.01 | –            | 0.10         | 0.32         | 0.13 | 175.8               |
| Color Flow with M Doppler                    | 1.01 | –            | 0.17         | 0.54         | 0.24 | 254.5               |

Table 17–10. ERB7 (Convex) Maximum Output Levels

### ERB7 (Linear)

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|
| B-Mode                                       | 0.70 | 0.15         | –            | –            | 0.24 | 14.7                |
| M-Mode                                       | 0.70 | –            | 0.10         | 0.12         | 0.10 | 31.4                |
| B/M-Mode                                     | 0.70 | 0.10         | 0.10         | 0.12         | 0.26 | 41.2                |
| Pulsed Wave Doppler                          | 0.95 | –            | 0.29         | 0.80         | 0.33 | 397.2               |
| B/Pulsed Wave Doppler                        | 0.95 | 0.10         | 0.29         | 0.80         | 0.49 | 407.0               |
| B/Color Flow Doppler                         | 0.91 | 0.95         | –            | –            | 1.40 | 87.2                |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.95 | 0.53         | 0.15         | 0.40         | 0.94 | 247.1               |
| M-Mode/Color M-Mode                          | 0.70 | –            | 0.51         | 1.41         | 0.73 | 425.9               |
| Color Flow with M Doppler                    | 0.70 | –            | 0.29         | 0.80         | 0.33 | 397.2               |

Table 17–11. ERB7 (Linear) Maximum Output Levels



**I739**

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|
| B-Mode                                       | 0.68 | 0.16         | –            | –            | 0.16 | 25.2                |
| M-Mode                                       | 0.68 | –            | <0.1         | <0.1         | <0.1 | 34.8                |
| B/M-Mode                                     | 0.68 | 0.11         | <0.1         | <0.1         | 0.13 | 51.6                |
| Pulsed Wave Doppler                          | 0.67 | –            | 0.36         | 0.96         | 0.71 | 476.4               |
| B/Pulsed Wave Doppler                        | 0.68 | 0.11         | 0.36         | 0.96         | 0.82 | 493.2               |
| B/Color Flow Doppler                         | 0.82 | 0.69         | –            | –            | 0.85 | 219.4               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.82 | 0.40         | 0.18         | 0.48         | 0.83 | 356.3               |
| M-Mode/Color M-Mode                          | 0.82 | –            | 0.59         | 0.89         | 0.74 | 336.3               |
| Color Flow with M Doppler                    | 0.82 | –            | 0.59         | 0.89         | 0.74 | 336.3               |

Table 17–12. I739 Maximum Acoustic Output Levels

**546L**

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> | I <sub>SPPA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|---------------------|
| B-Mode                                       | 1.13 | 0.94         | –            | –            | 1.32 | 105.9               | 808.6               |
| M-Mode                                       | 1.13 | –            | 0.11         | 0.38         | 0.16 | 149.6               | 808.6               |
| B/M-Mode                                     | 1.13 | 0.63         | 0.11         | 0.38         | 1.03 | 220.2               | 808.6               |
| Pulsed Wave Doppler                          | 1.26 | –            | 0.39         | 1.25         | 0.54 | 529.2               | 832.8               |
| B/Pulsed Wave Doppler                        | 1.26 | 0.63         | 0.39         | 1.25         | 1.41 | 559.8               | 832.8               |
| B/Color Flow Doppler                         | 1.42 | 2.88         | –            | –            | 3.94 | 351.5               | 808.6               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.42 | 1.75         | 0.20         | 0.62         | 2.68 | 475.6               | 832.8               |
| M-Mode/Color M-Mode                          | 1.13 | –            | 0.52         | 1.64         | 0.72 | 575.4               | 808.6               |
| Color Flow with M Doppler                    | 1.13 | –            | 0.52         | 1.64         | 0.72 | 575.4               | 808.6               |

Table 17–13. 546L Maximum Acoustic Output Levels



## FDA Acoustic Output Data

---

### 739L

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> | I <sub>SPPA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|---------------------|
| B-Mode                                       | 1.39 | 0.91         | –            | –            | 1.04 | 87.1                | 1117.2              |
| M-Mode                                       | 1.39 | –            | < 0.1        | 0.28         | 0.10 | 139.1               | 1117.2              |
| B/M-Mode                                     | 1.39 | 0.61         | < 0.1        | 0.28         | 0.80 | 197.2               | 1117.2              |
| Pulsed Wave Doppler                          | 1.67 | –            | 0.36         | 1.53         | 0.64 | 560.6               | 817.2               |
| B/Pulsed Wave Doppler                        | 1.67 | 0.61         | 0.36         | 1.53         | 1.33 | 618.7               | 1117.2              |
| B/Color Flow Doppler                         | 1.39 | 3.16         | –            | –            | 3.85 | 417.2               | 1117.2              |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.67 | 1.88         | 0.18         | 0.77         | 2.59 | 518.0               | 1117.2              |
| M-Mode/Color M-Mode                          | 1.39 | –            | 1.04         | 2.00         | 1.28 | 540.0               | 1117.2              |
| Color Flow with M Doppler                    | 1.39 | –            | 1.04         | 2.00         | 1.28 | 540.0               | 1117.2              |

Table 17–14. 739L Maximum Acoustic Output Levels

### L764

| Mode                                         | MI  | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC   | I <sub>SPTA.3</sub> |
|----------------------------------------------|-----|--------------|--------------|--------------|-------|---------------------|
| B-Mode                                       | 0.4 | 0.2          | –            | –            | 0.2   | 4.4                 |
| M-Mode                                       | 0.4 | –            | < 0.1        | < 0.1        | < 0.1 | 2.8                 |
| B/M-Mode                                     | 0.4 | 0.2          | < 0.1        | < 0.1        | < 0.1 | 7.1                 |
| Pulsed Wave Doppler                          | 0.4 | –            | 0.1          | 0.5          | 0.2   | 226.4               |
| B/Pulsed Wave Doppler                        | 0.4 | 0.2          | 0.1          | 0.5          | 0.4   | 230.7               |
| B/Color Flow Doppler                         | 0.4 | 0.5          | –            | –            | 0.5   | 46.1                |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.4 | 0.5          | 0.1          | 1.6          | 0.7   | 272.4               |
| M-Mode/Color M-Mode                          | 0.4 | –            | < 0.1        | < 0.1        | 0.3   | 44.5                |
| Color Flow with M Doppler                    | 0.4 | 0.5          | < 0.1        | < 0.1        | 0.5   | 48.8                |

Table 17–15. L764 Maximum Acoustic Output Levels

## LA39

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA,3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|
| B-Mode                                       | 1.52 | 0.38         | –            | –            | 0.63 | 110.8               |
| M-Mode                                       | 1.52 | –            | < 0.1        | 0.14         | 0.06 | 203.2               |
| B/M-Mode                                     | 1.52 | 0.25         | < 0.1        | 0.14         | 0.49 | 277.0               |
| Pulsed Wave Doppler                          | 1.36 | –            | 0.24         | 1.08         | 0.59 | 572.0               |
| B/Pulsed Wave Doppler                        | 1.52 | 0.25         | 0.24         | 1.08         | 1.01 | 645.8               |
| B/Color Flow Doppler                         | 1.52 | 0.93         | –            | –            | 2.10 | 389.6               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.52 | 0.59         | 0.12         | 0.54         | 1.56 | 517.7               |
| M-Mode/Color M-Mode                          | 1.52 | –            | 0.57         | 1.68         | 1.40 | 604.3               |
| Color Flow with M Doppler                    | 1.52 | –            | 0.24         | 1.08         | 0.59 | 572.0               |

Table 17–16. LA39 Maximum Acoustic Output Levels

## LD

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA,3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|
| B-Mode                                       | 0.63 | 0.21         | –            | –            | 0.49 | 4.4                 |
| M-Mode                                       | 0.63 | –            | <0.1         | <0.1         | <0.1 | 26.1                |
| B/M-Mode                                     | 0.63 | 0.14         | <0.1         | <0.1         | 0.38 | 29.1                |
| Pulsed Wave Doppler                          | 0.58 | –            | 0.33         | 2.32         | 1.70 | 680.6               |
| B/Pulsed Wave Doppler                        | 0.63 | 0.14         | 0.33         | 2.32         | 2.02 | 683.6               |
| B/Color Flow Doppler                         | 0.90 | 1.35         | –            | –            | 1.37 | 82.6                |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.90 | 0.74         | 0.16         | 1.16         | 1.69 | 383.1               |
| M-Mode/Color M-Mode                          | 0.79 | –            | 0.29         | 2.13         | 1.23 | 655.5               |
| Color Flow with M Doppler                    | 0.79 | –            | 0.29         | 2.13         | 1.23 | 655.5               |

Table 17–17. LD Maximum Acoustic Output



## FDA Acoustic Output Data

---

### P509

| Mode                                         | MI    | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC   | I <sub>SPTA,3</sub> |
|----------------------------------------------|-------|--------------|--------------|--------------|-------|---------------------|
| B-Mode                                       | 0.88  | 0.25         | –            | –            | 0.26  | 287.1               |
| M-Mode                                       | 0.88  | –            | < 0.1        | 0.18         | < 0.1 | 43.9                |
| B/M-Mode                                     | 0.88  | 0.16         | < 0.1        | 0.18         | 0.26  | 235.3               |
| Pulsed Wave Doppler                          | 0.81  | –            | 0.56         | 1.63         | 0.66  | 405.9               |
| B/Pulsed Wave Doppler                        | 0.88  | 0.16         | 0.56         | 1.63         | 0.83  | 597.4               |
| B/Color Flow Doppler                         | 0.88  | 0.63         | –            | –            | 0.72  | 470.7               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.88  | 0.40         | 0.28         | 0.82         | 0.77  | 534.0               |
| M-Mode/Color M-Mode                          | 0.88  | –            | 0.48         | 1.56         | 0.54  | 455.4               |
| Color Flow with M Doppler                    | 0.88  | –            | 0.48         | 1.56         | 0.54  | 455.4               |
| CWD                                          | < 0.1 | –            | 1.33         | 1.73         | 2.49  | 232.8               |

Table 17–18. L764 Maximum Acoustic Output Levels

### S220

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA,3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|
| B-Mode                                       | 0.51 | 0.44         | –            | –            | 0.64 | 201.50              |
| M-Mode                                       | 0.51 | –            | 0.07         | 0.20         | 0.10 | 53.80               |
| B/M-Mode                                     | 0.51 | 0.29         | 0.07         | 0.20         | 0.52 | 188.13              |
| Pulsed Wave Doppler                          | 0.58 | –            | 0.59         | 2.52         | 1.12 | 334.94              |
| B/Pulsed Wave Doppler                        | 0.58 | 0.29         | 0.59         | 2.52         | 1.55 | 469.27              |
| B/Color Flow Doppler                         | 0.83 | 1.12         | –            | –            | 1.35 | 462.63              |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.83 | 0.71         | 0.29         | 1.26         | 1.45 | 465.95              |
| M-Mode/Color M-Mode                          | 0.83 | –            | 0.53         | 2.71         | 0.93 | 496.88              |
| Color Flow with M Doppler                    | 0.83 | –            | 0.53         | 2.71         | 0.93 | 496.88              |
| CW Doppler                                   | 0.03 | –            | 0.84         | 1.16         | 2.17 | 183.66              |

Table 17–19. S220 Maximum Acoustic Output Levels

## S222

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> | I <sub>SPPA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|---------------------|
| B-Mode                                       | 1.44 | 0.85         | –            | –            | 1.97 | 382.4               | 587.8               |
| M-Mode                                       | 1.44 | –            | 0.17         | 1.75         | 0.65 | 446.8               | 587.8               |
| B/M-Mode                                     | 1.44 | 0.57         | 0.17         | 1.75         | 1.96 | 701.7               | 587.8               |
| Pulsed Wave Doppler                          | 0.70 | –            | 0.22         | 2.27         | 0.84 | 674.2               | 141.6               |
| B/Pulsed Wave Doppler                        | 1.44 | < 0.1        | 0.22         | 2.27         | 0.87 | 680.0               | 587.8               |
| B/Color Flow Doppler                         | 1.44 | 2.05         | –            | –            | 2.41 | 346.0               | 587.8               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.44 | 1.03         | 0.11         | 1.14         | 1.63 | 512.0               | 587.8               |
| M-Mode/Color M-Mode                          | 1.44 | –            | 0.21         | 4.01         | 0.86 | 578.2               | 587.8               |
| Color Flow with M Doppler                    | 1.44 | –            | 0.22         | 2.27         | 0.84 | 674.2               | 587.8               |
| CW Doppler                                   | 0.15 | –            | 2.20         | 2.25         | 4.32 | 406.6               | –                   |

Table 17–20. S222 Maximum Acoustic Output Levels

## S316

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|
| B-Mode                                       | 0.64 | 0.16         | –            | –            | 0.26 | 24.6                |
| M-Mode                                       | 0.64 | –            | <0.1         | 0.14         | <0.1 | 40.3                |
| B/M-Mode                                     | 0.64 | 0.11         | <0.1         | 0.14         | 0.25 | 56.7                |
| Pulsed Wave Doppler                          | 0.70 | –            | 0.11         | 0.88         | 0.45 | 247.8               |
| B/Pulsed Wave Doppler                        | 0.70 | 0.11         | 0.11         | 0.88         | 0.62 | 264.2               |
| B/Color Flow Doppler                         | 0.72 | 1.24         | –            | –            | 1.47 | 71.2                |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.72 | 0.67         | <0.1         | 0.44         | 1.05 | 167.7               |
| M-Mode/Color M-Mode                          | 0.64 | –            | 0.16         | 0.66         | 0.58 | 124.5               |
| Color Flow with M Doppler                    | 0.64 | –            | 0.16         | 0.66         | 0.58 | 124.5               |
| CW Doppler                                   | 0.02 | –            | 0.8          | 0.5          | 3.1  | 268                 |

Table 17–21. S316 Maximum Acoustic Output Levels



## FDA Acoustic Output Data

---

### S317

| Mode                                         | MI    | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> | I <sub>SPPA.3</sub> |
|----------------------------------------------|-------|--------------|--------------|--------------|------|---------------------|---------------------|
| B-Mode                                       | 1.09  | 1.34         | –            | –            | 2.19 | 314.3               | 452.8               |
| M-Mode                                       | 1.09  | –            | 0.11         | 0.91         | 0.36 | 238.2               | 452.8               |
| B/M-Mode                                     | 1.09  | 0.89         | 0.11         | 0.91         | 1.82 | 447.7               | 452.8               |
| Pulsed Wave Doppler                          | 0.57  | –            | 0.26         | 2.01         | 1.07 | 482.6               | 169.1               |
| B/Pulsed Wave Doppler                        | 1.09  | 0.89         | 0.26         | 2.01         | 2.52 | 692.2               | 452.8               |
| B/Color Flow Doppler                         | 1.09  | 2.67         | –            | –            | 3.66 | 287.2               | 452.8               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.09  | 1.78         | 0.13         | 1.00         | 3.09 | 489.7               | 452.8               |
| M-Mode/Color M-Mode                          | 1.09  | –            | 0.34         | 2.41         | 1.23 | 668.0               | 452.8               |
| Color Flow with M Doppler                    | 1.09  | –            | 0.34         | 2.41         | 1.23 | 668.0               | 452.8               |
| CW Doppler                                   | < 0.1 | –            | 1.20         | 3.07         | 2.37 | 482.5               | –                   |

Table 17–22. S317 Maximum Acoustic Output Levels

### S611

| Mode                                         | MI   | TIS/TIB Scan | TIS Non-Scan | TIB Non-Scan | TIC  | I <sub>SPTA.3</sub> |
|----------------------------------------------|------|--------------|--------------|--------------|------|---------------------|
| B-Mode                                       | 0.74 | 0.44         | –            | –            | 0.48 | 89.5                |
| M-Mode                                       | 0.74 | –            | <0.1         | 0.13         | <0.1 | 57.4                |
| B/M-Mode                                     | 0.74 | 0.29         | <0.1         | 0.13         | 0.37 | 117.1               |
| Pulsed Wave Doppler                          | 0.72 | –            | 0.96         | 1.57         | 0.94 | 584.4               |
| B/Pulsed Wave Doppler                        | 0.74 | 0.29         | 0.96         | 1.57         | 1.26 | 644.1               |
| B/Color Flow Doppler                         | 0.74 | 1.68         | –            | –            | 1.99 | 311.3               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.74 | 0.98         | 0.48         | 0.78         | 1.63 | 477.7               |
| M-Mode/Color M-Mode                          | 0.51 | –            | 1.38         | 1.53         | 1.67 | 710.6               |
| Color Flow with M Doppler                    | 0.51 | –            | 1.38         | 1.53         | 1.67 | 710.6               |
| CW Doppler                                   | 0.05 | –            | 1.21         | 2.85         | 2.41 | 579.6               |

Table 17–23. S611 Maximum Acoustic Output Levels



---

*FDA Acoustic Output Data*

---

**T739**

| Mode                                         | MI   | TIS/TIB<br>Scan | TIS<br>Non-Scan | TIB<br>Non-Scan | TIC  | I <sub>SPTA,3</sub> |
|----------------------------------------------|------|-----------------|-----------------|-----------------|------|---------------------|
| B-Mode                                       | 0.72 | 0.13            | –               | –               | 0.13 | 22.9                |
| M-Mode                                       | 0.72 | –               | <0.1            | <0.1            | <0.1 | 25.5                |
| B/M-Mode                                     | 0.72 | <0.1            | <0.1            | <0.1            | 0.11 | 40.8                |
| Pulsed Wave Doppler                          | 0.71 | –               | 0.28            | 0.98            | 0.56 | 474.1               |
| B/Pulsed Wave Doppler                        | 0.72 | <0.1            | 0.28            | 0.98            | 0.64 | 489.4               |
| B/Color Flow Doppler                         | 0.93 | 0.48            | –               | –               | 0.83 | 248.8               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.93 | 0.28            | 0.14            | 0.49            | 0.74 | 369.1               |
| M-Mode/Color M-Mode                          | 0.93 | –               | 0.39            | 0.89            | 0.74 | 307.5               |
| Color Flow with M Doppler                    | 0.93 | –               | 0.39            | 0.89            | 0.74 | 307.5               |

Table 17–24. T739 Maximum Acoustic Output Levels



## FDA Acoustic Output Data

### Summary of system control settings that yield maximum output

| Operating Mode                               | B510                                                                | C364                                                                 | C386                                                                 | C358                                                                  | C551                                                                |
|----------------------------------------------|---------------------------------------------------------------------|----------------------------------------------------------------------|----------------------------------------------------------------------|-----------------------------------------------------------------------|---------------------------------------------------------------------|
| B-Mode                                       | Depth=4                                                             | Depth=16                                                             | Depth=10                                                             | Depth=16                                                              | Depth=20                                                            |
| M-Mode                                       | Depth=4                                                             | Depth=16                                                             | Depth=10                                                             | Depth=16                                                              | Depth=20                                                            |
| B/M-Mode                                     | Depth=4                                                             | Depth=16                                                             | Depth=10                                                             | Depth=16                                                              | Depth=20                                                            |
| CWD-Mode                                     | Depth=5.4                                                           | Depth= na                                                            | Depth=na                                                             | Depth= na                                                             | Depth= na                                                           |
| Pulsed Wave Doppler                          | Vel<br>Scale= $\pm 0.33$<br>SV=2<br>Penet=On                        | Vel<br>Scale= $\pm 7.9$<br>SV=2<br>Penet=Off                         | Vel<br>Scale= $\pm 3.96$<br>SV=2<br>Penet=Off                        | Vel<br>Scale= $\pm 0.33$<br>SV=1<br>Penet=Off                         | Vel<br>Scale= $\pm 0.5$<br>SV=6<br>Penet=Off                        |
| B/Pulsed Wave Doppler                        | Vel<br>Scale= $\pm 0.33$<br>SV=2<br>Penet=On                        | Vel<br>Scale= $\pm 7.9$<br>SV=2<br>Penet=Off                         | Vel<br>Scale= $\pm 3.96$<br>SV=2<br>Penet=Off                        | Vel<br>Scale= $\pm 0.33$<br>SV=1<br>Penet=Off                         | Vel<br>Scale= $\pm 0.5$<br>SV=6<br>Penet=Off                        |
| B/Color Flow Doppler                         | Vel<br>Scale= $\pm 1.3$<br>Penet=Off<br>ROI=Min<br>PS=8             | Vel<br>Scale= $\pm 5.2$<br>Penet=Off<br>ROI=Min<br>PS=16             | Vel<br>Scale= $\pm 0.21$<br>Penet=On<br>ROI=Min<br>PS=16             | Vel<br>Scale= $\pm 0.14$<br>Penet=Off<br>ROI=Max<br>PS=16             | Vel<br>Scale= $\pm 3.42$<br>Penet=On<br>ROI=Min<br>PS=10            |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Vel<br>Scale= $\pm 1.3$<br>Penet=Off<br>ROI=Min<br>PS=8             | Vel<br>Scale= $\pm 5.2$<br>Penet=Off<br>ROI=Min<br>PS=16             | Vel<br>Scale= $\pm 0.21$<br>Penet=On<br>ROI=Min<br>PS=16             | Vel<br>Scale= $\pm 0.14$<br>Penet=Off<br>ROI=Max<br>PS=16             | Vel<br>Scale= $\pm 3.42$<br>Penet=On<br>ROI=Min<br>PS=10            |
| M-Mode/<br>Color M-Mode                      | Vel<br>Scale= $\pm 1.3$<br>Penet=Off<br>ROI=Min<br>Depth=4<br>PS=16 | Vel<br>Scale= $\pm 5.2$<br>Penet=Off<br>ROI=Min<br>Depth=10<br>PS=16 | Vel<br>Scale= $\pm 0.21$<br>Penet=On<br>ROI=Min<br>Depth=10<br>PS=16 | Vel<br>Scale= $\pm 0.74$<br>Penet=Off<br>ROI=Max<br>Depth=24<br>PS=16 | Vel<br>Scale= $\pm 3.42$<br>Penet=On<br>ROI=Min<br>Depth=5<br>PS=16 |
| Color Flow with<br>M Doppler                 | Vel<br>Scale= $\pm 1.3$<br>Penet=Off<br>ROI=Min<br>Depth=4<br>PS=16 | Vel<br>Scale= $\pm 5.2$<br>Penet=Off<br>ROI=Min<br>Depth=10<br>PS=16 | Vel<br>Scale= $\pm 0.21$<br>Penet=On<br>ROI=Min<br>Depth=10<br>PS=16 | Vel<br>Scale= $\pm 0.74$<br>Penet=Off<br>ROI=Max<br>Depth=24<br>PS=16 | Vel<br>Scale= $\pm 3.42$<br>Penet=On<br>ROI=Min<br>Depth=5<br>PS=16 |

Table 17–25. Maximum Output Control Settings (continued on next page)



Summary of system control settings that yield maximum output (cont'd)

| Operating Mode                               | C721                                                            | E721                                                                   | I739                                                            | 546L                                                                  | 739L                                                            |
|----------------------------------------------|-----------------------------------------------------------------|------------------------------------------------------------------------|-----------------------------------------------------------------|-----------------------------------------------------------------------|-----------------------------------------------------------------|
| B-Mode                                       | Depth=4                                                         | Depth=14                                                               | Depth=4                                                         | Depth=16                                                              | Depth=4                                                         |
| M-Mode                                       | Depth=4                                                         | Depth=14                                                               | Depth=4                                                         | Depth=16                                                              | Depth=4                                                         |
| B/M-Mode                                     | Depth=4                                                         | Depth=14                                                               | Depth=4                                                         | Depth=16                                                              | Depth=4                                                         |
| CWD-Mode                                     | Depth= na                                                       | Depth= na                                                              | Depth=na                                                        | Depth=na                                                              | Depth=na                                                        |
| Pulsed Wave Doppler                          | Vel<br>Scale= ± 3.0<br>SV=2<br>Penet=Off                        | Vel<br>Scale= ± 0.33<br>SV=1<br>Penet=Off                              | Vel<br>Scale= ± 9.1<br>SV=1<br>Penet=On                         | Vel<br>Scale= ± 0.33<br>SV=1<br>Penet=On                              | Vel<br>Scale= ± 0.33<br>SV=1<br>Penet=On                        |
| B/Pulsed Wave Doppler                        | Vel<br>Scale= ± 3.0<br>SV=2<br>Penet= Off                       | Vel<br>Scale= ± 0.33<br>SV=1<br>Penet=Off                              | Vel<br>Scale= ± 9.1<br>SV=1<br>Penet=On                         | Vel<br>Scale= ± 0.33<br>SV=1<br>Penet=On                              | Vel<br>Scale= ± 0.33<br>SV=1<br>Penet=On                        |
| B/Color Flow Doppler                         | Vel<br>Scale= ± 0.98<br>Penet=On<br>ROI=Max<br>PS=10            | Vel<br>Scale= ± 0.14<br>Penet=Off<br>ROI=Min –Max<br>PS=10             | Vel<br>Scale= ± 1.30<br>Penet=On<br>ROI=Min<br>PS=16            | Vel<br>Scale= ± 0.14<br>Penet=On<br>ROI=Min –Max<br>PS=10             | Vel<br>Scale= ± 0.86<br>Penet=On<br>ROI=Min<br>PS=16            |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Vel<br>Scale= ± 0.98<br>Penet=On<br>ROI=Max<br>PS=10            | Vel<br>Scale= ± 0.14<br>Penet=Off<br>ROI=Min –Max<br>PS=10             | Vel<br>Scale= ± 1.30<br>Penet=On<br>ROI=Min<br>PS=16            | Vel<br>Scale= ± 0.14<br>Penet=On<br>ROI=Min –Max<br>PS=10             | Vel<br>Scale= ± 0.86<br>Penet=On<br>ROI=Min<br>PS=16            |
| M-Mode/<br>Color M-Mode                      | Vel<br>Scale= ± 0.98<br>Penet=On<br>ROI=Max<br>Depth=6<br>PS=10 | Vel<br>Scale= ± 0.74<br>Penet=Off<br>ROI=Min –Max<br>Depth=14<br>PS=16 | Vel<br>Scale= ± 1.30<br>Penet=On<br>ROI=Min<br>Depth=4<br>PS=16 | Vel<br>Scale= ± 0.74<br>Penet=On<br>ROI=Min –Max<br>Depth=16<br>PS=16 | Vel<br>Scale= ± 0.85<br>Penet=On<br>ROI=Min<br>Depth=4<br>PS=16 |
| Color Flow with<br>M Doppler                 | Vel<br>Scale= ± 0.98<br>Penet=On<br>ROI=Max<br>Depth=6<br>PS=10 | Vel<br>Scale= ± 0.74<br>Penet=Off<br>ROI=Min –Max<br>Depth=14<br>PS=16 | Vel<br>Scale= ± 1.30<br>Penet=On<br>ROI=Min<br>Depth=4<br>PS=16 | Vel<br>Scale= ± 0.74<br>Penet=On<br>ROI=Min –Max<br>Depth=16<br>PS=16 | Vel<br>Scale= ± 0.85<br>Penet=On<br>ROI=Min<br>Depth=4<br>PS=16 |

Table 17–25. Maximum Output Control Settings (cont'd)



## FDA Acoustic Output Data

### Summary of system control settings that yield maximum output (cont'd)

| Operating Mode                               | L764                                                                | P509                                                                 | S220                                                                 | S222                                                                 | S316                                                                |
|----------------------------------------------|---------------------------------------------------------------------|----------------------------------------------------------------------|----------------------------------------------------------------------|----------------------------------------------------------------------|---------------------------------------------------------------------|
| B-Mode                                       | Depth=4                                                             | Depth=24                                                             | Depth=10                                                             | Depth=10                                                             | Depth=6                                                             |
| M-Mode                                       | Depth=4                                                             | Depth=24                                                             | Depth=10                                                             | Depth=10                                                             | Depth=6                                                             |
| B/M-Mode                                     | Depth=4                                                             | Depth=24                                                             | Depth=10                                                             | Depth=10                                                             | Depth=6                                                             |
| CWD-Mode                                     | Depth= na                                                           | Depth=11.9                                                           | Depth=10                                                             | Depth=10                                                             | Depth=6                                                             |
| Pulsed Wave Doppler                          | Vel<br>Scale= $\pm 2.0$<br>SV=1<br>Penet=On                         | Vel<br>Scale= $\pm 3.0$<br>SV=4<br>Penet=On                          | Vel<br>Scale= $\pm 2.0$<br>SV=6<br>Penet=Off                         | Vel<br>Scale= $\pm 0.86$<br>SV=3<br>Penet=Off                        | Vel<br>Scale= $\pm 2.28$<br>SV=2<br>Penet=Off                       |
| B/Pulsed Wave Doppler                        | Vel<br>Scale= $\pm 2.0$<br>SV=1<br>Penet=On                         | Vel<br>Scale= $\pm 3.0$<br>SV=4<br>Penet=On                          | Vel<br>Scale= $\pm 2.0$<br>SV=6<br>Penet=Off                         | Vel<br>Scale= $\pm 0.86$<br>SV=3<br>Penet=Off                        | Vel<br>Scale= $\pm 2.28$<br>SV=2<br>Penet=Off                       |
| B/Color Flow Doppler                         | Vel<br>Scale= $\pm 5.2$<br>Penet=Off<br>ROI=Min<br>PS=16            | Vel<br>Scale= $\pm 2.98$<br>Penet=On<br>ROI=Min<br>PS=14             | Vel<br>Scale= $\pm 3.4$<br>Penet=Off<br>ROI=Min<br>PS=14             | Vel<br>Scale= $\pm 2.98$<br>Penet=Off<br>ROI=Min<br>PS=16            | Vel<br>Scale= $\pm 2.6$<br>Penet=Off<br>ROI=Min<br>PS=14            |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Vel<br>Scale= $\pm 5.2$<br>Penet=Off<br>ROI=Min<br>PS=16            | Vel<br>Scale= $\pm 2.98$<br>Penet=On<br>ROI=Min<br>PS=14             | Vel<br>Scale= $\pm 3.4$<br>Penet=Off<br>ROI=Min<br>PS=14             | Vel<br>Scale= $\pm 2.98$<br>Penet=Off<br>ROI=Min<br>PS=16            | Vel<br>Scale= $\pm 2.6$<br>Penet=Off<br>ROI=Min<br>PS=14            |
| M-Mode/<br>Color M-Mode                      | Vel<br>Scale= $\pm 5.2$<br>Penet=Off<br>ROI=Min<br>Depth=4<br>PS=16 | Vel<br>Scale= $\pm 2.98$<br>Penet=On<br>ROI=Min<br>Depth=24<br>PS=14 | Vel<br>Scale= $\pm 3.4$<br>Penet=Off<br>ROI=Min<br>Depth=14<br>PS=14 | Vel<br>Scale= $\pm 2.98$<br>Penet=Off<br>ROI=Min<br>Depth=8<br>PS=16 | Vel<br>Scale= $\pm 2.6$<br>Penet=Off<br>ROI=Min<br>Depth=6<br>PS=14 |
| Color Flow with<br>M Doppler                 | Vel<br>Scale= $\pm 5.2$<br>Penet=Off<br>ROI=Min<br>Depth=4<br>PS=16 | Vel<br>Scale= $\pm 2.98$<br>Penet=On<br>ROI=Min<br>Depth=24<br>PS=14 | Vel<br>Scale= $\pm 3.4$<br>Penet=Off<br>ROI=Min<br>Depth=14<br>PS=14 | Vel<br>Scale= $\pm 2.98$<br>Penet=Off<br>ROI=Min<br>Depth=8<br>PS=16 | Vel<br>Scale= $\pm 2.6$<br>Penet=Off<br>ROI=Min<br>Depth=6<br>PS=14 |

Table 17–25. Maximum Output Control Settings (cont'd)

Summary of system control settings that yield maximum output (cont'd)

| Operating Mode                               | S317                                                           | S611                                                          | T739                                                           | LA39                                                            | LD                                                               |
|----------------------------------------------|----------------------------------------------------------------|---------------------------------------------------------------|----------------------------------------------------------------|-----------------------------------------------------------------|------------------------------------------------------------------|
| B-Mode                                       | Depth=10                                                       | Depth=4                                                       | Depth=4                                                        | Depth=4                                                         | Depth=24                                                         |
| M-Mode                                       | Depth=10                                                       | Depth=4                                                       | Depth=4                                                        | Depth=4                                                         | Depth=24                                                         |
| B/M-Mode                                     | Depth=10                                                       | Depth=4                                                       | Depth=4                                                        | Depth=4                                                         | Depth=24                                                         |
| CWD-Mode                                     | Depth=12.9                                                     | Depth=6.4                                                     | Depth=na                                                       | Depth=na                                                        | Depth=na                                                         |
| Pulsed Wave Doppler                          | Vel<br>Scale=± 0.33<br>SV=1<br>Penet=On                        | Vel<br>Scale=± 0.75<br>SV=5<br>Penet=Off                      | Vel Scale=± 9.1<br>SV=1<br>Penet=On                            | Vel<br>Scale=± 1.31<br>SV=1<br>Penet=Off                        | Vel<br>Scale=± 0.33<br>SV=1<br>Penet=Off                         |
| B/Pulsed Wave Doppler                        | Vel<br>Scale=± 0.33<br>SV=1<br>Penet=On                        | Vel<br>Scale=± 0.75<br>SV=5<br>Penet=Off                      | Vel Scale=± 9.1<br>SV=1<br>Penet=On                            | Vel<br>Scale=± 1.31<br>SV=1<br>Penet=Off                        | Vel<br>Scale=± 0.14<br>SV=1<br>Penet=Off                         |
| B/Color Flow Doppler                         | Vel<br>Scale=± 0.25<br>Penet=On<br>ROI=Min<br>PS=14            | Vel<br>Scale=± 2.98<br>Penet=On<br>ROI=Min<br>PS=8            | Vel<br>Scale=± 1.30<br>Penet=On<br>ROI=Min<br>PS=16            | Vel<br>Scale=± 0.33<br>Penet=Off<br>ROI=Min<br>PS=16            | Vel<br>Scale=± 0.14<br>Penet=Off<br>ROI=Max<br>PS=10             |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Vel<br>Scale=± 0.25<br>Penet=On<br>ROI=Min<br>PS=14            | Vel<br>Scale=± 2.98<br>Penet=On<br>ROI=Min<br>PS=8            | Vel<br>Scale=± 1.30<br>Penet=On<br>ROI=Min<br>PS=16            | Vel<br>Scale=± 0.33<br>Penet=Off<br>ROI=Min<br>PS=16            | Vel<br>Scale=± 0.74<br>Penet=Off<br>ROI=Max<br>PS=10             |
| M-Mode/<br>Color M-Mode                      | Vel<br>Scale=± 0.74<br>Penet=On<br>ROI=Min<br>Depth=8<br>PS=16 | Vel<br>Scale=± 2.98<br>Penet=On<br>ROI=Min<br>Depth=8<br>PS=8 | Vel<br>Scale=± 1.30<br>Penet=On<br>ROI=Min<br>Depth=4<br>PS=16 | Vel<br>Scale=± 0.74<br>Penet=Off<br>ROI=Min<br>Depth=4<br>PS=16 | Vel<br>Scale=± 0.74<br>Penet=Off<br>ROI=Max<br>Depth=24<br>PS=16 |
| Color Flow with<br>M Doppler                 | Vel<br>Scale=± 0.74<br>Penet=On<br>ROI=Min<br>Depth=8<br>PS=16 | Vel<br>Scale=± 2.98<br>Penet=On<br>ROI=Min<br>Depth=8<br>PS=8 | Vel<br>Scale=± 1.30<br>Penet=On<br>ROI=Min<br>Depth=4<br>PS=16 | Vel<br>Scale=± 0.74<br>Penet=Off<br>ROI=Min<br>Depth=4<br>PS=16 | Vel<br>Scale=± 0.74<br>Penet=Off<br>ROI=Max<br>Depth=24<br>PS=16 |

Table 17–25. Maximum Output Control Settings (cont'd)



**Summary of system control settings that yield maximum output (cont'd)**

| <b>Operating Mode</b>                        | <b>ERB7 (Convex)</b>                                                  | <b>ERB7 (Linear)</b>                                                |
|----------------------------------------------|-----------------------------------------------------------------------|---------------------------------------------------------------------|
| B-Mode                                       | Depth=14                                                              | Depth=14                                                            |
| M-Mode                                       | Depth=14                                                              | Depth=14                                                            |
| B/M-Mode                                     | Depth=14                                                              | Depth=14                                                            |
| CWD-Mode                                     | Depth= na                                                             | Depth=na                                                            |
| Pulsed Wave Doppler                          | Vel<br>Scale= $\pm 0.33$<br>SV=1<br>Penet=On                          | Vel<br>Scale= $\pm 0.33$<br>SV=1<br>Penet=On                        |
| B/Pulsed Wave Doppler                        | Vel<br>Scale= $\pm 0.33$<br>SV=1<br>Penet=On                          | Vel<br>Scale= $\pm 0.33$<br>SV=1<br>Penet=On                        |
| B/Color Flow Doppler                         | Vel<br>Scale= $\pm 0.14$<br>Penet=On<br>ROI=Max<br>PS=10              | Vel<br>Scale= $\pm 0.86$<br>Penet=On<br>ROI=Min<br>PS=16            |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Vel<br>Scale= $\pm 0.14$<br>Penet=On<br>ROI=Max<br>PS=10              | Vel<br>Scale= $\pm 0.86$<br>Penet=On<br>ROI=Min<br>PS=16            |
| M-Mode/<br>Color M-Mode                      | Vel<br>Scale= $\pm 0.74$<br>Penet=Off<br>ROI=Max<br>Depth=14<br>PS=16 | Vel<br>Scale= $\pm 0.85$<br>Penet=On<br>ROI=Min<br>Depth=4<br>PS=16 |
| Color Flow with<br>M Doppler                 | Vel<br>Scale= $\pm 0.74$<br>Penet=Off<br>ROI=Max<br>Depth=14<br>PS=16 | Vel<br>Scale= $\pm 0.85$<br>Penet=On<br>ROI=Min<br>Depth=4<br>PS=16 |

Table 17–25. Maximum Output Control Settings (cont'd)

## Maximum Thermal Indices

Indices for probe/mode combinations not shown are less than 1.0 for all control settings. Refer to the Key to Tables on 17–81.

### B510

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|--------------------------------------|-------|----------|--------------------------|---------------|
| Pulsed Wave Doppler                          | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.21    | 4.02  | 11.04 | 4.34                                 | —     | 1.76     | 1.00                     | 3.70          |
| B/Pulsed Wave Doppler                        | Scan          | 0.17    | 4.38  | 8.33  | 3.44                                 | —     | 1.76     | 1.00                     | 6.15          |
|                                              | Non-Scan      | 0.21    | 4.02  | 11.04 | 4.34                                 | —     | 1.76     | 1.00                     | 3.88          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.22    | 4.04  | 11.25 | 4.71                                 | —     | 1.76     | 1.00                     | 3.70          |
|                                              | Non-Scan      | 0.11    | 4.02  | 5.52  | 2.17                                 | —     | 1.76     | 1.00                     | 3.88          |
| Color Flow with<br>M Doppler                 | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.11    | 4.04  | 5.86  | 2.53                                 | —     | 1.76     | 1.00                     | 3.70          |
| CW Doppler                                   | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.36    | 4.01  | 18.76 | 7.67                                 | —     | 1.76     | 1.00                     | 5.40          |

Table 17–26. B510 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq@PII\ max}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|-------------------|
| Pulsed Wave Doppler                          | Scan          | —           | —     | —     | —        | —                | —                        | —                 |
|                                              | Non-Scan      | 0.32        | 4.02  | 11.04 | 3.36     | 0.52             | 1.00                     | —                 |
| B/Pulsed Wave Doppler                        | Scan          | 0.17        | 4.38  | 8.33  | 3.02     | 1.51             | 1.00                     | —                 |
|                                              | Non-Scan      | 0.32        | 4.02  | 11.04 | 3.36     | 0.52             | 1.00                     | —                 |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.22        | 4.04  | 11.25 | 3.07     | 0.61             | 1.00                     | —                 |
|                                              | Non-Scan      | 0.16        | 4.02  | 5.52  | 3.36     | 0.52             | 1.00                     | —                 |
| Color Flow with<br>M Doppler                 | Scan          | —           | —     | —     | —        | —                | —                        | —                 |
|                                              | Non-Scan      | 0.11        | 4.86  | 4.64  | 3.30     | 0.65             | 1.00                     | —                 |
| CW Doppler                                   | Scan          | —           | —     | —     | —        | —                | —                        | —                 |
|                                              | Non-Scan      | 0.67        | 4.01  | 18.76 | 3.23     | 0.56             | 1.00                     | —                 |

Table 17–27. B510 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### B510 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| Pulsed Wave Doppler                          | 0.24        | 4.02  | 11.04 | 1.00                            | 3.88          |
| B/Pulsed Wave Doppler                        | 0.41        | 4.02  | 19.37 | 1.00                            | 3.88          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | .036        | 4.02  | 16.78 | 1.00                            | 3.88          |
| Color Flow with M Doppler                    | 0.13        | 4.04  | 5.86  | 1.00                            | 3.70          |
| CW Doppler                                   | 0.40        | 4.01  | 18.76 | 1.00                            | 5.40          |

Table 17–28. B510 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | fc [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|----------|-----------------|----------|--------------|
| B–Mode                                    | Scan          | 0.42       | 0.90       | 4.03     | 4.51     | 2.86            | 2800     | 0.90         |
| M–Mode                                    | Non–Scan      | 0.42       | 0.90       | 4.03     | 4.51     | 2.86            | 187      | 0.90         |
| B/M–Mode                                  | Scan          | 0.42       | 0.90       | 4.03     | 4.51     | 2.86            | 2800     | 0.90         |
|                                           | Non–Scan      | 0.42       | 0.90       | 4.03     | 4.51     | 2.86            | 187      | 0.90         |
| Pulsed Wave Doppler                       | Non–Scan      | 0.50       | 1.01       | 3.36     | 4.02     | 0.98            | 660      | 1.01         |
| B/Pulsed Wave Doppler                     | Scan          | 0.42       | 0.90       | 4.03     | 4.51     | 2.86            | 2800     | 0.90         |
|                                           | Non–Scan      | 0.50       | 1.01       | 3.36     | 4.02     | 0.98            | 660      | 1.01         |
| B/Color Flow Mode                         | Scan          | 0.42       | 0.90       | 3.07     | 4.05     | 0.77            | 280      | 0.90         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.42       | 0.90       | 3.07     | 4.05     | 0.77            | 280      | 0.90         |
|                                           | Non–Scan      | 0.50       | 1.01       | 3.36     | 4.02     | 0.98            | 660      | 1.01         |
| M/M–CFM                                   | Non–Scan      | 0.42       | 0.90       | 3.07     | 4.05     | 0.77            | 280      | 0.90         |
| M–CFM/PD                                  | Non–Scan      | 0.42       | 0.90       | 3.07     | 4.05     | 0.77            | 280      | 0.90         |
| CWD Mode                                  | Non–Scan      | 0.02       | 0.11       | 3.23     | 4.01     | –               | –        | 0.11         |

Table 17–29. B510 Maximum Mechanical Index (MI)

C358

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$  | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|--------|--------------------------------------|-------|----------|--------------------------|---------------|
| B-Mode                                       | Scan          | 1.12    | 2.62  | 93.26  | 44.90                                | —     | 1.64     | 0.94                     | 3.40          |
| M-Mode                                       | Non-Scan      | 0.29    | 2.62  | 23.30  | 11.22                                | —     | 1.64     | 0.94                     | 3.40          |
| B/M-Mode                                     | Scan          | 0.75    | 2.62  | 62.17  | 29.93                                | —     | 1.64     | 0.94                     | 3.40          |
|                                              | Non-Scan      | 0.29    | 2.62  | 23.30  | 11.22                                | —     | 1.64     | 0.94                     | 3.40          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.32    | 3.32  | 99.40  | 20.38                                | —     | 2.74     | 2.63                     | 6.80          |
| B/Pulsed Wave Doppler                        | Scan          | 0.75    | 2.62  | 62.17  | 29.93                                | —     | 1.64     | 0.94                     | 3.40          |
|                                              | Non-Scan      | 0.32    | 3.32  | 99.40  | 20.38                                | —     | 2.74     | 2.63                     | 6.80          |
| B/Color Flow Doppler                         | Scan          | 1.94    | 3.31  | 138.40 | 46.59                                | —     | 2.74     | 2.63                     | 6.80          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 1.34    | 3.31  | 100.28 | 38.26                                | —     | 2.74     | 2.63                     | 6.80          |
|                                              | Non-Scan      | 0.16    | 3.32  | 49.70  | 10.19                                | —     | 2.74     | 2.63                     | 6.80          |
| M/M-CFM                                      | Non-Scan      | 0.80    | 3.31  | 187.58 | 43.50                                | —     | 2.74     | 2.63                     | 6.80          |
| M-CFM/PD                                     | Non-Scan      | 0.80    | 3.31  | 187.58 | 43.50                                | —     | 2.74     | 2.63                     | 6.80          |

Table 17–30. C358 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$  | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII-max$ |
|----------------------------------------------|---------------|-------------|-------|--------|----------|------------------|--------------------------|------------------|
| B-Mode                                       | Scan          | 1.12        | 2.62  | 93.26  | 3.92     | 0.98             | 0.94                     | —                |
| M-Mode                                       | Non-Scan      | 1.08        | 2.62  | 23.30  | 3.92     | 0.37             | 0.94                     | —                |
| B/M-Mode                                     | Scan          | 0.75        | 2.62  | 62.17  | 3.92     | 0.98             | 0.94                     | —                |
|                                              | Non-Scan      | 1.08        | 2.62  | 23.30  | 3.92     | 0.37             | 0.94                     | —                |
| Pulsed Wave Doppler                          | Non-Scan      | 2.05        | 3.32  | 99.40  | 7.05     | 0.51             | 2.63                     | —                |
| B/Pulsed Wave Doppler                        | Scan          | 0.75        | 2.62  | 62.17  | 3.92     | 0.98             | 0.94                     | —                |
|                                              | Non-Scan      | 2.05        | 3.32  | 99.40  | 7.05     | 0.51             | 2.63                     | —                |
| B/Color Flow Doppler                         | Scan          | 1.94        | 3.31  | 138.40 | 6.97     | 3.11             | 2.63                     | —                |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 1.34        | 3.31  | 100.28 | 6.97     | 3.11             | 2.63                     | —                |
|                                              | Non-Scan      | 1.03        | 3.32  | 49.70  | 7.05     | 0.51             | 2.63                     | —                |
| M/M-CFM                                      | Non-Scan      | 3.38        | 3.31  | 187.58 | 6.97     | 0.77             | 2.63                     | —                |
| M-CFM/PD                                     | Non-Scan      | 3.38        | 3.31  | 187.58 | 6.97     | 0.77             | 2.63                     | —                |

Table 17–31. C358 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### C358 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$  | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|--------|---------------------------------|---------------|
| B-Mode                                       | 2.13        | 2.62  | 93.26  | 0.94                            | 3.40          |
| M-Mode                                       | 0.53        | 2.62  | 23.30  | 0.94                            | 3.40          |
| B/M-Mode                                     | 1.95        | 2.62  | 85.47  | 0.94                            | 3.40          |
| Pulsed Wave Doppler                          | 1.36        | 3.32  | 99.40  | 2.63                            | 6.80          |
| B/Pulsed Wave Doppler                        | 2.78        | 3.32  | 161.57 | 2.63                            | 6.80          |
| B/Color Flow Doppler                         | 2.46        | 3.31  | 138.40 | 2.63                            | 6.80          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 2.62        | 3.32  | 149.99 | 2.63                            | 6.80          |
| M/M-CFM                                      | 2.78        | 3.31  | 187.58 | 2.63                            | 6.80          |
| M-CFM/PD                                     | 2.78        | 3.31  | 187.58 | 2.63                            | 6.80          |

Table 17–32. C358 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | fc [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|----------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 1.23       | 1.96       | 3.92     | 2.62     | 1.00            | 3687     | 1.96         |
| M-Mode                                    | Non-Scan      | 1.23       | 1.96       | 3.92     | 2.62     | 1.00            | 921      | 1.96         |
| B/M-Mode                                  | Scan          | 1.23       | 1.96       | 3.92     | 2.62     | 1.00            | 3987     | 1.96         |
|                                           | Non-Scan      | 1.23       | 1.96       | 3.92     | 2.62     | 1.00            | 921      | 1.96         |
| Pulsed Wave Doppler                       | Non-Scan      | 0.91       | 1.63       | 6.99     | 3.25     | 0.81            | 651      | 1.63         |
| B/Pulsed Wave Doppler                     | Scan          | 1.23       | 1.96       | 3.92     | 2.62     | 1.00            | 3687     | 1.96         |
|                                           | Non-Scan      | 0.91       | 1.63       | 6.99     | 3.25     | 0.81            | 651      | 1.63         |
| B/Color Flow Mode                         | Scan          | 1.23       | 1.96       | 6.96     | 3.26     | 0.86            | 280      | 1.96         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 1.23       | 1.96       | 6.96     | 3.26     | 0.86            | 280      | 1.96         |
|                                           | Non-Scan      | 0.91       | 1.63       | 6.99     | 3.25     | 0.81            | 651      | 1.63         |
| M/M-CFM                                   | Non-Scan      | 1.23       | 1.96       | 6.96     | 3.26     | 0.86            | 1480     | 1.96         |
| M-CFM/PD                                  | Non-Scan      | 1.23       | 1.96       | 6.96     | 3.26     | 0.86            | 1480     | 1.96         |
| CWD Mode                                  | Non-Scan      | –          | –          | –        | –        | –               | –        | –            |

Table 17–33. C358 Maximum Mechanical Index (MI)



C364

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|--------------------------------------|-------|----------|--------------------------|---------------|
| B-Mode                                       | Scan          | 0.27    | 3.36  | 16.97 | 5.86                                 | —     | 1.87     | 1.22                     | 4.20          |
| M-Mode                                       | Non-Scan      | 0.02    | 3.34  | 4.36  | 1.57                                 | —     | 1.87     | 1.22                     | 4.20          |
| B/M-Mode                                     | Scan          | 0.18    | 3.36  | 11.31 | 3.91                                 | —     | 1.87     | 1.22                     | 4.20          |
|                                              | Non-Scan      | 0.02    | 3.34  | 4.36  | 1.57                                 | —     | 1.87     | 1.22                     | 4.20          |
| Pulsed Wave Doppler                          | Non-Scan      | 1.06    | 3.29  | 66.92 | 42.30                                | —     | 1.21     | 0.51                     | 2.30          |
| B/Pulsed Wave Doppler                        | Scan          | 0.18    | 3.36  | 11.31 | 3.91                                 | —     | 1.87     | 1.22                     | 4.20          |
|                                              | Non-Scan      | 1.06    | 3.29  | 66.92 | 42.30                                | —     | 1.21     | 0.51                     | 2.30          |
| B/Color Flow Doppler                         | Scan          | 1.01    | 3.31  | 63.73 | 36.46                                | —     | 1.21     | 0.51                     | 2.30          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.60    | 3.31  | 37.52 | 20.18                                | —     | 1.21     | 0.51                     | 2.30          |
|                                              | Non-Scan      | 0.53    | 3.29  | 33.46 | 21.15                                | —     | 1.21     | 0.51                     | 2.30          |
| M/M-CFM                                      | Non-Scan      | 0.33    | 3.31  | 24.24 | 13.85                                | —     | 1.21     | 0.51                     | 2.30          |
| M-CFM/PD                                     | Non-Scan      | 1.06    | 3.31  | 24.24 | 13.85                                | —     | 1.21     | 0.51                     | 2.30          |

Table 17–34. C364 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII_{max}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| B-Mode                                       | Scan          | 0.27        | 3.36  | 16.97 | 4.61     | 0.98             | 1.22                     | —                  |
| M-Mode                                       | Non-Scan      | 0.20        | 3.34  | 4.36  | 4.57     | 0.32             | 1.22                     | —                  |
| B/M-Mode                                     | Scan          | 0.18        | 3.36  | 11.31 | 4.61     | 0.98             | 1.22                     | —                  |
|                                              | Non-Scan      | 0.20        | 3.34  | 4.36  | 4.57     | 0.32             | 1.22                     | —                  |
| Pulsed Wave Doppler                          | Non-Scan      | 3.23        | 3.29  | 66.92 | 1.91     | 0.38             | 0.51                     | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.18        | 3.36  | 11.31 | 4.61     | 0.98             | 1.22                     | —                  |
|                                              | Non-Scan      | 3.23        | 3.29  | 66.92 | 1.91     | 0.38             | 0.51                     | —                  |
| B/Color Flow Doppler                         | Scan          | 1.01        | 3.31  | 63.73 | 1.97     | 0.78             | 0.51                     | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.60        | 3.31  | 37.52 | 1.97     | 0.78             | 0.51                     | —                  |
|                                              | Non-Scan      | 1.61        | 3.29  | 33.46 | 1.91     | 0.38             | 0.51                     | —                  |
| M/M-CFM                                      | Non-Scan      | 1.92        | 3.31  | 24.24 | 1.97     | 0.21             | 0.51                     | —                  |
| M-CFM/PD                                     | Non-Scan      | 3.23        | 3.31  | 24.24 | 1.97     | 0.21             | 0.51                     | —                  |

Table 17–35. C364 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### C364 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| B-Mode                                       | 0.34        | 3.36  | 16.97 | 1.22                            | 4.20          |
| M-Mode                                       | 0.09        | 3.34  | 4.36  | 1.22                            | 4.20          |
| B/M-Mode                                     | 0.31        | 3.36  | 15.68 | 1.22                            | 4.20          |
| Pulsed Wave Doppler                          | 2.08        | 3.29  | 66.92 | 0.51                            | 2.30          |
| B/Pulsed Wave Doppler                        | 2.30        | 3.29  | 78.23 | 0.51                            | 2.30          |
| B/Color Flow Doppler                         | 1.85        | 3.31  | 63.73 | 0.51                            | 2.30          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 2.08        | 3.29  | 70.98 | 0.51                            | 2.30          |
| M/M-CFM                                      | 0.67        | 3.31  | 24.24 | 0.51                            | 2.30          |
| M-CFM/PD                                     | 2.08        | 3.31  | 24.24 | 0.51                            | 2.30          |

Table 17–36. C364 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | fc [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | Pri3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|----------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 0.68       | 1.22       | 4.57     | 3.34     | 0.39            | 2667     | 1.22         |
| M-Mode                                    | Non-Scan      | 0.68       | 1.22       | 4.57     | 3.34     | 0.39            | 885      | 1.22         |
| B/M-Mode                                  | Scan          | 0.68       | 1.22       | 4.57     | 3.34     | 0.39            | 2667     | 1.22         |
|                                           | Non-Scan      | 0.68       | 1.22       | 4.57     | 3.34     | 0.39            | 885      | 1.22         |
| Pulsed Wave Doppler                       | Non-Scan      | 0.95       | 1.67       | 2.01     | 3.31     | 0.45            | 650      | 1.67         |
| B/Pulsed Wave Doppler                     | Scan          | 0.68       | 1.22       | 4.57     | 3.34     | 0.39            | 2667     | 1.22         |
|                                           | Non-Scan      | 0.95       | 1.67       | 2.01     | 3.31     | 0.45            | 650      | 1.67         |
| B/Color Flow Mode                         | Scan          | 0.78       | 1.41       | 1.97     | 3.31     | 1.53            | 10400    | 1.41         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.78       | 1.41       | 1.97     | 3.31     | 1.53            | 10400    | 1.41         |
|                                           | Non-Scan      | 0.95       | 1.67       | 2.01     | 3.31     | 0.45            | 650      | 1.67         |
| M/M-CFM                                   | Non-Scan      | 0.68       | 1.22       | 1.97     | 3.31     | 1.47            | 10400    | 1.22         |
| M-CFM/PD                                  | Non-Scan      | 0.68       | 1.22       | 1.97     | 3.31     | 1.47            | 10400    | 1.22         |
| CWD Mode                                  | Non-Scan      | –          | –          | –        | –        | –               | –        | –            |

Table 17–37. C364 Maximum Mechanical Index (MI)

C386

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|--------------------------------------|-------|----------|--------------------------|---------------|
| Pulsed Wave Doppler                          | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.15    | 3.4   | 43.9  | 8.9                                  | —     | 2.7      | 2.6                      | 6.0           |
| B/Pulsed Wave Doppler                        | Scan          | 0.17    | 3.8   | 28.1  | 3.0                                  | —     | 2.5      | 2.2                      | 7.8           |
|                                              | Non-Scan      | 0.15    | 3.4   | 43.9  | 8.9                                  | —     | 2.7      | 2.6                      | 7.8           |
| B/Color Flow Doppler                         | Scan          | 1.39    | 2.6   | 128.7 | 39.6                                 | —     | 2.4      | 2.0                      | 6.0           |
|                                              | Non-Scan      | —       | —     | —     | —                                    | —     | —        | —                        | —             |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.78    | 2.6   | 78.4  | 21.3                                 | —     | 2.4      | 2.0                      | 6.0           |
|                                              | Non-Scan      | <0.1    | 3.4   | 21.9  | 4.5                                  | —     | 2.7      | 2.6                      | 7.8           |
| Color M-Mode                                 | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.45    | 2.6   | 100.6 | 36.6                                 | —     | 2.4      | 2.0                      | 6.0           |
| Color Flow with<br>M Doppler                 | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.45    | 2.6   | 100.6 | 36.6                                 | —     | 2.4      | 2.0                      | 6.0           |

Table 17–38. C386 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII_{max}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| Pulsed Wave Doppler                          | Scan          | —           | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 0.93        | 3.4   | 43.9  | 7.3      | 0.4              | 2.6                      | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.17        | 3.8   | 28.1  | 7.0      | 8.5              | 2.2                      | —                  |
|                                              | Non-Scan      | 0.93        | 3.4   | 43.9  | 7.3      | 0.4              | 2.6                      | —                  |
| B/Color Flow Doppler                         | Scan          | 1.39        | 2.6   | 128.7 | 6.1      | 1.4              | 2.0                      | —                  |
|                                              | Non-Scan      | —           | —     | —     | —        | —                | —                        | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.78        | 2.6   | 78.4  | 6.1      | 1.4              | 2.0                      | —                  |
|                                              | Non-Scan      | 0.46        | 3.4   | 21.9  | 7.3      | 0.4              | 2.6                      | —                  |
| Color M-Mode                                 | Scan          | —           | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 0.97        | 2.6   | 100.6 | 6.1      | 1.5              | 2.0                      | —                  |
| Color Flow with<br>M Doppler                 | Scan          | —           | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 0.97        | 2.6   | 100.6 | 6.1      | 1.5              | 2.0                      | —                  |

Table 17–39. C386 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### C386 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| Pulsed Wave Doppler                          | 0.60        | 3.4   | 43.9  | 2.6                             | 7.8           |
| B/Pulsed Wave Doppler                        | 0.95        | 3.4   | 71.9  | 2.6                             | 7.8           |
| B/Color Flow Doppler                         | 1.93        | 2.6   | 128.7 | 2.0                             | 6.0           |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.44        | 3.4   | 100.3 | 2.6                             | 7.8           |
| Color M-Mode                                 | 1.58        | 2.6   | 100.6 | 2.0                             | 6.0           |
| Color Flow with M Doppler                    | 1.58        | 2.6   | 100.6 | 2.0                             | 6.0           |

Table 17–40. C386 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | $f_c$ [MHz] | $P_d$ [ $\mu$ sec] | PRF [Hz] | PrI3 @PII <sub>max</sub> |
|-------------------------------------------|---------------|------------|------------|----------|-------------|--------------------|----------|--------------------------|
| B–Mode                                    | Scan          | 0.62       | 1.15       | 7.05     | 3.69        | 0.35               | 2667     | 1.15                     |
| M–Mode                                    | Non–Scan      | 0.62       | 1.15       | 7.05     | 3.69        | 0.35               | 886      | 1.15                     |
| B/M–Mode                                  | Scan          | 0.62       | 1.15       | 7.05     | 3.69        | 0.35               | 2667     | 1.15                     |
|                                           | Non–Scan      | 0.62       | 1.15       | 7.05     | 3.69        | 0.35               | 886      | 1.15                     |
| Pulsed Wave Doppler                       | Non–Scan      | 0.61       | 1.19       | 7.07     | 3.42        | 0.39               | 650      | 1.19                     |
| B/Pulsed Wave Doppler                     | Scan          | 0.62       | 1.15       | 7.05     | 3.69        | 0.35               | 2667     | 1.15                     |
|                                           | Non–Scan      | 0.61       | 1.19       | 7.07     | 3.42        | 0.39               | 650      | 1.19                     |
| B/Color Flow Mode                         | Scan          | 0.85       | 1.27       | 6.05     | 2.55        | 2.03               | 420      | 1.27                     |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.85       | 1.27       | 6.05     | 2.55        | 2.03               | 420      | 1.27                     |
|                                           | Non–Scan      | 0.61       | 1.19       | 7.07     | 3.42        | 0.39               | 650      | 1.19                     |
| M/M–CFM                                   | Non–Scan      | 0.85       | 1.18       | 6.05     | 2.55        | 2.03               | 420      | 1.18                     |
| M–CFM/PD                                  | Non–Scan      | 0.85       | 1.18       | 6.05     | 2.55        | 2.03               | 420      | 1.18                     |
| CWD Mode                                  | Non–Scan      | –          | –          | –        | –           | –                  | –        | –                        |

Table 17–41. C386 Maximum Mechanical Index (MI)

C551

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$  | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|--------|--------------------------------------|-------|----------|--------------------------|---------------|
| B-Mode                                       | Scan          | 0.64    | 3.70  | 37.06  | 14.91                                | —     | 1.51     | 0.79                     | 4.20          |
| M-Mode                                       | Non-Scan      | 0.26    | 4.29  | 12.72  | 4.06                                 | —     | 1.51     | 0.79                     | 4.20          |
| B/M-Mode                                     | Scan          | 0.43    | 3.70  | 24.71  | 9.94                                 | —     | 1.51     | 0.79                     | 4.20          |
|                                              | Non-Scan      | 0.26    | 4.29  | 12.72  | 4.06                                 | —     | 1.51     | 0.79                     | 4.20          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.54    | 4.97  | 101.60 | 22.85                                | —     | 2.08     | 1.51                     | 4.20          |
| B/Pulsed Wave Doppler                        | Scan          | 0.43    | 3.70  | 24.71  | 9.94                                 | —     | 1.51     | 0.79                     | 4.20          |
|                                              | Non-Scan      | 0.54    | 4.97  | 101.60 | 22.85                                | —     | 2.08     | 1.51                     | 4.20          |
| B/Color Flow Doppler                         | Scan          | 1.25    | 4.02  | 67.42  | 22.40                                | —     | 2.08     | 1.51                     | 4.20          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.84    | 4.02  | 46.06  | 16.17                                | —     | 2.08     | 1.51                     | 4.20          |
|                                              | Non-Scan      | 0.27    | 4.97  | 50.80  | 11.43                                | —     | 2.08     | 1.51                     | 4.20          |
| M/M-CFM                                      | Non-Scan      | 0.37    | 4.02  | 32.83  | 9.83                                 | —     | 2.08     | 1.51                     | 4.20          |
| M-CFM/PD                                     | Non-Scan      | 0.37    | 4.02  | 32.83  | 9.83                                 | —     | 2.08     | 1.51                     | 4.20          |

Table 17–42. C551 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$  | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII_{max}$ |
|----------------------------------------------|---------------|-------------|-------|--------|----------|------------------|--------------------------|--------------------|
| B-Mode                                       | Scan          | 0.64        | 3.70  | 37.06  | 3.81     | 1.05             | 0.79                     | —                  |
| M-Mode                                       | Non-Scan      | 0.30        | 4.29  | 12.72  | 3.79     | 0.61             | 0.79                     | —                  |
| B/M-Mode                                     | Scan          | 0.43        | 3.70  | 24.71  | 3.81     | 1.05             | 0.79                     | —                  |
|                                              | Non-Scan      | 0.30        | 4.29  | 12.72  | 3.79     | 0.61             | 0.79                     | —                  |
| Pulsed Wave Doppler                          | Non-Scan      | 2.38        | 4.97  | 101.60 | 4.41     | 0.47             | 1.51                     | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.43        | 3.70  | 24.71  | 3.81     | 1.05             | 0.79                     | —                  |
|                                              | Non-Scan      | 2.38        | 4.97  | 101.60 | 4.41     | 0.47             | 1.51                     | —                  |
| B/Color Flow Doppler                         | Scan          | 1.25        | 4.02  | 67.42  | 4.39     | 1.24             | 1.51                     | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.84        | 4.02  | 46.06  | 4.39     | 1.24             | 1.51                     | —                  |
|                                              | Non-Scan      | 1.19        | 4.97  | 50.80  | 4.41     | 0.47             | 1.51                     | —                  |
| M/M-CFM                                      | Non-Scan      | 1.40        | 4.02  | 32.83  | 4.39     | 0.26             | 1.51                     | —                  |
| M-CFM/PD                                     | Non-Scan      | 1.40        | 4.02  | 32.83  | 4.39     | 0.26             | 1.51                     | —                  |

Table 17–43. C551 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### C551 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$  | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|--------|---------------------------------|---------------|
| B-Mode                                       | 0.92        | 3.70  | 37.06  | 0.79                            | 4.20          |
| M-Mode                                       | 0.32        | 4.29  | 12.72  | 0.79                            | 4.20          |
| B/M-Mode                                     | 0.93        | 3.70  | 37.43  | 0.79                            | 4.20          |
| Pulsed Wave Doppler                          | 1.83        | 4.97  | 101.60 | 1.51                            | 4.20          |
| B/Pulsed Wave Doppler                        | 2.45        | 4.97  | 126.31 | 1.51                            | 4.20          |
| B/Color Flow Doppler                         | 1.38        | 4.02  | 67.42  | 1.51                            | 4.20          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.91        | 4.97  | 96.87  | 1.51                            | 4.20          |
| M/M-CFM                                      | 0.68        | 4.02  | 32.83  | 1.51                            | 4.20          |
| M-CFM/PD                                     | 0.68        | 4.02  | 32.83  | 1.51                            | 4.20          |

Table 17-44. C551 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | fc [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PILmax |
|-------------------------------------------|---------------|------------|------------|----------|----------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 0.76       | 1.52       | 3.79     | 4.29     | 0.18            | 3095     | 1.52         |
| M-Mode                                    | Non-Scan      | 0.76       | 1.52       | 3.79     | 4.29     | 0.18            | 1040     | 1.52         |
| B/M-Mode                                  | Scan          | 0.76       | 1.52       | 3.79     | 4.29     | 0.18            | 3095     | 1.52         |
|                                           | Non-Scan      | 0.76       | 1.52       | 3.79     | 4.29     | 0.18            | 1040     | 1.52         |
| Pulsed Wave Doppler                       | Non-Scan      | 0.98       | 2.13       | 4.19     | 4.01     | 0.66            | 651      | 2.13         |
| B/Pulsed Wave Doppler                     | Scan          | 0.76       | 1.52       | 3.79     | 4.29     | 0.18            | 3095     | 1.52         |
|                                           | Non-Scan      | 0.98       | 2.13       | 4.19     | 4.01     | 0.66            | 651      | 2.13         |
| B/Color Flow Mode                         | Scan          | 1.18       | 2.38       | 4.16     | 3.99     | 1.28            | 280      | 2.38         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 1.18       | 2.38       | 4.16     | 3.99     | 1.28            | 280      | 2.38         |
|                                           | Non-Scan      | 0.98       | 2.13       | 4.19     | 4.01     | 0.66            | 651      | 2.13         |
| M/M-CFM                                   | Non-Scan      | 0.76       | 1.52       | 4.16     | 3.99     | 1.30            | 1480     | 1.52         |
| M-CFM/PD                                  | Non-Scan      | 0.76       | 1.52       | 4.16     | 3.99     | 1.30            | 1480     | 1.52         |

Table 17-45. C551 Maximum Mechanical Index (MI)

**C721**

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|--------------------------------------|-------|----------|--------------------------|---------------|
| Pulsed Wave Doppler                          | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.31    | 6.6   | 9.8   | 1.9                                  | —     | 1.3      | 0.6                      | 3.0           |
| B/Pulsed Wave Doppler                        | Scan          | 0.23    | 5.3   | 9.2   | 2.6                                  | —     | 1.5      | 0.8                      | 3.0           |
|                                              | Non-Scan      | 0.31    | 6.6   | 9.8   | 1.9                                  | —     | 1.3      | 0.6                      | 3.0           |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.54    | 5.1   | 22.1  | 6.8                                  | —     | 1.3      | 0.6                      | 3.0           |
|                                              | Non-Scan      | 0.15    | 6.6   | 4.9   | 0.9                                  | —     | 1.3      | 0.6                      | 3.0           |
| Color Flow with<br>M Doppler                 | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.62    | 5.1   | 25.8  | 8.4                                  | —     | 1.3      | 0.6                      | 3.0           |

Table 17–46. C721 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII_{max}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| Pulsed Wave Doppler                          | Scan          | —           | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 0.24        | 6.6   | 9.8   | 3.6      | 0.6              | 0.6                      | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.23        | 5.3   | 9.2   | 3.5      | 1.8              | 0.8                      | —                  |
|                                              | Non-Scan      | 0.24        | 6.6   | 9.8   | 3.6      | 0.6              | 0.6                      | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.54        | 5.1   | 22.1  | 3.1      | 1.1              | 0.6                      | —                  |
|                                              | Non-Scan      | 0.12        | 6.6   | 4.9   | 3.6      | 0.6              | 0.6                      | —                  |
| Color Flow with<br>M Doppler                 | Scan          | —           | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 0.52        | 5.1   | 25.8  | 3.1      | 0.7              | 0.6                      | —                  |

Table 17–47. C721 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### C721 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| Pulsed Wave Doppler                          | 0.20        | 5.1   | 8.3   | 0.9                             | 3.0           |
| B/Pulsed Wave Doppler                        | 0.43        | 5.1   | 17.5  | 0.9                             | 3.0           |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.46        | 5.1   | 19.0  | 0.9                             | 3.0           |
| Color Flow with M Doppler                    | 0.27        | 5.2   | 11.3  | 0.9                             | 3.0           |

Table 17–48. C721 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | $f_c$ [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|-------------|-----------------|----------|--------------|
| B–Mode                                    | Scan          | 0.46       | 0.90       | 3.51     | 5.46        | 0.14            | 8658     | 0.90         |
| M–Mode                                    | Non–Scan      | 0.46       | 0.60       | 3.51     | 5.46        | 0.15            | 865      | 0.60         |
| B/M–Mode                                  | Scan          | 0.46       | 0.90       | 3.51     | 5.46        | 0.14            | 8658     | 0.90         |
|                                           | Non–Scan      | 0.46       | 0.60       | 3.51     | 5.46        | 0.15            | 865      | 0.60         |
| Pulsed Wave Doppler                       | Non–Scan      | 0.42       | 1.00       | 3.09     | 5.04        | 0.53            | 650      | 1.00         |
| B/Pulsed Wave Doppler                     | Scan          | 0.46       | 0.90       | 3.51     | 5.46        | 0.14            | 8658     | 0.90         |
|                                           | Non–Scan      | 0.42       | 1.00       | 3.09     | 5.04        | 0.53            | 650      | 1.00         |
| B/Color Flow Mode                         | Scan          | 0.46       | 0.90       | 3.06     | 5.04        | 0.55            | 280      | 0.90         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.46       | 0.90       | 3.06     | 5.04        | 0.55            | 280      | 0.90         |
|                                           | Non–Scan      | 0.42       | 1.00       | 3.09     | 5.04        | 0.53            | 650      | 1.00         |
| M/M–CFM                                   | Non–Scan      | 0.46       | 0.70       | 3.06     | 5.04        | 0.55            | 280      | 0.70         |
| M–CFM/PD                                  | Non–Scan      | 0.46       | 0.70       | 3.06     | 5.04        | 0.55            | 280      | 0.70         |
| CWD Mode                                  | Non–Scan      | –          | –          | –        | –           | –               | –        | –            |

Table 17–49. C721 Maximum Mechanical Index (MI)



E721

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|--------------------------------------|-------|----------|--------------------------|---------------|
| B-Mode                                       | Scan          | 0.24    | 4.46  | 11.47 | 3.61                                 | —     | 1.54     | 0.82                     | 2.70          |
| M-Mode                                       | Non-Scan      | 0.10    | 5.71  | 3.82  | 1.05                                 | —     | 1.54     | 0.82                     | 2.70          |
| B/M-Mode                                     | Scan          | 0.16    | 4.46  | 7.64  | 2.41                                 | —     | 1.54     | 0.82                     | 2.70          |
|                                              | Non-Scan      | 0.10    | 5.71  | 3.82  | 1.05                                 | —     | 1.54     | 0.82                     | 2.70          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.67    | 4.99  | 28.56 | 8.16                                 | —     | 1.54     | 0.82                     | 2.70          |
| B/Pulsed Wave Doppler                        | Scan          | 0.16    | 4.46  | 7.64  | 2.41                                 | —     | 1.54     | 0.82                     | 2.70          |
|                                              | Non-Scan      | 0.67    | 4.99  | 28.56 | 8.16                                 | —     | 1.54     | 0.82                     | 2.70          |
| B/CFM-Mode                                   | Scan          | 0.39    | 5.05  | 16.97 | 4.86                                 | —     | 1.49     | 0.77                     | 2.60          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.27    | 5.05  | 12.31 | 3.63                                 | —     | 1.49     | 0.77                     | 2.60          |
|                                              | Non-Scan      | 0.33    | 4.99  | 14.28 | 4.08                                 | —     | 1.54     | 0.82                     | 2.70          |
| M/M-CFM                                      | Non-Scan      | 0.73    | 5.05  | 30.09 | 7.95                                 | —     | 1.49     | 0.77                     | 2.60          |
| M-CFM/PD                                     | Non-Scan      | 0.73    | 5.05  | 30.09 | 7.95                                 | —     | 1.49     | 0.77                     | 2.60          |

Table 17–50. E721 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PIL_{max}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| B-Mode                                       | Scan          | 0.24        | 4.46  | 11.47 | 3.70     | 0.87             | 0.82                     | —                  |
| M-Mode                                       | Non-Scan      | 0.12        | 5.71  | 3.82  | 3.44     | 0.53             | 0.82                     | —                  |
| B/M-Mode                                     | Scan          | 0.16        | 4.46  | 7.64  | 3.70     | 0.87             | 0.82                     | —                  |
|                                              | Non-Scan      | 0.12        | 5.71  | 3.82  | 3.44     | 0.53             | 0.82                     | —                  |
| Pulsed Wave Doppler                          | Non-Scan      | 1.42        | 4.99  | 28.56 | 3.75     | 0.27             | 0.82                     | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.16        | 4.46  | 7.64  | 3.70     | 0.87             | 0.82                     | —                  |
|                                              | Non-Scan      | 1.42        | 4.99  | 28.56 | 3.75     | 0.27             | 0.82                     | —                  |
| B-CFM                                        | Scan          | 0.39        | 5.05  | 16.97 | 3.74     | 1.26             | 0.77                     | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.27        | 5.05  | 12.31 | 3.74     | 1.26             | 0.77                     | —                  |
|                                              | Non-Scan      | 0.71        | 4.99  | 14.28 | 3.75     | 0.27             | 0.82                     | —                  |
| M/M-CFM                                      | Non-Scan      | 0.68        | 5.05  | 30.09 | 3.74     | 0.59             | 0.77                     | —                  |
| M-CFM/PD                                     | Non-Scan      | 0.68        | 5.05  | 30.09 | 3.74     | 0.59             | 0.77                     | —                  |

Table 17–51. E721 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### E721 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| B-Mode                                       | 0.28        | 4.46  | 11.47 | 0.82                            | 2.70          |
| M-Mode                                       | < 0.1       | 5.71  | 3.82  | 0.82                            | 2.70          |
| B/M-Mode                                     | 0.28        | 4.46  | 11.46 | 0.82                            | 2.70          |
| Pulsed Wave Doppler                          | 0.70        | 4.99  | 28.56 | 0.82                            | 2.70          |
| B/Pulsed Wave Doppler                        | 0.88        | 4.99  | 36.21 | 0.82                            | 2.70          |
| B-CFM                                        | 0.42        | 5.05  | 16.97 | 0.77                            | 2.60          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.65        | 4.99  | 26.59 | 0.82                            | 2.70          |
| M/M-CFM                                      | 0.76        | 5.05  | 30.09 | 0.77                            | 2.60          |
| Color Flow with M Doppler                    | 0.76        | 5.05  | 30.09 | 0.77                            | 2.60          |

Table 17–52. E721 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | $f_c$ [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|-------------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 0.81       | 1.83       | 3.44     | 5.71        | 0.15            | 4075     | 1.83         |
| M-Mode                                    | Non-Scan      | 0.81       | 1.83       | 3.44     | 5.71        | 0.15            | 866      | 1.83         |
| B/M-Mode                                  | Scan          | 0.81       | 1.83       | 3.44     | 5.71        | 0.15            | 4075     | 1.83         |
|                                           | Non-Scan      | 0.81       | 1.83       | 3.44     | 5.71        | 0.15            | 866      | 1.83         |
| Pulsed Wave Doppler                       | Non-Scan      | 0.74       | 1.79       | 2.58     | 6.43        | 0.59            | 651      | 1.79         |
| B/Pulsed Wave Doppler                     | Scan          | 0.81       | 1.83       | 3.44     | 5.71        | 0.15            | 4075     | 1.83         |
|                                           | Non-Scan      | 0.74       | 1.79       | 2.58     | 6.43        | 0.59            | 651      | 1.79         |
| B/Color Flow Mode                         | Scan          | 0.81       | 1.83       | 2.63     | 6.57        | 0.77            | 280      | 1.83         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.81       | 1.83       | 2.63     | 6.57        | 0.77            | 280      | 1.83         |
|                                           | Non-Scan      | 0.74       | 1.79       | 2.58     | 6.43        | 0.59            | 651      | 1.79         |
| M/M-CFM                                   | Non-Scan      | 0.81       | 1.83       | 2.63     | 6.57        | 0.77            | 1480     | 1.83         |
| M-CFM/PD                                  | Non-Scan      | 0.81       | 1.83       | 2.63     | 6.57        | 0.77            | 1480     | 1.83         |

Table 17–53. E721 Maximum Mechanical Index (MI)

**ERB7 (Convex)**

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_{3.3}(z_1)$<br>$I_{TA,3.3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|--------------------------------------------|-------|----------|--------------------------|---------------|
| B-Mode                                       | Scan          | 0.22    | 4.49  | 10.16 | 5.78                                       | —     | 1.15     | 0.46                     | 1.90          |
| M-Mode                                       | Non-Scan      | 0.13    | 5.78  | 4.76  | 2.29                                       | —     | 1.15     | 0.46                     | 1.90          |
| B/M-Mode                                     | Scan          | 0.14    | 4.49  | 6.77  | 3.85                                       | —     | 1.15     | 0.46                     | 1.90          |
|                                              | Non-Scan      | 0.13    | 5.78  | 4.76  | 2.29                                       | —     | 1.15     | 0.46                     | 1.90          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.17    | 4.98  | 7.35  | 3.20                                       | —     | 1.15     | 0.46                     | 2.15          |
| B/Pulsed Wave Doppler                        | Scan          | 0.14    | 4.49  | 6.77  | 3.85                                       | —     | 1.15     | 0.46                     | 1.90          |
|                                              | Non-Scan      | 0.17    | 4.98  | 7.35  | 3.20                                       | —     | 1.15     | 0.46                     | 2.15          |
| B/CFM-Mode                                   | Scan          | 0.22    | 5.02  | 9.73  | 5.02                                       | —     | 1.15     | 0.46                     | 2.15          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.18    | 5.02  | 8.25  | 4.44                                       | —     | 1.15     | 0.46                     | 2.15          |
|                                              | Non-Scan      | < 0.1   | 4.98  | 3.68  | 1.60                                       | —     | 1.15     | 0.46                     | 2.15          |
| M/M-CFM                                      | Non-Scan      | 0.10    | 5.02  | 7.97  | 3.69                                       | —     | 1.15     | 0.46                     | 2.15          |
| M-CFM/PD                                     | Non-Scan      | 0.17    | 5.02  | 7.97  | 3.69                                       | —     | 1.15     | 0.46                     | 2.15          |

Table 17–54. ERB7 (Convex) Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PIL_{max}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| B-Mode                                       | Scan          | 0.22        | 4.49  | 10.16 | 1.95     | 0.79             | 0.46                     | —                  |
| M-Mode                                       | Non-Scan      | 0.18        | 5.78  | 4.76  | 1.87     | 0.47             | 0.46                     | —                  |
| B/M-Mode                                     | Scan          | 0.14        | 4.49  | 6.77  | 1.95     | 0.79             | 0.46                     | —                  |
|                                              | Non-Scan      | 0.18        | 5.78  | 4.76  | 1.87     | 0.47             | 0.46                     | —                  |
| Pulsed Wave Doppler                          | Non-Scan      | 0.54        | 4.98  | 7.35  | 2.28     | 0.21             | 0.46                     | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.14        | 4.49  | 6.77  | 1.95     | 0.79             | 0.46                     | —                  |
|                                              | Non-Scan      | 0.54        | 4.98  | 7.35  | 2.28     | 0.21             | 0.46                     | —                  |
| B-CFM                                        | Scan          | 0.22        | 5.02  | 9.73  | 2.20     | 1.33             | 0.46                     | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.18        | 5.02  | 8.25  | 2.20     | 1.33             | 0.46                     | —                  |
|                                              | Non-Scan      | 0.27        | 4.98  | 3.68  | 2.28     | 0.21             | 0.46                     | —                  |
| M/M-CFM                                      | Non-Scan      | 0.32        | 5.02  | 7.97  | 2.20     | 0.21             | 0.46                     | —                  |
| M-CFM/PD                                     | Non-Scan      | 0.54        | 5.02  | 7.97  | 2.20     | 0.21             | 0.46                     | —                  |

Table 17–55. ERB7 (Convex) Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### ERB7 (Convex) (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| B-Mode                                       | 0.33        | 4.49  | 10.16 | 0.46                            | 1.90          |
| M-Mode                                       | 0.16        | 5.78  | 4.76  | 0.46                            | 1.90          |
| B/M-Mode                                     | 0.38        | 4.49  | 11.53 | 0.46                            | 1.90          |
| Pulsed Wave Doppler                          | 0.24        | 4.98  | 7.35  | 0.46                            | 2.15          |
| B/Pulsed Wave Doppler                        | 0.46        | 4.98  | 14.13 | 0.46                            | 2.15          |
| B-CFM                                        | 0.32        | 5.02  | 9.73  | 0.46                            | 2.15          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.39        | 4.98  | 11.93 | 0.46                            | 2.15          |
| M/M-CFM                                      | 0.13        | 5.02  | 7.97  | 0.46                            | 2.15          |
| Color Flow with M Doppler                    | 0.24        | 5.02  | 7.97  | 0.46                            | 2.15          |

Table 17–56. ERB7 (Convex) Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | $f_c$ [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|-------------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 1.01       | 2.20       | 1.87     | 5.78        | 0.13            | 4075     | 2.20         |
| M-Mode                                    | Non-Scan      | 1.01       | 2.20       | 1.87     | 5.78        | 0.13            | 816      | 2.20         |
| B/M-Mode                                  | Scan          | 1.01       | 2.20       | 1.87     | 5.78        | 0.13            | 4075     | 2.20         |
|                                           | Non-Scan      | 1.01       | 2.20       | 1.87     | 5.78        | 0.13            | 816      | 2.20         |
| Pulsed Wave Doppler                       | Non-Scan      | 1.03       | 2.22       | 2.10     | 5.03        | 0.52            | 651      | 2.22         |
| B/Pulsed Wave Doppler                     | Scan          | 1.01       | 2.20       | 1.87     | 5.78        | 0.13            | 4075     | 2.20         |
|                                           | Non-Scan      | 1.03       | 2.22       | 2.10     | 5.03        | 0.52            | 651      | 2.22         |
| B/Color Flow Mode                         | Scan          | 1.01       | 2.20       | 2.64     | 5.03        | 0.99            | 280      | 2.20         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 1.01       | 2.20       | 2.64     | 5.03        | 0.99            | 280      | 2.20         |
|                                           | Non-Scan      | 1.03       | 2.22       | 2.10     | 5.03        | 0.52            | 651      | 2.22         |
| M/M-CFM                                   | Non-Scan      | 1.01       | 2.20       | 2.64     | 5.03        | 0.99            | 1480     | 2.20         |
| M-CFM/PD                                  | Non-Scan      | 1.01       | 2.20       | 2.64     | 5.03        | 0.99            | 1480     | 2.20         |

Table 17–57. ERB7 (Convex) Maximum Mechanical Index (MI)

ERB7 (Linear)

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_{3.3}(z_1)$<br>$I_{TA.3.3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|--------------------------------------------|-------|----------|--------------------------|---------------|
| B-Mode                                       | Scan          | 0.15    | 4.52  | 6.90  | 3.63                                       | —     | 1.07     | 0.40                     | 1.90          |
| M-Mode                                       | Non-Scan      | 0.10    | 6.06  | 3.59  | 1.10                                       | —     | 1.33     | 0.62                     | 2.90          |
| B/M-Mode                                     | Scan          | 0.10    | 4.52  | 4.60  | 2.42                                       | —     | 1.07     | 0.40                     | 1.90          |
|                                              | Non-Scan      | 0.10    | 6.06  | 3.59  | 1.10                                       | —     | 1.33     | 0.62                     | 2.90          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.29    | 6.62  | 9.27  | 3.54                                       | —     | 1.07     | 0.40                     | 1.90          |
| B/Pulsed Wave Doppler                        | Scan          | 0.10    | 4.52  | 4.60  | 2.42                                       | —     | 1.07     | 0.40                     | 1.90          |
|                                              | Non-Scan      | 0.29    | 6.62  | 9.27  | 3.54                                       | —     | 1.07     | 0.40                     | 1.90          |
| B/CFM-Mode                                   | Scan          | 0.95    | 5.07  | 39.69 | 19.29                                      | —     | 1.07     | 0.40                     | 1.90          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.53    | 5.07  | 22.15 | 10.86                                      | —     | 1.07     | 0.40                     | 1.90          |
|                                              | Non-Scan      | 0.15    | 6.62  | 4.64  | 1.77                                       | —     | 1.07     | 0.40                     | 1.90          |
| M/M-CFM                                      | Non-Scan      | 0.51    | 5.07  | 24.01 | 11.11                                      | —     | 1.07     | 0.40                     | 1.90          |
| M-CFM/PD                                     | Non-Scan      | 0.29    | 5.07  | 24.01 | 11.11                                      | —     | 1.07     | 0.40                     | 1.90          |

Table 17–58. ERB7 (Linear) Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PIL_{max}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| B-Mode                                       | Scan          | 0.15        | 4.52  | 6.90  | 2.00     | 0.86             | 0.40                     | —                  |
| M-Mode                                       | Non-Scan      | 0.12        | 6.06  | 3.59  | 2.74     | 0.48             | 0.62                     | —                  |
| B/M-Mode                                     | Scan          | 0.10        | 4.52  | 4.60  | 2.00     | 0.86             | 0.40                     | —                  |
|                                              | Non-Scan      | 0.12        | 6.06  | 3.59  | 2.74     | 0.48             | 0.62                     | —                  |
| Pulsed Wave Doppler                          | Non-Scan      | 0.80        | 6.62  | 9.27  | 2.08     | 0.17             | 0.40                     | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.10        | 4.52  | 4.60  | 2.00     | 0.86             | 0.40                     | —                  |
|                                              | Non-Scan      | 0.80        | 6.62  | 9.27  | 2.08     | 0.17             | 0.40                     | —                  |
| B-CFM                                        | Scan          | 0.95        | 5.07  | 39.69 | 1.96     | 0.92             | 0.40                     | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.53        | 5.07  | 22.15 | 1.96     | 0.92             | 0.40                     | —                  |
|                                              | Non-Scan      | 0.40        | 6.62  | 4.64  | 2.08     | 0.17             | 0.40                     | —                  |
| M/M-CFM                                      | Non-Scan      | 1.41        | 5.07  | 24.01 | 1.96     | 0.28             | 0.40                     | —                  |
| M-CFM/PD                                     | Non-Scan      | 0.80        | 5.07  | 24.01 | 1.96     | 0.28             | 0.40                     | —                  |

Table 17–59. ERB7 (Linear) Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### ERB7 (Linear) (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| B-Mode                                       | 0.24        | 4.52  | 6.90  | 0.40                            | 1.90          |
| M-Mode                                       | 0.10        | 6.06  | 3.59  | 0.62                            | 2.90          |
| B/M-Mode                                     | 0.26        | 4.52  | 8.19  | 0.40                            | 1.90          |
| Pulsed Wave Doppler                          | 0.33        | 6.62  | 9.27  | 0.40                            | 1.90          |
| B/Pulsed Wave Doppler                        | 0.49        | 6.62  | 13.88 | 0.40                            | 1.90          |
| B-CFM                                        | 1.40        | 5.07  | 39.69 | 0.40                            | 1.90          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.94        | 6.62  | 26.78 | 0.40                            | 1.90          |
| M/M-CFM                                      | 0.73        | 5.07  | 24.01 | 0.40                            | 1.90          |
| Color Flow with M Doppler                    | 0.33        | 5.07  | 24.01 | 0.40                            | 1.90          |

Table 17–60. ERB7 (Linear) Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | $f_c$ [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|-------------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 0.70       | 1.81       | 2.74     | 6.06        | 0.15            | 4075     | 1.81         |
| M-Mode                                    | Non-Scan      | 0.70       | 1.81       | 2.74     | 6.06        | 0.15            | 820      | 1.81         |
| B/M-Mode                                  | Scan          | 0.70       | 1.81       | 2.74     | 6.06        | 0.15            | 4075     | 1.81         |
|                                           | Non-Scan      | 0.70       | 1.81       | 2.74     | 6.06        | 0.15            | 820      | 1.81         |
| Pulsed Wave Doppler                       | Non-Scan      | 0.95       | 2.19       | 1.99     | 5.19        | 0.52            | 651      | 2.19         |
| B/Pulsed Wave Doppler                     | Scan          | 0.70       | 1.81       | 2.74     | 6.06        | 0.15            | 4075     | 1.81         |
|                                           | Non-Scan      | 0.95       | 2.19       | 1.99     | 5.19        | 0.52            | 651      | 2.19         |
| B/Color Flow Mode                         | Scan          | 0.91       | 2.03       | 1.96     | 5.07        | 0.99            | 1720     | 2.03         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.91       | 2.03       | 1.96     | 5.07        | 0.99            | 1720     | 2.03         |
|                                           | Non-Scan      | 0.95       | 2.19       | 1.99     | 5.19        | 0.52            | 651      | 2.19         |
| M/M-CFM                                   | Non-Scan      | 0.70       | 1.81       | 1.96     | 5.07        | 0.99            | 1700     | 1.81         |
| M-CFM/PD                                  | Non-Scan      | 0.70       | 1.81       | 1.96     | 5.07        | 0.99            | 1700     | 1.81         |

Table 17–61. ERB7 (Linear) Maximum Mechanical Index (MI)

I739

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|--------------------------------------|-------|----------|--------------------------|---------------|
| Pulsed Wave Doppler                          | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.36    | 5.1   | 15.0  | 7.0                                  | —     | 0.8      | 0.2                      | 1.6           |
| B/Pulsed Wave Doppler                        | Scan          | 0.11    | 7.0   | 3.2   | 0.8                                  | —     | 1.1      | 0.4                      | 2.1           |
|                                              | Non-Scan      | 0.36    | 5.1   | 15.0  | 7.0                                  | —     | 0.8      | 0.2                      | 1.6           |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.40    | 5.1   | 13.9  | 9.9                                  | —     | 0.8      | 0.2                      | 1.6           |
|                                              | Non-Scan      | 0.18    | 5.1   | 7.5   | 3.5                                  | —     | 0.8      | 0.2                      | 1.6           |
| Color Flow with<br>M Doppler                 | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.59    | 5.1   | 21.4  | 18.1                                 | —     | 0.8      | 0.2                      | 1.6           |

Table 17–62. I739 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII_{max}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| Pulsed Wave Doppler                          | Scan          | —           | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 0.96        | 5.1   | 15.0  | 1.7      | 0.3              | 0.2                      | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.11        | 7.0   | 3.2   | 2.1      | 0.8              | 0.4                      | —                  |
|                                              | Non-Scan      | 0.96        | 5.1   | 15.0  | 1.7      | 0.3              | 0.2                      | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.40        | 5.1   | 13.9  | 1.8      | 0.4              | 0.2                      | —                  |
|                                              | Non-Scan      | 0.48        | 5.1   | 7.5   | 1.7      | 0.3              | 0.2                      | —                  |
| Color Flow with<br>M Doppler                 | Scan          | —           | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 0.89        | 5.1   | 21.4  | 1.8      | 0.5              | 0.2                      | —                  |

Table 17–63. I739 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### I739 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| Pulsed Wave Doppler                          | 0.71        | 5.1   | 15.0  | 0.2                             | 1.6           |
| B/Pulsed Wave Doppler                        | 0.82        | 5.1   | 18.3  | 0.2                             | 1.6           |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.83        | 5.1   | 21.5  | 0.2                             | 1.6           |
| Color Flow with M Doppler                    | 0.74        | 5.1   | 21.4  | 0.2                             | 1.6           |

Table 17–64. I739 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | $f_c$ [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PII <sub>max</sub> |
|-------------------------------------------|---------------|------------|------------|----------|-------------|-----------------|----------|--------------------------|
| B–Mode                                    | Scan          | 0.68       | 1.73       | 2.09     | 6.99        | 0.20            | 8658     | 1.73                     |
| M–Mode                                    | Non–Scan      | 0.68       | 1.73       | 2.09     | 6.99        | 0.20            | 866      | 1.73                     |
| B/M–Mode                                  | Scan          | 0.68       | 1.73       | 2.09     | 6.99        | 0.20            | 8658     | 1.73                     |
|                                           | Non–Scan      | 0.68       | 1.73       | 2.09     | 6.99        | 0.20            | 866      | 1.73                     |
| Pulsed Wave Doppler                       | Non–Scan      | 0.67       | 1.65       | 2.07     | 6.64        | 2.21            | 650      | 1.65                     |
| B/Pulsed Wave Doppler                     | Scan          | 0.68       | 1.73       | 2.09     | 6.99        | 0.20            | 8658     | 1.73                     |
|                                           | Non–Scan      | 0.67       | 1.65       | 2.07     | 6.64        | 2.21            | 650      | 1.65                     |
| B/Color Flow Mode                         | Scan          | 0.82       | 1.78       | 1.76     | 5.11        | 0.53            | 280      | 1.78                     |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.82       | 1.78       | 1.76     | 5.11        | 0.53            | 280      | 1.78                     |
|                                           | Non–Scan      | 0.67       | 1.65       | 2.07     | 6.64        | 2.21            | 650      | 1.65                     |
| M/M–CFM                                   | Non–Scan      | 0.93       | 1.73       | 1.84     | 5.08        | 0.52            | 280      | 1.73                     |
| M–CFM/PD                                  | Non–Scan      | 0.93       | 1.73       | 1.84     | 5.08        | 0.52            | 280      | 1.73                     |
| CWD Mode                                  | Non–Scan      | –          | –          | –        | –           | –               | –        | –                        |

Table 17–65. I739 Maximum Mechanical Index (MI)



546L

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$  | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|--------|--------------------------------------|-------|----------|--------------------------|---------------|
| B-Mode                                       | Scan          | 0.94    | 4.45  | 45.09  | 19.70                                | —     | 1.28     | 0.58                     | 2.60          |
| M-Mode                                       | Non-Scan      | 0.11    | 4.42  | 5.37   | 2.58                                 | —     | 1.28     | 0.58                     | 2.60          |
| B/M-Mode                                     | Scan          | 0.63    | 4.45  | 30.06  | 13.13                                | —     | 1.28     | 0.58                     | 2.60          |
|                                              | Non-Scan      | 0.11    | 4.42  | 5.37   | 2.58                                 | —     | 1.28     | 0.58                     | 2.60          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.39    | 3.99  | 20.54  | 8.14                                 | —     | 1.44     | 0.72                     | 2.80          |
| B/Pulsed Wave Doppler                        | Scan          | 0.63    | 4.45  | 30.06  | 13.13                                | —     | 1.28     | 0.58                     | 2.60          |
|                                              | Non-Scan      | 0.39    | 3.99  | 20.54  | 8.14                                 | —     | 1.44     | 0.72                     | 2.80          |
| B/CFM-Mode                                   | Scan          | 2.88    | 4.04  | 147.33 | 63.96                                | —     | 1.44     | 0.72                     | 2.80          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 1.75    | 4.04  | 88.69  | 38.55                                | —     | 1.44     | 0.72                     | 2.80          |
|                                              | Non-Scan      | 0.20    | 3.99  | 10.27  | 4.07                                 | —     | 1.44     | 0.72                     | 2.80          |
| M/M-CFM                                      | Non-Scan      | 0.52    | 4.04  | 26.77  | 11.86                                | —     | 1.44     | 0.72                     | 2.80          |
| M-CFM/PD                                     | Non-Scan      | 0.52    | 4.04  | 26.77  | 11.86                                | —     | 1.44     | 0.72                     | 2.80          |

Table 17–66. 546L Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$  | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PIL_{max}$ |
|----------------------------------------------|---------------|-------------|-------|--------|----------|------------------|--------------------------|--------------------|
| B-Mode                                       | Scan          | 0.94        | 4.45  | 45.09  | 2.66     | 0.81             | 0.58                     | —                  |
| M-Mode                                       | Non-Scan      | 0.38        | 4.42  | 5.37   | 2.63     | 0.25             | 0.58                     | —                  |
| B/M-Mode                                     | Scan          | 0.63        | 4.45  | 30.06  | 2.66     | 0.81             | 0.58                     | —                  |
|                                              | Non-Scan      | 0.38        | 4.42  | 5.37   | 2.63     | 0.25             | 0.58                     | —                  |
| Pulsed Wave Doppler                          | Non-Scan      | 1.25        | 3.99  | 20.54  | 3.29     | 0.26             | 0.72                     | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.63        | 4.45  | 30.06  | 2.66     | 0.81             | 0.58                     | —                  |
|                                              | Non-Scan      | 1.25        | 3.99  | 20.54  | 3.29     | 0.26             | 0.72                     | —                  |
| B-CFM                                        | Scan          | 2.88        | 4.04  | 147.33 | 3.03     | 0.74             | 0.72                     | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 1.75        | 4.04  | 88.69  | 3.03     | 0.74             | 0.72                     | —                  |
|                                              | Non-Scan      | 0.62        | 3.99  | 10.27  | 3.29     | 0.26             | 0.72                     | —                  |
| M/M-CFM                                      | Non-Scan      | 1.64        | 4.04  | 26.77  | 3.03     | 0.26             | 0.72                     | —                  |
| M-CFM/PD                                     | Non-Scan      | 1.64        | 4.04  | 26.77  | 3.03     | 0.26             | 0.72                     | —                  |

Table 17–67. 546L Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### 546L (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$  | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|--------|---------------------------------|---------------|
| B-Mode                                       | 1.32        | 4.45  | 45.09  | 0.58                            | 2.60          |
| M-Mode                                       | 0.16        | 4.42  | 5.37   | 0.58                            | 2.60          |
| B/M-Mode                                     | 1.03        | 4.45  | 35.42  | 0.58                            | 2.60          |
| Pulsed Wave Doppler                          | 0.54        | 3.99  | 20.54  | 0.72                            | 2.80          |
| B/Pulsed Wave Doppler                        | 1.41        | 3.99  | 50.60  | 0.72                            | 2.80          |
| B-CFM                                        | 3.94        | 4.04  | 147.33 | 0.72                            | 2.80          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 2.68        | 3.99  | 98.96  | 0.72                            | 2.80          |
| M/M-CFM                                      | 0.72        | 4.04  | 26.77  | 0.72                            | 2.80          |
| Color Flow with M Doppler                    | 0.72        | 4.04  | 26.77  | 0.72                            | 2.80          |

Table 17-68. 546L Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | $f_c$ [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|-------------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 1.13       | 2.30       | 2.63     | 4.42        | 0.21            | 3687     | 2.30         |
| M-Mode                                    | Non-Scan      | 1.13       | 2.30       | 2.63     | 4.42        | 0.21            | 914      | 2.30         |
| B/M-Mode                                  | Scan          | 1.13       | 2.30       | 2.63     | 4.42        | 0.21            | 3687     | 2.30         |
|                                           | Non-Scan      | 1.13       | 2.30       | 2.63     | 4.42        | 0.21            | 914      | 2.30         |
| Pulsed Wave Doppler                       | Non-Scan      | 1.26       | 2.57       | 2.98     | 4.06        | 0.65            | 651      | 2.57         |
| B/Pulsed Wave Doppler                     | Scan          | 1.13       | 2.30       | 2.63     | 4.42        | 0.21            | 3687     | 2.30         |
|                                           | Non-Scan      | 1.26       | 2.57       | 2.98     | 4.06        | 0.65            | 651      | 2.57         |
| B/Color Flow Mode                         | Scan          | 1.42       | 2.93       | 2.98     | 4.02        | 1.24            | 280      | 2.93         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 1.42       | 2.93       | 2.98     | 4.02        | 1.24            | 280      | 2.93         |
|                                           | Non-Scan      | 1.26       | 2.57       | 2.98     | 4.06        | 0.65            | 651      | 2.57         |
| M/M-CFM                                   | Non-Scan      | 1.13       | 2.30       | 2.98     | 4.02        | 1.22            | 1480     | 2.30         |
| M-CFM/PD                                  | Non-Scan      | 1.13       | 2.30       | 2.98     | 4.02        | 1.22            | 1480     | 2.30         |

Table 17-69. 546L Maximum Mechanical Index (MI)

739L

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$  | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|--------|--------------------------------------|-------|----------|--------------------------|---------------|
| B-Mode                                       | Scan          | 0.91    | 6.18  | 31.07  | 14.71                                | —     | 1.12     | 0.44                     | 1.90          |
| M-Mode                                       | Non-Scan      | < 0.1   | 6.18  | 3.11   | 1.47                                 | —     | 1.12     | 0.44                     | 1.90          |
| B/M-Mode                                     | Scan          | 0.61    | 6.18  | 20.71  | 9.81                                 | —     | 1.12     | 0.44                     | 1.90          |
|                                              | Non-Scan      | < 0.1   | 6.18  | 3.11   | 1.47                                 | —     | 1.12     | 0.44                     | 1.90          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.36    | 5.01  | 14.96  | 10.22                                | —     | 0.88     | 0.27                     | 1.00          |
| B/Pulsed Wave Doppler                        | Scan          | 0.61    | 6.18  | 20.71  | 9.81                                 | —     | 1.12     | 0.44                     | 1.90          |
|                                              | Non-Scan      | 0.36    | 5.01  | 14.96  | 10.22                                | —     | 0.88     | 0.27                     | 1.00          |
| B/CFM-Mode                                   | Scan          | 3.16    | 5.02  | 127.39 | 59.47                                | —     | 1.27     | 0.56                     | 2.10          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 1.88    | 5.02  | 74.05  | 34.64                                | —     | 1.27     | 0.56                     | 2.10          |
|                                              | Non-Scan      | 0.18    | 5.01  | 7.48   | 5.11                                 | —     | 0.88     | 0.27                     | 1.00          |
| M/M-CFM                                      | Non-Scan      | 1.04    | 5.02  | 42.78  | 19.93                                | —     | 1.27     | 0.56                     | 2.10          |
| M-CFM/PD                                     | Non-Scan      | 1.04    | 5.02  | 42.78  | 19.93                                | —     | 1.27     | 0.56                     | 2.10          |

Table 17–70. 739L Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$  | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PIL_{max}$ |
|----------------------------------------------|---------------|-------------|-------|--------|----------|------------------|--------------------------|--------------------|
| B-Mode                                       | Scan          | 0.91        | 6.18  | 31.07  | 1.91     | 0.69             | 0.44                     | —                  |
| M-Mode                                       | Non-Scan      | 0.28        | 6.18  | 3.11   | 1.91     | 0.18             | 0.44                     | —                  |
| B/M-Mode                                     | Scan          | 0.61        | 6.18  | 20.71  | 1.91     | 0.69             | 0.44                     | —                  |
|                                              | Non-Scan      | 0.28        | 6.18  | 3.11   | 1.91     | 0.18             | 0.44                     | —                  |
| Pulsed Wave Doppler                          | Non-Scan      | 1.53        | 5.01  | 14.96  | 1.10     | 0.20             | 0.27                     | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.61        | 6.18  | 20.71  | 1.91     | 0.69             | 0.44                     | —                  |
|                                              | Non-Scan      | 1.53        | 5.01  | 14.96  | 1.10     | 0.20             | 0.27                     | —                  |
| B-CFM                                        | Scan          | 3.16        | 5.02  | 127.39 | 2.17     | 0.63             | 0.56                     | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 1.88        | 5.02  | 74.05  | 2.17     | 0.63             | 0.56                     | —                  |
|                                              | Non-Scan      | 0.77        | 5.01  | 7.48   | 1.10     | 0.20             | 0.27                     | —                  |
| M/M-CFM                                      | Non-Scan      | 2.00        | 5.02  | 42.78  | 2.17     | 0.36             | 0.56                     | —                  |
| M-CFM/PD                                     | Non-Scan      | 2.00        | 5.02  | 42.78  | 2.17     | 0.36             | 0.56                     | —                  |

Table 17–71. 739L Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### 739L (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$  | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|--------|---------------------------------|---------------|
| B-Mode                                       | 1.04        | 6.18  | 31.07  | 0.44                            | 1.90          |
| M-Mode                                       | 0.10        | 6.18  | 3.11   | 0.44                            | 1.90          |
| B/M-Mode                                     | 0.80        | 6.18  | 23.82  | 0.44                            | 1.90          |
| Pulsed Wave Doppler                          | 0.64        | 5.01  | 14.96  | 0.27                            | 1.00          |
| B/Pulsed Wave Doppler                        | 1.33        | 5.01  | 35.67  | 0.27                            | 1.00          |
| B-CFM                                        | 3.85        | 5.02  | 127.39 | 0.56                            | 2.10          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 2.59        | 5.01  | 81.53  | 0.27                            | 1.00          |
| M/M-CFM                                      | 1.28        | 5.02  | 42.78  | 0.56                            | 2.10          |
| Color Flow with M Doppler                    | 1.28        | 5.02  | 42.78  | 0.56                            | 2.10          |

Table 17-72. 739L Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | $f_c$ [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|-------------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 1.39       | 3.29       | 1.91     | 6.18        | 0.14            | 8658     | 3.29         |
| M-Mode                                    | Non-Scan      | 1.39       | 3.29       | 1.91     | 6.18        | 0.14            | 866      | 3.29         |
| B/M-Mode                                  | Scan          | 1.39       | 3.29       | 1.91     | 6.18        | 0.14            | 8658     | 3.29         |
|                                           | Non-Scan      | 1.39       | 3.29       | 1.91     | 6.18        | 0.14            | 866      | 3.29         |
| Pulsed Wave Doppler                       | Non-Scan      | 1.67       | 3.53       | 1.08     | 5.09        | 0.52            | 651      | 3.53         |
| B/Pulsed Wave Doppler                     | Scan          | 1.39       | 3.29       | 1.91     | 6.18        | 0.14            | 8658     | 3.29         |
|                                           | Non-Scan      | 1.67       | 3.53       | 1.08     | 5.09        | 0.52            | 651      | 3.53         |
| B/Color Flow Mode                         | Scan          | 1.39       | 3.29       | 2.17     | 5.02        | 1.01            | 1720     | 3.29         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 1.39       | 3.29       | 2.17     | 5.02        | 1.01            | 1720     | 3.29         |
|                                           | Non-Scan      | 1.67       | 3.53       | 1.08     | 5.09        | 0.52            | 651      | 3.53         |
| M/M-CFM                                   | Non-Scan      | 1.39       | 3.29       | 2.17     | 5.02        | 1.01            | 1700     | 3.29         |
| M-CFM/PD                                  | Non-Scan      | 1.39       | 3.29       | 2.17     | 5.02        | 1.01            | 1700     | 3.29         |

Table 17-73. 739L Maximum Mechanical Index (MI)

L764

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_3(z_1)$<br>$I_{TA.3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|--------------------------------------|-------|----------|--------------------------|---------------|
| Pulsed Wave Doppler                          | Scan          | NA      | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.2     | 5.16  | 5.9   | —                                    | —     | —        | 0.8x0.7                  | 2.88          |
| B/Pulsed Wave Doppler                        | Scan          | 0.2     | 6.12  | 6.1   | —                                    | —     | —        | 0.9x0.7                  | 3.48          |
|                                              | Non-Scan      | 0.2     | 5.16  | 5.9   | —                                    | —     | —        | 0.8x0.7                  | 2.88          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.3     | 6.12  | 10.0  | —                                    | —     | —        | 0.9x0.7                  | 3.48          |
|                                              | Non-Scan      | 0.2     | 5.16  | 5.9   | —                                    | —     | —        | 0.8x0.7                  | 2.88          |
| Color Flow with<br>M Doppler                 | Scan          | 0.3     | 6.45  | 10.1  | —                                    | —     | —        | 0.8x0.7                  | 2.88          |
|                                              | Non-Scan      | < 0.1   | 6.45  | 0.9   | —                                    | —     | —        | 0.8x0.7                  | 2.72          |

Table 17–74. L764 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII_{max}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| Pulsed Wave Doppler                          | Scan          | NA          | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 0.5         | 5.17  | 5.9   | 2.5      | 0.1              | 0.8x0.7                  | 0.1                |
| B/Pulsed Wave Doppler                        | Scan          | 0.1         | 6.10  | 6.1   | —        | —                | 0.9x0.7                  | —                  |
|                                              | Non-Scan      | 0.5         | 5.17  | 5.9   | 2.5      | 0.1              | 0.8x0.7                  | 0.1                |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.3         | 6.10  | 17.0  | —        | —                | 0.9x0.7                  | —                  |
|                                              | Non-Scan      | 0.5         | 5.17  | 5.9   | 2.5      | 0.1              | 0.8x0.7                  | 0.1                |
| Color Flow with<br>M Doppler                 | Scan          | 0.3         | 6.10  | 10.9  | —        | —                | 0.8x0.7                  | —                  |
|                                              | Non-Scan      | < 0.1       | 6.10  | 0.9   | 2.5      | 0.1              | 0.8x0.7                  | 0.1                |

Table 17–75. L764 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### L764 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| Pulsed Wave Doppler                          | 0.2         | 5.16  | 5.9   | 0.8x0.7                         | 2.9           |
| B/Pulsed Wave Doppler                        | 0.4         | 5.16  | 17.0  | 0.9x0.7                         | 2.9           |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.7         | 5.16  | 22.9  | 0.9x0.7                         | 2.9           |
| Color Flow with M Doppler                    | 0.3         | 6.30  | 6.9   | 0.8x0.7                         | 2.9           |

Table 17–76. L764 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | $f_c$ [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PII <sub>max</sub> |
|-------------------------------------------|---------------|------------|------------|----------|-------------|-----------------|----------|--------------------------|
| B–Mode                                    | Scan          | 0.89       | 2.25       | 2.75     | 6.36        | 0.20            | 870      | 2.25                     |
| M–Mode                                    | Non–Scan      | 0.83       | 1.49       | 2.75     | 6.36        | 0.20            | 870      | 1.49                     |
| B/M–Mode                                  | Scan          | 0.89       | 2.25       | 2.75     | 6.36        | 0.20            | 870      | 2.25                     |
|                                           | Non–Scan      | 0.83       | 1.49       | 2.75     | 6.36        | 0.20            | 870      | 1.49                     |
| Pulsed Wave Doppler                       | Non–Scan      | 0.59       | 0.94       | 2.47     | 5.06        | 0.96            | 2000     | 0.94                     |
| B/Pulsed Wave Doppler                     | Scan          | 0.89       | 2.25       | 2.75     | 6.36        | 0.20            | 870      | 2.25                     |
|                                           | Non–Scan      | 0.59       | 0.94       | 2.47     | 5.06        | 0.96            | 2000     | 0.94                     |
| B/Color Flow Mode                         | Scan          | 0.89       | 2.25       | 2.50     | 5.16        | 0.57            | 10400    | 2.25                     |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.89       | 2.25       | 2.50     | 5.16        | 0.57            | 10400    | 2.25                     |
|                                           | Non–Scan      | 0.59       | 0.94       | 2.47     | 5.06        | 0.96            | 2000     | 0.94                     |
| M/M–CFM                                   | Non–Scan      | 0.83       | 1.49       | 2.50     | 5.16        | 0.57            | 10400    | 1.49                     |
| M–CFM/PD                                  | Non–Scan      | 0.83       | 1.49       | 2.50     | 5.16        | 0.57            | 10400    | 1.49                     |
| CWD Mode                                  | Non–Scan      | –          | –          | –        | –           | –               | –        | –                        |

Table 17–77. L764 Maximum Mechanical Index (MI)

LA39

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|--------------------------------------|-------|----------|--------------------------|---------------|
| B-Mode                                       | Scan          | 0.38    | 7.46  | 10.40 | 6.21                                 | —     | 0.62     | 0.13                     | 0.90          |
| M-Mode                                       | Non-Scan      | < 0.1   | 7.46  | 1.04  | 0.62                                 | —     | 0.62     | 0.13                     | 0.90          |
| B/M-Mode                                     | Scan          | 0.25    | 7.46  | 6.93  | 4.14                                 | —     | 0.62     | 0.13                     | 0.90          |
|                                              | Non-Scan      | < 0.1   | 7.46  | 1.04  | 0.62                                 | —     | 0.62     | 0.13                     | 0.90          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.24    | 6.74  | 7.47  | 4.76                                 | —     | 0.58     | 0.08                     | 0.75          |
| B/Pulsed Wave Doppler                        | Scan          | 0.25    | 7.46  | 6.93  | 4.14                                 | —     | 0.62     | 0.13                     | 0.90          |
|                                              | Non-Scan      | 0.24    | 6.74  | 7.47  | 4.76                                 | —     | 0.58     | 0.08                     | 0.75          |
| B/Color Flow Doppler                         | Scan          | 0.93    | 6.72  | 28.31 | 18.30                                | —     | 0.48     | 0.08                     | 0.75          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.59    | 6.72  | 17.62 | 11.22                                | —     | 0.48     | 0.08                     | 0.75          |
|                                              | Non-Scan      | 0.12    | 6.74  | 3.73  | 2.38                                 | —     | 0.48     | 0.08                     | 0.75          |
| M/M-CFM                                      | Non-Scan      | 0.57    | 6.72  | 18.85 | 12.40                                | —     | 0.48     | 0.08                     | 0.75          |
| M-CFM/PD                                     | Non-Scan      | 0.24    | 6.72  | 18.85 | 12.40                                | —     | 0.48     | 0.08                     | 0.75          |

Table 17–78. LA39 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII_{max}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| B-Mode                                       | Scan          | 0.38        | 7.46  | 10.40 | 1.00     | 0.45             | 0.13                     | —                  |
| M-Mode                                       | Non-Scan      | 0.14        | 7.46  | 1.04  | 1.00     | 0.11             | 0.13                     | —                  |
| B/M-Mode                                     | Scan          | 0.25        | 7.46  | 6.93  | 1.00     | 0.45             | 0.13                     | —                  |
|                                              | Non-Scan      | 0.14        | 7.46  | 1.04  | 1.00     | 0.11             | 0.13                     | —                  |
| Pulsed Wave Doppler                          | Non-Scan      | 1.08        | 6.74  | 7.47  | 0.99     | 0.14             | 0.12                     | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.25        | 7.46  | 6.93  | 1.00     | 0.45             | 0.13                     | —                  |
|                                              | Non-Scan      | 1.08        | 6.74  | 7.47  | 0.99     | 0.14             | 0.12                     | —                  |
| B-CFM                                        | Scan          | 0.93        | 6.72  | 28.31 | 0.90     | 0.49             | 0.08                     | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.59        | 6.72  | 17.62 | 0.90     | 0.49             | 0.08                     | —                  |
|                                              | Non-Scan      | 0.54        | 6.74  | 3.73  | 0.99     | 0.14             | 0.12                     | —                  |
| M/M-CFM                                      | Non-Scan      | 1.68        | 6.72  | 18.85 | 0.90     | 0.20             | 0.08                     | —                  |
| Color Flow with M Doppler                    | Non-Scan      | 1.08        | 6.72  | 18.85 | 0.90     | 0.20             | 0.08                     | —                  |

Table 17–79. LA39 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### LA39 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| B-Mode                                       | 0.63        | 7.46  | 10.40 | 0.13                            | 0.90          |
| M-Mode                                       | < 0.1       | 7.46  | 1.04  | 0.13                            | 0.90          |
| B/M-Mode                                     | 0.49        | 7.46  | 7.98  | 0.13                            | 0.90          |
| Pulsed Wave Doppler                          | 0.48        | 6.74  | 7.47  | 0.12                            | 0.75          |
| B/Pulsed Wave Doppler                        | 0.90        | 6.74  | 14.40 | 0.12                            | 0.75          |
| B/Color Flow Doppler                         | 2.10        | 6.72  | 28.31 | 0.08                            | 0.75          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.50        | 6.74  | 21.36 | 0.12                            | 0.75          |
| M/M-CFM                                      | 1.40        | 6.72  | 18.85 | 0.08                            | 0.75          |
| Color Flow with M Doppler                    | 0.48        | 6.72  | 18.85 | 0.08                            | 0.75          |

Table 17–80. LA39 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | fc [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|----------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 1.52       | 3.97       | 1.00     | 7.46     | 0.13            | 8658     | 3.97         |
| M-Mode                                    | Non-Scan      | 1.52       | 3.97       | 1.00     | 7.46     | 0.13            | 867      | 3.97         |
| B/M-Mode                                  | Scan          | 1.52       | 3.97       | 1.00     | 7.46     | 0.13            | 8658     | 3.97         |
|                                           | Non-Scan      | 1.52       | 3.97       | 1.00     | 7.46     | 0.13            | 867      | 3.97         |
| Pulsed Wave Doppler                       | Non-Scan      | 1.36       | 3.63       | 0.95     | 6.69     | 0.57            | 651      | 3.63         |
| B/Pulsed Wave Doppler                     | Scan          | 1.52       | 3.97       | 1.00     | 7.46     | 0.13            | 8658     | 3.97         |
|                                           | Non-Scan      | 1.36       | 3.63       | 0.95     | 6.69     | 0.57            | 651      | 3.63         |
| B/Color Flow Mode                         | Scan          | 1.52       | 3.97       | 0.90     | 6.72     | 0.42            | 660      | 3.97         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 1.52       | 3.97       | 0.90     | 6.72     | 0.42            | 660      | 3.97         |
|                                           | Non-Scan      | 1.36       | 3.63       | 0.95     | 6.69     | 0.57            | 651      | 3.63         |
| M/M-CFM                                   | Non-Scan      | 1.52       | 3.97       | 0.90     | 6.72     | 0.42            | 1480     | 3.97         |
| M-CFM/PD                                  | Non-Scan      | 1.52       | 3.97       | 0.90     | 6.72     | 0.42            | 1480     | 3.97         |

Table 17–81. LA39 Maximum Mechanical Index (MI)



LD

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$  | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|--------|--------------------------------------|-------|----------|--------------------------|---------------|
| B-Mode                                       | Scan          | 0.21    | 3.50  | 38.80  | 4.43                                 | —     | 2.98     | 3.11                     | 8.00          |
| M-Mode                                       | Non-Scan      | <0.1    | 3.28  | 4.33   | 0.88                                 | —     | 2.98     | 3.11                     | 8.00          |
| B/M-Mode                                     | Scan          | 0.14    | 3.50  | 25.87  | 2.95                                 | —     | 2.98     | 3.11                     | 8.00          |
|                                              | Non-Scan      | <0.1    | 3.28  | 4.33   | 0.88                                 | —     | 2.98     | 3.11                     | 8.00          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.33    | 3.32  | 125.61 | 20.61                                | —     | 2.78     | 2.69                     | 8.20          |
| B/Pulsed Wave Doppler                        | Scan          | 0.14    | 3.50  | 25.87  | 2.95                                 | —     | 2.98     | 3.11                     | 8.00          |
|                                              | Non-Scan      | 0.33    | 3.32  | 125.61 | 20.61                                | —     | 2.78     | 2.69                     | 8.20          |
| B/CFM-Mode                                   | Scan          | 1.35    | 3.30  | 103.01 | 19.11                                | —     | 2.78     | 2.69                     | 8.20          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.74    | 3.30  | 64.44  | 11.03                                | —     | 2.78     | 2.69                     | 8.20          |
|                                              | Non-Scan      | 0.16    | 3.32  | 62.80  | 10.31                                | —     | 2.78     | 2.69                     | 8.20          |
| M/M-CFM                                      | Non-Scan      | 0.29    | 3.30  | 91.59  | 18.30                                | —     | 2.78     | 2.69                     | 8.20          |
| M-CFM/PD                                     | Non-Scan      | 0.29    | 3.30  | 91.59  | 18.30                                | —     | 2.78     | 2.69                     | 8.20          |

Table 17–82. LD Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$  | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PIL_{max}$ |
|----------------------------------------------|---------------|-------------|-------|--------|----------|------------------|--------------------------|--------------------|
| B-Mode                                       | Scan          | 0.21        | 3.50  | 38.80  | 6.64     | 4.07             | 3.11                     | —                  |
| M-Mode                                       | Non-Scan      | <0.1        | 3.28  | 4.33   | 7.07     | 0.48             | 3.11                     | —                  |
| B/M-Mode                                     | Scan          | 0.14        | 3.50  | 25.87  | 6.64     | 4.07             | 3.11                     | —                  |
|                                              | Non-Scan      | <0.1        | 3.28  | 4.33   | 7.07     | 0.48             | 3.11                     | —                  |
| Pulsed Wave Doppler                          | Non-Scan      | 2.32        | 3.32  | 125.61 | 8.08     | 0.57             | 2.69                     | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.14        | 3.50  | 25.87  | 6.64     | 4.07             | 3.11                     | —                  |
|                                              | Non-Scan      | 2.32        | 3.32  | 125.61 | 8.08     | 0.57             | 2.69                     | —                  |
| B-CFM                                        | Scan          | 1.35        | 3.30  | 103.01 | 7.31     | 1.45             | 2.69                     | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.74        | 3.30  | 64.44  | 7.31     | 1.45             | 2.69                     | —                  |
|                                              | Non-Scan      | 1.16        | 3.32  | 62.80  | 8.08     | 0.57             | 2.69                     | —                  |
| M/M-CFM                                      | Non-Scan      | 2.13        | 3.30  | 91.59  | 7.31     | 0.52             | 2.69                     | —                  |
| M-CFM/PD                                     | Non-Scan      | 2.13        | 3.30  | 91.59  | 7.31     | 0.52             | 2.69                     | —                  |

Table 17–83. LD Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### LD (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$  | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|--------|---------------------------------|---------------|
| B-Mode                                       | 0.49        | 3.50  | 38.80  | 3.11                            | 8.00          |
| M-Mode                                       | <0.1        | 3.28  | 4.33   | 3.11                            | 8.00          |
| B/M-Mode                                     | 0.38        | 3.50  | 30.20  | 3.11                            | 8.00          |
| Pulsed Wave Doppler                          | 1.70        | 3.32  | 125.61 | 2.69                            | 8.20          |
| B/Pulsed Wave Doppler                        | 2.02        | 3.32  | 151.47 | 2.69                            | 8.20          |
| B-CFM                                        | 1.37        | 3.30  | 103.01 | 2.69                            | 8.20          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.69        | 3.32  | 127.24 | 2.69                            | 8.20          |
| M/M-CFM                                      | 1.23        | 3.30  | 91.59  | 2.69                            | 8.20          |
| Color Flow with M Doppler                    | 1.23        | 3.30  | 91.59  | 2.69                            | 8.20          |

Table 17-84. LD Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | fc [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|----------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 0.63       | 1.05       | 7.07     | 3.28     | 0.21            | 2667     | 1.05         |
| M-Mode                                    | Non-Scan      | 0.63       | 1.05       | 7.07     | 3.28     | 0.21            | 888      | 1.05         |
| B/M-Mode                                  | Scan          | 0.63       | 1.05       | 7.07     | 3.28     | 0.21            | 2667     | 1.05         |
|                                           | Non-Scan      | 0.63       | 1.05       | 7.07     | 3.28     | 0.21            | 888      | 1.05         |
| Pulsed Wave Doppler                       | Non-Scan      | 0.58       | 0.99       | 7.56     | 3.22     | 0.80            | 651      | 0.99         |
| B/Pulsed Wave Doppler                     | Scan          | 0.63       | 1.05       | 7.07     | 3.28     | 0.21            | 2667     | 1.05         |
|                                           | Non-Scan      | 0.58       | 0.99       | 7.56     | 3.22     | 0.80            | 651      | 0.99         |
| B/Color Flow Mode                         | Scan          | 0.90       | 1.65       | 6.91     | 3.28     | 1.53            | 280      | 1.65         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.90       | 1.65       | 6.91     | 3.28     | 1.53            | 280      | 1.65         |
|                                           | Non-Scan      | 0.58       | 0.99       | 7.56     | 3.22     | 0.80            | 651      | 0.99         |
| M/M-CFM                                   | Non-Scan      | 0.79       | 1.05       | 6.91     | 3.28     | 1.53            | 1480     | 1.05         |
| M-CFM/PD                                  | Non-Scan      | 0.79       | 1.05       | 6.91     | 3.28     | 1.53            | 1480     | 1.05         |

Table 17-85. LD Maximum Mechanical Index (MI)

## FDA Acoustic Output Data

### P509

| Imaging Mode                                 | Scanning Mode | Max TIS | fc   | Wo    | Min of W.3(Z1) ITA.3(Z1) | Z1 | Zbp  | Dimen- sions of Aaprt | FL Focal Zone |
|----------------------------------------------|---------------|---------|------|-------|--------------------------|----|------|-----------------------|---------------|
| B-Mode                                       | Scan          | 0.25    | 4.93 | 11.03 | 5.37                     | –  | 1.57 | 0.86                  | 2.00          |
| M-Mode                                       | Non-Scan      | < 0.1   | 4.93 | 3.72  | 1.81                     | –  | 1.57 | 0.86                  | 2.00          |
| B/M-Mode                                     | Scan          | 0.16    | 4.93 | 7.35  | 3.58                     | –  | 1.57 | 0.86                  | 2.00          |
|                                              | Non-Scan      | < 0.1   | 4.93 | 3.72  | 1.81                     |    | 1.57 | 0.86                  | 2.00          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.56    | 4.02 | 29.55 | 16.28                    | –  | 1.68 | 0.99                  | 2.00          |
| B/Pulsed Wave Doppler                        | Scan          | 0.16    | 4.93 | 7.35  | 3.58                     | –  | 1.57 | 0.86                  | 2.00          |
|                                              | Non-Scan      | 0.56    | 4.02 | 29.55 | 16.28                    |    | 1.68 | 0.99                  | 2.00          |
| B/CFM-Mode                                   | Scan          | 0.63    | 4.05 | 31.52 | 17.76                    | –  | 1.68 | 0.99                  | 2.00          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.40    | 4.05 | 19.44 | 10.67                    | –  | 1.68 | 0.99                  | 2.00          |
|                                              | Non-Scan      | 0.28    | 4.02 | 14.77 | 8.14                     |    | 1.68 | 0.99                  | 2.00          |
| M/M-CFM                                      | Non-Scan      | 0.48    | 4.05 | 24.14 | 13.47                    | –  | 1.68 | 0.99                  | 2.00          |
| M-CFM/PD                                     | Non-Scan      | 0.48    | 4.05 | 24.14 | 13.47                    | –  | 1.68 | 0.99                  | 2.00          |
| CWD                                          | Non-Scan      | 1.33    | 4.00 | 69.43 | 36.62                    | –  | 1.05 | 0.38                  | 11.90         |

Table 17–86. P509 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                         | Scanning Mode | Max TIB | fc   | Wo    | Zsp  | deq(Zsp) | Dimen- sions of Aaprt | deq @Pllmax |
|--------------------------------------|---------------|---------|------|-------|------|----------|-----------------------|-------------|
| B-Mode                               | Scan          | 0.25    | 4.93 | 11.03 | 2.06 | 0.36     | 0.86                  | –           |
| M-Mode                               | Non-Scan      | 0.18    | 4.93 | 3.72  | 2.06 | 0.35     | 0.86                  | –           |
| B/M-Mode                             | Scan          | 0.16    | 4.93 | 7.35  | 2.06 | 0.36     | 0.86                  | –           |
|                                      | Non-Scan      | 0.18    | 4.93 | 3.72  | 2.06 | 0.35     | 0.86                  | –           |
| Pulsed Wave Doppler                  | Non-Scan      | 1.63    | 4.02 | 29.55 | 2.20 | 0.31     | 0.99                  | –           |
| B/Pulsed Wave Doppler                | Scan          | 0.16    | 4.93 | 7.35  | 2.06 | 0.36     | 0.86                  | –           |
|                                      | Non-Scan      | 1.63    | 4.02 | 29.55 | 2.20 | 0.31     | 0.99                  | –           |
| B-CFM                                | Scan          | 0.63    | 4.05 | 31.52 | 2.13 | 0.38     | 0.99                  | –           |
| B/Color Flow Doppler/<br>Pulsed Wave | Scan          | 0.40    | 4.05 | 19.44 | 2.13 | 0.38     | 0.99                  | –           |
|                                      | Non-Scan      | 0.82    | 4.02 | 14.77 | 2.20 | 0.31     | 0.99                  | –           |
| M/M-CFM                              | Non-Scan      | 1.56    | 4.05 | 24.14 | 2.13 | 0.33     | 0.99                  | –           |
| M-CFM/PD                             | Non-Scan      | 1.56    | 4.05 | 24.14 | 2.13 | 0.33     | 0.99                  | –           |
| CWD                                  | Non-Scan      | 1.73    | 4.00 | 69.43 | 2.42 | 0.83     | 0.38                  | 0.83        |

Table 17–87. P509 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### P509 (cont'd)

| Imaging Mode                                 | Max TIC | fc   | Wo    | Dimensions of Aaprt | FL Focal Zone |
|----------------------------------------------|---------|------|-------|---------------------|---------------|
| B–Mode                                       | 0.26    | 4.93 | 11.03 | 0.86                | 2.00          |
| M–Mode                                       | < 0.1   | 4.93 | 3.72  | 0.86                | 2.00          |
| B/M–Mode                                     | 0.26    | 4.93 | 11.08 | 0.86                | 2.00          |
| Pulsed Wave Doppler                          | 0.66    | 4.02 | 29.55 | 0.99                | 2.00          |
| B/Pulsed Wave Doppler                        | 0.83    | 4.02 | 36.9  | 0.99                | 2.00          |
| B/Color Flow Doppler                         | 0.72    | 4.05 | 31.52 | 0.99                | 2.00          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.77    | 4.02 | 34.21 | 0.99                | 2.00          |
| M/M–CFM                                      | 0.54    | 4.05 | 24.14 | 0.99                | 2.00          |
| M–CFM/PD                                     | 0.54    | 4.05 | 24.14 | 0.99                | 2.00          |
| CWD                                          | 2.49    | 4.00 | 69.43 | 0.38                | 11.90         |

Table 17–88. P509 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | fc [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PII <sub>max</sub> |
|-------------------------------------------|---------------|------------|------------|----------|----------|-----------------|----------|--------------------------|
| B–Mode                                    | Scan          | 0.42       | 0.96       | 4.19     | 4.93     | 0.28            | 2667     | 0.96                     |
| M–Mode                                    | Non–Scan      | 0.42       | 0.96       | 4.19     | 4.93     | 0.28            | 900      | 0.96                     |
| B/M–Mode                                  | Scan          | 0.42       | 0.96       | 4.19     | 4.93     | 0.28            | 2667     | 0.96                     |
|                                           | Non–Scan      | 0.42       | 0.96       | 4.19     | 4.93     | 0.28            | 900      | 0.96                     |
| Pulsed Wave Doppler                       | Non–Scan      | 0.62       | 1.25       | 3.61     | 4.23     | 0.43            | 650      | 1.25                     |
| B/Pulsed Wave Doppler                     | Scan          | 0.42       | 0.96       | 4.19     | 4.93     | 0.28            | 2667     | 0.96                     |
|                                           | Non–Scan      | 0.62       | 1.25       | 3.61     | 4.23     | 0.43            | 650      | 1.25                     |
| B/Color Flow Mode                         | Scan          | 0.42       | 0.96       | 4.53     | 5.00     | 0.91            | 5200     | 0.96                     |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.42       | 0.96       | 4.53     | 5.00     | 0.91            | 5200     | 0.96                     |
|                                           | Non–Scan      | 0.62       | 1.25       | 3.61     | 4.23     | 0.43            | 650      | 1.25                     |
| M/M–CFM                                   | Non–Scan      | 0.42       | 0.96       | 4.53     | 5.00     | 0.92            | 5200     | 0.96                     |
| M–CFM/PD                                  | Non–Scan      | 0.42       | 0.96       | 4.53     | 5.00     | 0.92            | 5200     | 0.96                     |
| CWD Mode                                  | Non–Scan      | 0.01       | 0.06       | 4.88     | 4.00     | –               | –        | 0.06                     |

Table 17–89. P509 Maximum Mechanical Index (MI)

S220

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$  | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|--------|--------------------------------------|-------|----------|--------------------------|---------------|
| Pulsed Wave Doppler                          | Scan          | —       | —     | —      | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.59    | 2.49  | 85.73  | 29.04                                | 2.87  | 2.87     | 2.87                     | 6.75          |
| B/Pulsed Wave Doppler                        | Scan          | 0.29    | 2.44  | 34.92  | 10.10                                | —     | 3.07     | 3.28                     | 9.60          |
|                                              | Non-Scan      | 0.59    | 2.49  | 85.73  | 29.04                                | 2.87  | 2.87     | 2.87                     | 6.75          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.71    | 2.27  | 72.80  | 23.15                                | 3.07  | 3.07     | 3.28                     | 6.75          |
|                                              | Non-Scan      | 0.29    | 2.49  | 42.86  | 14.52                                | 2.87  | 2.87     | 2.87                     | 6.75          |
| Color Flow with<br>M Doppler                 | Scan          | —       | —     | —      | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.53    | 2.27  | 75.76  | 26.09                                | 3.07  | 3.07     | 3.28                     | 6.75          |
| CW Doppler                                   | Scan          | —       | —     | —      | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.84    | 2.50  | 117.15 | 35.20                                | 4.5   | 2.03     | 1.43                     | 10.00         |

Table 17–90. S220 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$  | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq@PII_{max}}$ |
|----------------------------------------------|---------------|-------------|-------|--------|----------|------------------|--------------------------|--------------------|
| Pulsed Wave Doppler                          | Scan          | —           | —     | —      | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 2.52        | 2.22  | 66.01  | 5.95     | 0.77             | 2.87                     | 1.18               |
| B/Pulsed Wave Doppler                        | Scan          | 0.29        | 2.44  | 34.92  | 7.30     | 1.14             | 3.28                     | —                  |
|                                              | Non-Scan      | 2.52        | 2.22  | 66.01  | 5.95     | 0.77             | 2.87                     | 0.33               |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.71        | 2.27  | 72.80  | 6.78     | 1.16             | 2.87                     | 1.18               |
|                                              | Non-Scan      | 1.26        | 2.22  | 33.01  | 5.95     | 0.77             | 2.87                     | 0.33               |
| Color Flow with<br>M Doppler                 | Scan          | —           | —     | —      | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 2.71        | 2.27  | 75.76  | 6.78     | 0.79             | 3.28                     | 1.18               |
| CW Doppler                                   | Scan          | —           | —     | —      | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 1.61        | 2.50  | 117.15 | 6.97     | 0.97             | 1.23                     | 0.97               |

Table 17–91. S220 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### S220 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$  | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|--------|---------------------------------|---------------|
| Pulsed Wave Doppler                          | 1.12        | 2.49  | 85.73  | 2.87                            | 6.75          |
| B/Pulsed Wave Doppler                        | 1.55        | 2.49  | 120.65 | 2.87                            | 6.75          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.45        | 2.49  | 105.71 | 2.87                            | 6.75          |
| Color Flow with M Doppler                    | 0.93        | 2.43  | 55.86  | 2.87                            | 6.75          |
| CW Doppler                                   | 2.17        | 2.50  | 117.15 | 1.43                            | 10.00         |

Table 17–92. S220 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | fc [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | Pri3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|----------|-----------------|----------|--------------|
| B–Mode                                    | Scan          | 0.51       | 0.79       | 7.30     | 2.44     | 0.70            | 4500     | 0.79         |
| M–Mode                                    | Non–Scan      | 0.51       | 0.79       | 7.30     | 2.44     | 0.70            | 814      | 0.79         |
| B/M–Mode                                  | Scan          | 0.51       | 0.79       | 7.30     | 2.44     | 0.70            | 4500     | 0.79         |
|                                           | Non–Scan      | 0.51       | 0.79       | 7.30     | 2.44     | 0.70            | 814      | 0.79         |
| Pulsed Wave Doppler                       | Non–Scan      | 0.58       | 0.89       | 6.30     | 2.41     | 1.08            | 600      | 0.89         |
| B/Pulsed Wave Doppler                     | Scan          | 0.51       | 0.79       | 7.30     | 2.44     | 0.70            | 4500     | 0.79         |
|                                           | Non–Scan      | 0.58       | 0.89       | 6.30     | 2.41     | 1.08            | 600      | 0.89         |
| B/Color Flow Mode                         | Scan          | 0.83       | 1.29       | 6.14     | 2.42     | 1.07            | 200      | 1.29         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.83       | 1.29       | 6.14     | 2.42     | 1.07            | 200      | 1.29         |
|                                           | Non–Scan      | 0.58       | 0.89       | 6.30     | 2.41     | 1.08            | 600      | 0.89         |
| M/M–CFM                                   | Non–Scan      | 0.83       | 1.29       | 6.14     | 2.42     | 1.07            | 200      | 1.29         |
| M–CFM/PD                                  | Non–Scan      | 0.83       | 1.29       | 6.14     | 2.42     | 1.07            | 200      | 1.29         |
| CWD Mode                                  | Non–Scan      | 0.03       | 0.14       | 6.97     | 2.50     | –               | –        | 0.14         |

Table 17–93. S220 Maximum Mechanical Index (MI)

S222

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$  | Min of W.3(Z1) ITA.3(Z1) | Z1 | Zbp  | Dimensions of Aaprt | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|--------|--------------------------|----|------|---------------------|---------------|
| B-Mode                                       | Scan          | 0.85    | 1.99  | 102.61 | 53.83                    | —  | 1.95 | 1.33                | 4.60          |
| M-Mode                                       | Non-Scan      | 0.17    | 1.99  | 33.86  | 17.76                    | —  | 1.95 | 1.33                | 4.60          |
| B/M-Mode                                     | Scan          | 0.57    | 1.99  | 68.41  | 35.89                    | —  | 1.95 | 1.33                | 4.60          |
|                                              | Non-Scan      | 0.17    | 1.99  | 33.86  | 17.76                    | —  | 1.95 | 1.33                | 4.60          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.22    | 2.49  | 61.79  | 18.68                    | —  | 2.76 | 2.66                | 7.30          |
| B/Pulsed Wave Doppler                        | Scan          | < 0.1   | 1.99  | 68.41  | 35.89                    | —  | 1.95 | 1.33                | 4.60          |
|                                              | Non-Scan      | 0.22    | 2.49  | 61.79  | 18.68                    | —  | 2.76 | 2.66                | 7.30          |
| B/CFM-Mode                                   | Scan          | 2.05    | 2.49  | 241.43 | 84.42                    | —  | 2.76 | 2.66                | 7.30          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 1.03    | 2.49  | 154.92 | 60.16                    | —  | 2.76 | 2.66                | 7.30          |
|                                              | Non-Scan      | 0.11    | 2.49  | 30.89  | 9.34                     | —  | 2.76 | 2.66                | 7.30          |
| Color Flow with M Doppler                    | Scan          | 0.21    | 2.49  | 95.63  | 35.20                    | —  | 2.76 | 2.66                | 7.30          |
|                                              | Non-Scan      | 0.22    | 2.49  | 95.63  | 35.20                    | —  | 2.76 | 2.66                | 7.30          |
| CWD Doppler                                  | Non-Scan      | 2.20    | 2.51  | 194.47 | 40.43                    | —  | 1.69 | 1.00                | 10.00         |

Table 17–94. S222 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$  | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of Aaprt | $d_{eq@PII-max}$ |
|----------------------------------------------|---------------|-------------|-------|--------|----------|------------------|---------------------|------------------|
| B-Mode                                       | Scan          | 0.85        | 1.99  | 102.61 | 4.66     | 0.60             | 1.33                | —                |
| M-Mode                                       | Non-Scan      | 1.75        | 1.99  | 33.86  | 4.66     | 0.33             | 1.33                | —                |
| B/M-Mode                                     | Scan          | 0.57        | 1.99  | 68.41  | 4.66     | 0.60             | 1.33                | —                |
|                                              | Non-Scan      | 1.75        | 1.99  | 33.86  | 4.66     | 0.33             | 1.33                | —                |
| Pulsed Wave Doppler                          | Non-Scan      | 2.27        | 2.49  | 61.79  | 7.16     | 0.51             | 2.66                | —                |
| B/Pulsed Wave Doppler                        | Scan          | < 0.1       | 1.99  | 68.41  | 4.66     | 0.60             | 1.33                | —                |
|                                              | Non-Scan      | 2.27        | 2.49  | 61.79  | 7.16     | 0.51             | 2.66                | —                |
| B/CFM-Mode                                   | Scan          | 2.05        | 2.49  | 241.43 | 7.32     | 1.16             | 2.66                | —                |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 1.03        | 2.49  | 154.92 | 7.32     | 1.16             | 2.66                | —                |
|                                              | Non-Scan      | 1.14        | 2.49  | 30.89  | 7.16     | 0.51             | 2.66                | —                |
| M/M-CFM                                      | Non-Scan      | 4.01        | 2.49  | 95.63  | 7.32     | 0.42             | 2.66                | —                |
| Color Flow with M Doppler                    | Scan          | —           | —     | —      | —        | —                | —                   | —                |
|                                              | Non-Scan      | 2.27        | 2.49  | 95.63  | 7.32     | 0.42             | 2.66                | —                |
| CW Doppler                                   | Scan          | —           | —     | —      | —        | —                | —                   | —                |
|                                              | Non-Scan      | 2.25        | 2.51  | 194.47 | 7.61     | 1.33             | 1.00                | 1.33             |

Table 17–95. S222 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### S222 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$  | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|--------|---------------------------------|---------------|
| B-Mode                                       | 1.97        | 1.99  | 102.61 | 1.33                            | 4.60          |
| M-Mode                                       | 0.65        | 1.99  | 33.86  | 1.33                            | 4.60          |
| B/M-Mode                                     | 1.96        | 1.99  | 102.27 | 1.33                            | 4.60          |
| Pulsed Wave Doppler                          | 0.84        | 2.49  | 61.79  | 2.66                            | 7.30          |
| B/Pulsed Wave Doppler                        | 0.87        | 2.49  | 130.19 | 2.66                            | 7.30          |
| B/Color Flow Mode                            | 2.41        | 2.49  | 241.43 | 2.66                            | 7.30          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.63        | 2.49  | 185.81 | 2.66                            | 7.30          |
| M/M-CFM                                      | 0.86        | 2.49  | 95.63  | 2.66                            | 7.30          |
| M-CFM/PD                                     | 0.84        | 2.49  | 95.63  | 2.66                            | 7.30          |
| CW Doppler                                   | 4.32        | 2.51  | 194.47 | 1.00                            | 10.00         |

Table 17-96. S222 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | fc [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|----------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 1.44       | 1.94       | 4.66     | 1.99     | 0.83            | 2667     | 1.94         |
| M-Mode                                    | Non-Scan      | 1.44       | 1.94       | 4.66     | 1.99     | 0.83            | 880      | 1.94         |
| B/M-Mode                                  | Scan          | 1.44       | 1.94       | 4.66     | 1.99     | 0.83            | 2667     | 1.94         |
|                                           | Non-Scan      | 1.44       | 1.94       | 4.66     | 1.99     | 0.83            | 880      | 1.94         |
| Pulsed Wave Doppler                       | Non-Scan      | 0.70       | 0.99       | 7.46     | 2.51     | 0.54            | 651      | 0.99         |
| B/Pulsed Wave Doppler                     | Scan          | 1.44       | 1.94       | 4.66     | 1.99     | 0.83            | 2667     | 1.94         |
|                                           | Non-Scan      | 0.70       | 0.99       | 7.46     | 2.51     | 0.54            | 651      | 0.99         |
| B/Color Flow Mode                         | Scan          | 1.44       | 1.94       | 7.32     | 2.49     | 1.97            | 5960     | 1.94         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 1.44       | 1.94       | 7.32     | 2.49     | 1.97            | 5960     | 1.94         |
|                                           | Non-Scan      | 0.70       | 0.99       | 7.46     | 2.51     | 0.54            | 651      | 0.99         |
| M/M-CFM                                   | Non-Scan      | 1.44       | 1.94       | 7.32     | 2.49     | 1.98            | 5960     | 1.94         |
| M-CFM/PD                                  | Non-Scan      | 1.44       | 1.94       | 7.32     | 2.49     | 1.98            | 5960     | 1.94         |
| CWD Mode                                  | Non-Scan      | 0.15       | 0.17       | 7.61     | 2.51     | —               | —        | 0.17         |

Table 17-97. S222 Maximum Mechanical Index (MI)



**S316**

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_{3(z_1)}$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|----------------------------------------|-------|----------|--------------------------|---------------|
| B-Mode                                       | Scan          | 0.16    | 3.42  | 16.83 | 4.68                                   | —     | 2.44     | 2.08                     | 5.77          |
| M-Mode                                       | Non-Scan      | < 0.1   | 3.36  | 5.00  | 1.38                                   | —     | 2.44     | 2.08                     | 5.77          |
| B/M-Mode                                     | Scan          | 0.11    | 3.42  | 11.22 | 3.12                                   | —     | 2.44     | 2.08                     | 5.77          |
|                                              | Non-Scan      | < 0.1   | 3.36  | 5.00  | 1.38                                   | —     | 2.44     | 2.08                     | 5.77          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.11    | 3.00  | 27.98 | 7.71                                   | —     | 2.32     | 1.89                     | 5.00          |
| B/Pulsed Wave Doppler                        | Scan          | 0.11    | 3.42  | 11.22 | 3.12                                   | —     | 2.44     | 2.08                     | 5.77          |
|                                              | Non-Scan      | 0.11    | 3.00  | 27.98 | 7.71                                   | —     | 2.32     | 1.89                     | 5.00          |
| B/CFM-Mode                                   | Scan          | 1.24    | 2.98  | 91.88 | 28.66                                  | —     | 2.32     | 1.89                     | 5.00          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.67    | 2.98  | 51.55 | 15.89                                  | —     | 2.32     | 1.89                     | 5.00          |
|                                              | Non-Scan      | < 0.1   | 3.00  | 13.99 | 3.86                                   | —     | 2.32     | 1.89                     | 5.00          |
| M/M-CFM                                      | Non-Scan      | 0.16    | 2.98  | 36.42 | 11.42                                  | —     | 2.32     | 1.89                     | 5.00          |
| M-CFM/PD                                     | Non-Scan      | 0.16    | 2.98  | 36.42 | 11.42                                  | —     | 2.32     | 1.89                     | 5.00          |

Table 17–98. S316 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII_{max}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| B-Mode                                       | Scan          | 0.16        | 3.42  | 16.83 | 5.89     | 1.29             | 2.08                     | —                  |
| M-Mode                                       | Non-Scan      | 0.14        | 3.36  | 5.00  | 5.57     | 0.46             | 2.08                     | —                  |
| B/M-Mode                                     | Scan          | 0.11        | 3.42  | 11.22 | 5.89     | 1.29             | 2.08                     | —                  |
|                                              | Non-Scan      | 0.14        | 3.36  | 5.00  | 5.57     | 0.46             | 2.08                     | —                  |
| Pulsed Wave Doppler                          | Non-Scan      | 0.88        | 3.00  | 27.98 | 6.23     | 0.49             | 1.89                     | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.11        | 3.42  | 11.22 | 5.89     | 1.29             | 2.08                     | —                  |
|                                              | Non-Scan      | 0.88        | 3.00  | 27.98 | 6.23     | 0.49             | 1.89                     | —                  |
| B-CFM                                        | Scan          | 1.24        | 2.98  | 91.88 | 5.75     | 1.46             | 1.89                     | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.67        | 2.98  | 51.55 | 5.75     | 1.46             | 1.89                     | —                  |
|                                              | Non-Scan      | 0.44        | 3.00  | 13.99 | 6.23     | 0.49             | 1.89                     | —                  |
| M/M-CFM                                      | Non-Scan      | 0.66        | 2.98  | 36.42 | 5.75     | 0.79             | 1.89                     | —                  |
| M-CFM/PD                                     | Non-Scan      | 0.66        | 2.98  | 36.42 | 5.75     | 0.79             | 1.89                     | —                  |

Table 17–99. S316 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### S316 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| B-Mode                                       | 0.26        | 3.42  | 16.83 | 2.08                            | 5.77          |
| M-Mode                                       | < 0.1       | 3.36  | 5.00  | 2.08                            | 5.77          |
| B/M-Mode                                     | 0.25        | 3.42  | 16.22 | 2.08                            | 5.77          |
| Pulsed Wave Doppler                          | 0.45        | 3.00  | 27.98 | 1.89                            | 5.00          |
| B/Pulsed Wave Doppler                        | 0.62        | 3.00  | 39.20 | 1.89                            | 5.00          |
| B-CFM                                        | 1.47        | 2.98  | 91.88 | 1.89                            | 5.00          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.05        | 3.00  | 65.54 | 1.89                            | 5.00          |
| M/M-CFM                                      | 0.58        | 2.98  | 36.42 | 1.89                            | 5.00          |
| Color Flow with M Doppler                    | 0.58        | 2.98  | 36.42 | 1.89                            | 5.00          |

Table 17–100. S316 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | fc [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|----------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 0.64       | 1.09       | 5.57     | 3.36     | 0.25            | 2667     | 1.09         |
| M-Mode                                    | Non-Scan      | 0.64       | 1.09       | 5.57     | 3.36     | 0.25            | 890      | 1.09         |
| B/M-Mode                                  | Scan          | 0.64       | 1.09       | 5.57     | 3.36     | 0.25            | 2667     | 1.09         |
|                                           | Non-Scan      | 0.64       | 1.09       | 5.57     | 3.36     | 0.25            | 890      | 1.09         |
| Pulsed Wave Doppler                       | Non-Scan      | 0.70       | 1.24       | 5.77     | 2.98     | 0.89            | 650      | 1.24         |
| B/Pulsed Wave Doppler                     | Scan          | 0.64       | 1.09       | 5.57     | 3.36     | 0.25            | 2667     | 1.09         |
|                                           | Non-Scan      | 0.70       | 1.24       | 5.77     | 2.98     | 0.89            | 650      | 1.24         |
| B/Color Flow Mode                         | Scan          | 0.72       | 1.28       | 5.75     | 2.98     | 0.89            | 5200     | 1.28         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.72       | 1.28       | 5.75     | 2.98     | 0.89            | 5200     | 1.28         |
|                                           | Non-Scan      | 0.70       | 1.24       | 5.77     | 2.98     | 0.89            | 650      | 1.24         |
| M/M-CFM                                   | Non-Scan      | 0.64       | 1.09       | 5.75     | 2.98     | 0.89            | 5200     | 1.09         |
| M-CFM/PD                                  | Non-Scan      | 0.64       | 1.09       | 5.75     | 2.98     | 0.89            | 5200     | 1.09         |
| CWD Mode                                  | Non-Scan      | 0.06       | 0.08       | 5.85     | 2.51     | –               | –        | 0.08         |

Table 17–101. S316 Maximum Mechanical Index (MI)

S317

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$  | Min of $W_{3.3}(z_1)$<br>$I_{TA.3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|--------|------------------------------------------|-------|----------|--------------------------|---------------|
| B-Mode                                       | Scan          | 1.34    | 2.67  | 122.81 | 53.31                                    | —     | 2.10     | 1.55                     | 4.60          |
| M-Mode                                       | Non-Scan      | 0.11    | 2.67  | 20.36  | 8.84                                     | —     | 2.10     | 1.55                     | 4.60          |
| B/M-Mode                                     | Scan          | 0.89    | 2.67  | 81.87  | 35.54                                    | —     | 2.10     | 1.55                     | 4.60          |
|                                              | Non-Scan      | 0.11    | 2.67  | 20.36  | 8.84                                     | —     | 2.10     | 1.55                     | 4.60          |
| Pulsed Wave Doppler                          | Non-Scan      | 0.26    | 2.51  | 70.54  | 22.11                                    | —     | 2.48     | 2.15                     | 6.70          |
| B/Pulsed Wave Doppler                        | Scan          | 0.89    | 2.67  | 81.87  | 35.54                                    | —     | 2.10     | 1.55                     | 4.60          |
|                                              | Non-Scan      | 0.26    | 2.51  | 70.54  | 22.11                                    | —     | 2.48     | 2.15                     | 6.70          |
| B/CFM-Mode                                   | Scan          | 2.67    | 2.56  | 227.34 | 83.63                                    | —     | 2.48     | 2.15                     | 6.70          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 1.78    | 2.56  | 154.61 | 59.59                                    | —     | 2.48     | 2.15                     | 6.70          |
|                                              | Non-Scan      | 0.13    | 2.51  | 35.27  | 11.05                                    | —     | 2.48     | 2.15                     | 6.70          |
| M/M-CFM                                      | Non-Scan      | 0.34    | 2.56  | 77.86  | 27.68                                    | —     | 2.48     | 2.15                     | 6.70          |
| M-CFM/PD                                     | Non-Scan      | 0.34    | 2.56  | 77.86  | 27.68                                    | —     | 2.48     | 2.15                     | 6.70          |
| CWD                                          | Non-Scan      | 1.20    | 2.49  | 96.14  | 49.86                                    | —     | 1.52     | 0.81                     | 12.90         |

Table 17–102. S317 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$  | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII-max}$ |
|----------------------------------------------|---------------|-------------|-------|--------|----------|------------------|--------------------------|-------------------|
| B-Mode                                       | Scan          | 1.34        | 2.67  | 122.81 | 4.54     | 0.71             | 1.55                     | —                 |
| M-Mode                                       | Non-Scan      | 0.91        | 2.67  | 20.36  | 4.54     | 0.34             | 1.55                     | —                 |
| B/M-Mode                                     | Scan          | 0.89        | 2.67  | 81.87  | 4.54     | 0.71             | 1.55                     | —                 |
|                                              | Non-Scan      | 0.91        | 2.67  | 20.36  | 4.54     | 0.34             | 1.55                     | —                 |
| Pulsed Wave Doppler                          | Non-Scan      | 2.01        | 2.51  | 70.54  | 6.60     | 0.54             | 2.15                     | —                 |
| B/Pulsed Wave Doppler                        | Scan          | 0.89        | 2.67  | 81.87  | 4.54     | 0.71             | 1.55                     | —                 |
|                                              | Non-Scan      | 2.01        | 2.51  | 70.54  | 6.60     | 0.54             | 2.15                     | —                 |
| B-CFM                                        | Scan          | 2.67        | 2.56  | 227.34 | 6.32     | 1.75             | 2.15                     | —                 |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 1.78        | 2.56  | 154.61 | 6.32     | 1.75             | 2.15                     | —                 |
|                                              | Non-Scan      | 1.00        | 2.51  | 35.27  | 6.60     | 0.54             | 2.15                     | —                 |
| M/M-CFM                                      | Non-Scan      | 2.41        | 2.56  | 77.86  | 6.32     | 0.49             | 2.15                     | —                 |
| M-CFM/PD                                     | Non-Scan      | 2.41        | 2.56  | 77.86  | 6.32     | 0.49             | 2.15                     | —                 |
| CWD                                          | Non-Scan      | 3.07        | 2.49  | 96.14  | 3.99     | 0.54             | 0.81                     | 0.54              |

Table 17–103. S317 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### S317 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$  | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|--------|---------------------------------|---------------|
| B-Mode                                       | 2.19        | 2.67  | 122.81 | 1.55                            | 4.60          |
| M-Mode                                       | 0.36        | 2.67  | 20.36  | 1.55                            | 4.60          |
| B/M-Mode                                     | 1.82        | 2.67  | 102.23 | 1.55                            | 4.60          |
| Pulsed Wave Doppler                          | 1.07        | 2.51  | 70.54  | 2.15                            | 6.70          |
| B/Pulsed Wave Doppler                        | 2.52        | 2.51  | 152.42 | 2.15                            | 6.70          |
| B-CFM                                        | 3.66        | 2.56  | 227.34 | 2.15                            | 6.70          |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 3.09        | 2.51  | 189.88 | 2.15                            | 6.70          |
| M/M-CFM                                      | 1.23        | 2.56  | 77.86  | 2.15                            | 6.70          |
| Color Flow with M Doppler                    | 1.23        | 2.56  | 77.86  | 2.15                            | 6.70          |
| CWD                                          | 2.37        | 2.49  | 96.14  | 0.81                            | 12.90         |

Table 17–104. S317 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | $f_c$ [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|-------------|-----------------|----------|--------------|
| B-Mode                                    | Scan          | 1.09       | 1.86       | 4.54     | 2.67        | 0.63            | 5176     | 1.86         |
| M-Mode                                    | Non-Scan      | 1.09       | 1.86       | 4.54     | 2.67        | 0.63            | 858      | 1.86         |
| B/M-Mode                                  | Scan          | 1.09       | 1.86       | 4.54     | 2.67        | 0.63            | 5176     | 1.86         |
|                                           | Non-Scan      | 1.09       | 1.86       | 4.54     | 2.67        | 0.63            | 858      | 1.86         |
| Pulsed Wave Doppler                       | Non-Scan      | 0.57       | 1.07       | 6.43     | 2.65        | 0.64            | 651      | 1.07         |
| B/Pulsed Wave Doppler                     | Scan          | 1.09       | 1.86       | 4.54     | 2.67        | 0.63            | 5176     | 1.86         |
|                                           | Non-Scan      | 0.57       | 1.07       | 6.43     | 2.65        | 0.64            | 651      | 1.07         |
| B/Color Flow Mode                         | Scan          | 1.09       | 1.86       | 6.38     | 2.56        | 1.01            | 500      | 1.86         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 1.09       | 1.86       | 6.38     | 2.56        | 1.01            | 500      | 1.86         |
|                                           | Non-Scan      | 0.57       | 1.07       | 6.43     | 2.65        | 0.64            | 651      | 1.07         |
| M/M-CFM                                   | Non-Scan      | 1.09       | 1.86       | 6.38     | 2.56        | 1.01            | 1480     | 1.86         |
| M-CFM/PD                                  | Non-Scan      | 1.09       | 1.86       | 6.38     | 2.56        | 1.01            | 1480     | 1.86         |
| CWD Mode                                  | Non-Scan      | < 0.1      | < 0.1      | 3.99     | 2.49        | –               | –        | <0.1         |

Table 17–105. S317 Maximum Mechanical Index (MI)

**S611**

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_3(z_1)$<br>$I_{TA.3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|--------------------------------------|-------|----------|--------------------------|---------------|
| Pulsed Wave Doppler                          | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.96    | 5.0   | 40.5  | 8.7                                  | —     | 1.6      | 0.9                      | 7.5           |
| B/Pulsed Wave Doppler                        | Scan          | 0.29    | 4.4   | 13.9  | 5.1                                  | —     | 1.6      | 0.9                      | 3.7           |
|                                              | Non-Scan      | 0.96    | 5.0   | 40.5  | 8.7                                  | —     | 1.6      | 0.9                      | 4.5           |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.98    | 4.0   | 50.1  | 12.3                                 | —     | 1.6      | 0.9                      | 7.5           |
|                                              | Non-Scan      | 0.48    | 5.0   | 20.3  | 4.4                                  | —     | 1.6      | 0.9                      | 4.5           |
| Color Flow with<br>M Doppler                 | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 1.38    | 4.0   | 72.3  | 14.4                                 | —     | 1.6      | 0.9                      | 7.5           |
| CW Doppler                                   | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 1.21    | 4.0   | 63.3  | 33.0                                 | —     | 1.0      | 0.3                      | 6.4           |

Table 17–106. S611 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq@PII_{max}}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| Pulsed Wave Doppler                          | Scan          | —           | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 1.57        | 5.0   | 40.55 | 4.5      | 0.4              | 0.9                      | —                  |
| B/Pulsed Wave Doppler                        | Scan          | 0.29        | 4.4   | 13.92 | 3.7      | 0.7              | 0.9                      | —                  |
|                                              | Non-Scan      | 1.57        | 5.0   | 40.55 | 4.5      | 0.4              | 0.9                      | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.98        | 4.0   | 50.05 | 5.1      | 1.0              | 0.9                      | —                  |
|                                              | Non-Scan      | 0.78        | 5.0   | 20.27 | 4.5      | 0.4              | 0.9                      | —                  |
| Color Flow with<br>M Doppler                 | Scan          | —           | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 2.16        | 4.0   | 72.26 | 5.1      | 0.4              | 0.9                      | —                  |
| CW Doppler                                   | Scan          | —           | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 2.85        | 4.0   | 63.31 | 2.0      | 0.5              | 0.3                      | —                  |

Table 17–107. S611 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### S611 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| Pulsed Wave Doppler                          | 0.94        | 5.0   | 40.5  | 0.9                             | 4.5           |
| B/Pulsed Wave Doppler                        | 1.26        | 5.0   | 54.5  | 0.9                             | 4.5           |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 1.63        | 5.0   | 70.3  | 0.9                             | 4.5           |
| Color Flow with M Doppler                    | 1.67        | 4.0   | 72.3  | 0.9                             | 7.5           |
| CW Doppler                                   | 2.41        | 4.0   | 63.3  | 0.3                             | 6.4           |

Table 17–108. S611 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | fc [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | Pri3 @PIImax |
|-------------------------------------------|---------------|------------|------------|----------|----------|-----------------|----------|--------------|
| B–Mode                                    | Scan          | 0.74       | 1.43       | 3.71     | 4.43     | 0.42            | 8658     | 1.43         |
| M–Mode                                    | Non–Scan      | 0.74       | 1.43       | 3.71     | 4.43     | 0.42            | 854      | 1.43         |
| B/M–Mode                                  | Scan          | 0.74       | 1.43       | 3.71     | 4.43     | 0.42            | 8658     | 1.43         |
|                                           | Non–Scan      | 0.74       | 1.43       | 3.71     | 4.43     | 0.42            | 854      | 1.43         |
| Pulsed Wave Doppler                       | Non–Scan      | 0.72       | 1.47       | 4.47     | 4.98     | 3.00            | 1130     | 1.47         |
| B/Pulsed Wave Doppler                     | Scan          | 0.74       | 1.43       | 3.71     | 4.43     | 0.42            | 8658     | 1.43         |
|                                           | Non–Scan      | 0.72       | 1.47       | 4.47     | 4.98     | 3.00            | 1130     | 1.47         |
| B/Color Flow Mode                         | Scan          | 0.74       | 1.43       | 4.98     | 4.02     | 1.53            | 280      | 1.43         |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.74       | 1.43       | 4.98     | 4.02     | 1.53            | 280      | 1.43         |
|                                           | Non–Scan      | 0.72       | 1.47       | 4.47     | 4.98     | 3.00            | 1130     | 1.47         |
| M/M–CFM                                   | Non–Scan      | 0.74       | 1.43       | 4.98     | 4.02     | 1.53            | 280      | 1.43         |
| M–CFM/PD                                  | Non–Scan      | 0.74       | 1.43       | 4.98     | 4.02     | 1.53            | 280      | 1.43         |
| CWD Mode                                  | Non–Scan      | 0.05       | 0.17       | 1.98     | 4.00     | –               | –        | 0.17         |

Table 17–109. S611 Maximum Mechanical Index (MI)

T739

| Imaging Mode                                 | Scanning Mode | Max TIS | $f_c$ | $W_0$ | Min of $W_3(z_1)$<br>$I_{TA,3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|----------------------------------------------|---------------|---------|-------|-------|--------------------------------------|-------|----------|--------------------------|---------------|
| Pulsed Wave Doppler                          | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.28    | 5.1   | 11.8  | 5.5                                  | —     | 0.8      | 0.2                      | 1.6           |
| B/Pulsed Wave Doppler                        | Scan          | <0.1    | 7.0   | 2.6   | 1.0                                  | —     | 1.1      | 0.4                      | 2.1           |
|                                              | Non-Scan      | 0.28    | 5.1   | 11.8  | 5.5                                  | —     | 0.8      | 0.2                      | 1.6           |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.28    | 5.1   | 10.5  | 8.4                                  | —     | 0.8      | 0.2                      | 1.6           |
|                                              | Non-Scan      | 0.14    | 5.1   | 5.9   | 2.7                                  | —     | 0.8      | 0.2                      | 1.6           |
| Color Flow with<br>M Doppler                 | Scan          | —       | —     | —     | —                                    | —     | —        | —                        | —             |
|                                              | Non-Scan      | 0.39    | 5.1   | 15.7  | 14.8                                 | —     | 0.8      | 0.2                      | 1.6           |

Table 17–110. I739 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode                                 | Scanning Mode | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII_{max}$ |
|----------------------------------------------|---------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| Pulsed Wave Doppler                          | Scan          | —           | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 0.98        | 5.1   | 11.8  | 1.7      | 0.2              | 0.2                      | —                  |
| B/Pulsed Wave Doppler                        | Scan          | <0.1        | 7.0   | 2.6   | 2.1      | 0.9              | 0.4                      | —                  |
|                                              | Non-Scan      | 0.98        | 5.1   | 11.8  | 1.7      | 0.2              | 0.2                      | —                  |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | Scan          | 0.28        | 5.1   | 10.5  | 1.8      | 0.8              | 0.2                      | —                  |
|                                              | Non-Scan      | 0.49        | 5.1   | 5.9   | 1.7      | 0.2              | 0.2                      | —                  |
| Color Flow with<br>M Doppler                 | Scan          | —           | —     | —     | —        | —                | —                        | —                  |
|                                              | Non-Scan      | 0.89        | 5.1   | 15.7  | 1.8      | 0.4              | 0.2                      | —                  |

Table 17–111. I739 Maximum Bone Thermal Index (TIB)



## FDA Acoustic Output Data

### T739 (cont'd)

| Imaging Mode                                 | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{\text{aprt}}$ | FL Focal Zone |
|----------------------------------------------|-------------|-------|-------|---------------------------------|---------------|
| Pulsed Wave Doppler                          | 0.56        | 5.1   | 11.8  | 0.2                             | 1.6           |
| B/Pulsed Wave Doppler                        | 0.64        | 5.1   | 14.4  | 0.2                             | 1.6           |
| B/Color Flow Doppler/<br>Pulsed Wave Doppler | 0.74        | 5.1   | 16.3  | 0.2                             | 1.6           |
| Color Flow with M Doppler                    | 0.74        | 5.1   | 15.7  | 0.2                             | 1.6           |

Table 17–112. T739 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode                              | Scanning Mode | Maximum MI | Pr.3 [MPa] | Zsp [cm] | fc [MHz] | Pd [ $\mu$ sec] | PRF [Hz] | PrI3 @PII <sub>max</sub> |
|-------------------------------------------|---------------|------------|------------|----------|----------|-----------------|----------|--------------------------|
| B–Mode                                    | Scan          | 0.72       | 1.97       | 2.14     | 6.96     | 0.18            | 8658     | 1.97                     |
| M–Mode                                    | Non–Scan      | 0.72       | 1.98       | 2.14     | 6.96     | 0.18            | 866      | 1.98                     |
| B/M–Mode                                  | Scan          | 0.72       | 1.97       | 2.14     | 6.96     | 0.18            | 8658     | 1.97                     |
|                                           | Non–Scan      | 0.72       | 1.98       | 2.14     | 6.96     | 0.18            | 866      | 1.98                     |
| Pulsed Wave Doppler                       | Non–Scan      | 0.71       | 1.65       | 2.34     | 6.63     | 2.18            | 650      | 1.65                     |
| B/Pulsed Wave Doppler                     | Scan          | 0.72       | 1.97       | 2.14     | 6.96     | 0.18            | 8658     | 1.97                     |
|                                           | Non–Scan      | 0.71       | 1.65       | 2.34     | 6.63     | 2.18            | 650      | 1.65                     |
| B/Color Flow Mode                         | Scan          | 0.93       | 1.97       | 1.84     | 5.08     | 0.52            | 280      | 1.97                     |
| B/Color Flow Mode/<br>Pulsed Wave Doppler | Scan          | 0.93       | 1.97       | 1.84     | 5.08     | 0.52            | 280      | 1.97                     |
|                                           | Non–Scan      | 0.71       | 1.65       | 2.34     | 6.63     | 2.18            | 650      | 1.65                     |
| M/M–CFM                                   | Non–Scan      | 0.93       | 1.98       | 1.84     | 5.08     | 0.52            | 280      | 1.98                     |
| M–CFM/PD                                  | Non–Scan      | 0.93       | 1.98       | 1.84     | 5.08     | 0.52            | 280      | 1.98                     |
| CWD Mode                                  | Non–Scan      | –          | –          | –        | –        | –               | –        | –                        |

Table 17–113. T739 Maximum Mechanical Index (MI)



**CWD2**

| Imaging Mode | Scanning Mode    | Max TIS | $f_c$ | $W_0$ | Min of $W_3(z_1)$<br>$I_{TA.3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|--------------|------------------|---------|-------|-------|--------------------------------------|-------|----------|--------------------------|---------------|
| CW Doppler   | Scan<br>Non-Scan | 0.49    | 2.00  | 51.19 | 44.56                                | 1.59  | 1.59     | 0.88                     | 1.00          |

Table 17–114. CWD2 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode | Scanning Mode    | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII_{max}$ |
|--------------|------------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| CW Doppler   | Scan<br>Non-Scan | 3.30        | 2.00  | 51.19 | 1.00     | 7.22             | 0.88                     | 7.22               |

Table 17–115. CWD2 Maximum Bone Thermal Index (TIB)

| Imaging Mode | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|--------------|-------------|-------|-------|--------------------------|---------------|
| CW Doppler   | 1.21        | 2.00  | 3.48  | 0.88                     | 1.00          |

Table 17–116. CWD2 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode | Scanning Mode | Maximum MI | Pr.3 [MPa] | $Z_{sp}$ [cm] | $f_c$ [MHz] | $P_d$ [ $\mu$ sec] | PRF [Hz] | PrI3 @PII <sub>max</sub> |
|--------------|---------------|------------|------------|---------------|-------------|--------------------|----------|--------------------------|
| CWD Mode     | Non-Scan      | 0.10       | 0.16       | 1.00          | 2.01        | –                  | –        | 0.16                     |

Table 17–117. CWD2 Maximum Mechanical Index (MI)



## FDA Acoustic Output Data

---

### CWD5

| Imaging Mode | Scanning Mode    | Max TIS | $f_c$ | $W_0$ | Min of $W_{.3}(z_1)$<br>$I_{TA.3}(z_1)$ | $z_1$ | $z_{bp}$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|--------------|------------------|---------|-------|-------|-----------------------------------------|-------|----------|--------------------------|---------------|
| CW Doppler   | Scan<br>Non-Scan | 1.92    | 5.00  | 80.41 | 56.90                                   | 0.85  | 0.85     | 0.25                     | 1.00          |

Table 17–118. CWD5 Maximum Soft Tissue Thermal Index (TIS)

| Imaging Mode | Scanning Mode    | Maximum TIB | $f_c$ | $W_0$ | $z_{sp}$ | $d_{eq}(z_{sp})$ | Dimensions of $A_{aprt}$ | $d_{eq}@PII_{max}$ |
|--------------|------------------|-------------|-------|-------|----------|------------------|--------------------------|--------------------|
| CW Doppler   | Scan<br>Non-Scan | 3.55        | 5.00  | 80.41 | 1.00     | 2.39             | 0.25                     | 2.39               |

Table 17–119. CWD5 Maximum Bone Thermal Index (TIB)

| Imaging Mode | Maximum TIC | $f_c$ | $W_0$ | Dimensions of $A_{aprt}$ | FL Focal Zone |
|--------------|-------------|-------|-------|--------------------------|---------------|
| CW Doppler   | 3.56        | 5.00  | 80.41 | 0.25                     | 1.00          |

Table 17–120. CWD5 Maximum Cranial Bone Thermal Index (TIC)

| Imaging Mode | Scanning Mode | Maximum MI | Pr.3 [MPa] | $Z_{sp}$ [cm] | $f_c$ [MHz] | $P_d$ [ $\mu$ sec] | PRF [Hz] | PrI3 @PII <sub>max</sub> |
|--------------|---------------|------------|------------|---------------|-------------|--------------------|----------|--------------------------|
| CWD Mode     | Non-Scan      | 0.06       | 0.19       | 1.00          | 5.01        | —                  | —        | 0.19                     |

Table 17–121. CWD5 Maximum Mechanical Index (MI)

## Key to Tables

| Symbol                                            | Units              | Description                                                                                                                           |
|---------------------------------------------------|--------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| <b>Output Display Indices:</b>                    |                    |                                                                                                                                       |
| MI                                                | n/a                | Mechanical Index                                                                                                                      |
| TIS <sub>scan</sub>                               | n/a                | Soft Tissue Thermal Index (scanning mode)                                                                                             |
| TIS <sub>non-scan</sub>                           | n/a                | Soft Tissue Thermal Index (non-scanning mode)                                                                                         |
| TIB                                               | n/a                | Bone Thermal Index                                                                                                                    |
| TIC                                               | n/a                | Cranial Thermal Index                                                                                                                 |
| <b>Intensity, Power, and Pressure Parameters:</b> |                    |                                                                                                                                       |
| I <sub>SPTA,3</sub>                               | mW/cm <sup>2</sup> | Spatial Peak-Temporal Average Intensity derated to 0.3 dB/cm/MHz                                                                      |
| I <sub>SPPA,3</sub>                               | W/cm <sup>2</sup>  | Spatial Peak-Pulse Average Intensity derated to 0.3 dB/cm/MHz                                                                         |
| I <sub>TA,3</sub>                                 | mW/cm <sup>2</sup> | Time average Intensity, derated to 0.3 dB/cm/MHz                                                                                      |
| PII <sub>max</sub>                                | n/a                | Maximum free-field pulse intensity integral                                                                                           |
| W <sub>0</sub>                                    | mW                 | Ultrasonic Power                                                                                                                      |
| W <sub>.3(z<sub>1</sub>)</sub>                    | mW                 | Derated ultrasonic power at axial distance z <sub>1</sub>                                                                             |
| P <sub>r</sub>                                    | MPa                | Peak rarefactional pressure                                                                                                           |
| P <sub>r,3</sub>                                  | MPa                | Peak rarefactional pressure, derated to 0.3 dB/cm/MHz                                                                                 |
| <b>Transmit Parameters:</b>                       |                    |                                                                                                                                       |
| f <sub>c</sub>                                    | MHz                | Center Frequency                                                                                                                      |
| PD                                                | μs                 | Pulse Duration                                                                                                                        |
| PRF                                               | kHz                | Pulse Repetition Frequency                                                                                                            |
| <b>Distance/Depth and Area Parameters:</b>        |                    |                                                                                                                                       |
| d <sub>eq</sub> (z)                               | cm                 | Equivalent beam diameter as a function of z, = $[(4/\pi)(W_0/I_{TA}(z))]^{1/2}$                                                       |
| z <sub>sp</sub>                                   | cm                 | Axial distance to spatial peak                                                                                                        |
| z <sub>1</sub>                                    | cm                 | Axial distance to the location of max[min(W <sub>.3</sub> (z), I <sub>TA,3</sub> (z) × 1 cm <sup>2</sup> )] where z ≥ z <sub>bp</sub> |
| z <sub>bp</sub>                                   | cm                 | Break point depth = $1.69(A_{aprt})^{1/2}$                                                                                            |
| ROC                                               | cm                 | Radius of curvature                                                                                                                   |
| FL                                                | cm                 | Focal Length                                                                                                                          |
| A <sub>aprt</sub>                                 | cm <sup>2</sup>    | Area of Active Aperture                                                                                                               |

Table 17–122. Index Tables Key



## Measurement Precision and Uncertainty

Precision and measurement uncertainty of the system used to measure the acoustic output is provided in the following table. These values are determined in accordance with Section 6.4 of the *Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment*, AIUM/NEMA 1992 and the *Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment*, NEMA 1992.

|                         | Center Frequency | Acoustic Power | Peak Rarefactional Pressure | Acoustic Intensity |
|-------------------------|------------------|----------------|-----------------------------|--------------------|
| Precision               | 4.4%             | 1.5%           | 2.6%                        | 5.1%               |
| Measurement Uncertainty | 10.0%            | 3.4%           | 6.0%                        | 11.6%              |

Table 17–123. Measurement Precision and Uncertainty

## Acoustic Output Display Operation and Accuracy

The Acoustic Output display shows MI (Mechanical Index) while in B-Mode. When other modes are selected TIS, TIB or TIC are displayed. The Setup/Custom Display menu page 9 will determine which value (TIS, TIB or TIC) is displayed.

Display increments are 0.2 for values of one or less and 1.0 for values greater than one. Display accuracy for MI is –15% to +38%. Display accuracy for TI is –37% to +95%. Index values less than 0.4 are not displayed.

## Endnotes

- <sup>1</sup> Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment, AIUM/NEMA, 1992
- <sup>2</sup> Implementation of the Principle of As Low As Reasonably Achievable (ALARA) for Medical and Dental Personnel, National Council on Radiation Protection and Measurements (NCRP), Report No.107, December 31, 1990.
- <sup>3</sup> FDA Center for Devices and Radiological Health (CDRH), 510(k) Guidance for Diagnostic Ultrasound and Fetal Doppler Ultrasound Medical Devices, September 8, 1989 draft.
- <sup>4</sup> FDA Center for Devices and Radiological Health (CDRH), 510(k) Guidance for Diagnostic Ultrasound and Fetal Doppler Ultrasound Medical Devices, September 8, 1989 draft.
- <sup>5</sup> Rochelson BL, Schulman H, Fleischer A, Farmakides G, Bracero L, Ducey J, Winter D, Penny B: The clinical significance of Doppler umbilical artery velocimetry in the small for gestational age fetus. *Am J Obstet Gynecol* 156:1223-1226, 1987.
- <sup>6</sup> Soothill PW, Nicolaides KH, Bilardo CM, Campbell S: Relation of fetal hypoxia in growth retardation to mean blood velocity in the fetal aorta. *Lancet* 2:1118-1119, 1986.
- <sup>7</sup> Jouppila P, Kirkinen P: Increased vascular resistance in the descending aorta of the human fetus in hypoxia. *R J Obstet Gynecol* 91:853-856, 1984.
- <sup>8</sup> Reuwer PJHM, Rietman GW, Sijmons EA, et al: Intrauterine growth retardation: Prediction of perinatal distress by Doppler ultrasound. *Lancet* 2:415-418, 1987.
- <sup>9</sup> Copel JA, Grannum PA, Hobbins JC, Cunningham FG: Doppler Ultrasound in Obstetrics. In: Pritchard, MacDonald, GNT (eds) Williams Obstetrics Seventeenth Edition, Philadelphia, PA, Appleton & Lange, 1988, Supplement No. 16.
- <sup>10</sup> Rochelson B, Schulman H, Farmakides G, et al: The significance of absent end-diastolic velocity in umbilical artery velocity waveforms. *Am J Obstet Gynecol* 156:1213, 1987.



**Endnotes (cont'd)**

- <sup>11</sup> Woo JSK, Liang ST, Lo RLS: Significance of an absent or reversed end diastolic flow in Doppler umbilical waveforms. *J Ultrasound Med* 6:291, 1987.
- <sup>12</sup> Kleinman CS, Weinstein EM, Copel JA: Pulsed Doppler Analysis of human fetal blood flow. *Clin Diagnostic Ultrasound* 17:173-185, 1986.
- <sup>13</sup> Strasburger JF, Huhta JC, Carpenter RJ, et al: Doppler echocardiography in the diagnosis and management of persistent fetal arrhythmias. *J Am Coll Cardiol* 7:1386-1391, 1986.
- <sup>14</sup> Steinfeld L, Rappaport HL, Rossbach HC, Martinez E: Diagnosis of fetal arrhythmias using echocardiographic and Doppler techniques. *J Am Coll Cardiol* 9:1425-1433, 1986.
- <sup>15</sup> Reed KL, Sahn DJ, Marx GR, et al: Cardiac Doppler flows during fetal arrhythmias: Physiologic consequences. *Obstet Gynecol* 70:1-6, 1987.
- <sup>16</sup> Silverman NH, Kleinman CS, Rudolph AM, et al: Fetal atrioventricular valve insufficiency associated with non-immune hydrops: a two-dimensional echocardiographic and pulsed Doppler ultrasound study. *Circulation* 72:825-832, 1985.
- <sup>17</sup> FDA/CDRH, 510(k) Diagnostic Ultrasound Guidance Update of 1991, April 26, 1991 draft.
- <sup>18</sup> Biological Effects of Ultrasound: Mechanisms and Clinical Implications, NCRP Report No. 74, December 30, 1983.
- <sup>19</sup> Exposure Criteria for Medical Diagnostic Ultrasound: I. Criteria Based on Thermal Mechanisms, NCRP Report No. 113, June 1, 1992.
- <sup>20</sup> Bioeffects Considerations for the Safety of Diagnostic Ultrasound, *Journal of Ultrasound in Medicine*, AIUM, September 1988.
- <sup>21</sup> Geneva Report on Safety and Standardization in Medical Ultrasound, WFUMB, May 1990.

# IEC Acoustic Output Tables

Acoustical parameters represent the maximum values for a probe/mode combination; other parameters refer to the operating conditions which yield these maximum acoustic parameters.

## Key to Tables

| Parameter                | Unit               | Description                                                                                                                                                                                                                                        |
|--------------------------|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| $P_-$                    | (MPa)              | Peak-negative acoustic pressure in the plane perpendicular to the beam-alignment axis containing the maximum pulse-pressure-squared integral (or maximum mean square acoustic pressure for continuous wave systems) in the whole ultrasonic field. |
| $I_{spta}$               | mW/cm <sup>2</sup> | Spatial-peak temporal-average derived intensity in the whole ultrasonic field.                                                                                                                                                                     |
| $I_p$                    | mm                 | Distance from the transducer output face to the point of maximum pulse-pressure-squared integral (or maximum mean square acoustic pressure for continuous wave systems).                                                                           |
| $W_{pb6}$<br>(  )<br>(⊥) | mm                 | −6 dB pulse beam-width at the point of maximum pulse-pressure-squared integral (or maximum mean square acoustic pressure for continuous wave systems). These directions shall be parallel (  ) and perpendicular (⊥) to the reference direction.   |
| $pr$                     | KHz                | Pulse repetition rate for non-scanning modes.                                                                                                                                                                                                      |
| $sr$                     | Hz                 | Scan repetition rate (srr) for scanning modes.                                                                                                                                                                                                     |
| Output Beam Dimensions   | mm                 | Dimensions parallel and perpendicular to the reference direction.                                                                                                                                                                                  |
| $f_{awf}$                | MHz                | Arithmetic-mean acoustic-working frequency measured by a hydrophone placed at the point of maximum pulse-pressure-squared integral (or maximum mean square acoustic pressure for continuous wave systems).                                         |
| APF                      | %                  | Acoustic power-up fraction.                                                                                                                                                                                                                        |
| AIF                      | %                  | Acoustic initialization fraction.                                                                                                                                                                                                                  |

Table 17–124. Key to IEC Table Value



## IEC Acoustic Output Tables

---

### Key to Tables (cont'd)

| Parameter              | Unit               | Description                                                                                                                                                                  |
|------------------------|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Maximum Power          | mW                 | Maximum temporal-average power output. For scanning modes, this shall be the total power output of all the acoustic pulses.                                                  |
| $I_{ob}$               | mW/cm <sup>2</sup> | Output beam intensity.                                                                                                                                                       |
| Power-up Mode          |                    | In systems in which the user defines the power-up mode, this shall be stated as either "user defined" or "not applicable" (n/a).                                             |
| Initialization Mode    |                    | In systems in which the user defines the initialization mode, this shall be stated as either "user defined" or "not applicable" (n/a).                                       |
| Acoustic Output Freeze |                    | If the system has acoustic output freeze, then this shall be stated as "yes;" otherwise, it shall be stated as "no."                                                         |
| $I_{tt}$               | mm                 | Transducer to transducer-output-face distance.                                                                                                                               |
| $I_{ts}$               | mm                 | Typical value for the transducer stand-off distance. If the transducer assembly is normally used in contact with the patient, then this shall be specified "contact" system. |

Table 17–124. Key to IEC Table Value (cont'd)



C358 Probe

| Parameter                                                     | B            | M           | B+M         | D            | B+D          | B+rD         | B+rD<br>+D   | M+cM         | cM+D         |
|---------------------------------------------------------------|--------------|-------------|-------------|--------------|--------------|--------------|--------------|--------------|--------------|
| <b>P<sub>-</sub> (MPa)</b>                                    | 2.57         | 2.57        | 2.57        | 3.25         | 3.25         | 3.11         | 3.25         | 3.09         | 3.09         |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                   | 315.6        | 545.3       | 755.6       | 3354.6       | 3565.0       | 319.8        | 1942.4       | 2490.4       | 2490.4       |
| <b>I<sub>p</sub> (mm)</b>                                     | 39.2         | 39.2        | 39.2        | 70.5         | 70.5         | 69.7         | 69.7         | 69.7         | 69.7         |
| <b>W<sub>pb6</sub></b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>  | 2.4<br>3.6   | 2.4<br>3.6  | 2.4<br>3.6  | 1.4<br>3.5   | 1.4<br>3.5   | 1.3<br>3.7   | 1.3<br>3.7   | 1.3<br>3.7   | 1.3<br>3.7   |
| <b>prf (kHz)</b>                                              | –            | 0.92        | 0.92        | 18.12        | 18.12        | 2.26         | 2.26         | 2.26         | 2.26         |
| <b>srr (Hz)</b>                                               | 40           | –           | 40          | –            | 40           | 7            | 7            | –            | –            |
| <b>Output Beam Dim.</b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b> | 10.2<br>9.2  | 10.2<br>9.2 | 10.2<br>9.2 | 10.2<br>25.8 | 10.2<br>25.8 | 10.2<br>25.8 | 10.2<br>25.8 | 10.2<br>25.8 | 10.2<br>25.8 |
| <b>f<sub>awf</sub> (MHz)</b>                                  | 2.62         | 2.62        | 2.62        | 3.32         | 3.32         | 3.31         | 3.31         | 3.31         | 3.31         |
| <b>APF (%)</b><br><b>(Acoustic Power up Fraction)</b>         | na           | na          | na          | na           | na           | na           | na           | na           | na           |
| <b>AIF (%)</b><br><b>(Acoustic Initialization Fraction)</b>   | na           | na          | na          | na           | na           | na           | na           | na           | na           |
| <b>Maximum Power (mW)</b>                                     | 93.26        | 23.30       | 85.47       | 99.40        | 161.57       | 138.40       | 149.99       | 187.58       | 187.58       |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                     | 99.38        | 24.83       | 91.08       | 37.83        | 104.08       | 95.26        | 99.67        | 87.35        | 87.35        |
| <b>Power-up Mode</b>                                          | na           | na          | na          | na           | na           | na           | na           | na           | na           |
| <b>Initialization Mode</b>                                    | na           | na          | na          | na           | na           | na           | na           | na           | na           |
| <b>Acoustic Output Freeze</b>                                 | yes          | yes         | yes         | yes          | yes          | yes          | yes          | yes          | yes          |
| <b>I<sub>tt</sub> (mm)</b>                                    | 0.00         | 0.00        | 0.00        | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         |
| <b>I<sub>ts</sub> (mm)</b>                                    | con-<br>tact | contact     | contact     | contact      | contact      | contact      | contact      | contact      | contact      |

Table 17–125. C358 IEC Acoustic Output Information



## IEC Acoustic Output Tables

### C364 Probe

| Parameter                                                                       | B       | M       | B+M     | D       | B+D     | B+rD    | B+rD<br>+D | M+cM    | cM+D    |
|---------------------------------------------------------------------------------|---------|---------|---------|---------|---------|---------|------------|---------|---------|
| <b>P<sub>-</sub> (MPa)</b>                                                      | 1.89    | 1.89    | 1.89    | 1.90    | 1.90    | 1.89    | 1.90       | 1.89    | 1.89    |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                                     | 75.0    | 177.0   | 227.0   | 851.2   | 901.2   | 254.0   | 577.6      | 1240.6  | 851.2   |
| <b>I<sub>p</sub> (mm)</b>                                                       | 46.1    | 45.7    | 46.1    | 19.1    | 19.1    | 19.7    | 19.7       | 19.7    | 19.7    |
| <b>W<sub>pb6</sub></b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>                    | 1.9     | 1.8     | 1.9     | 1.9     | 1.9     | 1.9     | 1.9        | 1.9     | 1.9     |
|                                                                                 | 2.5     | 2.5     | 2.5     | 7.1     | 7.1     | 6.4     | 6.4        | 6.4     | 6.4     |
| <b>pr<sub>r</sub> (kHz)</b>                                                     | –       | 0.89    | 0.89    | 15.72   | 15.72   | 10.40   | 10.40      | 10.40   | 10.40   |
| <b>s<sub>rr</sub> (Hz)</b>                                                      | 70      | –       | 70      | –       | 70      | 9       | 9          | –       | –       |
| <b>Output Beam</b><br><b>Dim. (  ) (mm)</b><br><b>(⊥)(mm)</b>                   | 10.2    | 10.2    | 10.2    | 10.2    | 10.2    | 10.2    | 10.2       | 10.2    | 10.2    |
|                                                                                 | 12.0    | 12.0    | 12.0    | 5.0     | 5.0     | 5.0     | 5.0        | 5.0     | 5.0     |
| <b>f<sub>awf</sub> (MHz)</b>                                                    | 3.36    | 3.34    | 3.36    | 3.29    | 3.29    | 3.31    | 3.31       | 3.31    | 3.31    |
| <b>APF (%)</b><br><b>(Acoustic Power</b><br><b>up Fraction)</b>                 | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>AIF (%)</b><br><b>(Acoustic</b><br><b>Initialization</b><br><b>Fraction)</b> | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Maximum Power</b><br><b>(mW)</b>                                             | 16.97   | 4.36    | 15.68   | 66.92   | 78.23   | 63.73   | 70.98      | 24.24   | 24.24   |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                                       | 13.86   | 3.57    | 12.81   | 131.22  | 140.46  | 112.01  | 126.24     | 42.53   | 42.53   |
| <b>Power-up Mode</b>                                                            | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Initialization</b><br><b>Mode</b>                                            | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Acoustic Output</b><br><b>Freeze</b>                                         | yes     | yes     | yes     | yes     | yes     | yes     | yes        | yes     | yes     |
| <b>I<sub>tt</sub> (mm)</b>                                                      | 0.00    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00       | 0.00    | 0.00    |
| <b>I<sub>ts</sub> (mm)</b>                                                      | contact | contact | contact | contact | contact | contact | contact    | contact | contact |

Table 17–126. C364 IEC Acoustic Output Information

C551 Probe

| Parameter                                                                       | B       | M       | B+M     | D       | B+D     | B+rD    | B+rD<br>+D | M+cM    | cM+D    |
|---------------------------------------------------------------------------------|---------|---------|---------|---------|---------|---------|------------|---------|---------|
| <b>P<sub>-</sub> (MPa)</b>                                                      | 2.48    | 2.48    | 2.48    | 3.64    | 3.64    | 4.15    | 4.15       | 2.48    | 2.48    |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                                     | 132.9   | 172.5   | 261.1   | 3219.5  | 3308.1  | 329.4   | 1818.7     | 2029.3  | 2029.3  |
| <b>I<sub>p</sub> (mm)</b>                                                       | 38.1    | 37.9    | 38.1    | 44.1    | 44.1    | 43.9    | 43.9       | 43.9    | 43.9    |
| <b>W<sub>pb6</sub></b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>                    | 1.9     | 1.7     | 1.9     | 0.9     | 0.9     | 1.1     | 1.1        | 1.1     | 1.1     |
|                                                                                 | 2.9     | 3.7     | 2.9     | 2.9     | 2.9     | 2.2     | 2.2        | 2.2     | 2.2     |
| <b>pr<sub>r</sub> (kHz)</b>                                                     | –       | 1.04    | 1.04    | 0.98    | 0.98    | 6.84    | 6.84       | 6.84    | 6.84    |
| <b>s<sub>r</sub> (Hz)</b>                                                       | 40      | –       | 40      | –       | 40      | 30      | 30         | –       | –       |
| <b>Output Beam</b><br><b>Dim. (  ) (mm)</b><br><b>(⊥)(mm)</b>                   | 9.0     | 9.0     | 9.0     | 9.0     | 9.0     | 9.0     | 9.0        | 9.0     | 9.0     |
|                                                                                 | 8.8     | 8.8     | 8.8     | 16.8    | 16.8    | 16.8    | 16.8       | 16.8    | 16.8    |
| <b>f<sub>awf</sub> (MHz)</b>                                                    | 3.70    | 4.29    | 3.70    | 4.97    | 4.97    | 4.02    | 4.02       | 4.02    | 4.02    |
| <b>APF (%)</b><br><b>(Acoustic Power</b><br><b>up Fraction)</b>                 | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>AIF (%)</b><br><b>(Acoustic</b><br><b>Initialization</b><br><b>Fraction)</b> | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Maximum Power</b><br><b>(mW)</b>                                             | 37.06   | 12.72   | 37.43   | 101.60  | 126.31  | 67.42   | 96.87      | 32.83   | 32.83   |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                                       | 46.79   | 16.07   | 47.26   | 67.20   | 98.39   | 59.45   | 78.92      | 29.36   | 29.36   |
| <b>Power-up Mode</b>                                                            | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Initialization</b><br><b>Mode</b>                                            | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Acoustic Output</b><br><b>Freeze</b>                                         | yes     | yes     | yes     | yes     | yes     | yes     | yes        | yes     | yes     |
| <b>I<sub>tt</sub> (mm)</b>                                                      | 0.00    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00       | 0.00    | 0.00    |
| <b>I<sub>ts</sub> (mm)</b>                                                      | contact | contact | contact | contact | contact | contact | contact    | contact | contact |

Table 17–127. C551 IEC Acoustic Output Information



## IEC Acoustic Output Tables

### C386 Probe

| Parameter                                                     | B              | M              | B+M            | D              | B+D            | B+rD           | B+rD<br>+D     | M+cM           | cM+D           |
|---------------------------------------------------------------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|----------------|
| <b>P<sub>-</sub> (MPa)</b>                                    | 2.82           | 2.82           | 2.82           | 2.67           | 2.82           | 2.82           | 2.82           | 2.19           | 21.9           |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                   | 28.76          | 330.33         | 349.51         | 1938.3<br>0    | 1957.4<br>8    | 421.78         | 1189.6<br>3    | 523.13         | 523.13         |
| <b>I<sub>p</sub> (mm)</b>                                     | 70.30          | 70.30          | 70.30          | 73.20          | 73.20          | 60.50          | 60.50          | 60.50          | 60.50          |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)                | 2.11<br>1.64   | 2.11<br>1.64   | 2.11<br>1.64   | 2.46<br>1.75   | 2.46<br>1.75   | 2.77<br>2.19   | 2.77<br>2.19   | 2.77<br>2.19   | 2.77<br>2.19   |
| <b>pr<sub>r</sub> (kHz)</b>                                   | –              | 0.87           | 0.87           | 7.86           | 7.86           | 0.42           | 0.42           | 0.42           | 0.42           |
| <b>s<sub>r</sub> (Hz)</b>                                     | 29.00          | –              | 29.00          | –              | 29.00          | 7.00           | 7.00           | –              | –              |
| <b>Output Beam<br/>Dim.</b> (  ) (mm)<br>(⊥)(mm)              | 17.10<br>13.00 | 17.10<br>13.00 | 17.10<br>13.00 | 19.80<br>13.00 | 17.10<br>13.00 | 17.10<br>13.00 | 17.10<br>13.00 | 19.80<br>13.00 | 19.80<br>13.00 |
| <b>f<sub>awf</sub> (MHz)</b>                                  | 3.75           | 3.75           | 3.75           | 3.36           | 3.36           | 2.55           | 2.55           | 2.55           | 2.55           |
| <b>APF (%)<br/>(Acoustic Power<br/>up Fraction)</b>           | na             | na             | na             | na             | na             | na             | na             | na             | na             |
| <b>AIF (%)<br/>(Acoustic<br/>Initialization<br/>Fraction)</b> | na             | na             | na             | na             | na             | na             | na             | na             | na             |
| <b>Maximum Power<br/>(mW)</b>                                 | 42.09          | 6.68           | 34.74          | 43.86          | 71.92          | 128.68         | 100.30         | 100.62         | 100.62         |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                     | 18.96          | 3.01           | 15.65          | 17.19          | 29.83          | 63.23          | 46.53          | 53.60          | 53.60          |
| <b>Power-up Mode</b>                                          | na             | na             | na             | na             | na             | na             | na             | na             | na             |
| <b>Initialization<br/>Mode</b>                                | na             | na             | na             | na             | na             | na             | na             | na             | na             |
| <b>Acoustic Output<br/>Freeze</b>                             | yes            | yes            | yes            | yes            | yes            | yes            | yes            | yes            | yes            |
| <b>I<sub>tt</sub> (mm)</b>                                    | 0.00           | 0.00           | 0.00           | 0.00           | 0.00           | 0.00           | 0.00           | 0.00           | 0.00           |
| <b>I<sub>ts</sub> (mm)</b>                                    | contact        | contact        | contact        | contact        | contact        | contact        | contact        | contact        | contact        |

Table 17–128. C386 IEC Acoustic Output Information

C721 Probe

| Parameter                                                                       | B       | M       | B+M     | D       | B+D     | B+rD    | B+rD<br>+D | M+cM    | cM+D    |
|---------------------------------------------------------------------------------|---------|---------|---------|---------|---------|---------|------------|---------|---------|
| <b>P<sub>-</sub> (MPa)</b>                                                      | 1.70    | 1.15    | 1.70    | 1.66    | 1.70    | 1.94    | 1.94       | 1.20    | 1.20    |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                                     | 44.32   | 54.41   | 83.96   | 453.84  | 483.39  | 156.46  | 319.92     | 262.89  | 262.89  |
| <b>I<sub>p</sub> (mm)</b>                                                       | 35.20   | 35.20   | 35.20   | 36.10   | 36.10   | 31.00   | 31.00      | 31.00   | 31.00   |
| <b>W<sub>pb6</sub></b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>                    | 1.26    | 1.26    | 1.26    | 1.43    | 1.43    | 1.63    | 1.63       | 1.63    | 1.63    |
|                                                                                 | 1.79    | 1.79    | 1.79    | 1.69    | 1.69    | 2.15    | 2.15       | 2.15    | 2.15    |
| <b>pr<sub>r</sub> (kHz)</b>                                                     | –       | 0.90    | 0.90    | 5.95    | 5.95    | 1.96    | 1.96       | 1.96    | 1.96    |
| <b>s<sub>r</sub> (Hz)</b>                                                       | 34.00   | –       | 34.00   | –       | 34.00   | 9.00    | 9.00       | –       | –       |
| <b>Output Beam</b><br><b>Dim. (  ) (mm)</b><br><b>(⊥)(mm)</b>                   | 10.29   | 10.29   | 10.29   | 7.64    | 10.29   | 10.29   | 10.29      | 7.64    | 7.64    |
|                                                                                 | 8.00    | 8.00    | 8.00    | 8.00    | 8.00    | 8.00    | 8.00       | 8.00    | 8.00    |
| <b>f<sub>awf</sub> (MHz)</b>                                                    | 5.25    | 5.25    | 5.25    | 6.63    | 6.63    | 5.07    | 5.07       | 5.07    | 5.07    |
| <b>APF (%)</b><br><b>(Acoustic Power</b><br><b>up Fraction)</b>                 | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>AIF (%)</b><br><b>(Acoustic</b><br><b>Initialization</b><br><b>Fraction)</b> | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Maximum Power</b><br><b>(mW)</b>                                             | 13.83   | 1.38    | 10.60   | 9.80    | 19.02   | 10.98   | 15.00      | 25.80   | 25.80   |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                                       | 16.80   | 1.68    | 12.88   | 16.03   | 27.24   | 53.43   | 40.33      | 42.23   | 42.23   |
| <b>Power-up Mode</b>                                                            | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Initialization</b><br><b>Mode</b>                                            | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Acoustic Output</b><br><b>Freeze</b>                                         | yes     | yes     | yes     | yes     | yes     | yes     | yes        | yes     | yes     |
| <b>I<sub>tt</sub> (mm)</b>                                                      | 0.00    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00       | 0.00    | 0.00    |
| <b>I<sub>ts</sub> (mm)</b>                                                      | contact | contact | contact | contact | contact | contact | contact    | contact | contact |

Table 17–129. C721 IEC Acoustic Output Information



## IEC Acoustic Output Tables

### E721 Probe

| Parameter                                                  | B           | M           | B+M         | D           | B+D         | B+rD       | B+rD<br>+D | M+cM       | cM+D       |
|------------------------------------------------------------|-------------|-------------|-------------|-------------|-------------|------------|------------|------------|------------|
| <b>P<sub>-</sub> (MPa)</b>                                 | 3.27        | 3.27        | 3.27        | 2.99        | 3.27        | 3.27       | 3.27       | 3.27       | 3.27       |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                | 124.7       | 128.7       | 211.8       | 2189.7      | 2272.8      | 125.5      | 1199.2     | 564.4      | 564.4      |
| <b>I<sub>p</sub> (mm)</b>                                  | 37.0        | 34.4        | 37.0        | 37.5        | 37.5        | 37.4       | 37.4       | 37.4       | 37.4       |
| <b>W<sub>pb6</sub></b><br>(∥) (mm)<br>(⊥)(mm)              | 1.4<br>2.1  | 1.1<br>2.9  | 1.4<br>2.1  | 1.2<br>1.8  | 1.2<br>1.8  | 1.4<br>2.0 | 1.4<br>2.0 | 1.4<br>2.0 | 1.4<br>2.0 |
| <b>pr<sub>r</sub> (kHz)</b>                                | –           | 0.87        | 0.87        | 3.93        | 3.93        | 0.28       | 0.28       | 1.48       | 1.48       |
| <b>s<sub>r</sub> (Hz)</b>                                  | 180         | –           | 180         | –           | 180         | 30         | 30         | –          | –          |
| <b>Output Beam Dim.</b><br>(∥) (mm)<br>(⊥)(mm)             | 8.0<br>10.3 | 8.0<br>10.3 | 8.0<br>10.3 | 8.0<br>10.3 | 8.0<br>10.3 | 8.0<br>9.6 | 8.0<br>9.6 | 8.0<br>9.6 | 8.0<br>9.6 |
| <b>f<sub>awf</sub> (MHz)</b>                               | 4.46        | 5.71        | 4.46        | 4.99        | 4.99        | 5.05       | 5.05       | 5.05       | 5.05       |
| <b>APF (%)</b><br>(Acoustic Power<br>up Fraction)          | na          | na          | na          | na          | na          | na         | na         | na         | na         |
| <b>AIF (%)</b><br>(Acoustic<br>Initialization<br>Fraction) | na          | na          | na          | na          | na          | na         | na         | na         | na         |
| <b>Maximum Power<br/>(mW)</b>                              | 11.47       | 3.82        | 11.46       | 28.56       | 36.21       | 16.97      | 26.59      | 30.09      | 30.09      |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                  | 13.93       | 4.64        | 13.92       | 34.69       | 43.98       | 21.40      | 32.69      | 38.74      | 38.74      |
| <b>Power-up Mode</b>                                       | na          | na          | na          | na          | na          | na         | na         | na         | na         |
| <b>Initialization<br/>Mode</b>                             | na          | na          | na          | na          | na          | na         | na         | na         | na         |
| <b>Acoustic Output<br/>Freeze</b>                          | yes         | yes         | yes         | yes         | yes         | yes        | yes        | yes        | yes        |
| <b>I<sub>tt</sub> (mm)</b>                                 | 0.00        | 0.00        | 0.00        | 0.00        | 0.00        | 0.00       | 0.00       | 0.00       | 0.00       |
| <b>I<sub>ts</sub> (mm)</b>                                 | contact     | contact     | contact     | contact     | contact     | contact    | contact    | contact    | contact    |

Table 17–130. E721 IEC Acoustic Output Information

## IEC Acoustic Output Tables

### 546L Probe

| Parameter                                                   | B          | M          | B+M        | D           | B+D         | B+rD         | B+rD<br>+D  | M+cM        | cM+D        |
|-------------------------------------------------------------|------------|------------|------------|-------------|-------------|--------------|-------------|-------------|-------------|
| <b>P<sub>-</sub> (MPa)</b>                                  | 3.50       | 3.50       | 3.50       | 3.72        | 3.72        | 4.14         | 4.14        | 3.50        | 3.50        |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                 | 233.8      | 325.4      | 481.3      | 1337.3      | 1493.2      | 766.3        | 1129.8      | 1306.4      | 1306.4      |
| <b>I<sub>p</sub> (mm)</b>                                   | 26.6       | 26.3       | 26.6       | 32.9        | 32.9        | 30.3         | 30.3        | 30.3        | 30.3        |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)              | 1.0<br>1.7 | 1.1<br>1.8 | 1.0<br>1.7 | 1.0<br>2.1  | 1.0<br>2.1  | 1.0<br>2.1   | 1.0<br>2.1  | 1.0<br>2.1  | 1.0<br>2.1  |
| <b>pr<sub>r</sub> (kHz)</b>                                 | –          | 0.91       | 0.91       | 0.65        | 0.65        | 4.54         | 4.54        | 4.54        | 4.54        |
| <b>s<sub>r</sub> (Hz)</b>                                   | 66         | –          | 66         | –           | 66          | 16           | 16          | –           | –           |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)             | 6.0<br>9.6 | 6.0<br>9.6 | 6.0<br>9.6 | 6.0<br>12.0 | 6.0<br>12.0 | 6.0<br>12.0  | 6.0<br>12.0 | 6.0<br>12.0 | 6.0<br>12.0 |
| <b>f<sub>awf</sub> (MHz)</b>                                | 4.45       | 4.42       | 4.45       | 3.99        | 3.99        | 4.04         | 4.04        | 4.04        | 4.04        |
| <b>APF (%)</b><br><b>(Acoustic Power up Fraction)</b>       | na         | na         | na         | na          | na          | na           | na          | na          | na          |
| <b>AIF (%)</b><br><b>(Acoustic Initialization Fraction)</b> | na         | na         | na         | na          | na          | na           | na          | na          | na          |
| <b>Maximum Power (mW)</b>                                   | 45.09      | 5.37       | 35.42      | 20.54       | 50.60       | 147.33       | 98.96       | 26.77       | 26.77       |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                   | 78.27      | 9.32       | 61.50      | 28.53       | 80.71       | 215.06       | 147.88      | 39.05       | 39.05       |
| <b>Power-up Mode</b>                                        | na         | na         | na         | na          | na          | na           | na          | na          | na          |
| <b>Initialization Mode</b>                                  | na         | na         | na         | na          | na          | na           | na          | na          | na          |
| <b>Acoustic Output Freeze</b>                               | yes        | yes        | yes        | yes         | yes         | yes          | yes         | yes         | yes         |
| <b>I<sub>tt</sub> (mm)</b>                                  | 0.00       | 0.00       | 0.00       | 0.00        | 0.00        | 0.00         | 0.00        | 0.00        | 0.00        |
| <b>I<sub>ts</sub> (mm)</b>                                  | contact    | contact    | contact    | contact     | contact     | con-<br>tact | contact     | contact     | contact     |

Table 17–131. 546L IEC Acoustic Output Information



## IEC Acoustic Output Tables

### ERB7L Probe

| Parameter                                           | B          | M           | B+M        | D          | B+D        |
|-----------------------------------------------------|------------|-------------|------------|------------|------------|
| P <sub>-</sub> (MPa)                                | 3.04       | 3.04        | 3.04       | 3.00       | 3.04       |
| I <sub>spta</sub> (mW/cm <sup>2</sup> )             | 28.3       | 105.8       | 124.7      | 1064.4     | 1083.3     |
| I <sub>p</sub> (mm)                                 | 20.0       | 27.4        | 20.0       | 20.8       | 20.8       |
| W <sub>pb6</sub><br>(  ) (mm)<br>(⊥)(mm)            | 1.1<br>1.7 | 0.8<br>1.4  | 1.1<br>1.7 | 0.8<br>1.3 | 0.8<br>1.3 |
| pr <sub>r</sub> (kHz)                               | –          | 0.82        | 0.82       | 1.13       | 1.13       |
| s <sub>r</sub> (Hz)                                 | 82         | –           | 82         | –          | 82         |
| Output Beam Dim.<br>(  ) (mm)<br>(⊥)(mm)            | 5.5<br>7.2 | 5.5<br>11.2 | 5.5<br>7.2 | 5.5<br>7.2 | 5.5<br>7.2 |
| f <sub>awf</sub> (MHz)                              | 4.52       | 6.06        | 4.52       | 6.62       | 6.62       |
| APF (%)<br>(Acoustic Power<br>up Fraction)          | na         | na          | na         | na         | na         |
| AIF (%)<br>(Acoustic<br>Initialization<br>Fraction) | na         | na          | na         | na         | na         |
| Maximum Power (mW)                                  | 6.90       | 3.59        | 8.19       | 9.27       | 13.88      |
| I <sub>ob</sub> (mW/cm <sup>2</sup> )               | 17.43      | 5.82        | 17.45      | 23.42      | 35.04      |
| Power-up Mode                                       | na         | na          | na         | na         | na         |
| Initialization Mode                                 | na         | na          | na         | na         | na         |
| Acoustic Output<br>Freeze                           | yes        | yes         | yes        | yes        | yes        |
| I <sub>tt</sub> (mm)                                | 0.00       | 0.00        | 0.00       | 0.00       | 0.00       |
| I <sub>ts</sub> (mm)                                | contact    | contact     | contact    | contact    | contact    |

Table 17–132. ERB7L IEC Acoustic Output Information



ERB7L Probe (cont'd)

| Parameter                                                  | B+rD       | B+rD+D     | M+cM       | cM+D       |
|------------------------------------------------------------|------------|------------|------------|------------|
| <b>P<sub>-</sub> (MPa)</b>                                 | 3.04       | 3.04       | 3.04       | 3.04       |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                | 178.4      | 630.8      | 878.8      | 1064.4     |
| <b>I<sub>p</sub> (mm)</b>                                  | 19.6       | 19.6       | 19.6       | 19.6       |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)             | 0.9<br>1.3 | 0.9<br>1.3 | 0.9<br>1.3 | 0.9<br>1.3 |
| <b>prf (kHz)</b>                                           | 1.72       | 1.72       | 1.70       | 1.70       |
| <b>srr (Hz)</b>                                            | 12         | 12         | –          | –          |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)            | 5.5<br>7.2 | 5.5<br>7.2 | 5.5<br>7.2 | 5.5<br>7.2 |
| <b>f<sub>awf</sub> (MHz)</b>                               | 5.07       | 5.07       | 5.07       | 5.07       |
| <b>APF (%)</b><br>(Acoustic Power<br>up Fraction)          | na         | na         | na         | na         |
| <b>AIF (%)</b><br>(Acoustic<br>Initialization<br>Fraction) | na         | na         | na         | na         |
| <b>Maximum Power (mW)</b>                                  | 39.69      | 26.78      | 24.01      | 24.01      |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                  | 100.23     | 67.64      | 57.38      | 57.38      |
| <b>Power-up Mode</b>                                       | na         | na         | na         | na         |
| <b>Initialization Mode</b>                                 | na         | na         | na         | na         |
| <b>Acoustic Output<br/>Freeze</b>                          | yes        | yes        | yes        | yes        |
| <b>I<sub>tt</sub> (mm)</b>                                 | 0.00       | 0.00       | 0.00       | 0.00       |
| <b>I<sub>ts</sub> (mm)</b>                                 | contact    | contact    | contact    | contact    |

Table 17–133. ERB7L IEC Acoustic Output Information (cont'd)



## IEC Acoustic Output Tables

### ERB7C Probe

| Parameter                                           | B          | M          | B+M        | D          | B+D        |
|-----------------------------------------------------|------------|------------|------------|------------|------------|
| P <sub>-</sub> (MPa)                                | 2.98       | 2.98       | 2.98       | 3.06       | 3.06       |
| I <sub>spta</sub> (mW/cm <sup>2</sup> )             | 61.7       | 114.1      | 155.2      | 581.9      | 623.0      |
| I <sub>p</sub> (mm)                                 | 19.5       | 18.7       | 19.5       | 22.8       | 22.8       |
| W <sub>pb6</sub><br>(  ) (mm)<br>(⊥)(mm)            | 1.0<br>1.7 | 0.8<br>2.2 | 1.0<br>1.7 | 1.0<br>1.5 | 1.0<br>1.5 |
| prf (kHz)                                           | –          | 0.82       | 0.82       | 2.60       | 2.60       |
| srr (Hz)                                            | 180        | –          | 180        | –          | 180        |
| Output Beam Dim.<br>(  ) (mm)<br>(⊥)(mm)            | 5.5<br>8.4 | 5.5<br>8.4 | 5.5<br>8.4 | 5.5<br>8.4 | 5.5<br>8.4 |
| f <sub>awf</sub> (MHz)                              | 4.49       | 5.78       | 4.49       | 4.98       | 4.98       |
| APF (%)<br>(Acoustic Power<br>up Fraction)          | na         | na         | na         | na         | na         |
| AIF (%)<br>(Acoustic<br>Initialization<br>Fraction) | na         | na         | na         | na         | na         |
| Maximum Power (mW)                                  | 10.16      | 4.76       | 11.53      | 7.35       | 14.13      |
| I <sub>ob</sub> (mW/cm <sup>2</sup> )               | 22.06      | 10.34      | 25.05      | 15.97      | 30.68      |
| Power-up Mode                                       | na         | na         | na         | na         | na         |
| Initialization Mode                                 | na         | na         | na         | na         | na         |
| Acoustic Output<br>Freeze                           | yes        | yes        | yes        | yes        | yes        |
| I <sub>tt</sub> (mm)                                | 0.00       | 0.00       | 0.00       | 0.00       | 0.00       |
| I <sub>ts</sub> (mm)                                | contact    | contact    | contact    | contact    | contact    |

Table 17–134. ERB7C IEC Acoustic Output Information

ERB7C Probe (cont'd)

| Parameter                                                  | B+rD       | B+rD+D     | M+cM       | cM+D       |
|------------------------------------------------------------|------------|------------|------------|------------|
| <b>P<sub>-</sub> (MPa)</b>                                 | 2.98       | 3.06       | 2.98       | 2.98       |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                | 51.8       | 337.4      | 368.5      | 581.9      |
| <b>I<sub>p</sub> (mm)</b>                                  | 22.0       | 22.0       | 22.0       | 22.0       |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)             | 1.0<br>1.5 | 1.0<br>1.5 | 1.0<br>1.5 | 1.0<br>1.5 |
| <b>prr (kHz)</b>                                           | 5.20       | 5.20       | 5.20       | 5.20       |
| <b>srr (Hz)</b>                                            | 12         | 12         | –          | –          |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)            | 5.5<br>8.4 | 5.5<br>8.4 | 5.5<br>8.4 | 5.5<br>8.4 |
| <b>f<sub>awf</sub> (MHz)</b>                               | 5.02       | 5.02       | 5.02       | 5.02       |
| <b>APF (%)</b><br>(Acoustic Power<br>up Fraction)          | na         | na         | na         | na         |
| <b>AIF (%)</b><br>(Acoustic<br>Initialization<br>Fraction) | na         | na         | na         | na         |
| <b>Maximum Power (mW)</b>                                  | 9.73       | 11.93      | 7.97       | 7.97       |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                  | 21.13      | 25.91      | 17.31      | 17.31      |
| <b>Power-up Mode</b>                                       | na         | na         | na         | na         |
| <b>Initialization Mode</b>                                 | na         | na         | na         | na         |
| <b>Acoustic Output<br/>Freeze</b>                          | yes        | yes        | yes        | yes        |
| <b>I<sub>tt</sub> (mm)</b>                                 | 0.00       | 0.00       | 0.00       | 0.00       |
| <b>I<sub>ts</sub> (mm)</b>                                 | contact    | contact    | contact    | contact    |

Table 17–135. ERB7C IEC Acoustic Output Information (cont'd)



## IEC Acoustic Output Tables

### L764 (LH) Probe

| Parameter                                                                       | B       | M       | B+M     | D       | B+D     | B+rD    | B+rD<br>+D | M+cM    | cM+D    |
|---------------------------------------------------------------------------------|---------|---------|---------|---------|---------|---------|------------|---------|---------|
| <b>P<sub>-</sub> (MPa)</b>                                                      | 4.44    | 4.44    | 4.44    | 3.64    | 4.44    | 4.44    | 4.44       | 3.02    | 3.02    |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                                     | 40.96   | 94.91   | 122.22  | 749.20  | 776.51  | 185.91  | 481.21     | 1502.8  | 1502.8  |
| <b>I<sub>p</sub> (mm)</b>                                                       | 27.46   | 27.46   | 27.46   | 24.70   | 25.00   | 25.00   | 25.00      | 25.00   | 25.00   |
| <b>W<sub>pb6</sub></b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>                    | 1.36    | 1.36    | 1.36    | 1.17    | 1.17    | 0.94    | 0.94       | 0.94    | 0.94    |
|                                                                                 | 1.39    | 1.39    | 1.39    | 1.43    | 1.43    | 1.44    | 1.44       | 1.44    | 1.44    |
| <b>pr<sub>r</sub> (kHz)</b>                                                     | –       | 0.87    | 0.87    | 2.00    | 2.00    | 10.40   | 10.40      | 10.40   | 10.40   |
| <b>s<sub>r</sub> (Hz)</b>                                                       | 33.00   | –       | 33.00   | –       | 33.00   | 35.00   | 35.00      | –       | –       |
| <b>Output Beam</b><br><b>Dim. (  ) (mm)</b><br><b>(⊥)(mm)</b>                   | 1.36    | –       | 1.36    | –       | 1.36    | 0.94    | 0.94       | –       | –       |
|                                                                                 | 1.39    | –       | 1.39    | –       | 1.39    | 1.44    | 1.44       | –       | –       |
| <b>f<sub>awf</sub> (MHz)</b>                                                    | 6.36    | 6.36    | 6.36    | 5.06    | 5.06    | 5.16    | 5.16       | 5.16    | 5.16    |
| <b>APF (%)</b><br><b>(Acoustic Power</b><br><b>up Fraction)</b>                 | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>AIF (%)</b><br><b>(Acoustic</b><br><b>Initialization</b><br><b>Fraction)</b> | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Maximum Power</b><br><b>(mW)</b>                                             | 7.16    | 7.16    | 11.93   | 9.16    | 13.93   | 18.40   | 16.16      | 13.63   | 13.63   |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                                       | 11.36   | 11.36   | 18.93   | 14.53   | 22.10   | 29.21   | 25.66      | 21.64   | 21.64   |
| <b>Power-up Mode</b>                                                            | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Initialization</b><br><b>Mode</b>                                            | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Acoustic Output</b><br><b>Freeze</b>                                         | yes     | yes     | yes     | yes     | yes     | yes     | yes        | yes     | yes     |
| <b>I<sub>tt</sub> (mm)</b>                                                      | 0.00    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00       | 0.00    | 0.00    |
| <b>I<sub>ts</sub> (mm)</b>                                                      | contact | contact | contact | contact | contact | contact | contact    | contact | contact |

Table 17–136. L764 (LH) IEC Acoustic Output Information

## IEC Acoustic Output Tables

### 739L Probe

| Parameter                                            | B          | M          | B+M        | D          | B+D        | B+rD       | B+rD<br>+D | M+cM       | cM+D       |
|------------------------------------------------------|------------|------------|------------|------------|------------|------------|------------|------------|------------|
| <b>P<sub>-</sub> (MPa)</b>                           | 4.85       | 4.85       | 4.85       | 4.06       | 4.85       | 4.85       | 4.85       | 4.85       | 4.85       |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>          | 194.7      | 289.1      | 418.9      | 822.2      | 952.0      | 925.3      | 938.7      | 1147.8     | 1147.8     |
| <b>I<sub>p</sub> (mm)</b>                            | 19.1       | 19.1       | 19.1       | 11.0       | 11.0       | 21.7       | 21.7       | 21.7       | 21.7       |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)       | 0.8<br>1.5 | 0.8<br>1.5 | 0.8<br>1.5 | 0.8<br>3.2 | 0.8<br>3.2 | 0.8<br>1.5 | 0.8<br>1.5 | 0.8<br>1.5 | 0.8<br>1.5 |
| <b>prf (kHz)</b>                                     | –          | 0.87       | 0.87       | 0.65       | 0.65       | 1.72       | 1.72       | 1.70       | 1.70       |
| <b>srr (Hz)</b>                                      | 56         | –          | 56         | –          | 56         | 19         | 19         | –          | –          |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)      | 6.0<br>7.3 | 6.0<br>7.3 | 6.0<br>7.3 | 6.0<br>4.5 | 6.0<br>4.5 | 6.0<br>9.3 | 6.0<br>9.3 | 6.0<br>9.3 | 6.0<br>9.3 |
| <b>f<sub>awf</sub> (MHz)</b>                         | 6.18       | 6.18       | 6.18       | 5.01       | 5.01       | 5.02       | 5.02       | 5.02       | 5.02       |
| <b>APF (%)</b><br>(Acoustic Power up Fraction)       | na         | na         | na         | na         | na         | na         | na         | na         | na         |
| <b>AIF (%)</b><br>(Acoustic Initialization Fraction) | na         | na         | na         | na         | na         | na         | na         | na         | na         |
| <b>Maximum Power (mW)</b>                            | 31.07      | 3.11       | 23.82      | 14.96      | 35.67      | 127.39     | 81.53      | 42.78      | 42.78      |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>            | 70.79      | 7.08       | 54.28      | 55.76      | 102.95     | 237.41     | 170.18     | 77.82      | 7.82       |
| <b>Power-up Mode</b>                                 | na         | na         | na         | na         | na         | na         | na         | na         | na         |
| <b>Initialization Mode</b>                           | na         | na         | na         | na         | na         | na         | na         | na         | na         |
| <b>Acoustic Output Freeze</b>                        | yes        | yes        | yes        | yes        | yes        | yes        | yes        | yes        | yes        |
| <b>I<sub>tt</sub> (mm)</b>                           | 0.00       | 0.00       | 0.00       | 0.00       | 0.00       | 0.00       | 0.00       | 0.00       | 0.00       |
| <b>I<sub>ts</sub> (mm)</b>                           | contact    | contact    | contact    | contact    | contact    | contact    | contact    | contact    | contact    |

Table 17–137. 739L IEC Acoustic Output Information



## IEC Acoustic Output Tables

### LA39 Probe

| Parameter                                                     | B       | M       | B+M     | D       | B+D     | B+rD    | B+rD<br>+D | M+cM    | cM+D    |
|---------------------------------------------------------------|---------|---------|---------|---------|---------|---------|------------|---------|---------|
| <b>P<sub>-</sub> (MPa)</b>                                    | 5.14    | 5.14    | 5.14    | 4.52    | 5.14    | 5.14    | 5.14       | 5.14    | 5.14    |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                   | 190.1   | 338.4   | 465.2   | 981.5   | 1108.3  | 606.7   | 857.5      | 918.1   | 981.5   |
| <b>I<sub>p</sub> (mm)</b>                                     | 10.0    | 10.0    | 10.0    | 9.9     | 9.9     | 9.0     | 9.0        | 9.0     | 9.0     |
| <b>W<sub>pb6</sub></b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>  | 0.6     | 0.6     | 0.6     | 0.7     | 0.7     | 0.6     | 0.6        | 0.6     | 0.6     |
|                                                               | 0.8     | 0.8     | 0.8     | 0.9     | 0.9     | 0.9     | 0.9        | 0.9     | 0.9     |
| <b>pr<sub>r</sub> (kHz)</b>                                   | –       | 0.87    | 0.87    | 2.60    | 2.60    | 0.66    | 0.66       | 1.48    | 1.48    |
| <b>s<sub>r</sub> (Hz)</b>                                     | 57      | –       | 57      | –       | 57      | 19      | 19         | –       | –       |
| <b>Output Beam Dim.</b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b> | 3.3     | 3.3     | 3.3     | 3.3     | 3.3     | 3.3     | 3.3        | 3.3     | 3.3     |
|                                                               | 4.1     | 4.1     | 4.1     | 3.7     | 3.7     | 2.4     | 2.4        | 2.4     | 2.4     |
| <b>f<sub>awf</sub> (MHz)</b>                                  | 7.46    | 7.46    | 7.46    | 6.74    | 6.74    | 6.72    | 6.72       | 6.72    | 6.72    |
| <b>APF (%)</b><br><b>(Acoustic Power up Fraction)</b>         | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>AIF (%)</b><br><b>(Acoustic Initialization Fraction)</b>   | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Maximum Power (mW)</b>                                     | 10.40   | 1.04    | 7.98    | 7.47    | 14.40   | 28.31   | 21.36      | 18.85   | 18.85   |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                     | 78.76   | 7.89    | 60.39   | 62.81   | 115.31  | 322.25  | 218.78     | 232.57  | 232.57  |
| <b>Power-up Mode</b>                                          | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Initialization Mode</b>                                    | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Acoustic Output Freeze</b>                                 | yes     | yes     | yes     | yes     | yes     | yes     | yes        | yes     | yes     |
| <b>I<sub>tt</sub> (mm)</b>                                    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00       | 0.00    | 0.00    |
| <b>I<sub>ts</sub> (mm)</b>                                    | contact | contact | contact | contact | contact | contact | contact    | contact | contact |

Table 17–138. LA39 IEC Acoustic Output Information

LD Probe

| Parameter                                                  | B            | M            | B+M          | D            | B+D          |
|------------------------------------------------------------|--------------|--------------|--------------|--------------|--------------|
| <b>P<sub>-</sub> (MPa)</b>                                 | 2.23         | 2.23         | 2.23         | 2.20         | 2.23         |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                | 20.0         | 122.6        | 136.0        | 4626.0       | 4639.3       |
| <b>I<sub>p</sub> (mm)</b>                                  | 66.4         | 70.7         | 66.4         | 80.8         | 80.8         |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)             | 3.5<br>2.5   | 2.2<br>2.3   | 3.5<br>2.5   | 2.1<br>2.6   | 2.1<br>2.6   |
| <b>prf (kHz)</b>                                           | –            | 0.89         | 0.89         | 4.53         | 4.53         |
| <b>srr (Hz)</b>                                            | 26           | –            | 26           | –            | 26           |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)            | 13.8<br>22.5 | 13.8<br>22.5 | 13.8<br>22.5 | 13.8<br>19.5 | 13.8<br>19.5 |
| <b>f<sub>awf</sub> (MHz)</b>                               | 3.50         | 3.28         | 3.50         | 3.32         | 3.32         |
| <b>APF (%)</b><br>(Acoustic Power<br>up Fraction)          | na           | na           | na           | na           | na           |
| <b>AIF (%)</b><br>(Acoustic<br>Initialization<br>Fraction) | na           | na           | na           | na           | na           |
| <b>Maximum Power (mW)</b>                                  | 38.80        | 4.33         | 30.20        | 125.61       | 151.47       |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                  | 12.50        | 1.39         | 9.73         | 46.68        | 55.01        |
| <b>Power-up Mode</b>                                       | na           | na           | na           | na           | na           |
| <b>Initialization Mode</b>                                 | na           | na           | na           | na           | na           |
| <b>Acoustic Output<br/>Freeze</b>                          | yes          | yes          | yes          | yes          | yes          |
| <b>I<sub>tt</sub> (mm)</b>                                 | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         |
| <b>I<sub>ts</sub> (mm)</b>                                 | contact      | contact      | contact      | contact      | contact      |

Table 17–139. LD IEC Acoustic Output Information



## IEC Acoustic Output Tables

### LD Probe (cont'd)

| Parameter                                                                       | B+rD         | B+rD+D       | M+cM         | cM+D         |
|---------------------------------------------------------------------------------|--------------|--------------|--------------|--------------|
| <b>P<sub>-</sub> (MPa)</b>                                                      | 2.97         | 2.97         | 2.60         | 2.60         |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                                     | 361.7        | 2500.5       | 3057.6       | 3057.6       |
| <b>I<sub>p</sub> (mm)</b>                                                       | 73.1         | 73.1         | 73.1         | 73.1         |
| <b>w<sub>pb6</sub></b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>                    | 2.1<br>2.4   | 2.1<br>2.4   | 2.1<br>2.4   | 2.1<br>2.4   |
| <b>prf (kHz)</b>                                                                | 5.96         | 5.96         | 5.96         | 5.96         |
| <b>srr (Hz)</b>                                                                 | 16           | 16           | –            | –            |
| <b>Output Beam Dim.</b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>                   | 13.8<br>19.5 | 13.8<br>19.5 | 13.8<br>19.5 | 13.8<br>19.5 |
| <b>f<sub>awf</sub> (MHz)</b>                                                    | 3.30         | 3.30         | 3.30         | 3.30         |
| <b>APF (%)</b><br><b>(Acoustic Power</b><br><b>up Fraction)</b>                 | na           | na           | na           | na           |
| <b>AIF (%)</b><br><b>(Acoustic</b><br><b>Initialization</b><br><b>Fraction)</b> | na           | na           | na           | na           |
| <b>Maximum Power (mW)</b>                                                       | 103.01       | 127.24       | 91.59        | 91.59        |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                                       | 37.00        | 46.00        | 33.82        | 38.82        |
| <b>Power-up Mode</b>                                                            | na           | na           | na           | na           |
| <b>Initialization Mode</b>                                                      | na           | na           | na           | na           |
| <b>Acoustic Output</b><br><b>Freeze</b>                                         | yes          | yes          | yes          | yes          |
| <b>I<sub>tt</sub> (mm)</b>                                                      | 0.00         | 0.00         | 0.00         | 0.00         |
| <b>I<sub>ts</sub> (mm)</b>                                                      | contact      | contact      | contact      | contact      |

Table 17–140. LD IEC Acoustic Output Information (cont'd)



## IEC Acoustic Output Tables

### I739 Probe

| Parameter                                            | B            | M            | B+M          | D            | B+D          | B+rD         | B+rD<br>+D   | M+cM         | cM+D         |
|------------------------------------------------------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|
| <b>P<sub>-</sub> (MPa)</b>                           | 3.02         | 3.02         | 3.02         | 2.64         | 3.02         | 3.02         | 3.02         | 2.16         | 2.16         |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>          | 58.28        | 83.35        | 122.21       | 831.40       | 870.26       | 474.36       | 672.31       | 601.95       | 601.95       |
| <b>I<sub>p</sub> (mm)</b>                            | 20.90        | 20.90        | 20.90        | 16.60        | 16.60        | 17.60        | 17.60        | 17.60        | 17.60        |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)       | 1.11<br>0.84 | 1.11<br>0.84 | 1.11<br>0.84 | 1.35<br>1.50 | 1.35<br>1.50 | 1.35<br>1.50 | 1.35<br>1.50 | 1.35<br>1.50 | 1.35<br>1.50 |
| <b>pr<sub>r</sub> (kHz)</b>                          | –            | 0.87         | 0.87         | 18.12        | 18.12        | 6.84         | 6.84         | 6.84         | 6.84         |
| <b>s<sub>r</sub> (Hz)</b>                            | 57.00        | –            | 57.00        | –            | 57.00        | 10.00        | 10.00        | –            | –            |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)      | 7.30<br>6.00 | 7.30<br>6.00 | 7.30<br>6.00 | 3.70<br>6.00 | 7.30<br>6.00 | 7.30<br>6.00 | 7.30<br>6.00 | 3.70<br>6.00 | 3.70<br>6.00 |
| <b>f<sub>awf</sub> (MHz)</b>                         | 6.99         | 6.99         | 6.99         | 5.07         | 5.07         | 5.11         | 5.11         | 5.11         | 5.11         |
| <b>APF (%)</b><br>(Acoustic Power up Fraction)       | na           | na           | na           | na           | na           | na           | na           | na           | na           |
| <b>AIF (%)</b><br>(Acoustic Initialization Fraction) | na           | na           | na           | na           | na           | na           | na           | na           | na           |
| <b>Maximum Power (mW)</b>                            | 4.87         | 0.46         | 3.71         | 9.61         | 12.86        | 24.62        | 18.74        | 21.37        | 21.37        |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>            | 11.09        | 1.05         | 8.44         | 43.80        | 51.19        | 104.78       | 77.98        | 98.43        | 98.43        |
| <b>Power-up Mode</b>                                 | na           | na           | na           | na           | na           | na           | na           | na           | na           |
| <b>Initialization Mode</b>                           | na           | na           | na           | na           | na           | na           | na           | na           | na           |
| <b>Acoustic Output Freeze</b>                        | yes          | yes          | yes          | yes          | yes          | yes          | yes          | yes          | yes          |
| <b>I<sub>tt</sub> (mm)</b>                           | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         |
| <b>I<sub>ts</sub> (mm)</b>                           | contact      | contact      | contact      | contact      | contact      | contact      | contact      | contact      | contact      |

Table 17–141. I739 IEC Acoustic Output Information



## IEC Acoustic Output Tables

### T739 Probe

| Parameter                                                     | B            | M            | B+M          | D            | B+D          | B+rD         | B+rD<br>+D   | M+cM         | cM+D         |
|---------------------------------------------------------------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|
| <b>P<sub>-</sub> (MPa)</b>                                    | 3.45         | 3.45         | 3.45         | 2.82         | 3.45         | 3.45         | 3.45         | 2.65         | 2.65         |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                   | 54.48        | 65.96        | 102.28       | 784.40       | 820.72       | 527.22       | 673.97       | 767.06       | 767.06       |
| <b>I<sub>p</sub> (mm)</b>                                     | 21.40        | 21.40        | 21.40        | 17.30        | 17.30        | 18.40        | 18.40        | 18.40        | 18.40        |
| <b>W<sub>pb6</sub></b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>  | 1.29<br>0.85 | 1.29<br>0.85 | 1.29<br>0.85 | 1.57<br>1.55 | 1.57<br>1.55 | 1.57<br>1.55 | 1.57<br>1.55 | 1.57<br>1.55 | 1.57<br>1.55 |
| <b>pr<sub>r</sub> (kHz)</b>                                   | –            | 0.87         | 0.87         | 18.10        | 18.10        | 6.84         | 6.84         | 6.84         | 6.84         |
| <b>s<sub>r</sub> (Hz)</b>                                     | 57.00        | –            | 57.00        | –            | 57.00        | 10.00        | 10.00        | –            | –            |
| <b>Output Beam<br/>Dim. (  ) (mm)</b><br><b>(⊥)(mm)</b>       | 7.30<br>6.00 | 7.30<br>6.00 | 7.30<br>6.00 | 3.70<br>6.00 | 7.30<br>6.00 | 7.30<br>6.00 | 7.30<br>6.00 | 3.70<br>6.00 | 3.70<br>6.00 |
| <b>f<sub>awf</sub> (MHz)</b>                                  | 6.96         | 6.96         | 6.96         | 5.07         | 5.07         | 5.08         | 5.08         | 5.08         | 5.08         |
| <b>APF (%)<br/>(Acoustic Power<br/>up Fraction)</b>           | na           | na           | na           | na           | na           | na           | na           | na           | na           |
| <b>AIF (%)<br/>(Acoustic<br/>Initialization<br/>Fraction)</b> | na           | na           | na           | na           | na           | na           | na           | na           | na           |
| <b>Maximum Power<br/>(mW)</b>                                 | 3.91         | 0.54         | 3.15         | 11.78        | 14.39        | 18.30        | 16.34        | 15.69        | 15.69        |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                     | 8.92         | 1.22         | 7.16         | 53.68        | 59.62        | 77.44        | 68.53        | 72.72        | 72.72        |
| <b>Power-up Mode</b>                                          | na           | na           | na           | na           | na           | na           | na           | na           | na           |
| <b>Initialization<br/>Mode</b>                                | na           | na           | na           | na           | na           | na           | na           | na           | na           |
| <b>Acoustic Output<br/>Freeze</b>                             | yes          | yes          | yes          | yes          | yes          | yes          | yes          | yes          | yes          |
| <b>I<sub>tt</sub> (mm)</b>                                    | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         |
| <b>I<sub>ts</sub> (mm)</b>                                    | contact      | contact      | contact      | contact      | contact      | contact      | contact      | contact      | contact      |

Table 17–142. T739 IEC Acoustic Output Information

**S316 Probe**

| Parameter                                                                       | B       | M       | B+M     | D       | B+D     | B+rD    | B+rD<br>+D | M+cM    | cM+D    |
|---------------------------------------------------------------------------------|---------|---------|---------|---------|---------|---------|------------|---------|---------|
| <b>P<sub>-</sub> (MPa)</b>                                                      | 2.15    | 2.15    | 2.15    | 2.05    | 2.15    | 2.15    | 2.15       | 2.15    | 2.15    |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                                     | 92.8    | 148.7   | 210.6   | 851.8   | 913.7   | 257.2   | 585.5      | 436.3   | 436.3   |
| <b>I<sub>p</sub> (mm)</b>                                                       | 58.9    | 55.7    | 58.9    | 62.3    | 62.3    | 57.5    | 57.5       | 57.5    | 57.5    |
| <b>W<sub>pb6</sub></b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>                    | 2.0     | 1.7     | 2.0     | 2.3     | 2.3     | 1.9     | 1.9        | 1.9     | 1.9     |
|                                                                                 | 3.7     | 4.5     | 3.7     | 3.1     | 3.1     | 3.6     | 3.6        | 3.6     | 3.6     |
| <b>pr<sub>r</sub> (kHz)</b>                                                     | –       | 0.89    | 0.89    | 4.53    | 4.53    | 5.20    | 5.20       | 5.20    | 5.20    |
| <b>s<sub>r</sub> (Hz)</b>                                                       | 112     | –       | 112     | –       | 112     | 39      | 39         | –       | –       |
| <b>Output Beam</b><br><b>Dim. (  ) (mm)</b><br><b>(⊥)(mm)</b>                   | 13.0    | 13.0    | 13.0    | 13.0    | 13.0    | 13.0    | 13.0       | 13.0    | 13.0    |
|                                                                                 | 16.0    | 16.0    | 16.0    | 14.5    | 14.5    | 14.5    | 14.5       | 14.5    | 14.5    |
| <b>f<sub>awf</sub> (MHz)</b>                                                    | 3.42    | 3.36    | 3.42    | 3.00    | 3.00    | 2.98    | 2.98       | 2.98    | 2.98    |
| <b>APF (%)</b><br><b>(Acoustic Power</b><br><b>up Fraction)</b>                 | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>AIF (%)</b><br><b>(Acoustic</b><br><b>Initialization</b><br><b>Fraction)</b> | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Maximum Power</b><br><b>(mW)</b>                                             | 16.83   | 5.00    | 16.22   | 27.98   | 39.20   | 91.88   | 65.54      | 36.42   | 36.42   |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                                       | 8.09    | 2.40    | 7.80    | 14.84   | 20.24   | 48.18   | 34.21      | 19.07   | 19.07   |
| <b>Power-up Mode</b>                                                            | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Initialization</b><br><b>Mode</b>                                            | na      | na      | na      | na      | na      | na      | na         | na      | na      |
| <b>Acoustic Output</b><br><b>Freeze</b>                                         | yes     | yes     | yes     | yes     | yes     | yes     | yes        | yes     | yes     |
| <b>I<sub>tt</sub> (mm)</b>                                                      | 0.00    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00    | 0.00       | 0.00    | 0.00    |
| <b>I<sub>ts</sub> (mm)</b>                                                      | contact | contact | contact | contact | contact | contact | contact    | contact | contact |

Table 17–143. S316 IEC Acoustic Output Information



## IEC Acoustic Output Tables

### S317 Probe

| Parameter                                                  | B            | M            | B+M          | D            | B+D          |
|------------------------------------------------------------|--------------|--------------|--------------|--------------|--------------|
| <b>P<sub>-</sub> (MPa)</b>                                 | 2.69         | 2.69         | 2.69         | 1.96         | 2.69         |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                | 817.8        | 555.2        | 1100.4       | 1495.8       | 2041.0       |
| <b>I<sub>p</sub> (mm)</b>                                  | 45.4         | 45.4         | 45.4         | 66.0         | 66.0         |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)             | 2.0<br>3.6   | 2.0<br>3.6   | 2.0<br>3.6   | 2.3<br>3.9   | 2.3<br>3.9   |
| <b>pr<sub>r</sub> (kHz)</b>                                | –            | 0.86         | 0.86         | 1.71         | 1.71         |
| <b>s<sub>r</sub> (Hz)</b>                                  | 82           | –            | 82           | –            | 82           |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)            | 12.0<br>12.9 | 12.0<br>12.9 | 12.0<br>12.9 | 12.0<br>17.9 | 12.0<br>17.9 |
| <b>f<sub>awf</sub> (MHz)</b>                               | 2.67         | 2.67         | 2.67         | 2.51         | 2.51         |
| <b>APF (%)</b><br>(Acoustic Power<br>up Fraction)          | na           | na           | na           | na           | na           |
| <b>AIF (%)</b><br>(Acoustic<br>Initialization<br>Fraction) | na           | na           | na           | na           | na           |
| <b>Maximum Power (mW)</b>                                  | 122.81       | 20.36        | 102.23       | 70.54        | 152.42       |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                  | 79.46        | 13.17        | 66.14        | 32.81        | 85.78        |
| <b>Power-up Mode</b>                                       | na           | na           | na           | na           | na           |
| <b>Initialization Mode</b>                                 | na           | na           | na           | na           | na           |
| <b>Acoustic Output<br/>Freeze</b>                          | yes          | yes          | yes          | yes          | yes          |
| <b>I<sub>tt</sub> (mm)</b>                                 | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         |
| <b>I<sub>ts</sub> (mm)</b>                                 | contact      | contact      | contact      | contact      | contact      |

Table 17–144. S317 IEC Acoustic Output Information

**S317 Probe (cont'd)**

| Parameter                                                                       | B+rD         | B+rD+D       | M+cM         | cM+D         | cwD         |
|---------------------------------------------------------------------------------|--------------|--------------|--------------|--------------|-------------|
| <b>P<sub>-</sub> (MPa)</b>                                                      | 2.69         | 2.69         | 2.69         | 2.69         | 0.17        |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                                     | 787.4        | 1414.2       | 1941.9       | 1941.9       | 959.0       |
| <b>I<sub>p</sub> (mm)</b>                                                       | 63.2         | 63.2         | 63.2         | 63.2         | 39.9        |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)                                  | 2.2<br>3.7   | 2.2<br>3.7   | 2.2<br>3.7   | 2.2<br>3.7   | 4.6<br>3.6  |
| <b>pr<sub>r</sub> (kHz)</b>                                                     | 5.96         | 5.96         | 5.96         | 5.96         | —           |
| <b>s<sub>r</sub> (Hz)</b>                                                       | 43           | 43           | —            | —            | —           |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)                                 | 12.0<br>17.9 | 12.0<br>17.9 | 12.0<br>17.9 | 12.0<br>17.9 | 12.0<br>6.7 |
| <b>f<sub>awf</sub> (MHz)</b>                                                    | 2.56         | 2.56         | 2.56         | 2.56         | 2.49        |
| <b>APF (%)</b><br><b>(Acoustic Power</b><br><b>up Fraction)</b>                 | na           | na           | na           | na           | na          |
| <b>AIF (%)</b><br><b>(Acoustic</b><br><b>Initialization</b><br><b>Fraction)</b> | na           | na           | na           | na           | na          |
| <b>Maximum Power (mW)</b>                                                       | 227.34       | 189.88       | 77.86        | 77.86        | 96.14       |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                                       | 120.62       | 103.20       | 39.91        | 39.91        | 119.22      |
| <b>Power-up Mode</b>                                                            | na           | na           | na           | na           | na          |
| <b>Initialization Mode</b>                                                      | na           | na           | na           | na           | na          |
| <b>Acoustic Output</b><br><b>Freeze</b>                                         | yes          | yes          | yes          | yes          | yes         |
| <b>I<sub>tt</sub> (mm)</b>                                                      | 0.00         | 0.00         | 0.00         | 0.00         | 0.00        |
| <b>I<sub>ts</sub> (mm)</b>                                                      | contact      | contact      | contact      | contact      | contact     |

Table 17–144. S317 IEC Acoustic Output Information (cont'd)



## IEC Acoustic Output Tables

### S220 (W) Probe

| Parameter                                           | B            | M            | B+M          | D            | B+D          |
|-----------------------------------------------------|--------------|--------------|--------------|--------------|--------------|
| P <sub>-</sub> (MPa)                                | 2.74         | 2.74         | 2.74         | 2.54         | 2.74         |
| I <sub>spta</sub> (mW/cm <sup>2</sup> )             | 682.46       | 182.47       | 637.44       | 1518.40      | 1973.38      |
| I <sub>p</sub> (mm)                                 | 73.00        | 73.00        | 73.00        | 63.00        | 63.00        |
| W <sub>pb6</sub><br>(  ) (mm)<br>(⊥)(mm)            | 3.01<br>4.42 | 3.01<br>4.42 | 3.01<br>4.42 | 2.70<br>6.25 | 2.70<br>6.25 |
| prf (kHz)                                           | –            | 0.87         | 0.87         | 3.40         | 3.40         |
| srr (Hz)                                            | 85.00        | –            | 85.00        | –            | 85.00        |
| Output Beam Dim.<br>(  ) (mm)<br>(⊥)(mm)            | 3.01<br>4.42 | –<br>–       | 3.01<br>4.42 | –<br>–       | 3.01<br>4.42 |
| f <sub>awf</sub> (MHz)                              | 2.44         | 2.44         | 2.44         | 2.49         | 2.49         |
| APF (%)<br>(Acoustic Power<br>up Fraction)          | na           | na           | na           | na           | na           |
| AIF (%)<br>(Acoustic<br>Initialization<br>Fraction) | na           | na           | na           | na           | na           |
| Maximum Power (mW)                                  | 52.37        | 7.81         | 42.72        | 85.73        | 120.64       |
| I <sub>ob</sub> (mW/cm <sup>2</sup> )               | 15.98        | 2.38         | 13.04        | 2.90         | 13.56        |
| Power-up Mode                                       | na           | na           | na           | na           | na           |
| Initialization Mode                                 | na           | na           | na           | na           | na           |
| Acoustic Output<br>Freeze                           | yes          | yes          | yes          | yes          | yes          |
| I <sub>tt</sub> (mm)                                | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         |
| I <sub>ts</sub> (mm)                                | contact      | contact      | contact      | contact      | contact      |

Table 17–145. S220 (W) IEC Acoustic Output Information

## IEC Acoustic Output Tables

### S220 (W) Probe (cont'd)

| Parameter                                                  | B+rD         | B+rD+D       | M+cM         | cM+D         | cwD          |
|------------------------------------------------------------|--------------|--------------|--------------|--------------|--------------|
| <b>P<sub>-</sub> (MPa)</b>                                 | 3.56         | 3.56         | 3.56         | 3.56         | 0.21         |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                | 1410.08      | 1691.73      | 2157.20      | 2157.20      | 1340.20      |
| <b>I<sub>p</sub> (mm)</b>                                  | 67.70        | 67.70        | 67.75        | 67.75        | 68.50        |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)             | 3.12<br>4.15 | 3.12<br>4.15 | 3.12<br>4.15 | 3.12<br>4.15 | 0.06<br>0.06 |
| <b>pr<sub>r</sub> (kHz)</b>                                | 6.80         | 6.80         | 6.80         | 6.80         | —            |
| <b>s<sub>r</sub> (Hz)</b>                                  | 45.00        | 45.00        | —            | —            | —            |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)            | 3.12<br>4.15 | 3.12<br>4.15 | —<br>—       | —<br>—       | —<br>—       |
| <b>f<sub>awf</sub> (MHz)</b>                               | 2.27         | 2.27         | 2.27         | 2.27         | 2.48         |
| <b>APF (%)</b><br>(Acoustic Power<br>up Fraction)          | na           | na           | na           | na           | na           |
| <b>AIF (%)</b><br>(Acoustic<br>Initialization<br>Fraction) | na           | na           | na           | na           | na           |
| <b>Maximum Power (mW)</b>                                  | 110.67       | 115.66       | 75.76        | 75.76        | 103.44       |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                  | 33.77        | 23.66        | 23.12        | 23.12        | 84.18        |
| <b>Power-up Mode</b>                                       | na           | na           | na           | na           | na           |
| <b>Initialization Mode</b>                                 | na           | na           | na           | na           | na           |
| <b>Acoustic Output<br/>Freeze</b>                          | yes          | yes          | yes          | yes          | yes          |
| <b>I<sub>tt</sub> (mm)</b>                                 | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         |
| <b>I<sub>ts</sub> (mm)</b>                                 | contact      | contact      | contact      | contact      | contact      |

Table 17–145. S220 (W) IEC Acoustic Output Information (cont'd)



## IEC Acoustic Output Tables

### S222 Probe

| Parameter                                                                       | B            | M            | B+M          | D            | B+D          |
|---------------------------------------------------------------------------------|--------------|--------------|--------------|--------------|--------------|
| <b>P<sub>-</sub> (MPa)</b>                                                      | 2.66         | 2.66         | 2.66         | 1.77         | 2.66         |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                                     | 691.6        | 837.9        | 1298.9       | 1818.9       | 1829.3       |
| <b>I<sub>p</sub> (mm)</b>                                                       | 46.6         | 46.6         | 46.6         | 71.6         | 71.6         |
| <b>W<sub>pb6</sub></b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>                    | 2.5<br>3.4   | 2.5<br>3.4   | 2.5<br>3.4   | 3.4<br>2.2   | 3.4<br>2.2   |
| <b>pr<sub>r</sub> (kHz)</b>                                                     | –            | 0.88         | 0.88         | 2.49         | 2.49         |
| <b>s<sub>r</sub> (Hz)</b>                                                       | 42           | –            | 42           | –            | 42           |
| <b>Output Beam Dim.</b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>                   | 13.0<br>10.2 | 13.0<br>10.2 | 13.0<br>10.2 | 13.0<br>20.5 | 13.0<br>20.5 |
| <b>f<sub>awf</sub> (MHz)</b>                                                    | 1.99         | 1.99         | 1.99         | 2.49         | 2.49         |
| <b>APF (%)</b><br><b>(Acoustic Power</b><br><b>up Fraction)</b>                 | na           | na           | na           | na           | na           |
| <b>AIF (%)</b><br><b>(Acoustic</b><br><b>Initialization</b><br><b>Fraction)</b> | na           | na           | na           | na           | na           |
| <b>Maximum Power (mW)</b>                                                       | 102.61       | 33.86        | 102.27       | 61.79        | 130.19       |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                                       | 77.08        | 25.44        | 76.83        | 23.21        | 74.60        |
| <b>Power-up Mode</b>                                                            | na           | na           | na           | na           | na           |
| <b>Initialization Mode</b>                                                      | na           | na           | na           | na           | na           |
| <b>Acoustic Output</b><br><b>Freeze</b>                                         | yes          | yes          | yes          | yes          | yes          |
| <b>I<sub>tt</sub> (mm)</b>                                                      | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         |
| <b>I<sub>ts</sub> (mm)</b>                                                      | contact      | contact      | contact      | contact      | contact      |

Table 17–146. S222 IEC Acoustic Output Information



## IEC Acoustic Output Tables

### S222 Probe (cont'd)

| Parameter                                                                       | B+rD         | B+rD+D       | M+cM         | cM+D         | cwD         |
|---------------------------------------------------------------------------------|--------------|--------------|--------------|--------------|-------------|
| <b>P<sub>-</sub> (MPa)</b>                                                      | 2.66         | 2.66         | 2.66         | 2.66         | 0.17        |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                                     | 1113.2       | 1469.5       | 2044.1       | 1818.9       | 990.8       |
| <b>I<sub>p</sub> (mm)</b>                                                       | 73.2         | 73.2         | 73.2         | 73.2         | 76.1        |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)                                  | 2.1<br>3.4   | 2.1<br>3.4   | 2.1<br>3.4   | 2.1<br>3.4   | 6.3<br>4.1  |
| <b>prf (kHz)</b>                                                                | 5.96         | 5.96         | 5.96         | 5.96         | —           |
| <b>srr (Hz)</b>                                                                 | 22           | 22           | —            | —            | —           |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)                                 | 13.0<br>20.5 | 13.0<br>20.5 | 13.0<br>20.5 | 13.0<br>20.5 | 13.0<br>7.7 |
| <b>f<sub>awf</sub> (MHz)</b>                                                    | 2.49         | 2.49         | 2.49         | 2.49         | 2.51        |
| <b>APF (%)</b><br><b>(Acoustic Power</b><br><b>up Fraction)</b>                 | na           | na           | na           | na           | na          |
| <b>AIF (%)</b><br><b>(Acoustic</b><br><b>Initialization</b><br><b>Fraction)</b> | na           | na           | na           | na           | na          |
| <b>Maximum Power (mW)</b>                                                       | 241.43       | 185.81       | 95.63        | 95.63        | 194.47      |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                                       | 116.38       | 95.49        | 48.64        | 48.64        | 194.79      |
| <b>Power-up Mode</b>                                                            | na           | na           | na           | na           | na          |
| <b>Initialization Mode</b>                                                      | na           | na           | na           | na           | na          |
| <b>Acoustic Output</b><br><b>Freeze</b>                                         | yes          | yes          | yes          | yes          | yes         |
| <b>I<sub>tt</sub> (mm)</b>                                                      | 0.00         | 0.00         | 0.00         | 0.00         | 0.00        |
| <b>I<sub>ts</sub> (mm)</b>                                                      | contact      | contact      | contact      | contact      | contact     |

Table 17–146. S222 IEC Acoustic Output Information (cont'd)



## IEC Acoustic Output Tables

### S611 Probe

| Parameter                                                  | B             | M             | B+M           | D             | B+D           |
|------------------------------------------------------------|---------------|---------------|---------------|---------------|---------------|
| <b>P<sub>-</sub> (MPa)</b>                                 | 2.59          | 2.59          | 2.59          | 2.96          | 2.96          |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                | 295.29        | 199.40        | 396.26        | 2358.70       | 2555.56       |
| <b>I<sub>p</sub> (mm)</b>                                  | 37.10         | 37.10         | 37.10         | 45.00         | 45.00         |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)             | 1.14<br>2.22  | 1.14<br>2.22  | 1.14<br>2.22  | 1.33<br>1.94  | 1.33<br>1.94  |
| <b>prr (kHz)</b>                                           | –             | 0.85          | 0.85          | 1.49          | 1.49          |
| <b>srr (Hz)</b>                                            | 142.00        | –             | 142.00        | –             | 142.00        |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)            | 11.50<br>8.00 | 11.50<br>8.00 | 11.50<br>8.00 | 11.20<br>8.00 | 11.50<br>8.00 |
| <b>f<sub>awf</sub> (MHz)</b>                               | 4.43          | 4.43          | 4.43          | 4.98          | 4.98          |
| <b>APF (%)</b><br>(Acoustic Power<br>up Fraction)          | na            | na            | na            | na            | na            |
| <b>AIF (%)</b><br>(Acoustic<br>Initialization<br>Fraction) | na            | na            | na            | na            | na            |
| <b>Maximum Power (mW)</b>                                  | 20.88         | 2.25          | 16.17         | 40.55         | 54.47         |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                  | 22.69         | 2.50          | 17.63         | 44.56         | 59.68         |
| <b>Power-up Mode</b>                                       | na            | na            | na            | na            | na            |
| <b>Initialization Mode</b>                                 | na            | na            | na            | na            | na            |
| <b>Acoustic Output<br/>Freeze</b>                          | yes           | yes           | yes           | yes           | yes           |
| <b>I<sub>tt</sub> (mm)</b>                                 | 0.00          | 0.00          | 0.00          | 0.00          | 0.00          |
| <b>I<sub>ts</sub> (mm)</b>                                 | contact       | contact       | contact       | contact       | contact       |

Table 17–147. S611 IEC Acoustic Output Information

## IEC Acoustic Output Tables

### S611 Probe (cont'd)

| Parameter                                                  | B+rD          | B+rD+D        | M+cM          | cM+D          | cwD          |
|------------------------------------------------------------|---------------|---------------|---------------|---------------|--------------|
| <b>P<sub>-</sub> (MPa)</b>                                 | 2.59          | 2.96          | 2.10          | 2.10          | 0.17         |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                | 1092.16       | 1823.86       | 2994.90       | 2994.90       | 956.60       |
| <b>I<sub>p</sub> (mm)</b>                                  | 51.40         | 51.40         | 51.40         | 51.40         | 19.80        |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)             | 3.05<br>2.47  | 3.05<br>2.47  | 3.05<br>2.47  | 3.05<br>2.47  | 5.00<br>2.00 |
| <b>prf (kHz)</b>                                           | 5.96          | 5.96          | 5.96          | 5.96          | —            |
| <b>srr (Hz)</b>                                            | 66.00         | 66.00         | —             | —             | —            |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)            | 11.50<br>8.00 | 11.50<br>8.00 | 11.20<br>8.00 | 11.20<br>8.00 | 5.04<br>8.00 |
| <b>f<sub>awf</sub> (MHz)</b>                               | 4.02          | 4.02          | 4.02          | 4.02          | 4.00         |
| <b>APF (%)</b><br>(Acoustic Power<br>up Fraction)          | na            | na            | na            | na            | na           |
| <b>AIF (%)</b><br>(Acoustic<br>Initialization<br>Fraction) | na            | na            | na            | na            | na           |
| <b>Maximum Power (mW)</b>                                  | 86.18         | 70.32         | 72.26         | 72.26         | 63.31        |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                  | 93.67         | 76.68         | 81.04         | 81.04         | 186.21       |
| <b>Power-up Mode</b>                                       | na            | na            | na            | na            | na           |
| <b>Initialization Mode</b>                                 | na            | na            | na            | na            | na           |
| <b>Acoustic Output<br/>Freeze</b>                          | yes           | yes           | yes           | yes           | yes          |
| <b>I<sub>tt</sub> (mm)</b>                                 | 0.00          | 0.00          | 0.00          | 0.00          | 0.00         |
| <b>I<sub>ts</sub> (mm)</b>                                 | contact       | contact       | contact       | contact       | contact      |

Table 17–147. S611 IEC Acoustic Output Information (cont'd)



## IEC Acoustic Output Tables

### B510 Probe

| Parameter                                           | B            | M            | B+M          | D            | B+D          |
|-----------------------------------------------------|--------------|--------------|--------------|--------------|--------------|
| P <sub>-</sub> (MPa)                                | 3.22         | 3.22         | 3.22         | 2.57         | 3.22         |
| I <sub>spta</sub> (mW/cm <sup>2</sup> )             | 32.22        | 23.90        | 45.38        | 258.70       | 280.18       |
| I <sub>p</sub> (mm)                                 | 40.30        | 4.03         | 40.30        | 3.36         | 3.36         |
| W <sub>pb6</sub><br>(  ) (mm)<br>(⊥)(mm)            | 2.54<br>1.99 | 1.61<br>3.93 | 2.54<br>1.99 | 1.73<br>2.73 | 1.73<br>2.73 |
| prf (kHz)                                           | –            | 869.00       | 869.00       | 0.66         | 0.66         |
| srr (Hz)                                            | 26.00        | –            | 26.00        | –            | 26.00        |
| Output Beam Dim.<br>(  ) (mm)<br>(⊥)(mm)            | 2.54<br>1.99 | –<br>–       | 2.54<br>1.99 | –<br>–       | 2.54<br>1.99 |
| f <sub>awf</sub> (MHz)                              | 4.38         | 4.38         | 4.38         | 4.02         | 4.02         |
| APF (%)<br>(Acoustic Power<br>up Fraction)          | na           | na           | na           | na           | na           |
| AIF (%)<br>(Acoustic<br>Initialization<br>Fraction) | na           | na           | na           | na           | na           |
| Maximum Power (mW)                                  | 12.49        | 0.99         | 9.31         | 11.04        | 19.37        |
| I <sub>ob</sub> (mW/cm <sup>2</sup> )               | 12.49        | 0.99         | 9.31         | 11.04        | 19.37        |
| Power-up Mode                                       | na           | na           | na           | na           | na           |
| Initialization Mode                                 | na           | na           | na           | na           | na           |
| Acoustic Output<br>Freeze                           | yes          | yes          | yes          | yes          | yes          |
| I <sub>tt</sub> (mm)                                | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         |
| I <sub>ts</sub> (mm)                                | contact      | contact      | contact      | contact      | contact      |

Table 17–148. B510 IEC Acoustic Output Information

**B510 Probe (cont'd)**

| Parameter                                                  | B+rD         | B+rD+D       | M+cM         | cM+D         | cwD          |
|------------------------------------------------------------|--------------|--------------|--------------|--------------|--------------|
| <b>P<sub>-</sub> (MPa)</b>                                 | 3.22         | 3.22         | 1.77         | 1.77         | 0.11         |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                | 117.08       | 198.63       | 57.30        | 57.30        | 358.80       |
| <b>I<sub>p</sub> (mm)</b>                                  | 3.07         | 3.07         | 3.30         | 3.30         | 3.23         |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)             | 1.69<br>3.54 | 1.69<br>3.54 | 1.51<br>3.19 | 1.51<br>3.19 | 4.17<br>2.37 |
| <b>prf (kHz)</b>                                           | 0.28         | 0.28         | 2.60         | 2.60         | —            |
| <b>srr (Hz)</b>                                            | 20.00        | 20.00        | —            | —            | —            |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)            | 1.69<br>3.54 | 1.69<br>3.54 | —<br>—       | —<br>—       | —<br>—       |
| <b>f<sub>awf</sub> (MHz)</b>                               | 4.05         | 4.05         | 4.86         | 4.86         | 4.01         |
| <b>APF (%)</b><br>(Acoustic Power<br>up Fraction)          | na           | na           | na           | na           | na           |
| <b>AIF (%)</b><br>(Acoustic<br>Initialization<br>Fraction) | na           | na           | na           | na           | na           |
| <b>Maximum Power (mW)</b>                                  | 14.18        | 16.78        | 5.86         | 5.86         | 18.76        |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                  | 14.18        | 16.78        | 5.86         | 5.86         | 18.76        |
| <b>Power-up Mode</b>                                       | na           | na           | na           | na           | na           |
| <b>Initialization Mode</b>                                 | na           | na           | na           | na           | na           |
| <b>Acoustic Output<br/>Freeze</b>                          | yes          | yes          | yes          | yes          | yes          |
| <b>I<sub>tt</sub> (mm)</b>                                 | 0.00         | 0.00         | 0.00         | 0.00         | 0.00         |
| <b>I<sub>ts</sub> (mm)</b>                                 | contact      | contact      | contact      | contact      | contact      |

Table 17–148. B510 IEC Acoustic Output Information (cont'd)



## IEC Acoustic Output Tables

### P509 Probe

| Parameter                                                                       | B           | M           | B+M         | D           | B+D         |
|---------------------------------------------------------------------------------|-------------|-------------|-------------|-------------|-------------|
| <b>P<sub>-</sub> (MPa)</b>                                                      | 2.55        | 2.55        | 2.55        | 2.15        | 2.55        |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                                     | 544.9       | 90.9        | 454.1       | 624.9       | 988.2       |
| <b>I<sub>p</sub> (mm)</b>                                                       | 20.6        | 20.6        | 20.6        | 22.0        | 22.0        |
| <b>W<sub>pb6</sub></b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>                    | 5.4<br>0.9  | 5.4<br>0.9  | 5.4<br>0.9  | 5.4<br>1.0  | 5.4<br>1.0  |
| <b>prr (kHz)</b>                                                                | –           | 0.90        | 0.90        | 5.95        | 5.95        |
| <b>srr (Hz)</b>                                                                 | 43          | –           | 43          | –           | 43          |
| <b>Output Beam Dim.</b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>                   | 10.0<br>8.6 | 10.0<br>8.6 | 10.0<br>8.6 | 10.0<br>9.9 | 10.0<br>9.9 |
| <b>f<sub>awf</sub> (MHz)</b>                                                    | 4.93        | 4.93        | 4.93        | 4.02        | 40.2        |
| <b>APF (%)</b><br><b>(Acoustic Power</b><br><b>up Fraction)</b>                 | n/a         | n/a         | n/a         | n/a         | n/a         |
| <b>AIF (%)</b><br><b>(Acoustic</b><br><b>Initialization</b><br><b>Fraction)</b> | n/a         | n/a         | n/a         | n/a         | n/a         |
| <b>Maximum Power (mW)</b>                                                       | 11.03       | 3.72        | 11.08       | 29.55       | 36.9        |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                                       | 12.85       | 4.34        | 12.90       | 29.97       | 38.54       |
| <b>Power-up Mode</b>                                                            | n/a         | n/a         | n/a         | n/a         | n/a         |
| <b>Initialization Mode</b>                                                      | n/a         | n/a         | n/a         | n/a         | n/a         |
| <b>Acoustic Output</b><br><b>Freeze</b>                                         | yes         | yes         | yes         | yes         | yes         |
| <b>I<sub>tt</sub> (mm)</b>                                                      | 0.0         | 0.0         | 0.0         | 0.0         | 0.0         |
| <b>I<sub>ts</sub> (mm)</b>                                                      | contact     | contact     | contact     | contact     | contact     |

Table 17–149. P509 IEC Acoustic Output Information

## IEC Acoustic Output Tables

### P509 Probe (cont'd)

| Parameter                                                     | B+rD        | B+rD+D      | M+cM        | cM+D        | cwD         |
|---------------------------------------------------------------|-------------|-------------|-------------|-------------|-------------|
| <b>P<sub>-</sub> (MPa)</b>                                    | 2.55        | 2.55        | 2.55        | 2.55        | 0.11        |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                   | 963.5       | 975.8       | 819.8       | 819.8       | 401.9       |
| <b>I<sub>p</sub> (mm)</b>                                     | 21.3        | 21.3        | 21.3        | 21.3        | 24.2        |
| <b>W<sub>pb6</sub></b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b>  | 5.0<br>1.0  | 5.0<br>1.0  | 5.0<br>1.0  | 5.0<br>1.0  | 4.6<br>2.7  |
| <b>pr<sub>r</sub> (kHz)</b>                                   | 5.96        | 5.96        | 5.96        | 5.96        | —           |
| <b>srr (Hz)</b>                                               | 33          | 33          | —           | —           | —           |
| <b>Output Beam Dim.</b><br><b>(  ) (mm)</b><br><b>(⊥)(mm)</b> | 10.0<br>9.9 | 10.0<br>9.9 | 10.0<br>9.9 | 10.0<br>9.9 | 10.0<br>3.8 |
| <b>f<sub>awf</sub> (MHz)</b>                                  | 4.05        | 4.05        | 4.05        | 4.05        | 4.00        |
| <b>APF (%)</b><br><b>(Acoustic Power up Fraction)</b>         | n/a         | n/a         | n/a         | n/a         | n/a         |
| <b>AIF (%)</b><br><b>(Acoustic Initialization Fraction)</b>   | n/a         | n/a         | n/a         | n/a         | n/a         |
| <b>Maximum Power (mW)</b>                                     | 31.52       | 34.21       | 24.14       | 24.14       | 69.43       |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                     | 33.08       | 35.81       | 25.05       | 25.05       | 181.93      |
| <b>Power-up Mode</b>                                          | n/a         | n/a         | n/a         | n/a         | n/a         |
| <b>Initialization Mode</b>                                    | n/a         | n/a         | n/a         | n/a         | n/a         |
| <b>Acoustic Output Freeze</b>                                 | yes         | yes         | yes         | yes         | yes         |
| <b>I<sub>tt</sub> (mm)</b>                                    | 0.0         | 0.0         | 0.0         | 0.0         | 0.0         |
| <b>I<sub>ts</sub> (mm)</b>                                    | contact     | contact     | contact     | contact     | contact     |

Table 17–149. P509 IEC Acoustic Output Information (cont'd)



## IEC Acoustic Output Tables

---

### CWD2 Probe

| Parameter                                           | cwD          |
|-----------------------------------------------------|--------------|
| P <sub>-</sub> (MPa)                                | 0.16         |
| I <sub>spta</sub> (mW/cm <sup>2</sup> )             | 322.58       |
| I <sub>p</sub> (mm)                                 | 10.00        |
| W <sub>pb6</sub><br>(  ) (mm)<br>(⊥)(mm)            | 8.00<br>5.00 |
| prr (kHz)                                           | –            |
| srr (Hz)                                            | –            |
| Output Beam Dim.<br>(  ) (mm)<br>(⊥)(mm)            | –<br>–       |
| f <sub>awf</sub> (MHz)                              | 2.01         |
| APF (%)<br>(Acoustic Power<br>up Fraction)          | na           |
| AIF (%)<br>(Acoustic<br>Initialization<br>Fraction) | na           |
| Maximum Power (mW)                                  | 51.19        |
| I <sub>ob</sub> (mW/cm <sup>2</sup> )               | 57.97        |
| Power-up Mode                                       | na           |
| Initialization Mode                                 | na           |
| Acoustic Output<br>Freeze                           | yes          |
| I <sub>tt</sub> (mm)                                | 0.00         |
| I <sub>ts</sub> (mm)                                | contact      |

Table 17–150. CWD2 IEC Acoustic Output Information





---

*IEC Acoustic Output Tables*

---

**CWD5 Probe**

| Parameter                                                  | cwD          |
|------------------------------------------------------------|--------------|
| <b>P<sub>-</sub> (MPa)</b>                                 | 0.19         |
| <b>I<sub>spta</sub> (mW/cm<sup>2</sup>)</b>                | 686.50       |
| <b>I<sub>p</sub> (mm)</b>                                  | 10.00        |
| <b>W<sub>pb6</sub></b><br>(  ) (mm)<br>(⊥)(mm)             | 3.00<br>3.00 |
| <b>pr<sub>r</sub> (kHz)</b>                                | –            |
| <b>sr<sub>r</sub> (Hz)</b>                                 | –            |
| <b>Output Beam Dim.</b><br>(  ) (mm)<br>(⊥)(mm)            | –<br>–       |
| <b>f<sub>awf</sub> (MHz)</b>                               | 5.02         |
| <b>APF (%)</b><br>(Acoustic Power<br>up Fraction)          | na           |
| <b>AIF (%)</b><br>(Acoustic<br>Initialization<br>Fraction) | na           |
| <b>Maximum Power (mW)</b>                                  | 80.41        |
| <b>I<sub>ob</sub> (mW/cm<sup>2</sup>)</b>                  | 320.09       |
| <b>Power-up Mode</b>                                       | na           |
| <b>Initialization Mode</b>                                 | na           |
| <b>Acoustic Output<br/>Freeze</b>                          | yes          |
| <b>I<sub>tt</sub> (mm)</b>                                 | 0.00         |
| <b>I<sub>ts</sub> (mm)</b>                                 | contact      |

Table 17–151. CWD5 IEC Acoustic Output Information



## *IEC Acoustic Output Tables*

---

This page left blank intentionally.

# OB Tables

## List of OB Tables

|                                   |       |
|-----------------------------------|-------|
| Table 18–1. BD : Berkowitz .....  | 18–3  |
| Table 18–2. BPD : Campbell .....  | 18–3  |
| Table 18–3. BD : Campbell .....   | 18–4  |
| Table 18–4. CRL : Campbell .....  | 18–4  |
| Table 18–5. FL : Campbell .....   | 18–5  |
| Table 18–6. TAD : Eriksen .....   | 18–5  |
| Table 18–7. AC : Hadlock .....    | 18–6  |
| Table 18–8. BPD : Hadlock .....   | 18–6  |
| Table 18–9. CRL : Hadlock .....   | 18–7  |
| Table 18–10. FL : Hadlock .....   | 18–7  |
| Table 18–11. HC : Hadlock .....   | 18–8  |
| Table 18–12. AC : Hansmann .....  | 18–9  |
| Table 18–13. BPD : Hansmann ..... | 18–10 |
| Table 18–14. CRL : Hansmann ..... | 18–11 |
| Table 18–15. FL : Hansmann .....  | 18–12 |
| Table 18–16. GS : Hansmann .....  | 18–12 |
| Table 18–17. HC : Hansmann .....  | 18–13 |
| Table 18–18. OFD : Hansmann ..... | 18–14 |
| Table 18–19. TAD : Hansmann ..... | 18–15 |
| Table 18–20. ThD : Hansmann ..... | 18–15 |
| Table 18–21. GS : Hellman .....   | 18–16 |
| Table 18–22. AC : Jeanty .....    | 18–16 |
| Table 18–23. BD : Jeanty .....    | 18–17 |
| Table 18–24. BPD : Jeanty .....   | 18–17 |
| Table 18–25. CRL : Jeanty .....   | 18–18 |
| Table 18–26. FL : Jeanty .....    | 18–18 |
| Table 18–27. HC : Jeanty .....    | 18–19 |
| Table 18–28. BPD : Kurtz .....    | 18–19 |
| Table 18–29. CRL : Nelson .....   | 18–20 |
| Table 18–30. BPD : Osaka .....    | 18–20 |



## OB Tables List

---

### List of OB Tables

|                                                    |       |
|----------------------------------------------------|-------|
| Table 18–31. CRL : Osaka .....                     | 18–21 |
| Table 18–32. EFBW : Osaka .....                    | 18–21 |
| Table 18–33. FL : Osaka .....                      | 18–22 |
| Table 18–34. FTA : Osaka .....                     | 18–23 |
| Table 18–35. HL : Osaka .....                      | 18–23 |
| Table 18–36. BPD : Paris .....                     | 18–24 |
| Table 18–37. CRL : Paris .....                     | 18–24 |
| Table 18–38. FL : Paris .....                      | 18–25 |
| Table 18–39. Ft : Paris .....                      | 18–25 |
| Table 18–40. TAD : Paris .....                     | 18–26 |
| Table 18–41. BPD : Rempen .....                    | 18–26 |
| Table 18–42. CRL : Rempen .....                    | 18–27 |
| Table 18–43. GS : Rempen .....                     | 18–28 |
| Table 18–44. CRL : Robinson .....                  | 18–29 |
| Table 18–45. AC : Sostoa .....                     | 18–29 |
| Table 18–46. BD : Sostoa .....                     | 18–30 |
| Table 18–47. BPD : Sostoa .....                    | 18–30 |
| Table 18–48. FL : Sostoa .....                     | 18–31 |
| Table 18–49. HC : Sostoa .....                     | 18–31 |
| Table 18–50. OFD : Sostoa .....                    | 18–32 |
| Table 18–51. BPD : Tokyo .....                     | 18–32 |
| Table 18–52. CRL : Tokyo .....                     | 18–33 |
| Table 18–53. FL : Tokyo .....                      | 18–33 |
| Table 18–54. GS : Tokyo .....                      | 18–34 |
| Table 18–55. LV : Tokyo .....                      | 18–34 |
| Table 18–56. FL/HC .....                           | 18–35 |
| Table 18–57. HC/AC .....                           | 18–35 |
| Table 18–58. EFW : Tokyo Shinozuka .....           | 18–36 |
| Table 18–59. APTDxTTD (AxT): Tokyo Shinozuka ..... | 18–36 |
| Table 18–60. FL : Tokyo Shinozuka .....            | 18–37 |
| Table 18–61. AC : Tokyo Shinozuka .....            | 18–37 |
| Table 18–62. BPD : Tokyo Shinozuka .....           | 18–38 |
| Table 18–63. CRL : Tokyo Shinozuka .....           | 18–38 |
| Table 18–64. AC : Australia .....                  | 18–39 |
| Table 18–65. CRL : Australia .....                 | 18–39 |
| Table 18–66. BPD : Australia .....                 | 18–40 |

## OB Tables

BD : Berkowitz

Unit : BD (mm)  
Age (Day)  
SD (mm)

BPD : Campbell  
King's College Hosp.  
London  
(Am.J.obst.gynecol).  
Oct 1, 1982

Unit : BPD (mm)  
Age (Day)  
SD (mm)

| BD  | Age | SD | BD  | Age | SD |
|-----|-----|----|-----|-----|----|
| <13 | n/a | —  | 53  | 230 | 0  |
| 13  | 81  | 0  | 54  | 237 | 0  |
| 14  | 82  | 0  | 55  | 244 | 0  |
| 15  | 84  | 0  | 56  | 251 | 0  |
| 16  | 86  | 0  | 57  | 258 | 0  |
| 17  | 88  | 0  | 58  | 266 | 0  |
| 18  | 91  | 0  | 59  | 275 | 0  |
| 19  | 95  | 0  | >59 | n/a | —  |
| 20  | 98  | 0  |     |     |    |
| 21  | 102 | 0  |     |     |    |
| 22  | 105 | 0  |     |     |    |
| 23  | 109 | 0  |     |     |    |
| 24  | 110 | 0  |     |     |    |
| 25  | 112 | 0  |     |     |    |
| 26  | 116 | 0  |     |     |    |
| 27  | 120 | 0  |     |     |    |
| 28  | 124 | 0  |     |     |    |
| 29  | 128 | 0  |     |     |    |
| 30  | 130 | 0  |     |     |    |
| 31  | 132 | 0  |     |     |    |
| 32  | 135 | 0  |     |     |    |
| 33  | 138 | 0  |     |     |    |
| 34  | 142 | 0  |     |     |    |
| 35  | 145 | 0  |     |     |    |
| 36  | 149 | 0  |     |     |    |
| 37  | 152 | 0  |     |     |    |
| 38  | 155 | 0  |     |     |    |
| 39  | 159 | 0  |     |     |    |
| 40  | 162 | 0  |     |     |    |
| 41  | 166 | 0  |     |     |    |
| 42  | 169 | 0  |     |     |    |
| 43  | 173 | 0  |     |     |    |
| 44  | 179 | 0  |     |     |    |
| 45  | 185 | 0  |     |     |    |
| 46  | 191 | 0  |     |     |    |
| 47  | 197 | 0  |     |     |    |
| 48  | 202 | 0  |     |     |    |
| 49  | 207 | 0  |     |     |    |
| 50  | 212 | 0  |     |     |    |
| 51  | 217 | 0  |     |     |    |
| 52  | 223 | 0  |     |     |    |

Table 18–1. BD : Berkowitz

| BPD | Age | SD | BPD | Age | SD |
|-----|-----|----|-----|-----|----|
| <20 | n/a | —  | 60  | 173 | 4  |
| 20  | 80  | 3  | 61  | 175 | 4  |
| 21  | 82  | 3  | 62  | 177 | 4  |
| 22  | 84  | 3  | 63  | 180 | 4  |
| 23  | 86  | 3  | 64  | 182 | 4  |
| 24  | 89  | 3  | 65  | 184 | 4  |
| 25  | 91  | 3  | 66  | 187 | 5  |
| 26  | 93  | 3  | 67  | 189 | 5  |
| 27  | 96  | 3  | 68  | 192 | 5  |
| 28  | 98  | 3  | 69  | 195 | 5  |
| 29  | 100 | 3  | 70  | 198 | 5  |
| 30  | 103 | 3  | 71  | 201 | 5  |
| 31  | 105 | 3  | 72  | 204 | 5  |
| 32  | 107 | 3  | 73  | 207 | 5  |
| 33  | 110 | 3  | 74  | 210 | 5  |
| 34  | 112 | 3  | 75  | 213 | 5  |
| 35  | 114 | 3  | 76  | 216 | 5  |
| 36  | 117 | 3  | 77  | 219 | 5  |
| 37  | 119 | 3  | 78  | 222 | 5  |
| 38  | 121 | 3  | 79  | 225 | 5  |
| 39  | 124 | 3  | 80  | 228 | 5  |
| 40  | 126 | 3  | 81  | 231 | 5  |
| 41  | 128 | 3  | 82  | 234 | 5  |
| 42  | 131 | 3  | 83  | 237 | 5  |
| 43  | 133 | 4  | 84  | 240 | 5  |
| 44  | 135 | 4  | 85  | 243 | 5  |
| 45  | 138 | 4  | 86  | 246 | 5  |
| 46  | 140 | 4  | 87  | 249 | 5  |
| 47  | 142 | 4  | 88  | 252 | 5  |
| 48  | 145 | 4  | 89  | 255 | 5  |
| 49  | 147 | 4  | 90  | 258 | 5  |
| 50  | 149 | 4  | 91  | 261 | 5  |
| 51  | 152 | 4  | 92  | 264 | 5  |
| 52  | 154 | 4  | 93  | 267 | 5  |
| 53  | 156 | 4  | 94  | 270 | 5  |
| 54  | 159 | 4  | 95  | 273 | 5  |
| 55  | 161 | 4  | 96  | 276 | 5  |
| 56  | 163 | 4  | 97  | 279 | 5  |
| 57  | 166 | 4  | >97 | n/a | —  |
| 58  | 168 | 4  |     |     |    |
| 59  | 170 | 4  |     |     |    |

Table 18–2. BPD : Campbell



## OB Tables

BD : Campbell  
King's College Hosp.  
London  
(Am.J.obst.gynecol).  
Oct 1, 1982

Unit : BD (mm)  
Age (Day)  
SD (mm)

CRL : Campbell  
King's College Hosp.  
London  
(Am.J.obst.gynecol).  
Oct 1, 1982

Unit : CRL (mm)  
Age (Day)  
SD (mm)

| BD  | Age | SD | BD  | Age | SD |
|-----|-----|----|-----|-----|----|
| <14 | n/a | —  | 54  | 237 | 0  |
| 14  | 81  | 0  | 55  | 245 | 0  |
| 15  | 84  | 0  | 56  | 252 | 0  |
| 16  | 87  | 0  | 57  | 260 | 0  |
| 17  | 90  | 0  | 58  | 267 | 0  |
| 18  | 93  | 0  | 59  | 275 | 0  |
| 19  | 96  | 0  | >59 | n/a | —  |
| 20  | 99  | 0  |     |     |    |
| 21  | 102 | 0  |     |     |    |
| 22  | 105 | 0  |     |     |    |
| 23  | 108 | 0  |     |     |    |
| 24  | 111 | 0  |     |     |    |
| 25  | 114 | 0  |     |     |    |
| 26  | 117 | 0  |     |     |    |
| 27  | 120 | 0  |     |     |    |
| 28  | 123 | 0  |     |     |    |
| 29  | 126 | 0  |     |     |    |
| 30  | 129 | 0  |     |     |    |
| 31  | 132 | 0  |     |     |    |
| 32  | 135 | 0  |     |     |    |
| 33  | 138 | 0  |     |     |    |
| 34  | 141 | 0  |     |     |    |
| 35  | 145 | 0  |     |     |    |
| 36  | 148 | 0  |     |     |    |
| 37  | 151 | 0  |     |     |    |
| 38  | 154 | 0  |     |     |    |
| 39  | 157 | 0  |     |     |    |
| 40  | 160 | 0  |     |     |    |
| 41  | 163 | 0  |     |     |    |
| 42  | 166 | 0  |     |     |    |
| 43  | 170 | 0  |     |     |    |
| 44  | 175 | 0  |     |     |    |
| 45  | 181 | 0  |     |     |    |
| 46  | 186 | 0  |     |     |    |
| 47  | 192 | 0  |     |     |    |
| 48  | 198 | 0  |     |     |    |
| 49  | 203 | 0  |     |     |    |
| 50  | 209 | 0  |     |     |    |
| 51  | 215 | 0  |     |     |    |
| 52  | 222 | 0  |     |     |    |
| 53  | 230 | 0  |     |     |    |

Table 18–3. BD : Campbell

| CRL | Age | SD | CRL | Age | SD |
|-----|-----|----|-----|-----|----|
| <10 | n/a | —  | 50  | 79  | 7  |
| 10  | 49  | 5  | 51  | 80  | 7  |
| 11  | 49  | 5  | 52  | 81  | 7  |
| 12  | 50  | 5  | 53  | 81  | 7  |
| 13  | 51  | 5  | 54  | 82  | 7  |
| 14  | 52  | 5  | 55  | 82  | 7  |
| 15  | 53  | 5  | 56  | 83  | 7  |
| 16  | 54  | 5  | 57  | 83  | 7  |
| 17  | 55  | 5  | 58  | 84  | 7  |
| 18  | 56  | 5  | 59  | 85  | 7  |
| 19  | 57  | 5  | 60  | 85  | 7  |
| 20  | 58  | 5  | 61  | 86  | 7  |
| 21  | 59  | 5  | 62  | 86  | 7  |
| 22  | 59  | 6  | 63  | 87  | 7  |
| 23  | 60  | 6  | 64  | 88  | 7  |
| 24  | 61  | 6  | 65  | 88  | 7  |
| 25  | 62  | 6  | 66  | 89  | 7  |
| 26  | 63  | 6  | 67  | 89  | 7  |
| 27  | 64  | 6  | 68  | 90  | 7  |
| 28  | 65  | 6  | 69  | 90  | 7  |
| 29  | 66  | 6  | 70  | 91  | 7  |
| 30  | 67  | 6  | 71  | 92  | 7  |
| 31  | 68  | 6  | 72  | 92  | 7  |
| 32  | 69  | 6  | 73  | 93  | 7  |
| 33  | 70  | 7  | 74  | 93  | 7  |
| 34  | 70  | 7  | 75  | 94  | 7  |
| 35  | 71  | 7  | 76  | 95  | 7  |
| 36  | 71  | 7  | 77  | 95  | 7  |
| 37  | 72  | 7  | 78  | 95  | 7  |
| 38  | 72  | 7  | 79  | 96  | 7  |
| 39  | 73  | 7  | 80  | 96  | 7  |
| 40  | 74  | 7  | 81  | 96  | 7  |
| 41  | 74  | 7  | 82  | 97  | 7  |
| 42  | 75  | 7  | 83  | 97  | 7  |
| 43  | 75  | 7  | 84  | 97  | 7  |
| 44  | 76  | 7  | 85  | 98  | 7  |
| 45  | 76  | 7  | >85 | n/a | —  |
| 46  | 77  | 7  |     |     |    |
| 47  | 78  | 7  |     |     |    |
| 48  | 78  | 7  |     |     |    |
| 49  | 79  | 7  |     |     |    |

Table 18–4. CRL : Campbell

## OB Tables

FL : Campbell  
King's College Hosp.  
London  
(Am.J.obst.gynecol).  
Oct 1, 1982

Unit : FL (mm)  
Age (Day)  
SD (mm)

TAD : Eriksen

Unit : TAD (mm)  
Age (Day)  
SD (mm)

| FL  | Age | SD | FL  | Age | SD |
|-----|-----|----|-----|-----|----|
| <16 | n/a | —  | 56  | 203 | 5  |
| 16  | 100 | 4  | 57  | 205 | 5  |
| 17  | 102 | 4  | 58  | 208 | 5  |
| 18  | 104 | 4  | 59  | 211 | 5  |
| 19  | 106 | 4  | 60  | 214 | 5  |
| 20  | 109 | 4  | 61  | 218 | 5  |
| 21  | 111 | 4  | 62  | 222 | 5  |
| 22  | 113 | 4  | 63  | 226 | 5  |
| 23  | 116 | 4  | 64  | 230 | 5  |
| 24  | 118 | 4  | 65  | 234 | 5  |
| 25  | 120 | 4  | 66  | 238 | 5  |
| 26  | 123 | 4  | 67  | 242 | 5  |
| 27  | 125 | 4  | 68  | 246 | 5  |
| 28  | 127 | 4  | 69  | 250 | 5  |
| 29  | 130 | 4  | 70  | 255 | 5  |
| 30  | 132 | 4  | 71  | 260 | 5  |
| 31  | 134 | 4  | 72  | 265 | 5  |
| 32  | 137 | 4  | 73  | 270 | 5  |
| 33  | 139 | 4  | 74  | 275 | 5  |
| 34  | 141 | 4  | 75  | 280 | 5  |
| 35  | 144 | 4  | 76  | 285 | 5  |
| 36  | 146 | 4  | >76 | n/a | 5  |
| 37  | 148 | 4  |     |     |    |
| 38  | 151 | 5  |     |     |    |
| 39  | 154 | 5  |     |     |    |
| 40  | 157 | 5  |     |     |    |
| 41  | 160 | 5  |     |     |    |
| 42  | 163 | 5  |     |     |    |
| 43  | 166 | 5  |     |     |    |
| 44  | 169 | 5  |     |     |    |
| 45  | 172 | 5  |     |     |    |
| 46  | 175 | 5  |     |     |    |
| 47  | 178 | 5  |     |     |    |
| 48  | 181 | 5  |     |     |    |
| 40  | 184 | 5  |     |     |    |
| 50  | 186 | 5  |     |     |    |
| 51  | 189 | 5  |     |     |    |
| 52  | 192 | 5  |     |     |    |
| 53  | 194 | 5  |     |     |    |
| 54  | 197 | 5  |     |     |    |
| 55  | 200 | 5  |     |     |    |

Table 18–5. FL : Campbell

| TAD | Age | SD | TAD | Age | SD | TAD  | Age | SD |
|-----|-----|----|-----|-----|----|------|-----|----|
| <23 | n/a | —  | 63  | 171 | 0  | 104  | 261 | 0  |
| 23  | 91  | 0  | 64  | 173 | 0  | 105  | 263 | 0  |
| 24  | 93  | 0  | 65  | 175 | 0  | 106  | 266 | 0  |
| 25  | 95  | 0  | 66  | 177 | 0  | 107  | 268 | 0  |
| 26  | 97  | 0  | 67  | 179 | 0  | 108  | 270 | 0  |
| 27  | 99  | 0  | 68  | 182 | 0  | 109  | 273 | 0  |
| 28  | 101 | 0  | 69  | 184 | 0  | 110  | 275 | 0  |
| 29  | 103 | 0  | 70  | 186 | 0  | 111  | 277 | 0  |
| 30  | 105 | 0  | 71  | 188 | 0  | 112  | 280 | 0  |
| 31  | 107 | 0  | 72  | 190 | 0  | >112 | n/a | —  |
| 32  | 109 | 0  | 73  | 192 | 0  |      |     |    |
| 33  | 111 | 0  | 74  | 195 | 0  |      |     |    |
| 34  | 113 | 0  | 75  | 197 | 0  |      |     |    |
| 35  | 115 | 0  | 76  | 199 | 0  |      |     |    |
| 36  | 117 | 0  | 77  | 201 | 0  |      |     |    |
| 37  | 119 | 0  | 78  | 203 | 0  |      |     |    |
| 38  | 120 | 0  | 79  | 206 | 0  |      |     |    |
| 39  | 122 | 0  | 80  | 208 | 0  |      |     |    |
| 40  | 124 | 0  | 81  | 210 | 0  |      |     |    |
| 41  | 126 | 0  | 82  | 212 | 0  |      |     |    |
| 42  | 128 | 0  | 83  | 214 | 0  |      |     |    |
| 43  | 130 | 0  | 84  | 217 | 0  |      |     |    |
| 44  | 132 | 0  | 85  | 219 | 0  |      |     |    |
| 45  | 134 | 0  | 86  | 221 | 0  |      |     |    |
| 46  | 136 | 0  | 87  | 223 | 0  |      |     |    |
| 47  | 138 | 0  | 88  | 225 | 0  |      |     |    |
| 48  | 140 | 0  | 89  | 228 | 0  |      |     |    |
| 49  | 142 | 0  | 90  | 230 | 0  |      |     |    |
| 50  | 144 | 0  | 91  | 232 | 0  |      |     |    |
| 51  | 146 | 0  | 92  | 234 | 0  |      |     |    |
| 52  | 148 | 0  | 93  | 236 | 0  |      |     |    |
| 53  | 150 | 0  | 94  | 239 | 0  |      |     |    |
| 54  | 152 | 0  | 95  | 241 | 0  |      |     |    |
| 55  | 154 | 0  | 96  | 243 | 0  |      |     |    |
| 56  | 156 | 0  | 97  | 245 | 0  |      |     |    |
| 57  | 158 | 0  | 98  | 247 | 0  |      |     |    |
| 58  | 161 | 0  | 99  | 250 | 0  |      |     |    |
| 59  | 163 | 0  | 100 | 252 | 0  |      |     |    |
| 60  | 165 | 0  | 101 | 254 | 0  |      |     |    |
| 61  | 167 | 0  | 102 | 256 | 0  |      |     |    |
| 62  | 169 | 0  | 103 | 259 | 0  |      |     |    |

Table 18–6. TAD : Eriksen



## OB Tables

AC : Hadlock  
Hadlock, Radiology  
1984, Vol 152:497

Unit : AC (mm)  
Age (Week)  
2SD (Week)

BPD : Hadlock  
Hadlock, Radiology  
1984, Vol 152:497

Unit: BPD (mm)  
Age (Week)  
2SD (Week)

| AC  | Age  | 2SD   | AC   | Age  | 2SD   |
|-----|------|-------|------|------|-------|
| <50 | n/a  | —     | 230  | 27.4 | ± 2.2 |
| 50  | 12.0 | ± 1.7 | 235  | 27.8 | ± 2.2 |
| 55  | 12.4 | ± 1.7 | 240  | 28.3 | ± 2.2 |
| 60  | 12.8 | ± 1.7 | 245  | 28.7 | ± 2.2 |
| 65  | 13.2 | ± 1.7 | 250  | 29.2 | ± 2.2 |
| 70  | 13.6 | ± 1.7 | 255  | 29.7 | ± 2.2 |
| 75  | 14.0 | ± 1.7 | 258  | 30.0 | ± 2.2 |
| 80  | 14.4 | ± 1.7 | 259  | 30.1 | ± 3.0 |
| 85  | 14.8 | ± 1.7 | 260  | 30.2 | ± 3.0 |
| 90  | 15.2 | ± 1.7 | 265  | 30.6 | ± 3.0 |
| 95  | 15.6 | ± 1.7 | 270  | 31.1 | ± 3.0 |
| 100 | 16.0 | ± 1.7 | 275  | 31.6 | ± 3.0 |
| 105 | 16.4 | ± 1.7 | 280  | 32.0 | ± 3.0 |
| 110 | 16.9 | ± 1.7 | 285  | 32.5 | ± 3.0 |
| 115 | 17.3 | ± 1.7 | 290  | 33.0 | ± 3.0 |
| 120 | 17.7 | ± 1.7 | 295  | 33.5 | ± 3.0 |
| 123 | 17.9 | ± 1.7 | 300  | 34.0 | ± 3.0 |
| 124 | 18.0 | ± 2.1 | 305  | 34.5 | ± 3.0 |
| 125 | 18.1 | ± 2.1 | 310  | 34.9 | ± 3.0 |
| 130 | 18.5 | ± 2.1 | 315  | 35.4 | ± 3.0 |
| 135 | 19.0 | ± 2.1 | 320  | 35.9 | ± 3.0 |
| 140 | 19.4 | ± 2.1 | 321  | 36.0 | ± 3.1 |
| 145 | 19.8 | ± 2.1 | 325  | 36.4 | ± 3.1 |
| 150 | 20.2 | ± 2.1 | 330  | 36.9 | ± 3.1 |
| 155 | 20.7 | ± 2.1 | 335  | 37.4 | ± 3.1 |
| 160 | 21.1 | ± 2.1 | 340  | 37.9 | ± 3.1 |
| 165 | 21.5 | ± 2.1 | 345  | 38.4 | ± 3.1 |
| 170 | 22.0 | ± 2.1 | 350  | 38.9 | ± 3.1 |
| 175 | 22.4 | ± 2.1 | 355  | 39.4 | ± 3.1 |
| 180 | 22.9 | ± 2.1 | 360  | 39.9 | ± 3.1 |
| 185 | 23.3 | ± 2.1 | 365  | 40.4 | ± 3.1 |
| 190 | 23.7 | ± 2.1 | 370  | 40.9 | ± 3.1 |
| 192 | 23.9 | ± 2.1 | 375  | 41.4 | ± 3.1 |
| 193 | 24.0 | ± 2.2 | 380  | 42.0 | ± 3.1 |
| 195 | 24.2 | ± 2.2 | 385  | 42.5 | ± 3.1 |
| 200 | 24.6 | ± 2.2 | >385 | n/a  | —     |
| 205 | 25.1 | ± 2.2 |      |      |       |
| 210 | 25.5 | ± 2.2 |      |      |       |
| 215 | 26.0 | ± 2.2 |      |      |       |
| 220 | 26.4 | ± 2.2 |      |      |       |
| 225 | 26.9 | ± 2.2 |      |      |       |

Table 18–7. AC : Hadlock

| BPD | Age  | 2SD  | BPD | Age  | 2SD  | BPD  | Age  | 2SD  |
|-----|------|------|-----|------|------|------|------|------|
| <14 | n/a  | —    | 54  | 22.4 | ±1.7 | 95   | 38.7 | ±3.2 |
| 14  | 11.9 | ±1.2 | 55  | 22.8 | ±1.7 | 96   | 39.2 | ±3.2 |
| 15  | 12.1 | ±1.2 | 56  | 23.1 | ±1.7 | 97   | 39.7 | ±3.2 |
| 16  | 12.3 | ±1.2 | 57  | 23.4 | ±1.7 | 98   | 40.2 | ±3.2 |
| 17  | 12.5 | ±1.2 | 58  | 23.8 | ±1.7 | 99   | 40.6 | ±3.2 |
| 18  | 12.8 | ±1.2 | 59  | 24.1 | ±2.2 | 100  | 41.1 | ±3.2 |
| 19  | 13.0 | ±1.2 | 60  | 24.5 | ±2.2 | 101  | 41.6 | ±3.2 |
| 20  | 13.2 | ±1.2 | 61  | 24.8 | ±2.2 | 102  | 42.1 | ±3.2 |
| 21  | 13.4 | ±1.2 | 62  | 25.2 | ±2.2 | 103  | 42.6 | ±3.2 |
| 22  | 13.6 | ±1.2 | 63  | 25.5 | ±2.2 | >103 | n/a  | —    |
| 23  | 13.8 | ±1.2 | 64  | 25.9 | ±2.2 |      |      |      |
| 24  | 14.1 | ±1.2 | 65  | 26.3 | ±2.2 |      |      |      |
| 25  | 14.3 | ±1.2 | 66  | 26.6 | ±2.2 |      |      |      |
| 26  | 14.5 | ±1.2 | 67  | 27.0 | ±2.2 |      |      |      |
| 27  | 14.8 | ±1.2 | 68  | 27.4 | ±2.2 |      |      |      |
| 28  | 15.0 | ±1.2 | 69  | 27.7 | ±2.2 |      |      |      |
| 29  | 15.2 | ±1.2 | 70  | 28.1 | ±2.2 |      |      |      |
| 30  | 15.5 | ±1.2 | 71  | 28.5 | ±2.2 |      |      |      |
| 31  | 15.7 | ±1.2 | 72  | 28.9 | ±2.2 |      |      |      |
| 32  | 16.0 | ±1.2 | 73  | 29.3 | ±2.2 |      |      |      |
| 33  | 16.3 | ±1.2 | 74  | 29.7 | ±2.2 |      |      |      |
| 34  | 16.5 | ±1.2 | 75  | 30.1 | ±3.1 |      |      |      |
| 35  | 16.8 | ±1.2 | 76  | 30.5 | ±3.1 |      |      |      |
| 36  | 17.0 | ±1.2 | 77  | 30.9 | ±3.1 |      |      |      |
| 37  | 17.3 | ±1.2 | 78  | 31.3 | ±3.1 |      |      |      |
| 38  | 17.6 | ±1.2 | 79  | 31.7 | ±3.1 |      |      |      |
| 39  | 17.9 | ±1.2 | 80  | 32.1 | ±3.1 |      |      |      |
| 40  | 18.1 | ±1.7 | 81  | 32.5 | ±3.1 |      |      |      |
| 41  | 18.4 | ±1.7 | 82  | 33.0 | ±3.1 |      |      |      |
| 42  | 18.7 | ±1.7 | 83  | 33.4 | ±3.1 |      |      |      |
| 43  | 19.0 | ±1.7 | 84  | 33.8 | ±3.1 |      |      |      |
| 44  | 19.3 | ±1.7 | 85  | 34.2 | ±3.1 |      |      |      |
| 45  | 19.6 | ±1.7 | 86  | 34.7 | ±3.1 |      |      |      |
| 46  | 19.9 | ±1.7 | 87  | 35.1 | ±3.1 |      |      |      |
| 47  | 20.2 | ±1.7 | 88  | 35.6 | ±3.1 |      |      |      |
| 48  | 20.5 | ±1.7 | 89  | 36.0 | ±3.2 |      |      |      |
| 49  | 20.8 | ±1.7 | 90  | 36.5 | ±3.2 |      |      |      |
| 50  | 21.1 | ±1.7 | 91  | 36.9 | ±3.2 |      |      |      |
| 51  | 21.5 | ±1.7 | 92  | 37.4 | ±3.2 |      |      |      |
| 52  | 21.8 | ±1.7 | 93  | 37.8 | ±3.2 |      |      |      |
| 53  | 22.1 | ±1.7 | 94  | 38.3 | ±3.2 |      |      |      |

Table 18–8 BPD : Hadlock  
Variability of GA estimate by BPD at term is  
± 2SD (6Wks)



OB Tables

CRL : Hadlock  
Hadlock, Radiology  
1992, Vol. 182:501

Unit : CRL (mm)  
Age (Week)  
SD (Week)

FL : Hadlock  
Hadlock, Radiology  
1984, Vol 152:497

Unit : FL (mm)  
Age (Week)  
2SD (Week)

| CRL | Age  | SD    | CRL | Age  | SD    | CRL  | Age  | SD    |
|-----|------|-------|-----|------|-------|------|------|-------|
| <2  | n/a  | —     | 42  | 11.1 | ± 0.5 | 83   | 14.2 | ± 0.6 |
| 2   | 5.7  | ± 0.3 | 43  | 11.2 | ± 0.5 | 84   | 14.3 | ± 0.6 |
| 3   | 5.9  | ± 0.3 | 44  | 11.2 | ± 0.5 | 85   | 14.4 | ± 0.6 |
| 4   | 6.1  | ± 0.3 | 45  | 11.3 | ± 0.5 | 86   | 14.5 | ± 0.6 |
| 5   | 6.2  | ± 0.3 | 46  | 11.4 | ± 0.5 | 87   | 14.6 | ± 0.6 |
| 6   | 6.4  | ± 0.3 | 47  | 11.5 | ± 0.5 | 88   | 14.7 | ± 0.7 |
| 7   | 6.6  | ± 0.3 | 48  | 11.6 | ± 0.5 | 89   | 14.8 | ± 0.7 |
| 8   | 6.7  | ± 0.3 | 49  | 11.7 | ± 0.5 | 90   | 14.9 | ± 0.7 |
| 9   | 6.9  | ± 0.3 | 50  | 11.7 | ± 0.5 | 91   | 15.0 | ± 0.7 |
| 10  | 7.1  | ± 0.3 | 51  | 11.8 | ± 0.5 | 92   | 15.1 | ± 0.7 |
| 11  | 7.2  | ± 0.3 | 52  | 11.9 | ± 0.5 | 93   | 15.2 | ± 0.7 |
| 12  | 7.4  | ± 0.3 | 53  | 12.0 | ± 0.5 | 94   | 15.3 | ± 0.7 |
| 13  | 7.5  | ± 0.3 | 54  | 12.0 | ± 0.5 | 95   | 15.3 | ± 0.7 |
| 14  | 7.7  | ± 0.3 | 55  | 12.1 | ± 0.5 | 96   | 15.4 | ± 0.7 |
| 15  | 7.9  | ± 0.4 | 56  | 12.2 | ± 0.5 | 97   | 15.5 | ± 0.7 |
| 16  | 8.0  | ± 0.4 | 57  | 12.3 | ± 0.5 | 98   | 15.6 | ± 0.7 |
| 17  | 8.1  | ± 0.4 | 58  | 12.3 | ± 0.5 | 99   | 15.7 | ± 0.7 |
| 18  | 8.3  | ± 0.4 | 59  | 12.4 | ± 0.6 | 100  | 15.9 | ± 0.7 |
| 19  | 8.4  | ± 0.4 | 60  | 12.5 | ± 0.6 | 101  | 16.0 | ± 0.7 |
| 20  | 8.6  | ± 0.4 | 61  | 12.6 | ± 0.6 | 102  | 16.1 | ± 0.7 |
| 21  | 8.7  | ± 0.4 | 62  | 12.6 | ± 0.6 | 103  | 16.2 | ± 0.7 |
| 22  | 8.9  | ± 0.4 | 63  | 12.7 | ± 0.6 | 104  | 16.3 | ± 0.7 |
| 23  | 9.0  | ± 0.4 | 64  | 12.8 | ± 0.6 | 105  | 16.4 | ± 0.7 |
| 24  | 9.1  | ± 0.4 | 65  | 12.8 | ± 0.6 | 106  | 16.5 | ± 0.7 |
| 25  | 9.2  | ± 0.4 | 66  | 12.9 | ± 0.6 | 107  | 16.6 | ± 0.7 |
| 26  | 9.4  | ± 0.4 | 67  | 13.0 | ± 0.6 | 108  | 16.7 | ± 0.7 |
| 27  | 9.5  | ± 0.4 | 68  | 13.1 | ± 0.6 | 109  | 16.8 | ± 0.7 |
| 28  | 9.6  | ± 0.4 | 69  | 13.1 | ± 0.6 | 110  | 16.9 | ± 0.8 |
| 29  | 9.7  | ± 0.4 | 70  | 13.2 | ± 0.6 | 111  | 17.0 | ± 0.8 |
| 30  | 9.9  | ± 0.4 | 71  | 13.3 | ± 0.6 | 112  | 17.1 | ± 0.8 |
| 31  | 10.0 | ± 0.4 | 72  | 13.4 | ± 0.6 | 113  | 17.2 | ± 0.8 |
| 32  | 10.1 | ± 0.5 | 73  | 13.4 | ± 0.6 | 114  | 17.3 | ± 0.8 |
| 33  | 10.2 | ± 0.5 | 74  | 13.5 | ± 0.6 | 115  | 17.4 | ± 0.8 |
| 34  | 10.3 | ± 0.5 | 75  | 13.6 | ± 0.6 | 116  | 17.5 | ± 0.8 |
| 35  | 10.4 | ± 0.5 | 76  | 13.7 | ± 0.6 | 117  | 17.6 | ± 0.8 |
| 36  | 10.5 | ± 0.5 | 77  | 13.8 | ± 0.6 | 118  | 17.7 | ± 0.8 |
| 37  | 10.6 | ± 0.5 | 78  | 13.8 | ± 0.6 | 119  | 17.8 | ± 0.8 |
| 38  | 10.7 | ± 0.5 | 79  | 13.9 | ± 0.6 | 120  | 17.9 | ± 0.8 |
| 39  | 10.8 | ± 0.5 | 80  | 14.0 | ± 0.6 | 121  | 18.0 | ± 0.8 |
| 40  | 10.9 | ± 0.5 | 81  | 14.1 | ± 0.6 | >121 | n/a  | —     |
| 41  | 11.0 | ± 0.5 | 82  | 14.2 | ± 0.6 |      |      |       |

Table 18–9. CRL : Hadlock

| FL | Age  | 2SD  | FL  | Age  | 2SD  |
|----|------|------|-----|------|------|
| <6 | n/a  | —    | 46  | 25.3 | ±2.1 |
| 6  | 11.9 | ±1.4 | 47  | 25.7 | ±2.1 |
| 7  | 12.2 | ±1.4 | 48  | 26.1 | ±2.1 |
| 8  | 12.4 | ±1.4 | 49  | 26.5 | ±2.1 |
| 9  | 12.7 | ±1.4 | 50  | 26.9 | ±2.1 |
| 10 | 13.0 | ±1.4 | 51  | 27.3 | ±2.1 |
| 11 | 13.3 | ±1.4 | 52  | 27.7 | ±2.1 |
| 12 | 13.5 | ±1.4 | 53  | 28.2 | ±2.1 |
| 13 | 13.8 | ±1.4 | 54  | 28.6 | ±2.1 |
| 14 | 14.1 | ±1.4 | 55  | 29.0 | ±2.1 |
| 15 | 14.4 | ±1.4 | 56  | 29.5 | ±2.1 |
| 16 | 14.7 | ±1.4 | 57  | 29.9 | ±2.1 |
| 17 | 15.0 | ±1.4 | 58  | 30.3 | ±3.0 |
| 18 | 15.3 | ±1.4 | 59  | 30.8 | ±3.0 |
| 19 | 15.6 | ±1.4 | 60  | 31.2 | ±3.0 |
| 20 | 16.0 | ±1.4 | 61  | 31.7 | ±3.0 |
| 21 | 16.3 | ±1.4 | 62  | 32.1 | ±3.0 |
| 22 | 16.6 | ±1.4 | 63  | 32.6 | ±3.0 |
| 23 | 16.9 | ±1.4 | 64  | 33.1 | ±3.0 |
| 24 | 17.2 | ±1.4 | 65  | 33.5 | ±3.0 |
| 25 | 17.6 | ±1.4 | 66  | 34.0 | ±3.0 |
| 26 | 17.9 | ±1.4 | 67  | 34.5 | ±3.0 |
| 27 | 18.2 | ±1.8 | 68  | 34.9 | ±3.0 |
| 28 | 18.6 | ±1.8 | 69  | 35.4 | ±3.0 |
| 29 | 18.9 | ±1.8 | 70  | 35.9 | ±3.0 |
| 30 | 19.3 | ±1.8 | 71  | 36.4 | ±3.1 |
| 31 | 19.6 | ±1.8 | 72  | 36.9 | ±3.1 |
| 32 | 20.0 | ±1.8 | 73  | 37.4 | ±3.1 |
| 33 | 20.3 | ±1.8 | 74  | 37.9 | ±3.1 |
| 34 | 20.7 | ±1.8 | 75  | 38.4 | ±3.1 |
| 35 | 21.0 | ±1.8 | 76  | 38.9 | ±3.1 |
| 36 | 21.4 | ±1.8 | 77  | 39.4 | ±3.1 |
| 37 | 21.8 | ±1.8 | 78  | 39.9 | ±3.1 |
| 38 | 22.2 | ±1.8 | 79  | 40.4 | ±3.1 |
| 39 | 22.5 | ±1.8 | 80  | 40.9 | ±3.1 |
| 40 | 22.9 | ±1.8 | 81  | 41.4 | ±3.1 |
| 41 | 23.3 | ±1.8 | 82  | 42.0 | ±3.1 |
| 42 | 23.7 | ±1.8 | 83  | 42.5 | ±3.1 |
| 43 | 24.1 | ±2.1 | >83 | n/a  | —    |
| 44 | 24.5 | ±2.1 |     |      |      |
| 45 | 24.9 | ±2.1 |     |      |      |

Table 18–10. FL : Hadlock



## OB Tables

HC : Hadlock                      Unit :    HC (mm)  
Hadlock, Radiology              Age (Week)  
1984, Vol 152:497                2SD (Week)

| HC  | Age  | 2SD   | HC   | Age  | 2SD   |
|-----|------|-------|------|------|-------|
| <55 | n/a  | —     | 245  | 26.6 | ± 2.1 |
| 55  | 12.0 | ± 1.2 | 250  | 27.1 | ± 2.1 |
| 60  | 12.3 | ± 1.2 | 255  | 27.7 | ± 2.1 |
| 65  | 12.6 | ± 1.2 | 260  | 28.3 | ± 2.1 |
| 70  | 12.8 | ± 1.2 | 265  | 28.9 | ± 2.1 |
| 75  | 13.1 | ± 1.2 | 270  | 29.4 | ± 2.1 |
| 80  | 13.4 | ± 1.2 | 274  | 29.9 | ± 2.1 |
| 85  | 13.7 | ± 1.2 | 275  | 30.0 | ± 3.0 |
| 90  | 14.0 | ± 1.2 | 280  | 30.7 | ± 3.0 |
| 95  | 14.3 | ± 1.2 | 285  | 31.3 | ± 3.0 |
| 100 | 14.7 | ± 1.2 | 290  | 31.9 | ± 3.0 |
| 105 | 15.0 | ± 1.2 | 295  | 32.6 | ± 3.0 |
| 110 | 15.3 | ± 1.2 | 300  | 33.3 | ± 3.0 |
| 115 | 15.6 | ± 1.2 | 305  | 33.9 | ± 3.0 |
| 120 | 16.0 | ± 1.2 | 310  | 34.6 | ± 3.0 |
| 125 | 16.3 | ± 1.2 | 315  | 35.3 | ± 3.0 |
| 130 | 16.6 | ± 1.2 | 319  | 35.9 | ± 3.0 |
| 135 | 17.0 | ± 1.2 | 320  | 36.1 | ± 2.7 |
| 140 | 17.3 | ± 1.2 | 325  | 36.8 | ± 2.7 |
| 145 | 17.7 | ± 1.2 | 330  | 37.6 | ± 2.7 |
| 149 | 18.0 | ± 1.2 | 335  | 38.3 | ± 2.7 |
| 150 | 18.1 | ± 1.5 | 340  | 39.1 | ± 2.7 |
| 155 | 18.4 | ± 1.5 | 345  | 39.9 | ± 2.7 |
| 160 | 18.8 | ± 1.5 | 350  | 40.7 | ± 2.7 |
| 165 | 19.2 | ± 1.5 | 355  | 41.6 | ± 2.7 |
| 170 | 19.6 | ± 1.5 | 360  | 42.4 | ± 2.7 |
| 175 | 20.0 | ± 1.5 | >360 | n/a  | —     |
| 180 | 20.4 | ± 1.5 |      |      |       |
| 185 | 20.8 | ± 1.5 |      |      |       |
| 190 | 21.3 | ± 1.5 |      |      |       |
| 195 | 21.7 | ± 1.5 |      |      |       |
| 200 | 22.2 | ± 1.5 |      |      |       |
| 205 | 22.6 | ± 1.5 |      |      |       |
| 210 | 23.1 | ± 1.5 |      |      |       |
| 215 | 23.6 | ± 1.5 |      |      |       |
| 219 | 23.9 | ± 1.5 |      |      |       |
| 220 | 24.0 | ± 2.1 |      |      |       |
| 225 | 24.5 | ± 2.1 |      |      |       |
| 230 | 25.0 | ± 2.1 |      |      |       |
| 235 | 25.5 | ± 2.1 |      |      |       |
| 240 | 26.1 | ± 2.1 |      |      |       |

Table 18–11. HC : Hadlock

AC : Hansmann

Hansmann : M and AI : Geburtsh, u, Frauenheilk 39 : 656,1979

Unit : AC (mm)      Age (Weeks/Days)      SD (mm)

| AC  | Age   | SD | AC  | Age   | SD | AC  | Age   | SD | AC  | Age   | SD | AC  | Age   | SD | AC   | Age   | SD |
|-----|-------|----|-----|-------|----|-----|-------|----|-----|-------|----|-----|-------|----|------|-------|----|
| <53 | n/a   | —  | 97  | 15W1D | 0  | 142 | 19W5D | 0  | 187 | 24W5D | 0  | 232 | 29W2D | 0  | 277  | 34W1D | 0  |
| 53  | 11W1D | 0  | 98  | 15W2D | 0  | 143 | 19W6D | 0  | 188 | 24W5D | 0  | 233 | 29W3D | 0  | 278  | 34W2D | 0  |
| 54  | 11W2D | 0  | 99  | 15W2D | 0  | 144 | 20W0D | 0  | 189 | 24W6D | 0  | 234 | 29W3D | 0  | 279  | 34W2D | 0  |
| 55  | 11W2D | 0  | 100 | 15W3D | 0  | 145 | 20W1D | 0  | 190 | 25W0D | 0  | 235 | 29W4D | 0  | 280  | 34W3D | 0  |
| 56  | 11W3D | 0  | 101 | 15W4D | 0  | 146 | 20W2D | 0  | 191 | 25W1D | 0  | 236 | 29W4D | 0  | 281  | 34W3D | 0  |
| 57  | 11W3D | 0  | 102 | 15W5D | 0  | 147 | 20W2D | 0  | 192 | 25W2D | 0  | 237 | 29W5D | 0  | 282  | 34W4D | 0  |
| 58  | 11W4D | 0  | 103 | 15W5D | 0  | 148 | 20W3D | 0  | 193 | 25W2D | 0  | 238 | 29W6D | 0  | 283  | 34W4D | 0  |
| 59  | 11W4D | 0  | 104 | 15W6D | 0  | 149 | 20W3D | 0  | 194 | 25W3D | 0  | 239 | 30W0D | 0  | 284  | 34W5D | 0  |
| 60  | 11W5D | 0  | 105 | 16W0D | 0  | 150 | 20W4D | 0  | 195 | 25W4D | 0  | 240 | 30W1D | 0  | 285  | 34W6D | 0  |
| 61  | 11W6D | 0  | 106 | 16W0D | 0  | 151 | 20W4D | 0  | 196 | 25W4D | 0  | 241 | 30W2D | 0  | 286  | 35W0D | 0  |
| 62  | 12W0D | 0  | 107 | 16W1D | 0  | 152 | 20W5D | 0  | 197 | 25W5D | 0  | 242 | 30W3D | 0  | 287  | 35W1D | 0  |
| 63  | 12W1D | 0  | 108 | 16W2D | 0  | 153 | 20W6D | 0  | 198 | 25W5D | 0  | 243 | 30W3D | 0  | 288  | 35W2D | 0  |
| 64  | 12W2D | 0  | 109 | 16W3D | 0  | 154 | 21W0D | 0  | 199 | 25W6D | 0  | 244 | 30W4D | 0  | 289  | 35W3D | 0  |
| 65  | 12W2D | 0  | 110 | 16W3D | 0  | 155 | 21W1D | 0  | 200 | 26W0D | 0  | 245 | 30W5D | 0  | 290  | 35W3D | 0  |
| 66  | 12W3D | 0  | 111 | 16W4D | 0  | 156 | 21W2D | 0  | 201 | 26W0D | 0  | 246 | 30W6D | 0  | 291  | 35W4D | 0  |
| 67  | 12W3D | 0  | 112 | 16W5D | 0  | 157 | 21W2D | 0  | 202 | 26W1D | 0  | 247 | 30W6D | 0  | 292  | 35W5D | 0  |
| 68  | 12W4D | 0  | 113 | 16W6D | 0  | 158 | 21W3D | 0  | 203 | 26W2D | 0  | 248 | 31W0D | 0  | 293  | 35W6D | 0  |
| 69  | 12W5D | 0  | 114 | 16W6D | 0  | 159 | 21W3D | 0  | 204 | 26W3D | 0  | 249 | 31W1D | 0  | 294  | 35W6D | 0  |
| 70  | 12W5D | 0  | 115 | 17W0D | 0  | 160 | 21W4D | 0  | 205 | 26W3D | 0  | 250 | 31W2D | 0  | 295  | 36W0D | 0  |
| 71  | 12W6D | 0  | 116 | 17W1D | 0  | 161 | 21W4D | 0  | 206 | 26W4D | 0  | 251 | 31W3D | 0  | 296  | 36W1D | 0  |
| 72  | 12W6D | 0  | 117 | 17W2D | 0  | 162 | 21W5D | 0  | 207 | 26W5D | 0  | 252 | 31W3D | 0  | 297  | 36W2D | 0  |
| 73  | 13W0D | 0  | 118 | 17W2D | 0  | 163 | 21W6D | 0  | 208 | 26W6D | 0  | 253 | 31W4D | 0  | 298  | 36W2D | 0  |
| 74  | 13W0D | 0  | 119 | 17W3D | 0  | 164 | 22W0D | 0  | 209 | 26W6D | 0  | 254 | 31W5D | 0  | 299  | 36W3D | 0  |
| 75  | 13W1D | 0  | 120 | 17W3D | 0  | 165 | 22W1D | 0  | 210 | 27W0D | 0  | 255 | 31W6D | 0  | 300  | 36W3D | 0  |
| 76  | 13W2D | 0  | 121 | 17W4D | 0  | 166 | 22W2D | 0  | 211 | 27W1D | 0  | 256 | 31W6D | 0  | 301  | 36W4D | 0  |
| 77  | 13W2D | 0  | 122 | 17W4D | 0  | 167 | 22W3D | 0  | 212 | 27W2D | 0  | 257 | 32W0D | 0  | 302  | 36W4D | 0  |
| 78  | 13W3D | 0  | 123 | 17W5D | 0  | 168 | 22W4D | 0  | 213 | 27W2D | 0  | 258 | 32W1D | 0  | 303  | 36W5D | 0  |
| 79  | 13W3D | 0  | 124 | 17W6D | 0  | 169 | 22W5D | 0  | 214 | 27W3D | 0  | 259 | 32W2D | 0  | 304  | 36W6D | 0  |
| 80  | 13W4D | 0  | 125 | 18W0D | 0  | 170 | 22W5D | 0  | 215 | 27W4D | 0  | 260 | 32W2D | 0  | 305  | 37W0D | 0  |
| 81  | 13W4D | 0  | 126 | 18W1D | 0  | 171 | 22W6D | 0  | 216 | 27W4D | 0  | 261 | 32W3D | 0  | 306  | 37W1D | 0  |
| 82  | 13W5D | 0  | 127 | 18W2D | 0  | 172 | 23W0D | 0  | 217 | 27W5D | 0  | 262 | 32W3D | 0  | 307  | 37W2D | 0  |
| 83  | 13W6D | 0  | 128 | 18W3D | 0  | 173 | 23W1D | 0  | 218 | 27W5D | 0  | 263 | 32W4D | 0  | 308  | 37W3D | 0  |
| 84  | 14W0D | 0  | 129 | 18W3D | 0  | 174 | 23W2D | 0  | 219 | 27W6D | 0  | 264 | 32W4D | 0  | 309  | 37W3D | 0  |
| 85  | 14W1D | 0  | 130 | 18W4D | 0  | 175 | 23W2D | 0  | 220 | 28W0D | 0  | 265 | 32W5D | 0  | 310  | 37W4D | 0  |
| 86  | 14W2D | 0  | 131 | 18W5D | 0  | 176 | 23W3D | 0  | 221 | 28W0D | 0  | 266 | 32W6D | 0  | 311  | 37W5D | 0  |
| 87  | 14W2D | 0  | 132 | 18W6D | 0  | 177 | 23W3D | 0  | 222 | 28W1D | 0  | 267 | 33W0D | 0  | 312  | 37W6D | 0  |
| 88  | 14W3D | 0  | 133 | 18W6D | 0  | 178 | 23W4D | 0  | 223 | 28W2D | 0  | 268 | 33W1D | 0  | 313  | 37W6D | 0  |
| 89  | 14W3D | 0  | 134 | 19W0D | 0  | 179 | 23W4D | 0  | 224 | 28W3D | 0  | 269 | 33W2D | 0  | 314  | 38W0D | 0  |
| 90  | 14W4D | 0  | 135 | 19W1D | 0  | 180 | 23W5D | 0  | 225 | 28W4D | 0  | 270 | 33W3D | 0  | 315  | 38W1D | 0  |
| 91  | 14W5D | 0  | 136 | 19W2D | 0  | 181 | 23W6D | 0  | 226 | 28W5D | 0  | 271 | 33W3D | 0  | 316  | 38W2D | 0  |
| 92  | 14W5D | 0  | 137 | 19W2D | 0  | 182 | 24W0D | 0  | 227 | 28W5D | 0  | 272 | 33W4D | 0  | 317  | 38W4D | 0  |
| 93  | 14W6D | 0  | 138 | 19W3D | 0  | 183 | 24W1D | 0  | 228 | 28W6D | 0  | 273 | 33W5D | 0  | 318  | 38W5D | 0  |
| 94  | 14W6D | 0  | 139 | 19W3D | 0  | 184 | 24W2D | 0  | 229 | 29W0D | 0  | 274 | 33W6D | 0  | 319  | 39W0D | 0  |
| 95  | 15W0D | 0  | 140 | 19W4D | 0  | 185 | 24W3D | 0  | 230 | 29W1D | 0  | 275 | 33W6D | 0  | 320  | 39W1D | 0  |
| 96  | 15W0D | 0  | 141 | 19W4D | 0  | 186 | 24W4D | 0  | 231 | 29W2D | 0  | 276 | 34W0D | 0  | >320 | n/a   | —  |

Table 18–12. AC : Hansmann



## OB Tables

BPD : Hansmann

Unit : BPD (mm)

2SD = mm

Age (Week/sDays)

2SD (mm or day)

2SD = day

| BPD | Age   | 2SD | BPD  | Age   | 2SD | BPD | Age   | 2SD | BPD | Age   | 2SD | BPD  | Age   | 2SD |
|-----|-------|-----|------|-------|-----|-----|-------|-----|-----|-------|-----|------|-------|-----|
| <14 | n/a   | —   | 58   | 22W2D | 5   | <14 | n/a   | —   | 58  | 22W5D | 9   | 103  | 40W0D | 19  |
| 14  | 10W0D | 0   | 59   | 22W4D | 5   | 14  | 9W1D  | 7   | 59  | 23W0D | 10  | 104  | 40W1D | 19  |
| 15  | 10W1D | 0   | 60   | 22W6D | 5   | 15  | 9W3D  | 7   | 60  | 23W2D | 10  | 105  | 40W2D | 17  |
| 16  | 10W2D | 0   | 61   | 23W1D | 5   | 16  | 9W5D  | 7   | 61  | 23W4D | 10  | >105 | n/a   | —   |
| 17  | 10W5D | 0   | 62   | 23W4D | 5   | 17  | 10W0D | 7   | 62  | 24W0D | 10  |      |       |     |
| 18  | 10W6D | 0   | 63   | 23W6D | 5   | 18  | 10W2D | 7   | 63  | 24W2D | 10  |      |       |     |
| 19  | 11W1D | 0   | 64   | 24W1D | 6   | 19  | 10W4D | 7   | 64  | 24W4D | 10  |      |       |     |
| 20  | 11W3D | 0   | 65   | 24W4D | 6   | 20  | 10W6D | 7   | 65  | 24W6D | 10  |      |       |     |
| 21  | 11W5D | 0   | 66   | 24W6D | 6   | 21  | 11W1D | 7   | 66  | 25W1D | 11  |      |       |     |
| 22  | 12W0D | 0   | 67   | 25W1D | 6   | 22  | 11W3D | 7   | 67  | 25W3D | 12  |      |       |     |
| 23  | 12W2D | 0   | 68   | 25W3D | 6   | 23  | 11W5D | 7   | 68  | 25W6D | 10  |      |       |     |
| 24  | 12W4D | 5   | 69   | 25W5D | 6   | 24  | 12W0D | 7   | 69  | 26W1D | 10  |      |       |     |
| 25  | 12W6D | 5   | 70   | 26W1D | 6   | 25  | 12W2D | 7   | 70  | 26W3D | 10  |      |       |     |
| 26  | 13W1D | 5   | 71   | 26W3D | 6   | 26  | 12W4D | 7   | 71  | 26W5D | 12  |      |       |     |
| 27  | 13W2D | 5   | 72   | 26W6D | 6   | 27  | 12W6D | 7   | 72  | 27W1D | 11  |      |       |     |
| 28  | 13W4D | 4   | 73   | 27W1D | 6   | 28  | 13W1D | 7   | 73  | 27W3D | 13  |      |       |     |
| 29  | 13W6D | 4   | 74   | 27W3D | 6   | 29  | 13W3D | 8   | 74  | 27W6D | 12  |      |       |     |
| 30  | 14W1D | 4   | 75   | 27W6D | 6   | 30  | 13W5D | 7   | 75  | 28W1D | 12  |      |       |     |
| 31  | 14W3D | 4   | 76   | 28W1D | 6   | 31  | 14W0D | 8   | 76  | 28W4D | 13  |      |       |     |
| 32  | 14W4D | 4   | 77   | 28W4D | 6   | 32  | 14W2D | 8   | 77  | 28W6D | 13  |      |       |     |
| 33  | 14W6D | 4   | 78   | 28W6D | 6   | 33  | 14W4D | 9   | 78  | 29W2D | 15  |      |       |     |
| 34  | 15W2D | 4   | 79   | 29W2D | 6   | 34  | 15W0D | 9   | 79  | 29W5D | 16  |      |       |     |
| 35  | 15W4D | 4   | 80   | 29W5D | 6   | 35  | 15W2D | 8   | 80  | 30W0D | 15  |      |       |     |
| 36  | 15W6D | 4   | 81   | 30W0D | 6   | 36  | 15W4D | 9   | 81  | 30W3D | 15  |      |       |     |
| 37  | 16W1D | 4   | 82   | 30W3D | 6   | 37  | 16W0D | 8   | 82  | 31W0D | 15  |      |       |     |
| 38  | 16W3D | 4   | 83   | 30W5D | 6   | 38  | 16W2D | 9   | 83  | 31W2D | 16  |      |       |     |
| 39  | 16W5D | 4   | 84   | 31W2D | 6   | 39  | 16W4D | 9   | 84  | 31W6D | 17  |      |       |     |
| 40  | 17W0D | 4   | 85   | 31W5D | 6   | 40  | 17W0D | 9   | 85  | 32W2D | 17  |      |       |     |
| 41  | 17W2D | 4   | 86   | 32W1D | 6   | 41  | 17W2D | 9   | 86  | 32W5D | 18  |      |       |     |
| 42  | 17W4D | 4   | 87   | 32W4D | 6   | 42  | 17W4D | 9   | 87  | 33W2D | 20  |      |       |     |
| 43  | 17W6D | 4   | 88   | 33W0D | 7   | 43  | 17W6D | 9   | 88  | 33W5D | 19  |      |       |     |
| 44  | 18W1D | 4   | 89   | 33W3D | 7   | 44  | 18W1D | 9   | 89  | 34W2D | 19  |      |       |     |
| 45  | 18W3D | 4   | 90   | 33W6D | 7   | 45  | 18W4D | 9   | 90  | 34W5D | 19  |      |       |     |
| 46  | 18W5D | 4   | 91   | 34W3D | 7   | 46  | 18W6D | 9   | 91  | 35W1D | 25  |      |       |     |
| 47  | 19W0D | 4   | 92   | 34W6D | 7   | 47  | 19W1D | 10  | 92  | 35W6D | 24  |      |       |     |
| 48  | 19W2D | 5   | 93   | 35W3D | 7   | 48  | 19W3D | 10  | 93  | 36W5D | 21  |      |       |     |
| 49  | 19W4D | 5   | 94   | 36W0D | 7   | 49  | 19W5D | 10  | 94  | 37W3D | 19  |      |       |     |
| 50  | 19W6D | 5   | 95   | 36W3D | 7   | 50  | 20W0D | 10  | 95  | 38W3D | 22  |      |       |     |
| 51  | 20W1D | 5   | 96   | 37W1D | 7   | 51  | 20W3D | 10  | 96  | 38W6D | 25  |      |       |     |
| 52  | 20W3D | 5   | 97   | 37W6D | 7   | 52  | 20W5D | 10  | 97  | 39W0D | 22  |      |       |     |
| 53  | 20W6D | 5   | 98   | 38W4D | 7   | 53  | 21W0D | 11  | 98  | 39W2D | 20  |      |       |     |
| 54  | 21W1D | 5   | 99   | 39W3D | 7   | 54  | 21W3D | 10  | 99  | 39W3D | 22  |      |       |     |
| 55  | 21W2D | 5   | 100  | 40W3D | 7   | 55  | 21W5D | 10  | 100 | 39W4D | 20  |      |       |     |
| 56  | 21W4D | 5   | 101  | 41W3D | 7   | 56  | 22W0D | 9   | 101 | 39W5D | 20  |      |       |     |
| 57  | 21W6D | 5   | >101 | n/a   | —   | 57  | 22W2D | 9   | 102 | 39W6D | 19  |      |       |     |

Table 18–13. BPD : Hansmann—Known LMP (left)—Unknown LMP (right)

## OB Tables

CRL : Hansmann

Unit : BPD (mm)

2SD = mm

Age (Weeks/Days)

2SD (mm or day)

2SD = day

| CRL | Age   | 2SD | CRL | Age   | 2SD | CRL | Age   | 2SD | CRL  | Age   | 2SD | CRL | Age   | 2SD | CRL  | Age   | 2SD |
|-----|-------|-----|-----|-------|-----|-----|-------|-----|------|-------|-----|-----|-------|-----|------|-------|-----|
| <13 | n/a   | —   | 57  | 12W2D | 16  | 102 | 15W6D | 10  | 147  | 20W2D | 17  | <6  | n/a   | —   | 86   | 14W4D | 12  |
| 13  | 7W4D  | 0   | 58  | 12W2D | 16  | 103 | 15W6D | 10  | 148  | 20W2D | 17  | 6   | 6W1D  | 7   | 90   | 14W6D | 12  |
| 14  | 7W5D  | 0   | 59  | 12W3D | 16  | 104 | 16W0D | 10  | 149  | 20W3D | 17  | 7   | 6W2D  | 7   | 93   | 15W1D | 12  |
| 15  | 8W0D  | 0   | 60  | 12W3D | 16  | 105 | 16W1D | 10  | 150  | 20W4D | 17  | 8   | 6W4D  | 7   | 96   | 15W3D | 12  |
| 16  | 8W1D  | 0   | 61  | 12W4D | 15  | 106 | 16W2D | 10  | 151  | 20W4D | 0   | 9   | 6W6D  | 7   | 100  | 15W5D | 12  |
| 17  | 8W2D  | 0   | 62  | 12W4D | 15  | 107 | 16W2D | 10  | 152  | 20W5D | 0   | 10  | 7W0D  | 7   | 103  | 16W0D | 12  |
| 18  | 8W3D  | 0   | 63  | 12W5D | 15  | 108 | 16W3D | 10  | 153  | 20W5D | 0   | 11  | 7W2D  | 7   | 106  | 16W2D | 13  |
| 19  | 8W4D  | 7   | 64  | 12W5D | 15  | 109 | 16W3D | 10  | 154  | 20W6D | 0   | 12  | 7W3D  | 7   | 110  | 16W4D | 14  |
| 20  | 8W5D  | 7   | 65  | 12W6D | 15  | 110 | 16W4D | 10  | 155  | 21W0D | 0   | 13  | 7W4D  | 7   | 113  | 17W0D | 14  |
| 21  | 8W6D  | 8   | 66  | 12W6D | 15  | 111 | 16W4D | 11  | 156  | 21W0D | 0   | 14  | 7W6D  | 7   | 116  | 17W2D | 14  |
| 22  | 9W0D  | 8   | 67  | 13W0D | 15  | 112 | 16W5D | 11  | 157  | 21W1D | 0   | 15  | 8W0D  | 7   | 120  | 17W4D | 14  |
| 23  | 9W1D  | 10  | 68  | 13W1D | 15  | 113 | 16W5D | 11  | 158  | 21W1D | 0   | 16  | 8W2D  | 7   | 123  | 18W0D | 14  |
| 24  | 9W2D  | 10  | 69  | 13W1D | 15  | 114 | 16W6D | 11  | 159  | 21W2D | 0   | 17  | 8W3D  | 7   | 126  | 18W2D | 14  |
| 25  | 9W3D  | 11  | 70  | 13W2D | 15  | 115 | 17W0D | 11  | 160  | 21W3D | 0   | 18  | 8W4D  | 7   | 130  | 18W6D | 14  |
| 26  | 9W4D  | 11  | 71  | 13W3D | 15  | 116 | 17W1D | 12  | 161  | 21W3D | 0   | 19  | 8W5D  | 7   | 133  | 19W1D | 15  |
| 27  | 9W4D  | 11  | 72  | 13W3D | 15  | 117 | 17W2D | 12  | 162  | 21W4D | 0   | 20  | 8W6D  | 7   | 136  | 19W4D | 16  |
| 28  | 9W5D  | 11  | 73  | 13W4D | 15  | 118 | 17W2D | 12  | 163  | 21W4D | 0   | 21  | 9W0D  | 7   | 140  | 20W0D | 16  |
| 29  | 9W6D  | 11  | 74  | 13W4D | 15  | 119 | 17W3D | 12  | 164  | 21W5D | 0   | 22  | 9W1D  | 7   | 143  | 20W3D | 16  |
| 30  | 10W0D | 12  | 75  | 13W5D | 15  | 120 | 17W3D | 12  | 165  | 21W6D | 0   | 23  | 9W2D  | 7   | 146  | 20W6D | 16  |
| 31  | 10W0D | 12  | 76  | 13W5D | 15  | 121 | 17W4D | 13  | 166  | 21W6D | 0   | 24  | 9W3D  | 7   | 150  | 21W3D | 15  |
| 32  | 10W1D | 12  | 77  | 13W6D | 15  | 122 | 17W5D | 13  | 167  | 22W0D | 0   | 26  | 9W5D  | 7   | >150 | n/a   | —   |
| 33  | 10W2D | 12  | 78  | 13W6D | 15  | 123 | 17W5D | 13  | 168  | 22W0D | 0   | 28  | 10W0D | 8   |      |       |     |
| 34  | 10W3D | 12  | 79  | 14W0D | 15  | 124 | 17W6D | 13  | 169  | 22W1D | 0   | 30  | 10W2D | 8   |      |       |     |
| 35  | 10W3D | 13  | 80  | 14W0D | 15  | 125 | 18W0D | 13  | 170  | 22W1D | 0   | 32  | 10W3D | 8   |      |       |     |
| 36  | 10W4D | 13  | 81  | 14W1D | 13  | 126 | 18W1D | 14  | 171  | 22W2D | 0   | 34  | 10W5D | 8   |      |       |     |
| 37  | 10W5D | 13  | 82  | 14W1D | 13  | 127 | 18W1D | 14  | 172  | 22W2D | 0   | 36  | 10W6D | 8   |      |       |     |
| 38  | 10W5D | 13  | 83  | 14W2D | 13  | 128 | 18W2D | 14  | 173  | 22W3D | 0   | 38  | 11W1D | 8   |      |       |     |
| 39  | 10W6D | 13  | 84  | 14W2D | 13  | 129 | 18W2D | 14  | 174  | 22W3D | 0   | 40  | 11W2D | 8   |      |       |     |
| 40  | 10W6D | 13  | 85  | 14W3D | 13  | 130 | 18W3D | 15  | 175  | 22W4D | 0   | 42  | 11W3D | 8   |      |       |     |
| 41  | 11W0D | 14  | 86  | 14W3D | 13  | 131 | 18W4D | 15  | >175 | n/a   | —   | 44  | 11W4D | 9   |      |       |     |
| 42  | 11W1D | 14  | 87  | 14W4D | 13  | 132 | 18W4D | 15  |      |       |     | 46  | 11W6D | 9   |      |       |     |
| 43  | 11W1D | 14  | 88  | 14W4D | 13  | 133 | 18W5D | 15  |      |       |     | 48  | 12W0D | 9   |      |       |     |
| 44  | 11W2D | 14  | 89  | 14W5D | 13  | 134 | 18W6D | 15  |      |       |     | 50  | 12W1D | 9   |      |       |     |
| 45  | 11W2D | 14  | 90  | 14W6D | 13  | 135 | 19W0D | 15  |      |       |     | 52  | 12W2D | 9   |      |       |     |
| 46  | 11W3D | 14  | 91  | 14W6D | 12  | 136 | 19W1D | 15  |      |       |     | 54  | 12W3D | 9   |      |       |     |
| 47  | 11W3D | 15  | 92  | 15W0D | 12  | 137 | 19W1D | 15  |      |       |     | 56  | 12W4D | 9   |      |       |     |
| 48  | 11W4D | 15  | 93  | 15W0D | 12  | 138 | 19W2D | 15  |      |       |     | 58  | 12W5D | 9   |      |       |     |
| 49  | 11W4D | 15  | 94  | 15W1D | 12  | 139 | 19W3D | 15  |      |       |     | 60  | 12W6D | 9   |      |       |     |
| 50  | 11W5D | 15  | 95  | 15W2D | 12  | 140 | 19W4D | 15  |      |       |     | 63  | 13W0D | 10  |      |       |     |
| 51  | 11W5D | 15  | 96  | 15W3D | 11  | 141 | 19W4D | 16  |      |       |     | 66  | 13W2D | 10  |      |       |     |
| 52  | 11W6D | 15  | 97  | 15W3D | 11  | 142 | 19W5D | 16  |      |       |     | 70  | 13W3D | 10  |      |       |     |
| 53  | 11W6D | 15  | 98  | 15W4D | 11  | 143 | 19W5D | 16  |      |       |     | 73  | 13W5D | 11  |      |       |     |
| 54  | 12W0D | 15  | 99  | 15W4D | 11  | 144 | 19W6D | 16  |      |       |     | 76  | 13W6D | 11  |      |       |     |
| 55  | 12W1D | 16  | 100 | 15W5D | 11  | 145 | 20W0D | 16  |      |       |     | 80  | 14W1D | 11  |      |       |     |
| 56  | 12W1D | 16  | 101 | 15W5D | 10  | 146 | 20W1D | 17  |      |       |     | 83  | 14W2D | 12  |      |       |     |

Table 18–14. CRL : Hansmann—Known LMP (left)—Unknown LMP (right)



## OB Tables

FL : Hansmann  
Unit : FL (mm)  
Age (Weeks/Days)  
2SD (mm)

GS : Hansmann      Unit : GS (mm)  
Hansmann : M      Age (Day)  
and AI : Geburtsh, u,      SD (mm)  
Frauenheilk 39 :  
656,1979

| FL  | Age   | 2SD | FL  | Age   | 2SD |
|-----|-------|-----|-----|-------|-----|
| <12 | n/a   | —   | 52  | 27W5D | 5   |
| 12  | 13W4D | 0   | 53  | 28W1D | 5   |
| 13  | 13W6D | 0   | 54  | 28W4D | 5   |
| 14  | 14W1D | 0   | 55  | 29W0D | 6   |
| 15  | 14W3D | 0   | 56  | 29W3D | 6   |
| 16  | 14W5D | 5   | 57  | 29W6D | 6   |
| 17  | 15W1D | 5   | 58  | 30W2D | 6   |
| 18  | 15W2D | 4   | 59  | 30W5D | 5   |
| 19  | 15W5D | 4   | 60  | 31W2D | 5   |
| 20  | 16W0D | 4   | 61  | 31W5D | 5   |
| 21  | 16W2D | 4   | 62  | 32W1D | 5   |
| 22  | 16W4D | 4   | 63  | 32W5D | 5   |
| 23  | 16W6D | 4   | 64  | 33W1D | 6   |
| 24  | 17W2D | 4   | 65  | 33W5D | 6   |
| 25  | 17W4D | 4   | 66  | 34W1D | 6   |
| 26  | 17W6D | 4   | 67  | 34W5D | 6   |
| 27  | 18W2D | 4   | 68  | 35W1D | 6   |
| 28  | 18W4D | 4   | 69  | 35W5D | 6   |
| 29  | 18W6D | 4   | 70  | 36W1D | 6   |
| 30  | 19W2D | 4   | 71  | 36W5D | 6   |
| 31  | 19W4D | 4   | 72  | 37W2D | 6   |
| 32  | 20W0D | 4   | 73  | 37W6D | 6   |
| 33  | 20W3D | 4   | 74  | 38W3D | 7   |
| 34  | 20W5D | 4   | 75  | 39W0D | 7   |
| 35  | 21W1D | 5   | >75 | n/a   | —   |
| 36  | 21W3D | 5   |     |       |     |
| 37  | 21W6D | 5   |     |       |     |
| 38  | 22W1D | 5   |     |       |     |
| 39  | 22W4D | 5   |     |       |     |
| 40  | 22W6D | 5   |     |       |     |
| 41  | 23W2D | 5   |     |       |     |
| 42  | 23W5D | 5   |     |       |     |
| 43  | 24W0D | 5   |     |       |     |
| 44  | 24W3D | 5   |     |       |     |
| 45  | 24W6D | 5   |     |       |     |
| 46  | 25W2D | 5   |     |       |     |
| 47  | 25W4D | 5   |     |       |     |
| 48  | 26W0D | 5   |     |       |     |
| 49  | 26W3D | 5   |     |       |     |
| 50  | 26W6D | 5   |     |       |     |
| 51  | 27W3D | 5   |     |       |     |

Table 18–15. FL : Hansmann  
Known/Unknown LMP

| GS  | Age | SD | GS  | Age | SD |
|-----|-----|----|-----|-----|----|
| <10 | n/a | —  | 38  | 60  | 5  |
| 10  | 33  | 5  | 39  | 61  | 5  |
| 11  | 34  | 5  | 40  | 62  | 5  |
| 12  | 35  | 5  | 41  | 63  | 5  |
| 13  | 36  | 5  | 42  | 64  | 5  |
| 14  | 37  | 5  | 43  | 65  | 5  |
| 15  | 38  | 5  | 44  | 66  | 5  |
| 16  | 39  | 5  | 45  | 67  | 5  |
| 17  | 40  | 5  | 46  | 68  | 5  |
| 18  | 41  | 5  | 47  | 69  | 5  |
| 19  | 42  | 5  | 48  | 70  | 5  |
| 20  | 43  | 5  | 49  | 71  | 5  |
| 21  | 44  | 5  | 50  | 72  | 5  |
| 22  | 45  | 5  | 51  | 73  | 5  |
| 23  | 46  | 5  | 52  | 74  | 5  |
| 24  | 47  | 5  | 53  | 75  | 5  |
| 25  | 48  | 5  | 54  | 76  | 5  |
| 26  | 49  | 5  | 55  | 77  | 5  |
| 27  | 50  | 5  | 56  | 78  | 5  |
| 28  | 51  | 5  | 57  | 79  | 5  |
| 29  | 52  | 5  | 58  | 80  | 5  |
| 30  | 53  | 5  | 59  | 81  | 5  |
| 31  | 54  | 5  | 60  | 82  | 5  |
| 32  | 55  | 5  | 61  | 83  | 5  |
| 33  | 56  | 5  | 62  | 84  | 5  |
| 34  | 57  | 5  | 63  | 85  | 5  |
| 35  | 58  | 5  | 64  | 86  | 5  |
| 36  | 58  | 5  | 65  | 87  | 5  |
| 37  | 59  | 5  | >65 | n/a | —  |

Table 18–16. GS : Hansmann



## OB Tables

HC : Hansmann  
Unit : HC (mm)  
Age (Weeks/Days)  
2SD (mm)

| HC   | Age   | 2SD | HC   | Age   | 2SD |
|------|-------|-----|------|-------|-----|
| <105 | n/a   | —   | 305  | 30W5D | 19  |
| 105  | 14W0D | 0   | 310  | 31W2D | 19  |
| 110  | 14W3D | 0   | 315  | 32W1D | 19  |
| 115  | 14W6D | 14  | 320  | 32W5D | 20  |
| 120  | 15W3D | 14  | 325  | 33W3D | 20  |
| 125  | 15W5D | 14  | 330  | 34W2D | 20  |
| 130  | 16W1D | 14  | 335  | 35W1D | 20  |
| 135  | 16W4D | 14  | 340  | 36W2D | 20  |
| 140  | 17W0D | 14  | 345  | 37W6D | 20  |
| 145  | 17W3D | 15  | >345 | n/a   | —   |
| 150  | 17W6D | 15  |      |       |     |
| 155  | 18W1D | 16  |      |       |     |
| 160  | 18W4D | 16  |      |       |     |
| 165  | 19W0D | 16  |      |       |     |
| 170  | 19W0D | 16  |      |       |     |
| 175  | 19W3D | 16  |      |       |     |
| 180  | 19W5D | 16  |      |       |     |
| 185  | 20W1D | 17  |      |       |     |
| 190  | 20W4D | 17  |      |       |     |
| 195  | 21W0D | 17  |      |       |     |
| 200  | 21W2D | 17  |      |       |     |
| 205  | 21W5D | 17  |      |       |     |
| 210  | 22W1D | 17  |      |       |     |
| 215  | 22W4D | 17  |      |       |     |
| 220  | 23W0D | 17  |      |       |     |
| 225  | 23W3D | 17  |      |       |     |
| 230  | 23W5D | 18  |      |       |     |
| 235  | 24W1D | 18  |      |       |     |
| 240  | 24W4D | 18  |      |       |     |
| 245  | 25W0D | 18  |      |       |     |
| 250  | 25W3D | 18  |      |       |     |
| 255  | 25W6D | 18  |      |       |     |
| 260  | 26W2D | 18  |      |       |     |
| 265  | 26W5D | 18  |      |       |     |
| 270  | 27W1D | 18  |      |       |     |
| 275  | 27W4D | 19  |      |       |     |
| 280  | 28W1D | 19  |      |       |     |
| 285  | 28W5D | 19  |      |       |     |
| 290  | 29W1D | 19  |      |       |     |
| 295  | 29W5D | 19  |      |       |     |
| 300  | 30W2D | 19  |      |       |     |

Table 18–17. HC : Hansmann  
Known/Unknown LMP



## OB Tables

OFD : Hansmann

Unit : OFD (mm)

Age : (Weeks/Days)

2SD (mm)

| OFD | Age   | 2SD | OFD | Age   | 2SD | OFD  | Age   | 2SD |
|-----|-------|-----|-----|-------|-----|------|-------|-----|
| <34 | n/a   | —   | 74  | 23W1D | 7   | 115  | 39W0D | 8   |
| 34  | 13W5D | 0   | 75  | 23W2D | 7   | >115 | n/a   | —   |
| 35  | 14W0D | 0   | 76  | 23W4D | 7   |      |       |     |
| 36  | 14W2D | 0   | 77  | 23W6D | 7   |      |       |     |
| 37  | 14W4D | 5   | 78  | 24W1D | 7   |      |       |     |
| 38  | 14W4D | 5   | 79  | 24W2D | 7   |      |       |     |
| 39  | 14W6D | 5   | 80  | 24W4D | 7   |      |       |     |
| 40  | 15W1D | 5   | 81  | 24W6D | 7   |      |       |     |
| 41  | 15W3D | 5   | 82  | 25W1D | 7   |      |       |     |
| 42  | 15W5D | 5   | 83  | 25W2D | 7   |      |       |     |
| 43  | 16W0D | 5   | 84  | 25W4D | 7   |      |       |     |
| 44  | 16W1D | 5   | 85  | 25W6D | 7   |      |       |     |
| 45  | 16W3D | 5   | 86  | 26W1D | 7   |      |       |     |
| 46  | 16W4D | 5   | 87  | 26W3D | 7   |      |       |     |
| 47  | 16W6D | 5   | 88  | 26W5D | 7   |      |       |     |
| 48  | 17W1D | 5   | 89  | 27W0D | 7   |      |       |     |
| 49  | 17W3D | 5   | 90  | 27W2D | 7   |      |       |     |
| 50  | 17W4D | 5   | 91  | 27W4D | 8   |      |       |     |
| 51  | 17W6D | 5   | 92  | 27W6D | 8   |      |       |     |
| 52  | 18W1D | 5   | 93  | 28W1D | 8   |      |       |     |
| 53  | 18W2D | 5   | 94  | 28W3D | 8   |      |       |     |
| 54  | 18W4D | 5   | 95  | 28W5D | 8   |      |       |     |
| 55  | 18W6D | 5   | 96  | 29W0D | 8   |      |       |     |
| 56  | 19W0D | 6   | 97  | 29W3D | 8   |      |       |     |
| 57  | 19W2D | 6   | 98  | 29W5D | 8   |      |       |     |
| 58  | 19W3D | 6   | 99  | 30W0D | 8   |      |       |     |
| 59  | 19W5D | 6   | 100 | 30W3D | 8   |      |       |     |
| 60  | 20W0D | 6   | 101 | 30W5D | 8   |      |       |     |
| 61  | 20W1D | 6   | 102 | 31W1D | 8   |      |       |     |
| 62  | 20W2D | 6   | 103 | 31W4D | 8   |      |       |     |
| 63  | 20W4D | 6   | 104 | 32W0D | 8   |      |       |     |
| 64  | 20W6D | 6   | 105 | 32W3D | 8   |      |       |     |
| 65  | 21W0D | 6   | 106 | 32W6D | 8   |      |       |     |
| 66  | 21W2D | 6   | 107 | 33W3D | 8   |      |       |     |
| 67  | 21W4D | 6   | 108 | 33W6D | 8   |      |       |     |
| 68  | 21W5D | 6   | 109 | 34W3D | 8   |      |       |     |
| 69  | 22W0D | 6   | 110 | 35W0D | 8   |      |       |     |
| 70  | 22W1D | 7   | 111 | 35W4D | 8   |      |       |     |
| 71  | 22W3D | 7   | 112 | 36W2D | 8   |      |       |     |
| 72  | 22W4D | 7   | 113 | 37W0D | 8   |      |       |     |
| 73  | 22W6D | 7   | 114 | 38W0D | 8   |      |       |     |

Table 18–18. OFD : Hansmann  
Known/Unknown LMP



## OB Tables

TAD : Hansmann  
Hansmann : M  
and AI : Geburtsh, u,  
Frauenheilk

Unit : TAD (mm)  
Age (Day)  
SD (mm)

ThD : Hansmann  
Unit : THQ (mm)  
Age (Weeks/Days)  
2SD (mm)

| TAD | Age | SD | TAD | Age | SD | TAD  | Age | SD |
|-----|-----|----|-----|-----|----|------|-----|----|
| <20 | n/a | —  | 49  | 148 | 4  | 79   | 217 | 5  |
| 20  | 87  | 4  | 50  | 150 | 4  | 80   | 220 | 5  |
| 21  | 89  | 4  | 51  | 152 | 4  | 81   | 222 | 5  |
| 22  | 91  | 4  | 52  | 155 | 4  | 82   | 225 | 5  |
| 23  | 93  | 4  | 53  | 157 | 4  | 83   | 227 | 5  |
| 24  | 95  | 4  | 54  | 159 | 4  | 84   | 230 | 5  |
| 25  | 97  | 4  | 55  | 161 | 4  | 85   | 232 | 5  |
| 26  | 99  | 4  | 56  | 164 | 4  | 86   | 235 | 5  |
| 27  | 101 | 4  | 57  | 166 | 4  | 87   | 237 | 5  |
| 28  | 103 | 4  | 58  | 168 | 4  | 88   | 240 | 5  |
| 29  | 105 | 4  | 59  | 170 | 4  | 89   | 242 | 5  |
| 30  | 107 | 4  | 60  | 173 | 4  | 90   | 245 | 5  |
| 31  | 109 | 4  | 61  | 175 | 4  | 91   | 247 | 5  |
| 32  | 111 | 4  | 62  | 177 | 4  | 92   | 250 | 5  |
| 33  | 113 | 4  | 63  | 179 | 4  | 93   | 252 | 5  |
| 34  | 115 | 4  | 64  | 182 | 4  | 94   | 255 | 5  |
| 35  | 117 | 4  | 65  | 184 | 4  | 95   | 258 | 5  |
| 36  | 119 | 4  | 66  | 186 | 4  | 96   | 261 | 5  |
| 37  | 122 | 4  | 67  | 188 | 4  | 97   | 264 | 5  |
| 38  | 124 | 4  | 68  | 191 | 5  | 98   | 267 | 5  |
| 39  | 126 | 4  | 69  | 193 | 5  | 99   | 270 | 5  |
| 40  | 128 | 4  | 70  | 195 | 5  | 100  | 273 | 5  |
| 41  | 130 | 4  | 71  | 198 | 5  | 101  | 276 | 5  |
| 42  | 132 | 4  | 72  | 200 | 5  | 102  | 279 | 5  |
| 43  | 135 | 4  | 73  | 203 | 5  | 103  | 282 | 5  |
| 44  | 137 | 4  | 74  | 205 | 5  | >103 | n/a | —  |
| 45  | 139 | 4  | 75  | 208 | 5  |      |     |    |
| 46  | 141 | 4  | 76  | 210 | 5  |      |     |    |
| 47  | 143 | 4  | 77  | 212 | 5  |      |     |    |
| 48  | 146 | 4  | 78  | 215 | 5  |      |     |    |

Table 18–19. TAD : Hansmann

| ThD | Age   | zsd | ThD | Age   | zsd | ThD  | Age   | zsd |
|-----|-------|-----|-----|-------|-----|------|-------|-----|
| <20 | n/a   | —   | 60  | 24W6D | 7   | 101  | 39W3D | 12  |
| 20  | 12W4D | 0   | 61  | 25W1D | 7   | 102  | 39W6D | 14  |
| 21  | 12W6D | 0   | 62  | 25W3D | 7   | 103  | 40W2D | 14  |
| 22  | 13W1D | 0   | 63  | 25W5D | 7   | 104  | 40W5D | 14  |
| 23  | 13W3D | 4   | 64  | 26W1D | 7   | 105  | 41W2D | 14  |
| 24  | 13W4D | 4   | 65  | 26W3D | 7   | >105 | n/a   | —   |
| 25  | 13W6D | 4   | 66  | 26W5D | 7   |      |       |     |
| 26  | 14W1D | 4   | 67  | 27W0D | 7   |      |       |     |
| 27  | 14W3D | 4   | 68  | 27W3D | 8   |      |       |     |
| 28  | 14W6D | 4   | 69  | 27W5D | 8   |      |       |     |
| 29  | 15W1D | 4   | 70  | 28W0D | 8   |      |       |     |
| 30  | 15W2D | 4   | 71  | 28W3D | 8   |      |       |     |
| 31  | 15W4D | 4   | 72  | 28W5D | 8   |      |       |     |
| 32  | 15W6D | 4   | 73  | 29W1D | 8   |      |       |     |
| 33  | 16W2D | 4   | 74  | 29W3D | 8   |      |       |     |
| 34  | 16W4D | 4   | 75  | 29W5D | 8   |      |       |     |
| 35  | 16W5D | 4   | 76  | 30W1D | 8   |      |       |     |
| 36  | 17W1D | 5   | 77  | 30W3D | 8   |      |       |     |
| 37  | 17W3D | 5   | 78  | 30W5D | 8   |      |       |     |
| 38  | 17W5D | 5   | 79  | 31W1D | 8   |      |       |     |
| 39  | 18W1D | 5   | 80  | 31W3D | 8   |      |       |     |
| 40  | 18W3D | 5   | 81  | 31W5D | 8   |      |       |     |
| 41  | 18W5D | 5   | 82  | 32W1D | 9   |      |       |     |
| 42  | 19W0D | 5   | 83  | 32W4D | 9   |      |       |     |
| 43  | 19W3D | 5   | 84  | 32W6D | 9   |      |       |     |
| 44  | 19W5D | 5   | 85  | 33W1D | 9   |      |       |     |
| 45  | 19W6D | 5   | 86  | 33W4D | 9   |      |       |     |
| 46  | 20W2D | 5   | 87  | 33W6D | 9   |      |       |     |
| 47  | 20W4D | 6   | 88  | 34W2D | 9   |      |       |     |
| 48  | 20W6D | 6   | 89  | 34W4D | 9   |      |       |     |
| 49  | 21W1D | 6   | 90  | 35W0D | 9   |      |       |     |
| 50  | 21W4D | 6   | 91  | 35W3D | 10  |      |       |     |
| 51  | 21W6D | 6   | 92  | 35W5D | 10  |      |       |     |
| 52  | 22W1D | 6   | 93  | 36W1D | 10  |      |       |     |
| 53  | 22W4D | 6   | 94  | 36W3D | 10  |      |       |     |
| 54  | 22W6D | 6   | 95  | 36W6D | 10  |      |       |     |
| 55  | 23W1D | 6   | 96  | 37W1D | 10  |      |       |     |
| 56  | 23W3D | 6   | 97  | 37W4D | 10  |      |       |     |
| 57  | 23W6D | 7   | 98  | 38W1D | 11  |      |       |     |
| 58  | 24W1D | 7   | 99  | 38W4D | 11  |      |       |     |
| 59  | 24W3D | 7   | 100 | 38W6D | 11  |      |       |     |

Table 18–20. ThD : Hansmann  
Known/Unknown LMP



## OB Tables

GS : Hellman      Unit :      GS (mm)      AC : Jeanty      Unit :      AC (mm)  
 Age (Week)      Jeanty : Radiology      Age (Day)  
 (SD: Standard Deviation )      SD (Week)      143 : 513,1982      SD (mm)

| GS  | Age | SD    | GS  | Age  | SD    |
|-----|-----|-------|-----|------|-------|
| <10 | n/a | ——    | 36  | 8.8  | ± 1.0 |
| 10  | 5.0 | ± 1.0 | 37  | 8.9  | ± 1.0 |
| 11  | 5.2 | ± 1.0 | 38  | 9.0  | ± 1.0 |
| 12  | 5.3 | ± 1.0 | 39  | 9.2  | ± 1.0 |
| 13  | 5.5 | ± 1.0 | 40  | 9.3  | ± 1.0 |
| 14  | 5.6 | ± 1.0 | 41  | 9.5  | ± 1.0 |
| 15  | 5.8 | ± 1.0 | 42  | 9.6  | ± 1.0 |
| 16  | 5.9 | ± 1.0 | 43  | 9.7  | ± 1.0 |
| 17  | 6.0 | ± 1.0 | 44  | 9.9  | ± 1.0 |
| 18  | 6.2 | ± 1.0 | 45  | 10.0 | ± 1.0 |
| 19  | 6.3 | ± 1.0 | 46  | 10.2 | ± 1.0 |
| 20  | 6.5 | ± 1.0 | 47  | 10.3 | ± 1.0 |
| 21  | 6.6 | ± 1.0 | 48  | 10.5 | ± 1.0 |
| 22  | 6.8 | ± 1.0 | 49  | 10.6 | ± 1.0 |
| 23  | 6.9 | ± 1.0 | 50  | 10.7 | ± 1.0 |
| 24  | 7.0 | ± 1.0 | 51  | 10.9 | ± 1.0 |
| 25  | 7.2 | ± 1.0 | 52  | 11.0 | ± 1.0 |
| 26  | 7.3 | ± 1.0 | 53  | 11.2 | ± 1.0 |
| 27  | 7.5 | ± 1.0 | 54  | 11.3 | ± 1.0 |
| 28  | 7.6 | ± 1.0 | 55  | 11.5 | ± 1.0 |
| 29  | 7.8 | ± 1.0 | 56  | 11.6 | ± 1.0 |
| 30  | 7.9 | ± 1.0 | 57  | 11.7 | ± 1.0 |
| 31  | 8.0 | ± 1.0 | 58  | 11.9 | ± 1.0 |
| 32  | 8.2 | ± 1.0 | 59  | 12.0 | ± 1.0 |
| 33  | 8.3 | ± 1.0 | 60  | 12.2 | ± 1.0 |
| 34  | 8.5 | ± 1.0 | >60 | n/a  | ——    |
| 35  | 8.6 | ± 1.0 |     |      |       |

Table 18–21. GS : Hellman

| AC  | Age | SD | AC   | Age | SD |
|-----|-----|----|------|-----|----|
| <50 | n/a | —— | 185  | 169 | 22 |
| 50  | 79  | 22 | 190  | 172 | 22 |
| 55  | 82  | 22 | 195  | 176 | 22 |
| 60  | 85  | 22 | 200  | 179 | 22 |
| 65  | 89  | 22 | 205  | 182 | 22 |
| 70  | 92  | 22 | 210  | 186 | 22 |
| 75  | 95  | 22 | 215  | 189 | 22 |
| 80  | 99  | 22 | 220  | 192 | 22 |
| 85  | 102 | 22 | 225  | 196 | 22 |
| 90  | 105 | 22 | 230  | 199 | 22 |
| 95  | 109 | 22 | 235  | 203 | 22 |
| 100 | 112 | 22 | 240  | 206 | 22 |
| 105 | 115 | 22 | 245  | 210 | 22 |
| 110 | 119 | 22 | 250  | 214 | 22 |
| 115 | 122 | 22 | 255  | 218 | 22 |
| 120 | 125 | 22 | 260  | 222 | 22 |
| 125 | 129 | 22 | 265  | 226 | 22 |
| 130 | 132 | 22 | 270  | 230 | 22 |
| 135 | 135 | 22 | 275  | 234 | 22 |
| 140 | 139 | 22 | 280  | 239 | 22 |
| 145 | 142 | 22 | 285  | 244 | 22 |
| 150 | 146 | 22 | 290  | 249 | 22 |
| 155 | 149 | 22 | 295  | 254 | 22 |
| 160 | 152 | 22 | 300  | 259 | 22 |
| 165 | 156 | 22 | 305  | 265 | 22 |
| 170 | 159 | 22 | 310  | 272 | 22 |
| 175 | 162 | 22 | 315  | 279 | 22 |
| 180 | 166 | 22 | >315 | n/a | —— |

Table 18–22. AC : Jeanty

OB Tables

BD : Jeanty  
Jeanty : Radiology  
143 : 513,1982

Unit : BD (mm)  
Age (Day)  
SD (mm)

BPD : Jeanty  
Jeanty : Radiology  
143 : 513,1982

Unit : BPD (mm)  
Age (Day)  
SD (mm)

| BD  | Age | SD | BD  | Age | SD |
|-----|-----|----|-----|-----|----|
| <15 | n/a | —  | 41  | 181 | 0  |
| 15  | 73  | 0  | 42  | 185 | 0  |
| 16  | 77  | 0  | 43  | 189 | 0  |
| 17  | 81  | 0  | 44  | 193 | 0  |
| 18  | 85  | 0  | 45  | 197 | 0  |
| 19  | 89  | 0  | 46  | 201 | 0  |
| 20  | 93  | 0  | 47  | 206 | 0  |
| 21  | 97  | 0  | 48  | 210 | 0  |
| 22  | 102 | 0  | 49  | 214 | 0  |
| 23  | 106 | 0  | 50  | 218 | 0  |
| 24  | 110 | 0  | 51  | 222 | 0  |
| 25  | 114 | 0  | 52  | 226 | 0  |
| 26  | 118 | 0  | 53  | 231 | 0  |
| 27  | 122 | 0  | 54  | 235 | 0  |
| 28  | 127 | 0  | 55  | 239 | 0  |
| 29  | 131 | 0  | 56  | 243 | 0  |
| 30  | 135 | 0  | 57  | 247 | 0  |
| 31  | 139 | 0  | 58  | 251 | 0  |
| 32  | 143 | 0  | 59  | 256 | 0  |
| 33  | 147 | 0  | 60  | 260 | 0  |
| 34  | 152 | 0  | 61  | 264 | 0  |
| 35  | 156 | 0  | 62  | 268 | 0  |
| 36  | 160 | 0  | 63  | 272 | 0  |
| 37  | 164 | 0  | 64  | 276 | 0  |
| 38  | 168 | 0  | 65  | 281 | 0  |
| 39  | 172 | 0  | >65 | n/a | —  |
| 40  | 177 | 0  |     |     |    |

Table 18–23. BD : Jeanty

| BPD | Age | SD | BPD | Age | SD | BPD | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|
| <10 | n/a | —  | 39  | 119 | 3  | 69  | 187 | 5  |
| 10  | 64  | 3  | 40  | 121 | 3  | 70  | 189 | 5  |
| 11  | 65  | 3  | 41  | 123 | 3  | 71  | 192 | 5  |
| 12  | 67  | 3  | 42  | 125 | 3  | 72  | 195 | 5  |
| 13  | 69  | 3  | 43  | 127 | 3  | 73  | 198 | 5  |
| 14  | 71  | 3  | 44  | 129 | 3  | 74  | 201 | 5  |
| 15  | 73  | 3  | 45  | 131 | 3  | 75  | 204 | 5  |
| 16  | 75  | 3  | 46  | 133 | 3  | 76  | 206 | 5  |
| 17  | 77  | 3  | 47  | 135 | 3  | 77  | 209 | 5  |
| 18  | 79  | 3  | 48  | 137 | 3  | 78  | 212 | 5  |
| 19  | 81  | 3  | 49  | 139 | 4  | 79  | 215 | 5  |
| 20  | 83  | 3  | 50  | 141 | 4  | 80  | 218 | 5  |
| 21  | 85  | 3  | 51  | 143 | 4  | 81  | 221 | 5  |
| 22  | 87  | 3  | 52  | 146 | 4  | 82  | 224 | 5  |
| 23  | 89  | 3  | 53  | 148 | 4  | 83  | 227 | 5  |
| 24  | 90  | 3  | 54  | 151 | 4  | 84  | 230 | 5  |
| 25  | 92  | 3  | 55  | 153 | 4  | 85  | 234 | 5  |
| 26  | 94  | 3  | 56  | 155 | 4  | 86  | 237 | 5  |
| 27  | 96  | 3  | 57  | 158 | 4  | 87  | 240 | 5  |
| 28  | 98  | 3  | 58  | 160 | 4  | 88  | 244 | 5  |
| 29  | 100 | 3  | 59  | 163 | 4  | 89  | 247 | 5  |
| 30  | 102 | 3  | 60  | 165 | 4  | 90  | 251 | 5  |
| 31  | 104 | 3  | 61  | 167 | 4  | 91  | 254 | 5  |
| 32  | 106 | 3  | 62  | 170 | 4  | 92  | 257 | 5  |
| 33  | 108 | 3  | 63  | 172 | 4  | 93  | 261 | 5  |
| 34  | 110 | 3  | 64  | 175 | 4  | 94  | 264 | 5  |
| 35  | 112 | 3  | 65  | 177 | 4  | 95  | 268 | 5  |
| 36  | 114 | 3  | 66  | 179 | 4  | >95 | n/a | —  |
| 37  | 115 | 3  | 67  | 182 | 4  |     |     |    |
| 38  | 117 | 3  | 68  | 184 | 4  |     |     |    |

Table 18–24. BPD : Jeanty



## OB Tables

CRL : Jeanty  
Jeanty : Radiology  
143 : 513,1982

Unit : CRL (mm)  
Age (Day)  
SD (mm)

FL : Jeanty  
Jeanty : Radiology  
143 : 513,1982

Unit : CRL (mm)  
Age (Day)  
SD (mm)

| CRL | Age | SD | CRL | Age | SD |
|-----|-----|----|-----|-----|----|
| <5  | n/a | —  | 30  | 69  | 7  |
| 5   | 44  | 4  | 31  | 70  | 7  |
| 6   | 45  | 4  | 32  | 70  | 7  |
| 7   | 46  | 4  | 33  | 71  | 7  |
| 8   | 48  | 4  | 34  | 72  | 7  |
| 9   | 50  | 4  | 35  | 73  | 7  |
| 10  | 51  | 4  | 36  | 73  | 7  |
| 11  | 52  | 4  | 37  | 74  | 7  |
| 12  | 53  | 4  | 38  | 75  | 7  |
| 13  | 54  | 4  | 39  | 76  | 7  |
| 14  | 55  | 4  | 40  | 76  | 7  |
| 15  | 56  | 5  | 41  | 76  | 7  |
| 16  | 57  | 5  | 42  | 77  | 7  |
| 17  | 58  | 5  | 43  | 77  | 7  |
| 18  | 59  | 5  | 44  | 78  | 7  |
| 19  | 60  | 5  | 45  | 79  | 7  |
| 20  | 61  | 5  | 46  | 79  | 7  |
| 21  | 62  | 6  | 47  | 80  | 7  |
| 22  | 63  | 6  | 48  | 81  | 7  |
| 23  | 64  | 6  | 49  | 81  | 7  |
| 24  | 65  | 6  | 50  | 82  | 7  |
| 25  | 66  | 6  | 51  | 83  | 7  |
| 26  | 67  | 7  | 52  | 83  | 7  |
| 27  | 67  | 7  | 53  | 84  | 7  |
| 28  | 67  | 7  | 54  | 85  | 7  |
| 29  | 68  | 7  | >54 | n/a | —  |

Table 18–25. CRL : Jeanty

| FL  | Age | SD | FL | Age | SD | FL  | Age | SD |
|-----|-----|----|----|-----|----|-----|-----|----|
| <10 | n/a | —  | 34 | 146 | 5  | 59  | 215 | 5  |
| 10  | 88  | 4  | 35 | 149 | 5  | 60  | 218 | 5  |
| 11  | 90  | 4  | 36 | 152 | 5  | 61  | 221 | 5  |
| 12  | 92  | 4  | 37 | 154 | 5  | 62  | 224 | 5  |
| 13  | 95  | 4  | 38 | 157 | 5  | 63  | 227 | 5  |
| 14  | 97  | 4  | 39 | 160 | 5  | 64  | 230 | 5  |
| 15  | 100 | 4  | 40 | 162 | 5  | 65  | 233 | 5  |
| 16  | 102 | 4  | 41 | 165 | 5  | 66  | 236 | 5  |
| 17  | 104 | 4  | 42 | 168 | 5  | 67  | 239 | 5  |
| 18  | 107 | 4  | 43 | 170 | 5  | 68  | 242 | 5  |
| 19  | 109 | 4  | 44 | 173 | 5  | 69  | 245 | 5  |
| 20  | 112 | 4  | 45 | 176 | 5  | 70  | 248 | 5  |
| 21  | 114 | 4  | 46 | 178 | 5  | 71  | 251 | 5  |
| 22  | 116 | 4  | 47 | 181 | 5  | 72  | 254 | 5  |
| 23  | 119 | 4  | 48 | 184 | 5  | 73  | 257 | 5  |
| 24  | 121 | 4  | 49 | 186 | 5  | 74  | 260 | 5  |
| 25  | 124 | 4  | 50 | 189 | 5  | 75  | 263 | 5  |
| 26  | 126 | 4  | 51 | 192 | 5  | 76  | 266 | 5  |
| 27  | 128 | 4  | 52 | 194 | 5  | 77  | 270 | 5  |
| 28  | 131 | 4  | 53 | 197 | 5  | 78  | 274 | 5  |
| 29  | 133 | 4  | 54 | 200 | 5  | 79  | 277 | 5  |
| 30  | 136 | 5  | 55 | 203 | 5  | 80  | 280 | 5  |
| 31  | 138 | 5  | 56 | 206 | 5  | >80 | n/a | —  |
| 32  | 141 | 5  | 57 | 209 | 5  |     |     |    |
| 33  | 144 | 5  | 58 | 212 | 5  |     |     |    |

Table 18–26. FL : Jeanty

## OB Tables

HC : Jeanty  
Jeanty : Radiology  
143 : 513,1982

Unit : HC (mm)  
Age (Day)  
SD (mm)

BPD : Kurtz

Unit : BPD (mm)  
Age (Day)  
SD (mm)

| HC  | Age | SD | HC   | Age | SD |
|-----|-----|----|------|-----|----|
| <80 | n/a | —  | 225  | 169 | 22 |
| 80  | 93  | 15 | 230  | 173 | 22 |
| 85  | 95  | 16 | 235  | 176 | 22 |
| 90  | 98  | 18 | 240  | 180 | 22 |
| 95  | 100 | 18 | 245  | 184 | 22 |
| 100 | 102 | 18 | 250  | 188 | 22 |
| 105 | 104 | 18 | 255  | 192 | 22 |
| 110 | 107 | 18 | 260  | 196 | 22 |
| 115 | 109 | 18 | 265  | 200 | 22 |
| 120 | 111 | 18 | 270  | 204 | 22 |
| 125 | 114 | 18 | 275  | 208 | 22 |
| 130 | 116 | 18 | 280  | 213 | 22 |
| 135 | 119 | 19 | 285  | 217 | 22 |
| 140 | 121 | 19 | 290  | 221 | 22 |
| 145 | 124 | 19 | 295  | 225 | 22 |
| 150 | 126 | 19 | 300  | 230 | 22 |
| 155 | 129 | 19 | 305  | 235 | 22 |
| 160 | 131 | 20 | 310  | 240 | 22 |
| 165 | 134 | 20 | 315  | 245 | 22 |
| 170 | 136 | 20 | 320  | 250 | 22 |
| 175 | 139 | 20 | 325  | 255 | 22 |
| 180 | 141 | 20 | 330  | 260 | 22 |
| 185 | 144 | 21 | 335  | 265 | 22 |
| 190 | 146 | 21 | 340  | 270 | 22 |
| 195 | 149 | 21 | 345  | 275 | 22 |
| 200 | 151 | 21 | 350  | 281 | 22 |
| 205 | 154 | 22 | 355  | 286 | 22 |
| 210 | 157 | 22 | 360  | 292 | 22 |
| 215 | 161 | 22 | >360 | n/a | —  |
| 220 | 165 | 22 |      |     |    |

Table 18–27. HC : Jeanty

| BPD | Age | SD | BPD | Age | SD | BPD | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|
| <21 | n/a | —  | 47  | 139 | 4  | 74  | 204 | 5  |
| 21  | 84  | 4  | 48  | 141 | 4  | 75  | 207 | 5  |
| 22  | 87  | 4  | 49  | 143 | 4  | 76  | 210 | 5  |
| 23  | 91  | 4  | 50  | 145 | 4  | 77  | 213 | 5  |
| 24  | 93  | 4  | 51  | 147 | 4  | 78  | 216 | 5  |
| 25  | 95  | 4  | 52  | 149 | 4  | 79  | 219 | 5  |
| 26  | 97  | 4  | 53  | 151 | 4  | 80  | 222 | 5  |
| 27  | 99  | 4  | 54  | 153 | 4  | 81  | 225 | 5  |
| 28  | 101 | 4  | 55  | 155 | 5  | 82  | 229 | 5  |
| 29  | 103 | 4  | 56  | 157 | 5  | 83  | 232 | 5  |
| 30  | 105 | 4  | 57  | 160 | 5  | 84  | 235 | 5  |
| 31  | 107 | 4  | 58  | 162 | 5  | 85  | 238 | 5  |
| 32  | 109 | 4  | 59  | 165 | 5  | 86  | 241 | 5  |
| 33  | 111 | 4  | 60  | 168 | 5  | 87  | 244 | 5  |
| 34  | 113 | 4  | 61  | 170 | 5  | 88  | 248 | 5  |
| 35  | 115 | 4  | 62  | 173 | 5  | 89  | 252 | 5  |
| 36  | 117 | 4  | 63  | 175 | 5  | 90  | 257 | 5  |
| 37  | 119 | 4  | 64  | 178 | 5  | 91  | 262 | 5  |
| 38  | 121 | 4  | 65  | 181 | 5  | 92  | 267 | 5  |
| 39  | 123 | 4  | 66  | 183 | 5  | 93  | 272 | 5  |
| 40  | 125 | 4  | 67  | 186 | 5  | 94  | 276 | 5  |
| 41  | 127 | 4  | 68  | 188 | 5  | 95  | 280 | 5  |
| 42  | 129 | 4  | 69  | 191 | 5  | 96  | 284 | 5  |
| 43  | 131 | 4  | 70  | 194 | 5  | 97  | 288 | 5  |
| 44  | 133 | 4  | 71  | 196 | 5  | 98  | 293 | 5  |
| 45  | 135 | 4  | 72  | 199 | 5  | >98 | n/a | —  |
| 46  | 137 | 4  | 73  | 201 | 5  |     |     |    |

Table 18–28. BPD : Kurtz



## OB Tables

CRL : Nelson

Unit : CRL (mm)  
Age (Day)  
SD (mm)

BPD : Osaka

Unit : BPD (mm)  
Age (Day)  
SD (mm)

| CRL | Age | SD | CRL | Age | SD | CRL | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|
| <3  | n/a | —  | 27  | 67  | 5  | 52  | 82  | 7  |
| 3   | 53  | 4  | 28  | 67  | 5  | 53  | 82  | 7  |
| 4   | 53  | 4  | 29  | 68  | 5  | 54  | 83  | 7  |
| 5   | 54  | 4  | 30  | 69  | 5  | 55  | 84  | 7  |
| 6   | 54  | 4  | 31  | 69  | 6  | 56  | 84  | 7  |
| 7   | 55  | 4  | 32  | 70  | 6  | 57  | 85  | 7  |
| 8   | 55  | 4  | 33  | 70  | 6  | 58  | 85  | 7  |
| 9   | 56  | 4  | 34  | 71  | 6  | 59  | 86  | 7  |
| 10  | 57  | 4  | 35  | 72  | 6  | 60  | 87  | 7  |
| 11  | 57  | 4  | 36  | 72  | 6  | 61  | 87  | 7  |
| 12  | 58  | 4  | 37  | 73  | 6  | 62  | 88  | 7  |
| 13  | 58  | 4  | 38  | 73  | 6  | 63  | 88  | 7  |
| 14  | 59  | 4  | 39  | 74  | 6  | 64  | 89  | 7  |
| 15  | 60  | 4  | 40  | 75  | 6  | 65  | 90  | 7  |
| 16  | 60  | 4  | 41  | 75  | 6  | 66  | 90  | 7  |
| 17  | 61  | 5  | 42  | 76  | 6  | 67  | 91  | 7  |
| 18  | 61  | 5  | 43  | 76  | 6  | 68  | 91  | 7  |
| 19  | 62  | 5  | 44  | 77  | 6  | 69  | 92  | 7  |
| 20  | 63  | 5  | 45  | 78  | 7  | 70  | 93  | 7  |
| 21  | 63  | 5  | 46  | 78  | 7  | 71  | 93  | 7  |
| 22  | 64  | 5  | 47  | 79  | 7  | 72  | 94  | 7  |
| 23  | 64  | 5  | 48  | 79  | 7  | 73  | 95  | 7  |
| 24  | 65  | 5  | 49  | 80  | 7  | >73 | n/a | —  |
| 25  | 66  | 5  | 50  | 81  | 7  |     |     |    |
| 26  | 66  | 5  | 51  | 81  | 7  |     |     |    |

Table 18–29. CRL : Nelson

| BPD | Age | SD  | BPD | Age | SD  | BPD | Age | SD  |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| <13 | n/a | —   | 40  | 122 | 2.6 | 68  | 185 | 3.3 |
| 13  | 70  | 1.9 | 41  | 124 | 2.7 | 69  | 187 | 3.3 |
| 14  | 71  | 1.9 | 42  | 126 | 2.7 | 70  | 190 | 3.4 |
| 15  | 73  | 1.9 | 43  | 128 | 2.7 | 71  | 193 | 3.4 |
| 16  | 75  | 1.9 | 44  | 130 | 2.7 | 72  | 195 | 3.4 |
| 17  | 77  | 2.0 | 45  | 132 | 2.8 | 73  | 198 | 3.4 |
| 18  | 78  | 2.0 | 46  | 135 | 2.8 | 74  | 200 | 3.5 |
| 19  | 80  | 2.0 | 47  | 137 | 2.8 | 75  | 203 | 3.5 |
| 20  | 82  | 2.1 | 48  | 139 | 2.8 | 76  | 206 | 3.5 |
| 21  | 84  | 2.1 | 49  | 141 | 2.9 | 77  | 209 | 3.5 |
| 22  | 86  | 2.1 | 50  | 143 | 2.9 | 78  | 212 | 3.5 |
| 23  | 88  | 2.1 | 51  | 145 | 2.9 | 79  | 214 | 3.6 |
| 24  | 90  | 2.2 | 52  | 148 | 2.9 | 80  | 217 | 3.6 |
| 25  | 92  | 2.2 | 53  | 150 | 3.0 | 81  | 220 | 3.6 |
| 26  | 94  | 2.2 | 54  | 152 | 3.0 | 82  | 224 | 3.6 |
| 27  | 96  | 2.3 | 55  | 154 | 3.0 | 83  | 227 | 3.6 |
| 28  | 98  | 2.3 | 56  | 157 | 3.0 | 84  | 230 | 3.7 |
| 29  | 99  | 2.3 | 57  | 159 | 3.1 | 85  | 234 | 3.7 |
| 30  | 101 | 2.3 | 58  | 161 | 3.1 | 86  | 237 | 3.7 |
| 31  | 103 | 2.4 | 59  | 164 | 3.1 | 87  | 238 | 3.7 |
| 32  | 105 | 2.4 | 60  | 166 | 3.1 | 88  | 245 | 3.7 |
| 33  | 107 | 2.4 | 61  | 168 | 3.2 | 89  | 249 | 3.8 |
| 34  | 109 | 2.5 | 62  | 171 | 3.2 | 90  | 254 | 3.8 |
| 35  | 112 | 2.5 | 63  | 173 | 3.2 | 91  | 259 | 3.8 |
| 36  | 114 | 2.5 | 64  | 175 | 3.2 | 92  | 265 | 3.8 |
| 37  | 116 | 2.5 | 65  | 178 | 3.3 | 93  | 273 | 3.9 |
| 38  | 118 | 2.6 | 66  | 180 | 3.3 | 94  | 280 | 3.9 |
| 39  | 120 | 2.6 | 67  | 182 | 3.3 | >94 | n/a | —   |

Table 18–30. BPD : Osaka

OB Tables

CRL : Osaka

Unit : CRL (mm)  
Age (Day)  
SD (mm)

EFBW : Osaka  
(gram)

Unit : EFBW

Age (Day)  
SD (gram)

| CRL | Age | SD  | CRL | Age | SD  | CRL | Age | SD  |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| <9  | n/a | —   | 27  | 68  | 4.5 | 46  | 80  | 6.3 |
| 9   | 50  | 1.7 | 28  | 69  | 4.6 | 47  | 80  | 6.3 |
| 10  | 52  | 2.0 | 29  | 69  | 4.6 | 48  | 81  | 6.4 |
| 11  | 53  | 2.2 | 30  | 70  | 4.8 | 49  | 82  | 6.6 |
| 12  | 55  | 2.5 | 31  | 71  | 4.9 | 50  | 83  | 6.7 |
| 13  | 56  | 2.6 | 32  | 71  | 4.9 | 51  | 83  | 6.7 |
| 14  | 57  | 2.8 | 33  | 72  | 5.1 | 52  | 83  | 6.7 |
| 15  | 58  | 2.9 | 34  | 73  | 5.2 | 53  | 84  | 6.9 |
| 16  | 59  | 3.1 | 35  | 73  | 5.2 | 54  | 85  | 7.0 |
| 17  | 60  | 3.2 | 36  | 74  | 5.4 | 55  | 85  | 7.0 |
| 18  | 61  | 3.4 | 37  | 74  | 5.4 | 56  | 86  | 7.2 |
| 19  | 62  | 3.5 | 38  | 75  | 5.5 | 57  | 86  | 7.2 |
| 20  | 63  | 3.7 | 39  | 76  | 5.7 | 58  | 87  | 7.3 |
| 21  | 63  | 3.7 | 40  | 76  | 5.7 | 59  | 87  | 7.3 |
| 22  | 64  | 3.8 | 41  | 77  | 5.8 | 60  | 88  | 7.5 |
| 23  | 65  | 4.0 | 42  | 77  | 5.8 | 61  | 89  | 7.6 |
| 24  | 66  | 4.1 | 43  | 78  | 6.0 | 62  | 89  | 7.6 |
| 25  | 66  | 4.1 | 44  | 79  | 6.1 | 63  | 90  | 7.8 |
| 26  | 67  | 4.3 | 45  | 79  | 6.1 | >63 | n/a | —   |

Table 18–31. CRL : Osaka

| EFBW | AGE | SD | EFBW | AGE | SD  | EFBW | AGE | SD  |
|------|-----|----|------|-----|-----|------|-----|-----|
| <137 | n/a | —  | 430  | 148 | 61  | 780  | 171 | 102 |
| 137  | 112 | 29 | 440  | 149 | 63  | 800  | 173 | 106 |
| 140  | 113 | 29 | 450  | 149 | 63  | 820  | 174 | 108 |
| 150  | 115 | 29 | 460  | 150 | 65  | 840  | 175 | 110 |
| 160  | 116 | 30 | 470  | 151 | 66  | 860  | 176 | 112 |
| 170  | 118 | 30 | 480  | 152 | 68  | 880  | 177 | 114 |
| 180  | 120 | 31 | 490  | 153 | 69  | 900  | 178 | 116 |
| 190  | 121 | 32 | 500  | 153 | 69  | 920  | 179 | 118 |
| 200  | 123 | 33 | 510  | 154 | 71  | 940  | 180 | 120 |
| 210  | 124 | 34 | 520  | 155 | 73  | 960  | 181 | 123 |
| 220  | 126 | 35 | 530  | 155 | 73  | 980  | 182 | 125 |
| 230  | 127 | 36 | 540  | 156 | 74  | 1000 | 183 | 127 |
| 240  | 128 | 37 | 550  | 157 | 76  | 1020 | 185 | 131 |
| 250  | 130 | 39 | 560  | 157 | 76  | 1040 | 186 | 133 |
| 260  | 131 | 40 | 570  | 158 | 78  | 1060 | 187 | 135 |
| 270  | 132 | 41 | 580  | 159 | 80  | 1080 | 188 | 138 |
| 280  | 133 | 42 | 590  | 160 | 81  | 1100 | 189 | 140 |
| 290  | 134 | 43 | 600  | 160 | 81  | 1120 | 190 | 142 |
| 300  | 135 | 44 | 610  | 161 | 83  | 1140 | 191 | 144 |
| 310  | 136 | 45 | 620  | 162 | 85  | 1160 | 192 | 146 |
| 320  | 137 | 46 | 630  | 162 | 85  | 1180 | 193 | 149 |
| 330  | 138 | 48 | 640  | 163 | 87  | 1200 | 194 | 151 |
| 340  | 139 | 49 | 650  | 164 | 89  | 1220 | 195 | 153 |
| 350  | 140 | 50 | 660  | 164 | 89  | 1240 | 195 | 153 |
| 360  | 141 | 51 | 670  | 165 | 91  | 1260 | 196 | 155 |
| 370  | 142 | 53 | 680  | 165 | 91  | 1280 | 197 | 158 |
| 380  | 143 | 54 | 690  | 166 | 92  | 1300 | 198 | 160 |
| 390  | 144 | 56 | 700  | 167 | 94  | 1320 | 199 | 162 |
| 400  | 145 | 57 | 720  | 168 | 96  | 1340 | 200 | 164 |
| 410  | 146 | 58 | 740  | 169 | 98  | 1360 | 201 | 167 |
| 420  | 147 | 60 | 760  | 170 | 100 | 1380 | 202 | 169 |

Table 18–32. EFBW : Osaka



## OB Tables

EFBW : Osaka (cont'd)  
EFBW (gram)

Unit :

FL : Osaka

Unit : FL (mm)  
Age (Day)  
SD (mm)

Age (Day)  
SD (gram)

| EFBW | AGE | SD  | EFBW | AGE | SD  | EFBW  | AGE | SD  |
|------|-----|-----|------|-----|-----|-------|-----|-----|
| 1400 | 203 | 171 | 2020 | 229 | 234 | 2640  | 254 | 302 |
| 1420 | 203 | 171 | 2040 | 229 | 234 | 2660  | 254 | 302 |
| 1440 | 204 | 174 | 2060 | 230 | 237 | 2680  | 255 | 305 |
| 1460 | 205 | 176 | 2080 | 231 | 239 | 2700  | 256 | 308 |
| 1480 | 206 | 178 | 2100 | 232 | 242 | 2720  | 257 | 311 |
| 1500 | 207 | 181 | 2120 | 233 | 244 | 2740  | 258 | 314 |
| 1520 | 208 | 183 | 2140 | 233 | 244 | 2760  | 259 | 317 |
| 1540 | 209 | 185 | 2160 | 234 | 247 | 2780  | 259 | 317 |
| 1560 | 210 | 188 | 2180 | 235 | 250 | 2800  | 260 | 320 |
| 1580 | 210 | 188 | 2200 | 236 | 252 | 2820  | 261 | 323 |
| 1600 | 211 | 190 | 2220 | 236 | 252 | 2840  | 262 | 326 |
| 1620 | 212 | 192 | 2240 | 237 | 255 | 2860  | 263 | 329 |
| 1640 | 213 | 195 | 2260 | 238 | 257 | 2880  | 264 | 332 |
| 1660 | 214 | 197 | 2280 | 239 | 260 | 2900  | 265 | 335 |
| 1680 | 215 | 200 | 2300 | 240 | 263 | 2920  | 266 | 339 |
| 1700 | 216 | 202 | 2320 | 241 | 265 | 2940  | 266 | 339 |
| 1720 | 216 | 202 | 2340 | 241 | 265 | 2960  | 267 | 342 |
| 1740 | 217 | 204 | 2360 | 242 | 268 | 2980  | 268 | 345 |
| 1760 | 218 | 207 | 2380 | 243 | 271 | 3000  | 269 | 348 |
| 1780 | 219 | 209 | 2400 | 244 | 274 | 3020  | 270 | 352 |
| 1800 | 220 | 212 | 2420 | 245 | 276 | 3040  | 271 | 355 |
| 1820 | 220 | 212 | 2440 | 245 | 276 | 3060  | 272 | 358 |
| 1840 | 221 | 214 | 2460 | 246 | 279 | 3080  | 273 | 362 |
| 1860 | 222 | 217 | 2480 | 247 | 282 | 3100  | 274 | 365 |
| 1880 | 223 | 219 | 2500 | 248 | 285 | 3120  | 275 | 369 |
| 1900 | 224 | 222 | 2520 | 249 | 288 | 3140  | 276 | 372 |
| 1920 | 224 | 222 | 2540 | 249 | 288 | 3160  | 277 | 376 |
| 1940 | 225 | 224 | 2560 | 250 | 290 | 3180  | 278 | 379 |
| 1960 | 226 | 227 | 2580 | 251 | 293 | 3200  | 279 | 383 |
| 1980 | 227 | 229 | 2600 | 252 | 296 | 3220  | 280 | 387 |
| 2000 | 228 | 232 | 2620 | 253 | 299 | >3220 | n/a | —   |

Table 18–32. EFBW : Osaka (cont'd)

| FL | Age | SD  | FL | Age | SD  | FL  | Age | SD  |
|----|-----|-----|----|-----|-----|-----|-----|-----|
| <9 | n/a | —   | 30 | 140 | 2.4 | 52  | 202 | 2.6 |
| 9  | 91  | 2.1 | 31 | 142 | 2.4 | 53  | 205 | 2.8 |
| 10 | 93  | 2.1 | 32 | 145 | 2.4 | 54  | 209 | 2.8 |
| 11 | 95  | 2.1 | 33 | 147 | 2.4 | 55  | 212 | 2.8 |
| 12 | 97  | 2.2 | 34 | 150 | 2.4 | 56  | 216 | 2.8 |
| 13 | 99  | 2.2 | 35 | 152 | 2.5 | 57  | 220 | 2.8 |
| 14 | 102 | 2.2 | 36 | 155 | 2.5 | 58  | 223 | 2.9 |
| 15 | 104 | 2.2 | 37 | 158 | 2.5 | 59  | 227 | 2.9 |
| 16 | 106 | 2.2 | 38 | 162 | 2.5 | 60  | 230 | 2.9 |
| 17 | 108 | 2.2 | 39 | 163 | 2.5 | 61  | 235 | 2.9 |
| 18 | 110 | 2.2 | 40 | 166 | 2.5 | 62  | 239 | 2.9 |
| 19 | 113 | 2.2 | 41 | 169 | 2.6 | 63  | 242 | 3.0 |
| 20 | 115 | 2.3 | 42 | 172 | 2.6 | 64  | 247 | 3.0 |
| 21 | 118 | 2.3 | 43 | 175 | 2.6 | 65  | 250 | 3.0 |
| 22 | 120 | 2.3 | 44 | 178 | 2.6 | 66  | 255 | 3.0 |
| 23 | 122 | 2.3 | 45 | 181 | 2.6 | 67  | 258 | 3.0 |
| 24 | 125 | 2.3 | 46 | 184 | 2.6 | 68  | 260 | 3.1 |
| 25 | 127 | 2.3 | 47 | 186 | 2.6 | 69  | 269 | 3.1 |
| 26 | 130 | 2.3 | 48 | 190 | 2.7 | 70  | 274 | 3.1 |
| 27 | 132 | 2.3 | 49 | 193 | 2.7 | 71  | 279 | 3.2 |
| 28 | 135 | 2.4 | 50 | 196 | 2.7 | >71 | n/a | —   |
| 29 | 137 | 2.4 | 51 | 199 | 2.7 |     |     |     |

Table 18–33. FL : Osaka



# OB Tables

FTA : Osaka

Unit : FTA (mm<sup>2</sup>)  
Age (Day)  
SD (mm<sup>2</sup>)

HL : Osaka

Unit : HL (mm)  
Age (Day)  
SD (mm)

| FTA  | Age | SD  | FTA  | Age | SD  | FTA   | Age | SD   |
|------|-----|-----|------|-----|-----|-------|-----|------|
| <560 | n/a | —   | 3300 | 175 | 400 | 6200  | 231 | 710  |
| 560  | 98  | 120 | 3400 | 177 | 410 | 6300  | 233 | 720  |
| 600  | 100 | 120 | 3500 | 179 | 420 | 6400  | 235 | 730  |
| 700  | 103 | 130 | 3600 | 181 | 430 | 6500  | 237 | 750  |
| 800  | 108 | 150 | 3700 | 183 | 440 | 6600  | 238 | 750  |
| 900  | 113 | 160 | 3800 | 185 | 450 | 6700  | 240 | 760  |
| 1000 | 115 | 170 | 3900 | 187 | 460 | 6800  | 242 | 780  |
| 1100 | 117 | 170 | 4000 | 189 | 470 | 6900  | 244 | 790  |
| 1200 | 122 | 190 | 4100 | 191 | 480 | 7000  | 246 | 800  |
| 1300 | 125 | 200 | 4200 | 193 | 490 | 7100  | 248 | 820  |
| 1400 | 128 | 210 | 4300 | 195 | 500 | 7200  | 250 | 830  |
| 1500 | 130 | 220 | 4400 | 197 | 510 | 7300  | 252 | 840  |
| 1600 | 134 | 230 | 4500 | 199 | 520 | 7400  | 254 | 860  |
| 1700 | 137 | 240 | 4600 | 201 | 530 | 7500  | 256 | 870  |
| 1800 | 139 | 250 | 4700 | 203 | 540 | 7600  | 258 | 880  |
| 1900 | 142 | 260 | 4800 | 205 | 560 | 7700  | 260 | 900  |
| 2000 | 145 | 270 | 4900 | 207 | 570 | 7800  | 262 | 910  |
| 2100 | 147 | 280 | 5000 | 209 | 580 | 7900  | 264 | 930  |
| 2200 | 150 | 290 | 5100 | 211 | 590 | 8000  | 265 | 930  |
| 2300 | 152 | 300 | 5200 | 213 | 600 | 8100  | 268 | 960  |
| 2400 | 155 | 310 | 5300 | 215 | 610 | 8200  | 270 | 970  |
| 2500 | 157 | 330 | 5400 | 216 | 620 | 8300  | 273 | 990  |
| 2600 | 159 | 330 | 5500 | 218 | 630 | 8400  | 274 | 1000 |
| 2700 | 162 | 340 | 5600 | 220 | 640 | 8500  | 276 | 1010 |
| 2800 | 164 | 350 | 5700 | 222 | 650 | 8600  | 279 | 1040 |
| 2900 | 166 | 360 | 5800 | 224 | 670 | 8660  | 280 | 1040 |
| 3000 | 168 | 370 | 5900 | 226 | 680 | >8660 | n/a | —    |
| 3100 | 170 | 380 | 6000 | 227 | 680 |       |     |      |
| 3200 | 173 | 390 | 6100 | 229 | 700 |       |     |      |

Table 18–34. FTA : Osaka

| HL  | Age | SD  | HL  | Age | SD  |
|-----|-----|-----|-----|-----|-----|
| <10 | n/a | —   | 37  | 164 | 2.4 |
| 10  | 91  | 2.0 | 38  | 167 | 2.4 |
| 11  | 93  | 2.0 | 39  | 170 | 2.4 |
| 12  | 96  | 2.0 | 40  | 174 | 2.4 |
| 13  | 98  | 2.1 | 41  | 178 | 2.4 |
| 14  | 100 | 2.1 | 42  | 182 | 2.5 |
| 15  | 103 | 2.1 | 43  | 185 | 2.5 |
| 16  | 105 | 2.1 | 44  | 188 | 2.5 |
| 17  | 108 | 2.1 | 45  | 192 | 2.5 |
| 18  | 110 | 2.1 | 46  | 196 | 2.5 |
| 19  | 113 | 2.1 | 47  | 200 | 2.5 |
| 20  | 115 | 2.1 | 48  | 204 | 2.6 |
| 21  | 117 | 2.1 | 49  | 208 | 2.6 |
| 22  | 121 | 2.2 | 50  | 213 | 2.6 |
| 23  | 123 | 2.2 | 51  | 217 | 2.6 |
| 24  | 126 | 2.2 | 52  | 222 | 2.6 |
| 25  | 129 | 2.2 | 53  | 227 | 2.7 |
| 26  | 132 | 2.2 | 54  | 232 | 2.7 |
| 27  | 134 | 2.2 | 55  | 237 | 2.7 |
| 28  | 137 | 2.2 | 56  | 242 | 2.7 |
| 29  | 140 | 2.3 | 57  | 248 | 2.8 |
| 30  | 143 | 2.3 | 58  | 254 | 2.8 |
| 31  | 145 | 2.3 | 59  | 260 | 2.8 |
| 32  | 149 | 2.3 | 60  | 267 | 2.9 |
| 33  | 151 | 2.3 | 61  | 275 | 2.9 |
| 34  | 155 | 2.3 | 62  | 280 | 2.9 |
| 35  | 158 | 2.3 | >62 | n/a | —   |
| 36  | 161 | 2.4 |     |     |     |

Table 18–35. HL : Osaka



## OB Tables

BPD : Paris

Unit : BPD (mm)  
Age (Day)  
SD (mm)

CRL : Paris

Unit : CRL (mm)  
Age (Day)  
SD (mm)

| BPD | Age | SD | BPD | Age | SD | BPD | SD  | Age |
|-----|-----|----|-----|-----|----|-----|-----|-----|
| <13 | n/a | —  | 39  | 123 | 3  | 66  | 187 | 4   |
| 13  | 77  | 3  | 40  | 126 | 4  | 67  | 189 | 4   |
| 14  | 78  | 3  | 41  | 128 | 4  | 68  | 192 | 4   |
| 15  | 79  | 3  | 42  | 130 | 4  | 69  | 194 | 4   |
| 16  | 80  | 3  | 43  | 133 | 4  | 70  | 197 | 4   |
| 17  | 81  | 3  | 44  | 135 | 4  | 71  | 199 | 4   |
| 18  | 82  | 3  | 45  | 137 | 4  | 72  | 202 | 4   |
| 19  | 83  | 3  | 46  | 140 | 4  | 73  | 204 | 4   |
| 20  | 84  | 3  | 47  | 142 | 4  | 74  | 207 | 4   |
| 21  | 85  | 3  | 48  | 144 | 4  | 75  | 210 | 5   |
| 22  | 87  | 3  | 49  | 147 | 4  | 76  | 213 | 5   |
| 23  | 89  | 3  | 50  | 149 | 4  | 77  | 217 | 5   |
| 24  | 91  | 3  | 51  | 151 | 4  | 78  | 220 | 5   |
| 25  | 93  | 3  | 52  | 154 | 4  | 79  | 224 | 5   |
| 26  | 95  | 3  | 53  | 156 | 4  | 80  | 227 | 5   |
| 27  | 97  | 3  | 54  | 158 | 4  | 81  | 231 | 5   |
| 28  | 100 | 3  | 55  | 161 | 4  | 82  | 234 | 5   |
| 29  | 102 | 3  | 56  | 163 | 4  | 83  | 238 | 5   |
| 30  | 104 | 3  | 57  | 165 | 4  | 84  | 242 | 5   |
| 31  | 106 | 3  | 58  | 168 | 4  | 85  | 247 | 5   |
| 32  | 108 | 3  | 59  | 170 | 4  | 86  | 252 | 5   |
| 33  | 110 | 3  | 60  | 172 | 4  | 87  | 256 | 5   |
| 34  | 113 | 3  | 61  | 175 | 4  | 88  | 261 | 5   |
| 35  | 115 | 3  | 62  | 177 | 4  | 89  | 266 | 5   |
| 36  | 117 | 3  | 63  | 179 | 4  | 90  | 287 | 5   |
| 37  | 119 | 3  | 64  | 182 | 4  | >90 | n/a | —   |
| 38  | 121 | 3  | 65  | 184 | 4  |     |     |     |

Table 18–36. BPD : Paris

| CRL | Age | SD | CRL | Age | SD | CRL | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|
| <5  | n/a | —  | 32  | 70  | 7  | 60  | 86  | 7  |
| 5   | 42  | 4  | 33  | 70  | 7  | 61  | 87  | 7  |
| 6   | 43  | 4  | 34  | 71  | 7  | 62  | 87  | 7  |
| 7   | 44  | 4  | 35  | 71  | 7  | 63  | 88  | 7  |
| 8   | 46  | 4  | 36  | 72  | 7  | 64  | 88  | 7  |
| 9   | 47  | 4  | 37  | 73  | 7  | 65  | 89  | 7  |
| 10  | 49  | 4  | 38  | 73  | 7  | 66  | 89  | 7  |
| 11  | 50  | 4  | 39  | 74  | 7  | 67  | 90  | 7  |
| 12  | 51  | 4  | 40  | 74  | 7  | 68  | 90  | 7  |
| 13  | 52  | 4  | 41  | 75  | 7  | 69  | 91  | 7  |
| 14  | 53  | 4  | 42  | 76  | 7  | 70  | 91  | 7  |
| 15  | 54  | 4  | 43  | 76  | 7  | 71  | 91  | 7  |
| 16  | 55  | 5  | 44  | 77  | 7  | 72  | 92  | 7  |
| 17  | 56  | 5  | 45  | 77  | 7  | 73  | 92  | 7  |
| 18  | 57  | 5  | 46  | 78  | 7  | 74  | 93  | 7  |
| 19  | 58  | 6  | 47  | 79  | 7  | 75  | 93  | 7  |
| 20  | 59  | 6  | 48  | 79  | 7  | 76  | 94  | 7  |
| 21  | 60  | 6  | 49  | 80  | 7  | 77  | 94  | 7  |
| 22  | 61  | 6  | 50  | 80  | 7  | 78  | 94  | 7  |
| 23  | 63  | 6  | 51  | 81  | 7  | 79  | 95  | 7  |
| 24  | 63  | 7  | 52  | 82  | 7  | 80  | 95  | 7  |
| 25  | 64  | 7  | 53  | 82  | 7  | 81  | 96  | 7  |
| 26  | 65  | 7  | 54  | 83  | 7  | 82  | 96  | 7  |
| 27  | 66  | 7  | 55  | 84  | 7  | 83  | 97  | 7  |
| 28  | 66  | 7  | 56  | 84  | 7  | 84  | 97  | 7  |
| 29  | 67  | 7  | 57  | 85  | 7  | 85  | 98  | 7  |
| 30  | 68  | 7  | 58  | 85  | 7  | >85 | n/a | —  |
| 31  | 69  | 7  | 59  | 86  | 7  |     |     |    |

Table 18–37. CRL : Paris



## OB Tables

FL : Paris

Unit : FL (mm)  
Age (Day)  
SD (mm)

Ft : Paris

Unit : Ft (mm)  
Age (Day)  
SD (mm)

| FL  | Age | SD | FL | Age | SD | FL  | Age | SD |
|-----|-----|----|----|-----|----|-----|-----|----|
| <15 | n/a | —  | 36 | 150 | 5  | 58  | 213 | 5  |
| 15  | 98  | 4  | 37 | 153 | 5  | 59  | 217 | 5  |
| 16  | 100 | 4  | 38 | 156 | 5  | 60  | 219 | 5  |
| 17  | 102 | 4  | 39 | 159 | 5  | 61  | 221 | 5  |
| 18  | 105 | 4  | 40 | 161 | 5  | 62  | 224 | 5  |
| 19  | 107 | 4  | 41 | 164 | 5  | 63  | 231 | 5  |
| 20  | 109 | 4  | 42 | 167 | 5  | 64  | 234 | 5  |
| 21  | 112 | 4  | 43 | 170 | 5  | 65  | 238 | 5  |
| 22  | 114 | 4  | 44 | 172 | 5  | 66  | 241 | 5  |
| 23  | 116 | 4  | 45 | 175 | 5  | 67  | 245 | 5  |
| 24  | 119 | 4  | 46 | 178 | 5  | 68  | 248 | 5  |
| 25  | 121 | 4  | 47 | 181 | 5  | 69  | 252 | 5  |
| 26  | 123 | 4  | 48 | 183 | 5  | 70  | 255 | 5  |
| 27  | 126 | 5  | 49 | 186 | 5  | 71  | 259 | 5  |
| 28  | 128 | 5  | 50 | 189 | 5  | 72  | 262 | 5  |
| 29  | 131 | 5  | 51 | 192 | 5  | 73  | 266 | 5  |
| 30  | 134 | 5  | 52 | 194 | 5  | 74  | 271 | 5  |
| 31  | 137 | 5  | 53 | 197 | 5  | 75  | 276 | 5  |
| 32  | 139 | 5  | 54 | 200 | 5  | 76  | 281 | 5  |
| 33  | 142 | 5  | 55 | 203 | 5  | 77  | 287 | 5  |
| 34  | 145 | 5  | 56 | 206 | 5  | >77 | n/a | —  |
| 35  | 148 | 5  | 57 | 210 | 5  |     |     |    |

Table 18–38. FL : Paris

| FT  | Age | SD | FT | Age | SD | FT  | Age | SD |
|-----|-----|----|----|-----|----|-----|-----|----|
| <13 | n/a | —  | 34 | 144 | 4  | 56  | 199 | 4  |
| 13  | 91  | 2  | 35 | 147 | 4  | 57  | 202 | 4  |
| 14  | 94  | 2  | 36 | 149 | 4  | 58  | 205 | 4  |
| 15  | 97  | 2  | 37 | 151 | 4  | 59  | 208 | 4  |
| 16  | 100 | 2  | 38 | 154 | 4  | 60  | 211 | 4  |
| 17  | 103 | 3  | 39 | 156 | 4  | 61  | 215 | 4  |
| 18  | 106 | 3  | 40 | 158 | 4  | 62  | 218 | 4  |
| 19  | 109 | 3  | 41 | 161 | 4  | 63  | 221 | 4  |
| 20  | 112 | 4  | 42 | 163 | 4  | 64  | 224 | 4  |
| 21  | 114 | 4  | 43 | 165 | 4  | 65  | 227 | 4  |
| 22  | 116 | 4  | 44 | 168 | 4  | 66  | 231 | 5  |
| 23  | 119 | 4  | 45 | 170 | 4  | 67  | 234 | 5  |
| 24  | 121 | 4  | 46 | 173 | 4  | 68  | 238 | 5  |
| 25  | 123 | 4  | 47 | 175 | 4  | 69  | 242 | 5  |
| 26  | 126 | 4  | 48 | 178 | 4  | 70  | 246 | 5  |
| 27  | 128 | 4  | 49 | 180 | 4  | 71  | 250 | 5  |
| 28  | 130 | 4  | 50 | 183 | 4  | 72  | 254 | 5  |
| 29  | 133 | 4  | 51 | 185 | 4  | 73  | 258 | 5  |
| 30  | 135 | 4  | 52 | 188 | 4  | 74  | 262 | 5  |
| 31  | 137 | 4  | 53 | 190 | 4  | 75  | 266 | 6  |
| 32  | 140 | 4  | 54 | 193 | 4  | >75 | n/a | —  |
| 33  | 142 | 4  | 55 | 196 | 4  |     |     |    |

Table 18–39. Ft : Paris

Unit : GS (mm)  
2SD (mm or day)

2SD = day

Table 18–41. BPD : Rempen  
Known LMP (left)—Unknown LMP (right)

CRL : Rempen  
 Unit : GS (mm)      Age (Weeks/Days)      2SD (mm or day)  
 2SD = mm                      2SD = day

| CRL | Age   | 2SD | CRL | Age   | 2SD | CRL | Age   | 2SD | CRL | Age   | 2SD |
|-----|-------|-----|-----|-------|-----|-----|-------|-----|-----|-------|-----|
| <1  | n/a   | —   | 45  | 11W2D | 8   | <2  | n/a   | —   | 46  | 11W2D | 7   |
| 1   | 5W5D  | 8   | 46  | 11W2D | 8   | 2   | 6W0D  | 6   | 47  | 11W2D | 7   |
| 2   | 5W6D  | 8   | 47  | 11W3D | 8   | 3   | 6W1D  | 6   | 48  | 11W3D | 6   |
| 3   | 6W0D  | 8   | 48  | 11W4D | 8   | 4   | 6W2D  | 6   | 49  | 11W4D | 7   |
| 4   | 6W1D  | 8   | 49  | 11W4D | 8   | 5   | 6W3D  | 6   | 50  | 11W4D | 6   |
| 5   | 6W2D  | 8   | 50  | 11W5D | 8   | 6   | 6W4D  | 6   | 51  | 11W5D | 6   |
| 6   | 6W3D  | 8   | 51  | 11W6D | 8   | 7   | 6W5D  | 6   | 52  | 11W5D | 7   |
| 7   | 6W4D  | 8   | 52  | 12W0D | 8   | 8   | 6W6D  | 6   | 53  | 11W6D | 6   |
| 8   | 6W6D  | 8   | 53  | 12W0D | 8   | 9   | 7W0D  | 6   | 54  | 12W0D | 7   |
| 9   | 6W6D  | 8   | 54  | 12W1D | 8   | 10  | 7W1D  | 6   | 55  | 12W0D | 7   |
| 10  | 7W0D  | 8   | 55  | 12W2D | 8   | 11  | 7W2D  | 6   | 56  | 12W1D | 6   |
| 11  | 7W2D  | 8   | 56  | 12W2D | 8   | 12  | 7W3D  | 6   | 57  | 12W1D | 7   |
| 12  | 7W2D  | 8   | 57  | 12W3D | 8   | 13  | 7W4D  | 7   | 58  | 12W2D | 6   |
| 13  | 7W4D  | 8   | 58  | 12W3D | 8   | 14  | 7W5D  | 7   | 59  | 12W3D | 7   |
| 14  | 7W4D  | 8   | 59  | 12W4D | 8   | 15  | 7W6D  | 7   | 60  | 12W3D | 6   |
| 15  | 7W5D  | 8   | 60  | 12W5D | 8   | 16  | 7W6D  | 7   | 61  | 12W4D | 7   |
| 16  | 7W6D  | 8   | 61  | 12W5D | 8   | 17  | 8W0D  | 7   | 62  | 12W4D | 6   |
| 17  | 8W0D  | 8   | 62  | 12W6D | 8   | 18  | 8W1D  | 6   | 63  | 12W5D | 7   |
| 18  | 8W1D  | 8   | 63  | 13W0D | 8   | 19  | 8W2D  | 6   | 64  | 12W5D | 7   |
| 19  | 8W2D  | 8   | 64  | 13W0D | 8   | 20  | 8W3D  | 6   | 65  | 12W6D | 6   |
| 20  | 8W3D  | 8   | 65  | 13W1D | 8   | 21  | 8W4D  | 7   | 66  | 12W6D | 7   |
| 21  | 8W4D  | 8   | 66  | 13W2D | 8   | 22  | 8W5D  | 7   | 67  | 13W0D | 6   |
| 22  | 8W5D  | 8   | >66 | n/a   | —   | 23  | 8W5D  | 7   | 68  | 13W0D | 7   |
| 23  | 8W6D  | 8   |     |       |     | 24  | 8W6D  | 7   | 69  | 13W1D | 6   |
| 24  | 8W6D  | 8   |     |       |     | 25  | 9W0D  | 6   | 70  | 13W1D | 7   |
| 25  | 9W0D  | 8   |     |       |     | 26  | 9W1D  | 6   | 71  | 13W2D | 7   |
| 26  | 9W1D  | 8   |     |       |     | 27  | 9W2D  | 7   | 72  | 13W2D | 6   |
| 27  | 9W2D  | 8   |     |       |     | 28  | 9W3D  | 7   | 73  | 13W3D | 7   |
| 28  | 9W3D  | 8   |     |       |     | 29  | 9W3D  | 7   | 74  | 13W3D | 6   |
| 29  | 9W4D  | 8   |     |       |     | 30  | 9W4D  | 7   | 75  | 13W4D | 7   |
| 30  | 9W4D  | 8   |     |       |     | 31  | 9W5D  | 7   | 76  | 13W4D | 6   |
| 31  | 9W5D  | 8   |     |       |     | 32  | 9W6D  | 7   | 77  | 13W4D | 7   |
| 32  | 9W6D  | 8   |     |       |     | 33  | 9W6D  | 7   | 78  | 13W5D | 6   |
| 33  | 10W0D | 8   |     |       |     | 34  | 10W0D | 6   | >78 | n/a   | —   |
| 34  | 10W1D | 8   |     |       |     | 35  | 10W1D | 6   |     |       |     |
| 35  | 10W1D | 8   |     |       |     | 36  | 10W2D | 7   |     |       |     |
| 36  | 10W2D | 8   |     |       |     | 37  | 10W2D | 7   |     |       |     |
| 37  | 10W3D | 8   |     |       |     | 38  | 10W3D | 6   |     |       |     |
| 38  | 10W4D | 8   |     |       |     | 39  | 10W4D | 6   |     |       |     |
| 39  | 10W4D | 8   |     |       |     | 40  | 10W5D | 7   |     |       |     |
| 40  | 10W5D | 8   |     |       |     | 41  | 10W5D | 7   |     |       |     |
| 41  | 10W6D | 8   |     |       |     | 42  | 10W6D | 6   |     |       |     |
| 42  | 11W0D | 8   |     |       |     | 43  | 11W0D | 7   |     |       |     |
| 43  | 11W0D | 8   |     |       |     | 44  | 11W0D | 7   |     |       |     |
| 44  | 11W1D | 8   |     |       |     | 45  | 11W1D | 6   |     |       |     |

Table 18–42. CRL : Rempen—Known LMP (left)—Unknown LMP (right)



## OB Tables

GS : Rempen

Unit : GS (mm)

2SD = mm

Age (Weeks/Days)

2SD (mm or day)

2SD = day

| GS | Age   | 2SD | GS  | Age   | 2SD | GS | Age  | 2SD | GS  | Age   | 2SD |
|----|-------|-----|-----|-------|-----|----|------|-----|-----|-------|-----|
| <1 | n/a   | —   | 45  | 10W2D | 11  | <1 | n/a  | —   | 45  | 10W0D | 10  |
| 1  | 4W4D  | 11  | 46  | 10W3D | 11  | 1  | 4W5D | 10  | 46  | 10W1D | 10  |
| 2  | 4W5D  | 11  | 47  | 10W4D | 11  | 2  | 4W6D | 10  | 47  | 10W2D | 10  |
| 3  | 4W6D  | 11  | 48  | 10W6D | 11  | 3  | 5W0D | 10  | 48  | 10W3D | 10  |
| 4  | 5W0D  | 11  | 49  | 11W0D | 11  | 4  | 5W1D | 10  | 49  | 10W4D | 10  |
| 5  | 5W0D  | 11  | 50  | 11W1D | 11  | 5  | 5W2D | 10  | 50  | 10W5D | 10  |
| 6  | 5W1D  | 11  | 51  | 11W2D | 11  | 6  | 5W2D | 10  | 51  | 10W6D | 10  |
| 7  | 5W2D  | 11  | 52  | 11W4D | 11  | 7  | 5W3D | 10  | 52  | 11W0D | 10  |
| 8  | 5W3D  | 11  | 53  | 11W5D | 11  | 8  | 5W4D | 10  | 53  | 11W1D | 10  |
| 9  | 5W3D  | 11  | 54  | 12W0D | 11  | 9  | 5W5D | 10  | 54  | 11W2D | 10  |
| 10 | 5W4D  | 11  | 55  | 12W1D | 11  | 10 | 5W5D | 10  | 55  | 11W3D | 10  |
| 11 | 5W5D  | 11  | 56  | 12W2D | 11  | 11 | 5W6D | 10  | 56  | 11W4D | 10  |
| 12 | 5W6D  | 11  | 57  | 12W4D | 11  | 12 | 6W0D | 10  | 57  | 11W5D | 10  |
| 13 | 6W0D  | 11  | 58  | 12W5D | 11  | 13 | 6W1D | 10  | 58  | 11W6D | 10  |
| 14 | 6W0D  | 11  | 59  | 13W0D | 11  | 14 | 6W2D | 10  | 59  | 12W0D | 10  |
| 15 | 6W1D  | 11  | 60  | 13W1D | 11  | 15 | 6W2D | 10  | 60  | 12W1D | 10  |
| 16 | 6W2D  | 11  | >60 | n/a   | —   | 16 | 6W3D | 10  | 61  | 12W2D | 10  |
| 17 | 6W3D  | 11  |     |       |     | 17 | 6W4D | 10  | 62  | 12W3D | 10  |
| 18 | 6W4D  | 11  |     |       |     | 18 | 6W5D | 10  | 63  | 12W4D | 10  |
| 19 | 6W5D  | 11  |     |       |     | 19 | 6W6D | 10  | 64  | 12W5D | 10  |
| 20 | 6W6D  | 11  |     |       |     | 20 | 6W6D | 10  | 65  | 12W6D | 10  |
| 21 | 6W6D  | 11  |     |       |     | 21 | 7W0D | 10  | 66  | 13W0D | 10  |
| 22 | 7W0D  | 11  |     |       |     | 22 | 7W1D | 10  | 67  | 13W1D | 10  |
| 23 | 7W1D  | 11  |     |       |     | 23 | 7W2D | 10  | 68  | 13W2D | 10  |
| 24 | 7W2D  | 11  |     |       |     | 24 | 7W3D | 10  | 69  | 13W3D | 10  |
| 25 | 7W3D  | 11  |     |       |     | 25 | 7W4D | 10  | 70  | 13W4D | 10  |
| 26 | 7W4D  | 11  |     |       |     | 26 | 7W4D | 10  | 71  | 13W5D | 10  |
| 27 | 7W5D  | 11  |     |       |     | 27 | 7W5D | 10  | 72  | 14W0D | 10  |
| 28 | 7W6D  | 11  |     |       |     | 28 | 7W6D | 10  | 73  | 14W1D | 10  |
| 29 | 8W0D  | 11  |     |       |     | 29 | 8W0D | 10  | >73 | n/a   | —   |
| 30 | 8W0D  | 11  |     |       |     | 30 | 8W1D | 10  |     |       |     |
| 31 | 8W1D  | 11  |     |       |     | 31 | 8W2D | 10  |     |       |     |
| 32 | 8W2D  | 11  |     |       |     | 32 | 8W3D | 10  |     |       |     |
| 33 | 8W3D  | 11  |     |       |     | 33 | 8W3D | 10  |     |       |     |
| 34 | 8W4D  | 11  |     |       |     | 34 | 8W4D | 10  |     |       |     |
| 35 | 8W5D  | 11  |     |       |     | 35 | 8W5D | 10  |     |       |     |
| 36 | 8W6D  | 11  |     |       |     | 36 | 8W6D | 10  |     |       |     |
| 37 | 9W0D  | 11  |     |       |     | 37 | 9W0D | 10  |     |       |     |
| 38 | 9W1D  | 11  |     |       |     | 38 | 9W1D | 10  |     |       |     |
| 39 | 9W2D  | 11  |     |       |     | 39 | 9W2D | 10  |     |       |     |
| 40 | 9W4D  | 11  |     |       |     | 40 | 9W3D | 10  |     |       |     |
| 41 | 9W5D  | 11  |     |       |     | 41 | 9W4D | 10  |     |       |     |
| 42 | 9W6D  | 11  |     |       |     | 42 | 9W5D | 10  |     |       |     |
| 43 | 10W0D | 11  |     |       |     | 43 | 9W6D | 10  |     |       |     |
| 44 | 10W1D | 11  |     |       |     | 44 | 9W6D | 10  |     |       |     |

Table 18–43. GS : Rempen—Known LMP (left)—Unknown LMP (right)

## OB Tables

CRL : Robinson  
Robinson :  
Robinson and AI  
BrJGynecol,82  
702, 1975

Unit : CRL (mm)  
Age (Day)  
SD (mm)

AC : Sostoa  
Sostoa : Hospital  
de la Santa Cruz y San  
Pablo, serviejo de  
obst.ygynecol

Unit : AC (mm)  
Age (Day)  
SD (mm)

| CRL | Age | SD | CRL | Age | SD | CRL | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|
| <7  | n/a | —  | 32  | 69  | 7  | 58  | 85  | 7  |
| 7   | 45  | 4  | 33  | 70  | 7  | 59  | 85  | 7  |
| 8   | 46  | 4  | 34  | 70  | 7  | 60  | 86  | 7  |
| 9   | 47  | 4  | 35  | 71  | 7  | 61  | 86  | 7  |
| 10  | 48  | 4  | 36  | 72  | 7  | 62  | 87  | 7  |
| 11  | 50  | 4  | 37  | 72  | 7  | 63  | 88  | 7  |
| 12  | 52  | 4  | 38  | 73  | 7  | 64  | 89  | 7  |
| 13  | 53  | 4  | 39  | 74  | 7  | 65  | 90  | 7  |
| 14  | 54  | 4  | 40  | 74  | 7  | 66  | 90  | 7  |
| 15  | 55  | 4  | 41  | 75  | 7  | 67  | 90  | 7  |
| 16  | 56  | 4  | 42  | 75  | 7  | 68  | 91  | 7  |
| 17  | 57  | 4  | 43  | 76  | 7  | 69  | 91  | 7  |
| 18  | 58  | 4  | 44  | 77  | 7  | 70  | 91  | 7  |
| 19  | 59  | 4  | 45  | 77  | 7  | 71  | 92  | 7  |
| 20  | 60  | 4  | 46  | 78  | 7  | 72  | 92  | 7  |
| 21  | 60  | 4  | 47  | 79  | 7  | 73  | 93  | 7  |
| 22  | 61  | 4  | 48  | 79  | 7  | 74  | 93  | 7  |
| 23  | 62  | 4  | 49  | 80  | 7  | 75  | 93  | 7  |
| 24  | 63  | 5  | 50  | 81  | 7  | 76  | 94  | 7  |
| 25  | 64  | 5  | 51  | 82  | 7  | 77  | 94  | 7  |
| 26  | 64  | 5  | 52  | 83  | 7  | 78  | 95  | 7  |
| 27  | 65  | 5  | 53  | 83  | 7  | 79  | 95  | 7  |
| 28  | 66  | 6  | 54  | 83  | 7  | 80  | 96  | 7  |
| 29  | 67  | 6  | 55  | 84  | 7  | 81  | 97  | 7  |
| 30  | 68  | 6  | 56  | 84  | 7  | 82  | 98  | 7  |
| 31  | 69  | 7  | 57  | 84  | 7  | >82 | n/a | —  |

Table 18–44. CRL : Robinson

| AC  | Age | SD | AC   | Age | SD |
|-----|-----|----|------|-----|----|
| <70 | n/a | —  | 210  | 180 | 22 |
| 70  | 98  | 22 | 215  | 184 | 22 |
| 75  | 101 | 22 | 220  | 187 | 22 |
| 80  | 105 | 22 | 225  | 190 | 22 |
| 85  | 107 | 22 | 230  | 194 | 22 |
| 90  | 109 | 22 | 235  | 197 | 22 |
| 95  | 110 | 22 | 240  | 200 | 22 |
| 100 | 113 | 22 | 245  | 204 | 22 |
| 105 | 116 | 22 | 250  | 207 | 22 |
| 110 | 119 | 22 | 255  | 210 | 22 |
| 115 | 123 | 22 | 260  | 214 | 22 |
| 120 | 126 | 22 | 265  | 217 | 22 |
| 125 | 129 | 22 | 270  | 220 | 22 |
| 130 | 132 | 22 | 275  | 223 | 22 |
| 135 | 136 | 22 | 280  | 229 | 22 |
| 140 | 139 | 22 | 285  | 233 | 22 |
| 145 | 141 | 22 | 290  | 236 | 22 |
| 150 | 143 | 22 | 295  | 239 | 22 |
| 155 | 145 | 22 | 300  | 242 | 22 |
| 160 | 147 | 22 | 305  | 245 | 22 |
| 165 | 150 | 22 | 310  | 248 | 22 |
| 170 | 153 | 22 | 315  | 251 | 22 |
| 175 | 157 | 22 | 320  | 256 | 22 |
| 180 | 160 | 22 | 325  | 261 | 22 |
| 185 | 163 | 22 | 330  | 266 | 22 |
| 190 | 167 | 22 | 335  | 271 | 22 |
| 195 | 170 | 22 | 340  | 276 | 22 |
| 200 | 173 | 22 | 344  | 280 | 22 |
| 205 | 177 | 22 | >344 | n/a | —  |

Table 18–45. AC : Sostoa



## OB Tables

BD : Sostoa                      Unit :    BD (mm)  
 Sostoa : Hospital              Age (Day)  
 de la Santa Cruz y San        SD (mm)  
 Pablo, serviejo de  
 obst.ygynecol

BPD : Sostoa                    Unit :    BPD (mm)  
 Sostoa : Hospital              Age (Day)  
 de la Santa Cruz y San        SD (mm)  
 Pablo, serviejo de  
 obst.ygynecol

| BD  | Age | SD | BD  | Age | SD |
|-----|-----|----|-----|-----|----|
| <23 | n/a | —  | 42  | 189 | 0  |
| 23  | 119 | 0  | 43  | 192 | 0  |
| 24  | 122 | 0  | 44  | 196 | 0  |
| 25  | 126 | 0  | 45  | 199 | 0  |
| 26  | 128 | 0  | 46  | 203 | 0  |
| 27  | 130 | 0  | 47  | 206 | 0  |
| 28  | 133 | 0  | 48  | 210 | 0  |
| 29  | 136 | 0  | 49  | 213 | 0  |
| 30  | 140 | 0  | 50  | 217 | 0  |
| 31  | 143 | 0  | 51  | 224 | 0  |
| 32  | 147 | 0  | 52  | 231 | 0  |
| 33  | 150 | 0  | 53  | 238 | 0  |
| 34  | 154 | 0  | 54  | 245 | 0  |
| 35  | 157 | 0  | 55  | 252 | 0  |
| 36  | 161 | 0  | 56  | 255 | 0  |
| 37  | 164 | 0  | 57  | 259 | 0  |
| 38  | 168 | 0  | 58  | 266 | 0  |
| 39  | 175 | 0  | 59  | 273 | 0  |
| 40  | 182 | 0  | 60  | 280 | 0  |
| 41  | 185 | 0  | >60 | n/a | —  |

Table 18–46. BD : Sostoa

| BPD | Age | SD | BPD | Age | SD | BPD | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|
| <24 | n/a | —  | 46  | 147 | 5  | 69  | 198 | 5  |
| 24  | 98  | 4  | 47  | 149 | 5  | 70  | 200 | 5  |
| 25  | 100 | 4  | 48  | 151 | 5  | 71  | 203 | 5  |
| 26  | 102 | 4  | 49  | 153 | 5  | 72  | 206 | 5  |
| 27  | 104 | 4  | 50  | 155 | 5  | 73  | 209 | 5  |
| 28  | 106 | 4  | 51  | 158 | 5  | 74  | 212 | 5  |
| 29  | 108 | 4  | 52  | 160 | 5  | 75  | 215 | 5  |
| 30  | 110 | 4  | 53  | 162 | 5  | 76  | 218 | 5  |
| 31  | 112 | 4  | 54  | 164 | 5  | 77  | 221 | 5  |
| 32  | 114 | 4  | 55  | 167 | 5  | 78  | 224 | 5  |
| 33  | 116 | 4  | 56  | 169 | 5  | 79  | 227 | 5  |
| 34  | 119 | 4  | 57  | 171 | 5  | 80  | 230 | 5  |
| 35  | 121 | 4  | 58  | 173 | 5  | 81  | 233 | 5  |
| 36  | 123 | 4  | 59  | 176 | 5  | 82  | 236 | 5  |
| 37  | 125 | 4  | 60  | 178 | 5  | 83  | 239 | 5  |
| 38  | 127 | 4  | 61  | 180 | 5  | 84  | 242 | 5  |
| 39  | 129 | 4  | 62  | 182 | 5  | 85  | 245 | 5  |
| 40  | 131 | 4  | 63  | 185 | 5  | 86  | 252 | 5  |
| 41  | 133 | 4  | 64  | 187 | 5  | 87  | 259 | 5  |
| 42  | 135 | 4  | 65  | 189 | 5  | 88  | 266 | 5  |
| 43  | 137 | 4  | 66  | 191 | 5  | 89  | 273 | 5  |
| 44  | 140 | 5  | 67  | 194 | 5  | 90  | 280 | 5  |
| 45  | 143 | 5  | 68  | 196 | 5  | >90 | n/a | —  |

Table 18–47. BPD : Sostoa



## OB Tables

FL : Sostoa  
Sostoa : Hospital  
de la Santa Cruz y San  
Pablo, serviejo de  
obst.ygynecol

Unit : FL (mm)  
Age (Day)  
SD (mm)

HC : Sostoa

Unit : HC (mm)  
Age (Day)  
SD (mm)

| FL  | Age | SD | FL  | Age | SD |
|-----|-----|----|-----|-----|----|
| <10 | n/a | —  | 41  | 178 | 4  |
| 10  | 98  | 4  | 42  | 181 | 4  |
| 11  | 100 | 4  | 43  | 183 | 4  |
| 12  | 103 | 4  | 44  | 186 | 4  |
| 13  | 105 | 4  | 45  | 189 | 5  |
| 14  | 108 | 4  | 46  | 192 | 5  |
| 15  | 111 | 4  | 47  | 196 | 5  |
| 16  | 113 | 4  | 48  | 199 | 5  |
| 17  | 116 | 4  | 49  | 203 | 5  |
| 18  | 118 | 4  | 50  | 206 | 5  |
| 19  | 121 | 4  | 51  | 210 | 5  |
| 20  | 124 | 4  | 52  | 213 | 5  |
| 21  | 126 | 4  | 53  | 217 | 5  |
| 22  | 129 | 4  | 54  | 220 | 5  |
| 23  | 131 | 4  | 55  | 224 | 5  |
| 24  | 134 | 4  | 56  | 227 | 5  |
| 25  | 137 | 4  | 57  | 231 | 5  |
| 26  | 139 | 4  | 58  | 234 | 5  |
| 27  | 142 | 4  | 59  | 238 | 5  |
| 28  | 144 | 4  | 60  | 241 | 5  |
| 29  | 147 | 4  | 61  | 245 | 5  |
| 30  | 150 | 4  | 62  | 248 | 5  |
| 31  | 152 | 4  | 63  | 252 | 5  |
| 32  | 155 | 4  | 64  | 255 | 5  |
| 33  | 157 | 4  | 65  | 259 | 5  |
| 34  | 160 | 4  | 66  | 263 | 5  |
| 35  | 163 | 4  | 67  | 267 | 5  |
| 36  | 165 | 4  | 68  | 271 | 5  |
| 37  | 168 | 4  | 69  | 275 | 5  |
| 38  | 170 | 4  | 70  | 280 | 5  |
| 39  | 173 | 4  | >70 | n/a | —  |
| 40  | 176 | 4  |     |     |    |

Table 18–48. FL : Sostoa

| HC  | Age | SD | HC   | Age | SD |
|-----|-----|----|------|-----|----|
| <93 | n/a | —  | 220  | 171 | 22 |
| 93  | 98  | 15 | 225  | 174 | 22 |
| 95  | 99  | 15 | 230  | 177 | 22 |
| 100 | 104 | 17 | 235  | 180 | 22 |
| 105 | 106 | 18 | 240  | 183 | 22 |
| 110 | 109 | 19 | 245  | 186 | 22 |
| 115 | 111 | 19 | 250  | 189 | 22 |
| 120 | 113 | 20 | 255  | 192 | 22 |
| 125 | 116 | 21 | 260  | 195 | 22 |
| 130 | 118 | 21 | 265  | 198 | 22 |
| 135 | 121 | 22 | 270  | 201 | 22 |
| 140 | 124 | 22 | 275  | 204 | 22 |
| 145 | 127 | 22 | 280  | 207 | 22 |
| 150 | 130 | 22 | 285  | 210 | 22 |
| 155 | 133 | 22 | 290  | 215 | 22 |
| 160 | 136 | 22 | 295  | 219 | 22 |
| 165 | 139 | 22 | 300  | 224 | 22 |
| 170 | 142 | 22 | 305  | 228 | 22 |
| 175 | 145 | 22 | 310  | 233 | 22 |
| 180 | 148 | 22 | 315  | 240 | 22 |
| 185 | 151 | 22 | 320  | 247 | 22 |
| 190 | 154 | 22 | 325  | 254 | 22 |
| 195 | 157 | 22 | 330  | 261 | 22 |
| 200 | 160 | 22 | 335  | 268 | 22 |
| 205 | 163 | 22 | 340  | 275 | 22 |
| 210 | 165 | 22 | 343  | 280 | 22 |
| 215 | 168 | 22 | >343 | n/a | —  |

Table 18–49. HC : Sostoa



## OB Tables

OFD : Sostoa

Unit : OFD (mm)  
Age (Day)  
SD (mm)

BPD : Tokyo

Unit : BPD (mm)  
Age (Day)  
SD (Day)

| OFD | Age | SD | OFD | Age | SD | OFD  | Age | SD |
|-----|-----|----|-----|-----|----|------|-----|----|
| <28 | n/a | —  | 55  | 146 | 0  | 83   | 195 | 0  |
| 28  | 98  | 0  | 56  | 147 | 0  | 84   | 197 | 0  |
| 29  | 99  | 0  | 57  | 149 | 0  | 85   | 199 | 0  |
| 30  | 101 | 0  | 58  | 151 | 0  | 86   | 201 | 0  |
| 31  | 103 | 0  | 59  | 153 | 0  | 87   | 202 | 0  |
| 32  | 105 | 0  | 60  | 154 | 0  | 88   | 204 | 0  |
| 33  | 106 | 0  | 61  | 156 | 0  | 89   | 206 | 0  |
| 34  | 108 | 0  | 62  | 158 | 0  | 90   | 208 | 0  |
| 35  | 110 | 0  | 63  | 160 | 0  | 91   | 210 | 0  |
| 36  | 112 | 0  | 64  | 162 | 0  | 92   | 212 | 0  |
| 37  | 114 | 0  | 65  | 163 | 0  | 93   | 215 | 0  |
| 38  | 115 | 0  | 66  | 165 | 0  | 94   | 218 | 0  |
| 39  | 117 | 0  | 67  | 167 | 0  | 95   | 221 | 0  |
| 40  | 119 | 0  | 68  | 169 | 0  | 96   | 224 | 0  |
| 41  | 121 | 0  | 69  | 170 | 0  | 97   | 226 | 0  |
| 42  | 122 | 0  | 70  | 172 | 0  | 98   | 229 | 0  |
| 43  | 124 | 0  | 71  | 174 | 0  | 99   | 232 | 0  |
| 44  | 126 | 0  | 72  | 176 | 0  | 100  | 235 | 0  |
| 45  | 128 | 0  | 73  | 178 | 0  | 101  | 238 | 0  |
| 46  | 130 | 0  | 74  | 179 | 0  | 102  | 243 | 0  |
| 47  | 131 | 0  | 75  | 181 | 0  | 103  | 248 | 0  |
| 48  | 133 | 0  | 76  | 183 | 0  | 104  | 253 | 0  |
| 49  | 135 | 0  | 77  | 185 | 0  | 105  | 259 | 0  |
| 50  | 137 | 0  | 78  | 186 | 0  | 106  | 264 | 0  |
| 51  | 138 | 0  | 79  | 188 | 0  | 107  | 269 | 0  |
| 52  | 140 | 0  | 80  | 190 | 0  | 108  | 274 | 0  |
| 53  | 142 | 0  | 81  | 192 | 0  | 109  | 280 | 0  |
| 54  | 144 | 0  | 82  | 194 | 0  | >109 | n/a | —  |

Table 18–50. OFD : Sostoa

| BPD | Age | SD  | BPD | Age | SD  | BPD | Age | SD  |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| <20 | n/a | —   | 44  | 135 | ± 5 | 69  | 194 | ± 7 |
| 20  | 85  | ± 6 | 45  | 138 | ± 6 | 70  | 196 | ± 7 |
| 21  | 87  | ± 6 | 46  | 140 | ± 6 | 71  | 199 | ± 8 |
| 22  | 89  | ± 6 | 47  | 142 | ± 6 | 72  | 201 | ± 8 |
| 23  | 92  | ± 6 | 48  | 144 | ± 6 | 73  | 204 | ± 8 |
| 24  | 94  | ± 6 | 49  | 146 | ± 6 | 74  | 207 | ± 8 |
| 25  | 96  | ± 6 | 50  | 148 | ± 6 | 75  | 210 | ± 8 |
| 26  | 98  | ± 6 | 51  | 151 | ± 6 | 76  | 213 | ± 8 |
| 27  | 100 | ± 6 | 52  | 153 | ± 6 | 77  | 216 | ± 8 |
| 28  | 102 | ± 6 | 53  | 154 | ± 6 | 78  | 218 | ± 8 |
| 29  | 102 | ± 6 | 54  | 157 | ± 6 | 79  | 221 | ± 8 |
| 30  | 106 | ± 5 | 55  | 160 | ± 6 | 80  | 225 | ± 8 |
| 31  | 108 | ± 5 | 56  | 162 | ± 6 | 81  | 228 | ± 8 |
| 32  | 110 | ± 5 | 57  | 164 | ± 6 | 82  | 231 | ± 8 |
| 33  | 112 | ± 5 | 58  | 167 | ± 6 | 83  | 234 | ± 9 |
| 34  | 114 | ± 5 | 59  | 169 | ± 6 | 84  | 238 | ± 9 |
| 35  | 116 | ± 5 | 60  | 171 | ± 6 | 85  | 241 | ± 9 |
| 36  | 118 | ± 5 | 61  | 174 | ± 7 | 86  | 245 | ± 9 |
| 37  | 120 | ± 5 | 62  | 176 | ± 7 | 87  | 249 | ± 9 |
| 38  | 123 | ± 5 | 63  | 179 | ± 7 | 88  | 253 | ± 9 |
| 39  | 125 | ± 5 | 64  | 181 | ± 7 | 89  | 258 | ± 9 |
| 40  | 127 | ± 5 | 65  | 183 | ± 7 | 90  | 262 | ± 9 |
| 41  | 129 | ± 5 | 66  | 186 | ± 7 | >90 | n/a | —   |
| 42  | 131 | ± 5 | 67  | 188 | ± 7 |     |     |     |
| 43  | 133 | ± 5 | 68  | 191 | ± 7 |     |     |     |

Table 18–51. BPD : Tokyo



## OB Tables

CRL : Tokyo

Unit : CRL (mm)  
Age (Day)  
SD (Day)

FL : Tokyo

Unit : FL (mm)  
Age (Day)  
SD (Day)

| CRL | Age | SD   | CRL | Age | SD  |
|-----|-----|------|-----|-----|-----|
| <13 | n/a | —    | 32  | 73  | ± 7 |
| 13  | 55  | ± 8  | 33  | 74  | ± 7 |
| 14  | 56  | ± 9  | 34  | 74  | ± 7 |
| 15  | 57  | ± 10 | 35  | 75  | ± 7 |
| 16  | 58  | ± 8  | 36  | 76  | ± 7 |
| 17  | 59  | ± 9  | 37  | 77  | ± 7 |
| 18  | 60  | ± 10 | 38  | 78  | ± 7 |
| 19  | 61  | ± 8  | 39  | 78  | ± 7 |
| 20  | 62  | ± 9  | 40  | 79  | ± 7 |
| 21  | 63  | ± 7  | 41  | 80  | ± 7 |
| 22  | 64  | ± 7  | 42  | 81  | ± 7 |
| 23  | 65  | ± 7  | 43  | 81  | ± 7 |
| 24  | 66  | ± 7  | 44  | 82  | ± 7 |
| 25  | 67  | ± 7  | 45  | 83  | ± 7 |
| 26  | 68  | ± 7  | 46  | 84  | ± 7 |
| 27  | 68  | ± 7  | 47  | 84  | ± 7 |
| 28  | 69  | ± 7  | 48  | 85  | ± 7 |
| 29  | 70  | ± 7  | 49  | 86  | ± 7 |
| 30  | 71  | ± 7  | 50  | 86  | ± 7 |
| 31  | 72  | ± 7  | >50 | n/a | —   |

Table 18–52. CRL : Tokyo

| FL  | Age | SD  | FL  | Age | SD  |
|-----|-----|-----|-----|-----|-----|
| <33 | n/a | —   | 53  | 207 | ± 7 |
| 33  | 143 | ± 6 | 54  | 210 | ± 7 |
| 34  | 146 | ± 6 | 55  | 214 | ± 7 |
| 35  | 149 | ± 6 | 56  | 217 | ± 7 |
| 36  | 153 | ± 6 | 57  | 220 | ± 7 |
| 37  | 156 | ± 6 | 58  | 224 | ± 7 |
| 38  | 159 | ± 6 | 59  | 228 | ± 8 |
| 39  | 162 | ± 6 | 60  | 231 | ± 8 |
| 40  | 166 | ± 6 | 61  | 235 | ± 8 |
| 41  | 169 | ± 6 | 62  | 239 | ± 8 |
| 42  | 172 | ± 6 | 63  | 243 | ± 8 |
| 43  | 175 | ± 6 | 64  | 247 | ± 8 |
| 44  | 178 | ± 6 | 65  | 251 | ± 8 |
| 45  | 181 | ± 6 | 66  | 256 | ± 8 |
| 46  | 185 | ± 7 | 67  | 260 | ± 8 |
| 47  | 188 | ± 7 | 68  | 266 | ± 7 |
| 48  | 191 | ± 7 | 69  | 271 | ± 7 |
| 49  | 194 | ± 7 | 70  | 278 | ± 7 |
| 50  | 197 | ± 7 | 71  | 286 | ± 6 |
| 51  | 200 | ± 7 | >71 | n/a | —   |
| 52  | 204 | ± 7 |     |     |     |

Table 18–53. FL : Tokyo



## OB Tables

GS : Tokyo  
Tokyo University  
Method 1986,6 by  
Univ. of Tokyo

Unit : GS (mm)  
Age (Day)  
SD (Day)

LV : Tokyo

Unit : LV (mm)  
Age (Day)  
SD (Day)

| GS  | Age | SD   | GS  | Age | SD   |
|-----|-----|------|-----|-----|------|
| <12 | n/a | ——   | 32  | 55  | ± 0X |
| 12  | 31  | ± 7  | 33  | 56  | ± 0X |
| 13  | 32  | ± 7  | 34  | 57  | ± 0X |
| 14  | 33  | ± 7  | 35  | 58  | ± 0X |
| 15  | 34  | ± 7  | 36  | 59  | ± 0X |
| 16  | 36  | ± 7  | 37  | 60  | ± 0X |
| 17  | 37  | ± 7  | 38  | 61  | ± 0X |
| 18  | 38  | ± 7  | 39  | 62  | ± 0X |
| 19  | 40  | ± 7  | 40  | 63  | ± 0X |
| 20  | 41  | ± 7  | 41  | 64  | ± 0X |
| 21  | 42  | ± 7  | 42  | 65  | ± 0X |
| 22  | 43  | ± 7  | 43  | 65  | ± 0X |
| 23  | 44  | ± 7  | 44  | 66  | ± 0X |
| 24  | 46  | ± 7  | 45  | 67  | ± 0X |
| 25  | 47  | ± 7  | 46  | 68  | ± 0X |
| 26  | 48  | ± 8  | 47  | 69  | ± 0X |
| 27  | 49  | ± 9  | 48  | 70  | ± 0X |
| 28  | 50  | ± 10 | 49  | 71  | ± 0X |
| 29  | 51  | ± 0X | 50  | 72  | ± 0X |
| 30  | 52  | ± 0X | >50 | n/a | ——   |
| 31  | 53  | ± 0X | /   | /   | /    |

Table 18–54. GS : Tokyo

| LV  | Age | SD   | LV  | Age | SD   |
|-----|-----|------|-----|-----|------|
| <44 | n/a | ——   | 66  | 213 | ± 10 |
| 44  | 154 | ± 5  | 67  | 217 | ± 10 |
| 45  | 157 | ± 5  | 68  | 220 | ± 10 |
| 46  | 159 | ± 5  | 69  | 224 | ± 10 |
| 47  | 161 | ± 5  | 70  | 227 | ± 11 |
| 48  | 163 | ± 5  | 71  | 231 | ± 11 |
| 49  | 166 | ± 6  | 72  | 234 | ± 11 |
| 50  | 168 | ± 6  | 73  | 238 | ± 11 |
| 51  | 171 | ± 6  | 74  | 241 | ± 11 |
| 52  | 173 | ± 6  | 75  | 245 | ± 11 |
| 53  | 176 | ± 7  | 76  | 249 | ± 11 |
| 54  | 178 | ± 7  | 77  | 252 | ± 11 |
| 55  | 181 | ± 7  | 78  | 256 | ± 11 |
| 56  | 183 | ± 8  | 79  | 260 | ± 10 |
| 57  | 186 | ± 8  | 80  | 264 | ± 10 |
| 58  | 189 | ± 8  | 81  | 267 | ± 10 |
| 59  | 192 | ± 8  | 82  | 271 | ± 10 |
| 60  | 195 | ± 9  | 83  | 275 | ± 10 |
| 61  | 198 | ± 9  | 84  | 278 | ± 10 |
| 62  | 201 | ± 9  | 85  | 282 | ± 10 |
| 63  | 204 | ± 9  | 86  | 285 | ± 10 |
| 64  | 207 | ± 10 | >86 | n/a | ——   |
| 65  | 210 | ± 10 | /   | /   | /    |

Table 18–55. LV : Tokyo



---

*OB Tables*

---

FL/HC

HC/AC

GA (weeks)

GA (weeks)

| GA  | Min  | Max  | GA  | Min  | Max  |
|-----|------|------|-----|------|------|
| <15 | n/a  | ——   | 29  | 19.6 | 20.8 |
| 15  | 15.3 | 17.1 | 30  | 19.2 | 21.4 |
| 16  | 13.3 | 16.5 | 31  | 19.3 | 21.3 |
| 17  | 14.6 | 17.6 | 32  | 19.1 | 21.3 |
| 18  | 15.8 | 18.0 | 33  | 19.9 | 21.5 |
| 19  | 16.1 | 18.3 | 34  | 19.4 | 21.8 |
| 20  | 16.8 | 19.8 | 35  | 20.1 | 22.3 |
| 21  | 15.9 | 20.3 | 36  | 20.1 | 22.1 |
| 22  | 18.4 | 20.2 | 37  | 20.8 | 22.6 |
| 23  | 19.2 | 20.8 | 38  | 20.9 | 22.7 |
| 24  | 18.7 | 20.9 | 39  | 20.6 | 23.4 |
| 25  | 18.7 | 20.3 | 40  | 20.7 | 22.5 |
| 26  | 18.6 | 20.4 | 41  | 21.6 | 23.2 |
| 27  | 18.6 | 20.4 | 42  | 20.1 | 23.9 |
| 28  | 18.8 | 20.6 | >42 | n/a  | ——   |

Table 18–56. FL/HC

| GA  | Min  | Max  |
|-----|------|------|
| <13 | n/a  | ——   |
| 13  | 1.14 | 1.31 |
| 15  | 1.05 | 1.39 |
| 17  | 1.07 | 1.29 |
| 19  | 1.09 | 1.26 |
| 21  | 1.06 | 1.25 |
| 23  | 1.05 | 1.21 |
| 25  | 1.04 | 1.22 |
| 27  | 1.05 | 1.22 |
| 29  | 0.99 | 1.21 |
| 31  | 0.96 | 1.17 |
| 33  | 0.96 | 1.11 |
| 35  | 0.93 | 1.11 |
| 37  | 0.92 | 1.05 |
| 39  | 0.87 | 1.06 |
| 41  | 0.93 | 1.00 |
| >41 | n/a  | ——   |

Table 18–57. HC/AC



## OB Tables

EFW : Tokyo Shinozuka

Unit : EFW (grams)

Age (Day)

1SD (grams)

| EFW  | Age   | 1SD | EFW   | Age   | 1SD |
|------|-------|-----|-------|-------|-----|
| <250 | n/a   | —   | 2250  | 34W4D | 264 |
| 250  | 19W3D | 45  | 2300  | 34W6D | 269 |
| 300  | 20W0D | 51  | 2350  | 35W1D | 274 |
| 350  | 20W4D | 58  | 2400  | 35W3D | 279 |
| 400  | 21W2D | 66  | 2450  | 35W5D | 284 |
| 450  | 21W5D | 71  | 2500  | 35W7D | 290 |
| 500  | 22W2D | 78  | 2550  | 36W2D | 295 |
| 550  | 22W6D | 85  | 2600  | 36W4D | 301 |
| 600  | 23W2D | 90  | 2650  | 36W6D | 306 |
| 650  | 23W6D | 98  | 2700  | 37W2D | 314 |
| 700  | 24W2D | 103 | 2750  | 37W4D | 320 |
| 750  | 24W5D | 109 | 2800  | 37W6D | 325 |
| 800  | 25W2D | 116 | 2850  | 38W1D | 331 |
| 850  | 25W5D | 122 | 2900  | 38W4D | 340 |
| 900  | 26W1D | 128 | 2950  | 38W6D | 345 |
| 950  | 26W4D | 134 | 3000  | 39W2D | 354 |
| 1000 | 26W6D | 138 | >3000 | n/a   | —   |
| 1050 | 27W2D | 145 |       |       |     |
| 1100 | 27W5D | 151 |       |       |     |
| 1150 | 28W0D | 155 |       |       |     |
| 1200 | 28W3D | 162 |       |       |     |
| 1250 | 28W5D | 166 |       |       |     |
| 1300 | 29W1D | 173 |       |       |     |
| 1350 | 29W3D | 177 |       |       |     |
| 1400 | 29W5D | 181 |       |       |     |
| 1450 | 30W0D | 186 |       |       |     |
| 1500 | 30W2D | 191 |       |       |     |
| 1550 | 30W5D | 197 |       |       |     |
| 1600 | 31W0D | 202 |       |       |     |
| 1650 | 31W2D | 207 |       |       |     |
| 1700 | 31W4D | 211 |       |       |     |
| 1750 | 31W6D | 216 |       |       |     |
| 1800 | 32W1D | 221 |       |       |     |
| 1850 | 32W3D | 226 |       |       |     |
| 1900 | 32W5D | 231 |       |       |     |
| 1950 | 32W7D | 236 |       |       |     |
| 2000 | 33W1D | 238 |       |       |     |
| 2050 | 33W3D | 243 |       |       |     |
| 2100 | 33W5D | 248 |       |       |     |
| 2150 | 34W0D | 253 |       |       |     |
| 2200 | 34W2D | 258 |       |       |     |

Table 18–58. EFW : Tokyo Shinozuka

APTDxTTD (AxT) : Tokyo Shinozuka

Unit : APTDxTTD (mm)

Age (Day)

1SD (cm2)

| AxT | Age   | 1SD  | AxT | Age   | 1SD  |
|-----|-------|------|-----|-------|------|
| <10 | n/a   | —    | 90  | 39W2D | 12.0 |
| 10  | 16W1D | 2.5  | >90 | n/a   | —    |
| 12  | 17W0D | 2.7  |     |       |      |
| 14  | 17W6D | 2.9  |     |       |      |
| 16  | 18W4D | 3.1  |     |       |      |
| 18  | 19W3D | 3.4  |     |       |      |
| 20  | 20W1D | 3.6  |     |       |      |
| 22  | 20W6D | 3.8  |     |       |      |
| 24  | 21W4D | 4.0  |     |       |      |
| 26  | 22W2D | 4.3  |     |       |      |
| 28  | 22W6D | 4.4  |     |       |      |
| 30  | 23W4D | 4.7  |     |       |      |
| 32  | 24W1D | 4.9  |     |       |      |
| 34  | 24W5D | 5.1  |     |       |      |
| 36  | 25W2D | 5.3  |     |       |      |
| 38  | 25W6D | 5.5  |     |       |      |
| 40  | 26W3D | 5.7  |     |       |      |
| 42  | 27W0D | 6.0  |     |       |      |
| 44  | 27W3D | 6.1  |     |       |      |
| 46  | 28W0D | 6.4  |     |       |      |
| 48  | 28W4D | 6.6  |     |       |      |
| 50  | 29W0D | 6.8  |     |       |      |
| 52  | 29W3D | 7.0  |     |       |      |
| 54  | 30W0D | 7.2  |     |       |      |
| 56  | 30W3D | 7.4  |     |       |      |
| 58  | 31W0D | 7.7  |     |       |      |
| 60  | 31W3D | 7.9  |     |       |      |
| 62  | 31W6D | 8.1  |     |       |      |
| 64  | 32W3D | 8.4  |     |       |      |
| 66  | 32W6D | 8.6  |     |       |      |
| 68  | 33W3D | 8.8  |     |       |      |
| 70  | 33W6D | 9.1  |     |       |      |
| 72  | 34W2D | 9.3  |     |       |      |
| 74  | 34W6D | 9.6  |     |       |      |
| 76  | 35W3D | 9.9  |     |       |      |
| 78  | 35W6D | 10.1 |     |       |      |
| 80  | 36W3D | 10.2 |     |       |      |
| 82  | 37W0D | 10.7 |     |       |      |
| 84  | 37W4D | 11.0 |     |       |      |
| 86  | 38W1D | 11.3 |     |       |      |
| 88  | 38W5D | 11.7 |     |       |      |

Table 18–59. APTDxTTD (AxT): Tokyo Shinozuka

1SD (cm)

[illegible]

Table 18-61. AC : Tokyo Shinozuka



## OB Tables

BPD : Tokyo Shinozuka

Unit : (mm)

Age (Day)

1SD (mm)

| BPD | Age   | 1SD | BPD | Age   | 1SD |
|-----|-------|-----|-----|-------|-----|
| <13 | n/a   | —   | 53  | 22W1D | 3.0 |
| 13  | 10W1D | 2.3 | 54  | 22W3D | 3.0 |
| 14  | 10W3D | 2.3 | 55  | 22W5D | 3.0 |
| 15  | 10W5D | 2.3 | 56  | 23W1D | 3.0 |
| 16  | 11W0D | 2.3 | 57  | 23W3D | 3.0 |
| 17  | 11W2D | 2.4 | 58  | 23W5D | 3.1 |
| 18  | 11W4D | 2.4 | 59  | 24W1D | 3.1 |
| 19  | 11W6D | 2.4 | 60  | 24W3D | 3.1 |
| 20  | 12W1D | 2.4 | 61  | 24W5D | 3.1 |
| 21  | 12W3D | 2.4 | 62  | 25W1D | 3.1 |
| 22  | 12W6D | 2.4 | 63  | 25W3D | 3.1 |
| 23  | 13W1D | 2.5 | 64  | 25W5D | 3.2 |
| 24  | 13W3D | 2.5 | 65  | 26W1D | 3.2 |
| 25  | 13W5D | 2.5 | 66  | 26W3D | 3.2 |
| 26  | 14W0D | 2.5 | 67  | 26W6D | 3.2 |
| 27  | 14W2D | 2.5 | 68  | 27W2D | 3.3 |
| 28  | 14W4D | 2.5 | 69  | 27W4D | 3.3 |
| 29  | 14W6D | 2.6 | 70  | 28W0D | 3.3 |
| 30  | 15W1D | 2.6 | 71  | 28W3D | 3.3 |
| 31  | 15W3D | 2.6 | 72  | 28W5D | 3.3 |
| 32  | 15W5D | 2.6 | 73  | 29W1D | 3.4 |
| 33  | 16W0D | 2.6 | 74  | 29W4D | 3.4 |
| 34  | 16W2D | 2.6 | 75  | 30W0D | 3.4 |
| 35  | 16W4D | 2.7 | 76  | 30W3D | 3.4 |
| 36  | 16W6D | 2.7 | 77  | 30W6D | 3.4 |
| 37  | 17W1D | 2.7 | 78  | 31W2D | 3.5 |
| 38  | 17W4D | 2.7 | 79  | 31W5D | 3.5 |
| 39  | 17W6D | 2.7 | 80  | 32W1D | 3.5 |
| 40  | 18W1D | 2.7 | 81  | 32W5D | 3.6 |
| 41  | 18W3D | 2.8 | 82  | 33W1D | 3.6 |
| 42  | 18W5D | 2.8 | 83  | 33W5D | 3.6 |
| 43  | 19W0D | 2.8 | 84  | 34W2D | 3.6 |
| 44  | 19W2D | 2.8 | 85  | 34W6D | 3.7 |
| 45  | 19W4D | 2.8 | 86  | 35W3D | 3.7 |
| 46  | 20W0D | 2.8 | 87  | 36W0D | 3.7 |
| 47  | 20W2D | 2.9 | 88  | 36W5D | 3.8 |
| 48  | 20W4D | 2.9 | 89  | 37W4D | 3.8 |
| 49  | 20W6D | 2.9 | 90  | 38W3D | 3.9 |
| 50  | 21W1D | 2.9 | >90 | n/a   | —   |
| 51  | 21W3D | 2.9 |     |       |     |
| 52  | 21W6D | 2.9 |     |       |     |

Table 18–62. BPD : Tokyo Shinozuka

CRL : Tokyo Shinozuka

Unit : (mm)

Age (Day)

1SD (mm)

| CRL | Age   | 1SD | CRL | Age   | 1SD |
|-----|-------|-----|-----|-------|-----|
| <5  | n/a   | —   | 45  | 11W6D | 6.9 |
| 5   | 6W3D  | 1.1 | 46  | 11W6D | 6.9 |
| 6   | 6W4D  | 1.3 | 47  | 12W0D | 7.1 |
| 7   | 6W6D  | 1.6 | 48  | 12W1D | 7.2 |
| 8   | 7W0D  | 1.7 | 49  | 12W1D | 7.2 |
| 9   | 7W1D  | 1.9 | 50  | 12W2D | 7.4 |
| 10  | 7W2D  | 2.0 | >50 | n/a   | —   |
| 11  | 7W3D  | 2.2 |     |       |     |
| 12  | 7W4D  | 2.3 |     |       |     |
| 13  | 7W5D  | 2.5 |     |       |     |
| 14  | 7W6D  | 2.6 |     |       |     |
| 15  | 8W1D  | 2.9 |     |       |     |
| 16  | 8W2D  | 3.1 |     |       |     |
| 17  | 8W3D  | 3.3 |     |       |     |
| 18  | 8W4D  | 3.4 |     |       |     |
| 19  | 8W5D  | 3.6 |     |       |     |
| 20  | 8W6D  | 3.7 |     |       |     |
| 21  | 9W0D  | 3.9 |     |       |     |
| 22  | 9W1D  | 4.0 |     |       |     |
| 23  | 9W2D  | 4.2 |     |       |     |
| 24  | 9W3D  | 4.3 |     |       |     |
| 25  | 9W4D  | 4.5 |     |       |     |
| 26  | 9W4D  | 4.5 |     |       |     |
| 27  | 9W5D  | 4.6 |     |       |     |
| 28  | 9W6D  | 4.8 |     |       |     |
| 29  | 10W0D | 4.9 |     |       |     |
| 30  | 10W1D | 5.1 |     |       |     |
| 31  | 10W2D | 5.2 |     |       |     |
| 32  | 10W3D | 5.4 |     |       |     |
| 33  | 10W4D | 5.5 |     |       |     |
| 34  | 10W5D | 5.7 |     |       |     |
| 35  | 10W6D | 5.9 |     |       |     |
| 36  | 10W6D | 5.9 |     |       |     |
| 37  | 11W0D | 6.0 |     |       |     |
| 38  | 11W0D | 6.0 |     |       |     |
| 39  | 11W1D | 6.2 |     |       |     |
| 40  | 11W2D | 6.3 |     |       |     |
| 41  | 11W3D | 6.5 |     |       |     |
| 42  | 11W3D | 6.5 |     |       |     |
| 43  | 11W4D | 6.6 |     |       |     |
| 44  | 11W5D | 6.8 |     |       |     |

Table 18–63. CRL : Tokyo Shinozuka



AC : Australia  
Unit : (mm)  
Age (Day)  
±2SD (Day)

| AC   | Age | SD |
|------|-----|----|
| <35  | n/a | —  |
| 35   | 70  | 8  |
| 46   | 77  | 8  |
| 57   | 84  | 8  |
| 69   | 91  | 8  |
| 80   | 98  | 9  |
| 92   | 105 | 9  |
| 103  | 112 | 9  |
| 114  | 119 | 9  |
| 126  | 126 | 10 |
| 137  | 133 | 10 |
| 149  | 140 | 10 |
| 160  | 147 | 10 |
| 171  | 154 | 10 |
| 183  | 161 | 10 |
| 194  | 168 | 12 |
| 206  | 175 | 12 |
| 217  | 182 | 12 |
| 228  | 189 | 14 |
| 240  | 196 | 14 |
| 251  | 203 | 14 |
| 263  | 210 | 14 |
| 274  | 217 | 14 |
| 285  | 224 | 16 |
| 297  | 231 | 16 |
| 308  | 238 | 18 |
| 320  | 245 | 18 |
| 331  | 252 | 18 |
| 342  | 259 | 18 |
| 354  | 266 | 20 |
| 365  | 273 | 20 |
| 377  | 280 | 20 |
| >377 | n/a | —  |
|      |     |    |
|      |     |    |
|      |     |    |
|      |     |    |
|      |     |    |
|      |     |    |
|      |     |    |
|      |     |    |
|      |     |    |

Table 18–64. AC : Australia

CRL : Australia  
Unit : (mm)  
Age (Day)  
±2SD (Day)

| CRL | Age | SD | CRL | Age | SD |
|-----|-----|----|-----|-----|----|
| <2  | n/a | —  | 51  | 82  | *  |
| 2   | 42  | *  | 53  | 83  | *  |
| 3   | 43  | *  | 55  | 84  | *  |
| 4   | 44  | *  | 57  | 85  | *  |
| 5   | 45  | *  | 58  | 86  | *  |
| 6   | 46  | *  | 60  | 87  | *  |
| 7   | 47  | *  | 62  | 88  | *  |
| 8   | 48  | *  | 64  | 89  | *  |
| 9   | 49  | *  | 66  | 90  | *  |
| 10  | 50  | *  | 68  | 91  | *  |
| 11  | 51  | *  | 70  | 92  | *  |
| 12  | 52  | *  | 72  | 93  | *  |
| 13  | 53  | *  | 74  | 94  | *  |
| 14  | 54  | *  | 76  | 95  | *  |
| 15  | 55  | *  | 78  | 96  | *  |
| 16  | 56  | *  | 80  | 97  | *  |
| 17  | 57  | *  | 82  | 98  | *  |
| 18  | 58  | *  | >82 | n/a | —  |
| 19  | 59  | *  |     |     |    |
| 20  | 60  | *  |     |     |    |
| 22  | 61  | *  |     |     |    |
| 23  | 62  | *  |     |     |    |
| 24  | 63  | *  |     |     |    |
| 25  | 64  | *  |     |     |    |
| 26  | 65  | *  |     |     |    |
| 27  | 66  | *  |     |     |    |
| 29  | 67  | *  |     |     |    |
| 30  | 68  | *  |     |     |    |
| 31  | 69  | *  |     |     |    |
| 33  | 70  | *  |     |     |    |
| 34  | 71  | *  |     |     |    |
| 36  | 72  | *  |     |     |    |
| 37  | 73  | *  |     |     |    |
| 38  | 74  | *  |     |     |    |
| 40  | 75  | *  |     |     |    |
| 41  | 76  | *  |     |     |    |
| 43  | 77  | *  |     |     |    |
| 45  | 78  | *  |     |     |    |
| 46  | 79  | *  |     |     |    |
| 48  | 80  | *  |     |     |    |
| 50  | 81  | *  |     |     |    |

Table 18–65. CRL : Australia \* : No Data



## OB Tables

---

BPD : Australia

Unit : (mm)

Age (Day)

±2SD (Day)

| BPD | Age | 2SD | BPD | Age | 2SD |
|-----|-----|-----|-----|-----|-----|
| <20 | n/a | —   | 60  | 169 | 13  |
| 20  | 84  | 4   | 61  | 171 | 13  |
| 21  | 86  | 4   | 62  | 173 | 13  |
| 22  | 88  | 4   | 63  | 176 | 14  |
| 23  | 90  | 4   | 64  | 178 | 14  |
| 24  | 92  | 5   | 65  | 181 | 14  |
| 25  | 94  | 5   | 66  | 183 | 14  |
| 26  | 95  | 5   | 67  | 186 | 15  |
| 27  | 97  | 5   | 68  | 188 | 15  |
| 28  | 99  | 5   | 69  | 191 | 15  |
| 29  | 101 | 6   | 70  | 193 | 15  |
| 30  | 103 | 6   | 71  | 196 | 16  |
| 31  | 105 | 6   | 72  | 199 | 16  |
| 32  | 107 | 6   | 73  | 201 | 16  |
| 33  | 109 | 7   | 74  | 204 | 16  |
| 34  | 111 | 7   | 75  | 206 | 17  |
| 35  | 113 | 7   | 76  | 209 | 17  |
| 36  | 115 | 7   | 77  | 212 | 17  |
| 37  | 117 | 8   | 78  | 214 | 17  |
| 38  | 119 | 8   | 79  | 217 | 17  |
| 39  | 121 | 8   | 80  | 220 | 18  |
| 40  | 123 | 8   | 81  | 222 | 18  |
| 41  | 126 | 9   | 82  | 225 | 18  |
| 42  | 128 | 9   | 83  | 228 | 18  |
| 43  | 130 | 9   | 84  | 231 | 19  |
| 44  | 132 | 9   | 85  | 234 | *   |
| 45  | 134 | 9   | 86  | 237 | *   |
| 46  | 136 | 10  | 87  | 240 | *   |
| 47  | 139 | 10  | 88  | 244 | *   |
| 48  | 141 | 10  | 89  | 247 | *   |
| 49  | 143 | 10  | 90  | 251 | *   |
| 50  | 145 | 11  | 91  | 255 | *   |
| 51  | 147 | 11  | 92  | 259 | *   |
| 52  | 149 | 11  | 93  | 264 | *   |
| 53  | 152 | 11  | 94  | 270 | *   |
| 54  | 154 | 12  | 95  | 276 | *   |
| 55  | 157 | 12  | 96  | 284 | *   |
| 56  | 159 | 12  | 97  | 292 | *   |
| 57  | 161 | 12  | 98  | 301 | *   |
| 58  | 164 | 13  | >98 | n/a | —   |
| 59  | 166 | 13  |     |     |     |

Table 18–66. BPD : Australia \* : No Data



# VCR Operating Instructions

## Operating Manuals

To use the LOGIQ™ 500 ultrasound system and video cassette recorder (VCR) safely and properly, read this appendix and the SVO-9500MD VCR Operating Manual.

## Recording

When making important recordings, always make a trial recording in advance to ensure normal video and audio recording.

GE Medical Systems is not responsible for compensation for a recording failure resulting from a problem in the VCR or video tape during its use.

## Cassette tapes

Use S-VHS video cassette tapes.

## VTR-PB function

To use a search function in the LOGIQ™ 500 system, the optional VTR-PB function is required. For details of the VTR-PB function, contact the nearest GE Distributor, Affiliate or Sales Representative.



## **Introduction of VCR Features**

### **Remote Control Function**

Allows the LOGIQ™ 500 to remote control the following VCR functions:

1. Ext Video (for switching between the scan image and video image modes)
2. Play (playback)
3. Stop
4. Pause
5. Record
6. FF (fast forward)
7. REW (rewind)
8. FWD SEARCH (fast-forward search: forward speed is variable from 1/30 time to 10 times the normal playback speed)
9. REV SEARCH (reverse search: reverse speed is variable from 1/30 time to 10 times the normal playback speed)
10. FREEZE
11. FRAME SHIFT (frame forward/frame reverse)

### **Frame-forward search**

When VCR image playback is paused, each time FRAME SHIFT is pressed, the VCR will advance or rewind the paused playback image one frame.

### **Variable-speed search function**

When VCR image playback is paused, holding down FRAME SHIFT allows the playback speed to vary from 1/30 time to 10 times the normal playback rate.



## *VCR Operating Instructions*

---

### **Tape search function**

Allows the search of a tape number based on a patient name or scan date. The VCR allows a maximum of 2,600 tapes to be registered or searched.

The optional VTR–PB function is required.

### **Image search function**

Allows the advanced search of examination images based on a patient name, ID or scan date. The VCR allows a maximum of 39–patient examination images in each video tape to be searched.

The optional VTR–PB function is required.

### **Saving of data to be searched**

Writing the data to be searched to a video tape and saving it on a hard drive allows another LOGIQ™ 500 to perform an image search as well.

The optional VTR–PB function is required.

### **Playback with measurements**

Images recorded by the LOGIQ™ 500 or LOGIQ™ 700 allow a variety of frozen-image measurements and calculations including the measurement of distance, area, circumference, and other factors.

### **Image processing**

Allows the writing of comments or display of a body pattern on VCR playback images.



## *VCR Operating Instructions*

---

### **Safety**

### **Operating Precautions**

#### **Applicable VCR**

The Sony SVO-9500 MD is the only VCR recommended by GE as a LOGIQ™ 500-connectable recording device.

The Sony SVO-9500 MD requires an interface/metal fitting kit for installation on the CRT or console shelf. Commercially available VCRs do not have a dedicated interface kit; therefore, the LOGIQ™ 500 cannot support their remote control.

For more information, contact a GE Distributor, Affiliate or Sales Representative.

#### **Installation**

The VCR must be properly installed in an approved configuration using the dedicated VCR metal fittings to prevent accidents such as the VCR falling over, and to prevent electromagnetic interference. Metal-fittings must be installed by a qualified service engineer.

#### **Power supply**

To power the VCR, plug the power cable supplied with the VCR into the power outlet for peripherals on the rear panel of the LOGIQ™ 500. Do not use a power cable other than the power cable supplied with the VCR. Do not plug in the power cable to an outlet other than the power outlet for peripherals.



## *VCR Operating Instructions*

---

### **Connections**

Always use the cables supplied with the VCR to connect the VCR to the LOGIQ™ 500. Do not use any other cables. When the VCR is connected to the LOGIQ™ 500, do not connect any VCR options, such as the remote control unit, indicator unit and microphone.

### **Prohibition of modifications**

The video cassette recorder, metal-fittings, power cable and connection cables must not be modified in any way. Only the specified connections should be made.

### **Prevention of Electromagnetic Interference**

The LOGIQ™ 500 and the recommended peripherals, which are approved by GEMS, will meet Voluntary Control for Interference electronics Medical equipment (VCIM) standards established by Electronic Industries Association of Japan (EIAJ) for interference waves generated from a medical electronic device, provided that they are properly installed and connected using the specified connection cables.

Using the equipment under any other conditions may interfere with the receiving capability of radios, television sets and other equipment.

The LOGIQ™ 500 and VCR must be properly operated and used in accordance with the LOGIQ™ 500 User Manual and VCR Operating Instructions.

### **Operation of the VCR**

For the detailed operation of the VCR, see the SVO-9500MD Operating Manual provided for the VCR.

## Setting Up the VCR

### VCR Part Names and Functions/Settings

#### Part Names and Functions of the VCR Front Panel




Figure 19-1. VCR Front Panel

1. POWER switch. Setting fixed to ON. The POWER switch must always be in the ON position (recessed).
2. Cassette inlet.
3. EJECT button. Press this button to eject a video cassette tape from the VCR.



**NOTE:** The EJECT function cannot be remote controlled from the LOGIQ™ 500. To eject a video cassette tape from the recorder, press <EJECT> on the VCR front panel.

4. STOP button (remote control only).
5. REW (rewind) button/LED (light emitting diode) (remote control only). The REW LED lights up during rewinding.





## VCR Operating Instructions

---

### Part Names and Functions of the VCR Front Panel (cont'd)

6. PLAY (playback) button/LED (remote control only).  
The PLAY LED lights up during playback.
7. F FWD (fast forward) button/LED (remote control only).  
The F FWD LED lights up during fast-forward operation.
8. PAUSE button/LED (remote control only).  
The PAUSE LED lights up during pause.
9. REC (record) button/LED (remote control only).  
The REC LED lights up during recording.
10. MONITOR (Need not be set).
11. AUDIO (Setting fixed to Hi-Fi).
12. Hi-Fi (Setting fixed to ON).
13. TRACKING control (Setting fixed to FIXED).
14. REMOTE terminal (Disabled).
15. FOOT SWITCH terminal (Disabled).
16. MIC (microphone) terminal (Disabled).
17. REV SEARCH (reverse search) button/LED (remote control only). The REV SEARCH LED lights up during a reverse search.
18. FWD SEARCH (forward search) button/LED (remote control only). The FWD SEARCH LED lights up during a forward search.
19. BLANK SEARCH button (Disabled).
20. FREEZE button (remote control only).  
The FREEZE LED lights up during freezing.
21. PROGRAM switch (Setting fixed to OFF).
22. MARK IN A, B buttons (Disabled).



**NOTE:** *If the setting of a switch is fixed, do not change it.*

*Functions in 4, 5, 6, 7, 8, 9, 17, 18, and 20 above are the only functions available by remote control from the LOGIQ™ 500.*

Part Names and Functions of the VCR Front Panel (cont'd)




Figure 19–2. VCR Front Panel Indicators

23. Indicator block

| Indicators                                                                          | Description                                                                                                                     |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------|
|  | Lights up when a video cassette tape is in the VCR.                                                                             |
| AUTO OFF                                                                            | Lights up if there is condensation in the VCR. If this indicator lights up, no cassette tape can be inserted into the recorder. |
| S-VHS                                                                               | Lights up when a video cassette tape is recorded or played back in the S-VHS mode.                                              |
| MARK                                                                                | Lights up during an index marker recording.                                                                                     |

24. TAPE REMAIN (remaining volume of tape) indicator. Indicates the percentage of tape volume remaining. If almost 100% of a video cassette tape is recordable, all five indicators will light up.

25. CTL/MENU switch (Setting fixed to CTL).

26. INDEX LED.  
Lights up during a patient image search or index search.

27. Time counter display. Indicates tape–running time in hours, minutes, and seconds.

28. RESET button (Disabled).



**NOTE:** If the setting of a switch is fixed, do not change it.

### Part Names and Functions of the VCR Rear Panel




Figure 19–3. VCR Rear Panel

1. AUDIO IN (audio input) CH–1/L and CH–2/R terminals. Connect these terminals to the Video Out Audio terminals on the LOGIQ™ 500.
2. VIDEO IN (video input) terminals. Connect these terminals to the Video Out S–Video terminals on the LOGIQ™ 500.
3. 75–W terminator switch (Setting fixed to ON).
4. AUDIO OUT (audio output) CH–1/L and CH–2/R terminals. Connect these terminals to the Video In Audio terminals on the LOGIQ™ 500.
5. MONITOR OUT (monitor output) terminal (Disabled).
6. VIDEO OUT (video output) terminals. Connect these terminals to the Video In S–Video terminals on the LOGIQ™ 500.

**Part Names and Functions of the VCR Rear Panel** (cont'd)

7. AC IN (AC power input) terminal. Connect this power inlet to a power outlet for peripherals on the rear panel of the LOGIQ™ 500 using the cable supplied with the VCR.
8. Ground terminal (Disabled).
9. INDICATE (indicator) terminal (Disabled).




Figure 19-4. VCR Rear Panel (RS-232C Terminal & Dip Switch)

10. RS-232C (remote control) terminal. Connect this terminal to the VCR/Color Printer Port A or B on the LOGIQ™ 500 rear panel. The proper port must be selected in the Set Up/System Parameters page 5 menu.
11. DIP switch (Settings fixed). When the VCR is connected to the LOGIQ™ 500 for use, do not change the settings of the DIP switch.

Settings:

- |      |            |                     |
|------|------------|---------------------|
| SW1: | OFF (down) |                     |
| SW2: | OFF (down) |                     |
| SW3: | OFF (down) |                     |
| SW4: | OFF (down) | Remote/local switch |
| SW5: | ON (up)    |                     |
| SW6: | OFF (down) |                     |



**NOTE:** If the VCR is disconnected from the LOGIQ™ 500, set SW4 to ON (up).



## VCR Operating Instructions

### Connecting the VCR to the LOGIQ™ 500

Connect the cables from the VCR to the LOGIQ™ 500 by referring to the VCR Part Names and Functions/Settings in the previous section and the following connection diagram.




Figure 19-5. VCR Connection Diagram

## Setting up the LOGIQ™ 500

After the VCR is connected to the LOGIQ™ 500, the LOGIQ™ 500 must also be set up.

Select the top menu Set Up from the software menu and select System Parameters.

On page 5 of the System Parameters menu, set Port A or Port B to VCR.

Then select SAVE to exit the menu.



**NOTE:** After changing the setting, turn OFF the system power switch, then restart the system.

For reference, see *Customizing Your System* in Volume 2 of the LOGIQ™ 500 User Manual.

```

16/06/98 17:28:25
739L
USER PARAM  SYSTEM PARAM  PRESET PROG  CUSTOM DISP
PAGE       : 5/ 7        PRIOR   NEXT   CONTENTS
COMMAND    : SAVE      RESET   DELETE  EXIT

***** System Setup (Recording) *****
Report Video Inverse to Printer : Off On
Record1 B&W : Image Archive
Record1 Color: Image Archive
Record2 B&W : Color Printer
Record2 Color: Color Printer
Report      : Black and White Printer
External Video Signal: S-Video RGB Col-Camp
Color Printer Type   : UP3030 UP1850 UP1850P
Color Printer Memory : Single Quad-Frm Quad-Frms
Color Printer Signal : RGB Video S-Video
Plug1-BW             : B&W-Prnt B&W-Pala Col-Pala
Plug2-CLR            : B&W-Prnt B&W-Pala Col-Pala
Port MIC              : Off
Port A                : Off
Port B                : Off
Maskline Record      : Off On
B&W Printer Exposure Pulse Length [msec] : 100
B&W Printer Exposure Min. Interval [sec] : 6.0
B&W Image Archive     : B&W Color
Color Image Archive   : B&W Color
'TRACKBALL' and 'SET' are available.
  
```

Figure 19–6. Set Up/System Parameters Menu Page 5



## VCR Operating Instructions

---

### Starting the VCR

#### VCR Start-Up



*NOTE: If an alphanumeric character (such as y and n) needs to be typed in during the VCR operation, the system must not be in the Blue Shift mode.*

*If the **Blue Shift** key is lit up, press **Blue Shift** to cancel the Blue Shift mode.*

After all connections and setups are completed, the VCR can be started and will be recognized. The operation flow from the moment the LOGIQ™ 500 Power/Stand-by switch is turned on to the moment the VCR is set is summarized as follows:

#### Checking Switches



Check that the LOGIQ™ 500 Power/Stand-by switch is in the Off/Stand-by position and that the VCR Power switch is in the ON (recessed) position.


At this time, the VCR power will still be off.

If the switch is lit, the VCR power cable has been connected to an outlet other than the power outlet for peripherals on the rear panel of the LOGIQ™ 500.



*NOTE: For safety reasons, always connect the VCR power cable to the power outlet for peripherals on the rear panel of the LOGIQ™ 500.*

## Power On



Turn the LOGIQ™ 500 system power on.

As soon as the system power is turned on, power will be fed simultaneously to the VCR which will light up the VCR Power switch.


Check to ensure the power switch on the VCR panel is lit. If not, check the following in order:

1. Check that the Power switch is in the ON (recessed) position.
2. Check that the power cable is properly connected.



*NOTE: If power is still not being fed to the VCR after checking 1 and 2, the VCR may be defective. Contact an authorized service personnel.*

## Status Icon Display



When the VCR has been correctly connected and set up, the VCR status icon EJECT (▲) will appear at the bottom left corner of the monitor upon completion of start-up.

If the VCR status icon has not appeared, check that the port setting in the LOGIQ™ 500 System Parameters menu page 5 has been set to VCR.

If the stop status icon (■) has appeared instead of the VCR status icon EJECT, a video tape is in the VCR. So that the system will recognize the tape ID, eject the tape from the recorder.






## VCR Operating Instructions

---

### Status Icon Display (cont'd)




The Ext Video key on the control panel will be partially lit.

If Ext Video is fully lit, the VCR video mode has been selected. Press **Ext Video** to switch to the scan image mode.

If Ext Video is not lit, the video tape or VCR has not been recognized. If a video tape has already been inserted, eject it from the recorder and check the tape.

When using a new Tape, proceed to Registering a New Tape on page 19–25.

### Tape insertion



Insert a registered video tape into the VCR.

When a tape has been inserted, the VCR automatically rewinds it up to the tape header so that the system can read and collate the tape ID.

When rewinding is completed, the VCR will play it back to read the tape ID. When the tape ID has been verified, the following message appears at the bottom of the monitor.

System ID = #####, Tape ID = #####



After confirming the ID, the VCR automatically fast forwards the tape to the position of the tape when it was first inserted.

*NOTE: The time required for ID recognition and collation will differ depending on the tape position upon insertion. It usually takes at least one and half minutes. Do not press any key until ID collation is completed.*



## *VCR Operating Instructions*

---

### **Checking Pause/Record**

When collation is completed, the Pause/Record key on the control panel is partially lit.

If the Pause/Record key is not lit, the tape has not been recognized. The message:

“Cannot read Tape ID. Register the new tape?(y/n)”

appears at the bottom of the monitor. The inserted tape may be a new tape or unregistered tape.

If it is a new tape, type in “y” and proceed to Registering a New Tape. See 19–25 for details. In all other cases, type in “n” and try to read the tape again.

### **Start-Up Final Check**

If all the steps have been successfully completed up to this point, the monitor and control panel status will be as follows:

1. Ext Video on the control panel is partially lit.
2. Pause/Record on the control panel is partially lit .
3. The monitor shows the stop VCR status icon (■) and the VCR counter at the bottom left corner of the screen.

If these conditions have not yet been established, verify the procedure in Starting the VCR again.



## VCR Operating Instructions

---

### VCR Status Icon

The VCR status icon will indicate the operation status of the VCR. The icons are displayed at the bottom left corner of the monitor screen.

| VCR STATUS ICON                                                                                  | DESCRIPTION                            |
|--------------------------------------------------------------------------------------------------|----------------------------------------|
|  Tape Eject     | A video tape has not yet been inserted |
|  Stop           | Operation has stopped                  |
|  Record         | Recording in progress                  |
|  Record & Pause | Recording has been paused              |
|  Fast Forward   | Fast forwarding                        |
|  Rewind         | Rewinding                              |

During the display of a VCR playback image, the icon (●) indicating a recorded image appears on the monitor.

### VCR Counter Display

The VCR counter display next to the status icon shows the following information:




Figure 19–7. VCR Counter Display

## VCR Operations

The keys used for remote control of the VCR are as follows:




Figure 19–8. VCR Remote Control (Control Panel/Keyboard)



**NOTE:** A video tape is ejected from the recorder using the **EJECT** key on the front panel of the VCR.




## VCR Operating Instructions

---

### Control Panel

#### Ext Video



#### **Ext Video**

Allows the selection of the video (VCR playback image) mode or scan image mode.

When the video mode is activated, the image displayed on the monitor is switched from scan image to VCR image. The Ext Video key is fully lit in the video mode.

Under the condition in which the video mode can be activated, Ext Video is partially lit. If the video mode cannot be activated, Ext Video is not lit.

To return to the scan image mode, press **Exit Video** again.

Pressing **Exit Video** when the VCR is recording stops the recording.

Pressing **Exit Video** when the VCR is playing back a tape stops the playback and selects the scan image mode.

*NOTE: A video tape must be inserted into or ejected from the recorder in the scan image mode.*

*The Cine function cannot be activated in the video mode.*

#### **Mic**

#### Mic



This is a microphone-activating switch.

Press this key to simultaneously record anything spoken in the vicinity of the console during a recording.

The microphone is located on the bottom right side of the monitor.

*NOTE: Doppler audio will be recorded even if the Mic key is not activated.*



Control Panel (cont'd)

**Stop / Play**



**Stop/Play**

Used to begin VCR playback or stop VCR operation.

Pressing the **Stop/Play** key while the VCR status icon is in a stop condition (■) begins playback.

Pressing **Stop/Play** key while the VCR is operating stops the current VCR operation, changing the VCR status icon to ■.

Stopping playback in the video mode (changing to ■) eliminates the image from the monitor.



*NOTE: In the scan image mode, pressing **Stop/Play** to begin playback while the VCR status icon is ■ selects the external video mode automatically.*

**Control Z**

If the heading function (service software) is turned OFF, the VCR tape counter can be reset by pressing **Control Z**.

If the heading function is turned ON to enable other VCR controls, **Control Z** will NOT reset the counter



## VCR Operating Instructions

---

### Control Panel (cont'd)

#### Pause/Record

##### **Pause / Record**



Used to begin a recording or pause operation.

In the scan image mode, pressing the **Pause/Record** key while the VCR status icon is ■ starts the recording process.

When recording starts, the tape position, patient name, patient ID, and date are recorded in the system as data to be searched.



*NOTE: Recording is not possible in a system other than the system on which the tape was registered. (The Pause/Record key does not light up for any system other than the system on which the tape was registered.)*

- Pressing **Pause/Record** while the VCR is recording pauses the recording.
- Pressing **Pause/Record** when the VCR recording is paused cancels the pause and resumes the recording.



*NOTE: When a recording resumes, the tape position (index) and patient information will be re-recorded in the system.*


*A single tape can record up to 39 indexes. It is recommended that recording not be paused (or stopped) unnecessarily to ensure effective utilization of a video tape.*

Pressing **Pause/Record** while the VCR is operating pauses the recording or playback.

Fast forward or rewind cannot be paused. (Fast-forward search/Reverse search can be paused.)

Pressing **Pause/Record** when the VCR operation is paused cancels the pause and resumes the VCR operation.

Control Panel (cont'd)



Freeze

**Freeze**

Enables or disables the freezing of a VCR playback image.

In the scan image mode, pressing the **Freeze** key stops image acquisition.

To perform VCR playback with measurements, press **Freeze** while a VCR image is being played back or VCR image playback is paused. This freezes the VCR image.

*NOTE: When a VCR image is frozen, Cine playback is not available.*



Blue

**Blue Shift (Blue Shift mode)**

To use the VCR control keys on the keyboard, press the **Blue Shift** key to activate the Blue Shift mode (LED is lit up).

*NOTE: To register a new tape or enter patient information, the Blue Shift mode must be deactivated.*



FF ►►

**FF (Fast forward/fast-forward search)**

Pressing the **FF** key while the VCR is stopped or rewinding will activate fast forward.

A tape may also be fast forwarded in the scan image mode.

To stop fast forward, press **Stop/Play**.

Pressing **FF** while the VCR is playing, reverse searching or paused will activate a fast-forward search.

To cancel the fast-forward search and return to normal playback, press **FF** again.





## VCR Operating Instructions

---

### Control Panel (cont'd)



◀◀ Rev

#### **Rev (Reverse/reverse search)**

Pressing the **Rev** key while the VCR is stopped or fast forwarding will begin reverse.

Tape reverse is also possible in the scan image mode.

To stop reverse, press **Stop/Play**.

Pressing **Rev** while the VCR is playing, fast-forward searching or paused will activate a reverse search.

To cancel the reverse search and return to normal playback, press **Rev** again.

#### **Frame Shift (frame forward/frame reverse)**



◀ 1 1 ▶  
Frame Shift

Pressing a **Frame Shift** key once when VCR playback is paused or fast-forward search/reverse search is paused will begin frame forward/frame reverse.

Holding down **Frame Shift** while VCR playback is paused or fast-forward search/reverse search is paused changes the playback speed of fast-forward search/reverse search.

Releasing the Frame Shift will return the image to the pause status.

#### **Record End Function (Rec-End)**

With the heading function turned ON (service software), the Rec-End function moves the tape to the position at the end of the last recording.


If the heading function is turned OFF, recording can be started at any position on the tape.

## Registering a New Tape

To perform the image search functions described later with the LOGIQ™ 500 system, the video tape to be recorded must be registered in the system in advance.

To register a new video tape:

### Power on




Turn the LOGIQ™ 500 system power on.

When the VCR has been correctly connected and set up, the VCR status icon EJECT (▲) appears at the bottom left corner of the monitor upon completion of system start-up.

If this status icon does not appear, conduct a check in accordance with Starting the VCR on page 19–14.

If the stop status icon (■) has appeared, there is a video tape in the VCR. Press <EJECT> to eject the tape from the recorder.

### Checking the scan mode



Check if the system is in the scan image mode.

Check that the screen displayed on the monitor is a normal live image and that Ext Video on the control panel is partially lit.

If Ext Video is fully lit, press Ext Video to select the scan image mode.


If Ext Video is not lit, conduct a check in accordance with Starting the VCR on page 19–14.



## VCR Operating Instructions

---

### Inserting a tape



Insert a new S-VHS tape into the VCR.

The system automatically begins to play it back and attempts to read tape information to see if the inserted tape has been previously registered.

In the case of a new tape, the system takes about 30 seconds to complete the check and display the following confirmation message, with a prompt to initiate registration procedures.



“Cannot read Tape ID. Register the new tape? (y/n)”

Type in “y” at the new tape registration confirmation message.

This causes the system ID number and the new tape's tape ID number to appear on the screen.

“System ID = #####, New Tape ID = #####”



*NOTE: New registration of a video tape takes about one and half minutes. The monitor displays nothing during this period. Wait until registration is completed.*

### Recording ID

This tape ID is the tape number that will be used for the tape search function. Always write it down on the label to be attached to the tape.



## VCR Operating Instructions

---

### New Tape Registration Complete

Upon completing a new registration, the screen will display the VCR status icon (■) and VCR counter.



*NOTE: The tape header contains the system ID and tape ID registration number recorded. Thus, no data is recorded on the tape header.*

When all the steps up to this point have been successfully completed, the monitor and control panel status is as follows:

1. Ext Video on the control panel is partially lit.
2. Pause/Record on the control panel is partially lit.
3. The monitor shows the stop VCR status icon (■) and VCR counter on the bottom left corner of the screen.

If these conditions have not yet been established, conduct a check following the procedure in Registering a New Tape on page 19–25.



### Recording/Playback/Image Search

In order to record, playback and image search, the Starting the VCR procedures must be completed. The monitor shows the stop VCR status icon (■) and both the Ext Video and Pause/Record keys on the control panel being partially lit.



*NOTE: Recording is available only on the system if the tape has been previously registered on that system.*

### Recording Procedure

#### Registering a new patient

When Ext Video is partially lit (the scan image mode), press the **New Patient** key to enter new patient information.

The image search function, described later, enables tape search/image search based on a patient name/patient ID and scan date. Check the patient name/patient ID before starting each patient examination.

Start the examination and scan the patient. Scanned images are displayed on the monitor.

#### Start recording

Press **Pause/Record** to begin recording.

When recording begins:

- The Pause/Record key is fully lit.
- The VCR status icon changes to recording (●).
- The VCR counter is updated to indicate the recording time.

The patient information and tape position acquired at the beginning of a recording will be stored on the system hard drive.



## VCR Operating Instructions

---

### Start recording (cont'd)

If a paused recording restarts, the patient information and tape position will also be stored in the system hard drive. Therefore, each time the tape is paused/restarted or stopped/restarted, the index heading markers are generated with the current patient information.

Tape searches and image searches are performed based on this information.



*NOTE: A single video tape allows up to 39 index head search markers to be recorded. More than 39 index markers can be recorded, but the image search function can handle only up to 39 index markers.*

### Recording audio

Doppler audio can be recorded in stereo or monaural, while doing Pulsed or Color Doppler studies. The stereo/monaural selection is found in the Set Up/Preset program menu page two.



Mic  
p

If there is a need to add audio to the video tape while recording, press the **Mic** key. (Blue Shift key must be enabled.)

Anything spoken in the vicinity of the console will be recorded. The microphone is located on the bottom right side of the monitor.

### Recording Complete

Press the **Stop/Play** key to complete recording

OR

press the **Pause/Record** key to pause recording.

When recording is stopped, the VCR status icon changes to stop (■). For pause, the status icon changes to paused recording (●■).

In both cases, press **Pause/Record** to restart recording.



## VCR Operating Instructions

---

### Ejecting the tape

When the recording is completed, stop the VCR (status icon changes to ■), then press <EJECT> to eject the tape from the recorder. Attach a label marked with the tape ID to the tape for storage. Rewinding the tape before ejecting will allow quicker ID identification when the tape is reinserted.

### Playback Procedure

#### VCR playback

To playback VCR images, press the **Stop/Play** key.

When playback begins, the system automatically switches to the video mode, displaying a VCR playback image on the monitor.

The Ext Video key is fully lit in the video mode.

If nothing is recorded on the tape and/or during playback of non-recorded areas, the monitor displays nothing. Moreover, pressing **Stop/Play** during playback to stop VCR playback will make the screen display nothing.

#### Audio volume control

Use the **Audio Volume** control, in the Doppler/CFM section of the keyboard, to control the VCR playback volume.



#### Freezing a VCR image

During the playback of VCR images, a VCR image can be frozen by pressing **Freeze**.

When an image is frozen, playback with measurements is available.



## *VCR Operating Instructions*

---

### **VCR playback (cont'd)**

#### **Blue Shift mode**

Pressing the **Blue Shift** key activates the Blue Shift mode.

Remote control of the VCR using the keyboard (FF/Rev/ Frame Shift) is available only in the Blue Shift mode (Blue Shift key is fully lit).

#### **Frame-forward search and variable speed search**

While a playback image is paused, each time Frame Shift is pressed, the paused playback image advances or rewinds one frame.

Holding down Frame Shift changes the search (playback) speed continuously from frame-shift speed to 10 times faster than playback speed.

If Frame Shift is released, the playback image returns to the pause status.

### **Stopping playback**

Press **Ext Video** to return to the scan image mode.

When the system returns to the scan image mode, the monitor displays a scan image and the VCR playback automatically stops.

Ext Video is partially lit.





## Playback with Measurements

Measurements of a frozen VCR playback image are possible.

Press **Freeze** to stop image acquisition of a VCR playback image, then press **Measurement** to activate the measurement function.

The following measurements of playback images are not available.

1. Echo-level measurement
2. CFM-point velocity measurement
3. Biopsy guide measurement
4. TAMAX Auto (Auto trace disabled)

To measure the TAMAX value or to take measurements using the TAMAX value, the TAMAX value must be obtained using manual trace.

Any measurements other than those above, can be taken in the same way as measurements of normal scan images.

Writing comments to a frozen image and display of a body pattern are available.

Measurement results can also be reported in the same way as in normal scan.



*NOTE: When a VCR playback image is frozen, Cine review is not available.*



## *VCR Operating Instructions*

---

### **Precautions for playback with measurements**

For an examination using a scan date (such as determining the growth of a fetus), note that the scan date recorded on a VCR image differs from the current date.

The measurement results of playback images can be output in report form. The patient information and examination category typed in using the New Patient function cannot be modified even when the VCR is started.

Thus, press **New Patient** or **ID/Name** (depending on the patient information/examination category of the VCR playback images) to correct the patient information output in report form.

If the message:

“Please play VCR image and freeze at the other frame”

appears when a measurement is being taken, this indicates that measurement scale data has not been read in that frame. In this case, advance the VCR image a few frames, then freeze the VCR playback image and attempt the measurement again.



## Advanced Search

### Image search

#### Activating Image Search

Pressing **Ctrl** and **I** with Ext Video On enables an image search based on patient information and the scan date.

The search conditions can be entered by typing in a patient name, patient ID, scan date, or all three.

(An attempt to make a search with no information specified for the search conditions, causes all the patient information on the tape to appear in a list of the relevant patients.)



*NOTE: Image search requires the VTR–PB function option.*

*The image search function cannot be activated in the scan image mode.*

#### Entering the search conditions

Use the **Trackball** or press **Return** to move the cursor to either the PT NAME, PT ID or DATE field, depending on the type of search being done.

For a search based on a patient name, enter the patient name using the cursor in the PT NAME field of the image search menu.

For a search based on a patient ID, enter the patient ID with the cursor in the PT ID field of the image search menu.

For a search based on a scan date, enter a scan date with the cursor in the DATE field of the image search menu.

Pressing **Return** when the cursor is to the right of the appropriate field begins a search.



*NOTE: A search based on a combination of patient name, patient ID, and scan date is also possible.*



Image search (cont'd)

**Displaying the relevant patients**

The system searches the data on the system hard drive and displays a list of the relevant patients recorded on the tape.

[ IMAGE SEARCH MENU ]

PT NAME : 12345678901234567890123456789

PT ID : 12345678901234

DATE (MM / DD / YY) : 12 / 31 / 93

NUMBER : ☐

|    | NAME                     | ID | DATE         |
|----|--------------------------|----|--------------|
| 1  | NNNNNNNNNN ... NNNNNNNNN |    | ## / ## / ## |
| 2  | NNNNNNNNNN ... NNNNNNNNN |    | ## / ## / ## |
| 3  | NNNNNNNNNN ... NNNNNNNNN |    | ## / ## / ## |
| 4  | NNNNNNNNNN ... NNNNNNNNN |    | ## / ## / ## |
|    | .....                    |    |              |
|    | .....                    |    |              |
|    | .....                    |    |              |
| 10 | NNNNNNNNNN ... NNNNNNNNN |    | ## / ## / ## |

Ctrl+N : Next Page      Ctrl+P : Previous Page

Ctrl+Q : Quit          Ctrl+C : Cancel

Operator Message Area

Figure 19–9. List of relevant patients Screen



*NOTE: The list displays the relevant patients in the order of the latest scan dates.*

To return from the list screen to the Search Conditions Input Screen, press **Ctrl** and **Q** simultaneously.



## VCR Operating Instructions

---

### Image search (cont'd)

#### **Selecting the relevant patient**

Select the relevant patient from the list, enter the number on the list and press **Return**.

If the number of the relevant patients exceeds 10, press **Ctrl** and **N** to display the next page of the patient list. To return to the previous page of the patient list, press **Ctrl** and **P**.

If no relevant patient exists, the message

“Requested information doesn’t exist on the tape”

appears and the system returns to Image Search Menu.

#### **Selecting the patient**

If the patient is selected by entering the number or if only one patient is present, the system displays the patient name, patient ID, and scan date of the relevant patient. It then begins a search.

#### **Displaying the images**

When the search is completed, the message:

“Search was completed.”

appears and the VCR image playback of the relevant patient begins if a search was activated while the VCR was playing.

## Image search (cont'd)

**Hints**

- In a search using a patient name, upper-case and lower-case letters are not differentiated.
- The system lists all search candidates whose names include the entered word.

Example: Entering Yama as the patient name to conduct a search displays all names containing those letters such as YAMADA, yamamonto, and Nakayama as candidates in the list.

- To exit the search menu, press **Ctrl** and **C** to cancel the search function.
- Pressing **Ctrl** and **Q** returns to the previous entry screen.

If a mistake is made in data entry, press **Back Space** the required number of times to erase any wrong characters, then re-enter the correct characters.



### Tape Search

Pressing **Ctrl** and **T** with Ext Video On enables a tape search to be initiated based on patient information and scan date.

The search conditions are entered by typing in a patient name, patient ID, scan date, or all three.

(An attempt to conduct a search when no information is specified for the search conditions causes all patient information in the hard drive to appear in the list of relevant patients.)



*NOTE: A tape search requires the optional VTR-PB function.*

*The tape search function cannot be activated in the scan image mode.*

#### Entering the search conditions

Use the **Trackball** or press **Return** to move the cursor.

1. For a search based on a patient name, enter the patient name with the cursor in the PT NAME field of the tape search menu.
2. For a search based on a patient ID, enter the patient ID with the cursor in the PT ID field of the tape search menu.
3. For a search based on a scan date, enter a scan date with the cursor in the DATE field of the tape search menu.

Pressing **Return** when the cursor is to the right of the DATE field begins a search.



*NOTE: A search based on a combination of patient name, patient ID, and scan date is also possible.*



Tape Search (cont'd)

**Display of the relevant patients**

The system searches the data on the system hard drive and displays a list of the relevant patients recorded on the tape.

[ TAPE SEARCH MENU ]

PT NAME : 12345678901234567890123456789

PT ID : 12345678901234

DATE (MM / DD / YY) : 12 / 31 / 93

NUMBER :

|    | NAME                       | ID |
|----|----------------------------|----|
| 1  | NNNNNNNNNN . . . NNNNNNNNN |    |
| 2  | NNNNNNNNNN . . . NNNNNNNNN |    |
| 3  | NNNNNNNNNN . . . NNNNNNNNN |    |
| 4  | NNNNNNNNNN . . . NNNNNNNNN |    |
|    | .....                      |    |
|    | .....                      |    |
|    | .....                      |    |
| 10 | NNNNNNNNNN . . . NNNNNNNNN |    |

Ctrl+N : Next Page      Ctrl+P Previous Page

Ctrl+Q : Quit          Ctrl+C Cancel

Operator Message Area

Figure 19–10. List of Relevant Patients Screen



*NOTE: The list displays the relevant patients in the order of the latest scan dates.*

**Re-entering the search conditions**

To return from the list screen to Tape Search Menu, press **Ctrl** and **Q** simultaneously.





**Tape Search (cont'd)**

**Selecting the relevant patient**

Select the relevant patient from the list, enter the number on the list and press **Return**.

If the number of the relevant patients exceeds 10, press **Ctrl** and **N** to display the next page of the patient list. To return to the previous page of the patient list, press **Ctrl** and **P**.

If no relevant patient exists, the message

“Requested information doesn’t exist on the tape”

appears and the system returns to the Tape Search Menu.

**Selecting the patient**

If the patient is selected by entering the number or if only one patient is present, the system displays the patient name/patient ID/scan date of the relevant patient and begins a search.

## Tape Search (cont'd)

### Displaying the search results

When the search is completed, the system displays the message

“Search was completed.”

as well as the system ID and tape ID of the relevant tape.

| [ TAPE SEARCH MENU ]  |                                 |
|-----------------------|---------------------------------|
| PT NAME               | : 12345678901234567890123456789 |
| PT ID                 | : 12345678901234                |
| DATE (MM / DD / YY)   | : 12 / 31 / 93                  |
|                       |                                 |
| SYSTEM ID<br>XXXXX    | TAPE ID<br>XXXX                 |
|                       |                                 |
| Ctrl + Q : Quit       | Ctrl + C : Cancel               |
| Operator Message Area |                                 |

Figure 19–11. Definite Tape ID



## Hints

- In a search using a patient name, upper-case and lower-case letters are not differentiated.
  - The system lists all search candidates whose names include the entered word.  
  
Example: Entering Yama as the patient name to conduct a search displays all names containing those letters such as YAMADA, yamamonto, and Nakayama as candidates in the list.
  - To exit the search menu in mid-course, press **Ctrl** and **C** to cancel the search function.
  - Pressing **Ctrl** and **Q** returns to the previous entry screen.
- If a mistake is made in data entry, press **Back Space** the required number of times to erase any wrong characters, then re-enter the correct characters.



### Saving/Reading the Data to be Searched

To conduct a search for a video tape recorded on another LOGIQ™ 500, the index, patient information, and scan date acquired when a recording is made are stored in the system hard drive whenever recording begins. Also, an index marker is recorded on the tape.

When an image search is initiated, the system conducts a search of data in the system hard drive to determine the index number of the relevant image, and transfers it to the VCR. The VCR searches for the index marker on the tape in accordance with the transferred index number. Thus, an image search requires that the data being searched be stored in the hard disk drive.



*NOTE: To conduct an image search on a system other than the system on which the relevant tape was registered, the data to be searched needs to be copied onto the system in which an image search is conducted.*

In the LOGIQ™ 500 system, data in the system hard drive can be saved on a video tape, then a system on which an image search is to be conducted reads the data on the tape to copy the data to be searched.

### Saving the data to be searched

To save the data in the hard drive on a video tape, press **Ctrl** and **V** simultaneously with Ext Video On. This causes the following message to appear:

“Save patient information on VCR Tape? (y/n)”.

To save the data on the tape, type in “Y”. The system automatically rewinds the tape up to the tape header and records the data in the hard drive onto the tape in bar code. When the data has been saved onto the tape, the VCR stop (changing the status icon to ■).

At this point, press **Ext Video** to place the system in the scan image mode and then eject the tape from the recorder.



## VCR Operating Instructions

---

### Reading the data to be searched

To store the data on a video tape in the system, press **Ctrl** and **S** simultaneously with Ext Video On. This causes the following message to appear:

“Save patient information to System? (y/n)”.

To cancel this function, press “N”

To read the data into the system, press “Y”. The system automatically rewinds the tape up to the tape header, reads the information on the tape and records it on the system hard drive. When saving on the hard drive is complete, the VCR stops (changing the status icon to ■).

At this point, press **Ctrl** and **I** simultaneously to initiate an image search or **Play** to begin playback.



*NOTE: Recording is available only on the system on which the relevant tape was registered.*

*The data to be searched on a video tape cannot be written back into the system on which the tape was registered.*

The original data on the system hard drive cannot be rewritten by rewriting the information (copied) on the tape.

If a read fails or if an attempt is made to read the information from a video tape on which no data to be searched has been saved, the VCR plays for about 2 minutes, then rewind the tape to its initial point before stopping.



## Troubleshooting

### Inspection Sequence

If trouble is suspected, conduct a check in the following order.

1. An operation guide/error message appears at the bottom of the system monitor.  
Refer to Operation Guide/Error Messages on 19–44.
2. Check the items in A List of Problems and How To Trouble Shoot Them on 19–47.
3. Check the VCR switch settings.  
Refer to Setting Up the VCR on 19–11.

If the problem remains after conducting the above, contact a qualified service representative.

### Operation Guide/Error Messages

The table below shows each message and the appropriate action to take. The message column lists the operation guide/error messages displayed at the bottom of the monitor during VCR operation.

| No. | Message                                           | Possible Cause                                                                              | Possible Corrective Action                                                                                                                                                                                                                                                                |
|-----|---------------------------------------------------|---------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| (1) | Cannot read Tape ID. Register the new tape? (y/n) | 1) An unregistered new tape was inserted into the VCR.<br>2) The tape ID could not be read. | Type "y" to register a new tape. Refer to Registering a New Tape on 19–25.<br>Type "n" to proceed to message (08) to re-attempt reading the tape drive.                                                                                                                                   |
| (2) | Check VCR.                                        | There is no signal input or response from the VCR.                                          | 1) Check the VCR power cord and all the cable connections.<br>2) Check the settings of the LOGIQ™ 500 System Parameters menu page 5. Refer to Customizing Your System in LOGIQ™ 500 Advanced Reference Manual.<br>3) Check the VCR switch settings. Refer to Setting Up the VCR on 19–11. |

Table 19–1. Operation Guide/Error Messages

Operation Guide/Error Messages (cont'd)

| No.  | Message                                                    | Possible Cause                                                                          | Possible Corrective Action                                                                                                                                                                                                                                                                                                                                                                                                    |
|------|------------------------------------------------------------|-----------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| (3)  | Check VCR. No Cas-<br>sette.                               | The VCR has no tape.                                                                    | Insert a tape into the VCR.                                                                                                                                                                                                                                                                                                                                                                                                   |
| (4)  | Please play VCR<br>image and freeze at<br>the other frame. | The image information on a tape<br>cannot be read.                                      | Advance the tape a few frames, then<br>freeze the image.                                                                                                                                                                                                                                                                                                                                                                      |
| (5)  | Protected. Check<br>Tape.                                  | This tape has been write pro-<br>tected.                                                | Use an un-protected tape.                                                                                                                                                                                                                                                                                                                                                                                                     |
| (6)  | Requested information<br>doesn't exist on the<br>tape.     | The patient image is not on the<br>tape.                                                | 1) Verify search conditions and re-<br>attempt, if necessary.<br>2) When an inappropriate tape is<br>used, cancel image search and<br>press <b>Ctrl</b> and <b>T</b> to activate a<br>tape search. Find the tape on<br>which the relevant patient de-<br>tected by the tape search has<br>been recorded. Insert the rele-<br>vant tape into the VCR and re-<br>attempt an image search. Refer<br>to Advanced Search on 19–34. |
| (7)  | Requested information<br>doesn't exist in Sys-<br>tem.     | The patient information is not on<br>the hard drive.                                    | 1) Verify search conditions and re-<br>attempt, if necessary.<br>2) When an inappropriate system is<br>used, cancel tape search and re-<br>attempt the search on the correct<br>system. Refer to Tape Search on<br>19–38.                                                                                                                                                                                                     |
| (8)  | Searching Tape ID and<br>System ID                         | After Tape insertion system<br>searches                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                               |
| (9)  | Not registered on this<br>system. Can't record.            | Tape ID and system ID do not<br>match                                                   |                                                                                                                                                                                                                                                                                                                                                                                                                               |
| (10) | Rewinding to begining<br>for Patient Info.                 | Displayed during rewind after #9<br>is displayed                                        |                                                                                                                                                                                                                                                                                                                                                                                                                               |
| (11) | Reading Patient In-<br>formation                           | System reading info from begin-<br>ing of tape                                          |                                                                                                                                                                                                                                                                                                                                                                                                                               |
| (12) | Returning to tape inser-<br>tion position                  | Tape not recognized so system<br>fast forwards to tape position<br>when it was inserted |                                                                                                                                                                                                                                                                                                                                                                                                                               |

Table 19–1. Operation Guide/Error Messages (cont'd)



## VCR Operating Instructions

---

### Operation Guide/Error Messages (cont'd)

| No.  | Message                                     | Possible Cause                                                                                                                                               | Possible Corrective Action                                                                                                                         |
|------|---------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| (13) | Recording is not available                  | Tape is not registered to the system and all attempts to check have been completed                                                                           | Use new tape. Find proper system tape.                                                                                                             |
| (14) | Re-attempt to read Tape ID? (y/n)           | The tape ID cannot be read.                                                                                                                                  | y) Type "y" to re-attempt reading the tape ID.<br>n) Type "n" to cancel reading the tape ID.                                                       |
| (15) | Save patient information on VCR Tape? (y/n) | The operator pressed <b>Ctrl</b> and <b>V</b> to save patient information on the VCR tape.                                                                   | y) Type "y" to save information on the tape.<br>n) Type "n" to cancel the save.<br>Refer to Saving/Reading the Data to be Searched on 19–42.       |
| (16) | Save patient information on System? (y/n)   | The operator pressed Control and S to save patient information on the hard drive.                                                                            | y) Type "y" to save information on the hard drive.<br>n) Type "n" to cancel the save.<br>Refer to Saving/Reading the Data to be Searched on 19–42. |
| (17) | Searching                                   | The VCR is searching the tape. Wait for the search to be completed.                                                                                          |                                                                                                                                                    |
| (18) | Search was completed.                       | The VCR search is completed.                                                                                                                                 |                                                                                                                                                    |
| (19) | System ID = #####, Tape ID = #####          | The tape was mounted.                                                                                                                                        |                                                                                                                                                    |
| (20) | System ID = #####, New Tape ID = #####      | A new tape was registered.                                                                                                                                   | NOTE: Enter the tape ID on the label to be attached to the tape.                                                                                   |
| (21) | This function is not available.             | The function selected is not available for some reason (e.g., a function disabling another function from being activated simultaneously has been initiated). |                                                                                                                                                    |

Table 19–1. Operation Guide/Error Messages (cont'd)

## VCR Operating Instructions

### A List of Problems and How To Troubleshoot Them



| Problem                                                                                                                        | Possible Cause                                                         | Possible Corrective Action                                                                                                                                                                                                                                                                 |
|--------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>The screen displays nothing.</p>           | <p>Unregistered area is being played back in the video image mode.</p> | <p>If a tape advances to an unregistered area as a result of playback or fast-forward search, nothing will appear on the screen and the VCR counter will stop running.</p> <p>Even under this condition, the PLAY or FWD LED on the VCR is lit. Rewind the tape up to a recorded area.</p> |
| <p>The screen displays nothing.</p>           | <p>The VCR has stopped in the video image mode.</p>                    | <p>Check Ext Video. If it is fully lit, the system is in the video mode. Press Ext Video to select the scan image mode, or Press Play to begin VCR playback.</p>                                                                                                                           |
| <p>VCR status icon doesn't appear.</p>       | <p>The VCR has not been recognized.</p>                                | <p>Check the VCR connections and LOGIQ™ 500 settings, then restart the system following the procedure in Starting the VCR on 19–14.</p>                                                                                                                                                    |
| <p>The VCR status icon doesn't appear.</p>  | <p>New tape is being registered.</p>                                   | <p>Registration of a new tape takes about one and half minutes. During this time, the correct status icon doesn't appear.</p> <p>Refer to Registering a New Tape on 19–25.</p>                                                                                                             |

Table 19–2. Operation Guide/Troubleshooting Problems





## VCR Operating Instructions

### A List of Problems and How To Troubleshoot Them (cont'd)



| Problem                                                                                                                                             | Possible Cause                                            | Possible Corrective Action                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-----------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>The VCR status icon doesn't appear.</p>                         | <p>System ID and tape ID are being checked/collated.</p>  | <p>When a tape is inserted into the VCR, the system automatically re-winds it to the tape header to conduct tape ID check/collation. After checking/collating, the VCR automatically fast forwards it to the initial tape position when the tape was inserted. During this period, the correct status icon doesn't appear on the screen. Wait for the VCR to stop. When ID checking/collation has been completed, the message</p> <p>System ID = #####,<br/>Tape ID = #####<br/>appears.</p> |
| <p>ID collation is complete, but Ext Video doesn't light up.</p>  | <p>A VCR-related key was pressed during ID collation.</p> | <p>Eject the tape from the VCR, turn the system power off, and re-start the system.</p>                                                                                                                                                                                                                                                                                                                                                                                                      |
| <p>Tape cannot be inserted into the VCR.</p>                     | <p>The VCR already has a tape.</p>                        | <p>Check the VCR for a tape.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <p>Tape cannot be inserted into the VCR.</p>                     | <p>Condensation has occurred in the VCR.</p>              | <p>If the VCR indicator AUTO OFF is lit, condensation has occurred in the VCR. Wait for the AUTO OFF LED to go off.</p>                                                                                                                                                                                                                                                                                                                                                                      |

Table 19–2. Operation Guide/Troubleshooting Problems (cont'd)

If the system has a problem or suffers a failure, contact a qualified service representative or a sales representative for assistance.

## Helpful hints



### Hints

The following hints can help when using the VCR:

- S-VHS video tapes must be used to take full advantage of the S-VHS features.
- Images taped on an S-VHS unit **must** be played back on an S-VHS unit.
- VHS recorded images can be played on an S-VHS unit.
- S-VHS recorded tapes will **not** playback on standard VHS recorders.
- VHS tapes **can** be used in an S-VHS unit; however, they will **not** work in S-VHS Mode.
- Break out tab to disable **Record**.
- Tape over hole to record again.
- For keyboard entries like "Ctrl, I" or to answer "y/n" questions, the Blue Shift key must be off.
- To use the VCR control keys, the Blue Shift key must be on.

## CAUTION



The system can keep track of patient and tape information. The system can search a tape for patient images. However, the system CANNOT:

- Stop or indicate when it comes to the end of the last study on the tape.
- Prevent recording over previous studies if record is pressed while positioned in the middle of a study.



## *VCR Operating Instructions*

---

This page left blank intentionally.

# Storage/Print Option

## Overview

Digital Imaging and COmmunications in Medicine (DICOM) provides an interface between the LOGIQ™ 500 and imaging/recording devices on a network.

DICOM implementation changes how peripherals are set up. This section explains the use of the LOGIQ™ 500 with a DICOM network.

## DICOM Presets

Presets relating to DICOM operation can be found in Set Up/System Parameters page 7. Refer to 14–63 for details.

## Menu Access

| Archive         | DICOM | Auto Seq | CINE |
|-----------------|-------|----------|------|
| Printer Setup   |       |          |      |
| Network Utility |       |          |      |

Figure 20–1. DICOM Sub-Menu Page 2

To access the Network Configuration Host Verification and Image Transferring Queue Status menus described in the following pages, select Network Utility from the DICOM sub-menu. For DICOM printer setup, select Printer Setup from the DICOM sub-menu.



Patient Entry Menu

“Accession No:” is an entry added to each Exam Category Patient Entry Menu screen for DICOM. It is used to track all studies for a patient during a single admission or visit.

Network Configuration

The Network Configuration Menu is used to set up network name, ID addresses for the LOGIQ™ 500 router(s), storage device(s) and printer(s).

Before filling in the Network Configuration Menu, the following information should be obtained from the Network System Administrator.

The form is titled "NETWORK CONFIGURATION" and is divided into several sections. The first section contains fields for "LOCAL HOST NAME", "LOCAL IP ADDRESS", "NETMASK", and "AE TITLE", with a "PORT" field to the right. The second section is titled "ROUTING INFORMATION" and contains fields for "DESTINATION" and "GATEWAY", with a "DEFAULT" field below "DESTINATION". The third section is titled "DICOM APPLICATION INFORMATION" and contains a table with columns for "SERVICE", "DEVICE", "AE TITLE", "IP ADDR", and "PORT". The table has 7 rows, with the first row labeled "1." through "7." and the last row labeled "SCHEDULE". Below the table are three buttons: "VERIFY MENU", "QUEUE MENU", and "EXIT". At the bottom of the form is a "Message Area" with a hatched background.

Figure 20–2. Network Configuration Worksheet



**NOTE:** It is a good idea to keep this network information available for future reference.



## Storage/Print Option

---

### LOGIQ™ 500 Information

- Local Host Name: This is the network's name for the LOGIQ™ 500. Maximum of seven characters.
- Local IP Address: This is the LOGIQ™ 500's Internet Protocol (IP) address. Every device on the network has a unique IP address.
- Netmask: This is an IP address filter that eliminates communication/messages from the network devices that are of no interest to the LOGIQ™ 500.
- AE Title: This is the name assigned to the LOGIQ™ 500 by the Network System Administrator. A name that a remote DICOM application will recognize. Application refers to a service provided by the network like verify, store, print, etc.
- Port: This is the LOGIQ™ 500 port number (0–99999).

A diagram of a network configuration form titled "[NETWORK CONFIGURATION]". It contains five rows of input fields: "LOCAL HOST NAME" (7 characters), "LOCAL IP ADDRESS" (4 octets), "NETMASK" (4 octets), "AE TITLE" (16 characters), and "PORT" (5 digits).

Figure 20–3. LOGIQ™ 500 Information

### Routing Information

- Default Gateway IP Address: This is the IP address of a gateway (doorway) so that data can be moved from one network to another.
- Router IP Address: Used if applicable. Most networks will not have a router.

A diagram of a routing information form titled "[ROUTING INFORMATION]". It contains two main sections: "DESTINATION" and "GATEWAY". Each section has a "DEFAULT" row and a "ROUTER" row. The "DESTINATION" section has a "NET/HOST" field with a rocker switch. The "GATEWAY" section has a "NET/HOST" field with a rocker switch. The form is divided into four numbered regions: 1 (Default Gateway IP Address), 2 (Destination: DICOM Server IP Address), 3 (Gateway: Gateway IP Address of DICOM server), and 4 (NET/HOST: Select "NET" or "Host" by the Ellipse rocker switch).

- 1 Default Gateway IP Address
- 2 Destination: DICOM Server IP Address.
- 3 Gateway: Gateway IP Address of DICOM server (listed in 1).
- 4 NET/HOST: Select "NET" or "Host" by the Ellipse rocker switch

Figure 20–4. Routing Information



DICOM Application Information

- Service: This defines the service provided by the host device. The services available are STORE–SC, STORE–US, PRINT and NULL (blank). They can be selected by the **Ellipse** rocker switch.
- Device: The device(s) network name as provided by the Network System Administrator. Maximum seven characters.
- The sending image type: select G = Grayscale or C = Color.
- AE Title: The device(s) application entity title, as provided by the Network System Administrator.
- IP Addr: The device(s) IP address as provided by the Network System Administrator.
- Port: The device(s) port information as provided by the Network System Administrator (0–99999).

| DICOM APPLICATION INFORMATION |        |  |          |  |         |      |
|-------------------------------|--------|--|----------|--|---------|------|
| SERVICE                       | DEVICE |  | AE TITLE |  | IP ADDR | PORT |
| 1.                            |        |  |          |  |         |      |
| 2.                            |        |  |          |  |         |      |
| 3.                            |        |  |          |  |         |      |
| 4.                            |        |  |          |  |         |      |
| 5.                            |        |  |          |  |         |      |
| 6.                            |        |  |          |  |         |      |
| 7.                            |        |  |          |  |         |      |
| – SCHEDULE                    |        |  |          |  |         |      |

Figure 20–5. DICOM Application Information

Menu Navigation

- Verify Menu: Displays the Host Verification Menu.
- Queue Menu: Displays the Image Transfer Queue Status Menu.
- Exit: Exit DICOM Network Utility.

It is important to note that modifications to the Network Configuration Menu will not be effective until the LOGIQ™ 500 is rebooted. (Turn the power off and back on again). Existing parameter will continue to be used until the system is rebooted.



## Printer Setup

The DICOM print function supports a maximum of 16 image boxes (4x4) on a single sheet of film.

| [ PRINTER SETUP MENU ]    |                                                 |           |
|---------------------------|-------------------------------------------------|-----------|
| NETWORK PRINTER           |                                                 |           |
| TEST                      | 1x1 LAND REP DFT DFT DFT DFT 8x10IN DFT MED COL |           |
| PAGE SETUP                | FORMAT[ColxRow]                                 | 1x1       |
|                           | ORIENTATION                                     | LANDSCAPE |
|                           | MAGNIFICATION                                   | REPLICATE |
|                           | EMPTY                                           | DEFAULT   |
|                           | BORDER                                          | DEFAULT   |
|                           | TRIM                                            | DEFAULT   |
|                           | MIN DENSITY                                     | 10        |
|                           | MAX DENSITY                                     | 340       |
| MEDIUM SETUP              | TYPE                                            | DEFAULT   |
|                           | PAGE SIZE                                       | 81NX101N  |
|                           | DESTINATION                                     | DEFAULT   |
| PRIORITY                  |                                                 | MED       |
| CONFIGURATION INFORMATION |                                                 |           |
| TIME OUT [MIN.]           | 0                                               |           |
|                           |                                                 | [ EXIT ]  |

Figure 20–6. Printer Setup Menu



## Printer Setup (cont'd)

All DICOM printer specific configurations are set through the Printer Setup Menu. The following parameters are selectable via the Printer Setup Menu:

| Parameter                        | Values                                                                        |
|----------------------------------|-------------------------------------------------------------------------------|
| <b>Page Setup</b>                |                                                                               |
| Format                           | 1x1, 2x2, 2x3, 3x2, 3x3, 3x4, 4x3, 3x5, 5x3, 4x4                              |
| Orientation                      | Landscape/Portrait                                                            |
| Magnification                    | Replicate/Bilinear/Cubic/None                                                 |
| Empty                            | Black/White / 50% /Default                                                    |
| Border                           | White/Black / 50% / Default                                                   |
| Trim                             | Yes / No / Default                                                            |
| Min Density                      | 0–500                                                                         |
| Max Density                      | 0–500                                                                         |
| <b>Medium Setup</b>              |                                                                               |
| Type                             | Paper/Blue film/Clear/Film/Default                                            |
| Page Size                        | 8"x10", 10"x12", 10"x14", 11"x14", 14"x14", 14"x17", 24cm x 24cm, 24cm x 30cm |
| Destination                      | Magazine/Processor/Default                                                    |
| <b>Priority</b>                  | Low / Med / High                                                              |
| <b>Configuration Information</b> | Type A DICOM Printer Device Information                                       |
| <b>Time Out</b>                  | Default = 0 minutes (0–99)                                                    |

The printer AE Title entered in the Network Configuration menu will be displayed on the Printer Setup Menu. Maximum number of printers that can be configured is seven (7).

A Force Print selection located on the menu enables the user to print the job to the printer. With the Force Print selection, the user has no control over the layout of the images on the film. The images collected to that point in time will be printed.



## Hints

The printer setup parameters cannot be changed if a job is pending in the queue for that device.



## Host Verification

When the network is configured properly in the Network Configuration Menu, the DICOM (Network) function will be available after the system is turned on.

Before sending an image, it is recommended that communication between the LOGIQ™ 500 and other host devices be verified. The Host Verification Menu provides this function.

The Host Verification Menu screen is divided into several sections. At the top, the title '[HOST VERIFICATION MENU]' is centered. Below it are four input fields: 'LOCAL HOST NAME', 'LOCAL IP ADDRESS', 'NETMASK', and 'AE TITLE'. In the center, there is a table titled '[HOST INFORMATION]' with four columns: 'DEVICE', 'AE TITLE', 'VERIFY?', and 'RESULT'. The table contains eight rows, with the first row showing '1.' in the device column, 'Y' in the verify column, and 'SUCCESSFUL' in the result column. Below the table is a 'PING' button. At the bottom, there are four buttons: '[EXECUTE]', '[NETCONF.MENU]', '[QUEUE.MENU]', and '[EXIT]'. A 'Message Area' is located at the very bottom of the screen.

| DEVICE | AE TITLE | VERIFY? | RESULT     |
|--------|----------|---------|------------|
| 1.     |          | Y       | SUCCESSFUL |
| 2.     |          | Y       | FAILED     |
| 3.     |          | N       | SUCCESSFUL |
| 4.     |          |         |            |
| 5.     |          |         |            |
| 6.     |          |         |            |
| 7.     |          |         |            |
| 8.     |          |         |            |

Figure 20–7. Host Verification Menu

Most information in the Host Verification Menu comes from the Network Configuration Menu. The only field that can be edited is “Verify?” under [Host Information].

Input a “Y” into this verify field to execute a DICOM verify or an “N” or “No Input” to ignore a specific device.

Trackball to the [EXECUTE] command and press **Set**. The result field will display the results of the DICOM verify as “Successful” or “Failed”.



## **Host Verification (cont'd)**

When the result is successful, it shows that DICOM communication is available between the LOGIQ™ 500 and host server. If the verify failed, the LOGIQ™ 500 cannot communicate with the remote server.

Failure could be due to many reasons. First, check that the Network Configuration Menu parameters are correct. Second, check with the Network System Administrator on the status of other network devices.

## **Ping**

Ping can be used to verify an IP Address on the network. Enter the desired IP Address and set Verify to 'y' (yes). Trackball to the [EXECUTE] command and press **Set**.

The result, sending/receiving packet number and average speed of sending/receiving packet success, will be displayed.

## **Menu Navigation**

- Netconf Menu: Display the Network Configuration Menu.
- Queue Menu: Display the Image Transfer Queue Status Menu.
- Exit: Exit the Network Utility Menu.



## Image Transfer

### CAUTION



If special characters are entered in the New Patient Menu (using Blue Shift or Red Shift keys), the image will not transfer to the DICOM device.

The Service and Device Name of each DICOM Server (1–7) that is registered in the Network Configuration Menu will appear on the VFD display when the DICOM Top Menu is selected.




Figure 20–8. DICOM Top Menu

To transfer the image displayed on the screen, press either the top or bottom of the appropriate rocker switch under the DICOM server selections.

During the transfer process, the image is first stored to the hard disk drive (HDD) as a temporary file. Until this operation is complete, keyboard operation is prohibited.

A message is displayed at the bottom of the screen (in the case of 3x5 format in the Printer Setup Menu):

“Print Queued Image Number: 1/15”

After the temporary fifteen are stored on the HDD, DICOM transfer is executed as a background process when all images are collected for that job.

## Image Transfer (cont'd)

If the "Auto Deletion of Transferred Queue" setting on System Parameters Page 7 is Yes, when the DICOM image transfer is complete, the temporary file is deleted from the HDD. At the same time, the message "Image Transferring is complete" will be displayed. DICOM image transfer is successful.



**NOTE:** If the system does not have enough information related to the DICOM server selected, an error message will be displayed and the image will not be sent.



### Hints

The print job will be sent to the printer if:

1. A new patient is entered on the New Patient Menu.
2. A time out has occurred.
3. Power is switched off.

## Image Transfer Status

| IMAGE TRANSFERRING QUEUE STATUS MENU |                 |                 |      |              |  |
|--------------------------------------|-----------------|-----------------|------|--------------|--|
| RETRY COUNT :                        |                 |                 |      | PAGE : 1 / 1 |  |
| RETRY INTERVAL TIME (MIN) :          |                 |                 |      |              |  |
| PT ID                                | IMAGE DATA/TIME | DEVICE          | TYPE | STATUS       |  |
| 1.                                   |                 |                 |      |              |  |
| 2.                                   |                 |                 |      |              |  |
| 3.                                   |                 |                 |      |              |  |
| 4.                                   |                 |                 |      |              |  |
| 5.                                   |                 |                 |      |              |  |
| [ RETRY ]                            |                 | [ DELETE ]      |      |              |  |
| [ NETCONF MENU ]                     |                 | [ VERIFY MENU ] |      | [ EXIT ]     |  |
| Message Area                         |                 |                 |      |              |  |

Figure 20–9. Image Transferring Queue Status Menu



### Image Transfer Status (cont'd)

The Image Transferring Queue Status Menu shows the jobs that are in progress. It shows Patient ID, Image Data and Time, Device, Type and Status.

It is important to note that each line represents a job and not a single image. With the 4x4 page format, a job could be as many as 16 images.

The Image Transferring Queue Status Menu will read one of four "Status" possibilities:

|       |                                           |
|-------|-------------------------------------------|
| FAIL  | Printer returned a failed status          |
| TRANS | Image transfer is in progress             |
| WAIT  | Print job is in queue, waiting to be sent |
| COMP  | Print job is complete.                    |



*NOTE: If the "Auto Deletion of Transferred Queue" setting on System Parameters Page 7 is No, COMP will be displayed.*

If for some reason the image failed to transfer, the message "Image Transferring is Failed" is displayed.

The temporary file on the Hard Disk Drive (HDD) is not automatically deleted. The system will try to transfer the file again depending on the parameters set in the Image Transferring Queue Status Menu for:

- Retry Count: Number of times to try again (0 to 9999).
- Retry Interval Time (min): The time in minutes between trying again (1 to 9999).



*NOTE: With this Retry function, COMP files will not try to transfer.*



## Image Transfer Status (cont'd)

### Manual Retry

If all automatic retries have failed to transfer the image, the image file can be selected manually by:

- **Trackball** to the file desired and press **Set**.
- **Trackball** to the command [RETRY] and press **Set**.
- The DICOM image transfer attempt will be made again.

The temporary image file remains on the HDD until image transfer is successful or it is manually deleted using the delete command in the Image Transferring Queue Status Menu.

If a job has failed or a print device is inoperative, the user can:

- **Trackball** to the job and press **Set**.
- Use the **Ellipse** rocker switch to change the destination device. At that time, select the same type (color or b/w) device as the previous selection.
- **Trackball** to retry and press **Set**.

The job will be sent to the newly selected device.

### Force Print

Force Print can be used to execute the print function immediately. For example, if the page setup is 3x3 and only 6 images have been printed, Force Print can be used to print the page with only 6 images instead of 9.

### Power Down

When the system power is turned off, the image transfer that is in progress will be completed before shut down. The following message is displayed.

“Execute all dicom jobs before shut down. Y/N?”

If yes, all queues that have a “wait” status will be executed before shut down. If no or a time out occurs, the system shuts down immediately. The status of all queues will be changed to “fail” and remain on the Queue Menu.



## Storage/Print Option

---

### Image Transfer Status (cont'd)

#### Maximum Temporary

##### Files

The maximum number of temporary images that can be stored on the HDD for DICOM transfer is 200.



*NOTE: It is important to remember that each line in the menu represents a job and not an image. With the 4x4 page format, a job could contain as many as 16 images.*

If there are 200 temporary images stored on the HDD and another DICOM image transfer is attempted, the message “The image storage request que is full” is displayed. No temporary file is written on the HDD and the transfer process will not be executed.

Images in the queue must be deleted before additional DICOM transfer attempts can be made.



# Worklist Option

## Overview

The DICOM Worklist or Schedule function allows the user to select a patient from the supplied list and have the vital patient information (Patient Name [18 characters], Patient ID [10 characters], Date of Birth and Reference status) displayed.

## New Patient Selection

The preset “Menu Selection at New Patient” can be set to “Patient\_Entry” or “Schedule”.

- Patient Entry—Each time the New Patient key is pressed, the Patient Entry Menu is displayed first.
- Schedule—Each time the New Patient key is pressed, the DICOM Worklist (schedule) is displayed first.

Despite the selected preset, the user can open the Schedule Menu from the New Patient Menu or the New Patient Menu from the Schedule Menu.

## Data Transferred

Data transferred from the schedule server to the Patient Menu will be:

- Patient Name
- Patient ID
- Accession Number
- Date of Birth

The schedule Menu shows a patient data list that includes:

Procedure start Date/Time, Patient Name, Patient ID, Date of Birth and reference Status.

The Schedule Menu is displayed as shown in Figure 20–10.

```

[ SCHEDULE ]
UPDATED AT : PAGE
DATE TIME PATIENT NAME DOB PATIENT ID $

PATIENT MENU SORT NAME DAYS FORWARD : 7 (0-7)
REF STATUS SORT DATE DAYS BACKWARD : 7 (0-7)
EXIT DELETE WLST QUERY INTERVAL : 0 MIN
DELETE ALL
UPDATE WLST

```

Figure 20–10 Worklist (Schedule Menu)

## Change Pages

In order to move forward or backward through the available pages of the worklist, move the cursor to either page number in the upper right corner and press **Set**.



## Command Selections

**PATIENT MENU:** Highlight a patient in the worklist. Press SET to display the Patient Menu input data screen with worklist data added.

**REF STATUS:** Change the status of the patient to “R” Referenced or “N” Not Referenced.

**EXIT:** Exit schedule menu. Return to scanning.

**SORT NAME:** Sort the worklist by Patient Name.

**SORT DATE:** Sort the worklist by Date.

**DELETE WLIST:** Highlight a single patient record and delete it.

**DELETE ALL:** Delete all patient data from the system.

**UPDATE WLIST:** Updates the internal schedule database from a remote server database.

**DAYS FORWARD:** Specify the range for the number of days forward for which patient database information will be displayed.

**DAYS BACKWARD:** Specify the range for the number of days backward (previous) for which patient database information will be displayed.



***NOTE:** If only the current days patient data base information is to be displayed on the worklist, type “0” for the Days Forward and Days Backward settings.*

**QUERY INTERVAL:** Specify the amount of time between automatic database queries and updates.



### *Worklist Option*

---

This page left blank intentionally.

# Index

## Symbols

- % Stenosis, 8–3, 10–23, 11–5
- % Stenosis Calculation [EC, UP], Preset Program menu page 4, 14–72

## Numbers

- 3D Aspect Ratio [EC, UP, P, R], Custom Display menu page 9, 14–29
- 3D B Background Level (%) [EC, UP], Custom Display menu page 18, 14–45
- 3D B Frame Line Level (%) [EC, UP], Custom Display menu page 18, 14–45
- 3D B Gradient Threshold (%) [EC, UP], Custom Display menu page 18, 14–45
- 3D B Percent of Source Frames [EC, UP], Custom Display menu page 18, 14–45
- 3D B Projection Technique [EC, UP, R], Custom Display menu page 18, 14–45
- 3D CFM/PDI Frame Line Level (%) [EC, UP], Custom Display menu page 18, 14–46
- 3D CFM/PDI Gradient Threshold (%) [EC, UP], Custom Display menu page 18, 14–46
- 3D CFM/PDI Percent of Source Frames (%) [EC, UP], Custom Display menu page 18, 14–46
- 3D CFM/PDI Projection Technique [EC, UP, R], Custom Display menu page 18, 14–46
- 3D CFM/PDI Shade Boundary Threshold (%) [EC, UP], Custom Display menu page 18, 14–46
- 3D CFM/PDI Shade Object Threshold (%) [EC, UP], Custom Display menu page 18, 14–46
- 3D Mode
  - B–Mode, 5–20
  - Color Flow, 5–39
- 3D Number of Frames in Build [EC, UP], Custom Display menu page 18, 14–45
- 3D Scan Line Resolution [EC, UP], Custom Display menu page 18, 14–45
- 3D–Surface Mode, 5–79
  - Overview, 5–79
  - Standard Procedure, 5–80
- 3D Swing Angle (deg) [EC, UP, P, R], Custom Display menu page 9, 14–29
- 3D Swing Center (%) [EC, UP, P, R], Custom Display menu page 9, 14–29
- 3D Swing Speed [EC, UP, R], Custom Display menu page 18, 14–45
- 3D View [EC, UP, P, R], Custom Display menu page 9, 14–29
- 3DvieW Mode, 5–70
  - Activating, 5–72
  - CFM Map, 5–72
  - Functionality, 5–71
  - MIP Image Rendering, 5–75
  - Mode Changes, 5–71
  - Overview, 5–70
  - Techniques, 5–77
- 546L, 15–24
- 739L, 15–25

## A

- A/B Ratio, 9–23
- Abdomen/Small Parts, Measurements
  - A/B Ratio, 8–3, 10–29
  - Max PG, 8–3
  - Mean PG, 8–3
  - Pulsatility Index (PI), 8–3, 10–29
  - Resistance Index (RI), 8–3, 10–29
  - S/D Ratio, 8–3, 10–29
- Accessories
  - Connector Panel, 3–7
  - Placing an order, 16–38
  - Requesting a catalog, 16–38
- Acclimation Time, 3–16
- ACE
  - CFM/PDI Enhancement Option, 5–38
  - Color Flow, 5–38



## Index

---

- Acoustic Output, 5–11, 17–1
  - 546L, 17–23, 17–55, 17–56, 17–66
  - 739L, 17–24, 17–57
  - Adjusting, 5–8
  - B510, 17–18, 17–35, 17–36
  - Best Practices, 2–14
  - C358, 17–19, 17–37
  - C364, 17–19, 17–39
  - C386, 17–20, 17–41
  - C551, 17–20, 17–43
  - C721, 17–21, 17–45
  - CWD2, 17–18, 17–79
  - CWD5, 17–19, 17–80
  - Display, 5–7
  - E721, 17–21, 17–47
  - ERB7 (Convex), 17–22, 17–49
  - ERB7 (Linear), 17–22, 17–51
  - Fetal Exposure, 9–5
  - General Warning, 9–5
  - Hazard, 5–11
  - I739, 17–23, 17–53
  - Key to probe tables, 17–81
  - key to probe tables, 17–85
  - L764, 17–24, 17–59
  - LA39, 17–25, 17–61
  - LD, 17–25, 17–63
  - P509, 17–26, 17–65
  - Precision, 17–82
  - Prudent Use, 9–5
  - S220, 17–26, 17–67
  - S222, 17–27, 17–69
  - S316, 17–27, 17–71
  - S317, 17–28, 17–73
  - S611, 17–28, 17–75
  - Safety, 2–14
  - T739, 17–29, 17–77
  - Tables, 17–18
- Acoustic Output Hazard, 2–5
- Acoustic Power B/M (%) [EC], Custom Display menu page 1, 14–14
- Acoustic Power D/CFM/PDI (%) [EC], Custom Display menu page 1, 14–14
- Acoustic Power Knob Effect [EC], Custom Display menu page 10, 14–30
- Acoustic Zoom, 6–2
- Active Measurement Cursor, System Parameters menu page 3, 14–54
- Add 1 Week to EDD [EC], Preset Program menu page 3, 14–69
- Advanced Cardiac Measurement Option, Sub-Menu, 3–45
- Advanced Vascular, 11–9
  - Menu Selections, 11–9
  - Report Page Layout, 11–13
  - Venous Comments, 11–16
- AFI, OB Measurements, 9–20
- Air Filter
  - Cleaning, 16–18
  - Locating, 16–17
  - Removing, 16–18
- ALARA, 5–11, 17–9
- Amniotic Fluid Index (AFI), 9–20
- Angle
  - Calculation, 8–3, 10–22
  - Correction (Theta), Doppler, 5–51
- Annotation, 6–15
  - Blue Shift, 6–20
  - Clear, 6–16
  - CTRL+M, 6–16
  - Cursor Home, 6–16
  - Editing, 6–21
  - Image, 6–19
  - Library, 6–17
  - Red Shift, 6–20
  - Set, 6–16
  - Shift+TAB, 6–17
  - Special Keys, 6–20
  - Tab, 6–17
- Annotation Library 1–24 [EC, UP], Preset Program pages 8 & 9, 14–90
- Aortic Valve, 10–86, 10–91
- Apical 2 Chamber, 10–83
- Apical 4 Chamber, 10–81
- Application Presets, Defining a User Preset, 14–93
- Archive
  - Disk Verification, 13–27
  - Image Recall Process, 13–32
  - Image Selection for Recall, 13–32
  - Media Format, 13–26
  - Media Search, 13–33
  - MO Eject, 13–35
  - MOD, 13–23
  - Patient Search, 13–28
  - Preset Parameters, 13–24
  - Storing Images, 13–27
- Arrow Mark Display [EC, UP], Preset Program menu page 4, 14–72
- Audio
  - Speakers, 3–26
  - VCR Recording, 19–29
  - Volume Control, 19–30
- Audio Tone Equalizer [EC], Custom Display menu page 10, 14–32
- Audio Volume, Doppler, 5–50
- Auto B Melt at Unfreeze [EC, UP], Custom Display menu page 12, 14–34
- Auto Deletion of Transferred Queue, System Parameters menu page 7, 14–63

Auto Focus Control [EC], Custom Display menu page 10, 14–30

Auto Sequence  
Programming, 10–66  
Sub–Menu, 3–43  
Top Menu, 3–43

Aux  
Gain, 10–36  
Positn, 10–36  
Wave, 10–35

AUX Auto Display [EC], Custom Display menu page 17, 14–44

AUX Gain [EC, R], Custom Display menu page 17, 14–44

AUX Position [EC, R], Custom Display menu page 17, 14–44

Average Activity [EC], Preset Program menu page 3, 14–69

Average Number [EC], Preset Program menu page 3, 14–69

Averaging Number for Doppler Realtime Calc [EC, UP], Preset Program menu page 4, 14–71

## **B**

B Angle/Width (deg. mm) [EC, UP, P, R], Custom Display menu page 2, 14–17

B Auto Caliper Control [EC, UP], Preset Program menu page 4, 14–72

B Color Auto Start [EC, R], Custom Display menu page 11, 14–32

B Color Map [EC, UP], Custom Display menu page 11, 14–32

B Combi Focus Width [EC, UP, P], Custom Display menu page 2, 14–17

B Dynamic Range [EC, UP, P, R], Custom Display menu page 2, 14–15

B Edge Enhance [EC, UP, P, R], Custom Display menu page 2, 14–15

B Focus Combi Number [EC, UP, P, R], Custom Display menu page 2, 14–17

B Focus Dual Marker [EC, UP], Custom Display menu page 11, 14–32

B Format B Dual [EC], Custom Display menu page 11, 14–33

B Format B Single [EC], Custom Display menu page 11, 14–33

B Frame Average [EC, UP, P], Custom Display menu page 2, 14–15

B Gain [EC, UP, P, R], Custom Display menu page 2, 14–15

B Gray Scale Map [EC, UP, P, R], Custom Display menu page 2, 14–16

B Image Frequency [EC, UP, P, R], Custom Display menu page 2, 14–15

B Image Softener [EC, UP, P, R], Custom Display menu page 2, 14–15

B Pause, Doppler, 5–53

B Rejection [EC, UP, P, R], Custom Display menu page 2, 14–17

B Scale Mark [EC], Custom Display menu page 11, 14–33

B Softener Level [EC, UP, P], Custom Display menu page 2, 14–15

B Tag Auto Start [EC, R], Custom Display menu page 11, 14–32

B Tag Center (%) [EC, R], Custom Display menu page 11, 14–33

B Tag Width (%) [EC, R], Custom Display menu page 11, 14–33

B Target Frame Rate [EC, UP, P], Custom Display menu page 2, 14–17

B Video Inverse [EC], Custom Display menu page 11, 14–33

B&W Image Archive, System Parameters menu page 5, 14–61

B&W Printer Exposure Min Interval [sec], System Parameters menu page 5, 14–61

B&W Printer Exposure Pulse Length [msec], System Parameters menu page 5, 14–61

B–Mode  
3D Mode, 5–20  
Acoustic Output, 5–11  
B/M Gain, 5–13  
Biopsy Guidezone(s), 5–22  
Circumference/Area, 7–10  
Color, 5–20  
Color Tag, 5–21  
Colorized, 5–20  
Colorized Gray Scale, Color Tag, 5–21  
Depth, 5–13  
Display, 5–4  
Display Format(Dual), 5–15  
Distance Measurement, 7–7  
Dual Format, 5–15  
Dynamic Range, 5–16  
Echo Level Measurement, 7–11  
Edge Enhance, 5–23  
Focus Number, 5–17  
Focus Position, 5–17  
Frame Averaging, 5–18



## Index

---

- B-Mode (cont'd)
    - Gray Scale, 5-20
    - Gray Scale Map, 5-16
    - Image Rotation, 5-22
    - Image Softener, 5-20
    - Imaging Freq, 5-18
    - Optimizing the Image, 5-10
    - Rejection, 5-23
    - Reverse, 5-15
    - Scan Area Position, 5-14
    - Scan Area Size, 5-14
    - Sub-Menu, 3-37
    - TGC, 5-11
    - Top Menu, 3-37
    - With Doppler, 5-49
  - B-Mode Analysis
    - Apical 2 Chamber, 10-83
    - Apical 4 Chamber, 10-81
    - Parasternal Long Axis, 10-77
    - Parasternal Short Axis
      - Aortic Valve, 10-78
      - Mitral Valve, 10-79
      - Papillary Muscles, 10-80
  - B/M Gain, 3-47
    - B-Mode, 5-13
    - M-Mode, 5-65
  - B/W Video Page Printer, 13-6
  - B510, 15-27
  - Backup, User Data, 14-104
  - Baseline Shift
    - Color Flow, 5-30
    - Doppler, 5-53
  - Bi Plane Ellipsoid Method, 10-17, 10-74
  - Bioeffects
    - Acoustic output, 17-6
    - ALARA, 17-9
    - Concerns surrounding diagnostic ultrasound, 17-2
    - Endnotes, 17-83
    - Mechanical, 17-4
    - Operator awareness, 17-5
    - Thermal, 17-3
    - Tissue characteristics, 17-5
    - Training and user assistance, 17-17
  - Biological Hazard, 2-3, 2-6
  - Biopsy
    - Depth Cursor, 5-23, 15-32
    - Guidezone(s), 5-22
    - Procedures
      - Accessories & supplies, 15-31
      - Guide attachment, 15-36
      - Probes, 15-48
      - Special concerns, 15-30
  - Blue Shift, 6-20
  - Body Pattern Copy to Active Side, System Parameters menu page 4, 14-55
  - Body Pattern Package [EC, UP], Preset Program menu page 1, 14-66
  - Body Pattern Package 1-8, System Parameters menu page 4, 14-55
  - Body Pattern Probe Change Erase, System Parameters menu page 4, 14-55
  - Body Pattern Probe Mark Preset [EC, UP], Preset Program menu page 1, 14-66
  - Body Pattern Unfreeze Erase, System Parameters menu page 4, 14-55
  - Body Patterns, 3-48, 6-22
    - Illustrations, Patterns Available, 6-23
    - Package Selections, 14-66
    - Programming Packages 1-8, 14-55
    - Sub-Menu, 3-44
  - BSA Type, System Parameters menu page 2, 14-52
  - Bullet Method, 10-11, 10-69
- ## C
- C358, 15-21
  - C364, 15-21
  - C386, 15-22
  - C551, 15-22
  - C721, 15-22
  - Calc Dir, Doppler, 5-60, 5-61
  - Calculation, Accuracy, 16-6
  - Calculation Formulas
    - Biplane Ellipsoid Method, 10-18
    - Bullet Method, 10-12
    - General, 8-11
    - GYN, 9-71
    - Left Ventricular, 10-7
      - Bi Plane Ellipsoid Method, 10-74
      - BP-DISC Method, 10-71
      - Bullet Method, 10-69
      - Cubed Method, 10-67
      - Gibson Method, 10-75
      - Modified Simpson's Rule Method, 10-73
      - Single Plane Ellipsoid Method, 10-73
      - SP-DISC Method, 10-69
      - Teichholz Method, 10-68
    - Modified Simpson's Rule Method, 10-14
    - Single Plane Ellipsoid Method, 10-16
    - Teichholz Method, 10-9
    - Vascular, 11-8
      - Advanced, 11-18



- Calculations
  - Cardiology, 10–19
  - Erasing, 7–6
  - Vascular, 11–3
- Capture, Color Flow, 5–35
- Capture Frame, 6–14
- Cardiac Calculation, Preset Program menu page 7, 14–88
- Cardiac Calculation Formulas
  - Aortic Valve, 10–86, 10–91
  - Apical 2 Chamber, 10–83
  - Apical 4 Chamber, 10–81
  - Left/Right Ventricle, 10–84
  - Mitral Valve, 10–85, 10–89
  - Parasternal Long Axis, 10–77
  - Parasternal Short Axis
    - Aortic Valve, 10–78
    - Mitral Valve, 10–79
    - Papillary Muscles, 10–80
  - Pulmonic Valve, 10–87, 10–93
  - Tricuspid Valve, 10–88, 10–95
- Cardiac Measurement Choices [EC, UP], Preset Program menu page 5 & 6, 14–73
- Cardiac Output, 8–4
  - Automatic Calculation, 8–4
- Cardiology, 10–1
  - Bi Plane Ellipsoid Method, 10–17
  - BSA Calculation Methods, 10–6
  - Bullet Method, 10–11
  - Calculations, 10–19
  - Cubed Method, 10–7
  - ECG, 10–31
  - Measurements
    - % Stenosis, 10–23
    - Angle, 10–22
    - ET, 10–25
    - Max PG, 10–26
    - Mean PG, 10–27
    - MVA, 10–24
    - PHT, 10–24
    - Volume, 10–19
  - Modified Simpson's Rule Method, 10–13
  - Overview, 10–3
  - Single Plane Ellipsoid Method, 10–15
  - Summary Report, Recording, 10–6
  - Teichholz Method, 10–9
- Cardiology Calculations
  - Trace Auto, 10–28
  - Transf Calcs, 10–30
- Care and Maintenance, Cleaning the system
  - Air Filters, 16–17
  - Cabinet, 16–12
  - Foot Switch, 16–15
  - Monitor, 16–13
  - Multi Image Camera, 16–15
  - Operator control panel, 16–14
  - Page Printer, 16–16
  - VCR, 16–15
- Caution, 2–2, 16–39
- CFM ACE [EC, UP, P, R], Custom Display menu page 7, 14–25
- CFM B Angle/Width (deg/mm) [EC, UP, P, R], Custom Display menu page 6, 14–24
- CFM B Frame Average [EC, UP, P, R], Custom Display menu page 6, 14–24
- CFM B Gain [EC, UP, P, R], Custom Display menu page 6, 14–23
- CFM B High Resolution [EC, UP, P, R], Custom Display menu page 6, 14–24
- CFM B Wall Filter [EC, UP, P, R], Custom Display menu page 6, 14–24
- CFM Base Line (%) [EC, R], Custom Display menu page 15, 14–41
- CFM Capture Interval (s) [EC, R], Custom Display menu page 15, 14–40
- CFM Color Map 1 [EC, UP], Custom Display menu page 15, 14–41
- CFM Color Map 2 [EC, UP], Custom Display menu page 15, 14–41
- CFM Color Map 3 [EC, UP], Custom Display menu page 15, 14–41
- CFM Color Map 4 [EC, UP], Custom Display menu page 15, 14–41
- CFM Color Map 5 [EC, UP], Custom Display menu page 15, 14–41
- CFM Color Map 6 [EC, UP], Custom Display menu page 15, 14–41
- CFM Color Threshold [EC, UP, P, R], Custom Display menu page 6, 14–22
- CFM M Gain (Delta from B) [EC, UP, P, R], Custom Display menu page 7, 14–25
- CFM M Wall Filter [EC, UP, P, R], Custom Display menu page 7, 14–25
- CFM MR–Flow [EC, UP, P, R], Custom Display menu page 7, 14–25
- CFM MR–Flow Combi Number [EC, UP, P], Custom Display menu page 7, 14–25
- CFM Noise Blanker [EC, UP, P, R], Custom Display menu page 7, 14–25
- CFM Packet Size Map [EC, UP, P], Custom Display menu page 6, 14–23
- CFM Packet Size MpDtl [EC, UP, P], Custom Display menu page 6, 14–23
- CFM Packet Size SrvyDtl [EC, UP, P], Custom Display menu page 6, 14–23



## Index

---

- CFM Packet Size Survey [EC, UP, P], Custom Display menu page 6, 14–23
- CFM Persistence [EC, UP, P, R], Custom Display menu page 7, 14–25
- CFM Tag Auto Start [EC, R], Custom Display menu page 15, 14–41
- CFM Tag Center (cm/s) [EC, R], Custom Display menu page 15, 14–41
- CFM Tag Width (cm/s) [EC, R], Custom Display menu page 15, 14–41
- CFM Target Frame Rate [EC, UP, P], Custom Display menu page 6, 14–23
- CFM Target Frame Rate Dtl [EC, UP, P], Custom Display menu page 6, 14–23
- CFM Threshold Turbulence [EC, UP, P, R], Custom Display menu page 6, 14–22
- CFM Threshold Velocity [EC, UP, P, R], Custom Display menu page 6, 14–22
- CFM Velo Dynamic Range [EC, UP], Custom Display menu page 15, 14–41
- CFM Velocity (cm/s) [EC, UP, P, R]
  - Custom Display menu page 5, 14–6
  - Custom Display menu page 6, 14–23
- CFM Wall Echo Canceller [EC, UP, P, R], Custom Display menu page 6, 14–22
- CFM/PDI Enhancement, 5–38
- CFM/PDI Focus Depth (%) [EC, UP, P], Custom Display menu page 6, 14–24
- CFM/PDI Initial Mode [EC, R], Custom Display menu page 16, 14–42
- CFM/PDI Initial Sequence [EC], Custom Display menu page 15, 14–41
- CFM/PDI Sequence [EC, UP], Custom Display menu page 15, 14–41
- CFM/PDI Shrink Angle/Width (deg, mm) [EC, UP, P, R], Custom Display menu page 4, 14–21
- CFM/PDI Shrink in Triplex [EC, UP, P, R], Custom Display menu page 4, 14–21
- CFM/PDI Spatical Filter [EC, UP, P], Custom Display menu page 6, 14–23
- CFM/PDI Start Depth (mm) [EC, UP, P, R], Custom Display menu page 6, 14–24
- CFM/PDI Vertical Size (mm) [EC, UP, P, R], Custom Display menu page 6, 14–24
- CFM/PWD Ratio in Triplex [EC, UP, P, R], Custom Display menu page 4, 14–21
- Change Probe and Category, System Parameters menu page 1, 14–49
- CINE, CINE Capture, 6–13
  - Cine, 6–11
    - Accessing, 6–9
    - Capture Frame, 6–14
    - Cine Gauge, 6–13
    - Cine Gauge/Image Tracking, 6–14
    - Exiting, 6–14
    - Functionality, 6–9
    - Helpful Hints, 6–14
    - Introduction, 6–8
    - Loop Speed, 6–12
    - Memory, 6–8
    - Multipl, 6–12
    - Side Change, 6–13
    - Sub-Menu, 3–43
    - Top Menu, 3–43
    - Using Cine Loop, 6–10
  - CINE Capture, 6–13
  - Cine Gauge Auto Display [EC], Custom Display menu page 10, 14–31
  - Cine Scroll, 3–47
  - Cine Speed [EC, R], Custom Display menu page 10, 14–31
  - Circ/Area Method [EC], Preset Program menu page 3, 14–68
  - Circuit breaker, 3–21
  - Class I Equipment, 2–9
  - Cleaning
    - Foot Switch, 16–15
    - Monitor, 16–13
    - Multi Image Camera, 16–15
    - Operator control panel, 16–14
    - Page Printer, 16–16
    - System cabinet, 16–12
    - VCR, 16–15
  - Clear, 3–49
  - Clinical instructions for fetal use, 17–10
    - Fetal heart, 17–12
    - Umbilical artery, 17–10
    - Uterine artery, 17–11
  - Col Applcat, System Parameters menu page 2, 14–53
  - Col Aux, System Parameters menu page 2, 14–53
  - Col Back Gnd, System Parameters menu page 2, 14–52
  - Col Comment, System Parameters menu page 2, 14–53
  - Col Dop Trac, System Parameters menu page 2, 14–53
  - Col Imag Crs, System Parameters menu page 2, 14–53

- Col Mes Crs, System Parameters menu page 2, 14–53
- Col Moving, System Parameters menu page 2, 14–53
- Col Scale, System Parameters menu page 2, 14–53
- Col Sys Inf, System Parameters menu page 2, 14–53
- Color
  - B–Mode, 5–20
  - Doppler, 5–58
  - Gray Scale
    - B–Mode, 5–20
    - Color Tag, 5–21
  - Basics, 5–25
  - Doppler, 5–58
    - Color Tag, 5–58
    - Tag Position, 5–58
  - M–Mode, 5–68
  - M–Mode, 5–68
- Color Flow
  - 3D Mode, 5–39
  - ACE, 5–38
  - Activating, 5–26
  - Baseline Shift, 5–30
  - Capture, 5–35
  - Color Window, 5–26
  - Common Controls, 5–28
  - Control Layout, 5–28
  - Diag Mode, 5–33
  - Display, 5–27
  - Exiting, 5–26
  - Frame Average, 5–34
  - Gain, 5–29
  - High Resolution, 5–35
  - M–Mode, 5–26
  - Maps, 5–32
  - Mode
    - Distance, 7–12
    - Echo Level (gray scale), 7–12
    - Trace, 7–12
    - Velocity Point, 7–13
  - MR–Flow, 5–39
  - MTI Filter, 5–33
  - Noise Blanker, 5–38
  - Packet Size, 5–36
  - Penetration (Penet.), 5–34
  - Persistence, 5–39
  - Select, 5–29
  - Slant Scan, 5–32
  - Spatial Filter, 5–36
  - Spectrum Invert, 5–29
  - Sub–Menu, 3–38
  - Tag Position, 5–37
  - Threshold, 5–35
  - Top Menu, 3–38
  - Velocity Scale, 5–30
  - Velocity Tag, 5–37
  - W.E. Cancel, 5–37
  - Window, Linear Probe, 5–32
  - Window Size, 5–31
- Color Image Archive, System Parameters menu page 5, 14–61
- Color Printer, Troubleshooting, 16–28
- Color Printer Memory, System Parameters menu page 5, 14–59
- Color Printer Signal, System Parameters menu page 5, 14–60
- Color Printer Type, System Parameters menu page 5, 14–59
- Color Tag
  - B–Mode, 5–21
  - Color Flow, 5–37
  - Doppler, 5–58
- Color Threshold, Color Flow, 5–35
- Color Video Page Printer, 13–9
- Comment, 3–48, 6–15
  - Displaying, 6–18
  - Editing, 6–17, 6–21
  - Entering, 6–17
  - Scripts, 6–18
  - Sub–Menu, 3–44
- Comment Clear Key Function [EC], Preset Program menu page 1, 14–66
- Comment Erase at Display Change [EC], Preset Program menu page 1, 14–66
- Comment Home Position Horizontal [EC, R], Preset Program menu page 1, 14–66
- Comment Home Position Vertical [EC, R], Preset Program menu page 1, 14–66
- Comments, Venous Page, 11–16
- Console Labels, 2–19
- Console Overview, 3–4
- Contraindications, Doppler, 1–6
- Control
  - Acoustic Output, 5–11
  - B/M Gain, 3–47
  - Cine Scroll, 3–47
  - Depth, 3–46
  - Rotation, 3–49
  - TGC, 5–11
  - Zoom Size, 3–49
- Control Panel
  - Replacing key caps, 16–16
  - Replacing key lamps, 16–16
- Control Z, VCR Counter Reset, 19–21
- Controls, 3–31
- Counter, VCR Display, 19–18



## Index

---

- Coupling Gels, 15–19
- CUA Selection Reset on New Patient, System Parameters menu page 1, 14–50
- CUA/AUA for HADLOCK, System Parameters menu page 3, 14–54
- Cubed Method, 10–7, 10–67
- Cursor
  - Annotation, 6–15
  - Biopsy Depth, 5–23, 15–32
  - Color Flow Doppler, 5–32
  - Doppler, 5–50
- Cursors, 7–4
- Custom Display Presets, 14–7
  - Contents, 14–12
- Customizing Your System, 14–1
- CWD
  - Sub Menu, 3–40
  - Top Menu, 3–40
- CWD Base Line (%) [EC, R], Custom Display menu page 15, 14–40
- CWD Dynamic Range [EC, UP, P, R], Custom Display menu page 5, 14–22
- CWD Gain [EC, UP, P, R], Custom Display menu page 5, 14–22
- CWD Gray Scale Map [EC, UP, P, R], Custom Display menu page 5, 14–22
- CWD Non-aliasing mode [EC, R], Custom Display menu page 15, 14–40
- CWD Rejection [EC, UP, P, R], Custom Display menu page 5, 14–22
- CWD Velocity (cm/s) [EC, UP, P, R], Custom Display menu page 5, 14–22
- CWD Wall Filter (cm/s) [EC, UP, P, R], Custom Display menu page 5, 14–22
- CWD2, 15–29
- CWD5, 15–29
- D**
- D Audio to Speaker [EC], Custom Display menu page 14, 14–39
- D Audio Volume (%) [EC, R], Custom Display menu page 14, 14–39
- D Auto Caliper Control [EC, UP], Preset Program menu page 4, 14–72
- D Color Auto Start [EC, R], Custom Display menu page 13, 14–37
- D Color Map [EC, UP], Custom Display menu page 13, 14–37
- D Format (Side/Side, Full) [EC], Custom Display menu page 13, 14–37
- D Image Processing [EC], Custom Display menu page 14, 14–39
- D Mode Full B+ Cursor [EC], Custom Display menu page 13, 14–38
- D Mode Full D [EC], Custom Display menu page 13, 14–38
- D Mode Side/Side B Large [EC], Custom Display menu page 13, 14–38
- D Mode Side/Side B Small [EC], Custom Display menu page 13, 14–38
- D Mode Top/Bottom B Large [EC], Custom Display menu page 13, 14–38
- D Mode Top/Bottom B Mid [EC], Custom Display menu page 13, 14–38
- D Mode Top/Bottom B Small [EC], Custom Display menu page 13, 14–38
- D Realtime Auto Trace [EC, R], Custom Display menu page 13, 14–36
- D Realtime Trace Method [EC, R], Custom Display menu page 13, 14–37
- D Tag Auto Start [EC, R], Custom Display menu page 13, 14–37
- D Tag Center (%) [EC, R], Custom Display menu page 13, 14–37
- D Tag Width (%) [EC, R], Custom Display menu page 13, 14–37
- D Time Resolution [EC], Custom Display menu page 14, 14–39
- D Type for Wall Filter [EC], Custom Display menu page 13, 14–36
- D Video Inverse [EC], Custom Display menu page 14, 14–39
- D/CFM Scale Mark [EC], Custom Display menu page 13, 14–36
- Danger, 2–2, 2–5
- Data Management Center (DMC), 9–61
- Date Format, System Parameters, System Parameters menu page 1, 14–49
- Defibrillator Caution, 2–19
- Delay Trigger Source, System Parameters menu page 2, 14–51
- Depth, 3–46
  - B-Mode, 5–13

- Device Labels, 2–7
- Diag Category, System Parameters, System Parameters menu page 1, 14–49
- Diag Mode, Color Flow, 5–33
- Diagnostics, Operator, 16–34
  - Accessing, 16–35
- Diastole Velocity for PI [EC, UP], Preset Program menu page 4, 14–71
- Diastole Velocity for RI [EC, UP], Preset Program menu page 4, 14–71
- Diastole Velocity for S/D (D/S) [EC, UP], Preset Program menu page 4, 14–71
- Diastole/Systole Determination [EC], Preset Program menu page 3, 14–70
- Diastolic/Systolic Ratio (D/S), 9–23
- DICOM, 20–1, 20–15
  - Application Information, 20–5
  - Host Verification, 20–8
  - Image Transfer, 20–10
  - Network Configuration, 20–3
  - Network Verification, 20–8
  - Routing Information, 20–4
- Disk Drive, Floppy Disk, 3–6
- Display
  - A Output, 5–6
  - Acoustic Output %, 5–5
  - Angle, 5–46
  - B Gray Scale Map Number, 5–5
  - B–Mode, 5–4
  - B–Mode Format, 5–4
  - B/W Threshold Marker, 5–27
  - Body Pattern, 5–6
  - CFM Area Cursor, 5–27
  - CG, 5–27
  - Color Flow, 5–27
  - Color Scale, 5–5, 5–27
  - D. Velocity, 5–46
  - Date, 5–5
  - Depth, 5–6
  - DF, 5–46
  - DG, 5–46
  - Doppler, 5–45
    - Dual Spectrum, 5–63
    - Spectrum, 5–62
  - Dual B–Mode Display Format, 5–86
  - Dual Format, 5–87
  - Dyn, 5–6
  - FR/Cine, 5–6
  - Gain, 5–6
  - Gray Scale, 5–5
  - Hospital Name, 5–5
  - HPRF, 5–46
  - HR, 5–6
  - ID, 5–5
  - Image Rotation, 5–22
  - LMP or GA, 5–5
  - M–Mode
    - Dual Spectrum, 5–63
    - Spectrum, 5–62
  - Measurements, 5–6
  - MTI Filter, 5–27
  - Operator Messages, 5–6
  - Patient Name, 5–5
  - PRF, 5–46
  - Probe
    - Name, 5–5
    - Orientation, 5–5
  - Sample Volume Length (SVL), 5–46
  - Scale Marker, 5–9
  - Scroll Depth, 5–9
  - Soft Menu, 3–34
  - SPEC INV, 5–46
  - Spectrum Ref. Value, 5–46
  - TGC Curve, 5–9
  - Time, 5–5
  - Time Scale, 5–46
  - Top/Bottom B Large Preset, 5–85
  - Top/Bottom B Mid Preset, 5–84
  - Top/Bottom B Small Preset, 5–85
  - TV Counter, 5–6
  - Units (V or F), 5–27
  - Velocity Scale, 5–27
  - VTR Counter, 5–6
  - VTR Status, 5–6
  - Wall Filter (WV), 5–46
  - Zoom, 5–9
    - with Reference, 5–86
- Display Dop Frequency [EC], Custom Display menu page 14, 14–39
- Display Dop PRF [EC], Custom Display menu page 14, 14–39
- Display EDD with GA Value [EC], Preset Program menu page 3, 14–69
- Display Format(Dual), B–Mode, 5–15
- Display Height 1–3 [EC], Preset Program menu page 10, 14–90
- Display MaxPG with Velocity [EC, UP], Preset Program menu page 4, 14–72
- Display Messages, 16–26
- Display Probe Name, System Parameters menu page 1, 14–51
- Display Spectrm Ref. Value [EC], Custom Display menu page 14, 14–39
- Display TGC Curve [EC], Custom Display menu page 10, 14–32
- Display Unit Age 1–4 [EC], Preset Program menu page 10, 14–90
- Display Unit Weight 1–2 [EC], Preset Program menu page 10, 14–90
- Displayed TI [EC], Custom Display menu page 10, 14–30



## Index

---

### Doppler

- Activating CWD, 5–47
- Activating PWD, 5–47
- Angle Correction, 5–51
- Audio Volume, 5–50
- B Pause, 5–53
- B–Mode Controls, 5–49
- Baseline Shift, 5–53
- Calc Dir, 5–60, 5–61
- Color, 5–58
- Color Flow
  - Activating, 5–26
  - Common Controls, 5–28
  - Control Layout, 5–28
  - M–Mode, 5–26
- Color Gain, 5–29
- Color Tag, 5–58
- Color Window, 5–26
- Colorized Gray Scale, 5–58
  - Color Tag, 5–58
  - Tag Position, 5–58
- Continuous Wave
  - Definition, 5–44
  - Non–Imaging, 5–44
  - Steerable, 5–44
- Contraindications, 1–6
- Control Layout, 5–48
- CWD, Sub–Menus, 5–47
- Display, 5–45, 5–62
- Dual Display, 5–63
- Dynamic Range, 5–54
- Fetal Exams
  - Contraindications, 9–4
  - Indications, 9–4
- Fetal use, 1–5
- High PRF, 5–43
- HPRF, 5–59
- M/D Cursor, 5–50
- Optimization, 5–48
- Peak Velocity, 7–14
- Penetration, 5–57
- Pulsed Wave
  - Definition, 5–42
  - Frequencies Used, 5–43
- PWD, Sub–Menu, 5–47
- RealTime Trace, 5–60
- Rejection, 5–59
- Sample Volume Length, 5–56
- Slant Scan, 5–54
- Spectral Gain, 5–51
- Spectrum, Gain, 5–51
- Spectrum Invert, 5–52
- Sweep Speed, 5–57
- Tag Position, 5–58
- TAMAX/TAMIN/TAMEAN/TAMODE, 7–15
- Velocity Scale, 5–52
- Wall Filter, 5–56

### Doppler Analysis

- Aortic Valve, 10–91
- Mitral Valve, 10–89

- Pulmonic Valve, 10–93

- Tricuspid Valve, 10–95

- Doppler Channel [EC, R], Custom Display menu
  - page 13, 14–37

- Doppler Realtime Calculation (EC, UP), Preset Program menu
  - page 7, 14–88

- Dual, 10–34

- Dual Format, B–Mode, 5–15

### Dynamic Range

- B–Mode, 5–16
- Doppler, 5–54
- M–Mode, 5–66

## E

- E721, 15–23

### ECG, 10–31

- Cine Gauge/Image Tracking, 6–14, 10–36
- Gain, 10–36
- Lead Placement, 10–33
- Positn, 10–36
- Soft Menu, 10–32
- Sub–Menu, 3–42
- Top Menu, 3–42
- Wave, 10–35

- ECG Audio Beat Sound [EC, R], Custom Display menu
  - page 17, 14–43

- ECG Audio Beat Sound Tone [EC], Custom Display menu
  - page 17, 14–43

- ECG Auto Display [EC], Custom Display menu
  - page 17, 14–43

- ECG Gain [EC, R], Custom Display menu
  - page 17, 14–43

- ECG Heart Rate Display [EC], Custom Display menu
  - page 17, 14–43

- ECG Position [EC, R], Custom Display menu
  - page 17, 14–43

- ECG R Delay 1 (ms) [EC, UP, R], Custom Display menu
  - page 17, 14–43

- ECG R Delay 1 to 2 (ms) [EC, UP, R], Custom Display menu
  - page 17, 14–43

- ECG Sync Mark Display [EC], Custom Display menu
  - page 17, 14–44

- Echo Level Method [EC], Preset Program menu
  - page 3, 14–68

### Edge Enhance

- B–Mode, 5–23
- M–Mode, 5–67

- EFBW Form (Europe), System Parameters menu
  - page 2, 14–52

EFBW Form (Tokyo), System Parameters menu  
page 2, 14–52

Electrical Hazard, 2–3, 2–4, 2–6

Ellipse, 3–49, 7–3

EMC (Electromagnetic Compatibility), 2–10

End Frame, 6–10

End Systole Period from Prior Edge [EC, UP],  
Preset Program menu page 3, 14–70

Endometrium Thickness, 9–68

Enter Depth [EC, UP, P, R], Custom Display menu  
page 1, 14–13

Equipment Safety, 2–5

ET (Ejection Time), 10–25

Exam  
Application Preset Selection, 4–7  
Application Presets, Defining a User Preset,  
14–93  
Beginning, 4–2

Explosion Hazard, 2–3, 2–5

External Video, 3–50

External Video Signal, System Parameters menu  
page 5, 14–59

## F

Fetal Doppler, 9–4

Efficacy, 17–13  
Cardiac anomalies, 17–13  
IUGR, 17–13  
Summary, 17–14  
Suggested guidelines, 17–14  
Methodology, 17–15  
Variance studies, 17–16

Fetal Growth Trend  
Data List Management, 9–51  
Growth Trending, 9–44  
List ID Management, 9–45  
Storing Patient Information, 9–39

Fetal Heart Rate, 9–23

Fetal Trend Management, 9–39  
Multigestation option, 9–59

Fetus Number, 9–53

Floppy Disk Drive, 3–6

Flow Volume (FV), 8–7

Focus  
Number, B–Mode, 5–17

Position, B–Mode, 5–17

Focus Depth [EC, UP, P, R], Custom Display menu  
page 1, 14–13

Foot Switch, 3–9  
Cleaning, 16–15

Frame Averaging  
B–Mode, 5–18  
Color Flow, 5–34

Frame Rate Display [EC], Custom Display menu  
page 10, 14–31

Freeze, 3–47, 6–6  
Foot Switch, 6–7

## G

Gain, Color Doppler, 5–29

General Calculations  
% Stenosis, 8–3  
Angle, 8–3  
Cardiac Output (CO), 8–4  
CFM Mode, Velocity Point. *See General Calculations*  
Flow Volume (FV), 8–7  
Heart Rate (HR), 8–6, 8–8, 10–30  
Radiology/Abdomen. *See General Calculations*  
Small Parts. *See General Calculations*  
Stroke Volume Ratio (SV), 8–5  
Trace Auto, 8–3, 8–8, 11–5  
Transf Calcs, 8–3  
Urology. *See General Calculations*  
Volume, 8–2

General Measurements, 7–1

Gestational Sac (GS), 9–19

Gibson Method, 10–75

Gray Scale  
Colorized  
B–Mode, 5–20  
Color Tag, 5–21  
Basics, 5–25  
Doppler, 5–58  
Color Tag, 5–58  
Tag Position, 5–58  
M–Mode, 5–68  
M–Mode, 5–66  
Map, B–Mode, 5–16

Ground Point, 2–20

GS, OB Measurements, 9–19

Guidezone(s), Biopsy, B–Mode, 5–22





## Index

---

### GYN

- Calculation Formulas, 9–71
- Exam Preparation, 9–3
- Measurements, 9–67, 9–68
  - Endometrium Thickness, 9–68
  - Ovarian Height, 9–67
  - Ovarian Length, 9–67
  - Ovarian Width, 9–67
  - RI, 9–69
  - Uterine Height, 9–68
  - Uterine Length, 9–68
  - Uterine Width, 9–68
- Summary Report, Layout, 9–70

## H

### Hazard

- Acoustic Output, 2–3
- Symbols, 2–3

Heart Rate (HR), 8–6, 8–8, 9–23, 10–30, 11–5

Heart Rate Method [EC], Preset Program menu  
page 3, 14–68

High Resolution, Color Flow, 5–35

Hip Orientation [EC, UP], Preset Program menu  
page 4, 14–71

Hospital Name, System Parameters menu page 1,  
14–49

Host Verification, 20–8

HPRF, Doppler, 5–59

HR, 9–23, 11–5

## I

I739, 15–24

ID/Name, 4–6

ID/Name Prohibition After Measurement, System  
Parameters menu page 2, 14–52

Illumination, key, 3–32

### Image

- Annotation, 6–15, 6–19
- Freeze, 6–6
- Memory, 13–2

### Printing

- B/W Video Page Printer, 13–6
- Color Video Page Printer, 13–9
- Multi-Image Camera (MIC), 13–17
- Recall, 13–4
- Recording, 13–2
- Rotation, B-Mode, 5–22
- Softener, B-Mode, 5–20

### Image Archive

- Disk Verification, 13–27
- Image Recall Process, 13–32
- Image Selection for Recall, 13–32
- Media Format, 13–26
- Media Search, 13–33
- MO Eject, 13–35
- MOD, 13–23
- Patient Search, 13–28
- Preset Parameters, 13–24
- Storing Images, 13–27

Image Archive Compression [EC], Preset Program  
menu page 2, 14–67

### Image Archive Option

- Sub-Menu, 3–42
- Top Menu, 3–42

Image Back Color [EC, UP], Custom Display menu  
page 10, 14–31

Image Gamma (%) [EC, UP], Custom Display menu  
page 10, 14–30

Image Memory, 3–46

Image Mode [EC, R], Custom Display menu page  
10, 14–30

Image Recall, 3–46, 13–4  
Sub-Menu, 3–45

Image Reverse [EC, UP, P, R], Custom Display  
menu page 1, 14–13

Image Rotation [EC, UP, P, R], Custom Display  
menu page 1, 14–13

Image Transfer, 20–10

Imaging Freq, B-Mode, 5–18

Information, Requesting, 16–38

Insite Access Enable, System Parameters menu  
page 1, 14–50

Introduction, Who to Contact, 1–7



## K

### Key

- Body Pattern, 3–48
- Clear, 3–49
- Comment, 3–48
- Ellipse, 3–49, 7–3
- External Video, 3–50
- Freeze, 3–47
- Illumination, 3–32
- Image Memory, 3–46
- Image Recall, 3–46
- M/D Cursor, 3–48
- Measurement, 3–48, 7–3
- Mic, 3–50
- Mode, 3–46
- New Patient, 3–46
- Pause, 3–50
- Play, 3–50
- Probe Controls, 3–29
- Record, 3–50
- Record 1, 3–47
- Record 2, 3–47
- Reverse, 3–46
- Scan Area, 3–48
- Set, 3–49
- Soft Menu, 3–34
- Stop, 3–50
- Zoom, 3–49

- Keyboard, 3–31, 3–33
  - Back Space, 3–33
  - Blue Shift, 3–33
  - Caps Lock, 3–33
  - Control, 3–33
  - Cursor Home, 3–33
  - Escape, 3–33
  - Red Shift, 3–33
  - Return, 3–33
  - Tab, 3–33

- Keyboard LED Bright, System Parameters menu page 1, 14–50

- Keyboard Menu Bright, System Parameters menu page 1, 14–50

- Keyboard Tab, System Parameters menu page 1, 14–50

## L

- L764, 15–25

- LA39, 15–25

- Labels, Icon Description, 2–7

- Language, System Parameters menu page 1, 14–49

- Laser Camera, 13–22

- Leakage Current, 2–9

- Left Ventricular, Calculation Formulas, 10–7

- Bi Plane Ellipsoid Method, 10–74

- BP-DISC Method, 10–71

- Bullet Method, 10–69

- Cubed Method, 10–67

- Gibson Method, 10–75

- Modified Simpson's Rule Method, 10–73

- Single Plane Ellipsoid Method, 10–73

- SP-DISC Method, 10–69

- Teichholz Method, 10–68

- Left/Right Ventricle, 10–84

- Length Unit [EC], Preset Program menu page 3, 14–68

- Library, Annotation, 6–17

- Entering/Editing, 6–17

- Line Printer

- Basic Operation, 13–21

- Error Messages, 16–29

- Setup, 13–21

- Locatn (Location), 9–22

- LOGIQ 500

- Back view, 3–5

- Features, 1–6

- Front view, 3–5

- General Indications for Use, 1–5

- Left side view, 3–4

- Moving, 3–11, 3–12

- Relocating, 3–10

- Right side view, 3–4

- Transporting, 3–13

- Wheels, Setting the Lock, 3–14

- Loop, Cine, 6–10

- Loop Speed, 6–12

- LT Bifurc, 11–3

- LT CCA, 11–3

- LT ECA, 11–3

- LT ICA, 11–3

- LT ICA/CCA, 11–4

- LV BP-DISC Method, 10–71

- LV SP-DISC Method, 10–69



## Index

---

### M

- M Auto Caliper Control [EC, UP], Preset Program menu page 4, 14–72
- M Color Auto Start [EC, R], Custom Display menu page 12, 14–34
- M Color Map [EC, UP], Custom Display menu page 12, 14–34
- M Dynamic Range [EC, UP, P, R], Custom Display menu page 3, 14–18
- M Edge Enhance [EC, UP, P, R], Custom Display menu page 3, 14–18
- M Format (Side/Side, Full) [EC], Custom Display menu page 12, 14–35
- M Gain (Delta from B) [EC, UP, P, R], Custom Display menu page 3, 14–18
- M Gray Scale Map [EC, UP, P, R], Custom Display menu page 3, 14–19
- M Image Frequency [EC, UP, P, R], Custom Display menu page 3, 14–19
- M Image Softener [EC, UP, P, R], Custom Display menu page 3, 14–18
- M Mode FULL B + Cursor [EC], Custom Display menu page 12, 14–35
- M Mode Full M [EC], Custom Display menu page 12, 14–36
- M Mode Side/Side B Large [EC], Custom Display menu page 12, 14–35
- M Mode Side/Side B Small [EC], Custom Display menu page 12, 14–35
- M Mode Simul PD Activity [EC], Custom Display menu page 12, 14–34
- M Mode Top/Bottom B Large [EC], Custom Display menu page 12, 14–35
- M Mode Top/Bottom B Mid [EC], Custom Display menu page 12, 14–35
- M Mode Top/Bottom B Small [EC], Custom Display menu page 12, 14–36
- M Peak Hold [EC, UP, P, R], Custom Display menu page 3, 14–19
- M Rejection [EC, UP, P, R], Custom Display menu page 3, 14–19
- M Scale Mark [EC], Custom Display menu page 12, 14–36
- M Softener Level [EC, UP, P, R], Custom Display menu page 3, 14–18
- M Tag Auto Start [EC, UP], Custom Display menu page 12, 14–34
- M Tag Center (%) [EC, R], Custom Display menu page 12, 14–34
- M Tag Width (%) [EC, R], Custom Display menu page 12, 14–35
- M Video Inverse [EC], Custom Display menu page 12, 14–36
- M–Mode, 5–64
  - B/M Gain, 5–65
  - Color, 5–68
  - Color Flow, 5–26
  - Colorized Gray Scale, 5–68
  - Depth Difference, 7–25
  - Display, 5–62
  - Dual Display, 5–63
  - Dynamic Range, 5–66
  - Edge Enhance, 5–67
  - Gray Scale, 5–66
  - Introduction, 5–62
  - M/D Cursor, 5–65
  - Optimizing, 5–64
  - Rejection, 5–67
  - Sub Menu, 3–41
  - Sweep Speed, 5–68
  - Time, 7–18, 7–24
  - Tissue Depth, 7–23
  - Top Menu, 3–41
  - Typical Exam, 5–62
  - Zoom, 5–65, 6–4
- M–Mode Analysis
  - Aortic Valve, 10–86
  - Left/Right Ventricle, 10–84
  - Mitral Valve, 10–85
  - Pulmonic Valve, 10–87
  - Tricuspid Valve, 10–88
- M/D Cursor, 3–48
  - Doppler, 5–50
  - M–Mode, 5–65
- Maintenance, 16–11
  - Cleaning the Air Filters, 16–17
  - User, 16–1
- Mask Image [EC], Preset Program menu page 3, 14–69
- Mask Image Characters, System Parameters menu page 2, 14–51
- Maskline Record, System Parameters menu page 5, 14–60
- Max PG, 8–3, 10–26
  - Automatic Calculation, 10–26
  - Cardiology Calculations, 10–26
  - General Calculations, 8–3
- Max Velocity (Advanced Vascular) [EC], Preset Program menu page 3, 14–69
- Mean PG, 8–3, 10–27
  - Abdomen/Small Parts Calculations, 8–3
  - Cardiology Calculations, 10–27

- Measurement, 3–48, 7–3
  - Accuracy, 16–5
  - Cursors, 7–4
  - Key, 7–5
- Measurement Clear Operation [EC], Preset Program menu page 3, 14–70
- Measurement Menu B [EC, UP], Preset Program menu page 7, 14–85
- Measurement Menu CFM [EC, UP], Preset Program menu page 7, 14–85
- Measurement Menu D [EC, UP], Preset Program menu page 7, 14–85
- Measurement Menu M [EC, UP], Preset Program menu page 7, 14–85
- Measurement Pressed, Sub-Menu, 3–44
- Measurements
  - Abdomen/Small Parts
    - A/B Ratio, 8–3, 10–29
    - Pulsatility Index (PI), 8–3, 10–29
    - Resistance Index (RI), 8–3, 10–29
    - S/D Ratio, 8–3, 10–29
  - B-Mode
    - Circumference/Area, 7–10
    - Distance, 7–7
    - Echo Level, 7–11
  - Bi Plane Ellipsoid Method, 10–17
  - Bullet Method, 10–11
  - Calculations, 7–6
  - Cardiology
    - % Stenosis, 10–23
    - Angle, 10–22
    - ET, 10–25
    - MVA, 10–24
    - PHT, 10–24
    - Volume, 10–19
  - Color Flow
    - Distance, 7–12
    - Echo Level (gray scale), 7–12
    - Trace, 7–12
    - Velocity Point. *See General Calculations*
  - Controls, 7–3
  - Cubed Method, 10–7
  - Doppler
    - Peak Velocity, 7–14
    - TAMAX/TAMIN/TAMEAN/TAMODE, 7–15
  - General, 7–1
  - GYN, 9–67, 9–68
    - Endometrium Thickness, 9–68
    - Ovarian Height, 9–67
    - Ovarian Length, 9–67
    - Ovarian Width, 9–67
    - RI, 9–69
    - Uterine Height, 9–68
    - Uterine Length, 9–68
    - Uterine Width, 9–68
  - M-Mode
    - Depth Difference, 7–25
    - Time, 7–18, 7–24
  - Tissue Depth, 7–23
  - Modified Simpson's Rule Method, 10–13
  - OB
    - A/B Ratio, 9–23
    - AFI, 9–20
    - D/S Ratio, 9–23
    - GS, 9–19
    - Hints, 9–23
    - HR, 9–23
    - Locatn, 9–22
    - Pulsatility Index (PI), 9–23
    - Resistance Index (RI), 9–23
    - S/D Ratio, 9–23
    - Trace Auto, 9–23
  - Radiology/Abdomen. *See General Calculations*
  - Single Plane Ellipsoid Method, 10–15
  - Small Parts. *See General Calculations*
  - Teichholz Method, 10–9
  - Urology. *See General Calculations*
  - Vascular, 11–3
    - A/B Ratio, 11–5
    - Heart Rate, 11–5
    - Pulsatility Index (PI), 11–5
    - Resistance Index (RI), 11–5
    - RT(LT) Bifurc, 11–3
    - RT(LT) CCA, 11–3
    - RT(LT) ECA, 11–3
    - RT(LT) ICA, 11–3
    - RT(LT) ICA/CCA, 11–4
    - S/D Ratio, 11–5
  - Volume, 10–19
- Media Selection for Fetal Trend [EC], Preset Program menu page 3, 14–70
- Memory
  - Cine, 6–8
  - Image, 13–2
    - Storage Space, 13–2
  - Image Recall, 13–4
- Menu Selection at New Patient, System Parameters menu page 7, 14–63
- Mic, 3–50, 19–29
- MIC (Multi-Image Camera), 13–17
- Microphone. *See Mic*
- MIP Image Rendering, 5–75
- Mitral Valve, 10–85, 10–89
- Mode, Controls, 3–46
- Modified Simpson's Rule Method, 10–13, 10–73
- Monitor
  - Adjustment
    - Brightness, 3–23
    - Contrast, 3–23
    - Position, 3–22
  - Cleaning, 16–13
  - Labels, 2–16
- Moving Hazard, 2–3



## Index

---

MR–Flow, Color Flow, 5–39

MTI Filter, Color Flow, 5–33

Multi Image Camera

Cleaning, 16–15

Troubleshooting, 16–29

Multigestation

Change Number of Fetuses, 9–55

Distinguishing each Fetus, 9–54

Fetal Trend Management, 9–59

Measurements/Calcs, 9–54

OB Graph, 9–57

OB option, 9–53

Patient Entry Menu, 9–53

Report Page Layout, 9–56

MVA (Mitral Valve Area), 10–24

## N

Needle Guide Type [EC, UP, P, R], Custom Display  
menu page 1, 14–14, 15–33

Network Configuration, 20–8

Network Verification, 20–8

New Patient, 4–3

Patient Entry Menu, 4–3

New Patient Y/N Ask [EC], Preset Program menu  
page 10, 14–90

Noise Blanker, Color Flow, 5–38

Non–Ionizing Radiation, 2–3

## O

OB, 9–24

Anatomical Survey, 9–33

Editing, 9–34

User Programmed Features, 9–34

Data Management Center (DMC), 9–61

Exam Preparation, 9–3

Fetal Trend Management. *See Fetal Growth  
Trend*

Formulas

European Version, 9–14

Osaka University Method, 9–10

Other Available Formulas, 9–16

Tokyo University Method, 9–8

USA Version, 9–11

Graph, 9–35

Changing Selection, 9–38

Selection, 9–36

Measurements

A/B Ratio, 9–23

AFI, 9–20

D/S Ratio, 9–23

GS, 9–19

Hints, 9–23

HR, 9–23

Locatn, 9–22

Pulsatility Index (PI), 9–23

Resistance Index (RI), 9–23

S/D Ratio, 9–23

Soft Menu, 9–7

Trace Auto, 9–23

Summary Report, 9–24

Editing, 9–31

European Version Layout, 9–27

Recording, 9–32

USA Version Layout, 9–26

Table Data, 18–1

OB Calculations, Transf Calcs, 9–23, 11–5

OB Format Selection, 9–6

OB Graph, Multigestation, 9–57

OB Measurements, Procedures, Amniotic Fluid  
Index, 9–20

OB Report Format, System Parameters menu page 2,  
14–52

OFD, Measurement from HC, 9–26

Operation Error Beep, System Parameters menu  
page 1, 14–50

Operator intervention, 17–7

Optimizing, B–Mode, 5–10

Option, Advanced Vascular, 11–9

Options, Advanced Cardiac Calculations (AMCAL),  
10–37

Ov–H, 9–67

Ov–L, 9–67

Ov–W, 9–67

Ovarian

Height, 9–67

Length, 9–67

Width, 9–67

## P

P509, 15–28

Packet Size

Affect on frame rate, 5–36  
Color Flow, 5–36

Page Printer, Cleaning, 16–16

Parasternal Long Axis, 10–77

Parasternal Short Axis

Aortic Valve, 10–78  
Mitral Valve, 10–79  
Papillary Muscles, 10–80

Password, Protection at Power Up, 3–18

Password 1–8, System Parameters menu page 6, 14–61

Password Ask, System Parameters menu page 6, 14–61

Patient Entry Menu, 4–3

Patient Name Format, System Parameters menu page 7, 14–63

Patient Safety, 2–4

Pause, 3–50

PCG

Gain, 10–36  
Positn, 10–36  
Wave, 10–35

PCG Auto Display [EC], Custom Display menu page 17, 14–44

PCG Filter, System Parameters menu page 2, 14–51

PCG Gain [EC, R], Custom Display menu page 17, 14–44

PCG Position [EC, R], Custom Display menu page 17, 14–44

PDI, 5–41

Optimizing, 5–41

PDI ACE [EC, UP, P, R], Custom Display menu page 8, 14–27

PDI B Frame Average [EC, UP, P, R], Custom Display menu page 8, 14–27

PDI B Gain [EC, UP, P, R], Custom Display menu page 8, 14–27

PDI B High Resolution [EC, UP, P, R], Custom Display menu page 8, 14–27

PDI B Wall Filter [EC, UP, P, R], Custom Display menu page 8, 14–27

PDI B/W Display [EC], Custom Display menu page 16, 14–42

PDI Capture Interval (s) [EC, R], Custom Display menu page 16, 14–42

PDI Color Threshold [EC, UP, P, R], Custom Display menu page 8, 14–26

PDI Default Map [EC], Custom Display menu page 16, 14–42

PDI Dynamic Range [EC, UP], Custom Display menu page 16, 14–42

PDI M Gain (Delta from B) [EC, UP, P, R], Custom Display menu page 8, 14–27

PDI M Wall Filter [EC, UP, P, R], Custom Display menu page 8, 14–27

PDI MR–Flow [EC, UP, P, R], Custom Display menu page 8, 14–28

PDI MR–Flow Combi Number [EC, UP, P], Custom Display menu page 8, 14–28

PDI Noise Blanker [EC, UP, P, R], Custom Display menu page 8, 14–27

PDI Packet Size Map [EC, UP, P], Custom Display menu page 8, 14–26

PDI Packet Size MpDtl [EC, UP, P], Custom Display menu page 8, 14–26

PDI Packet Size SrvyDtl [EC, UP, P], Custom Display menu page 8, 14–26

PDI Packet Size Survey [EC, UP, P], Custom Display menu page 8, 14–26

PDI Persistence [EC, UP, P, R], Custom Display menu page 8, 14–28

PDI Tag Auto Start [EC, R], Custom Display menu page 16, 14–42

PDI Tag Center (cm/s) [EC, R], Custom Display menu page 16, 14–42

PDI Tag Width (cm/s) [EC, R], Custom Display menu page 16, 14–42

PDI Target Frame Rate [EC, UP, P], Custom Display menu page 8, 14–26

PDI Target Frame Rate Dtl [EC, UP, P], Custom Display menu page 8, 14–26

PDI Threshold Power [EC, UP, P, R], Custom Display menu page 8, 14–26

PDI Velocity (cm/s) [EC, UP, P, R], Custom Display menu page 8, 14–26

PDI/PWD Ratio in Triplex [EC, UP, P, R], Custom Display menu page 4, 14–21

Penet., Color Flow, 5–34

Penetration, Doppler, 5–57

Peripherals

Black/White Video Printer, 13–7  
Color Video Printer, 13–10  
Connector Panel, 3–7  
Multi–Image Camera, 13–18  
S–VHS Video Cassette Recorder, 13–13



## Index

---

- Persistence, Color Flow, 5–39
- Personnel Safety, 2–5
- Phantoms, 16–46
- PHT (Pressure Half Time), 10–24
- Physio Sweep Speed, 10–32
- Physio Sweep Speed on B [EC, UP, P, R], Custom Display menu page 3, 14–18
- Play, 3–50
- Plug 1–BW, System Parameters menu page 5, 14–60
- Plug 2–CLR, System Parameters menu page 5, 14–60
- Port A, System Parameters menu page 5, 14–60
- Port B, System Parameters menu page 5, 14–60
- Port MIC, System Parameters menu page 5, 14–60
- Post Processing
  - B–Mode
    - Color, 5–20
    - Color Tag, 5–21
    - Gray Scale Map, 5–16
  - Color Flow
    - Color Tag, 5–37
    - Color Threshold, 5–35
    - Map, 5–32
    - Packet Size, 5–36
    - Spatial Filter, 5–36
    - Tag Position, 5–37
    - W.E. Cancel, 5–37
  - Doppler
    - Color, 5–58
    - Color Tag, 5–58
    - Tag Position, 5–58
  - Image Rotation, 5–22
  - M–Mode
    - Color, 5–68
    - Gray Scale Map, 5–66
    - Rejection, 5–67
  - Rejection
    - B–Mode, 5–23
    - Doppler, 5–59
- Power, 3–15
  - Circuit Breaker, 3–21
  - Connection, 3–15
  - Cord, 3–10
  - Off, 3–20
  - On/Off Standby, 3–16
- Power Doppler Imaging, 5–41
- Power on Memory Test, System Parameters menu page 1, 14–51
- Power on Status, System Parameters menu page 1, 14–51
- Power Up Sequence, 3–17
- PPSA Coefficient (1) [EC, UP], Preset Program menu page 4, 14–72
- PPSA Coefficient (2) [EC, UP], Preset Program menu page 4, 14–72
- Pre–Processing
  - B–Mode
    - Dynamic Range, 5–16
    - Edge Enhance, 5–23
    - Frame Averaging, 5–18
    - Image Softener, 5–20
  - Color Flow
    - Capture, 5–35
    - Diag Mode, 5–33
    - Frame Average, 5–34
    - High Resolution, 5–35
    - Penetration (Penet.), 5–34
  - Doppler, Dynamic Range, 5–54
  - M–Mode
    - Dynamic Range, 5–66
    - Edge Enhance, 5–67
- Precaution Levels, 2–2
- Preset
  - Sub–Menu, 3–41
  - Top Menu, 3–41
- Preset Parameters, 14–1, 14–7, 14–47, 14–64
  - Overview, 14–5
  - Save Values, 14–91
  - Top Menu, Defining a User Preset, 14–93
- Preset Program, 14–64
- Presets
  - Changing a Value, 14–10
  - Command lines, 14–8
  - Custom Display Contents, 14–12
  - How to Customize, 14–7
- Printing, Images, 13–6
- Probe, 3–27, 15–1
  - 546L, 15–24
  - 739L, 15–25
  - Activating, 3–29
  - Application, 15–6
  - B510, 15–27
  - Biopsy Guidezone(s), 5–22
  - Body Pattern Marker, 6–22
  - C358, 15–21
  - C364, 15–21
  - C386, 15–22
  - C551, 15–22
  - C721, 15–22
  - Care and Maintenance, 15–10
  - Cleaning, 15–13
  - Connecting, 3–27
  - Controls, 3–29
  - Coupling Gels, 15–19
  - Curved Array, Biopsy guidezone, 15–21
  - Curved Array (Convex), 15–21
  - CWD2, 15–29

Probe (cont'd)  
 CWD5, 15–29  
 Deactivating, 3–30  
 Description, 15–20  
 Disconnecting, 3–30  
 Disinfecting, 15–16  
 Dual Frequency Operation  
   Imaging Freq, 5–18  
   Penetration, 5–57  
 E721, 15–23  
 Environmental Requirements, 15–10  
 I739, 15–24  
 Immersion Levels, 15–15  
 L764, 15–25  
 LA39, 15–25  
 Labelling, 15–3  
 LD, 15–25  
 Linear  
   Color Flow Window, 5–32  
   Doppler Cursor, 5–54  
 Linear Array, Biopsy guidezone, 15–24  
 Naming Conventions, 15–20  
 Orientation Marker, 6–22  
 P509, 15–28  
 Planned Maintenance, 15–19  
 S220, 15–26  
 S222, 15–26  
 S316, 15–27  
 S317, 15–27  
 S611, 15–27  
 Safety, 15–11  
   Coupling gels  
     Applying, 15–19  
     Precautions, 15–19  
   Sector, biopsy guidezone, 15–26  
   Special Handling Instructions, 15–12  
   Specifications, 15–7  
   Storing, 3–30  
   T739, 15–25  
 Probe Direction Mark Display [EC], Custom Display  
   menu page 10, 14–31  
 Probe Name Menu, 3–43  
 PSA Calculation, System Parameters menu page 3,  
   14–54  
 Pulmonic Valve, 10–87, 10–93  
 Pulsatility Index (PI), 8–3, 9–23, 10–29, 11–5  
 Pulsed Wave Doppler, Introduction, 5–42  
 PWD  
   *See also Doppler*  
   Sub-Menu, 3–39  
   Top Menu, 3–39

PWD Alternative Scan [EC, UP, P, R], Custom  
   Display menu page 4, 14–21  
 PWD Base Line (%) [EC, R], Custom Display menu  
   page 15, 14–40  
 PWD Dynamic Range [EC, UP, P, R], Custom  
   Display menu page 4, 14–20  
 PWD Gain [EC, UP, P, R], Custom Display menu  
   page 4, 14–20  
 PWD Gray Scale Map [EC, UP, P, R], Custom  
   Display menu page 4, 14–20  
 PWD HPRF [EC, UP], Custom Display menu page 15,  
   14–40  
 PWD Rejection [EC, UP, P, R], Custom Display  
   menu page 4, 14–20  
 PWD Sample Volume (mm) [EC, UP, P, R], Custom  
   Display menu page 4, 14–20  
 PWD Velocity (cm/s) [EC, UP, P, R], Custom Display  
   menu page 4, 14–20  
 PWD Wall Filter (cm/s) [EC, UP, P, R], Custom  
   Display menu page 4, 14–20  
 PWD/CFM/PDI Penetration [EC, UP, P, R], Custom  
   Display menu page 4, 14–20  
 PWD/CFM/PDI Slant Auto Start [EC, UP, P, R],  
   Custom Display menu page 4, 14–20

## Q

Quality Assurance, 16–44  
 Baselines, 16–48  
 Checklist, 16–60  
 Frequency of tests, 16–45  
 Introduction, 16–44  
 Periodic checks, 16–48  
 Phantoms, 16–46  
 Record keeping, 16–59  
 Results, 16–49  
 System setup, 16–50  
 Test descriptions, 16–50  
   Axial distance measurement, 16–51  
   Axial resolution, 16–53  
   Contrast resolution, 16–57  
   Functional resolution, 16–56  
   Gray scale photography, 16–58  
   Lateral distance measurement, 16–52  
   Lateral resolution, 16–54  
   Penetration, 16–55  
 Typical tests, 16–45





## Index

---

### R

R Delay, 10–35

Radiology/Abdomen  
Calculations. *See General Calculations*  
Measurements. *See General Calculations*

RDelay Time of End Systole Prior Edge [EC, UP],  
Preset Program menu page 3, 14–70

RealTime Trace, Doppler, 5–60

Record, 3–50

Record 1, 3–47, 13–6, 13–9

Record 1 B&W, System Parameters menu page 5,  
14–59

Record 1 Color, System Parameters menu page 5,  
14–59

Record 2, 3–47, 13–6, 13–9

Record 2 B&W, System Parameters menu page 5,  
14–59

Record 2 Color, System Parameters menu page 5,  
14–59

Recording  
Audio, 19–29  
Image, 13–2

Red Shift, 6–20

Ref Scan, 10–35

Regulatory Labels  
Americas Systems, 2–22  
European Systems, 2–21

Rejection  
B–Mode, 5–23  
Doppler, 5–59  
M–Mode, 5–67

Report, System Parameters menu page 5, 14–59

Report Cursor Blink, System Parameters menu  
page 1, 14–49

Report Cursor Type, System Parameters menu page 1,  
14–49

Report Format [EC], Preset Program menu page 7,  
14–74

Report Video Inverse to Printer, System Parameters  
menu page 5, 14–59

Resistance Index (RI), 8–3, 9–23, 9–69, 10–29,  
11–5

Reverse, 3–46  
B–Mode, 5–15

Review Loop, 6–11

RI, 9–69

RI Calculation to bidirectional Flow [EC, UP], Preset  
Program menu page 4, 14–71

Rotation, 3–49  
B–Mode Image, 5–22

Routing Information, 20–4

RT Bifurc, 11–3

RT CCA, 11–3

RT ECA, 11–3

RT ICA, 11–3

RT ICA/CCA, 11–4

### S

S220, 15–26

S222, 15–26

S316, 15–27

S317, 15–27

S611, 15–27

Safety, 2–1

Acoustic Output, 2–14

BF Equipment, 2–9

CF Equipment, 2–9

Class I Equipment, 2–9

EMC, 2–10

Equipment and personnel safety

Explosion hazard, 2–5

Smoke and fire hazard, 2–6

Hazard Symbols, Icon Description, 2–3

Labels, Icon Description, 2–7

Leakage Current, 2–9

Patient

Diagnostic Information, 2–4

Identification, 2–4

Mechanical Hazards, 2–4

Patient safety

Acoustic Output hazard, 2–5

Diagnostic information, 2–4

Electrical hazard, 2–4

Mechanical hazards, 2–4

Patient identification, 2–4

Training, ALARA, 2–5

Precaution Levels, Icon Description, 2–2

Safety hazards, defined, 2–3–2–4

Safety icons, defined, 2–2

Warning Label, Location, 2–16

Sample Volume Length, Doppler, 5–56

Save Values, 14–91



- Scan Area, 3–48
  - Color Flow Window Size, 5–31
  - Position, B–Mode, 5–14
  - Size, B–Mode, 5–14
- Select, Color Flow, 5–29
- Sequence 1–8 [EC], Preset Program menu page 7, 14–81
- Service, Requesting, 16–38
- Set, 3–49
- Set–Up
  - Sub–Menu, 3–42
  - Top Menu, 3–42
- Setup Blue Back Start, System Parameters menu page 2, 14–51
- Single, 10–34
- Single Plane Ellipsoid Method, 10–15, 10–73
- Site Requirements, 3–2
- Slant Scan
  - Color Flow, 5–32
  - Doppler, 5–54
- Small Parts
  - Calculations. *See General Calculations*
  - Measurements. *See General Calculations*
- Smoke & Fire Hazard, 2–3, 2–6
- Smoothing, B–Mode, Image Softener, 5–20
- Soft Key, 3–34
- Soft Menu
  - B–Mode
    - 3D Mode, 5–20
    - Biopsy Guidezone(s), 5–22
    - Color, 5–20
    - Color Tag, 5–21
    - Dynamic Range, 5–16
    - Edge Enhance, 5–23
    - Focus Number, 5–17
    - Focus Position, 5–17
    - Frame Averaging, 5–18
    - Gray Scale Map, 5–16
    - Image Rotation, 5–22
    - Image Softener, 5–20
    - Imaging Freq, 5–18
    - Rejection, 5–23
  - Color Flow
    - 3D Mode, 5–39
    - ACE, 5–38
    - Color Capture, 5–35
    - Color Flow Map, 5–32
    - Color Tag, 5–37
    - Color Threshold, 5–35
    - Diag Mode, 5–33
    - Frame Average, 5–34
    - High Resolution, 5–35
    - MR–Flow, 5–39
    - MTI Filter, 5–33
    - Noise Blanker, 5–38
    - Packet Size, 5–36
    - Penet., 5–34
    - Persistence, 5–39
    - Slant Scan, 5–32
    - Spatial Filter, 5–36
    - Tag Position, 5–37
    - W.E. Cancel, 5–37
  - Control Panel, 3–34
  - Doppler
    - Calc Dir, 5–60, 5–61
    - Color, 5–58
    - Color Tag, 5–58
    - Dynamic Range, 5–54
    - HPRF, 5–59
    - Penetration, 5–57
    - RealTime Trace, 5–60
    - Rejection, 5–59
    - Sample Volume Length, 5–56
    - Slant Scan, 5–54
    - Sweep Speed, 5–57
    - Tag Position, 5–58
    - Wall Filter, 5–56
  - M–Mode
    - Color, 5–68
    - Dynamic Range, 5–66
    - Edge Enhance, 5–67
    - Gray Scale Map, 5–66
    - Rejection, 5–67
    - Sweep Speed, 5–68
- Software Option
  - Cardiology, 10–1
  - OB, 9–1
  - Vascular, 11–1
- Sort Criteria for Schedule List, System Parameters menu page 7, 14–63
- Spatial Average, Affect on frame rate, 5–36
- Spatial Filter, Color Flow, 5–36
- Speakers, 3–26
- Spec. Inv in Triplex [EC], Custom Display menu page 14, 14–39
- Spectral Gain, Doppler, 5–51
- Spectrum
  - Doppler, Display, 5–62
  - Gain, Doppler, 5–51
  - Invert, 5–52
    - Color Flow, 5–29
  - M–Mode, Display, 5–62
- Standby Time (min), System Parameters menu page 1, 14–50
- Standoff Setting [EC, UP, P, R], Custom Display menu page 1, 14–14
- Start Frame, 6–10



## Index

---

Stenosis Ratio. *See* % Stenosis

Stepper Increment [EC, UP], Preset Program menu  
page 4, 14–71

Stop, 3–50

Storage areas, 3–6

Stroke Volume, 8–5

Automatic Calculation, 8–6

Style Fun, Style Name, Page# Sep

E721, 15–23

Probe, E721, 15–23

Sub–Menu, 3–34, 3–36

Advanced Cardiac Measurement Option, 3–45

Auto Sequence, 3–43

B–Mode, 3–37

Body Patterns, 3–44

Cine, 3–43

Color Flow, 3–38

Comment, 3–44

CWD, 3–40

ECG, 3–42

Image Archive Option, 3–42

Image Recall, 3–45

M–Mode, 3–41

Measurement Pressed, 3–44

Preset, 3–41

PWD, 3–39

Setup, 3–42

Sub–Menu Select, 3–36

Summary Report

GYN, Layout, 9–70

OB, 9–24

Vascular

Displaying, 11–6

Editing, 11–7

Sweep Speed

Doppler, 5–57

M–Mode, 5–68

Sync Selectn, 10–35

System

Care, 16–11

Specifications, 16–3

System Error Erase Interval, System Parameters  
menu page 1, 14–51

System Parameters, Preset Parameters, 14–47

Systolic/Diastolic Ratio (S/D), 8–3, 9–23, 10–29,  
11–5

## T

T739, 15–25

TA Max Auto [EC], Preset Program menu page 3,  
14–69

Tag Position

Color Flow, 5–37

Doppler, 5–58

Teichholz Method, 10–9, 10–68

TGC, 5–11

TGC Curve, Disable, 5–9

TGC Depth Remap [EC], Custom Display menu  
page 10, 14–31

Time Adjustment, 14–3, 14–3

Time Format, System Parameters, System  
Parameters menu page 1, 14–49

Time Unit [EC], Preset Program menu page 3, 14–68

Timeline, Common Controls, 5–64

Timeline Format [EC], Custom Display menu page  
12, 14–34

Timeline Sweep Method [EC], Custom Display menu  
page 12, 14–34

Timeline Sweep Speed [EC, UP, P, R], Custom  
Display menu page 3, 14–18

Top Menu, 3–34

Auto Sequence, 3–43

Automatic, 3–36

B–Mode, 3–37

Cine, 3–43

Color Flow, 3–38

CWD, 3–40

Disable, 3–36

ECG, 3–42

Image Archive Option, 3–42

Mode Default Menu, 3–35

Preset, 3–41

PWD, 3–39

Set–Up, 3–42

Top Menu Select, 3–35

Trace Auto, 8–3, 8–8, 10–28, 11–5  
Measurement Method, 9–23

Trackball, 3–49

Trackball Function at Freeze [EC], Preset Program  
menu page 1, 14–66

Transf Calcs, 8–3, 9–23, 10–30, 11–5

Tricuspid Valve, 10–88, 10–95

Troubleshooting, 16–19

Color printer, 16–28

Display Messages, 16–26

Multi Image Camera, 16–29

Operation Error Message, 16–28

System Error Message, 16–27

VCR, 16–29

Warning Message, 16–33

Type BF Equipment, 2–9

Type CF Equipment, 2–9

## U

Urology

Mechanical Stepping Device, Mechanical Stepper/Needle Placement Guide, 12–10

Stepper Volume Calculation, 12–7

Stepper Volume Formula, 12–9

Use A/N Keys for User Define, System Parameters menu page 3, 14–54

User

Data Backup, 14–104

Define Function, Programming, 14–98

Programmed Calculations, Tables, 9–22

Tables, 9–22

User Define Key Names and Locks, System Parameters menu page 3, 14–54

User ID 1–8, System Parameters menu page 6, 14–61

User Sequence 1–8, Programming, 14–81

User Table Editor, Preset Program menu page 7, 14–77

Ut–H, 9–68

Ut–L, 9–68

Ut–W, 9–68

Uterine

Height, 9–68

Length, 9–68

Width, 9–68

Uterine Cavity, 9–20

## V

Vascular, 11–1

Advanced, Calculation Formulas, 11–18

Advanced Option, 11–9

Calculation Formulas, 11–8

Calculations, 11–3

Carotid Artery Measurements, 11–3

Exam Preparation, 11–2

Measurements, 11–3

A/B Ratio, 11–5

Heart Rate, 11–5

Pulsatility Index (PI), 11–5

Resistance Index (RI), 11–5

RT(LT) Bifurc, 11–3

RT(LT) CCA, 11–3

RT(LT) ECA, 11–3

RT(LT) ICA, 11–3

RT(LT) ICA/CCA, 11–4

S/D Ratio, 11–5

Summary Report

Displaying, 11–6

Editing, 11–7

Recording, 11–7

Vascular Calculations, % Stenosis, 11–5

VCR, 13–12

Basic Recording, Audio, 19–29

Cleaning, 16–15

Counter Display, 19–18

Counter Reset, Control Z, 19–21

VCR Color Intensity [EC, R], Preset Program menu page 2, 14–67

VCR Operating Instructions, 19–1

VCR Play Audio [EC], Preset Program menu page 2, 14–67

VCR Play Audio Volume (%) [EC, R], Preset Program menu page 2, 14–67

VCR Rec Audio Doppler [EC], Preset Program menu page 2, 14–67

VCRPB Calib Method [EC, R], Preset Program menu page 3, 14–68

Velocity Ratio (A/B), 8–3, 9–23, 10–29, 11–5

Velocity Scale

Color Flow, 5–30

Doppler, 5–52

Velocity Unit [EC], Preset Program menu page 3, 14–68

Venous, Comments Page, 11–16

Video Cassette Recorder. *See* VCR

Volume, 8–2, 10–19

Doppler Audio, 5–50

VTR. *See* VCR



## Index

---

### W

W.E. Cancel, Color Flow, 5–37  
Wall Filter, Doppler, 5–56  
Warning, 2–2, 2–5, 13–5  
    Label Locations, 2–16  
Warranties  
    Duration, 16–7  
    Scope, 16–7  
Wheels, 3–14  
Word Wrap, 6–19

### Z

Zoom, 3–49  
    Display Format, 5–86  
    Introduction, 6–2  
    M-Mode, 5–65, 6–4  
    Methods  
        Acoustic Zoom, 6–2  
        Display, 6–3  
    Multi-Image, 6–5  
    Size, 3–49  
Zoom Auto Start [EC, UP, P], Custom Display menu  
    page 1, 14–13  
Zoom Factor [EC, UP, P, R], Custom Display menu  
    page 1, 14–13  
Zoom Reference [EC, UP, P], Custom Display menu  
    page 1, 14–13  
Zoom Reference Update [EC, UP], Custom Display  
    menu page 10, 14–31